

Bio SB 
BIOSCIENCE FOR THE WORLD



Products for Molecular Pathology

Volume 5

Bio SB, Inc. Mission Statement



At Bio SB our passion is providing biomedical laboratories with the tools to improve the diagnosis, prognosis and therapy prediction that benefit patients worldwide.

Bio SB is proud to develop, manufacture and distribute unique products for Immunohistochemistry (IHC), Immunocytochemistry (ICC), Immunofluorescence (IF), Fluorescent in situ Hybridization (FISH) and Chromogenic in situ Hybridization (CISH) technologies that meet the highest international standards for applications in Molecular Pathology. With a network of over 60 distributors that represent Bio SB in the Americas, Europe, Africa, Asia, Australia and the Middle East, Bio SB can ensure that our products can reach clinicians, pathologists, researchers, and patients worldwide.

Bio SB has a talented, dedicated, and vibrant team of professionals that excel in the R&D, Production and Marketing of immunochemicals and nucleic acids for Cancer Diagnostics and Research, Microbiology, Immunology and Immunotherapy. Bio SB is a CE/IVDD and ISO 13485:2016 certified company that manufactures and develops products in accordance with USA FDA QSR 21 CFR Part 820 cGMP. These guidelines enable us to produce IVD products that meet the highest *in vitro* diagnostic standards.

Our mission is to develop and supply high-quality products, while providing outstanding technical service to our customers. We develop rabbit and mouse monoclonal antibodies and detection systems for Immunohistochemistry, Immunocytochemistry and Immunofluorescence using proprietary technologies. Our ImmunoDetector, PolyDetector, PolyDetector Plus and TintoFast Mohs PolyDetector systems are unique, high-sensitivity detection systems already in use in laboratories worldwide.

The robotic automated TintoStainer Plus system enables Bio SB to offer the most reliable, flexible, and open automated system for use in IHC, ICC and IF applications. The Bio SB TintoStainer Plus is an automated system for use in laboratories that need to implement a wide variety of protocols yet need to maintain low operating costs. The Bio SB TintoStainer Plus is a fast, flexible, open and affordable solution which allows for fast implementation of laboratory protocols.

Bio SB continues to innovate a new line of innovative, sensitive, specific, and cost effective IVD Mouse and Rabbit Monoclonal Antibodies to help laboratories worldwide improve patient care. This allows Bio SB to have one of the largest selections of Rabbit Monoclonal antibodies for use in Diagnostic IHC.

Our strategic alliance with ZytoVision of Germany has added a complete range of products for Chromogenic and Fluorescent in situ hybridization. These technologies are designed for the evaluation and detection of genetic aberrations such as translocations, deletions, amplifications, chromosomal aneuploidy and gene amplifications. We also offer outstanding fast detection and discrimination of human pathogenic viruses in formalin-fixed, paraffin-embedded tissue sections, cell samples, blood or bone marrow smears, and metaphase chromosome spreads.

As industry research, genetic sequencing, and explosive growth in genomics pushes molecular pathology forward, Bio SB will continue to offer a large portfolio of products at a competitive price to Laboratories, Clinicians, Researchers and patients worldwide.

Dr. Alfonso Heras
President & CEO

Products for Molecular Pathology

Antibodies for IHC

YAP1- MMab
 PD-L1 - RMab
 ARID1A - MMab
 p16 - RMab
 HEG1 - MMab
 FGFR-3- MMab
 HER-2 neu - RMab

Ancillaries for IHC

TintoDeparaffinator Citrate
 Tween 20
 ChromoProtector
 FluoroMounter
 TintoHematoxylin Automation

CISH & FISH

ZytoLight FISH
 ZytoLight FlexISH
 ZytoDot CISH
 ZytoFast CISH
 Ancillaries & Equipment

Detection Systems

ImmunoDetector HRP with DAB
 PolyDetector HRP with AEC
 PolyDetector Plus HRP with DAB
 AmpliDetector Plus FITC
 Mohs PolyDetector with DAB
 Mohs PolyDetector Plus with DAB

Microarrays & Control Slides

11-Core NH-TMA
 7-Core NH Lymphoid TMA
 31-Core Multi Cancer CLMA
 7-Core EGFR Cell Line
 H.Pylori ID-Array (2-Core)

Equipment

TintoDetector
 TintoDetector Mini
 TintoStainer Plus
 TintoRetriever Pressure Cooker



BIO SB manufactures and develops products in accordance with FDA QSR 21 CFR Part 820 cGMP and ISO 13485:2016. These guidelines enable us to produce an IVD product that meets the highest in vitro diagnostic standards.

United States

Phone: 1-800-561-1145

International

Phone: 1-805-692-2768

Global Distribution

www.biosb.com/distributors

Table of Contents

New Products5

Antibodies by Application13

I. Immunohistochemistry23

Antibodies for Immunohistochemistry24

Fast Mohs Immunohistochemistry 218

TintoFast Antibodies 219

Fast Mohs Detection Systems..... 224

Immunofluorescence Autoimmunity..... 229

Direct Immunofluorescence Antibodies..... 230

InDirect Immunofluorescence Antibodies..... 235

InDirect Immunofluorescence Detection Systems..... 247

Detection Systems for IHC..... 239

PolyDetector AP Detection Systems 240

PolyDetector HRP Detection Systems 241

PolyDetector Plus HRP Detection Systems 242

ImmunoDetector AP Detection Systems 243

ImmunoDetector HRP Detection Systems..... 244

Complete Multiplex Detection Systems for IHC 247

PIN MultiDetector AP & HRP Kit..... 248

PNI Carcinoma MultiDetector AP & HRP..... 250

Cervical Cancer MultiDetector AP & HRP 252

Substrate-Chromogen Systems..... 255

HRP Substrate-Chromogen Detection Systems..... 256

PolyDetector Liquid AEC HRP Red Ready-To-Use 256

PolyDetector Liquid DAB HRP Brown Kit 256

PolyDetector HRP Black Kit 256

PolyDetector HRP Green Kit 256

PolyDetector HRP Yellow Kit 256

PolyDetector HRP Murberry Kit 256

PolyDetector HRP Fuchsia Kit..... 256

PolyDetector HRP Blue Plus Kit 256

AP Substrate-Chromogen Detection Systems 257

PolyDetector Alk Blue Ready-To-Use 257

PolyDetector Alk Magenta..... 257

PolyDetector ALK Scarlet..... 257

Environmentally Friendly Ancillaries for IHC 259

TintoDeparaffinator Citrate 20X 260

TintoDeparaffinator EDTA 20X 260

TintoDeparaffinator Hot Rinse 20X..... 260

XyGreen PermaMouncer..... 260

ChromoProtector 261

ImmunoRetriever 20X with Citrate..... 261

ImmunoRetriever 20X with EDTA 261

ImmunoDNA Protein Block/Antibody Diluent 261

Immuno/DNA Washer..... 261

ImmunoDNA Digestor 262

Hematoxylin Counterstainer..... 263

Nuclear Fast Red Counterstainer..... 263

Methyl Green Counterstainer..... 263

Hematoxylin Automation 263

Negative Controls..... 263

Mounting Media 263



| | |
|---|-----|
| Aqua Mounter | 263 |
| PermaMounter..... | 263 |
| FluoroMounter..... | 263 |
| FluoroMounter with DAPI | 264 |
| Tween 20 | 265 |
| CytoLayer Slide Treatment Solution..... | 265 |
| PolyDetector Peroxidase Blocker | 265 |
| PolyDetector AP Blocker..... | 265 |
| ImmunoDetector Biotin Blocker | 265 |
| ImmunoDNA Background Blocker..... | 265 |

| | |
|---|------------|
| Tissue and Cell Line Microarrays | 265 |
| Tissue Microarrays | 266 |
| Cell Line Microarrays | 267 |
| Infectious Disease Cell Line Microarrays..... | 268 |

II. FISH & CISH Products

271

| | |
|--|-----|
| FISH-ZytoLight..... | 272 |
| FISH-ZytoLight Probes by Application | 273 |
| ZytoLight FISH probes | 278 |
| ZytoLight FISH Kits..... | 279 |
| FlexISH | 280 |
| FlexISH Probes..... | 280 |
| FlexISH Kits | 280 |
| CISH-ZytoDot..... | 282 |
| CISH-ZytoDot Probes by Application..... | 285 |
| ZytoDot CISH probes..... | 287 |
| ZytoDot CISH Kits..... | 287 |
| CISH-ZytoFast..... | 288 |
| CISH-ZytoFast Probes by Application | 288 |
| ZytoFast Probes..... | 288 |
| ZytoFast Detection Systems | 288 |
| Ancillaries for FISH & CISH..... | 289 |

III. Equipment and Slides

291

| | |
|--|-----|
| TintoStainer Plus..... | 294 |
| TintoDetector..... | 296 |
| TintoDetector Mini..... | 302 |
| TintoRetriever Digital Pressure Cooker with Thermometer .. | 305 |
| Plastic Staining Dish..... | 307 |
| Slide Holder | 307 |
| TintoRetriever Slide Thermometers | 307 |
| Hydrophilic Plus Slides for Molecular Pathology..... | 309 |
| TintoDetector Cap Gap Plus Slides..... | 309 |

| | |
|--|------------|
| Molecular Pathology Technical Reference | 312 |
| Immunophenotyping of Cancer for IHC | 354 |
| ZytoLight Interpretation Guides..... | 364 |
| ZytoLight Lung Cancer Probe Pool | 370 |
| ZytoLight Sarcoma Probe Pool | 371 |
| Ordering Information..... | 372 |
| Distributors | 373 |




A large circular frame containing a microscopic image of tissue, likely a histological section. The tissue shows various cell types and structures, including what appears to be a glandular or epithelial structure. The frame is surrounded by several overlapping, semi-circular bands of different colors and textures, suggesting layers of tissue or different types of biological structures. The background is a light purple and blue gradient with a geometric, low-poly pattern.


Bio SB
BIOSCIENCE FOR THE WORLD




New Antibodies for Molecular Pathology




C-Met, RMab



MDM2, RMab




PMS2, RMab



CD10, RMab

New Antibodies Optimized for Automated IHC

Bio SB is proud to announce our new line of antibodies for immunohistochemistry. As well as developing new biomarkers for immunohistochemistry applications, we have developed new clones for currently existing antibodies that are specifically designed for use in automated IHC systems. Along with increased sensitivity and overall clearer picture, our highly concentrated, high quality antibodies have a much higher dilution factor compared to our competitors, allowing you to run more tests per mL of concentrated antibody!




C-Met / HGFR, RMAb (RBT-C-Met/HGFR)

C-Met is deregulated in many types of human malignancies, including cancers of kidney, liver, stomach, breast, and brain. Normally, only stem cells and progenitor cells express MET, which allows these cells to grow invasively in order to generate new tissues in an embryo or regenerate damaged tissues in an adult. However, cancer stem cells are thought to hijack the ability to express MET, and thus become the cause of cancer persistence and spread to other sites in the body (metastasis).


Application: Kidney & Urothelial Cancer, Liver Cancer, Breast Cancer, Neural & Neuroendocrine Cancer, Colon & Gastrointestinal Cancer, Lung Cancer

c-Myc, RMAb (RBT-CMYC)



c-Myc (RBT-CMYC) Clone Comparison



c-Myc (EP121) on Tonsil Tissue



c-Myc (RBT-CMYC) on Tonsil Tissue



Calcitonin, RMAb (RBT-Calcitonin)

Immunohistochemical staining with Calcitonin antibody has proven to be an effective way of demonstrating the existence of Calcitonin-producing cells in the thyroid. C-cell Hyperplasia and Medullary Thyroid Carcinomas stain positive for Calcitonin. Studies of Calcitonin have resulted in the identification of a wide spectrum of C-cell proliferative abnormalities.


Application: Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Head & Neck Cancer, Cytopathology

New Antibodies Optimized for Automated IHC

Calretinin, RMab (RBT-CALB2)

Calretinin antibody has been shown to be useful in differentiating Mesothelioma from Adenocarcinomas of the lung and other sources. It is also useful in differentiating adrenal-cortical neoplasms from Pheochromocytomas.


Application: Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Head & Neck Cancer, Cytopathology



CD4, MMab (BSB-179)


CD4 antigen is involved in the recognition of Type II Major Histocompatibility Complex antigens (MHC-II). CD4 is also the receptor for Human Immunodeficiency Virus (HIV). It is present on most T-helper cells and normal thymocytes.

Application: Melanoma & Skin Cancer, Lymphoma




CD10, RMab (RBT-CD10)


CD10 (RBT-CD10) Clone Comparison



CD10 (EP195) on Prostate Tissue




CD10 (RBT-CD10) on Prostate Tissue




CD138, RMab (RBT-CD138)

CD138/syndecan-1 is a useful marker for labeling normal and neoplastic plasma cells and Plasmacytoid Lymphomas. It is a selective marker for B-cell Lymphoblastic Leukemia and Lymphoplasmacytoid Leukemia. It is lost from the apoptotic myeloma cells, and thus, is a useful marker for viable Myeloma cells. Various forms of Hodgkin's Disease have also shown positive staining with this antibody.

Application: Hematopoietic, Lymphoma, Rejection & Autoimmunity



New Antibodies Optimized for Automated IHC



CDX2, RMAb (RBT-CDX2)


The CDX2 protein is expressed in Primary and Metastatic Colorectal Carcinomas and has also been demonstrated in the intestinal metaplasia of the stomach and intestinal-type gastric cancer. It is not expressed in the normal gastric mucosa. Anti-CDX2 antibody has been useful in distinguishing the gastrointestinal origin of Metastatic Adenocarcinomas and carcinoids. A high percentage of Mucinous Carcinomas of the Ovary also stain positively with this antibody, as well as Carcinomas from the upper gastrointestinal tract.

Application: Colon & GI Cancer

Clusterin / Apolipoprotein J, RMAb (RM437)

Clusterin is expressed in a wide variety of hematopoietic and non-hematopoietic tumors. Overexpression of Clusterin is associated with poor prognosis in breast cancer and chemosensitivity in cervical cancer.

Application: Rejection & Autoimmunity



Cytokeratin 5 & 6, RMAb (RM431)


Cytokeratin 5 & 6 antibody is positively seen in nearly 100% of Malignant Mesotheliomas and are rarely seen in Lung Adenocarcinomas. CK 5 and 6 can positively be seen in undifferentiated Large-cell Carcinoma as well as Squamous Carcinoma. Fewer than 10% of Carcinomas of the breast, colon, and prostate stain positively for this marker. Cytokeratin 5 & 6 antibody has also been used successfully as a myoepithelial cell marker in the prostate to determine malignancy.

Application: Mesothelioma, Lung Cancer, Prostate Cancer, Breast Cancer, Melanoma & Skin Cancer

EZH2, RMAb (RBT-EZH2)

EZH2 may inhibit apoptosis in ectopic Breast Cancer, overexpression is correlated with poor clinical outcome in Prostate and Gastric Cancer. EZH2 overexpression is positively correlated with epithelial-mesenchymal transition, tumor size, lymphatic invasion and TNM stage, and poor disease-free survival and overall survival of patients. EZH2 regulates ADAR1 expression, and EZH2 and BET inhibitors show synergistic inhibition in pancreatic cancer. EZH2 expression is significantly higher in BAP1-mutant renal clear cell carcinoma patients with progressed stage, grade, nodal invasion, and metastasis.


Application: Breast Cancer, Gastric Cancer




New Antibodies Optimized for Automated IHC

GLUT1, RMAb (RBT-GLUT1)


GLUT1 (RBT-GLUT1) Clone Comparison



GLUT1 (EP141) on Prostate Adenocarcinoma Tissue



GLUT1 (RBT-GLUT1) on Prostate Adenocarcinoma Tissue




LEF1, RMAb (RBT-LEF1)

LEF1 is highly overexpressed and associated with disease progression and poor prognosis in B-cell chronic lymphocytic leukemia. Strong nuclear expression of LEF1 has been observed in the majority of chronic lymphocytic leukemia/small lymphocytic lymphoma cases and LEF1 is not detected in other small B cell lymphomas. Gene expression profiling revealed overexpression of LEF1 in chronic lymphocytic leukemia (CLL)/small lymphocytic lymphoma (SLL). LEF1 immunostaining has been detected in all neoplastic cells of CLL/SLL cases. LEF1 was identified in 50% of high grade follicular lymphoma and 38% of diffuse large B-cell lymphoma, but not in mantle cell lymphoma or marginal zone lymphoma.


Application: Leukemia & Histiocytic, Lymphoma, Colon & Gastrointestinal Cancer, Brain Cancer

MDM2, RMAb (RBT-MDM2)


MDM2 (RBT-MDM2) Clone Comparison



MDM2 (BSB-64) on Testis Tissue




MDM2 (RBT-MDM2) on Testis Tissue



Myeloperoxidase / MPO, MAb (BSB-180)

Myeloperoxidase detects granulocytes and monocytes in blood and precursors of granulocytes in the bone marrow. This antibody can detect myeloid cell populations of the bone marrow as well as in other sites.



Application: Leukemia & Histiocytic




New Antibodies Optimized for Automated IHC

MyoD1, RMAb (RBT-MYOD1)

MyoD1 (RBT-MYOD1) Clone Comparison





MyoD1 (EP212) on Rhabdomyosarcoma Tissue




MyoD1 (RBT-MYOD1) on Rhabdomyosarcoma Tissue

OCT4 / POU5F1, RMAb (RBT-PCT4)

OCT4 (RBT-OCT4) Clone Comparison



OCT4 (EP141) on Testis Tissue




OCT4 (RBT-OCT4) on Testis Tissue

OLIG2, RMAb (RBT-OLIG2)



OLIG2 is universally expressed in glioblastoma and other diffuse gliomas (astrocytomas, oligodendrogliomas and oligoastrocytomas), and is a useful positive diagnostic marker of these brain tumors.

Application: Neural and Neuroendocrine Cancer




PMS2, RMAb (RBT-PMS2)

PMS2 (RBT-PMS2) Clone Comparison



PMS2 (EP51) on Prostate Adenocarcinoma




PMS2 (RBT-PMS2) on Prostate Adenocarcinoma

New Antibodies Optimized for Automated IHC

S100P, RMAb (RBT-S100P)

S100P antibody is an early development marker of pancreatic carcinogenesis and can be used as a marker for pancreatic ductal adenocarcinoma. It may also serve as a predictor of distant metastasis and poor survival in non-small cell lung carcinomas.


Application: Pancreatic Cancer, Lung Cancer



Synaptophysin, RMAb (RBT-Synaptophysin)

Synaptophysin reacts with neuroendocrine cells of human adrenal medulla, carotid body, skin, pituitary, thyroid, lung, pancreas and gastrointestinal mucosa. Positive staining is seen in neurons of the brain, spinal cord, retina, and Paneth's cells in the gastrointestinal tract and gastric parietal cells. This antibody identifies normal neuroendocrine cells and neuroendocrine neoplasms. Synaptophysin is an independent broad-range marker of neural and neuroendocrine differentiation.


Application: Lung Cancer, Colon & Gastrointestinal Cancer, Gall Bladder & Pancreatic Cancer, Sarcoma & Soft Tissue, Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Undifferentiated Tumor, Carcinomas Of Unknown Primary Site



TFF3, MAb (BSB-181)

TFF3 has been identified as an indicator of Barrett's Esophagus, a premalignant condition indicating predisposition to esophageal adenocarcinoma. In cases of acid reflux, TFF3 and p53 can be used to help identify likely cases of dysplasia. Overexpression of TFF3 has been found to promote proliferation and invasion in Cervical Cancer cells, through the regulation of E-Cadherin. Expression of TFF3 has also been found to be decreased in Colorectal Cancer, and may be a prognostic indicator with less TFF3 expression indicating higher pathologic stages of the tumor.


Application: Colon & Gi Cancer, Cervical Cancer



VEGFA, RMAb (RBT-VEGFA)

VEGFA has been implicated with poor prognosis in breast cancer. The overexpression of VEGFA may be an early step in the process of metastasis, a step involved in the "angiogenic" switch. Although VEGFA has been correlated with poor survival, its exact mechanism of action in the progression of tumors remains unclear. VEGFA is also released in rheumatoid arthritis in response to TNF-alpha, increasing endothelial permeability and swelling and also stimulating angiogenesis (formation of capillaries). Once released, VEGFA may elicit several responses. It may cause a cell to survive, move, or further differentiate.

Application: Angiosarcoma, Angioma



New Antibodies Optimized for Automated IHC


| Antibody | Clone | Application |
|----------------------------|-------------------|---|
| c-Met/HGFR | RBT-c-Met/HGFR | Kidney & Urothelial Cancer, Liver Cancer, Breast Cancer, Neural & Neuroendocrine Cancer, Colon & Gastrointestinal Cancer, Lung Cancer |
| c-Myc | RBT-CMYC | Leukemia & Histiocytic, Lymphoma, Prostate Cancer |
| Calcitonin | RBT-Calcitonin | Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Head & Neck Cancer, Cytopathology |
| Calretinin | RBT-CALB2 | Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Head & Neck Cancer, Cytopathology |
| CD10 | RBT-CD10 | Hodgkin's & Non-Hodgkin Lymphoma, Lymphoma, Kidney & Urothelial Cancer, Liver Cancer, Gall Bladder & Pancreatic Cancer, Endometrial & Genital Cancer, Breast Cancer |
| CD138 | RBT-CD138 | Hematopoietic, Lymphoma, Rejection & Autoimmunity |
| CD4 | BSB-179 | Melanoma & Skin Cancer, Lymphoma |
| CDX2 | RBT-CDX2 | Colon & GI Cancer |
| Clusterin/Apolipoprotein J | RM437 | Rejection & Autoimmunity |
| Cytokeratin 5 & 6 | RM341 | Mesothelioma, Lung Cancer, Prostate Cancer, Breast Cancer, Melanoma & Skin Cancer |
| EZH2 | RBT-EZH2 | Breast Cancer, Gastric Cancer |
| GLUT1 | RBT-GLUT1 | Colon & GI Cancer, Breast Cancer, Mesothelioma |
| LEF1 | RBT-LEF1 | Leukemia & Histiocytic, Lymphoma, Colon & Gastrointestinal Cancer, Brain Cancer |
| MDM2 | RBT-MDM2 | Sarcoma, Testicular Cancer |
| Myeloperoxidase/MPO | BSB-180 | Leukemia & Histiocytic |
| MyoD1 | RBT-MYOD1 | Rhabdomyosarcoma |
| OCT4/POU5F1 | RBT-OCT4 | Ovarian Cancer, Testicular Cancer, Germ Cell Tumor |
| OLIG2 | RBT-OLIG2 | Neural and Neuroendocrine Cancer |
| PMS2 | RBT-PMS2 | Colon & GI Cancer |
| S100P | RBT-S100P | Pancreatic Cancer, Lung Cancer |
| Synaptophysin | RBT-Synaptophysin | Lung Cancer, Colon & Gastrointestinal Cancer, Gall Bladder & Pancreatic Cancer, Sarcoma & Soft Tissue, Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer |
| TFF3 | BSB-181 | Colon & GI Cancer, Cervical Cancer |
| VEGFA | RBT-VEGFA | Angiosarcoma, Angioma |

Bio SB performs R&D, production, distribution and marketing of unique products for IHC, ICC, FISH, CISH and PCR technologies that meet the highest international standards for applications in Molecular Pathology, Cancer Diagnostics & Research, Microbiology, Immunology and Genetics.

BIO SB manufactures and develops products in accordance with FDA QSR 21 CFR Part 820 cGMP, CE/IVD and ISO 13485:2016 standards. These guidelines enable us to produce an IVD product that meets the highest in vitro diagnostic standards.

Antibodies by Application

All Bio SB Antibodies Organized by Application



Alzheimer's Disease

| | |
|----------------------------|---------|
| APOE4, MMab | BSB-170 |
| pTau Phospho/Thr 181, MMab | BSB-176 |

Breast Cancer Antibodies


| | |
|--|------------|
| Actin, Smooth Muscle, MMab | BSB-15 |
| Annexin VII, RMab | EP367 |
| ARID1A, RMab | EP303 |
| ATM, RMab | EP327 |
| Aurora B, RMab | RM278 |
| Bax, RMab | E63 |
| Bcl-x, RMab | EP94 |
| Beta-Catenin, MMab | 14 |
| CA 15-3, MMab | DF3 |
| Cathepsin K, MMab | BSB-172 |
| CA-125, RMab | EP48 |
| Caldesmon, MMab | BSB-19 |
| Calponin, MMab | BSB-20 |
| CD142 / TF, MMab | BSB-143 |
| CD146 / MUC18 / Mel-CAM, MMab | BSB-122 |
| CDK2, RMab | RBT-CDK2 |
| CITED1, MMab | BSB-177 |
| Claudin-1, RPab | Polyclonal |
| Claudin-7, RMab | EP399 |
| Collagen IV, RMab | RBT-COL4 |
| CXCL12 / SDF-1, MMab | BSB-165 |
| CXCR4 / CD184 / Fusin, RMab | EP394 |
| Cyclin D3, MMab | BSB-171 |
| Cytokeratin 14, MMab | LL002 |
| Cytokeratin 35BH11, MMab | 35betaH11 |
| Cytokeratin 8, RMab | EP17 |
| Cytokeratin LMW/AE1, MMab | AE1 |
| E-Cadherin, RMab | EP6 |
| EGFR Phospho, RMab | EP11 |
| EMA, MMab | E29 |
| Estrogen Receptor, MMab | BSB-1 |
| Estrogen Receptor, RMab | RBT11 |
| Estrogen Receptor, RMab | RM292 |
| Estrogen Receptor Alpha, RMab | EP1 |
| EZH, RMab | RBT-EZH2 |
| FOXA1/HNF-3A, RPab | Polyclonal |
| GAB1 / GRB2-associated binding protein 1, MMab | BSB-155 |

Breast Cancer Antibodies


| | |
|--|-----------|
| GATA3, MMab | L50-823 |
| GATA3, RMab | EP368 |
| GCDFP-15, MMab | 23A3 |
| HER-2 neu, RMab | EP3 |
| HER-2 neu, MMab | BSB-3 |
| HER-2 neu, RMab | RBT-HER2 |
| HER-2 neu Phospho, RMab | EP123 |
| HER-3, RMab | RBT-HER3 |
| HLA-DR Alpha Chain, RMab | EP96 |
| Ki-67, RMab | EP5 |
| Ki-67, RMab | RM360 |
| Laminin-R / RPSA, MMab | BSB-144 |
| Mammaglobin, RMab | EP249 |
| Maspin, MMab | BSB-92 |
| MCM3, RMab | EP202 |
| MDR-1, MMab | JSB-1 |
| MDR-1, RMab | EP271 |
| MGMT / AGAT, RMab | EP337 |
| MUC1, RMab | EP85 |
| Myosin Smooth Muscle Heavy Chain, MMab | BSB-17 |
| Nestin, RMab | EP287 |
| NGFR, MMab | BSB-18 |
| NKX3.1, RMab | RM430 |
| P120 Catenin, RMab | EP66 |
| p27, MMab | SX53G8 |
| p53, MMab | DO7 |
| p63, RMab | EP174 |
| p63, MMab | 4A4 |
| p120, RMab | EP66 |
| PELP1, RMab | RBT-PELP1 |
| PRAME, RMab | RBT-PRAME |
| Progesterone Receptor, RMab | EP2 |

Breast Cancer Antibodies

| | |
|-----------------------------------|-------------|
| Progesterone Receptor, MMab | BSB2 |
| Progesterone Receptor RBT22, RMab | RBT22 |
| PTEN, RMab | RM265 |
| PTEN, RMab | RBT-PTEN |
| TRPS1, RMab | EP392 |
| TAG-72, MMab | BSB-21 |
| TDP-43 / TARDBP, MMab | BSB-166 |
| Topoisomerase II alpha, RMab | RBT-Topo 2a |




HER-2



Progesterone Receptor

Carcinoma of Unknown Primary Site


| | |
|---------------------------------------|---------|
| Cytokeratin 7, RMab | RM284 |
| Cytokeratin 20, MMab | Ks20.8 |
| Cytokeratin Cocktail, AE1 & AE3, MMab | AE1&AE3 |
| Cytokeratin MNF116, MMab | MNF116 |
| Cytokeratin Oscar, MMab | OSCAR |
| Napsin A, RMab | EP205 |
| S100 Beta, RMab | EP32 |




Estrogen Receptor



GATA-3





CK MNF116







Ki-67

Cervical Cancer Antibodies

| | |
|-------------------------|------------|
| Cyclin B1, RMab | RBT-B1 |
| Cyclin B1, RMab | RM281 |
| Cyclin E1, RMab | EP126 |
| Cytokeratin 17, MMab | BSB-33 |
| Cytokeratin 17, RMab | EP98 |
| FGFR-3, MMab | BSB-150 |
| FOXL2, RPab | Polyclonal |
| HIF-1a, RMab | EP118 |
| HPV, MMab | BSB-66 |
| HPV16, MMab | CAMVIR-1 |
| HSP-27, MMab | G3.1 |
| HSP70, RMab | RM432 |
| MCM2, RMab | RBT-MCM2 |
| MCM5, RMab | RBT-MCM5 |
| p14 ARF/CDKN2A, RMab | RBT-p14 |
| p16, MMab | 16P04,JC2 |
| p16, RMab | RBT-p16 |
| p16, RMab | RM267 |
| Stathmin, RMab | EP247 |
| Topoisomerase IIa, RMab | RBT-Topo2a |
| VEGFA, RMab | RBT-VEGFA |

**HPV****p16****Endometrial Cancer**



| | |
|----------------------|------------|
| bcl-2, RMab | EP36 |
| hPL, RPab | Polyclonal |
| Inhibin Alpha, MMab | R1 |
| Inhibin Alpha, RMab | EP378 |
| LIN28, RMab | EP150 |
| p21, MMab | DCS-60.2 |
| PELP1, RMab | RBT-PELP1 |
| Retinoblastoma, MMab | 1F8 |

**PELP1****MUC5AC****MSH2****PMS2****Colon & G.I. Antibodies**

| | |
|-----------------------------------|-------------------|
| ATM, RMab | EP327 |
| Beta-Catenin, RMab | RM276 |
| BRAF V600E, RMab | RM8 |
| CA 19-9, MMab | 121SLE |
| Cadherin 17/ Li-Cadherin, RMab | EP86 |
| CD44, MMab | BSB-12 |
| CD45RA, MMab | 4KB5 |
| CD73 / NT5E, RMab | RM431 |
| CDX2, RMab | RBT-CDX2 |
| CDX2, RMab | EP25 |
| CEA, MMab | BSB-13 |
| CEACAM7, MMab | BSB-168 |
| Claudin-1, RPab | Rabbit Polyclonal |
| COX-2, RMab | RBT-COX2 |
| COX-2, RMab | EP293 |
| Cytokeratin 8 & 18, RMab | B22.1/B23.1 |
| Cytokeratin 18, RMab | EP30 |
| Cytokeratin 19, MMab | BSB-34 |
| Cytokeratin 20, RMab | EP23 |
| Cytokeratin LMW CAM5.2, MMab | CAM5.2 |
| Gastrin, RPab | Polyclonal |
| GLUT1, RMab | RBT- GLUT1 |
| ICOS / CD278, RMab | RM417 |
| IgG4, MMab | BSB-96 |
| Lamin-B1, RMab | RBT-LMNB1 |
| MLH1, MMab | G168-728 |
| MMP-9, RMab | EP127 |
| MSH2, MMab | BSB-147 |
| MSH2, RMab | RBT-MSH2 |
| MSH6, MMab | 44 |
| MSH6, RMab | EP49 |
| MUC1, MMab | BSB-44 |
| MUC2, MMab | BSB-45 |
| MUC2, RMab | EP187 |
| MUC5AC, MMab | CLH2 |
| MUC6, MMab | CLH5 |

Colon & G.I. Antibodies

| | |
|---------------------------------|----------|
| NPM1/B23, MMab | BSB-124 |
| p14 ARF/CDKN2A, RMab | RBT-p14 |
| PMS2, RMab | RBT-PMS2 |
| PMS2, RMab | EP51 |
| SATB2, RMab | EP281 |
| SOX-9, RMab | EP317 |
| TIM-3 / HAVCR2 / CD366, MMab | BSB-163 |
| TFF3, MMab | BSB-181 |
| Villin, MMab | CWWB1 |
| Villin, RMab | EP163 |

**Insulin****Osteonectin /
SPARC****Gall Bladder & Pancreatic
Cancer Antibodies**

| | |
|-------------------------|-----------|
| DDR1, MMab | BSB-173 |
| Collagen IV, RMab | RBT-COL4 |
| Collagen IV, MMab | CIV22 |
| Glucagon, RMab | EP74 |
| IgG4, RMab | EP138 |
| IMP-3/1GF2BP3, RMab | EP286 |
| Insulin, MMab | BSB-42 |
| Insulin, RMab | EP125 |
| Islet-1/ISL1, RMab | EP283 |
| LRG1, MMab | BSB-174 |
| MUC4, MMab | 8G7 |
| MUC4, RMab | EP256 |
| Osteonectin/SPARC, MMab | BSB-93 |
| PAX-6, RMab | EP341 |
| PDX1, RMab | EP139 |
| S100A8, RMab | EP90 |
| S100P, RMab | RBT-S100P |
| S100P, RMab | EP186 |
| SMAD4/DPC4, MMab | BSB-63 |
| SMAD4/DPC4, RMab | RBT-SMAD4 |
| Somatostatin, RMab | EP130 |

Germ Cell Tumor Antibodies

| | |
|-------------------------------------|---------|
| Alpha-Fetoprotein, MMab | BSB-23 |
| LIN28, RMab | EP150 |
| MAGEC2, RMab | EP405 |
| Nanog, RMab | EP225 |
| p57, MMab | Kp10 |
| PLAP, MMab | BSB-47 |
| SALL4, RMab | EP299 |
| SF-1 / Steroidogenic Factor 1, MMab | BSB-149 |
| SOX-2, RMab | EP103 |

GIST Antibodies

| | |
|-------------|----------|
| CD117, RMab | EP10 |
| CD117, RMab | RM359 |
| DOG-1, RMab | RBT-DOG1 |
| DOG-1, RMab | EP332 |

Hematopoietic Antibodies


| | |
|--------------------------------|-------|
| CD41/Integrin alpha IIB, RMab | EP178 |
| Clusterin/Apoloprotein J, RMab | RM437 |
| Clusterin/Apoloprotein J, RMab | EP181 |
| Tryptase, MMab | G3 |
| Tryptase, RMab | EP259 |

Immunotherapy Antibodies


| | |
|------------------------------|------------|
| B7H3 / CD276, RMab | RBT-B7H3 |
| CTLA-4 / CD152, MMab | BSB-88 |
| CTLA-4 / CD152, RMab | RBT-CTLA4 |
| FOXP3, RPab | Polyclonal |
| ICOS / CD178, RMab | RM417 |
| IFN-a, MMab | BSB-158 |
| IFN-y, MMab | BSB-161 |
| LAG-3/CD223, RMAB | EP294 |
| OX-40/CD134, MMab | BSB-90 |
| PD-1/CD279, MMab | NAT-105 |
| PD-1/CD279, RMab | EP239 |
| PD-L1/CD274, RMab | RBT-PDL1 |
| PD-L1/CD274, RMab | 28-8 |
| TIGIT, MMab | BSB-152 |
| TIM-3 / HAVCR2 / CD366, MMab | BSB-163 |

Hodgkin's & NHB Lymphoma Antibodies

| | |
|-------------|----------|
| bcl-6, MMab | BSB-26 |
| bcl-6, RMab | EP278 |
| bcl-6, RMab | RBT-bcl6 |
| BOB.1, RMab | RBT-BOB1 |
| CD3, RMab | RBT-CD3 |
| CD10, RMab | EP195 |
| CD15, MMab | BSB-119 |
| CD19, MMab | BSB-97 |
| CD20, MMab | L26 |
| CD23, MMab | 1B12 |
| CD23, RMab | EP75 |
| CD30, MMab | Ber-H2 |
| CD57, MMab | BSB-10 |




bcl-6




CD10

Hodgkin's & NHB Lymphoma Antibodies

| | |
|-----------------------|------------|
| CD137 / TNFRSF9, MMab | BSB-159 |
| Fascin, MMab | BSB-36 |
| FOXO1, RMab | EP290 |
| FOXP3, RMab | EP340 |
| Granzyme B, RPab | Polyclonal |
| MUM-1, RMab | EP190 |
| OCT-2, RMab | EP115 |
| PAX-5, RMab | RBT-PAX5 |
| PD-L1/CD274, RMab | 28-8 |
| TIA-1, RMab | RBT-TIA1 |



SARS-CoV-2 & CD163 (Dual IHC)



HPV

Infectious Disease Antibodies


| | |
|---|----------------------|
| Aspergillus, RMab | RBT-A.fumi |
| ACE2, MMab | BSB-135 |
| Adenovirus, MMab | 20/11 & 2/6 |
| CD147, MMab | BSB-137 |
| Candida Albicans, RMab | RBT- C.alb |
| Cytomegalovirus, MMab | 8B1.2, 1G5.2 & 2D4.2 |
| Epstein Barr Virus, LMP-1, MMab | CS1-4 |
| Factor H / Complement Factor H, MMab | BSB-164 |
| Helicobacter pylori, MMab | BSB-37 |
| Helicobacter pylori, RPab | Polyclonal |
| Helicobacter pylori, RMab | EP279 |
| Hepatitis B Virus Core Antigen, RPab | Polyclonal |
| Hepatitis B Virus Surface Antigen, MMab | A10F1 |
| Herpes Simplex Virus 1, RPab | Polyclonal |
| Herpes Simplex Virus 2, RPab | Polyclonal |
| Herpes Simplex Virus I, MMab | 10A3 |
| Herpes Simplex Virus I & II, MMab | 10A3/BSB-116 |

Infectious Disease Antibodies


| | |
|-----------------------------------|--|
| Herpes Simplex Virus I & II, RPab | Polyclonal |
| Herpes Simplex Virus II, MMab | BSB-116 |
| HHV-8, MMab | 13B10 |
| HHV-8, RMab | RBT-HHV8 |
| HPV, MMab | BSB-66 |
| HPV16, MMab | CAMVIR-1 |
| IFN-a, MMab | BSB-158 |
| IFN-y, MMab | BSB-161 |
| Mycobacterium Tuberculosis | Polyclonal |
| Parvovirus, MMab | R92F6 |
| Pneumocystis Jirovecii, MMab | 3F6 |
| SARS-CoV-2, MMab | BSB-134 |
| SV40, MMab | Pab101 |
| TMPRSS2, MMab | BSB-136 |
| Toxoplasma gondii, RPab | Polyclonal |
| Treponema Pallidum, RPab | Polyclonal |
| Varicella Zoster Virus, MMab | SG1, SG1-1, SG2-2E6, SG3, SG4, NCP-1 & IE-62 |

Kidney Cancer & Urothelial Antibodies

| | |
|----------------------------|----------------|
| Amyloid A, RMab | EP335 |
| Cadherin-6, RMab | EP217 |
| Carbonic Anhydrase 9, RMab | EP161 |
| CD147, MMab | BSB-137 |
| c-Met/HGFR, RMab | RBT-c-Met/HGFR |
| FGFR-3, MMab | BSB-150 |
| Ksp-Cadherin, MMab | 4H6/F9 |
| Parvalbumin, RMab | EP300 |
| PAX-2, RMab | EP235 |
| PAX-8, RMab | ZR-1 |
| PLA2R1, MMab | BSB-129 |
| Renal Cell Carcinoma, MMab | PN-15 |
| S100A1, RMab | EP184 |
| S100A6, RMab | EP313 |
| SDHB, MMab | BSB-131 |
| Smoothelin, MMab | R4A |
| STAR, RMab | EP226 |
| Survivin, RMab | EP119 |
| TFE3, RMab | EP285 |
| Thrombomodulin, RMab | EP175 |
| Uroplakin III, RPab | Polyclonal |
| Uroplakin III, RMab | EP321 |
| WT1, MMab | 6F-H2 |



PAX-2




Renal Cell Carcinoma

Leukemia/Histocytic Antibodies


| | |
|-------------------------------|------------|
| A-1-Antichymotrypsin, RMab | EP384 |
| Alpha1-Antichymotrypsin, RPab | Polyclonal |
| Alpha-1-Antitrypsin, RPab | Polyclonal |
| Annexin A1, MMab | BSB-95 |
| c-Myc, RMab | RBT-CMYC |
| c-Myc, MMab | 9E10 |
| c-Myc, RMab | EP121 |
| Caspase-3, RMab | RM250 |
| CD1a, RMab | EP80 |
| CD5, RMab | RM314 |
| CD7, MMab | LP15 |
| CD8, RMab | EP334 |
| CD11b, RMab | EP45 |
| CD11c, RMab | EP157 |
| CD13, MMab | 38C12 |
| CD14, MMab | 7 |
| CD16, RMab | EP364 |
| CD25, MMab | 4C9 |
| CD33, RMab | RBT-CD33 |
| CD34, MMab | QEnd/10 |
| CD34, RMab | EP88 |
| CD38, MMab | SPC32 |
| CD42b, RMab | EP409 |
| CD61, MMab | 2f2 |
| CD61, RMab | EP65 |
| CD68, MMab | BSB-8 |
| CD68, MMab | KP1 |
| CD71, MMab | 10F11 |
| CD99, MMab | BSB-9 |
| CD103/ITGAE, RMab | EP206 |
| CD163, MMab | 10D6 |
| Factor XIIIa, RMab | EP292 |
| Glycophorin A, MMab | GA-R2 |
| Hemaglobin A, RMab | EP124 |
| IDH R132H, RMab | RBT-IDH1 |

Leukemia/Histocytic Antibodies


| | |
|------------------------------|------------|
| IgA, MMab | BSB-39 |
| IgG, MMab | BSB-40 |
| IgM, MMab | BSB-41 |
| Langerin, MMab | 12D6 |
| Langerin/CD207, RMab | EP349 |
| LEF1, RMab | RBT-LEF1 |
| LEF-1, RMab | EP310 |
| Lysozyme, MMab | Polyclonal |
| Lysozyme, RMab | EP134 |
| Macrophage/HAM-56, MMab | HAM-56 |
| MNDA, MMab | BSB-157 |
| Musashi 2, RMab | RM422 |
| Myeloperoxidase, RPab | Polyclonal |
| Myeloperoxidase, RMab | EP151 |
| Myeloperoxidase MPO, MMab | BSB-180 |
| NPM1/B23, MMab | BSB-124 |
| NPM1 / B23, MMab | BSB-124 |
| PDGFR-B, RMab | RBT-PDGFRB |
| PU.1, RMab | EP18 |
| Spectrin, MMab | RBC2/3D5 |
| T-Bet/TBX-21, RMab | EP263 |
| TCL-1, RMab | EP105 |
| TCR Alpha, MMab | BSB-126 |
| TCR Beta, MMab | BSB-117 |
| TCR Delta, MMab | BSB-127 |
| TdT, MMab | Polyclonal |
| TdT, RMab | RBT-TdT |
| TdT, RMab | EP266 |
| TIA-1, MMab | TIA-1 |
| TIM-3 / HAVCR2 / CD366, MMab | BSB-163 |
| TRAcP, MMab | 9C5 |
| Zap-70, MMab | 2F3.2 |




Arginase-1



Glypican-3



LEF-1



TCR Delta

Liver Cancer Antibodies

| | |
|---|------------|
| Alpha-Fetoprotein, RMab | EP209 |
| Arginase-1, RMab | EP261 |
| CEA, RPab | Polyclonal |
| Glutamine Synthase, MMab | GS-6 |
| Glypican-3, MMab | 1G12 |
| HEG1, MMab | SKM9-2 |
| Hepatocyte Specific Antigen / Hep Par 1, MMab | OCH1E5 |
| Prealbumin / Transthyretin, MMab | BSB-125 |

Lung Cancer Antibodies

| | |
|--|-------------|
| ALK, RMab | EP302 |
| BAP1, MMab | BSB-109 |
| BG8 Lewis Y, MMab | F3 |
| BRG-1 / SMARCA4, MMab | BSB-154 |
| c-Met / HGFR, RMab | EP1454Y |
| Calretinin, RMab | RBT-CALB2 |
| CD142 / TF / Coagulation Factor III . Thromboplastin, MMab | BSB-143 |
| CD73 / NT5E, RMab | RM431 |
| CXCL12 / SDF-1, MMab | BSB-165 |
| CXCR5 / CD185, RPab | Polyclonal |
| Cytokeratin 6, RMab | EP67 |
| Cytokeratin 7, MMab | OV-TL 12/30 |

Lymphoma Antibodies

| | |
|-----------------------|----------|
| ALK-1, RMab | RBT-ALK1 |
| bcl-2, MMab | BSB-5 |
| bcl-6, RMab | RBT-bcl6 |
| bcl-10, MMab | BSB-22 |
| c-Myc, RMab | EP121 |
| CD2, MMab | AB75 |
| CD3 Epsilon, RMab | RBT-CD3e |
| CD4, RMab | RBT-CD4 |
| CD5, RMab | RBT-CD5 |
| CD5, RMab | RM324 |
| CD6, MMab | BSB-54 |
| CD7, RMab | EP132 |
| CD8, MMab | C8/144B |
| CD10, MMab | 56C6 |
| CD14, RMab | EP128 |
| CD21, RMab | EP64 |
| CD25, RMab | RBT-CD25 |
| CD35, MMab | BSB-132 |
| CD38, RMab | EP135 |
| CD43, MMab | MT1 |
| CD45R, MMab | MB1 |
| CD45RO, MMab | UCHL-1 |
| CD74, MMab | N2 |
| CD75, MMab | LN1 |
| CD79a, MMab | JCB117 |
| CD123 / IL-3Ra, MMab | BSB-59 |
| CD137 / TNFSFR9, MMab | BSB-159 |

Lung Cancer Antibodies


| | |
|----------------------------|-----------|
| Cytokeratin LMW/ AE1, MMab | AE1 |
| Desmoglein-3, RMab | EP306 |
| EGFR, MMab | 31G7 |
| ERCC1, RMab | EP219 |
| GLUT-1, RMab | EP141 |
| INSM1, MMab | BSB-123 |
| INSM1, RMab | RBT-INSM1 |
| Lamin-B1, RMab | RBT-LMNB1 |
| Musashi 2, RMab | RM422 |
| Napsin A, MMab | BSB-112 |
| pan-TRK, RMab | RM423 |
| pan-TRK, RMab | RBT-TRK |
| PGP 9.5, MMab | BSB-46 |
| Pygopus 2 / Pygo 2, MMab | BSB-156 |

Lymphoma Antibodies


| | |
|---|--------------|
| CD138, MMab | B-A38 |
| CD138, RMab | EP201 |
| CXCR5 / CD185, RPab | Polyclonal |
| Cyclin D1, RMab | RBT-14 |
| Cyclin D1, RMab | RM241 |
| FOXP1, RMab | EP137 |
| Granzyme B, RMab | EP230 |
| hGAL / GCET, RMab | EP316 |
| Histone H3 Phospho, RPab | Polyclonal |
| Histone H3 Phospho, RMab | EP233 |
| IgD, MMab | Polyclonal |
| IgD, RMab | EP173 |
| Kappa Light Chains,MMab | BSB-58 |
| Lambda, MMab | BSB-16 |
| LMO2, RMab | RBT-LMO2 |
| MNDA, MMab | BSB-157 |
| NPM1 / B23, MMab | BSB-124 |
| PD-1 / CD279, MMab | NAT-105 |
| PD-1/ CD279, RMab | EP239 |
| PD-L1 / CD271, RMab | 28-8 |
| PU.1, RMab | EP18 |
| SOX-11, MMab | CL0142 |
| SOX-11, MMab | BSB-167 |
| TCR Alpha, MMab | BSB-126 |
| TCR Delta, MMab | BSB-127 |
| Tri-Methyl-Histone H3(Lys27)/H3K27Me3, RMab | RBT-H3K27Me3 |

Lung Cancer Antibodies


| | |
|-----------------------------------|-------------------|
| ROS-1, RMab | EP282 |
| S100A9, RMab | EP185 |
| Somatostaton Receptor 2, RMab | EP149 |
| SOX-2, RMab | EP103 |
| SOX-2, RMab | RM427 |
| Surfactant Protein D / SP-D, MMab | BSB-162 |
| Synaptophysin, RMab | RBT-Synaptophysin |
| Thymidylate Synthase / TS, MMab | BSB-160 |
| TTF-1, MMab | 8G7G3/1 |
| YAP1, MMab | BSB-146 |




ALK



Desmoglein



FOXO1



PD-L1

Melanoma & Skin Cancer Antibodies

| | |
|---------------------------------|----------------------|
| Adipophilin/ADRP, MMab | BSB-91 |
| CD4, MMab | BSB-179 |
| CD63, MMab | NKi/C3 |
| Cytokeratin 5 & 6, RMab | EP24 & EP67 |
| Cytokeratin 10, RMab | EP97 |
| Cytokeratin 14, RMab | RM328 |
| Cytokeratin MNF116, MMab | MNF116 |
| HSP70, RMab | RM432 |
| MART-1 / Melan-A, M2-7C10, MMab | M2-7C10 |
| MART-1 / Melan-A, A103, MMab | A103 |
| Melanoma / KBA.62, MMab | KBA6.2 |
| Melanoma / PNL2, MMab | PNL2 |
| Melanoma Cocktail, MMab | HMB-45, A103 & BSB-6 |
| Melanosome, MMab | HMB-45 |

Melanoma & Skin Cancer Antibodies


| | |
|-----------------------|-----------|
| MiTF, MMab | C5/D5 |
| MLH1, MMab | G168-728 |
| MTAP, RMab | RBT-MTAP |
| NRAS (Q61R), RMab | RBT-NRAS |
| PRAME, RMab | RBT-PRAME |
| S-100, MMab | 4C4.9 |
| S100B, RMab | EP32 |
| SOX-10, RMab | EP268 |
| SOX-10, MMab | BSB-62 |
| TDP-43 / TARDBP, MMab | BSB-166 |
| TIGIT, MMab | BSB-152 |
| TORC3/CRTC3, MMab | BSB-175 |
| Tyrosinase, MMab | BSB-6 |

Mesothelioma Antibodies


| | |
|---|----------|
| BAP1, MMab | BSB-109 |
| Calretinin, RMab | EP1798 |
| Calretinin, RMab | RM324 |
| Caveolin 1, RMab | EP353 |
| Claudin-4, RMab | EP417 |
| Cytokeratin 5, RMab | EP24 |
| Cytokeratin 5 & 6, RMab | RM341 |
| Cytokeratin 5 & 6, MMab | D5/16B4 |
| EpCAM / Epithelial Specific Antigen, MMab | Ber-EP4 |
| EpCAM / Epithelial Specific Antigen, MMab | MOC-31 |
| HEG1, RMab | SKM9-2 |
| Mesothelial Cell, MMab | HBME-1 |
| Mesothelin, RMab | EP140 |
| MTAP, RMab | RBT-MTAP |

Ovarian Cancer Antibodies


| | |
|-------------------------------------|----------|
| ARID1A, RMab | EP303 |
| CA-125, MMab | OC125 |
| hCG, MMab | BSB-38 |
| HE4, RMab | EP370 |
| hPL, RMab | EP241 |
| LIN28, RMab | EP150 |
| Nanog, RMab | EP225 |
| OCT4/POU5F1, RMab | RBT-OCT4 |
| SF-1 / Steroidogenic Factor 1, MMab | BSB-149 |
| STAR, RMab | EP226 |
| Surfactant Protein D / SP-D, MMab | BSB-162 |




PRAME



Melanosome



BAP1



MTAP

Pituitary Antibodies

| | |
|---------------------------|------------|
| ACTH, MMab | BSB-25 |
| FSH, RMab | EP257 |
| FSH, MMab | BSB-55 |
| GH, MMab | BSB-99 |
| GH, RPab | Polyclonal |
| HMG2A, RMab | EP398 |
| Luteinizing Hormone, MMab | BSB-53 |
| Prolactin, MMab | PRL02 |
| SDHB, MMab | BSB-131 |
| STAR, RMab | EP226 |
| TSH, MMab | BSB-56 |

Neural & Neuroendocrine Cancer Antibodies

| | |
|---------------------------------|------------|
| IDH1 R132H, RMab | RBT-IDH1 |
| INI-1, RMab | RBT-INI1 |
| INSM1, MMab | BSB-123 |
| INSM1, RMab | RBT-INSM1 |
| MGMT / AGAT, RMab | EP337 |
| Myelin Basic Protein, RPab | Polyclonal |
| Myelin Basic Protein, RMab | EP207 |
| NeuN, MMab | A60 |
| NeuN, RMab | RBT-NeuN |
| Neurofilament, MMab | 2F11 |
| Neurofilament, RMab | EP79 |
| NSE, MMab | BSB-94 |
| OLIG2, RMab | RBT-OLIG2 |
| OLIG2, RMab | EP112 |
| pan-TRK, RMab | RM423 |
| panTRK, RMab | RBT-TRK |
| PAX-7, RMab | BSB-145 |
| PHOX2B, RMab | EP312 |
| Prealbumin/ Transthyretin, MMab | BSB-125 |

Neural & Neuroendocrine Cancer Antibodies


| | |
|-----------------------------|----------------|
| Alpha Synuclein, MMab | BSB-114 |
| Amyloid Beta, RMab | RBT-A4 |
| ATRX, MMab | BSB-108 |
| ATRX, RMab | RBT-ATRX |
| Brachyury, RMab | RBT-TBXT |
| Calcitonin, RMab | RBT-Calcitonin |
| Calcitonin, RMab | EP92 |
| CD56, MMab | 123C3.D5 |
| Chromogranin A, MMab | LK2H10 |
| CXCR4 / CD184 / Fusin, RMab | EP394 |
| GFAP, MMab | G-A-5 |
| GFAP, RMab | RM246 |
| GH, RMab | EP267 |
| Glucagon, MMab | BSB-111 |
| IDH1 R132H, MMab | IHC132 |

Neural & Neuroendocrine Cancer Antibodies


| | |
|-------------------------------|------------|
| Prolactin, RMab | EP193 |
| Pygopus 2 / Pygo 2, MMab | BSB-156 |
| SDHB, MMab | BSB-131 |
| Somatostatin, MMab | BSB-113 |
| Somatostatin Receptor 2, RMab | EP149 |
| SOX-2, RMab | RM427 |
| Synaptophysin, RPab | Polyclonal |
| Synaptophysin, RMab | EP158 |
| Tau, MMab | BSB-115 |
| TSH, RMab | EP254 |
| ZFAND3, MMab | BSB-169 |

Prostate Cancer Antibodies

| | |
|--------------------------------|-----------------|
| ALDH1A1, RMab | EP168 |
| AMACR / P504S, RMab | 13H4 |
| Androgen Receptor,MMab | BSB-4 |
| Annexin VII, RMab | EP367 |
| Aurora B, RMab | RM278 |
| Caspase-3, RMab | RM250 |
| Cytokeratin 5, 6 & ERG, RMab | EP24/EP67/EP111 |
| Cytokeratin HMW / 34BE12, MMab | 34BE12 |
| FOXA1 / HNF-3A, RMab | EP277 |
| NKX3.1, RMab | EP356 |
| NKX3.1, RMab | RM430 |
| p40, RMab | ZR8 |
| p63, RMab | EP174 |
| Prostein / P501S, RMab | EP381 |
| PSA, MMab | BSB-7 |
| PSA, RMab | RBT-PSA |
| PSAP, MMab | PASE/4LJ |
| PSAP, RMab | EP53 |
| PSMA, RMab | EP192 |
| PSP94 / MSMB, RMab | EP203 |



PSA, RMab



CDK4, RMab

Sarcoma & Soft Tissue Antibodies

| | |
|------------------------------|--------------|
| Actin, Muscle Specific, MMab | HHF35 |
| BCOR, MMab | BSB-128 |
| CD13, RMab | EP117 |
| CD35, RMab | EP197 |
| CD99, MMab | BSB-9 |
| CDK4, RMab | EP180 |
| Desmin, MMab | D33 |
| Desmin, RMab | EP15 |
| Fumarate Hydratase, MMab | BSB-151 |
| HHV-8, MMab | 13B10 |
| HHV-8, RMab | RBT-HHV8 |
| INI-1, MMab | 25 |
| INI-1, RMab | RBT-INI1 |
| MDM2, RMab | RBT-MDM2 |
| MDM2, MMab | BSB-64 |
| MyoD1, RMab | RBT-MyoD1 |
| MyoD1, RMab | EP212 |
| Myogenin, MMab | F5D |
| Myoglobin, MMab | BSB-104 |
| Myoglobin, RMab | EP87 |
| NKX2.2, RMab | EP336 |
| NUT / NUTM1, RMab | RBT-NUTM1 |
| PAX-7, MMab | BSB-145 |
| SS18-SSX, RMab | RBT-SS18-SSX |
| STAT6, RMab | EP325 |
| TLE1, RMab | RBT-TLE1 |
| TLE-1, MMab | 1F5 |
| TLE-1, MMab | BSB-142 |
| Vimentin, RMab | EP21 |

Rejection & Autoimmunity Antibodies


| | |
|----------------|------------|
| C3d, RPab | Polyclonal |
| C4d, RMab | EP272 |
| IL-1a, MMab | BSB-138 |
| IL-1b, MMab | BSB-139 |
| IL-6, MMab | BSB-140 |
| iNOS, RMab | RBT-iNOS |
| Perforin, MMab | 5B10 |
| PLA2R1, MMab | BSB-129 |
| TNFa-IP2, MMab | BSB-141 |

Thyroid & Parathyroid Antibodies


| | |
|--------------------------|------------|
| Calcitonin, RPab | Polyclonal |
| Cytokeratin 19, RMab | RM364 |
| Galectin-3, MMab | 9C4 |
| HMGA2, RMab | EP398 |
| panTRK, RMab | RBT-TRK |
| Parafibromin, MMab | BSB-50 |
| Parathyroid (PTH), MMab | BSB-24 |
| Thyroglobulin, MMab | BSB-49 |
| Thyroglobulin, RMab | EP250 |
| Thyroid Peroxidase, RMab | EP159 |
| Trop-2 / EGP-1, MMab | BSB-148 |

Undifferentiated Tumor Antibodies


| | |
|---------------------------|---------------|
| CD45, MMab | 2B11 & PD7/26 |
| Cytokeratin AE1/AE3, MMab | AE1 & AE3 |
| PLAP, RMab | EP194 |
| S-100, MMab | 4C4.9 |
| SOX-10, MMab | EP298 |
| SOX, RMab | BSB-62 |
| Synaptophysin, RPab | Polyclonal |
| Vimentin, MMab | V9 |



STAT6



MDM2







Antibodies for Immunohistochemistry

Over 650 IVD Antibodies for Immunohistochemistry
Higher Affinity | Higher Sensitivity | Higher Specificity

298 Rabbit Monoclonal Antibodies
270 Mouse Monoclonal Antibodies
63 Rabbit Polyclonal Antibodies



Alpha-1-Antichymotrypsin, RMAb




IHC of Alpha-1-Antichymotrypsin on a FFPE Tonsil Tissue

Alpha-1-Antichymotrypsin Glycoprotein found in alpha (1)-globulin fraction in human serum. It inhibits chymotrypsin-like proteinases in vivo and has cytotoxic killer-cell activity in vitro. The protein also has a role as an acute-phase protein and is active in the control of immunologic and inflammatory processes, and as a tumor marker. It is a member of the serpin superfamily.

Alpha-1-Antichymotrypsin antibody reacts with histiocytes and histiocytic neoplasms. Its major application is defining the presence of Alpha-1-Antichymotrypsin in histiocytes and tumors derived from them. In eosinophilic granuloma and malignant histiocytosis, the reaction for this marker is heterogeneous in intensity and distribution. In fibrous histiocytomas, under certain circumstances, a diffuse homogeneous reaction may be observed.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP384
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Thymus, Breast
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Alpha-1-Antichymotrypsin, RPAb




IHC of A-1-Antichymotrypsin on a FFPE Tonsil Tissue

Alpha-1-Antichymotrypsin Glycoprotein found in alpha (1)-globulin fraction in human serum. It inhibits chymotrypsin-like proteinases in vivo and has cytotoxic killer-cell activity in vitro. The protein also has a role as an acute-phase protein and is active in the control of immunologic and inflammatory processes, and as a tumor marker. It is a member of the serpin superfamily.

Alpha-1-Antichymotrypsin antibody reacts with histiocytes and histiocytic neoplasms. Its major application is defining the presence of Alpha-1-Antichymotrypsin in histiocytes and tumors derived from them. In eosinophilic granuloma and malignant histiocytosis, the reaction for this marker is heterogeneous in intensity and distribution. In fibrous histiocytomas, under certain circumstances, a diffuse homogeneous reaction may be observed.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Thymus, Spleen, Liver, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Alpha-1-Antitrypsin, RPAb



IHC of A-1-Antitrypsin on a FFPE Tonsil Tissue

Alpha-1-Antitrypsin (A1AT) is a glycoprotein generally known as serum trypsin inhibitor. Alpha-1-Antitrypsin is also referred to as alpha-1 proteinase inhibitor (A1PI) because it is a serine protease inhibitor (serpin), inhibiting a wide variety of proteases. It protects tissues from enzymes of inflammatory cells, especially elastase, and has a reference range in blood of 1.5 -3.5 gram/liter (in the U.S. the reference range is generally expressed as mg/dL or micromoles), but the concentration can rise many fold upon acute inflammation. In its absence, elastase is free to break down elastin, which contributes to the elasticity of the lungs, resulting in respiratory complications such as emphysema, or COPD (Chronic Obstructive Pulmonary Disease) in adults and cirrhosis in adults or children.

Alpha-1-Antitrypsin is considered to be very useful in the study of inherited AAT deficiency, benign and Malignant Hepatic Tumors and Yolk-Sac Carcinomas. Positive staining for A-1-Antitrypsin may also be used in detection of benign and malignant lesions of a histiocytic nature. Sensitivity and specificity of the results have made this antibody a useful tool in the screening of patients with Cryptogenic Cirrhosis or other forms of liver disease with portal fibrosis of unknown etiology.


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3497 | Tinto Predilute | 3.0 ml |
| BSB 3498 | Tinto Predilute | 7.0 ml |
| BSB 3499 | Tinto Predilute | 15.0 ml |
| BSB 3500 | Concentrate | 0.1 ml |
| BSB 3501 | Concentrate | 0.5 ml |
| BSB 3502 | Concentrate | 1.0 ml |
| BSB 3503 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5001 | Tinto Predilute | 3.0 ml |
| BSB 5002 | Tinto Predilute | 7.0 ml |
| BSB 5003 | Tinto Predilute | 15.0 ml |
| BSB 5004 | Concentrate | 0.1 ml |
| BSB 5005 | Concentrate | 0.5 ml |
| BSB 5006 | Concentrate | 1.0 ml |
| BSB 5007 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5008 | Tinto Predilute | 3.0 ml |
| BSB 5009 | Tinto Predilute | 7.0 ml |
| BSB 5010 | Tinto Predilute | 15.0 ml |
| BSB 5011 | Concentrate | 0.1 ml |
| BSB 5012 | Concentrate | 0.5 ml |
| BSB 5013 | Concentrate | 1.0 ml |
| BSB 5014 | Control Slides | 5 |

ACE2, MAb



IHC of ACE2 on a FFPE Infected Testis Tissue

Angiotensin Converting Enzyme 2 is a dipeptidyl carboxydipeptidase active in the renin-angiotensin pathway, which helps regulate cardiovascular and renal functions. ACE2 is a secreted protein that cleaves angiotensin I (AngI) into angiotensin 1-9, and angiotensin II (AngII) into vasodilator angiotensin 1-7. ACE2 is expressed in the heart and kidney, where it may counteract vasoconstriction by inactivating Ang II. It is also found in the GI tract, lungs, and testis, in endothelial cells and less in vascular smooth muscle cells. ACE2 can contribute to tumor inhibition by inactivating AngII, which has been found to participate in tumor proliferation, angiogenesis, and metastasis.

ACE2 has shown inhibitory effects on lung, prostate, breast, and liver cancer through various signaling mechanisms, and its expression may be correlated with immune cell penetration into the tissue. ACE2 also functions as the main receptor for the spike glycoprotein of human coronaviruses HCoV-NL63, SARS-CoV and SARS-CoV-2.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-135


ISOTYPE: IgG1

CONTROL: Kidney, Testis, Brain, Colon, Fallopian tube

LOCALIZATION: Membranous, Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Rat

ACTH, MAb



IHC of ACTH on a FFPE Pituitary Tissue

Adrenocorticotropic Hormone (ACTH or corticotropin) is a polypeptide hormone synthesized from POMC, (Pro-opiomelanocortin) and secreted from corticotropes in the anterior lobe of the pituitary gland in response to Corticotropin-releasing Hormone (CRH) released by the hypothalamus. It consists of 39 amino acids.

ACTH is a useful marker in the classification of pituitary tumors and the study of pituitary disease. It reacts with ACTH-producing cells (corticotrophs), as well as other tumors (e.g., some Small-Cell Carcinomas present in lung tissue) causing Paraneoplastic Syndromes by secreting ACTH.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-25


ISOTYPE: IgG1/K

CONTROL: Normal Pituitary

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Actin Muscle Specific, MAb



IHC of Actin, Muscle Specific on a FFPE Skeletal Muscle Tissue

Actin is a globular-structural, 345 kDa protein that polymerizes in a helical fashion to form an actin filament (or microfilament). Actin filaments provide mechanical support for the cell, determine the cell shape, enable cell movements (through lamellipodia, filopodia, or pseudopodia); and participate in certain cell junctions, in cytoplasmic streaming and in contraction of the cell during cytokinesis. In muscle cells they play an essential role, along with myosin, in muscle contraction. In the cytosol, actin is predominantly bound to ATP, but can also bind to ADP.

This antibody recognizes actin of skeletal, cardiac, and smooth-muscle cells. It is not reactive with other mesenchymal cells except for myoepithelium. Muscle-Specific Actin recognizes alpha and gamma isotypes of all muscle groups. Non-muscle cells such as vascular endothelial cells and connective tissues are nonreactive. Neoplastic cells of non-muscle-derived tissue such as Carcinomas, Melanomas and Lymphomas are negative. This antibody is useful in the identification of rhabdoid cellular elements.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: HHf-35

ISOTYPE: IgG1/K

CONTROL: Skeletal Muscle, Appendix, Prostate

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Dog, Cat, Mouse, Rat, Monkey, Rabbit, Chicken, Shrew

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3702-3 | Tinto Predilute | 3.0 ml |
| BSB-3702-7 | Tinto Predilute | 7.0 ml |
| BSB-3702-15 | Tinto Predilute | 15.0 ml |
| BSB-3702-01 | Concentrate | 0.1 ml |
| BSB-3702-05 | Concentrate | 0.5 ml |
| BSB-3702-1 | Concentrate | 1.0 ml |
| BSB-3702-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5015 | Tinto Predilute | 3.0 ml |
| BSB 5016 | Tinto Predilute | 7.0 ml |
| BSB 5017 | Tinto Predilute | 15.0 ml |
| BSB 5018 | Concentrate | 0.1 ml |
| BSB 5019 | Concentrate | 0.5 ml |
| BSB 5020 | Concentrate | 1.0 ml |
| BSB 5021 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5022 | Tinto Predilute | 3.0 ml |
| BSB 5023 | Tinto Predilute | 7.0 ml |
| BSB 5024 | Tinto Predilute | 15.0 ml |
| BSB 5025 | Concentrate | 0.1 ml |
| BSB 5026 | Concentrate | 0.5 ml |
| BSB 5027 | Concentrate | 1.0 ml |
| BSB 5028 | Control Slides | 5 |

Actin Smooth Muscle, MAb




IHC of Actin, Smooth Muscle on a FFPE Appendix Tissue

Adenovirus, MAb



IHC of Adenovirus on a FFPE Infected Liver Tissue

Adipophilin/ADRP, MAb



IHC of Adipophilin on a FFPE Squamous Cell Carcinoma Tissue

Actin is a major component of the cytoskeleton and is present in every cell type. Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells. In vertebrates 3 main groups of actin isoforms (alpha, beta and gamma) have been identified. The alpha actins are found in muscle tissues and are a major constituent of the contractile apparatus. The beta and gamma actins coexist in most cell types as components of the cytoskeleton and as mediators of internal cell motility.

Smooth-Muscle Actin antibody does not stain cardiac or skeletal muscle; however, it will stain myofibroblasts and myoepithelial cells. This antibody could be used together with Muscle-Specific Actin to distinguish Leiomyosarcoma from Rhabdomyosarcoma. In most cases of Rhabdomyosarcoma, this antibody gives negative results whereas M.S. Actin is positive in the rhabdomyoblasts. Leiomyosarcomas are positive with both M.S. Actin and S. M. Actin antibodies.

Adenoviruses belong to the family Adenoviridae. They infect both humans and animals. Adenovirus was first isolated in human adenoids (tonsils), from which the name is derived. Adenoviruses are classified as group 1 under the Baltimore classification scheme. They are medium-sized (60-90 nm), non-enveloped icosahedral viruses containing double-stranded DNA.

Adipose differentiation-related protein, also known as perilipin 2 (PLIN2), ADRP or adipophilin, is a protein which in humans is encoded by the ADFP gene. Adipocyte differentiation-related protein is associated with the globule surface membrane material. This protein is a major constituent of the globule surface. Increase in mRNA levels is one of the earliest indications of adipocyte differentiation. Adipophilin occurs in a wide range of cultured cell lines, including fibroblasts and endothelial and epithelial cells. In tissues, however, expression of adipophilin is restricted to certain cell types, such as lactating mammary epithelial cells, adrenal cortex cells, Sertoli and Leydig cells of the male reproductive system, and steatosis or fatty change hepatocytes in alcoholic liver cirrhosis.

Adipophilin expression in various sebaceous lesions and other cutaneous tumors with a clear cell histology that may mimic sebaceous differentiation. Adipophilin can be valuable in an immunohistochemical panel when evaluating cutaneous lesions with clear cell histology as it identifies intracytoplasmic lipid vesicles in sebaceous and xanthomatous lesions. In periocular lesions, it is effective in helping to exclude basal cell carcinoma and squamous cell carcinoma when sebaceous carcinoma is under consideration. It is especially helpful in identifying intracytoplasmic lipid vesicles in poorly differentiated sebaceous carcinomas in challenging cases such as small periocular biopsy specimens.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-15

ISOTYPE: IgG2a/K

CONTROL: Appendix, Uterus

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Rat, Rabbit, Cat

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 20/11 & 2/6

ISOTYPE: IgG1/K

CONTROL: Adenovirus Infected Tissue

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-91

ISOTYPE: IgG1

CONTROL: Adrenal, SCC, TCC and Sebaceous Neoplasms

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5029 | Tinto Predilute | 3.0 ml |
| BSB 5030 | Tinto Predilute | 7.0 ml |
| BSB 5031 | Tinto Predilute | 15.0 ml |
| BSB 5032 | Concentrate | 0.1 ml |
| BSB 5033 | Concentrate | 0.5 ml |
| BSB 5034 | Concentrate | 1.0 ml |
| BSB 5035 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5036 | Tinto Predilute | 3.0 ml |
| BSB 5037 | Tinto Predilute | 7.0 ml |
| BSB 5038 | Tinto Predilute | 15.0 ml |
| BSB 5039 | Concentrate | 0.1 ml |
| BSB 5040 | Concentrate | 0.5 ml |
| BSB 5041 | Concentrate | 1.0 ml |
| BSB 5042 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3246 | Tinto Predilute | 3.0 ml |
| BSB 3247 | Tinto Predilute | 7.0 ml |
| BSB 3248 | Tinto Predilute | 15.0 ml |
| BSB 3249 | Concentrate | 0.1 ml |
| BSB 3250 | Concentrate | 0.5 ml |
| BSB 3251 | Concentrate | 1.0 ml |
| BSB 3252 | Control Slides | 5 |

Albumin, RPab



IHC and IF of Albumin on a FFPE Lupus Erythematosus (IHC) and a FFPE Kidney Tissue (IF)

The albumins are a family of globular proteins, the most common of which are the serum albumins. Albumins are commonly found in blood plasma and differ from other blood proteins in that they are not glycosylated. Albumin functions primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume. Mutations in this gene on chromosome 4 result in various anomalous proteins.

Low albumin (hypoalbuminemia) may be caused by liver disease, nephrotic syndrome, burns, protein-losing enteropathy, malabsorption, malnutrition, late pregnancy, artefact, genetic variations and malignancy. High albumin (hyperalbuminemia) is almost always caused by dehydration. In some cases of retinol (Vitamin A) deficiency, the albumin level can be elevated to high-normal values.

It has been reported in systemic lupus erythematosus (SLE) patients an increased prevalence of IgG autoantibodies against human serum albumin (anti-HSA IgG) that are associated with SLE disease activity.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Salivary Gland, Kidney, Tonsil, Lupus Erythematosus

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ALDH1A1, RMAb



IHC of ALDH1A1 on a FFPE Pancreas Tissue

Aldehyde dehydrogenase 1 family member A1, also known as ALDH1A1 is an enzyme that in humans is encoded by the Aldh1A1 gene on chromosome 9. ALDH1A1 belongs to the aldehyde dehydrogenase family, which plays a role in alcohol metabolism. There are two major aldehyde dehydrogenase isozymes in the liver; cytosolic and mitochondrial. The ALDH1A1 gene encodes the cytosolic isozyme, and catalyzes the oxidation of retinaldehyde to retinoic acid. Studies have indicated that ALDH1A1 may also be involved in the regulation of the metabolic responses to high-fat diet.

ALDH1A1 has been a well established marker of hematopoietic stem cells and progenitor cells. Recent studies also show that ALDH1A1 is an important cancer stem marker associated with tumor progression in cancers of the breast, prostate and lung. This antibody labels epithelial cells of the stomach, liver, kidney and thyroid, neural cells and stromal cells including endothelial cells. In tumors, it stains stromal cells as well as tumor cells in many types of cancers. ALDH1A1 can be used with CD34 to aid in the differentiation between solitary fibrous tumors, hemangiopericytoma, meningioma and synovial sarcomas.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP168


ISOTYPE: IgG

CONTROL: Kidney, Liver, Testis, Colon Cancer, Breast Cancer

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Predicted: Mouse, Rat

ALK-1/CD246, RMAb



IHC of ALK-1/CD246 on a FFPE Anaplastic Large Cell Lymphoma Tissue

Anaplastic Lymphoma Kinase (ALK) was originally discovered as a NPM (Nucleophosmin)-ALK fusion protein. The ALK gene is on chromosome 2. Upon translocation between chromosome 2 and chromosome 5 t(2;5), the ALK gene fuses with the NPM gene. The chimeric product (NPM ALK) resulting from t(2;5) translocation is a protein of 80 kDa with the N terminal portion of NPM linked to the complete intracellular portion of ALK.

This antibody recognizes a human p80 protein, identified as a hybrid of the Anaplastic Lymphoma Kinase (ALK) gene and the Nucleophosmin (NPM) gene resulting from the t(2;5)(p23;q35) translocation found in a third of Large-Cell Lymphomas. ALK-1 is detected in 60% of Anaplastic Large-Cell Lymphomas and has proven to indicate a better prognosis in the ALK-1 (+) group.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-ALK-1

ISOTYPE: IgG


CONTROL: Anaplastic Large Cell Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 3012 | Tinto Predilute | 3.0 ml | BSB 2440 | Tinto Predilute | 3.0 ml | BSB 5043 | Tinto Predilute | 3.0 ml |
| BSB 3013 | Tinto Predilute | 7.0 ml | BSB 2441 | Tinto Predilute | 7.0 ml | BSB 5044 | Tinto Predilute | 7.0 ml |
| BSB 3014 | Tinto Predilute | 15.0 ml | BSB 2442 | Tinto Predilute | 15.0 ml | BSB 5045 | Tinto Predilute | 15.0 ml |
| BSB 3015 | Concentrate | 0.1 ml | BSB 2443 | Concentrate | 0.1 ml | BSB 5046 | Concentrate | 0.1 ml |
| BSB 3016 | Concentrate | 0.5 ml | BSB 2444 | Concentrate | 0.5 ml | BSB 5047 | Concentrate | 0.5 ml |
| BSB 3017 | Concentrate | 1.0 ml | BSB 2445 | Concentrate | 1.0 ml | BSB 5048 | Concentrate | 1.0 ml |
| BSB 3018 | Control Slides | 5 | BSB 2446 | Control Slides | 5 | BSB 5049 | Control Slides | 5 |

ALK-1/CD246, RMAb




IHC of ALK on a FFPE Anaplastic Large Cell Lymphoma Tissue

Anaplastic Lymphoma Kinase (ALK) was originally discovered as a NPM (Nucleophosmin)-ALK fusion protein. The ALK gene is on chromosome 2. Upon translocation between chromosome 2 and chromosome 5 t(2;5), the ALK gene fuses with the NPM gene. The chimeric product (NPM ALK) resulting from t(2;5) translocation is a protein of 80 kDa with the N terminal portion of NPM linked to the complete intracellular portion of ALK.

This antibody recognizes a human p80 protein, identified as a hybrid of the Anaplastic Lymphoma Kinase (ALK) gene and the Nucleophosmin (NPM) gene resulting from the t(2;5)(p23;q35) translocation found in a third of Large-Cell Lymphomas. ALK-1 is detected in 60% of Anaplastic Large-Cell Lymphomas and has proven to indicate a better prognosis in the ALK-1 (+) group.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP302
ISOTYPE: IgG
CONTROL: Anaplastic Large Cell Lymphoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Predicted: Mouse, Rat

Alpha Synuclein, MMAb



IHC of Alpha Synuclein on a FFPE Astrocytoma Tissue


Alpha-synuclein is a 140 amino acids protein encoded by the SNCA gene. It is predominantly expressed in the neocortex, hippocampus, substantia nigra, thalamus, and cerebellum with smaller amounts are found in the heart, muscles, and other tissues. In the brain, alpha-synuclein is found mainly in presynaptic terminals.

An alpha-synuclein fragment, known as the non-Abeta component (NAC) of Alzheimer's disease amyloid, originally found in an amyloid-enriched fraction, was shown to be a fragment of its precursor protein, NACP. Alpha-synuclein aggregates to form insoluble fibrils in pathological conditions characterized by Lewy bodies, such as Parkinson's disease, dementia with Lewy bodies and multiple system atrophy. These disorders are known as synucleinopathies. Occasionally, Lewy bodies contain tau protein; however, alpha-synuclein and tau constitute two distinctive subsets of filaments in the same inclusion bodies.

Alpha-synuclein pathology is also found in both sporadic and familial cases with Alzheimer's disease. In rare cases of familial forms of Parkinson's disease, there is a mutation in the gene coding for alpha-synuclein. Genomic duplication and triplication of the gene appear to be a rare cause of Parkinson's disease in other lineages, although more common than point mutations. Hence certain mutations of alpha-synuclein may cause it to form amyloid-like fibrils that contribute to Parkinson's disease.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-114
ISOTYPE: IgG2a/K
CONTROL: Brain, Breast, Skin, Tonsil, Bone Marrow, Alzheimer's Disease, Parkinson's Disease
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

Alpha-Fetoprotein, MMAb



IHC of Alpha-Fetoprotein on a FFPE Liver Tissue

Alpha-fetoprotein (AFP) is a protein which in humans is encoded by the AFP gene. This gene encodes alpha-fetoprotein, a major plasma protein produced by the yolk sac and the liver during fetal life. This protein is thought to be the fetal counterpart of serum albumin, and the alpha-fetoprotein and albumin genes are present in tandem on chromosome 4.

Positive staining with this antibody is seen in hepatocytes of fetal liver and hepatoma. Since only traces of AFP are found in adult serum, elevated levels suggest either a benign or malignant lesion of the liver, a Yolk-Sac Carcinoma, or one of a few other tumors. In conjunction with elevated serum levels, AFP has been immunohistochemically demonstrated in Yolk-Sac Carcinomas in gonadal and extragonadal sites of hepatic malignancies and a few other neoplasms.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-23
ISOTYPE: IgG2a/K
CONTROL: Fetal Liver, Hepatocellular Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Canine

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2796 | Tinto Predilute | 3.0 ml |
| BSB 2797 | Tinto Predilute | 7.0 ml |
| BSB 2698 | Tinto Predilute | 15.0 ml |
| BSB 2799 | Concentrate | 0.1 ml |
| BSB 2800 | Concentrate | 0.5 ml |
| BSB 2801 | Concentrate | 1.0 ml |
| BSB 2802 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3286 | Tinto Predilute | 3.0 ml |
| BSB 3287 | Tinto Predilute | 7.0 ml |
| BSB 3288 | Tinto Predilute | 15.0 ml |
| BSB 3289 | Concentrate | 0.1 ml |
| BSB 3290 | Concentrate | 0.5 ml |
| BSB 3291 | Concentrate | 1.0 ml |
| BSB 3292 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5050 | Tinto Predilute | 3.0 ml |
| BSB 5051 | Tinto Predilute | 7.0 ml |
| BSB 5052 | Tinto Predilute | 15.0 ml |
| BSB 5053 | Concentrate | 0.1 ml |
| BSB 5054 | Concentrate | 0.5 ml |
| BSB 5055 | Concentrate | 1.0 ml |
| BSB 5056 | Control Slides | 5 |

Alpha-Fetoprotein, RMAb




IHC of Alpha-Fetoprotein on a FFPE Yolk-Sac Carcinoma Tissue

Alpha-fetoprotein (AFP) is a protein which in humans is encoded by the AFP gene. This gene encodes alpha-fetoprotein, a major plasma protein produced by the yolk sac and the liver during fetal life. This protein is thought to be the fetal counterpart of serum albumin, and the alpha-fetoprotein and albumin genes are present in tandem on chromosome 4.

Positive staining with this antibody is seen in hepatocytes of fetal liver and hepatoma. Since only traces of AFP are found in adult serum, elevated levels suggest either a benign or malignant lesion of the liver, a Yolk-Sac Carcinoma, or one of a few other tumors. In conjunction with elevated serum levels, AFP has been immunohistochemically demonstrated in Yolk-Sac Carcinomas in gonadal and extragonadal sites of hepatic malignancies and a few other neoplasms.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP209
ISOTYPE: IgG
CONTROL: Fetal Liver, Hepatocellular Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

AMACR/P504S, RMAb




IHC of AMACR/P504S on a FFPE Prostatic Adenocarcinoma Tissue

AMACR (P504S) is an acronym for the protein alpha-methylacyl CoA racemase that helps to metabolize certain fatty acids within the body. AMACR has been recently described as a prostate cancer-specific gene that encodes a protein involved in the beta-oxidation of branched chain fatty acids. Expression of AMACR protein is found in Prostatic Adenocarcinoma but not in benign prostatic tissue. It stains premalignant lesions of the prostate: High-Grade Prostatic Intraepithelial Neoplasia (PIN) and Atypical Adenomatous Hyperplasia. Several studies have suggested that AMACR can be used as a prostate cancer biomarker.

High expression of AMACR (P504S) protein is usually found in Prostatic Adenocarcinoma but not in benign prostatic tissue by immunohistochemical staining in paraffin-embedded tissues. Using AMACR as a positive marker along with basal-cell staining (34βE12 or p63) as a negative marker could help to confirm the diagnosis of small foci of Prostate Carcinoma on needle biopsies.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: 13H4
ISOTYPE: IgG
CONTROL: Kidney, Liver, Salivary Gland, Prostate Lesions, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Amyloid A, RMAb



IHC and IF of Amyloid A on a FFPE Kidney Tissue


Serum amyloid A (SAA) proteins are a family of apolipoproteins associated with high-density lipoprotein (HDL) in plasma. Different isoforms of SAA are expressed constitutively (constitutive SAAs) at different levels or in response to inflammatory stimuli (acute phase SAAs). These proteins are produced predominantly by the liver. The conservation of these proteins throughout invertebrates and vertebrates suggests that SAAs play a highly essential role in all animals. Acute-phase serum amyloid A proteins (A-SAAs) are secreted during the acute phase of inflammation. A-SAAs are implicated in several chronic inflammatory diseases, such as amyloidosis, atherosclerosis, and rheumatoid arthritis.

Amyloidosis is a disease characterized by the abnormal build-up of amyloid, abnormal non-branching fibrillary β-pleated sheet proteins that are insoluble and highly resistant to proteolytic degradation that result in localized or systemic organ dysfunction. AA amyloidosis is associated with a variety of chronic inflammatory conditions and infections, derived from SAA. Immunohistochemical staining using a panel of antibodies including κ and λ Ig light chains, Amyloid A, and Transthyretin can aid in recognizing most forms of amyloid. The Amyloid A immunostaining detects tissue deposition of serum Amyloid A protein, an acute phase reactive protein. It is positive in AA Amyloidosis and familial Mediterranean fever. SAA concentrations have been reported to be a marker of poor prognosis, elevated in patients with advanced stages of cancer and those with malignant disease.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP335
ISOTYPE: IgG
CONTROL: Kidney, Amyloidosis
LOCALIZATION: Cytoplasmic, Extracellular
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2384 | Tinto Predilute | 3.0 ml | BSB 5057 | Tinto Predilute | 3.0 ml | BSB 2803 | Tinto Predilute | 3.0 ml |
| BSB 2385 | Tinto Predilute | 7.0 ml | BSB 5058 | Tinto Predilute | 7.0 ml | BSB 2804 | Tinto Predilute | 7.0 ml |
| BSB 2386 | Tinto Predilute | 15.0 ml | BSB 5059 | Tinto Predilute | 15.0 ml | BSB 2805 | Tinto Predilute | 15.0 ml |
| BSB 2387 | Concentrate | 0.1 ml | BSB 5060 | Concentrate | 0.1 ml | BSB 2806 | Concentrate | 0.1 ml |
| BSB 2388 | Concentrate | 0.5 ml | BSB 5061 | Concentrate | 0.5 ml | BSB 2807 | Concentrate | 0.5 ml |
| BSB 2389 | Concentrate | 1.0 ml | BSB 5062 | Concentrate | 1.0 ml | BSB 2808 | Concentrate | 1.0 ml |
| BSB 2390 | Control Slides | 5 | BSB 5063 | Control Slides | 5 | BSB 2809 | Control Slides | 5 |

Amyloid Beta, RMAb




IHC of Amyloid Beta on a FFPE Astrocytoma Tissue

Brain Amyloid beta is elevated in patients with sporadic Alzheimer's disease and is the main component of amyloid plaques. Similar plaques appear in some variants of Lewy body dementia and in inclusion body myositis, while Amyloid beta can also form the aggregates that coat cerebral blood vessels in cerebral amyloid angiopathy. The plaques are composed of a tangle of regularly ordered fibrillar aggregates called amyloid fibers, a protein fold shared by other peptides such as the prions associated with protein misfolding diseases.

The mechanism by which Amyloid beta may damage and kill neurons is by generating reactive oxygen species during the process of its self-aggregation. It has been reported that amyloid beta production follows a circadian rhythm, rising when an animal or a person is awake and falling during sleep. The wakefulness-promoting neuroprotein orexin has been shown to be necessary for the circadian rhythm of amyloid beta production. This is consistent with recent findings that chronic sleep deprivation is associated with early onset Alzheimer's disease.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-A4
ISOTYPE: IgG
CONTROL: Testis, Kidney, Pancreas, Salivary Gland, Alzheimer's Disease
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

Androgen Receptor, MMAb




IHC of Androgen Receptor on a FFPE Prostate Tissue

The androgen receptor (AR) is a type of nuclear receptor which is activated by binding of either of the androgenic hormones testosterone or dihydrotestosterone. The main function of the androgen receptor is as a DNA-binding transcription factor which regulates gene expression. However, the androgen receptor has additional functions independent of DNA binding. The AR signaling pathway plays a key role in development and function of male reproductive organs, including the prostate and epididymis. AR also plays a role in nonreproductive organs, such as muscle, hair follicles, and the brain.

This antibody reacts with the androgen receptor and also with the newly-described A form of the receptor. This antibody does not cross-react with estrogen, progesterone or glucocorticoid receptors. Abnormalities in the AR-signaling pathway have been linked to a number of diseases, including Prostate Cancer, Kennedy's Disease and male infertility.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-4
ISOTYPE: IgG1
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Annexin A1, MMAb



IHC of Annexin A1 on a FFPE Hairy Cell Leukemia Tissue

The protein Annexin A1 is encoded by the ANXA1 gene, which is upregulated in hairy cell leukemia. The NF-κB signal transduction pathway is exploited by cancerous cells to proliferate and avoid apoptosis. Annexin A1 inhibits that pathway by binding to the p65 subunit, thus making Annexin A1 of particular interest for use as a potential anti-cancer drug. It may also contain tumor suppressive and protective characteristics, which have been evidenced by its ability to protect against DNA damage induced by heat in breast cancer cells.

Annexin A1 is strongly expressed on the cell membrane and occasionally in the cytoplasm of tumor cells in 97% of samples from patients with hairy cell leukemia. By contrast, B-cell lymphomas other than hairy cell leukemia are ANXA1 negative. Thus, ANXA1 is a molecule specific to hairy cell leukemia that can be used to differentiate this disease from other B-cell lymphomas.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-95
ISOTYPE: IgG2b
CONTROL: Liver, Tonsil, Spleen, Thymus, Lung, Colon, Hairy Cell Leukemia
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Dog

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3441 | Tinto Predilute | 3.0 ml |
| BSB 3442 | Tinto Predilute | 7.0 ml |
| BSB 3443 | Tinto Predilute | 15.0 ml |
| BSB 3444 | Concentrate | 0.1 ml |
| BSB 3445 | Concentrate | 0.5 ml |
| BSB 3446 | Concentrate | 1.0 ml |
| BSB 3447 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6071 | Tinto Predilute | 3.0 ml |
| BSB 6072 | Tinto Predilute | 7.0 ml |
| BSB 6073 | Tinto Predilute | 15.0 ml |
| BSB 6074 | Concentrate | 0.1 ml |
| BSB 6075 | Concentrate | 0.5 ml |
| BSB 6076 | Concentrate | 1.0 ml |
| BSB 6077 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6359 | Tinto Predilute | 3.0 ml |
| BSB 6360 | Tinto Predilute | 7.0 ml |
| BSB 6361 | Tinto Predilute | 15.0 ml |
| BSB 6362 | Concentrate | 0.1 ml |
| BSB 6363 | Concentrate | 0.5 ml |
| BSB 6364 | Concentrate | 1.0 ml |
| BSB 6365 | Control Slides | 5 |

Annexin VII, RMAb



IHC of Annexin VII on a FFPE Prostate Adenocarcinoma Tissue

Annexin VII (or Annexin A7) is a Ca²⁺-dependent and phospholipid-binding protein. It is encoded by the ANX7 gene, which is located in human chromosome 10q21. It is a member of the annexin protein family, which regulates various endomembrane processes including vesicle fusion, segregation, and compartmentalization and is implicated in plasma membrane repair mechanisms.

Studies suggest that ANX7 is a tumor suppressor gene and it was found that Annexin VII is associated with several types of cancers, such as Prostate, Breast, Liver, and Gastric cancer. IHC analysis showed that a low or the loss of expression of ANX7 correlates with Prostate cancer progression. On the contrary, high levels of ANX7 expression were found in Metastatic Breast Tumors and that HER2-Negative patients suffer increased risk of death when ANX7 expression levels are elevated. Another study found that ANX7 predominantly shows low expression in other cancer types, such as Colon Adenocarcinoma or Bladder Transitional Cell Carcinoma. A study identified the prognostic impact of ANXA7 in Prostate Cancer using tissue microarrays, identified ANXA7 as a new prognostic factor and indicated a bimodal correlation to tumor progression.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP367


ISOTYPE: IgG

CONTROL: Breast, Colon, Kidney, Tonsil, Prostate, Testis, Transitional Cell Carcinoma

LOCALIZATION: Nuclear, Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

Arginase-1, RMAb



IHC of Arginase-1 on a FFPE Hepatocellular Carcinoma Tissue

Arginase is the catalyst for the fifth and final step in the urea cycle, which is a series of biochemical reactions in mammals during which the body disposes of harmful ammonia. Arginase works to convert L-arginine into L-ornithine and urea. Arginase-1 is located primarily in the cytoplasm of the liver. Arginase consists of three tetramers, and the enzyme requires a two-molecule metal cluster of manganese in order to maintain proper function. These Mn²⁺ ions coordinate with water, orienting and stabilizing the molecule and allowing water to act as a nucleophile and attack L-arginine, hydrolyzing it into ornithine and urea.

Arginase-1 is abundantly expressed in the liver and it represents a sensitive and specific marker of benign and malignant hepatocytes. In sections of normal liver, anti-Arginase-1 produces strong, diffuse cytoplasmic reactivity in all hepatocytes throughout the lobule. In a small percentage of cases, patchy nuclear reactivity is also evident in hepatocytes along with the strong cytoplasmic reactivity. Hepatocellular carcinoma usually shows higher protein expression of ARG1 than normal liver cells.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP261


ISOTYPE: IgG

CONTROL: Liver, Hepatocellular Carcinoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ARID1A, RMAb



IHC of ARID1A on a FFPE Colon Adenocarcinoma Tissue

Genes encoding subunits of SWI/SNF chromatin remodeling complexes are collectively mutated in 20% of all human cancers. ARID1A is the SWI/SNF subunit gene that is most frequently mutated, at variable frequencies across molecular and histological subtypes of cancer.

ARID1A is a tumour suppressor gene frequently mutated in Clear Cell and Endometrioid Carcinomas of the Ovary and Endometrium. Loss of ARID1A function as shown by loss of expression, presumably due to mutations, is an early molecular event, occurring before malignant transformation, in the development of the majority of Ovarian Clear Cell and Endometrioid Carcinomas arising in Endometriomas. Mutations were identified in 6/17 (35%) Ovarian Clear Cell Carcinomas, 5/8 (63%) Ovarian Endometrioid Carcinomas, 2/5 (40%) Endometrial Carcinomas, and 1/7 (14%) Carcinosarcomas. Some studies have demonstrated that ARID1A has a critical tumor suppressor role in the Colon, and that its inactivation leads to the development of Colon Cancers via a mechanism that is distinct from previously established genetic models. ARID1A inactivating mutations are present at a high frequency in advanced endocrine-resistant ER+ Breast Cancer. ARID1A may play an important role in and serves as a valuable prognostic marker in Gastric Cancer.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP303

ISOTYPE: IgG

CONTROL: Breast, Fallopian Tube, Prostate, Testis, Transitional Cell Carcinoma,


T Cell Lymphoblastic Lymphoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|-------------|-----------------|---------|
| BSB-3709-3 | Tinto Predilute | 3.0 ml | BSB 2447 | Tinto Predilute | 3.0 ml | BSB-3749-3 | Tinto Predilute | 3.0 ml |
| BSB-3709-7 | Tinto Predilute | 7.0 ml | BSB 2448 | Tinto Predilute | 7.0 ml | BSB-3749-7 | Tinto Predilute | 7.0 ml |
| BSB-3709-15 | Tinto Predilute | 15.0 ml | BSB 2449 | Tinto Predilute | 15.0 ml | BSB-3749-15 | Tinto Predilute | 15.0 ml |
| BSB-3709-01 | Concentrate | 0.1 ml | BSB 2450 | Concentrate | 0.1 ml | BSB-3749-01 | Concentrate | 0.1 ml |
| BSB-3709-05 | Concentrate | 0.5 ml | BSB 2451 | Concentrate | 0.5 ml | BSB-3749-05 | Concentrate | 0.5 ml |
| BSB-3709-1 | Concentrate | 1.0 ml | BSB 2452 | Concentrate | 1.0 ml | BSB-3749-1 | Concentrate | 1.0 ml |
| BSB-3709-CS | Control Slides | 5 | BSB 2453 | Control Slides | 5 | BSB-3749-CS | Control Slides | 5 |

ATM, RMAb




IHC of ATM on a FFPE Chronic Lymphocytic Lymphoma tissue

ATRX, MMAb



IHC of ATRX on a FFPE Astrocytoma Tissue

ATRX, RMAb



IHC of ATRX on a FFPE Astrocytoma tissue

Ataxia telangiectasia mutated (ATM) protein is a serine/threonine kinase that belongs to the family of phosphatidylinositol-3 kinase (PI3K)-related protein kinases. ATM is a key checkpoint protein of the DNA damage response pathway. When DNA double-strand breaks occur, ATM activates different cascades, resulting in activation of cell cycle checkpoints, cell cycle arrests, and apoptosis.

Mutations in the ATM gene lead to an increased risk of several cancer types, which include Leukemias, Lymphomas, Colorectal Cancer, Pancreatic Cancer and Adenocarcinoma of the Stomach. Studies show that Gastric Cancer tissues were found to be ATM-negative by IHC staining, indicating that ATM can be used as a biomarker. It has also been found that low expression of the ATM protein contributes to more aggressive progression and poor clinical outcome of Breast Cancer. Results of another study indicate loss of ATM protein expression is associated with development of Lung Cancer.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP327
ISOTYPE: IgG
CONTROL: Fallopian Tube, Brain, Colon, Breast, Testis, Tonsil, Transitional Cell Carcinoma, T Cell Lymphoblastic
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Mutation/loss of ATRX expression has been described in anaplastic gliomas. A study explored the role of ATRX status in the molecular classification of anaplastic gliomas and its impact on survival. Loss of ATRX expression was detected in 45 % of anaplastic astrocytomas (AA), 27 % of anaplastic oligoastrocytomas (AOA) and 10 % of anaplastic oligodendrogliomas (AO). Survival analysis showed a marked separation of IDH mutant astrocytic tumors into two groups based on ATRX status: tumors with ATRX loss had a significantly better prognosis.

Another study conducted an immunohistochemical analysis of ATRX expression in adult diffuse gliomas and found ATRX immunoreactivity of tumor cells was either almost totally absent or completely retained in all cases. There was perfect concordance between the IHC results and ATRX mutation status. ATRX loss was observed in 54.5, 30.8 and 0.0% of grades II/III astrocytomas, oligoastrocytomas and oligodendrogliomas, respectively, and 12.7% of glioblastomas. Another recent study analyzed the use of ATRX, IDH and 1p/19q codeletion in a series of astrocytomas, oligodendrogliomas, oligoastrocytomas and glioblastomas and presented an algorithm based on stepwise analysis with initial immunohistochemistry for ATRX and IDH1-R132H followed by 1p/19q analysis then by IDH sequencing, which reduces the number of molecular analyses and which has a far better association with patient outcome.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-108
ISOTYPE: IgG2a/K
CONTROL: Breast, Tonsil, Testis, Salivary Gland, Placenta, TCC, Astrocytoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

α -thalassemia/mental retardation syndrome X-linked (ATRX) gene is located on chromosome Xq21.1. Germline mutations of ATRX have been found to cause the complex genetic disorder called Alpha-Thalassemia mental retardation syndrome. Somatic mutations, deletions, and altered ATRX expression levels were found to be prevalent in several cancer types. A study reported the loss of ATRX expression was found to be a prognostic marker for chromosome instability in pancreatic neuroendocrine tumors.

Mutation/loss of ATRX expression has been described in anaplastic gliomas. A study explored the role of ATRX status in the molecular classification of anaplastic gliomas and its impact on survival. Survival analysis showed a marked separation of IDH mutant astrocytic tumors into two groups based on ATRX status: tumors with ATRX loss had a significantly better prognosis. Another recent study analyzed the use of ATRX, IDH and 1p/19q codeletion in a series of astrocytomas, oligodendrogliomas, oligoastrocytomas and glioblastomas and presented an algorithm based on stepwise analysis with initial immunohistochemistry for ATRX and IDH1-R132H followed by 1p/19q analysis, then by IDH sequencing, which reduces the number of molecular analyses and has a far better association with patient outcome.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-ATRX
ISOTYPE: IgG
CONTROL: Breast, Colon, Fallopian Tube, Brain, Tonsil, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3710-3 | Tinto Predilute | 3.0 ml |
| BSB-3710-7 | Tinto Predilute | 7.0 ml |
| BSB-3710-15 | Tinto Predilute | 15.0 ml |
| BSB-3710-01 | Concentrate | 0.1 ml |
| BSB-3710-05 | Concentrate | 0.5 ml |
| BSB-3710-1 | Concentrate | 1.0 ml |
| BSB-3710-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3293 | Tinto Predilute | 3.0 ml |
| BSB 3294 | Tinto Predilute | 7.0 ml |
| BSB 3295 | Tinto Predilute | 15.0 ml |
| BSB 3296 | Concentrate | 0.1 ml |
| BSB 3297 | Concentrate | 0.5 ml |
| BSB 3298 | Concentrate | 1.0 ml |
| BSB 3299 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3711-3 | Tinto Predilute | 3.0 ml |
| BSB-3711-7 | Tinto Predilute | 7.0 ml |
| BSB-3711-15 | Tinto Predilute | 15.0 ml |
| BSB-3711-01 | Concentrate | 0.1 ml |
| BSB-3711-05 | Concentrate | 0.5 ml |
| BSB-3711-1 | Concentrate | 1.0 ml |
| BSB-3711-CS | Control Slides | 5 |

Aurora B, RMab



IHC of Aurora B on a FFPE Lung Adenocarcinoma tissue

Aurora B kinase, a 39.3 kDa sized serine-threonine kinase, is a member of the Aurora family of mitotic kinases. The enzymatic activity of Aurora B kinase prevents stable kinetochore-microtubule attachments in early mitosis and promotes stabilization of the attachments in later mitosis due to low activity. The gene for Aurora B kinase (AURKB) is located on chromosome 17p13.1.

Abnormal expression of Aurora B kinase has been found in many cancer types, such as Non-Small Cell Lung Carcinoma, Mesothelioma, Glioblastoma, Oral Cancer and Hepatocellular Carcinoma. In Prostate and Colorectal Cancer, Aurora B expression directly correlates with the progression of cancer. In a study investigating the expression of Aurora B kinase in Thyroid Carcinoma, abundant expression of Aurora B kinase was detected in Anaplastic Carcinomas via IHC. A higher expression of Aurora B kinase in Anaplastic Thyroid Cancer than in differentiated Thyroid Cancer therefore suggests its use as a prognostic marker. In another immunohistochemical analysis, Aurora B kinase showed nuclear overexpression in high Gleason-grade Prostate Cancer, compared to low and intermediate grade cases. Additionally, there are numerous studies that indicate the use of Aurora B kinase as a prognostic marker in Ovarian Cancer. In conclusion, elevated Aurora-B expression in Breast Cancer patients contributes to chemoresistance and predicts poor prognosis.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM278


ISOTYPE: IgG

CONTROL: Tonsil, Colon, Stomach, Skin, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

B7H3/CD276, RMab



IHC of B7-H3 on a FFPE Lymphoblastic Lymphoma Tissue

B7-H3, also known as CD276, is a human protein encoded by the CD276 gene. The protein encoded by this gene belongs to the immunoglobulin superfamily, and thought to participate in the regulation of T-cell-mediated immune response. Studies show that while the transcript of this gene is ubiquitously expressed in normal tissues and solid tumors, the protein is preferentially expressed only in tumor tissues, such as melanoma, prostate cancer, and pancreatic cancer. Expression of B7-H3 protein can be induced on dendritic cells (DCs) and monocytes by inflammatory cytokines. Soluble B7-H3 protein binds a putative counter-receptor on activated T cells that is distinct from CD28, cytotoxic T lymphocyte antigen 4 (CTLA-4), inducible costimulator (ICOS) and PD-1. B7-H3 costimulates proliferation of both CD4+ and CD8+ T cells, enhances the induction of cytotoxic T cells and selectively stimulates interferon gamma (IFN-gamma) production in the presence of T cell receptor signaling.

Recently, B7-H3 expression has been reported in several human cancers indicating an additional function of B7-H3 as a regulator of antitumor immunity. B7H3 has been shown in recent years to be of clinical significance in different types of cancer. It has been demonstrated that knockdown of B7-H3 in melanoma and breast cancer cells results in both increased chemosensitivity and decreased metastatic potential, which has been observed in both in vitro and in vivo experiments.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-B7H3


ISOTYPE: IgG

CONTROL: Testis, Adrenal, Tonsil, Breast, Fallopian Tube, Breast Carcinoma, Prostate Carcinoma, Ovarian Carcinoma

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human, Mouse

BAP1, MAb



IHC of BAP1 on a FFPE Mesothelioma Tissue

In cancer, BAP1 can function both as a Tumor suppressor and as a metastasis suppressor. Exome sequencing identified inactivating mutations in BAP1 in 47% of Uveal melanomas, and BAP1 mutation have been found to be strongly associated with metastasis. The atypical melanocytic lesions resemble Spitz nevi and have been characterized as "atypical Spitz tumors" (ASTs), although they have a unique histology and exhibit both BRAF and BAP1 mutations.

BAP1 mutations have been identified in aggressive Mesotheliomas with similar mutations as seen in melanomas. Sequencing studies have been used to identify germline mutations in BAP1 in families with genetic predispositions to mesothelioma and melanocytic skin tumors. Mutations in the tumor suppressor gene BAP1 occur in approximately 15% of clear cell renal cell carcinoma cases. Sequencing efforts demonstrated worse outcomes in patients with BAP1 mutated clear cell renal cell carcinoma. Immunohistochemistry for BAP1 is a prognostic biomarker to predict poor oncologic outcomes and adverse clinicopathological features in patients with non-metastatic clear cell renal cell carcinoma. BAP1 assessment using immunohistochemistry on needle biopsy may benefit preoperative risk stratification and guide treatment planning.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-109

ISOTYPE: IgG1


CONTROL: Testis, TCC, Mesothelioma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB-3712-3 | Tinto Predilute | 3.0 ml | BSB 2810 | Tinto Predilute | 3.0 ml | BSB 3300 | Tinto Predilute | 3.0 ml |
| BSB-3712-7 | Tinto Predilute | 7.0 ml | BSB 2811 | Tinto Predilute | 7.0 ml | BSB 3301 | Tinto Predilute | 7.0 ml |
| BSB-3712-15 | Tinto Predilute | 15.0 ml | BSB 2812 | Tinto Predilute | 15.0 ml | BSB 3302 | Tinto Predilute | 15.0 ml |
| BSB-3712-01 | Concentrate | 0.1 ml | BSB 2813 | Concentrate | 0.1 ml | BSB 3303 | Concentrate | 0.1 ml |
| BSB-3712-05 | Concentrate | 0.5 ml | BSB 2814 | Concentrate | 0.5 ml | BSB 3304 | Concentrate | 0.5 ml |
| BSB-3712-1 | Concentrate | 1.0 ml | BSB 2815 | Concentrate | 1.0 ml | BSB 3305 | Concentrate | 1.0 ml |
| BSB-3712-CS | Control Slides | 5 | BSB 2816 | Control Slides | 5 | BSB 3306 | Control Slides | 5 |

BAX, RMAb




IHC of BAX on a FFPE Hodgkin's Lymphoma Tissue

BCA-225, MAb



IHC of BCA-225 on a FFPE Breast Carcinoma Tissue

bcl-2, MAb



IHC of bcl-2 on a FFPE Follicular Lymphoma Tissue

Bax is a protein of the bcl-2 gene family. It promotes apoptosis by competing with bcl-2 proper. The Bax gene contains a small promoter element that complements a binding domain on the multi-faceted p53 tumor suppressor. Wild-type p53 has been demonstrated to upregulate the transcription of a chimeric reporter plasmid, utilizing the consensus promoter sequence of Bax approx. 50-fold over mutant p53. Mutations in this consensus sequence eliminate transcription of the reporter gene. Thus, it is likely that p53 promotes Bax's apoptotic faculties in vivo as a primary transcription factor.

Bax exerts a pro-apoptotic rather than an anti-apoptotic effect on cells. Bax targets mitochondrial membranes, inducing mitochondrial damage and cell death in a caspase-independent manner. Bad plays a critical role in the Bax-mediated apoptosis pathway by dimerizing with BclL, causing the displacement of Bax. The displacement of Bax allows apoptosis to proceed.

This antibody recognizes a human breast carcinoma-associated glycoprotein, BCA-225 (220-225 kDa). This protein differs in size and distribution from other Breast Carcinoma antigens. Unlike other carcinoma antibodies against Breast Carcinoma antigens, this antibody does not react with benign or malignant colonic tissues. Since this antigen is localized in malignancies of Breast Carcinomas and Carcinoma of the Uterine Cervix, it can be effectively used to identify metastatic Breast Carcinoma lesions.

Strong intracytoplasmic staining is seen in primary and metastatic Breast Carcinoma tissue, as well as in Cervical Carcinomas. Apical staining is seen in normal kidney, lung, Fallopian tube, liver, skin (eccrine sweat glands) and uterus. Similar staining patterns are observed in lung, ovarian, and endometrial cancers. Carcinomas of the colon, stomach, prostate, urinary bladder, liver, pancreas, thyroid, and parotid are negative, as are Sarcomas and Lymphoid Cancers.

bcl-2 is an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of bcl-2, such as in the case of translocation of bcl-2 to Ig heavy-chain loci, is thought to be the cause of Follicular Lymphoma.

Anti-bcl-2 has shown consistent negative reaction on reactive germinal centers and positive staining of neoplastic follicles in Follicular Lymphoma. Consequently, this antibody is valuable when distinguishing between reactive and neoplastic follicular proliferation in lymph node biopsies. This antibody may also be used in distinguishing between those Follicular Lymphomas that express bcl-2 protein and the small number in which the neoplastic cells are bcl-2-negative. Anti-bcl-2 has been used as a predictive biomarker for recurrence of Cancer of the Breast and Non-Small-Cell Carcinoma of the Lung.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: E63

ISOTYPE: IgG

CONTROL: Breast, Tonsil, Cervix, Hodgkin's Lymphoma

LOCALIZATION: Cell Membranous, Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: Cu-18

ISOTYPE: IgG1/K

CONTROL: Breast, Lung, Uterus, Cervical Carcinoma, Breast Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-5 (BCL/A4)

ISOTYPE: IgG1/K

CONTROL: Tonsil, Lymph Node


LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6078 | Tinto Predilute | 3.0 ml |
| BSB 6079 | Tinto Predilute | 7.0 ml |
| BSB 6080 | Tinto Predilute | 15.0 ml |
| BSB 6081 | Concentrate | 0.1 ml |
| BSB 6082 | Concentrate | 0.5 ml |
| BSB 6083 | Concentrate | 1.0 ml |
| BSB 6084 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5064 | Tinto Predilute | 3.0 ml |
| BSB 5065 | Tinto Predilute | 7.0 ml |
| BSB 5066 | Tinto Predilute | 15.0 ml |
| BSB 5067 | Concentrate | 0.1 ml |
| BSB 5068 | Concentrate | 0.5 ml |
| BSB 5069 | Concentrate | 1.0 ml |
| BSB 5070 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5071 | Tinto Predilute | 3.0 ml |
| BSB 5072 | Tinto Predilute | 7.0 ml |
| BSB 5073 | Tinto Predilute | 15.0 ml |
| BSB 5074 | Concentrate | 0.1 ml |
| BSB 5075 | Concentrate | 0.5 ml |
| BSB 5076 | Concentrate | 1.0 ml |
| BSB 5077 | Control Slides | 5 |

bcl-2, RMab

IHC of bcl-2 on a FFPE Tonsil Tissue

bcl-2 is an integral outer mitochondrial membrane protein that regulates apoptosis of some cells such as lymphocytes. Constitutive expression of bcl-2, such as in the case of translocation of bcl-2 to Ig heavy-chain loci, is thought to be the cause of Follicular Lymphoma.


Anti-bcl-2 has shown consistent negative reaction on reactive germinal centers and positive staining of neoplastic follicles in Follicular Lymphoma. Consequently, this antibody is valuable when distinguishing between reactive and neoplastic follicular proliferation in lymph node biopsies. This antibody may also be used in distinguishing between those Follicular Lymphomas that express bcl-2 protein and the small number in which the neoplastic cells are bcl-2-negative. Anti-bcl-2 has been used as a predictive biomarker for recurrence of Cancer of the Breast and Non-Small-Cell Carcinoma of the Lung.

bcl-6, MAb

IHC of bcl-6 on a FFPE Colon Tissue

bcl-6 is a transcriptional regulator gene which codes for a 706-amino-acid nuclear zinc finger protein. Antibodies to this protein stain the germinal center cells in lymphoid follicles, follicular cells and interfollicular cells in Follicular Lymphoma, Diffuse Large B-Cell Lymphomas, Burkitt's Lymphoma, and the majority of the Reed-Sternberg cells in Nodular Lymphocyte-Predominant Hodgkin's Disease.

bcl-6 is also useful in identifying neoplastic cells in cases of nodular Lymphocyte-Predominant Hodgkin's Disease. In contrast, anti-bcl-6 rarely stains Mantle-Cell Lymphoma and MALT Lymphoma. bcl-6 expression is seen in approximately 45% of CD30+ Anaplastic Large-Cell Lymphomas but is consistently absent in other peripheral T-cell Lymphomas.

bcl-6, RMab

IHC of bcl-6 on a FFPE Tonsil Tissue


bcl-6 is a transcriptional regulator gene which codes for a 706-amino-acid nuclear zinc finger protein. Antibodies to this protein stain the germinal center cells in lymphoid follicles, follicular cells and interfollicular cells in Follicular Lymphoma, Diffuse Large B-Cell Lymphomas, Burkitt's Lymphoma, and the majority of the Reed-Sternberg cells in Nodular Lymphocyte-Predominant Hodgkin's Disease.

bcl-6 is also useful in identifying neoplastic cells in cases of nodular Lymphocyte-Predominant Hodgkin's Disease. In contrast, anti-bcl-6 rarely stains Mantle-Cell Lymphoma and MALT Lymphoma. bcl-6 expression is seen in approximately 45% of CD30+ Anaplastic Large-Cell Lymphomas but is consistently absent in other peripheral T-cell Lymphomas.

ANTIBODY TYPE: Rabbit Monoclonal**CLONE:** EP36**ISOTYPE:** IgG**CONTROL:** Tonsil, Lymph Node, Breast, Placenta, Fallopian Tube**LOCALIZATION:** Cytoplasmic, Membranous**SPECIES REACTIVITY:** Human, Mouse**ANTIBODY TYPE:** Mouse Monoclonal**CLONE:** BSB-26**ISOTYPE:** IgG1**CONTROL:** Tonsil, Lymph Node, Thymus, Skin, Breast, Brain, Follicular Lymphoma**LOCALIZATION:** Nuclear**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Rabbit Monoclonal**CLONE:** RBT-bcl6**ISOTYPE:** IgG**CONTROL:** Tonsil, Lymph Node, Thymus, Skin, Breast, Brain, Follicular Lymphoma**LOCALIZATION:** Nuclear**SPECIES REACTIVITY:** Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6541 | Tinto Predilute | 3.0 ml | BSB 3434 | Tinto Predilute | 3.0 ml | BSB 5078 | Tinto Predilute | 3.0 ml |
| BSB 6542 | Tinto Predilute | 7.0 ml | BSB 3435 | Tinto Predilute | 7.0 ml | BSB 5079 | Tinto Predilute | 7.0 ml |
| BSB 6543 | Tinto Predilute | 15.0 ml | BSB 3436 | Tinto Predilute | 15.0 ml | BSB 5080 | Tinto Predilute | 15.0 ml |
| BSB 6544 | Concentrate | 0.1 ml | BSB 3437 | Concentrate | 0.1 ml | BSB 5081 | Concentrate | 0.1 ml |
| BSB 6545 | Concentrate | 0.5 ml | BSB 3438 | Concentrate | 0.5 ml | BSB 5082 | Concentrate | 0.5 ml |
| BSB 6546 | Concentrate | 1.0 ml | BSB 3439 | Concentrate | 1.0 ml | BSB 5083 | Concentrate | 1.0 ml |
| BSB 6547 | Control Slides | 5 | BSB 3440 | Control Slides | 5 | BSB 5084 | Control Slides | 5 |

bcl-6, RMab




IHC of bcl-6 on a FFPE Tonsil Tissue

bcl-10, MMab



IHC of bcl-10 on a FFPE Tonsil Tissue

bcl-X, RMab



IHC of bcl-X on a FFPE Hodgkin's Lymphoma Tissue

bcl-6 is a transcriptional regulator gene which codes for a 706-amino-acid nuclear zinc finger protein. Antibodies to this protein stain the germinal center cells in lymphoid follicles, follicular cells and interfollicular cells in Follicular Lymphoma, Diffuse Large B-Cell Lymphomas, Burkitt's Lymphoma, and the majority of the Reed-Sternberg cells in Nodular Lymphocyte-Predominant Hodgkin's Disease.

bcl-6 is also useful in identifying neoplastic cells in cases of nodular Lymphocyte-Predominant Hodgkin's Disease. In contrast, anti-bcl-6 rarely stains Mantle-Cell Lymphoma and MALT Lymphoma. bcl-6 expression is seen in approximately 45% of CD30+ Anaplastic Large-Cell Lymphomas but is consistently absent in other peripheral T-cell Lymphomas.

bcl-10 (also known as B-cell lymphoma/leukemia 10) is a 233 amino acid protein encoded by the BCL10 gene on Chromosome 1. The protein encoded by this gene contains a caspase recruitment domain, and has been shown to induce apoptosis and activate NF-kappaB. This protein is reported to interact with caspase recruitment domains (CARD) containing proteins including CARD9, 10, 11 and 14. These proteins function as upstream regulators in NF-kappaB signaling.

Studies have shown that bcl-10 plays a critical role in the development of mucosa associated lymphoma tissue (MALT) lymphoma, and can be utilized in the classification of lymphomas. The antibody labels subpopulations of normal B and T cells. In MALT lymphomas with the t(1;14) translocation, the antibody strongly labeled the nuclei and cytoplasm, while 55% of MALT lymphomas lacking this translocation exhibited the same labeling pattern, although at a much lower level.

bcl-X, or bcl-2-like 1 protein, a member of the bcl-2 protein family, inhibits cell death, or apoptosis and functions as a regulator of apoptosis. bcl-X has two isoforms: bcl-XL (Long), a 241-amino acid protein; and bcl-XS (Short), a 178-amino acid protein lacking a 63-amino acid domain that is well conserved among members of the bcl-2 family.

bcl-X is typically present in the cytosol in association with the mitochondrial membrane. bcl-XL forms heterodimers with various proteins, including Bax, Bak and bcl-2. It has been found that heterodimerization with Bax does not seem to be required for anti-apoptotic activity.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP278
ISOTYPE: IgG
CONTROL: Tonsil, Follicular Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-22
ISOTYPE: IgG1/K
CONTROL: Tonsil, Kidney, Cervix, Bladder TCC, MALT Lymphomas
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP94
ISOTYPE: IgG
CONTROL: Kidney, Tonsil, Cervix, Hodgkin's Lymphoma
LOCALIZATION: Cytoplasmic Cell/Nuclear-Membrane
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2817 | Tinto Predilute | 3.0 ml |
| BSB 2818 | Tinto Predilute | 7.0 ml |
| BSB 2819 | Tinto Predilute | 15.0 ml |
| BSB 2820 | Concentrate | 0.1 ml |
| BSB 2821 | Concentrate | 0.5 ml |
| BSB 2822 | Concentrate | 1.0 ml |
| BSB 2823 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2698 | Tinto Predilute | 3.0 ml |
| BSB 2699 | Tinto Predilute | 7.0 ml |
| BSB 2700 | Tinto Predilute | 15.0 ml |
| BSB 2701 | Concentrate | 0.1 ml |
| BSB 2702 | Concentrate | 0.5 ml |
| BSB 2703 | Concentrate | 1.0 ml |
| BSB 2704 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6085 | Tinto Predilute | 3.0 ml |
| BSB 6086 | Tinto Predilute | 7.0 ml |
| BSB 6087 | Tinto Predilute | 15.0 ml |
| BSB 6088 | Concentrate | 0.1 ml |
| BSB 6089 | Concentrate | 0.5 ml |
| BSB 6090 | Concentrate | 1.0 ml |
| BSB 6091 | Control Slides | 5 |

BCOR, MAb



IHC of BCOR on a FFPE Cervix Tissue

BCOR is a gene that encodes for an epigenetic regulator involved in the specification of cell differentiation and body structure development. Various BCOR aberrations, represent driver elements of various sarcomas such as Clear Cell Sarcoma of the Kidney, Primitive Mesenchymal Myxoid Tumor of infancy, small round blue cell sarcoma, endometrial stromal sarcoma and histologically heterogeneous CNS neoplasms group with similar genomic methylation patterns known as CNS-HGNET-BCOR. Furthermore, other BCOR alterations (often loss of function mutations) recur in a large variety of Mesenchymal, Epithelial, Neural and Hematological Tumors, suggesting a central role in cancer evolution.


Recent studies have found the IHC of BCOR to be a highly sensitive marker for SBRCTs and CCSKs with BCOR abnormalities and YWHAE rearrangements and can be used as a useful diagnostic marker in these various molecular subsets.

It was concluded that BCOR immunohistochemical staining is a highly sensitive marker for YWHAE-NUTM2 high-grade Endometrial Stromal Sarcomas with both classic and unusual morphology and identifies a subset of high-grade Endometrial Stromal Sarcomas with BCOR alterations, including BCOR rearrangement and internal tandem duplication

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-128
ISOTYPE: IgG1/K
CONTROL: Testis, Cervix, Prostate, TCC, Angiosarcoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-2370-3 | Tinto Predilute | 3.0 ml |
| BSB-2370-7 | Tinto Predilute | 7.0 ml |
| BSB-2370-15 | Tinto Predilute | 15.0 ml |
| BSB-2370-01 | Concentrate | 0.1 ml |
| BSB-2370-05 | Concentrate | 0.5 ml |
| BSB-2370-1 | Concentrate | 1.0 ml |
| BSB-2370-CS | Control Slides | 5 |

Beta-Catenin, MAb



IHC of Beta-Catenin on a FFPE Breast Tissue


Beta-Catenin is a subunit of the Cadherin protein complex. Cadherins are a type of protein normally expressed on the surface of certain cells. Specifically, Beta Catenin is a 92 kDa protein normally found in the cytoplasm of the cell in the sub-membranous location. This protein is associated with E-Cadherin and may be essential for the function of E-Cadherin.

Mutations in the Beta-Catenin gene result in the nuclear accumulation of this protein. Nuclear accumulation of this protein has been demonstrated in Fibromatosis lesions of the breast and abdomen, and therefore is useful in differentiating this lesion from other spindle-cell lesions that may occur in these locations.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 14
ISOTYPE: IgG1
CONTROL: Fibromatosis of the Breast and Abdomen
LOCALIZATION: Cytoplasmic, Membranous, Nuclear
SPECIES REACTIVITY: Human, Dog, Mouse, Rat, Chicken

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5085 | Tinto Predilute | 3.0 ml |
| BSB 5086 | Tinto Predilute | 7.0 ml |
| BSB 5087 | Tinto Predilute | 15.0 ml |
| BSB 5088 | Concentrate | 0.1 ml |
| BSB 5089 | Concentrate | 0.5 ml |
| BSB 5090 | Concentrate | 1.0 ml |
| BSB 5091 | Control Slides | 5 |

Beta-Catenin, RMab



IHC of Beta-Catenin on a FFPE Salivary Gland Tissue


Beta-Catenin is a subunit of the Cadherin protein complex. Cadherins are a type of protein normally expressed on the surface of certain cells. Specifically, Beta Catenin is a 92 kDa protein normally found in the cytoplasm of the cell in the sub-membranous location. This protein is associated with E-Cadherin and may be essential for the function of E-Cadherin.

Mutations in the Beta-Catenin gene result in the nuclear accumulation of this protein. Nuclear accumulation of this protein has been demonstrated in Fibromatosis lesions of the breast and abdomen, and therefore is useful in differentiating this lesion from other spindle-cell lesions that may occur in these locations.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM276
ISOTYPE: IgG
CONTROL: Fibromatosis of the Breast & Abdomen, Breast, Adbomen, Colon, Testis, Pancreas
LOCALIZATION: Cytoplasmic, Membranous, Nuclear
SPECIES REACTIVITY: Human, Predicted: Mouse, Rat, Sheep, Hamster, Cow, Macaque Monkey, African Green Monkey


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3756-3 | Tinto Predilute | 3.0 ml |
| BSB-3756-7 | Tinto Predilute | 7.0 ml |
| BSB-3756-15 | Tinto Predilute | 15.0 ml |
| BSB-3756-01 | Concentrate | 0.1 ml |
| BSB-3756-05 | Concentrate | 0.5 ml |
| BSB-3756-1 | Concentrate | 1.0 ml |
| BSB-3756-CS | Control Slides | 5 |

BG8 LewisY, MMab




IHC of BG8 LewisY on a FFPE Adenocarcinoma Tissue

BOB.1, RMab



IHC of BOB.1 on a FFPE Tonsil Tissue

Brachyury, RMab



IHC of Brachyury on a FFPE Chordoma Tissue

Blood group antigens have been examined as potential discriminators between Pulmonary Adenocarcinoma (PACA) and Epithelioid Mesothelioma (EM). LewisY is the only one of these that appears to have some merit. BG8 is raised from the SK-LU-3 lung cancer line and is able to distinguish between PACA and EM. Studies of 231 cases of PACA and 197 cases of EM have shown that sensitivity and specificity for PACA were both 93%. It has been reported that sensitivity of nonmesothelial antigens for Adenocarcinoma is organ dependent, with BG8 Lewis performing at 98% in the breast cancer group, and 100% in the lung cancer group. The specificity of the nonmesothelial (non-EM) antigens for adenocarcinoma was 98% for BG8.

It has been concluded using logical regression analysis that a three-antibody immunohistochemical panel including Calretinin, BG8, and MOC-31 would provide 96% sensitivity and specificity for distinguishing EM from Adenocarcinoma in a variety of sources (lung, ovary, breast, stomach).

The BOB-1 protein is a co-activator that interacts with Oct1 and/or Oct2 transcription factors, and is critical in germinal center formation and immunoglobulin production. The strongest expression of BOB-1 is found in the germinal center, mantle-zone B cells, and plasma cells. Because BOB-1/OBF.1 are germinal center derived, L&H cells in Nodular Lymphocyte Predominant Hodgkin Lymphoma are consistently immunoreactive for BOB-1. Conversely, the Hodgkin/Reed-Sternberg cells in classical Hodgkin Lymphoma either do not express both or express only one of the two proteins.

In Diffuse Large B-cell Lymphomas, the highest expression levels for BOB-1/OBF.1 are reported in Follicular Center Lymphomas, Diffuse Large B-cell Lymphomas, and Burkitt Lymphomas. B-CLL, MALT-type, and Mantle Cell Lymphomas score negative or display a heterogenous/weaker activity. The strong nuclear expression of BOB-1 and Oct-2 by Germinal Center Derived Lymphomas makes these antibodies a novel class of broad spectrum B-lineage immunohistochemical markers in the differential diagnosis of Lymphomas, specifically between Primary Mediastinal B-cell Lymphoma from classical Hodgkin Disease.

Expression of the brachyury gene has been identified as a definitive diagnostic marker of chordoma, a malignant tumor that arises from remnant notochordal cells lodged in the vertebrae. Furthermore, germ line duplication of brachyury confers major susceptibility to chordoma. The chromosomal region on 6q27 containing the brachyury gene was gained in 6 of 21 chordomas (29%), and none of the 21 chordomas analyzed showed deletions that could have affected this gene.

Brachyury is an important factor in promoting the epithelial-mesenchymal transition (EMT). Cells that over-express brachyury have down-regulated expression of the adhesion molecule E-cadherin, which allows them to undergo EMT. Overexpression of brachyury has been linked to Hepatocellular carcinoma. While brachyury is promoting EMT, it can also induce metastasis of HCC cells. Brachyury expression is a prognostic biomarker for HCC, and the gene may be a target for cancer treatments in the future. Additionally, overexpression of brachyury may play a part in EMT associated with benign disease such as renal fibrosis.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: F3
ISOTYPE: IgM
CONTROL: Placenta, Tonsil, Pancreas, Cervix, Lung Adenocarcinoma, Bladder TCC
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBTBOB.1
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Nuclear, Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-TBXT
ISOTYPE: IgG
CONTROL: Testis, Chordoma, Chondroma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6373 | Tinto Predilute | 3.0 ml |
| BSB 6374 | Tinto Predilute | 7.0 ml |
| BSB 6375 | Tinto Predilute | 15.0 ml |
| BSB 6376 | Concentrate | 0.1 ml |
| BSB 6377 | Concentrate | 0.5 ml |
| BSB 6378 | Concentrate | 1.0 ml |
| BSB 6379 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3700-3 | Tinto Predilute | 3.0 ml |
| BSB-3700-7 | Tinto Predilute | 7.0 ml |
| BSB-3700-15 | Tinto Predilute | 15.0 ml |
| BSB-3700-01 | Concentrate | 0.1 ml |
| BSB-3700-05 | Concentrate | 0.5 ml |
| BSB-3700-1 | Concentrate | 1.0 ml |
| BSB-3700-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3490 | Tinto Predilute | 3.0 ml |
| BSB 3491 | Tinto Predilute | 7.0 ml |
| BSB 3492 | Tinto Predilute | 15.0 ml |
| BSB 3493 | Concentrate | 0.1 ml |
| BSB 3494 | Concentrate | 0.5 ml |
| BSB 3495 | Concentrate | 1.0 ml |
| BSB 3496 | Control Slides | 5 |

BRAF V600E, RMAb



IHC of BRAF-V600E on a FFPE Papillary Thyroid Carcinoma Tissue


BRAF is a human gene that makes a protein called B-Raf, which is more formally known as serine/threonine-protein kinase B-Raf. The B-Raf protein is involved in sending signals inside cells, which are involved in directing cell growth. Mutations in the BRAF gene can cause disease in two ways. First, mutations can be inherited and cause birth defects. Second, mutations can appear later in life and cause cancer, as an oncogene.

Mutations in this gene have been found in cancers, including non-Hodgkin lymphoma, colorectal cancer, malignant melanoma, papillary thyroid carcinoma, non-small-cell lung carcinoma, and adenocarcinoma of the lung. The frequency of BRAF mutations varies widely in human cancers, from more than 80% in melanomas and nevi, to as little as 0–18% in other tumors, such as 1–3% in lung cancers and 5% in colorectal cancer. In 90% of the cases, thymine is substituted with adenine at nucleotide 1799. This leads to valine (V) being substituted for by glutamate (E) at codon 600 (referred to as V600E) in the activation segment that has been found in human cancers. This mutation has been widely observed in papillary thyroid carcinoma, colorectal cancer, melanoma and non-small-cell lung cancer. BRAF-V600E mutation are present in 57% of Langerhans cell histiocytosis patients. The V600E mutation is a likely driver mutation in 100% of cases of hairy cell leukemia. High frequency of BRAF V600E mutations have been detected in ameloblastoma, a benign but locally infiltrative odontogenic neoplasm.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM8
ISOTYPE: IgG
CONTROL: BRAF 600E Mutated Melanoma Papillary Thyroid Cancer
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2824 | Tinto Predilute | 3.0 ml |
| BSB 2825 | Tinto Predilute | 7.0 ml |
| BSB 2826 | Tinto Predilute | 15.0 ml |
| BSB 2827 | Concentrate | 0.1 ml |
| BSB 2828 | Concentrate | 0.5 ml |
| BSB 2829 | Concentrate | 1.0 ml |
| BSB 2830 | Control Slides | 5 |

BRG-1/SMARCA4, MAb



IHC of BRG-1/SMARCA4 on a FFPE Lung Adenocarcinoma Tissue


Brahma-related gene-1 (BRG-1) protein is encoded by the gene SMARCA4, which is localized on chromosome 19. BRG-1 is the core catalytic ATPase subunit of the SWI/SNF complex. SWI/SNFs are a member of the family of ATP-dependent chromatin-remodeling complexes and the function of SWI/SNFs is to facilitate the transcriptional activation or repression of target genes. BRG-1 is essential for DNA repair, differentiation, and organ development.

There are several studies that suggest the involvement of BRG-1 in different cancer types. It was found that the loss of BRG-1 expression occurs in a portion of tested cancer types, including Breast, Colon, Head/Neck, Ovarian, Liver and Renal Cell Cancer. On the contrary, overexpression of BRG-1 was found in Breast, Colorectal, and Prostate Cancer, as well as Melanoma and Neuroblastoma. There is no established cutoff for determining high versus low expression; staining of the surrounding normal tissue has been used as a median value relative to which BRG-1/SMARCA4 expression may be considered increased or decreased. An IHC analysis of BRG-1 in Non-Small Cell Lung Cancer (NSCLC) patients revealed that the survival rate of BRG-1 negative patients was 0%, when compared to BRG-1 positive patients, indicating the prognostic value of BRG-1 as a biomarker.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-154
ISOTYPE: IgG1
CONTROL: Colon, Kidney, Prostate, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3714-3 | Tinto Predilute | 3.0 ml |
| BSB-3714-7 | Tinto Predilute | 7.0 ml |
| BSB-3714-15 | Tinto Predilute | 15.0 ml |
| BSB-3714-01 | Concentrate | 0.1 ml |
| BSB-3714-05 | Concentrate | 0.5 ml |
| BSB-3714-1 | Concentrate | 1.0 ml |
| BSB-3714-CS | Control Slides | 5 |

C-Met/HGFR, RMAb



IHC of c-Met on a FFPE Papillary Thyroid Carcinoma Tissue


c-Met is a proto-oncogene that encodes hepatocyte growth factor receptor (HGFR). The HGFR protein possesses tyrosinase-kinase activity. MET is a membrane receptor that is essential for embryonic development and wound healing, with its only known ligand being hepatocyte growth factor (HGF). Met is normally expressed by cells of epithelial origin, while expression of HGF is restricted to cells of mesenchymal origin. Upon HGF stimulation, MET induces several biological responses that collectively give rise to a program known as invasive growth.

MET is deregulated in many types of human malignancies, including cancers of kidney, liver, stomach, breast, and brain. Normally, only stem cells and progenitor cells express MET, which allows these cells to grow invasively in order to generate new tissues in an embryo or regenerate damaged tissues in an adult. However, cancer stem cells are thought to hijack the ability to express MET, and thus become the cause of cancer persistence and spread to other sites in the body (metastasis).

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP1454Y
ISOTYPE: IgG
CONTROL: Breast, Tonsil, Cervix, Papillary Thyroid Carcinoma, Colon Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6583 | Tinto Predilute | 3.0 ml |
| BSB 6584 | Tinto Predilute | 7.0 ml |
| BSB 6585 | Tinto Predilute | 15.0 ml |
| BSB 6586 | Concentrate | 0.1 ml |
| BSB 6587 | Concentrate | 0.5 ml |
| BSB 6588 | Concentrate | 1.0 ml |
| BSB 6589 | Control Slides | 5 |

C-Myc, MMab




IHC of c-Myc on a FFPE Colon Carcinoma Tissue

C-Myc, RMAb



IHC of c-Myc on a FFPE Burkitt's Lymphoma Tissue

C1q, RPAb



IHC and IF of C1q on a FFPE Lupus Erythematosus Tissue (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

Oncogene-encoded proteins c-Myc, n-Myc, and l-Myc function in cell proliferation, differentiation and neoplastic disease. A mutated version of Myc is found in many cancers, which causes Myc to be constitutively expressed. This leads to the unregulated expression of many genes, some of which are involved in cell proliferation, and results in the formation of cancer. c-Myc is a transcription factor and is a proto-oncogene that is the focal point in cell cycle regulation, metabolism, apoptosis, differentiation, cell adhesion, and tumorigenesis.

A common human translocation involving Myc is t(8;14), which is critical to the development of most cases of Burkitt's Lymphoma. Malfunctions in Myc have also been found in carcinoma of the cervix, colon, breast, lung, and stomach.

Oncogene-encoded proteins c-Myc, n-Myc, and l-Myc function in cell proliferation, differentiation and neoplastic disease. A mutated version of Myc is found in many cancers, which causes Myc to be constitutively expressed. This leads to the unregulated expression of many genes, some of which are involved in cell proliferation, and results in the formation of cancer. c-Myc is a transcription factor and is a proto-oncogene that is the focal point in cell cycle regulation, metabolism, apoptosis, differentiation, cell adhesion, and tumorigenesis.

A common human translocation involving Myc is t(8;14), which is critical to the development of most cases of Burkitt's Lymphoma. Malfunctions in Myc have also been found in carcinoma of the cervix, colon, breast, lung, and stomach.

The Complement Component 1q (C1q) is a protein complex involved in the complement system, which is part of the innate immune system. C1q together with C1r and C1s form the C1 complex. Antibodies of the adaptive immune system can bind antigen, forming an antigen-antibody complex. When C1q binds antigen antibody complexes, the C1 complex becomes activated. Activation of the C1 complex initiates the classical complement pathway of the complement system.

C1q nephropathy is a rare glomerular disease with characteristic mesangial C1q deposition noted on IHC or IF microscopy. It is histologically defined and poorly understood. Light microscopic features are heterogeneous and comprise minimal change disease (MCD), focal segmental glomerulosclerosis (FSGS), and proliferative glomerulonephritis. Clinical presentation is also diverse, and ranges from asymptomatic hematuria or proteinuria to frank nephritic or nephrotic syndrome in both children and adults. Hypertension and renal insufficiency at the time of diagnosis are common findings. Lupus nephritis is an inflammation of the kidneys caused by

Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 9E10
ISOTYPE: IgG1
CONTROL: Burkitt Lymphoma, Lung Cancer, Prostate Cancer, Breast Carcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP121
ISOTYPE: IgG
CONTROL: Burkitt Lymphoma, Lung Cancer, Prostate Cancer
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Kidney, Cervix, Spleen, Lupus Erythematosus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6863 | Tinto Predilute | 3.0 ml |
| BSB 6864 | Tinto Predilute | 7.0 ml |
| BSB 6865 | Tinto Predilute | 15.0 ml |
| BSB 6866 | Concentrate | 0.1 ml |
| BSB 6867 | Concentrate | 0.5 ml |
| BSB 6868 | Concentrate | 1.0 ml |
| BSB 6869 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6576 | Tinto Predilute | 3.0 ml |
| BSB 6577 | Tinto Predilute | 7.0 ml |
| BSB 6578 | Tinto Predilute | 15.0 ml |
| BSB 6579 | Concentrate | 0.1 ml |
| BSB 6580 | Concentrate | 0.5 ml |
| BSB 6581 | Concentrate | 1.0 ml |
| BSB 6582 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3019 | Tinto Predilute | 3.0 ml |
| BSB 3020 | Tinto Predilute | 7.0 ml |
| BSB 3021 | Tinto Predilute | 15.0 ml |
| BSB 3022 | Concentrate | 0.1 ml |
| BSB 3023 | Concentrate | 0.5 ml |
| BSB 3024 | Concentrate | 1.0 ml |
| BSB 3025 | Control Slides | 5 |


C3c, RPab



IHC and IF of C3c on a FFPE Lupus Erythematosus (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

Complement component 3, often simply called C3, is a protein of the immune system. It plays a central role in the complement system and contributes to innate immunity. C3 glomerulopathy was recently coined to describe renal biopsy appearances characterized by the presence of glomerular deposits composed predominantly of C3 in the absence of significant amounts of Ig. The presence of C3 in the absence of Ig suggests activation of complement by antibody-independent pathways, typically the alternative pathway, and many patients with this type of renal lesion have evidence of genetic or acquired alternative pathway dysregulation. C3 glomerulopathy has been further divided into dense deposit disease (DDD) and C3 glomerulonephritis (C3GN) based on electron microscopy (EM) appearances. The underlying genetic defect has been identified in some hereditary forms of C3GN such as CFHR5 nephropathy. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

C3d, RPab




IHC and IF of C3d on a FFPE Rejected Kidney Transplant (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

Complement component 3, or C3, is a protein of the immune system that plays a central role in the complement system and contributes to innate immunity. Its activation is required for both classical and alternative complement activation pathways. C3d deposition in the renal transplant PTCs (peritubular capillaries) is indicative of AR (acute rejection) with subsequent high probability of graft loss.

Anti-C3d combined with anti-C4d can be utilized as a tool for diagnosis of AR and warrant prompt and aggressive anti-rejection treatment. C3d is also a helpful adjunct in the diagnosis of bullous pemphigoid (BP) and perhaps pemphigus vulgaris (PV), especially in the cases in which only formalin-fixed, paraffin embedded tissue is available for analysis.

C4c, RPab



IHC and IF of C4c on a FFPE Kidney (IHC) and on a Frozen Glomerulonephritis Tissue (IF)

Complement component 4 (C4), in humans, is a protein involved in the intricate complement system, originating from the human leukocyte antigen (HLA) system. It serves several critical functions in immunity, tolerance, and autoimmunity with the other numerous components.

Low serum complement activity or low protein concentrations of complement C4 are found on Systemic Lupus Erythematosus (SLE) and it is often associated with Congenital C4 deficiency. Complete deficiencies of complement components are among the strongest genetic risk factors for SLE or lupus-like disease, across HLA haplotypes and racial backgrounds.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Placenta, Kidney, Fallopian Tube, Lupus Erythematosus

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Rejected Kidney Transplant

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Testis, Kidney, Pancreas, Salivary Gland, Colon

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3026 | Tinto Predilute | 3.0 ml |
| BSB 3027 | Tinto Predilute | 7.0 ml |
| BSB 3028 | Tinto Predilute | 15.0 ml |
| BSB 3029 | Concentrate | 0.1 ml |
| BSB 3030 | Concentrate | 0.5 ml |
| BSB 3031 | Concentrate | 1.0 ml |
| BSB 3032 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6387 | Tinto Predilute | 3.0 ml |
| BSB 6388 | Tinto Predilute | 7.0 ml |
| BSB 6389 | Tinto Predilute | 15.0 ml |
| BSB 6390 | Concentrate | 0.1 ml |
| BSB 6391 | Concentrate | 0.5 ml |
| BSB 6392 | Concentrate | 1.0 ml |
| BSB 6393 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3033 | Tinto Predilute | 3.0 ml |
| BSB 3034 | Tinto Predilute | 7.0 ml |
| BSB 3035 | Tinto Predilute | 15.0 ml |
| BSB 3036 | Concentrate | 0.1 ml |
| BSB 3037 | Concentrate | 0.5 ml |
| BSB 3038 | Concentrate | 1.0 ml |
| BSB 3039 | Control Slides | 5 |

C4d, RMAb




IHC and IF of C4d on a FFPE Kidney Tissue (IHC) and on a Frozen Kidney Rejection Tissue (IF)

Complement component 4, or C4, plays a central role in the complement system. C4d is the final proteolytic remnant of deposited C4b on endothelium and remains covalently attached to endothelium for little more than a week. It is easily detectable by Immunohistochemistry.

Anti-C4d combined with anti-C3d can be utilized as a tool for diagnosis of AR (Acute Rejection) and warrant prompt and aggressive anti-rejection treatment. C4d can be detected in peritubular capillaries in both chronic renal allograft rejection as well as hyperacute rejection, acute vascular rejection, acute cellular rejection, and borderline rejection. It has been shown to be a significant predictor of transplant kidney graft survival and is an aid in treating acute rejection.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP272
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Kidney Transplant Rejection
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Rat, Mouse

CA-125, MAb




IHC of CA-125 on a FFPE Ovarian Carcinoma Tissue

CA-125 reacts with malignant ovarian epithelial cells. CA-125 also reacts with antigens in Seminal Vesicle Carcinoma and Anaplastic Lymphoma.

In adult tissues, CA-125 is found in epithelial cells of Fallopian tube, endometrium and endocervix, pancreas, colon, gall bladder, stomach, kidney, apocrine sweat gland, and mammary gland. It is also found in mesothelial cell lining of pleura, pericardium and peritoneum. It is found in ovarian tumors of serous, endometrioid or clear-cell types and Adenocarcinomas of Mullerian type.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: OC125
ISOTYPE: IgG1/K
CONTROL: Colon, Pancreas, Epithelioid Mesothelioma, Ovarian Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

CA-125, RMAb



IHC of CA-125 on a FFPE Ovarian Carcinoma Tissue

CA-125 reacts with malignant ovarian epithelial cells. CA-125 also reacts with antigens in Seminal Vesicle Carcinoma and Anaplastic Lymphoma.


In adult tissues, CA-125 is found in epithelial cells of Fallopian tube, endometrium and endocervix, pancreas, colon, gall bladder, stomach, kidney, apocrine sweat gland, and mammary gland. It is also found in mesothelial cell lining of pleura, pericardium and peritoneum. It is found in ovarian tumors of serous, endometrioid or clear-cell types and Adenocarcinomas of Mullerian type.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP48
ISOTYPE: IgG
CONTROL: Colon, Pancreas, Epithelioid Mesothelioma, Ovarian Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2831 | Tinto Predilute | 3.0 ml |
| BSB 2832 | Tinto Predilute | 7.0 ml |
| BSB 2833 | Tinto Predilute | 15.0 ml |
| BSB 2834 | Concentrate | 0.1 ml |
| BSB 2835 | Concentrate | 0.5 ml |
| BSB 2836 | Concentrate | 1.0 ml |
| BSB 2837 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5092 | Tinto Predilute | 3.0 ml |
| BSB 5093 | Tinto Predilute | 7.0 ml |
| BSB 5094 | Tinto Predilute | 15.0 ml |
| BSB 5095 | Concentrate | 0.1 ml |
| BSB 5096 | Concentrate | 0.5 ml |
| BSB 5097 | Concentrate | 1.0 ml |
| BSB 5098 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6401 | Tinto Predilute | 3.0 ml |
| BSB 6402 | Tinto Predilute | 7.0 ml |
| BSB 6403 | Tinto Predilute | 15.0 ml |
| BSB 6404 | Concentrate | 0.1 ml |
| BSB 6405 | Concentrate | 0.5 ml |
| BSB 6406 | Concentrate | 1.0 ml |
| BSB 6407 | Control Slides | 5 |

CA15-3, MAb

IHC of CA15-3 on a FFPE Breast Tissue

This antibody has been used for evaluating the primary site of a metastatic carcinoma of unknown origin and distinguishing between benign and malignant lesions. It is believed that CA15-3 reacts primarily with the DF3-antigen, a 300 kDa mucin-like glycoprotein present on the apical border of secretory mammary epithelial cells.


CA15-3 has been detected with immunohistochemistry in a wide spectrum of carcinomas, including Breast Carcinomas (ductal and lobular), Sarcomas (Synovial Sarcoma and Malignant Fibrous Histiocytomas), and Lung Carcinomas. CA15-3 can be used as a supplementary marker for epithelial differentiation. CA15-3 does not stain Melanomas or Ewing's Sarcomas. Approximately 30% of Hepatocellular Carcinomas are positive for CA15-3.

CA19-9, MAb

IHC of CA19-9 on a FFPE Salivary Gland Tissue

CA19-9 (carbohydrate antigen 19-9 or sialylated Lewis (a) antigen) is a blood test from the tumor marker category. It was discovered in patients with Colon Cancer and Pancreatic Cancer in 1981. Increased levels of CA19-9 are also found in non-malignant conditions, such as Mirizzi's Syndrome and diseases of the bile duct and liver. The main use of CA19-9 is to determine whether a pancreatic tumor is secreting it; if that is the case, then the levels should fall when the tumor is treated, and they may rise again if the disease recurs.

CA19-9 antigen is highly expressed in Gastrointestinal (gastric, pancreatic, and colonic) Adenocarcinomas and salivary gland Mucoepidermoid Carcinomas. CA19-9 is usually not reactive with breast, kidney, and prostate Carcinomas, but is reactive with sialylated Lea-active pentasaccharide (sialylated lacto-N-fucopentaose II), which is enzymatically synthesized by sialylation of Type 1 carbohydrate chains.

Cadherin-6, RMAb

IHC of Cadherin-6 on a FFPE Kidney Tissue

Cadherin-6 is a member of the cadherin superfamily. Cadherins are membrane glycoproteins that mediate homophilic cell-cell adhesion and play critical roles in cell differentiation and morphogenesis. It is a type 1 cadherin and may play a role in kidney development as well as endometrium and placenta formation.

Cadherin-6 is highly expressed in kidney and the central nervous system. It has been found to be related to fetal kidney development and has been identified as a major cadherin in renal proximal tubules where conventional renal cell carcinoma originates. The expression of Cadherin-6 is associated with tumor progression in renal cell carcinoma.


ANTIBODY TYPE: Mouse Monoclonal**CLONE:** DF3**ISOTYPE:** IgG1**CONTROL:** Kidney, Breast, Pancreas, Lymph Node, Cervix, Salivary Gland, Bladder TCC**LOCALIZATION:** Cytoplasmic, Membranous**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Mouse Monoclonal**CLONE:** 121SLE**ISOTYPE:** IgM**CONTROL:** Colon, Cervix, Pancreas, Practretatic Cancer, Breast, Breast Carcinoma, Colon Carcinoma, Transitional Cell Carcinoma, Ovarian Carcinoma, Thyroid Carcinoma**LOCALIZATION:** Cytoplasmic**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Rabbit Monoclonal**CLONE:** EP217**ISOTYPE:** IgG**CONTROL:** Kidney, Renal Cell Carcinoma**LOCALIZATION:** Membranous**SPECIES REACTIVITY:** Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5099 | Tinto Predilute | 3.0 ml |
| BSB 5100 | Tinto Predilute | 7.0 ml |
| BSB 5101 | Tinto Predilute | 15.0 ml |
| BSB 5102 | Concentrate | 0.1 ml |
| BSB 5103 | Concentrate | 0.5 ml |
| BSB 5104 | Concentrate | 1.0 ml |
| BSB 5105 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5106 | Tinto Predilute | 3.0 ml |
| BSB 5107 | Tinto Predilute | 7.0 ml |
| BSB 5108 | Tinto Predilute | 15.0 ml |
| BSB 5109 | Concentrate | 0.1 ml |
| BSB 5110 | Concentrate | 0.5 ml |
| BSB 5111 | Concentrate | 1.0 ml |
| BSB 5112 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2391 | Tinto Predilute | 3.0 ml |
| BSB 2392 | Tinto Predilute | 7.0 ml |
| BSB 2393 | Tinto Predilute | 15.0 ml |
| BSB 2394 | Concentrate | 0.1 ml |
| BSB 2395 | Concentrate | 0.5 ml |
| BSB 2396 | Concentrate | 1.0 ml |
| BSB 2397 | Control Slides | 5 |

Cadherin17/LI-Cadherin, RMab




IHC of LI-Cadherin on a FFPE Colon Carcinoma Tissue

Cadherin 17 also known as LI-Cadherin, is part of the cadherin superfamily and is a calcium-dependent, membrane-associated glycoprotein. Cadherins are responsible for mediating cell-cell adhesion and are important for the structural integrity of epithelia. Cadherin 17 consists of an extracellular region containing 7 cadherin domains, and a transmembrane region but lacking the conserved cytoplasmic domain. It is a component of the gastrointestinal tract and pancreatic ducts, acting as an intestinal proton-dependent peptide transporter in the first step in oral absorption of many medically important peptide-based drugs. It may play a role in the morphological organization of liver and intestine.

In normal tissues, the Cadherin 17 antibody labels epithelial cells in the gastrointestinal tract and pancreatic duct, but not in kidney, liver and other tissues. In tumors, Cadherin 17 is expressed on adenocarcinoma of the digestive system, including liver cancer. It is a sensitive marker for the identification of gastric intestinal metaplasia and well differentiated adenocarcinomas.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP86
ISOTYPE: IgG
CONTROL: Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Calcitonin, RPaB




IHC of Calcitonin on a FFPE Thyroid Tissue

Calcitonin is a 32-amino acid polypeptide hormone that is produced in humans primarily by C-cells located in the thyroid, and in many other animals in the ultimobranchial gland. It acts to reduce blood calcium (Ca²⁺), opposing the effects of parathyroid hormone (PTH). It has been found in fish, reptiles, birds, and mammals. Its importance in humans has not been as well established as in other animals.

Immunohistochemical staining with Calcitonin antibody has proven to be an effective way of demonstrating the existence of Calcitonin-producing cells in the thyroid. C-cell Hyperplasia and Medullary Thyroid Carcinomas stain positive for Calcitonin. Studies of Calcitonin have resulted in the identification of a wide spectrum of C-cell proliferative abnormalities.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: N/A
CONTROL: Thyroid, Thyroid Medullary Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

Calcitonin, RMab



IHC of Calcitonin on a FFPE Thyroid Tissue

Calcitonin is a 32-amino acid polypeptide hormone that is produced in humans primarily by C-cells located in the thyroid, and in many other animals in the ultimobranchial gland. It acts to reduce blood calcium (Ca²⁺), opposing the effects of parathyroid hormone (PTH). It has been found in fish, reptiles, birds, and mammals. Its importance in humans has not been as well established as in other animals.

Immunohistochemical staining with Calcitonin antibody has proven to be an effective way of demonstrating the existence of Calcitonin-producing cells in the thyroid. C-cell Hyperplasia and Medullary Thyroid Carcinomas stain positive for Calcitonin. Studies of Calcitonin have resulted in the identification of a wide spectrum of C-cell proliferative abnormalities.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP92
ISOTYPE: IgG
CONTROL: Thyroid, Thyroid Medullary Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6849 | Tinto Predilute | 3.0 ml |
| BSB 6850 | Tinto Predilute | 7.0 ml |
| BSB 6851 | Tinto Predilute | 15.0 ml |
| BSB 6852 | Concentrate | 0.1 ml |
| BSB 6853 | Concentrate | 0.5 ml |
| BSB 6854 | Concentrate | 1.0 ml |
| BSB 6855 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5113 | Tinto Predilute | 3.0 ml |
| BSB 5114 | Tinto Predilute | 7.0 ml |
| BSB 5115 | Tinto Predilute | 15.0 ml |
| BSB 5116 | Concentrate | 0.1 ml |
| BSB 5117 | Concentrate | 0.5 ml |
| BSB 5118 | Concentrate | 1.0 ml |
| BSB 5119 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6408 | Tinto Predilute | 3.0 ml |
| BSB 6409 | Tinto Predilute | 7.0 ml |
| BSB 6410 | Tinto Predilute | 15.0 ml |
| BSB 6411 | Concentrate | 0.1 ml |
| BSB 6412 | Concentrate | 0.5 ml |
| BSB 6413 | Concentrate | 1.0 ml |
| BSB 6414 | Control Slides | 5 |

Caldesmon, MAb




IHC of Caldesmon on a FFPE Appendix Tissue

Caldesmon 1, also known as CALD1, is a human gene. Caldesmon is a calmodulin-binding protein. Like Calponin, Caldesmon tonically inhibits the ATPase activity of myosin in smooth muscle. This gene encodes a Calmodulin and actin-binding protein that play an essential role in the regulation of smooth muscle and nonmuscle contraction.

Two closely-related variants of human Caldesmon have been identified. The h-Caldesmon variant (120–150 kD) is predominantly expressed in smooth muscle, whereas l-Caldesmon (70–80 kD) is found in non-muscle tissue and cells. Neither of the two variants has been detected in skeletal muscle. Anti-Caldesmon recognizes only the h-Caldesmon variant. Anti-Caldesmon antibody labels smooth muscle and tumors of smooth muscle, myofibroblastic, and myoepithelial differentiation. Anti-Caldesmon has also been used to differentiate Epithelioid Mesothelioma from Serous Papillary Carcinoma of the ovary.

Calponin, MAb




IHC of Calponin on a FFPE Leiomyoma Tissue

Calponin is a 34 kDa polypeptide that interacts with actin, tropomyosin, and calmodulin. It is involved in smooth-muscle contraction mechanisms and is restricted exclusively to smooth-muscle tissue. Calponin is a calcium-binding protein. Calponin tonically inhibits the ATPase activity of myosin in smooth muscle. Phosphorylation of calponin by a protein kinase (which is dependent upon calcium binding to calmodulin) releases the calponin's inhibition of the smooth-muscle ATPase.

Calponin has been found to be useful in differentiating benign sclerosing lesions of the breast from Carcinoma. Calponin positivity has also been noted in Malignant Myoepithelioma and Pleomorphic Adenoma of Salivary Gland origin, as well as in Angiomatoid Malignant Fibrous Histiocytoma.

Calretinin, RMAb



IHC of Calretinin on a FFPE Mesothelioma Tissue

Calretinin is a vitamin D-dependent calcium-binding protein involved in calcium signaling. It is expressed in the central and peripheral nervous system and in many normal and pathological tissues. It stains Mesothelioma and can be used to help differentiate lung tumors. Calretinin is also considered an important diagnostic tool in the differential diagnosis of cystic and solid Ameloblastic Tumors.

Anti-calretinin has been shown to be useful in differentiating Mesothelioma from Adenocarcinomas of the lung and other sources. It is also useful in differentiating adrenal-cortical neoplasms from Pheochromocytomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-19

ISOTYPE: IgG1/K

CONTROL: Appendix, Uterus, Leiomyoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Rabbit, Pig

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-20

ISOTYPE: IgG1/K

CONTROL: Appendix, Uterus, Breast Ducts,

Leiomyoma, Prostate, Colon, Breast, Skin

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP1798

ISOTYPE: IgG

CONTROL: Brain, Testis, Colon, Benign


Mesotheliomal Cells, Malignant Mesothelioma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6099 | Tinto Predilute | 3.0 ml | BSB 5120 | Tinto Predilute | 3.0 ml | BSB 5127 | Tinto Predilute | 3.0 ml |
| BSB 6100 | Tinto Predilute | 7.0 ml | BSB 5121 | Tinto Predilute | 7.0 ml | BSB 5128 | Tinto Predilute | 7.0 ml |
| BSB 6101 | Tinto Predilute | 15.0 ml | BSB 5122 | Tinto Predilute | 15.0 ml | BSB 5129 | Tinto Predilute | 15.0 ml |
| BSB 6102 | Concentrate | 0.1 ml | BSB 5123 | Concentrate | 0.1 ml | BSB 5130 | Concentrate | 0.1 ml |
| BSB 6103 | Concentrate | 0.5 ml | BSB 5124 | Concentrate | 0.5 ml | BSB 5131 | Concentrate | 0.5 ml |
| BSB 6104 | Concentrate | 1.0 ml | BSB 5125 | Concentrate | 1.0 ml | BSB 5132 | Concentrate | 1.0 ml |
| BSB 6105 | Control Slides | 5 | BSB 5126 | Control Slides | 5 | BSB 5133 | Control Slides | 5 |

Calretinin, RMAb




IHC of Calretinin on a FFPE Mesothelioma Tissue

Calretinin is a vitamin D-dependent calcium-binding protein involved in calcium signaling. It is expressed in the central and peripheral nervous system and in many normal and pathological tissues. It stains Mesothelioma and can be used to help differentiate lung tumors. Calretinin is also considered an important diagnostic tool in the differential diagnosis of cystic and solid Ameloblastic Tumors.

Anti-calretinin has been shown to be useful in differentiating Mesothelioma from Adenocarcinomas of the lung and other sources. It is also useful in differentiating adrenal-cortical neoplasms from Pheochromocytomas.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM324
ISOTYPE: IgG
CONTROL: Brain, Testis, Colon, Benign Mesothelioma Cells, Malignant Mesothelioma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

Carbonic Anhydrase 9, RMAb




IHC of Carbonic Anhydrase 9 on a FFPE Kidney Tissue

Carbonic anhydrases (CAs) are a large family of zinc metalloenzymes that catalyze the reversible hydration of carbon dioxide. They participate in a variety of biological processes, including respiration, calcification, acid-base balance, bone resorption, and the formation of aqueous humor, cerebrospinal fluid, saliva, and gastric acid. They show extensive diversity in tissue distribution and in their subcellular localization.

CA9 is a transmembrane protein and the only tumor-associated CA isoenzyme known. It is expressed in all clear-cell renal cell carcinoma, but is not detected in normal kidney or most other normal tissues. It may be involved in cell proliferation and transformation. CA9 is considered to be one of the best cellular biomarkers of hypoxic regions in many solid tumors.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP161
ISOTYPE: IgG
CONTROL: Stomach, Gallbladder, Kidney Carcinoma, Cervix Carcinoma, Lung Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

Caspase-3, RMAb



IHC of Caspase-3 on a FFPE Small Lymphocytic Lymphoma Tissue

Caspase-3 is a cysteine-aspartic acid protease encoded by the gene CASP3 on chromosome 4q35.1. Caspase-3 belongs to the caspase family, which is a family of enzymes crucial to mediating apoptosis. Once activated, Caspase-3 degrades intracellular proteins as well as functional proteins and induces cell death.

Dysregulated apoptosis is a typical characteristic of human cancer. Abnormal Caspase-3 expression has been directly associated with acute Myelogenous Leukemia, where Caspase-3 overexpression has been detected. On the contrary, decreased expression of Caspase-3 has been found in Prostate Cancer. Another study has found a significant association between high Caspase-3 levels and increased death rates in Breast Cancer using IHC, highlighting the prognostic potential of Caspase-3. In another study, IHC analysis of Gastric, Ovarian, Cervical, and Colorectal Cancer demonstrated that patients with high expression of cleaved Caspase-3 had a significantly shorter overall survival time compared with those with low cleaved Caspase-3 expression. Caspase-3 is the predominant caspase involved in amyloid- β precursor protein (APP) cleavage, consistent with its marked elevation in dying neurons of Alzheimer's disease brains and colocalization of its APP cleavage product with A β in senile plaques.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM250
ISOTYPE: IgG
CONTROL: Colon, Tonsil, Testis, Fallopian Tube, Stomach, Transitional Cell Carcinoma
LOCALIZATION: Nuclear, Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3757-3 | Tinto Predilute | 3.0 ml |
| BSB-3757-7 | Tinto Predilute | 7.0 ml |
| BSB-3757-15 | Tinto Predilute | 15.0 ml |
| BSB-3757-01 | Concentrate | 0.1 ml |
| BSB-3757-05 | Concentrate | 0.5 ml |
| BSB-3757-1 | Concentrate | 1.0 ml |
| BSB-3757-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6415 | Tinto Predilute | 3.0 ml |
| BSB 6416 | Tinto Predilute | 7.0 ml |
| BSB 6417 | Tinto Predilute | 15.0 ml |
| BSB 6418 | Concentrate | 0.1 ml |
| BSB 6419 | Concentrate | 0.5 ml |
| BSB 6420 | Concentrate | 1.0 ml |
| BSB 6421 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3715-3 | Tinto Predilute | 3.0 ml |
| BSB-3715-7 | Tinto Predilute | 7.0 ml |
| BSB-3715-15 | Tinto Predilute | 15.0 ml |
| BSB-3715-01 | Concentrate | 0.1 ml |
| BSB-3715-05 | Concentrate | 0.5 ml |
| BSB-3715-1 | Concentrate | 1.0 ml |
| BSB-3715-CS | Control Slides | 5 |

Caveolin-1, RMAb



IHC of Caveolin-1 on a FFPE Placenta Tissue

CAV-1 is expressed at different levels in different tissues, with the highest in adipocytes, endothelial cells, fibroblasts, and mesothelial cells. CAV-1 is useful in assisting in the identification of epithelioid mesothelioma. CAV-1 IHC expression has been found in 100% epithelioid mesotheliomas, whereas only 7.5% of the lung adenocarcinomas were positive for CAV-1. Staining in most mesotheliomas has been reported as being strong and diffuse when compared with the weak, focal staining (no more than 1% of the tumor cells) seen in the lung adenocarcinomas and therefore CAV-1 is considered a very useful marker to help to differentiate these two malignancies. CAV-1 has been found to be comparable to other mesothelioma markers such as calretinin and podoplanin that are commonly used to assist in the differentiation between epithelioid mesotheliomas and lung adenocarcinomas. CAV-1 has also been found to be useful in the identification of Ewing sarcoma/PNET with expression in 96% cases of Ewing sarcoma/ PNET. CAV-1 is very useful in the differentiation of epithelioid mesothelioma from lung adenocarcinoma and identification of Ewing sarcoma/PNET.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP353


ISOTYPE: IgG

CONTROL: Placenta, Liver, Kidney, Spleen, Lung, Mesothelioma, Ewing's Sarcoma, RCC

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

CD1a, RMAb



IHC of CD1a on a FFPE Thymus Tissue

CD1 proteins have been demonstrated to restrict T-cell response to non-peptide lipid and glycolipid antigens. At least five CD1 genes (CD1a, b, c, d, and e) have been identified. CD1a belongs to a family of glycoproteins expressed on the surface of various human antigen-presenting cells. In particular, CD1a is a protein of 43 to 49 kDa, and has been shown to be expressed on dendritic cells and cortical thymocytes. Langerhans cells in the skin and some epithelia also express this protein. This antigen is expressed in cells comprising Langerhans Cell Histiocytosis and Langerhans Cell Sarcoma.

Anti-CD1a has been used to differentiate various cutaneous Lymphomas (T-cell) from B-cell Lymphomas and Pseudolymphomas. CD1a is also expressed by some malignancies of T-cell lineage and in Histiocytosis X.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP80


ISOTYPE: IgG

CONTROL: Skin, Thymus, Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

CD2, MMAb



IHC of CD2 on a FFPE T-Cell Lymphoma Tissue

CD2 is a cell-adhesion molecule found on the surface of T-cells and natural killer (NK) cells. It has also been called T-cell surface antigen T11/Leu-5, LFA-2, LFA-3 receptor, erythrocyte receptor and rosette receptor. Due to its structural characteristics, CD2 is a member of the immunoglobulin superfamily; it possesses two immunoglobulin-like domains in its extracellular portion. It interacts with other adhesion molecules, such as lymphocyte function-associated antigen-3 (LFA-3/CD58) in humans, or CD48 in rodents, which are expressed on the surfaces of other cells. In addition to its adhesive properties, CD2 also acts as a co-stimulatory molecule on T and NK cells.

CD2 is a surface antigen of the human T-lymphocyte lineage that is expressed on all peripheral blood T-cells. It is one of the earliest T-cell markers, being present on more than 95% of thymocytes; it is also found on some natural killer cells but not on B-lymphocytes. CD2 is implicated in the triggering of T-cells; the cytoplasmic domain is implicated in the signaling function. It is useful for the identification of Lymphomas and Leukemias of T-cell origin. As with other pan-T cell antigens, CD2 may be aberrantly deleted in some neoplastic T-cell populations, especially Peripheral T-cell Lymphomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: AB75

ISOTYPE: IgG1/K

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3314 | Tinto Predilute | 3.0 ml |
| BSB 3315 | Tinto Predilute | 7.0 ml |
| BSB 3316 | Tinto Predilute | 15.0 ml |
| BSB 3317 | Concentrate | 0.1 ml |
| BSB 3318 | Concentrate | 0.5 ml |
| BSB 3319 | Concentrate | 1.0 ml |
| BSB 3320 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5134 | Tinto Predilute | 3.0 ml |
| BSB 5135 | Tinto Predilute | 7.0 ml |
| BSB 5136 | Tinto Predilute | 15.0 ml |
| BSB 5137 | Concentrate | 0.1 ml |
| BSB 5138 | Concentrate | 0.5 ml |
| BSB 5139 | Concentrate | 1.0 ml |
| BSB 5140 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6205 | Tinto Predilute | 3.0 ml |
| BSB 6206 | Tinto Predilute | 7.0 ml |
| BSB 6207 | Tinto Predilute | 15.0 ml |
| BSB 6208 | Concentrate | 0.1 ml |
| BSB 6209 | Concentrate | 0.5 ml |
| BSB 6210 | Concentrate | 1.0 ml |
| BSB 6211 | Control Slides | 5 |

CD3, RMAb




IHC of CD3 on a FFPE Colon Tissue

CD3 Epsilon, RMAb



IHC of CD3 Epsilon on a FFPE Colon Tissue

CD4, RMAb



IHC of CD4 on a FFPE Tonsil Tissue

The CD3 antigen is a protein complex composed of three distinct chains (CD3 γ , CD3 δ and CD3 ϵ) that associate with T-cell receptors and the ζ -chain to generate an activation signal in T-lymphocytes. The TCR, ζ -chain and CD3 molecules together, comprise the TCR complex. The CD3 γ , CD3 δ , and CD3 ϵ chains are highly-related cell surface proteins of the immunoglobulin superfamily containing a single extracellular immunoglobulin domain. The intracellular tails of the CD3 molecules contain a single conserved motif known as an immunoreceptor tyrosine-based activation motif (or ITAM for short), which is essential for the signaling capacity of the TCR. Phosphorylation of the ITAM on CD3 renders the CD3 chain capable of binding the enzyme ZAP70 (zeta-associated protein), a kinase important in the signaling cascade of the T-cell.

CD3 has been considered the best all-around T-cell marker. This antibody reacts with an antigen present in early thymocytes. The positive staining of this marker may represent a sign of early commitment to the T-cell lineage.

The CD3 antigen is a protein complex composed of three distinct chains (CD3 γ , CD3 δ and CD3 ϵ) that associate with T-cell receptors and the ζ -chain to generate an activation signal in T-lymphocytes. The TCR, ζ -chain and CD3 molecules together, comprise the TCR complex. The CD3 γ , CD3 δ and CD3 ϵ chains are highly-related cell surface proteins of the immunoglobulin superfamily containing a single extracellular immunoglobulin domain. The intracellular tails of the CD3 molecules contain a single conserved motif known as an immunoreceptor tyrosine-based activation motif (or ITAM for short), which is essential for the signaling capacity of the TCR. Phosphorylation of the ITAM on CD3 renders the CD3 chain capable of binding the enzyme ZAP70 (zeta-associated protein), a kinase important in the signaling cascade of the T-cell.

CD3 has been considered the best all-around T-cell marker. This antibody reacts with an antigen present in early thymocytes. The positive staining of this marker may represent a sign of early commitment to the T-cell lineage.

CD4 is a glycoprotein expressed on the surface of T-helper cells, regulatory T-cells, monocytes, macrophages, and dendritic cells. On T-cells, CD4 is the co-receptor for the T-cell receptor (TCR). It amplifies the signal generated by the TCR by recruiting the tyrosine kinase that is essential for activating many molecules involved in the signaling cascade of an activated T-cell.

CD4 antigen is involved in the recognition of Type II Major Histocompatibility Complex antigens (MHC-II). CD4 is also the receptor for Human Immunodeficiency Virus (HIV). It is present on most T-helper cells and normal thymocytes.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-CD3

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Liver, Testis, Kidney, Colon, Spleen, Thymus, Lymphoblastic Lymphoma

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-CD3 ϵ

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-CD4


ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6422 | Tinto Predilute | 3.0 ml | BSB 5141 | Tinto Predilute | 3.0 ml | BSB 5148 | Tinto Predilute | 3.0 ml |
| BSB 6423 | Tinto Predilute | 7.0 ml | BSB 5142 | Tinto Predilute | 7.0 ml | BSB 5149 | Tinto Predilute | 7.0 ml |
| BSB 6424 | Tinto Predilute | 15.0 ml | BSB 5143 | Tinto Predilute | 15.0 ml | BSB 5150 | Tinto Predilute | 15.0 ml |
| BSB 6425 | Concentrate | 0.1 ml | BSB 5144 | Concentrate | 0.1 ml | BSB 5151 | Concentrate | 0.1 ml |
| BSB 6426 | Concentrate | 0.5 ml | BSB 5145 | Concentrate | 0.5 ml | BSB 5152 | Concentrate | 0.5 ml |
| BSB 6427 | Concentrate | 1.0 ml | BSB 5146 | Concentrate | 1.0 ml | BSB 5153 | Concentrate | 1.0 ml |
| BSB 6428 | Control Slides | 5 | BSB 5147 | Control Slides | 5 | BSB 5154 | Control Slides | 5 |

CD5, RMab

IHC of CD5 on a FFPE Thymus Tissue

CD5 is a glycoprotein monomer with an MW of 67 kDa belonging to the scavenger receptor cysteine-rich (SRCR) family of extracellular domain-like structures. It possesses a large cytoplasmic domain suitable for signal transduction.


CD5 is a T-cell marker that also reacts with a range of neoplastic B-cells, e.g., B-cell Chronic Lymphocytic Leukemia (B-CLL), B-cell Small Lymphocytic Lymphoma (BSLL), and Mantle Cell Lymphoma. CD5 is expressed in T-lymphocyte subsets and is modulated during cellular activation; however, it does not react with granulocytes or monocytes.

CD5, RMab

IHC of CD5 on a FFPE Tonsil Tissue

CD5 is a glycoprotein monomer with an MW of 67 kDa belonging to the scavenger receptor cysteine-rich (SRCR) family of extracellular domain-like structures. It possesses a large cytoplasmic domain suitable for signal transduction.

CD5 is a T-cell marker that also reacts with a range of neoplastic B-cells, e.g., B-cell Chronic Lymphocytic Leukemia (B-CLL), B-cell Small Lymphocytic Lymphoma (BSLL), and Mantle Cell Lymphoma. CD5 is expressed in T-lymphocyte subsets and is modulated during cellular activation; however, it does not react with granulocytes or monocytes.

CD6, MMab

IHC of CD6 on a FFPE Thymus Tissue

CD6 (also known as Cluster of Differentiation 6) is a human protein encoded by the CD6 gene on Chromosome 11. CD6 protein is found on the outer membrane of T-lymphocytes.

The encoded CD6 protein contains three scavenger receptor cysteine-rich (SRCR) domains and a binding site for an activated leukocyte cell adhesion molecule. Studies have shown CD6 regulates T-cell responses through activation-dependent recruitment of the positive regulator SLP-76. The gene product is important for continuation of T cell activation. This gene may be associated with susceptibility to multiple sclerosis.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-CD5
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM314
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen, Thymus
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-54
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node, Prostate, Colon, Spleen
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|-------------|-----------------|---------|----------|-----------------|---------|
| BSB 5155 | Tinto Predilute | 3.0 ml | BSB-3759-3 | Tinto Predilute | 3.0 ml | BSB 2705 | Tinto Predilute | 3.0 ml |
| BSB 5156 | Tinto Predilute | 7.0 ml | BSB-3759-7 | Tinto Predilute | 7.0 ml | BSB 2706 | Tinto Predilute | 7.0 ml |
| BSB 5157 | Tinto Predilute | 15.0 ml | BSB-3759-15 | Tinto Predilute | 15.0 ml | BSB 2707 | Tinto Predilute | 15.0 ml |
| BSB 5158 | Concentrate | 0.1 ml | BSB-3759-01 | Concentrate | 0.1 ml | BSB 2708 | Concentrate | 0.1 ml |
| BSB 5159 | Concentrate | 0.5 ml | BSB-3759-05 | Concentrate | 0.5 ml | BSB 2709 | Concentrate | 0.5 ml |
| BSB 5160 | Concentrate | 1.0 ml | BSB-3759-1 | Concentrate | 1.0 ml | BSB 2710 | Concentrate | 1.0 ml |
| BSB 5161 | Control Slides | 5 | BSB-3759-CS | Control Slides | 5 | BSB 2711 | Control Slides | 5 |

CD7, MMab




IHC of CD7 on a T-Cell Tonsil Tissue

CD7, R Mab



IHC of CD7 on a FFPE Tonsil Tissue

CD8, MMab



IHC of CD8 on a FFPE Tonsil Tissue

CD7 is a 40 kDa transmembrane, single-chain glycoprotein, which is a member of the immunoglobulin gene superfamily. It is expressed in the majority of immature and mature T-lymphocytes, and T-cell Leukemia. It is also found in natural killer cells, a small subpopulation of normal B-cells and in malignant B-cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

CD7 is a consistently-expressed T-cell antigen in Lymphoblastic Lymphomas and Leukemias; therefore, it is a useful marker in the identification of such neoplastic proliferations. CD7 is expressed in the majority of mature peripheral T-cells, the majority of post-thymic T-cells, NK cells, some myeloid cells, T-cell Acute Lymphoblastic Leukemia/Lymphoma, Acute Myelogenous Leukemia and Chronic Myelogenous Leukemia. Interestingly, CD7 is conspicuously absent in adult T-cell Leukemia/Lymphoma and is not expressed in Sezary cells.

CD7 is a 40 kDa transmembrane, single-chain glycoprotein, which is a member of the immunoglobulin gene superfamily. It is expressed in the majority of immature and mature T-lymphocytes and T-cell Leukemia. CD7 is also found in natural killer cells, a small subpopulation of normal B-cells and in malignant B-cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interactions during early lymphoid development.

CD7 is a consistently-expressed T-cell antigen in Lymphoblastic Lymphomas and Leukemias; therefore, it is a useful marker in the identification of such neoplastic proliferations. CD7 is expressed in the majority of mature peripheral T-cells, the majority of post-thymic T-cells, NK cells, some myeloid cells, T-cell Acute Lymphoblastic Leukemia/Lymphoma, Acute Myelogenous Leukemia and Chronic Myelogenous Leukemia. Interestingly, CD7 is conspicuously absent in adult T-cell Leukemia/Lymphoma and is not expressed in Sezary cells.

CD8 is a transmembrane glycoprotein that serves as a co-receptor for the T-cell receptor (TCR). Like the TCR, CD8 binds to a major histocompatibility complex (MHC) molecule that is specific for the Class I MHC protein. To function, CD8 forms a dimer, consisting of a pair of CD8 chains. The most common form of CD8 is composed of a CD8- α and CD8- β chain, both members of the immunoglobulin superfamily with an immunoglobulin variable (IgV)-like extracellular domain connected to the membrane by a thin stalk, and an intracellular tail.

CD8 is a T-cell marker for the detection of cytotoxic/suppressor cells of blood lymphocytes. CD8 is also detected on NK cells, most thymocytes, a subpopulation of null cells and bone marrow cells. This antibody is used to distinguish between reactive and neoplastic T-cells.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: LP15

ISOTYPE: IgG2b

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP132

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Colon, Liver, Spleen, Bone Marrow, Lymphoblastic Lymphoma

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: C8/144B

ISOTYPE: IgG1/K

CONTROL: Tonsil, Lymph Node


LOCALIZATION: Membranous

SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5162 | Tinto Predilute | 3.0 ml |
| BSB 5163 | Tinto Predilute | 7.0 ml |
| BSB 5164 | Tinto Predilute | 15.0 ml |
| BSB 5165 | Concentrate | 0.1 ml |
| BSB 5166 | Concentrate | 0.5 ml |
| BSB 5167 | Concentrate | 1.0 ml |
| BSB 5168 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2321 | Tinto Predilute | 3.0 ml |
| BSB 2322 | Tinto Predilute | 7.0 ml |
| BSB 2323 | Tinto Predilute | 15.0 ml |
| BSB 2324 | Concentrate | 0.1 ml |
| BSB 2325 | Concentrate | 0.5 ml |
| BSB 2326 | Concentrate | 1.0 ml |
| BSB 2327 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5169 | Tinto Predilute | 3.0 ml |
| BSB 5170 | Tinto Predilute | 7.0 ml |
| BSB 5171 | Tinto Predilute | 15.0 ml |
| BSB 5172 | Concentrate | 0.1 ml |
| BSB 5173 | Concentrate | 0.5 ml |
| BSB 5174 | Concentrate | 1.0 ml |
| BSB 5175 | Control Slides | 5 |

CD8, RMab

IHC of CD8 on a FFPE Spleen Tissue

CD8 is a transmembrane glycoprotein that serves as a co-receptor for the T-cell receptor (TCR). Like the TCR, CD8 binds to a major histocompatibility complex (MHC) molecule that is specific for the Class I MHC protein. To function, CD8 forms a dimer, consisting of a pair of CD8 chains. The most common form of CD8 is composed of a CD8- α and CD8- β chain, both members of the immunoglobulin superfamily with an immunoglobulin variable (IgV)-like extracellular domain connected to the membrane by a thin stalk, and an intracellular tail.


CD8 is a T-cell marker for the detection of cytotoxic/suppressor cells of blood lymphocytes. CD8 is also detected on NK cells, most thymocytes, a subpopulation of null cells and bone marrow cells. This antibody is used to distinguish between reactive and neoplastic T-cells.

CD10, MAb

IHC of CD10 on a FFPE Kidney Tissue

CD10, also known as neutral endopeptidase (NEP), Nephilysin, and common Acute Lymphoblastic Leukemia antigen (CALLA), is a zinc-dependent metalloprotease enzyme that degrades a number of small secreted peptides, most notably the amyloid beta peptide whose abnormal misfolding and aggregation in neural tissue has been implicated as a cause of Alzheimer's Disease.

CD10 is a useful marker for the characterization of childhood Leukemia and B-cell Lymphomas. This antibody reacts with the antigens of Lymphoblastic, Burkitt's, and Follicular Lymphomas, and Chronic Myelocytic Leukemia. Also, CD10 detects the antigen of glomerular epithelial cells and the brush border of the proximal tubules. This characteristic may be helpful in interpreting renal ontogenesis in conjunction with other markers. Other non-lymphoid cells that are reactive with CD10 are breast myoepithelial cells, bile canaliculi, neutrophils, a small population of bone marrow cells, fetal small intestine epithelium, and normal fibroblasts.

CD10, RMab

IHC of CD10 on a FFPE Tonsil Tissue

CD10, also known as neutral endopeptidase (NEP), Nephilysin, and common Acute Lymphoblastic Leukemia antigen (CALLA), is a zinc-dependent metalloprotease enzyme that degrades a number of small secreted peptides, most notably the amyloid beta peptide whose abnormal misfolding and aggregation in neural tissue has been implicated as a cause of Alzheimer's Disease.

CD10 is a useful marker for the characterization of childhood Leukemia and B-cell Lymphomas. This antibody reacts with the antigens of Lymphoblastic, Burkitt's, and Follicular Lymphomas, and Chronic Myelocytic Leukemia. Also, CD10 detects the antigen of glomerular epithelial cells and the brush border of the proximal tubules. This characteristic may be helpful in interpreting renal ontogenesis in conjunction with other markers. Other non-lymphoid cells that are reactive with CD10 are breast myoepithelial cells, bile canaliculi, neutrophils, a small population of bone marrow cells, fetal small intestine epithelium, and normal fibroblasts.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP334
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Liver, Colon
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 56C6
ISOTYPE: IgG1
CONTROL: Kidney, Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Dog, Cat, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP195
ISOTYPE: IgG
CONTROL: Kidney, Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2845 | Tinto Predilute | 3.0 ml | BSB 5176 | Tinto Predilute | 3.0 ml | BSB 6429 | Tinto Predilute | 3.0 ml |
| BSB 2846 | Tinto Predilute | 7.0 ml | BSB 5177 | Tinto Predilute | 7.0 ml | BSB 6430 | Tinto Predilute | 7.0 ml |
| BSB 2847 | Tinto Predilute | 15.0 ml | BSB 5178 | Tinto Predilute | 15.0 ml | BSB 6431 | Tinto Predilute | 15.0 ml |
| BSB 2848 | Concentrate | 0.1 ml | BSB 5179 | Concentrate | 0.1 ml | BSB 6432 | Concentrate | 0.1 ml |
| BSB 2849 | Concentrate | 0.5 ml | BSB 5180 | Concentrate | 0.5 ml | BSB 6433 | Concentrate | 0.5 ml |
| BSB 2850 | Concentrate | 1.0 ml | BSB 5181 | Concentrate | 1.0 ml | BSB 6434 | Concentrate | 1.0 ml |
| BSB 2851 | Control Slides | 5 | BSB 5182 | Control Slides | 5 | BSB 6435 | Control Slides | 5 |

CD11b, RMab




IHC of CD11b on a FFPE Spleen Tissue

CD11c, RMab



IHC of CD11c on a FFPE Spleen Tissue

CD13, MMab



IHC of CD13 on a FFPE Tonsil Tissue

Integrin alpha M (ITGAM) is one protein subunit that forms the heterodimeric integrin alpha-M beta-2 molecule, also known as macrophage-1 antigen (Mac-1) or complement receptor 3 (CR3). ITGAM is also known as CR3A and cluster of differentiation molecule 11b (CD11b). The $\alpha\beta 2$ molecule is expressed on the surface of many leukocytes involved in the innate immune system, including monocytes, granulocytes, macrophages, and natural killer cells. It mediates inflammation by regulating leukocyte adhesion and migration and has been implicated in several immune processes such as phagocytosis, cell-mediated cytotoxicity, chemotaxis and cellular activation. The ITGAM subunit of $\alpha\beta 2$ is directly involved in causing the adhesion and spreading of cells but cannot mediate cellular migration without the presence of the Beta-2 (CD18) subunit.

CD11b has been used as a common myeloid marker and is expressed in about 50% of acute myeloid leukemia (AML). In combination with CD117, CD11b is helpful in differentiating acute promyelocytic leukemia (CD11b negative) from recovering benign myeloid proliferation (CD11b positive, CD117 negative).

CD11c (1TGAX) is a member of the leukointegrin family and shares the same beta subunit with other members of the leukocyte adhesion molecule family, which include CD11a (LFA-1), CD11b (MAC-1), and CD11d (ITGAD), but has a unique alpha chain. CD11c is a type I transmembrane protein found at high levels on most human dendritic cells, but also on monocytes, macrophages, neutrophils, and some B cells that induces cellular activation and helps trigger neutrophil respiratory burst.

CD11c is expressed in Hairy Cell leukemias, Acute Nonlymphocytic Leukemias, and some B-cell Chronic Lymphocytic Leukemias.

CD13 (also known as aminopeptidase-N) is expressed on the majority of peripheral blood monocytes and granulocytes. It is also expressed by the majority of acute myeloid leukemias, chronic myeloid leukemias in myeloid blast crisis, a smaller percentage of lymphoid leukemias and myeloid cell lines. CD13 is absent from normal lymphocytes, platelets and erythrocytes. CD13 is also present on fibroblasts, endothelial cells, epithelial cells from renal proximal tubules and intestinal brush border, bone marrow stromal cells, osteoclasts, and cells forming bile canaliculi.

Anti-Human CD13 recognizes the human CD13 antigen expressed on the majority of peripheral blood monocytes and granulocytes and on endothelial cells. CD13 plays a role in biologically active peptide metabolism, in the control of growth and differentiation, in phagocytosis and in bactericidal/tumoricidal activities. CD13 also serves as a receptor for human coronaviruses (HCV).

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP45

ISOTYPE: IgG

CONTROL: Spleen, Leukemia

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP157

ISOTYPE: IgG

CONTROL: Bone Marrow, Spleen, Tonsil, Colon, Liver, Hairy Cell Leukemia

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 38C12

ISOTYPE: IgG1

CONTROL: Spleen, Tonsil, Prostate, Liver

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6436 | Tinto Predilute | 3.0 ml |
| BSB 6437 | Tinto Predilute | 7.0 ml |
| BSB 6438 | Tinto Predilute | 15.0 ml |
| BSB 6439 | Concentrate | 0.1 ml |
| BSB 6440 | Concentrate | 0.5 ml |
| BSB 6441 | Concentrate | 1.0 ml |
| BSB 6442 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6443 | Tinto Predilute | 3.0 ml |
| BSB 6444 | Tinto Predilute | 7.0 ml |
| BSB 6445 | Tinto Predilute | 15.0 ml |
| BSB 6446 | Concentrate | 0.1 ml |
| BSB 6447 | Concentrate | 0.5 ml |
| BSB 6448 | Concentrate | 1.0 ml |
| BSB 6449 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6324 | Tinto Predilute | 3.0 ml |
| BSB 6325 | Tinto Predilute | 7.0 ml |
| BSB 6326 | Tinto Predilute | 15.0 ml |
| BSB 6327 | Concentrate | 0.1 ml |
| BSB 6328 | Concentrate | 0.5 ml |
| BSB 6329 | Concentrate | 1.0 ml |
| BSB 6330 | Control Slides | 5 |

CD13, RMab




IHC of CD13 on a FFPE Tonsil Tissue

CD14, MAb



IHC of CD14 on a FFPE Placenta Tissue

CD14, RMab



IHC of CD14 on a FFPE Placenta Tissue

CD13 (also known as aminopeptidase-N) is expressed on the majority of peripheral blood monocytes and granulocytes. It is also expressed by the majority of acute myeloid leukemias, chronic myeloid leukemias in myeloid blast crisis, a smaller percentage of lymphoid leukemias and myeloid cell lines. CD13 is absent in normal lymphocytes, platelets and erythrocytes. CD13 is also present on fibroblasts, endothelial cells, epithelial cells from renal proximal tubules and intestinal brush border, bone marrow stromal cells, osteoclasts, and cells forming bile canaliculi.

Anti-Human CD13 recognizes the human CD13 antigen expressed on the majority of peripheral blood monocytes and granulocytes and on endothelial cells. CD13 plays a role in biologically active peptide metabolism, in the control of growth and differentiation, in phagocytosis and in bactericidal/tumoricidal activities. CD13 also serves as a receptor for human coronaviruses (HCV).

CD14 is a human gene. The protein encoded by this gene is a component of the innate immune system. CD14 exists in two forms. It is either anchored into the membrane by a glycosylphosphatidylinositol tail (mCD14), or it appears in a soluble form (sCD14). CD14 acts as a co-receptor (along with the Toll-like receptor TLR 4 and MD-2) for the detection of bacterial lipopolysaccharide (LPS). CD14 can only bind LPS in the presence of lipopolysaccharide-binding protein (LBP). Although LPS is considered its main ligand, CD14 also recognizes other pathogen associated molecular patterns.

CD14 is expressed mainly by macrophages and (at 10 times lesser extent) by neutrophil granulocytes. A soluble form sCD14 is secreted by the liver and monocytes and is sufficient in low concentrations to confer LPS-responsiveness to cells which otherwise do not express CD14. sCD14 is also present in human milk where it is believed to regulate microbial growth in the infant gut. Increased sCD14 levels are associated with inflammatory infectious diseases and high mortality in gram-negative shock. CD14 also appears to be involved in clearance of gram-negative bacteria via its high affinity binding to LPS-LBP complexes.

CD14 is a component of the innate immune system. CD14 exists in two forms: mCD14 (anchored by a glycosylphosphatidylinositol tail) or sCD14 (soluble). CD14 acts as a co-receptor (along with the Toll-like receptor TLR 4 and MD-2) for the detection of bacterial lipopolysaccharide (LPS). CD14 can only bind LPS in the presence of lipopolysaccharide-binding protein (LBP). Although LPS is considered its main ligand, CD14 also recognizes other pathogen associated molecular patterns.

CD14 is expressed mainly by macrophages and (at 10 times lesser extent) by neutrophil granulocytes. A soluble form sCD14 is secreted by the liver and monocytes and is sufficient in low concentrations to confer LPS-responsiveness to cells which otherwise do not express CD14. sCD14 is also present in human milk where it is believed to regulate microbial growth in the infant gut. Increased sCD14 levels are associated with inflammatory infectious diseases and high mortality in gram-negative shock. CD14 also appears to be involved in clearance of gram-negative bacteria via its high affinity binding to LPS-LBP complexes.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP117
ISOTYPE: IgG
CONTROL: Spleen, Tonsil, Prostate, Liver
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Rat

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 7
ISOTYPE: Ig2Ga
CONTROL: Placenta, Tonsil, Spleen, Diffuse Large B-cell Lymphoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP128
ISOTYPE: IgG
CONTROL: Placenta, Tonsil, Spleen, Diffuse Large B-cell Lymphoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6450 | Tinto Predilute | 3.0 ml |
| BSB 6451 | Tinto Predilute | 7.0 ml |
| BSB 6452 | Tinto Predilute | 15.0 ml |
| BSB 6453 | Concentrate | 0.1 ml |
| BSB 6454 | Concentrate | 0.5 ml |
| BSB 6455 | Concentrate | 1.0 ml |
| BSB 6456 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6310 | Tinto Predilute | 3.0 ml |
| BSB 6311 | Tinto Predilute | 7.0 ml |
| BSB 6312 | Tinto Predilute | 15.0 ml |
| BSB 6313 | Concentrate | 0.1 ml |
| BSB 6314 | Concentrate | 0.5 ml |
| BSB 6315 | Concentrate | 1.0 ml |
| BSB 6316 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6457 | Tinto Predilute | 3.0 ml |
| BSB 6458 | Tinto Predilute | 7.0 ml |
| BSB 6459 | Tinto Predilute | 15.0 ml |
| BSB 6460 | Concentrate | 0.1 ml |
| BSB 6461 | Concentrate | 0.5 ml |
| BSB 6462 | Concentrate | 1.0 ml |
| BSB 6463 | Control Slides | 5 |

CD15, MAb



IHC of CD15 on a FFPE Hodgkin's Lymphoma Tissue

CD15 is a phosphatidylinositol-anchored transmembrane protein found on neutrophils and which may be involved in phagocytosis. It is expressed in patients with Hodgkin's Disease, some B-cell Chronic Lymphocytic Leukemias, Acute Lymphoblastic Leukemias, and most Acute Non-Lymphocytic Leukemias. It is also called Lewis x.

A positive reaction for CD15 combined with a negative reaction for CD45 and other B and T-lineage markers provides support for Reed-Sternberg cells found in Hodgkin's disease. Also, this antibody does not detect Mesotheliomas, making it a more frequently used antibody to distinguish Epithelial Mesothelioma from Adenocarcinoma.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-119


ISOTYPE: IgM

CONTROL: Tonsil, Lymph Node, Hodgkin's Lymphoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

CD16, RMAb



IHC of CD16 on a FFPE Bone Marrow Tissue

CD16 is a low affinity Fc receptor, found on the surface of natural killer cells, neutrophil polymorphonuclear leukocytes, monocytes and macrophages. CD16 has been identified as Fc receptors FcγRIIIa (CD16a) and FcγRIIIb (CD16b). These receptors bind to the Fc portion of IgG antibodies which then activates NK cells for antibody-dependent cell-mediated cytotoxicity. A lack of CD16 in a given population of neutrophils may indicate prematurity, as could be caused by a left shift due

to neutrophilic leukocytosis induced by tissue necrosis or bacterial infection.

The IHC of CD16 is useful in the differential diagnosis of hepatosplenic gamma delta T-cell lymphoma and gamma delta T-cell large granular lymphocyte leukemia from other peripheral T-cell lymphomas, such as mucosal and cutaneous gamma delta T-cell lymphoma. A significant decrease can be seen in the number of granulocytes expressing CD16 in chronic myelomonocytic leukemia compared to chronic myelogenous leukemia and control bone marrow biopsy, probably related to dysgranulopoiesis. It has also been demonstrated that colorectal carcinoma patients with high CD16+ cell infiltration is associated with improved overall survival after adjusting for known prognostic factors and this association was independent from CD8+ lymphocyte infiltration and presence of metastases.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP364


ISOTYPE: IgG

CONTROL: Placenta, Liver, Breast, Spleen, Thymus, Lung

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

CD19, MAb



IHC of CD19 on a FFPE Colon Tissue

CD19 is a human protein encoded by the CD19 gene. CD19 is expressed on follicular dendritic cells and B-cells; it is present on B-cells from earliest recognizable B-lineage cells during development to B-cell blasts, but is lost on maturation to plasma cells. In normal lymphoid tissue, CD19 is observed in germinal centers (on both B-cells and follicular dendritic cells), in mantle-zone cells, and in scattered cells in the interfollicular areas, with an overall immunoreactivity pattern similar to that of CD20 and CD22. However, in contrast to CD20, CD19 is also expressed in pre-B-cells.

CD19 positivity is seen in the vast majority of B-cell neoplasms (B-Lymphoblastic Lymphoma, Small Lymphocytic Lymphoma/CLL, Mantle Cell Lymphoma, Follicular Lymphoma, Burkitt's Lymphoma, Marginal Zone Lymphoma, Diffuse Large B-cell Lymphoma, T-cell-rich B-cell Lymphoma, Lymphoblastic Lymphoma, Hairy Cell Leukemia), and commonly at a lower intensity than normal B-cell elements. Plasma cell neoplasms are consistently negative, as are T-cell neoplasms. CD19 expression is not seen in Reed-Sternberg cells of classic Hodgkin's Disease.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-97

ISOTYPE: IgG1

CONTROL: Tonsil, Lymph Node, Spleen, Colon

LOCALIZATION: Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5183 | Tinto Predilute | 3.0 ml |
| BSB 5184 | Tinto Predilute | 7.0 ml |
| BSB 5185 | Tinto Predilute | 15.0 ml |
| BSB 5186 | Concentrate | 0.1 ml |
| BSB 5187 | Concentrate | 0.5 ml |
| BSB 5188 | Concentrate | 1.0 ml |
| BSB 5189 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3321 | Tinto Predilute | 3.0 ml |
| BSB 3322 | Tinto Predilute | 7.0 ml |
| BSB 3323 | Tinto Predilute | 15.0 ml |
| BSB 3324 | Concentrate | 0.1 ml |
| BSB 3325 | Concentrate | 0.5 ml |
| BSB 3326 | Concentrate | 1.0 ml |
| BSB 3327 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6226 | Tinto Predilute | 3.0 ml |
| BSB 6227 | Tinto Predilute | 7.0 ml |
| BSB 6228 | Tinto Predilute | 15.0 ml |
| BSB 6229 | Concentrate | 0.1 ml |
| BSB 6230 | Concentrate | 0.5 ml |
| BSB 6231 | Concentrate | 1.0 ml |
| BSB 6232 | Control Slides | 5 |

CD20, MMab




IHC of CD20 on a FFPE Colon Tissue

CD20 is a transmembrane, non-glycosylated protein expressed on B-cell precursors and mature B-cells, but lost following differentiation into plasma cells. This antibody does not cross-react with non-hematopoietic neoplasms. CD20 (B-cell Pan) reacts with a membrane antigen present in B-cells.

This antibody strongly recognizes Reed-Sternberg cells predominant in Hodgkin's disease. Since no staining of histiocytes or plasma cells has been observed and CD20 has not been detected in T-cell malignancies, it is a very strong marker of B-cell Lymphomas. B-cell Panmarker recognizes a formalin-resistant intracytoplasmic antigen.

CD21, RMAb




IHC of CD21 on a FFPE Tonsil Tissue

CD21, also known as CR2, complement component (3d/Epstein Barr virus) receptor 2, is an integral membrane glycoprotein of molecular weight 140 kDa, involved in the complement system. CD21 binds to C3d. B-cells have CR2 receptors on their surfaces, allowing the complement system to play a role in B-cell activation and maturation. Complement component receptor-2 (CR2) is the membrane protein on B-lymphocytes to which the Epstein-Barr virus (EBV) binds during infection of these cells.

Anti-CD21 is useful in the identification of follicular dendritic cell matrixes found in normal lymph nodes and tonsillar tissue. This antibody also labels Follicular Dendritic Cell Tumor/Sarcomas. The antigen is absent on T-lymphocytes, monocytes, and granulocytes.

CD23, MMab



IHC of CD23 on a FFPE Tonsil Tissue

CD23, also known as Fc epsilon RII, is the "low affinity" receptor for IgE, an antibody isotype involved in allergy and (arguably) resistance to parasites, and is important in regulation of IgE levels. Unlike many of the antibody receptors, CD23 is a C-type lectin. It is found on mature B-cells, activated macrophages, eosinophils, follicular dendritic cells and platelets.

This is a B-cell antibody that is useful for differentiating between B-CLL and B-SLL's that are CD23-positive from Mantle-cell Lymphomas and Small-Cleaved Lymphomas that are CD23-negative. This antibody reacts with the antigen that is found on a subpopulation of peripheral blood cells, B-lymphocytes and on EBV-transformed B-lymphoblastoid cell lines.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: L26
ISOTYPE: IgG2A/k
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human, Dog, Cat, Mouse

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP64
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human, Mouse


ANTIBODY TYPE: Mouse Monoclonal
CLONE: 1B12
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5190 | Tinto Predilute | 3.0 ml |
| BSB 5191 | Tinto Predilute | 7.0 ml |
| BSB 5192 | Tinto Predilute | 15.0 ml |
| BSB 5193 | Concentrate | 0.1 ml |
| BSB 5194 | Concentrate | 0.5 ml |
| BSB 5195 | Concentrate | 1.0 ml |
| BSB 5196 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5197 | Tinto Predilute | 3.0 ml |
| BSB 5198 | Tinto Predilute | 7.0 ml |
| BSB 5199 | Tinto Predilute | 15.0 ml |
| BSB 5200 | Concentrate | 0.1 ml |
| BSB 5201 | Concentrate | 0.5 ml |
| BSB 5202 | Concentrate | 1.0 ml |
| BSB 5203 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5204 | Tinto Predilute | 3.0 ml |
| BSB 5205 | Tinto Predilute | 7.0 ml |
| BSB 5206 | Tinto Predilute | 15.0 ml |
| BSB 5207 | Concentrate | 0.1 ml |
| BSB 5208 | Concentrate | 0.5 ml |
| BSB 5209 | Concentrate | 1.0 ml |
| BSB 5210 | Control Slides | 5 |

CD23, RMAb




IHC of CD23 on a FFPE Lymphoma Tissue

CD25, MMAb



IHC of CD25 on a FFPE Thymus Tissue

CD25, RMAb



IHC of CD25 on a FFPE Pituitary Tissue

CD23, also known as Fc epsilon RII, is the “low affinity” receptor for IgE, an antibody isotype involved in allergy and (arguably) resistance to parasites, and is important in regulation of IgE levels. Unlike many of the antibody receptors, CD23 is a C-type lectin. It is found on mature B-cells, activated macrophages, eosinophils, follicular dendritic cells and platelets.

This is a B-cell antibody that is useful for differentiating between B-CLL and B-SLLs that are CD23-positive from Mantle-cell Lymphomas and Small-Cleaved Lymphomas that are CD23- negative. This antibody reacts with the antigen that is found on a subpopulation of peripheral blood cells, B-lymphocytes and on EBV-transformed B-lymphoblastoid cell lines.

CD25 is the alpha chain of the IL-2 receptor. It is a Type I transmembrane protein present on activated T-cells, activated B-cells, some thymocytes, myeloid precursors, and oligodendrocytes that associates with CD122 to form a heterodimer that can act as a high-affinity receptor for IL-2. It is expressed in most B-cell neoplasms, some Acute Non-lymphocytic Leukemias, and Neuroblastomas.

Expression of CD25 is a reliable diagnostic tool for distinguishing neoplastic mast-cell aggregates from reactive proliferations, and has, therefore, recently become a minor criterion for the diagnosis of Systemic Mastocytosis (SM). Anti-CD25 antibodies have also been useful in identifying mast cells in skin biopsies in the setting of Urticaria Pigmentosa, which is predictive of Systemic Mastocytosis. Quantitation of regulatory T-cells (Treg) in the setting of hepatocellular carcinoma has been used as an independent predictive factor for tumor recurrence after hepatic resection for HCC. Also, the percentage of tumor-infiltrating CD25+FOXP3+ regulatory T-cells among tumor cells, inside tumor parenchyma and at its periphery are significantly higher in recurrent Cutaneous Melanoma than in Non-recurrent Melanoma.

CD25 is the alpha chain of the IL-2 receptor. It is a type I transmembrane protein present on activated T cells, activated B cells, some thymocytes, myeloid precursors, and oligodendrocytes that associates with CD122 to form a heterodimer that can act as a high-affinity receptor for IL-2. Studies have shown that a large proportion of resting memory T cells constitutively express CD25.


CD25 is expressed in most B-cell neoplasms, some acute nonlymphocytic leukemias, neuroblastomas, and tumor infiltrating lymphocytes. Its soluble form, called sIL-2R may be elevated in these diseases and is occasionally used to track disease progression. CD25 is also utilized in cases of mastocytosis.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP75
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 4C9
ISOTYPE: IgG2B
CONTROL: Tonsil, Small Bowel, Colon, Spleen, Mastocytosis, Hodgkin's Lymphoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-CD25
ISOTYPE: IgG
CONTROL: Tonsil, Small Bowel, Colon, Spleen, Mastocytosis, Hodgkin's Lymphoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6471 | Tinto Predilute | 3.0 ml | BSB 6317 | Tinto Predilute | 3.0 ml | BSB 2454 | Tinto Predilute | 3.0 ml |
| BSB 6472 | Tinto Predilute | 7.0 ml | BSB 6318 | Tinto Predilute | 7.0 ml | BSB 2455 | Tinto Predilute | 7.0 ml |
| BSB 6473 | Tinto Predilute | 15.0 ml | BSB 6319 | Tinto Predilute | 15.0 ml | BSB 2456 | Tinto Predilute | 15.0 ml |
| BSB 6474 | Concentrate | 0.1 ml | BSB 6320 | Concentrate | 0.1 ml | BSB 2457 | Concentrate | 0.1 ml |
| BSB 6475 | Concentrate | 0.5 ml | BSB 6321 | Concentrate | 0.5 ml | BSB 2458 | Concentrate | 0.5 ml |
| BSB 6476 | Concentrate | 1.0 ml | BSB 6322 | Concentrate | 1.0 ml | BSB 2459 | Concentrate | 1.0 ml |
| BSB 6477 | Control Slides | 5 | BSB 6323 | Control Slides | 5 | BSB 2460 | Control Slides | 5 |

CD30, MMab

IHC of CD30 on a FFPE Hodgkin's Lymphoma Tissue

CD30 is a transmembrane cytokine receptor belonging to the tumor necrosis factor (TNF) receptor superfamily. Mature CD30 has a molecular mass of 120 kDa and is derived from a 90 kDa precursor protein.


CD30 antibody detects an epitope which is expressed by Reed-Sternberg cells in Hodgkin's Disease, the majority of Anaplastic Large-cell Lymphomas, and in Embryonal Carcinomas and Seminomas. This antibody also stains plasma cells intensely in paraffin-embedded tissue.

CD31, MMab

IHC of CD31 on a FFPE Placenta Tissue

CD31 is also called PECAM-1 for platelet endothelial cell-adhesion molecule. It plays a key role in removing aged neutrophils from the body. CD-31 is normally found on stem cells, endothelial cells, platelets, macrophages and Kupffer cells, granulocytes, T/NK cells, lymphocytes, megakaryocytes, fibroblasts, osteoclasts and neutrophils. CD-31 is also expressed in certain tumors, including Epithelioid Hemangioendothelioma, Epithelioid Sarcoma-like Hemangioendothelioma, other vascular tumors, Histiocytic malignancies, and Plasmacytomas. It is rarely found in some sarcomas and carcinomas. CD-31 and macrophages play a key role in tissue regeneration.

CD31 is widely used to identify the vascular origin of neoplasms, as it is a highly specific and sensitive marker for vascular endothelial cells.

CD33, R Mab

IHC of CD33 on a FFPE Colon Tissue

CD33 (gp62 or siglec-3) is a glycosylated transmembrane protein that is a member of the sialic acid-binding immunoglobulin-like lectin (siglec) family. The genomic locus of this protein has been mapped to chromosome 19q13.1-3.5. The function of CD33 is not known, but it may have a role in cell-to-cell adhesion. In maturing granulocytic cells, there is progressive down-regulation of CD33 from the blast stage to mature neutrophils. However, in monocytes and macrophages/histiocytes, strong expression of CD33 is maintained throughout maturation.

Detection of CD33 using monoclonal antibodies has been a critical component in immunophenotyping acute leukemias, particularly Acute Myeloid Leukemias.


ANTIBODY TYPE: Mouse Monoclonal**CLONE:** Ber-H2**ISOTYPE:** IgG1/K**CONTROL:** Tonsil, Lymph Node, Hodgkin's Lymphoma**LOCALIZATION:** Membranous**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Mouse Monoclonal**CLONE:** 1A10**ISOTYPE:** IgG1/K**CONTROL:** Tonsil, Placenta, Appendix, Spleen, Kidney**LOCALIZATION:** Cytoplasmic, Membranous**SPECIES REACTIVITY:** Human, Dog, Rat, Mouse**ANTIBODY TYPE:** Rabbit Monoclonal**CLONE:** RBT-CD33**ISOTYPE:** IgG**CONTROL:** Placenta, Myometrium, Lung, Colon, Spleen, Lymph Node, Tonsil, Acute Myeloid Leukemia**LOCALIZATION:** Membranous**SPECIES REACTIVITY:** Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5211 | Tinto Predilute | 3.0 ml |
| BSB 5212 | Tinto Predilute | 7.0 ml |
| BSB 5213 | Tinto Predilute | 15.0 ml |
| BSB 5214 | Concentrate | 0.1 ml |
| BSB 5215 | Concentrate | 0.5 ml |
| BSB 5216 | Concentrate | 1.0 ml |
| BSB 5217 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5218 | Tinto Predilute | 3.0 ml |
| BSB 5219 | Tinto Predilute | 7.0 ml |
| BSB 5220 | Tinto Predilute | 15.0 ml |
| BSB 5221 | Concentrate | 0.1 ml |
| BSB 5222 | Concentrate | 0.5 ml |
| BSB 5223 | Concentrate | 1.0 ml |
| BSB 5224 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3448 | Tinto Predilute | 3.0 ml |
| BSB 3449 | Tinto Predilute | 7.0 ml |
| BSB 3450 | Tinto Predilute | 15.0 ml |
| BSB 3451 | Concentrate | 0.1 ml |
| BSB 3452 | Concentrate | 0.5 ml |
| BSB 3453 | Concentrate | 1.0 ml |
| BSB 3454 | Control Slides | 5 |

CD34, MAb




IHC of CD34 on a FFPE Kidney Tissue

CD34, RMAb



IHC of CD34 on a FFPE Dermatofibrosarcoma Protuberans

CD35, MAb



IHC of CD35 on a FFPE Tonsil Tissue

CD34 functions as a cell-cell adhesion factor and cell-surface glycoprotein. It may also mediate the attachment of stem cells to bone marrow extracellular matrix or directly to stromal cells. Cells expressing CD34 are normally found in the umbilical cord and bone marrow as hematopoietic cells, and in vascular endothelium. In addition to stem cell recognition, CD34 is expressed by vascular endothelium; it appears that proliferating endothelial cells express this molecule in greater amounts than resting cells. In comparison to factor VIII R Antigen, CD34 stains are stronger and appear to be more sensitive in nature.

In tumors, CD34 is found in Alveolar Soft Part Sarcoma, pre B-ALL (positive in 75%), AML (40%), AMLM7 (most), Dermatofibrosarcoma Protuberans, Gastrointestinal Stromal Tumors, Giant Cell Fibroblastoma, Granulocytic Sarcoma, Kaposi's Sarcoma, Liposarcoma, Malignant Fibrous Histiocytoma, Malignant Peripheral Nerve Sheath tumors, Meningeal Hemangiopericytomas, Meningiomas, Neurofibromas, Schwannomas, and Papillary Thyroid Carcinoma. A negative CD34 may exclude Ewing's Sarcoma/PNET, Myofibrosarcoma of the breast, and Inflammatory Myofibroblastic tumors of the stomach.

CD34 functions as a cell-cell adhesion factor and cell-surface glycoprotein. It may also mediate the attachment of stem cells to bone marrow extracellular matrix or directly to stromal cells. Cells expressing CD34 are normally found in the umbilical cord and bone marrow as hematopoietic cells, and in vascular endothelium. In addition to stem cell recognition, CD34 is expressed by vascular endothelium; it appears that proliferating endothelial cells express this molecule in greater amounts than resting cells. In comparison to factor VIII R Antigen, CD34 stains are stronger and appear to be more sensitive in nature.

In tumors, CD34 is found in Alveolar Soft Part Sarcoma, pre B-ALL (positive in 75%), AML(40%), AMLM7 (most), Dermatofibrosarcoma Protuberans, Gastrointestinal Stromal Tumors, Giant Cell Fibroblastoma, Granulocytic Sarcoma, Kaposi's Sarcoma, Liposarcoma, Malignant Fibrous Histiocytoma, Malignant Peripheral Nerve Sheath tumors, Meningeal Hemangiopericytomas, Meningiomas, Neurofibromas, Schwannomas, and Papillary Thyroid Carcinoma. A negative CD34 may exclude Ewing's Sarcoma/PNET, Myofibrosarcoma of the breast, and Inflammatory Myofibroblastic tumors of the stomach.

CD35 (erythrocyte complement receptor 1 or CR1, also known as C3b/C4b receptor and immune adherence receptor) serves as the main system for processing and clearance of complement-opsonized immune complexes. The number of CR1 molecules decreases with aging of erythrocytes in normal individuals and is also decreased in pathological conditions such as Systemic Lupus Erythematosus (SLE), HIV infection, some Hemolytic Anemias and other conditions featuring immune complexes.

Anti-CD35 is considered a mature B-cell marker, which labels follicular dendritic reticulum cells and tumors derived from such cells such as Follicular Dendritic Cell Tumor/Sarcoma. CD35 antigen is found in erythrocytes, B-cells, and a subset of T-cells, monocytes, as well as in eosinophils and neutrophils.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: QBEnd/10
ISOTYPE: IgG1
CONTROL: Tonsil, Placenta, Appendix
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Dog

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP88
ISOTYPE: IgG
CONTROL: Tonsil, Placenta, Appendix
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse, Rat, Sheep, Dog, Pig, Loxodonta Africana


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-132
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5225 | Tinto Predilute | 3.0 ml |
| BSB 5226 | Tinto Predilute | 7.0 ml |
| BSB 5227 | Tinto Predilute | 15.0 ml |
| BSB 5228 | Concentrate | 0.1 ml |
| BSB 5229 | Concentrate | 0.5 ml |
| BSB 5230 | Concentrate | 1.0 ml |
| BSB 5231 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6485 | Tinto Predilute | 3.0 ml |
| BSB 6486 | Tinto Predilute | 7.0 ml |
| BSB 6487 | Tinto Predilute | 15.0 ml |
| BSB 6488 | Concentrate | 0.1 ml |
| BSB 6489 | Concentrate | 0.5 ml |
| BSB 6490 | Concentrate | 1.0 ml |
| BSB 6491 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-5237-3 | Tinto Predilute | 3.0 ml |
| BSB-5237-7 | Tinto Predilute | 7.0 ml |
| BSB-5237-15 | Tinto Predilute | 15.0 ml |
| BSB-5237-01 | Concentrate | 0.1 ml |
| BSB-5237-05 | Concentrate | 0.5 ml |
| BSB-5237-1 | Concentrate | 1.0 ml |
| BSB-5237-CS | Control Slides | 5 |

CD35, RMAb




IHC of CD35 on a FFPE Tonsil Tissue

CD38, MMAb



IHC of CD38 on a FFPE Tonsil Tissue

CD38, RMAb



IHC of CD38 on a FFPE Tonsil Tissue

CD35 (erythrocyte complement receptor 1 or CR1, also known as C3b/C4b receptor and immune adherence receptor) serves as the main system for processing and clearance of complement-opsonized immune complexes. The number of CR1 molecules decreases with aging of erythrocytes in normal individuals and is also decreased in pathological conditions such as Systemic Lupus Erythematosus (SLE), HIV infection, some Hemolytic Anemias and other conditions featuring immune complexes.

Anti-CD35 is considered a mature B-cell marker, which labels follicular dendritic reticulum cells and tumors derived from such cells such as Follicular Dendritic Cell Tumor/Sarcoma. CD35 antigen is found in erythrocytes, B-cells, and a subset of T-cells, monocytes, as well as in eosinophils and neutrophils.

CD38 is a glycoprotein found on the surface of many immune cells (white blood cells), including CD4+, CD8+, B and natural killer cells. It is a marker of cell activation. The CD38 protein has been connected to HIV infection, Leukemias, Myelomas, solid tumors, Type II Diabetes Mellitus and bone metabolism, as well as some genetically-determined conditions. It has also been used as a prognostic marker in Leukemia. CD38 is highly expressed on thymocytes. It is also expressed by early cells of B and T lineages, NK cells, plasma cells, monocytes and macrophages, and may be detected on cells from Multiple Myeloma, ALL (B and T) and some AML.

Monoclonal antibodies to CD38 have been shown to be useful in subtyping of Lymphomas and Leukemias, inhibition of B-lymphopoiesis, detection of plasma cells, protection of B-cells from apoptosis, and as a marker for activated B and T-cell proliferation.

CD38 is a glycoprotein found on the surface of many immune cells (white blood cells), including CD4+, CD8+, B and natural killer cells. It is a marker of cell activation. The CD38 protein has been linked to HIV infection, Leukemias, Myelomas, solid tumors, Type II Diabetes Mellitus and bone metabolism, as well as some genetically-determined conditions. It has also been used as a prognostic marker in Leukemia. CD38 is highly expressed on thymocytes. It is also expressed by early cells of B and T lineages, NK cells, plasma cells, monocytes and macrophages, and may be detected on cells from Multiple Myeloma, ALL (B and T) and some AML.

Monoclonal antibodies to CD38 have been shown to be useful in subtyping of Lymphomas and Leukemias, inhibition of B-lymphopoiesis, detection of plasma cells, protection of B-cells from apoptosis, and as a marker for activated B and T-cell proliferation.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP197
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: SPC32
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node, Spleen, Prostate, Salivary Gland
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human, Rabbit


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP135
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Breast
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6492 | Tinto Predilute | 3.0 ml |
| BSB 6493 | Tinto Predilute | 7.0 ml |
| BSB 6494 | Tinto Predilute | 15.0 ml |
| BSB 6495 | Concentrate | 0.1 ml |
| BSB 6496 | Concentrate | 0.5 ml |
| BSB 6497 | Concentrate | 1.0 ml |
| BSB 6498 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6198 | Tinto Predilute | 3.0 ml |
| BSB 6199 | Tinto Predilute | 7.0 ml |
| BSB 6200 | Tinto Predilute | 15.0 ml |
| BSB 6201 | Concentrate | 0.1 ml |
| BSB 6202 | Concentrate | 0.5 ml |
| BSB 6203 | Concentrate | 1.0 ml |
| BSB 6204 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6499 | Tinto Predilute | 3.0 ml |
| BSB 6500 | Tinto Predilute | 7.0 ml |
| BSB 6501 | Tinto Predilute | 15.0 ml |
| BSB 6502 | Concentrate | 0.1 ml |
| BSB 6503 | Concentrate | 0.5 ml |
| BSB 6504 | Concentrate | 1.0 ml |
| BSB 6505 | Control Slides | 5 |

CD41/Integrin alpha 11b, RMab




IHC of CD41/Integrin alpha 11b on a FFPE Bone Marrow Tissue

CD42b, RMab



IHC of CD42b on a FFPE Bone Marrow Tissue

CD43, MMab



IHC of CD43 on a FFPE Tonsil Tissue

ITGA2B encodes CD41, or integrin alpha 11b. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. Alpha chain 11b undergoes post-translational cleavage to yield disulfide-linked light and heavy chains that join with beta 3 to form a fibrinogen receptor expressed in platelets that plays a crucial role in coagulation. Mutations that interfere with this role result in thrombasthenia. In addition to adhesion, integrins are known to participate in cell-surface mediated signalling.

CD41 expression has been found on platelets, megakaryocytes, and immature hematopoietic progenitors.

CD42b Platelet glycoprotein Ib alpha chain also known as glycoprotein Ib (platelet), alpha polypeptide is a protein that in humans is encoded by the GP1BA gene. The Gp Ib functions as a receptor for von Willebrand factor (VWF), which includes noncovalent association of the alpha and beta subunits with platelet glycoprotein IX and platelet glycoprotein V. The binding of the GP Ib-IX-V complex to VWF facilitates initial platelet adhesion to vascular subendothelium after vascular injury, and also initiates signaling events within the platelet that lead to enhanced platelet activation, thrombosis, and hemostasis. This gene encodes the alpha subunit. Several polymorphisms and mutations have been described in this gene, some of which are the cause of Bernard-Soulier syndromes and platelet-type von Willebrand disease (1).

CD42b is a marker of megakaryocytes and platelets and used in the Diagnosis of AML-M7 to distinguish AML-M7 (CD42b+) from acute myeloid leukemia with myelofibrosis (usually CD42b -). Acute megakaryoblastic leukemia has a significantly higher proportion of blasts expressing CD42b (2). Patients with Histoplasma capsulatum infections have been noted with unusual IHC positivity of the platelet associated marker CD42b/GP1b expressed on the surface of this fungal organisms. Evaluation of additional cases demonstrated that a majority of histoplasmosis cases (83%) showed positive staining with CD42b/GP1b, comparable to GMS stain results (3).

CD43 (leukosialin, sialophorin, or leukocyte sialoglycoprotein) is one of the major glycoproteins expressed in all thymocytes and T-cells. It plays a role in the physiochemical properties of the T-cell surface and in lectin binding. During T-cell activation, CD43 is actively removed from the T-cell antigen-presenting cell contact site, suggesting a negative regulatory role in adaptive immune response.

This antibody has been found useful in identification and classification of T-cell malignancies and low grade B-cell Lymphomas. CD43 expression is seen in some cases of B-cell Lymphocytic Lymphoma and Centrocytic Lymphoma. When used in combination with CD45 and CD20, effective immunophenotyping of the majority of Lymphomas can be obtained. Co-staining of a lymphoid infiltrate with CD20 and CD3 argues against a reactive process and favors Lymphoma.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP178
ISOTYPE: IgG
CONTROL: Spleen, Bone Marrow
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP409
ISOTYPE: IgG
CONTROL: Bone Marrow, Spleen
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: MT1
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6506 | Tinto Predilute | 3.0 ml | BSB 3504 | Tinto Predilute | 3.0 ml | BSB 5239 | Tinto Predilute | 3.0 ml |
| BSB 6507 | Tinto Predilute | 7.0 ml | BSB 3505 | Tinto Predilute | 7.0 ml | BSB 5240 | Tinto Predilute | 7.0 ml |
| BSB 6508 | Tinto Predilute | 15.0 ml | BSB 3506 | Tinto Predilute | 15.0 ml | BSB 5241 | Tinto Predilute | 15.0 ml |
| BSB 6509 | Concentrate | 0.1 ml | BSB 3507 | Concentrate | 0.1 ml | BSB 5242 | Concentrate | 0.1 ml |
| BSB 6510 | Concentrate | 0.5 ml | BSB 3508 | Concentrate | 0.5 ml | BSB 5243 | Concentrate | 0.5 ml |
| BSB 6511 | Concentrate | 1.0 ml | BSB 3509 | Concentrate | 1.0 ml | BSB 5244 | Concentrate | 1.0 ml |
| BSB 6512 | Control Slides | 5 | BSB 3510 | Control Slides | 5 | BSB 5245 | Control Slides | 5 |

CD44, MAb



IHC of CD44 on a FFPE Breast Fibroadenoma Tissue

The CD44 protein is a cell-surface glycoprotein involved in cell-cell interactions, cell adhesion and migration. CD44 is also known as Homing-cell adhesion molecule (H-CAM) and Phagocytic glycoprotein-1 (PgP-1). A specialized sialofucosylated glycoform of CD44 called HCELL is found natively on human hematopoietic stem cells and functions as a "bone-homing receptor", directing migration of human hematopoietic stem cells and mesenchymal stem cells to bone marrow.

This protein participates in a wide variety of cellular functions including lymphocyte activation, recirculation and homing, hematopoiesis, and tumor metastasis. Transcripts for this gene undergo complex alternative splicing that results in many functionally distinct isoforms; however, the full-length nature of some of these variants has not been determined. Splice variants of CD44 on Colon Cancer cells display the HCELL glycoform, which mediates binding to vascular E-selectin under hemodynamic flow conditions, a critical step in Colon Cancer metastasis. In addition, variations in CD44 are reported as cell surface markers for some breast and prostate cancer stem cells and have been seen as an indicator of increased survival time in Epithelial Ovarian Cancer patients.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-12


ISOTYPE: IgG2A

CONTROL: Urothelium, Tonsil, Kidney, Breast, Liver, Skin, Prostate, Thymus, Spleen, Lymph Node, Esophageal Carcinoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

CD45, MAb



IHC of CD45 on a FFPE Tonsil Tissue

The CD45 antigen is a protein which was originally called Leukocyte Common Antigen. It is a Type I transmembrane protein which is in various forms present on all differentiated hematopoietic cells except erythrocytes and assists in the activation of those cells (a form of co-stimulation). It is expressed in Lymphomas, B-cell Chronic Lymphocytic Leukemia, Hairy Cell Leukemia, and Acute Non-Lymphocytic Leukemia.

CD45 is a monoclonal antibody that is routinely used to aid in the differential diagnosis of undifferentiated neoplasms, whenever malignant Lymphoma is suspected by the morphological or clinical data. It is a highly specific antibody; thus, a positive result is highly indicative of lymphoid or myeloid origin. Certain types of lymphoid neoplasms may lack CD45 (Hodgkin's Disease, some T-cell Lymphomas and some Leukemias) so its absence does not rule out a hematolymphoid tumor. This antibody is exclusively expressed by cells of hematopoietic lineage and is present in most benign and malignant lymphocytes, erythrocytes and plasma cell precursors.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 2B11 & PD7/26


ISOTYPE: IgG1/K

CONTROL: Tonsil, Lymph Node, Thymus, Spleen

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

CD45R, MAb



IHC of CD45R on a FFPE Tonsil Tissue

CD45R contains an extracellular domain, a single transmembrane segment and two tandem intracytoplasmic catalytic domains. It is specifically expressed in hematopoietic cells and has been shown to be an essential regulator of T and B-cell antigen-receptor signaling. It functions through either direct interaction with components of the antigen receptor complexes, or by activating various Src family kinases required for the antigen-receptor signaling. CD45R also suppresses JAK kinases, and thus functions as a regulator of cytokine-receptor signaling.

CD45R represents a restricted form of the CD45 family, which primarily recognizes only cells of B lineage from proB-cell through mature B lymphocytes and, prior to the availability of anti-CD19 MAbs, was commonly used as a pan B-cell marker. It also reacts with certain activated T-cells, as well as non-MHC restricted lytically active lymphokine-activated killer (LAK) cells. MB1 antibody stains preferentially B-cells and their neoplasms but is less specific, as it will also react with some T-cell Lymphomas and Non-lymphoid Tumors. The antigen for this antibody is in the membrane of all B-cells with the exception of plasma cells and some mature T-cells.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: MB1

ISOTYPE: IgG1

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6233 | Tinto Predilute | 3.0 ml |
| BSB 6234 | Tinto Predilute | 7.0 ml |
| BSB 6235 | Tinto Predilute | 15.0 ml |
| BSB 6236 | Concentrate | 0.1 ml |
| BSB 6237 | Concentrate | 0.5 ml |
| BSB 6238 | Concentrate | 1.0 ml |
| BSB 6239 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5246 | Tinto Predilute | 3.0 ml |
| BSB 5247 | Tinto Predilute | 7.0 ml |
| BSB 5248 | Tinto Predilute | 15.0 ml |
| BSB 5249 | Concentrate | 0.1 ml |
| BSB 5250 | Concentrate | 0.5 ml |
| BSB 5251 | Concentrate | 1.0 ml |
| BSB 5252 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6254 | Tinto Predilute | 3.0 ml |
| BSB 6255 | Tinto Predilute | 7.0 ml |
| BSB 6256 | Tinto Predilute | 15.0 ml |
| BSB 6257 | Concentrate | 0.1 ml |
| BSB 6258 | Concentrate | 0.5 ml |
| BSB 6259 | Concentrate | 1.0 ml |
| BSB 6260 | Control Slides | 5 |

CD45RA, MMab




IHC of CD45RA on a FFPE Tonsil Tissue

CD45RO, MMab



IHC of CD45RO on a FFPE Colon Tissue

CD56, MMab



IHC of CD56 on a FFPE Neuroblastoma Tissue

CD45 is a complex molecule and is comprised of different glycoproteins ranging from 180-240 kDa. Expression of CD45 is found on all hemopoietic cells. Detection of the different isoforms can distinguish between different cell forms (e.g., naive T-cells and memory T-cells). CD45RA is an isoform of the CD45 complex and has restricted expression between different subtypes of lymphoid cells.

CD45RA antibody reacts with mature, non-activated T and B-cells. CD45RA is also reactive with medullary thymocytes, mantle-zone lymphocytes in follicles of lymph nodes, spleen and lymphocytes of the paracortex. CD45RA shows no reactivity with cortical thymocytes, immature T-cells or activated B-cells in germinal centers.

The CD45 family consists of multiple members that are all products of a single complex gene. Three isoforms of CD45 exist: on B-lymphocytes, where the protein is called B220 (its molecular mass is 220 kDa); on naive T-lymphocytes, where it is called CD45 RA, and on activated and memory T-lymphocytes, where it is called CD45 RO. CD45RO is a single-chain, transmembraneous glycoprotein which represents the low molecular weight isoform of the Leukocyte Common Antigen (LCA). It is expressed on most thymocytes, about 45% of peripheral blood T-cells, virtually all T-cells in skin reactive infiltrates, and the majority of T-cell malignancies. It is also found on a subset of B-cells and on exceptional B-cell Lymphomas.

CD45RO (T-Cell, Pan) antibody reacts with thymocytes and activated T-cells, but only on a subpopulation of resting T-cells. This antibody shows no reactivity with B-cells, making it a good marker for T-cell tumors to be phenotyped. In addition, granulocytes and monocytes are also labeled with this antibody. T-Cell, Pan has been designated as CD45RO at The International Leukocyte Typing Workshop.

CD56 or Neural-Cell Adhesion Molecule (NCAM) is a homophilic binding glycoprotein expressed on the surface of neurons, glia and skeletal muscle. CD56 has been implicated in cell-cell adhesion, neurite outgrowth, synaptic plasticity, and learning and memory.

Normal cells that stain positively for CD56 include NK cells, activated T-cells, brain and cerebellum, and neuroendocrine tissues. Tumors that are CD56-positive are Myeloma, Myeloid Leukemia, Neuroendocrine tumors, Wilm's Tumor, Adult Neuroblastoma, NK/T cell Lymphomas, Pancreatic Acinar-cell Carcinoma, Pheochromocytoma, and Small-cell Lung Carcinoma. It is also expressed on some mesodermally-derived tumors (Rhabdomyosarcoma). Ewing's Sarcoma/PNET is CD56-negative.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 4KB5
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: UCHL-1
ISOTYPE: IgG2a/K
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human, Mouse, Rat, Non-human Primate

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 123C3.D5
ISOTYPE: IgG1/K
CONTROL: Pancreas, Tonsil, Neuroblastoma, Brain, Thyroid, Prostate, Colon, Lung SCC
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5253 | Tinto Predilute | 3.0 ml |
| BSB 5254 | Tinto Predilute | 7.0 ml |
| BSB 5255 | Tinto Predilute | 15.0 ml |
| BSB 5256 | Concentrate | 0.1 ml |
| BSB 5257 | Concentrate | 0.5 ml |
| BSB 5258 | Concentrate | 1.0 ml |
| BSB 5259 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5260 | Tinto Predilute | 3.0 ml |
| BSB 5261 | Tinto Predilute | 7.0 ml |
| BSB 5262 | Tinto Predilute | 15.0 ml |
| BSB 5263 | Concentrate | 0.1 ml |
| BSB 5264 | Concentrate | 0.5 ml |
| BSB 5265 | Concentrate | 1.0 ml |
| BSB 5266 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5267 | Tinto Predilute | 3.0 ml |
| BSB 5268 | Tinto Predilute | 7.0 ml |
| BSB 5269 | Tinto Predilute | 15.0 ml |
| BSB 5270 | Concentrate | 0.1 ml |
| BSB 5271 | Concentrate | 0.5 ml |
| BSB 5272 | Concentrate | 1.0 ml |
| BSB 5273 | Control Slides | 5 |

CD57, MAb

IHC of CD57 on a FFPE Tonsil Tissue

CD57 (NK-1) recognizes an oligosaccharide (MW 100-110 kDa) antigenic determinant on myeloid cells and on a variety of polypeptides, lipids and chondroitin sulfate proteoglycans. This surface antigen is associated with myelin-associated glycoprotein (MAG). The CD57 antigen is present on 15-20% of normal peripheral blood mononuclear cells. It is expressed on a subset of natural killer cells (60%) and on a subset of T-lymphocytes. This carbohydrate is also present on N-CAM in the nervous system.


Follicular Center-cell Lymphomas often contain many NK cells within the neoplastic follicles. NK-1 reportedly also reacts with a variety of cell types in non-lymphoid tissues. NK-1 stains neuroendocrine cells and their tumors, including Carcinoid Tumor and Medulloblastomas. NK-1 also reacts with a variety of cell types in non-lymphoid tissues, including Neurofibroma, Ganglioneuroma, and Prostate Carcinoma.

CD61, MAb

IHC of CD61 on a FFPE Bone Marrow Tissue

CD61 is a glycoprotein found on megakaryocytes (bone marrow cells), platelets and their precursors. CD61 antigen plays a role in platelet aggregation and also as a receptor for fibrinogen, fibronectin, von Willebrand factor and vitronectin.

CD61 labels the IIIa subunit of the noncovalently-linked glycoprotein heterodimer IIb/IIIa complex present on human platelets and their precursors. This antibody is useful in identifying megakaryoblastic differentiation as seen in Megakaryoblastic Leukemia.

CD61, RMAb

IHC of CD61 on an FFPE Kidney Tissue

CD61 is a glycoprotein found on megakaryocytes (bone marrow cells), platelets and their precursors. CD61 antigen plays a role in platelet aggregation and also as a receptor for fibrinogen, fibronectin, von Willebrand factor and vitronectin.

CD61 labels the IIIa subunit of the noncovalently-linked glycoprotein heterodimer IIb/IIIa complex present on human platelets and their precursors. This antibody is useful in identifying megakaryoblastic differentiation as seen in Megakaryoblastic Leukemia.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-10

ISOTYPE: IgM/K

CONTROL: Tonsil, Lymph Node, Spleen, Prostate, Breast, Brain

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 2F2

ISOTYPE: IgG1/K

CONTROL: Bone Marrow

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP65

ISOTYPE: IgG

CONTROL: Brain, Kidney, Testis, Bone Marrow

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5274 | Tinto Predilute | 3.0 ml |
| BSB 5275 | Tinto Predilute | 7.0 ml |
| BSB 5276 | Tinto Predilute | 15.0 ml |
| BSB 5277 | Concentrate | 0.1 ml |
| BSB 5278 | Concentrate | 0.5 ml |
| BSB 5279 | Concentrate | 1.0 ml |
| BSB 5280 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5281 | Tinto Predilute | 3.0 ml |
| BSB 5282 | Tinto Predilute | 7.0 ml |
| BSB 5283 | Tinto Predilute | 15.0 ml |
| BSB 5284 | Concentrate | 0.1 ml |
| BSB 5285 | Concentrate | 0.5 ml |
| BSB 5286 | Concentrate | 1.0 ml |
| BSB 5287 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3511 | Tinto Predilute | 3.0 ml |
| BSB 3512 | Tinto Predilute | 7.0 ml |
| BSB 3513 | Tinto Predilute | 15.0 ml |
| BSB 3514 | Concentrate | 0.1 ml |
| BSB 3515 | Concentrate | 0.5 ml |
| BSB 3516 | Concentrate | 1.0 ml |
| BSB 3517 | Control Slides | 5 |

CD63, MAb




IHC of CD63 on a FFPE Melanoma Tissue

CD68, MAb



IHC of CD68 on a FFPE Tonsil Tissue

CD68, MAb



IHC of CD68 on a FFPE Lymphoblastic Lymphoma Tissue

The protein encoded by CD63 gene is a member of the transmembrane-4 superfamily, also known as the tetraspanin family, and mediates signal-transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is a cell-surface glycoprotein that is known to complex with integrins. It may function as a blood-platelet activation marker. Deficiency of this protein is associated with Hermansky-Pudlak Syndrome. This gene has been associated with tumor progression. CD63 is a good marker for flow-cytometric quantification of in vitro-activated basophils for diagnosis of IgE-mediated allergy. The test is commonly designated as a basophil activation test.

Anti-CD63 reacts with a 53 kDa protein. The antigen was originally designated as a lysosomal membrane protein characterized as an activation-dependent platelet surface antigen. In fact, the CD63 antigen has a diverse distribution on the surface and in the cytoplasm of many cell types including lymphoid, myeloid and endothelial cells and Melanoma. It is weakly expressed in granulocytes, B and T-cells. It has been quite useful in identifying Malignant Melanoma. CD63 is thought to be associated with the early stages of Melanoma tumor progression (in regulation of motility and adhesion of Melanoma cells).

The CD68 antigen is a heavily glycosylated transmembrane protein of 87-115 kDa which is specifically expressed by tissue macrophages, Langerhans cells and, at low levels, by dendritic cells. CD68 could play a role in phagocytic activities of tissue macrophages, both in intracellular lysosomal metabolism and extracellular cell-cell and cell-pathogen interactions.

CD68 marks cells of monocyte/macrophage lineage. This antibody is capable of staining monocytes, Kupffer cells, osteoclasts, granulocytes and their precursors; Lymphomas are negative or show a few granules. This antibody may be useful for the identification of Myelomonocytic and Histiocytic Tumors. CD68 may help to distinguish Malignant Fibrous Histiocytoma from other Pleomorphic Sarcomas. However, since CD68 detects a formalin-resistant epitope that may be associated with lysosomal granules, other lysosome-rich cells may also produce positive results.

The CD68 antigen is a heavily glycosylated transmembrane protein of 87-115 kDa which is specifically expressed by tissue macrophages, Langerhans cells and, at low levels, by dendritic cells. CD68 could play a role in phagocytic activities of tissue macrophages, both in intracellular lysosomal metabolism and extracellular cell-cell and cell-pathogen interactions.

CD68 marks cells of monocyte/macrophage lineage. This antibody is capable of staining monocytes, Kupffer cells, osteoclasts, granulocytes and their precursors; Lymphomas are negative or show a few granules. This antibody may be useful for the identification of Myelomonocytic and Histiocytic Tumors. CD68 may help to distinguish Malignant Fibrous Histiocytoma from other Pleomorphic Sarcomas. However, since CD68 detects a formalin-resistant epitope that may be associated with lysosomal granules, other lysosome-rich cells may also produce positive results.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: NK1/C3
ISOTYPE: IgG1/K
CONTROL: Skin, Malignant Melanoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-8
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse, Rabbit, Rat

ANTIBODY TYPE: Mouse Monoclonal
CLONE: KP-1
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Hamster, Mouse, Non-Human Primate, Porcine, Rabbit, Rat, Cat, Monkey

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6296 | Tinto Predilute | 3.0 ml |
| BSB 6297 | Tinto Predilute | 7.0 ml |
| BSB 6298 | Tinto Predilute | 15.0 ml |
| BSB 6299 | Concentrate | 0.1 ml |
| BSB 6300 | Concentrate | 0.5 ml |
| BSB 6301 | Concentrate | 1.0 ml |
| BSB 6302 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5288 | Tinto Predilute | 3.0 ml |
| BSB 5289 | Tinto Predilute | 7.0 ml |
| BSB 5290 | Tinto Predilute | 15.0 ml |
| BSB 5291 | Concentrate | 0.1 ml |
| BSB 5292 | Concentrate | 0.5 ml |
| BSB 5293 | Concentrate | 1.0 ml |
| BSB 5294 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2712 | Tinto Predilute | 3.0 ml |
| BSB 2713 | Tinto Predilute | 7.0 ml |
| BSB 2714 | Tinto Predilute | 15.0 ml |
| BSB 2715 | Concentrate | 0.1 ml |
| BSB 2716 | Concentrate | 0.5 ml |
| BSB 2717 | Concentrate | 1.0 ml |
| BSB 2718 | Control Slides | 5 |

CD71, MAb

IHC of CD71 on a FFPE Bone Marrow Tissue

CD71, also known as Transferrin Receptor Protein 1 (TfR1) is a protein encoded by the TFRC gene. CD71 is required for iron delivery from transferrin to cells. It is most highly expressed on placental syncytiotrophoblasts, myocytes, basal keratinocytes, hepatocytes, endocrine pancreas, spermatocytes, and erythroid precursors. The level of transferrin receptor expression is highest in the early erythroid precursors through intermediate normoblast phase, after which expression decreases through the reticulocyte phase.


The high level of CD71 within erythroid precursors makes it an excellent marker for erythroid components within bone marrow biopsy specimens without interference from mature erythrocytes. It may also be used in the determination of erythroid leukemia, benign erythroid proliferative disorders, and myelodysplastic syndrome.

CD73/NT5E, RMAb

IHC of CD73/NT5E on a FFPE Papillary Thyroid Carcinoma Tissue

Ecto-5'-nucleotidase (NT5E), or also known as CD73, encoded by the NT5E gene on the human chromosome 6q14.3, is a dimer of two identical 70 kDa subunits. As a plasma membrane protein and expressed in various cell types, CD73 hydrolyzes extracellular adenosine monophosphate (AMP) into adenosine and inorganic phosphate. It performs homeostatic functions and plays an important role in immunity and inflammation.

In a study that analyzed several publications conducting IHC analysis, CD73 overexpression was found to be associated with several types of human cancers including Bladder, Brain, Invasive Lobular Breast, Esophageal, Gastric, Pancreatic, Rectal mucinous, Renal Cell, Lung Large Cell, oral cavity Squamous Cell, Melanoma, and Lung Adenocarcinoma Cancers. Studies also found that lymph node metastases correlated with high CD73 expression. One study found that CD73 expression is frequently higher in NSCLC tissue, compared to normal tissue, which indicates the prognostic value of CD73 in the tumorigenesis of NSCLC. Another study revealed the prognostic potential of CD73 as a biomarker in Breast Cancer, since CD73 is associated with cancer cell invasion, migration and lymph node metastasis. Other studies showed a positive correlation between CD73 expression with Colorectal Cancer malignancy, and metastatic Melanomas.

CD74, MAb

IHC of CD74 on a FFPE Tonsil Tissue

CD74, also known as the MHC Class II-associated invariant chain (II), is a Type II transmembrane protein which binds to the peptide-binding groove of newly-synthesized MHC class II alpha/beta heterodimers and prevents their premature association with endogenous polypeptides. CD74 is expressed primarily by antigen-presenting cells such as B-lymphocytes (from before the pre-B-cell stage to before the plasma-cell stage), macrophages and monocytes, together with many epithelial cells.

CD74 stains predominantly germinal-center lymphocytes and B-cell lymphomas but rarely T-cell lymphomas. It stains the cell membrane but a paranuclear globular labeling is also noted. CD74 is useful in differentiating Atypical Fibroxanthoma from Malignant Fibrous Histiocytoma, as well as Small-cell Lung Carcinoma from Non-small cell Lung Carcinomas.


ANTIBODY TYPE: Mouse Monoclonal**CLONE:** 10F11**ISOTYPE:** IgG2b**CONTROL:** Bone Marrow, Placenta, Adrenal, Tonsil, Skin**LOCALIZATION:** Cytoplasmic, Membranous**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Rabbit Monoclonal**CLONE:** RM431**ISOTYPE:** IgG**CONTROL:** Placenta, Adrenal Gland, Liver, Testis, Transitional Cell Carcinoma, Ovarian Serous Carcinoma**LOCALIZATION:** Membranous, Cytoplasmic**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Mouse Monoclonal**CLONE:** LN2**ISOTYPE:** IgG1/K**CONTROL:** Tonsil, Lymph Node, Thymus**LOCALIZATION:** Cytoplasmic, Membranous**SPECIES REACTIVITY:** Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6520 | Tinto Predilute | 3.0 ml |
| BSB 6521 | Tinto Predilute | 7.0 ml |
| BSB 6522 | Tinto Predilute | 15.0 ml |
| BSB 6523 | Concentrate | 0.1 ml |
| BSB 6524 | Concentrate | 0.5 ml |
| BSB 6525 | Concentrate | 1.0 ml |
| BSB 6526 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3716-3 | Tinto Predilute | 3.0 ml |
| BSB-3716-7 | Tinto Predilute | 7.0 ml |
| BSB-3716-15 | Tinto Predilute | 15.0 ml |
| BSB-3716-01 | Concentrate | 0.1 ml |
| BSB-3716-05 | Concentrate | 0.5 ml |
| BSB-3716-1 | Concentrate | 1.0 ml |
| BSB-3716-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5295 | Tinto Predilute | 3.0 ml |
| BSB 5296 | Tinto Predilute | 7.0 ml |
| BSB 5297 | Tinto Predilute | 15.0 ml |
| BSB 5298 | Concentrate | 0.1 ml |
| BSB 5299 | Concentrate | 0.5 ml |
| BSB 5300 | Concentrate | 1.0 ml |
| BSB 5301 | Control Slides | 5 |

CD75, MMab




IHC of CD75 on a FFPE Liver Tissue

CD79a, MMab



IHC of CD79a on a FFPE Colon Tissue

CD99, MMab



IHC of CD99 on a FFPE Thymus Tissue

CD75, also known as CDw75 or Beta-galactoside alpha-2,6-sialyltransferase 1 (ST6GAL-1) is an enzyme that in humans is encoded by the ST6GAL1 gene. The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein, which is normally found in the Golgi but which can be proteolytically processed to a soluble form, is involved in the generation of the cell-surface carbohydrate determinants and differentiation antigens HB-6, CDw75, and CD76.

CD75 is an antibody which labels germinal centre B-cells as well as epithelial cells, however is absent in plasma cells. It reacts with RBC precursors of bone marrow, ductal and ciliated epithelial cells of kidney, breast, prostate, pancreas, lung, and with glioblastomas, astrocytomas, and Reed Sternberg cells in lymphocyte predominant Hodgkin's disease. It is shown to be a helpful antibody for ascribing a B-cell phenotype in known lymphoid tissues.

CD79a is non-covalently associated with membrane-bound immunoglobulins on B-cells to constitute the B-cell Ag receptor. CD79a first appears at pre B-cell stage and persists until the plasma-cell stage, where it is found as an intracellular component. CD79a is found in the majority of Acute Leukemias of precursor B-cell type, in B-cell lines, B-cell Lymphomas, and in some Myelomas.

CD79a is a B-cell marker that is generally used to complement CD20. This antibody will stain many of the same Lymphomas as CD20, but also stains more B-precursor Lymphoid Leukemias than CD20. CD79a also stains more cases of Plasma-cell Myeloma and occasionally some types of endothelial cells as well. CD79a will stain many cases of Acute Promyelocytic Leukemia (FAB-M3), but only rarely stains other types of Myeloid Leukemia.

CD99, also known as MIC-2 or single-chain Type-1 glycoprotein, is a human protein encoded by the CD99 gene. The protein has a MW of 32 kD. It is expressed on all leukocytes but highest on thymocytes, and is believed to augment T-cell adhesion and apoptosis of double-positive T-cells. It also participates in migration and activation.

The CD99 antigen is found on the cell membrane of Ewing's Sarcoma and Primitive Peripheral Neuroectodermal Tumors (PNET). It is also present on a variety of other cell types including bone marrow, lymph nodes, spleen, cortical thymocytes, granulosa cells of the ovary, beta cells, CNS ependymal cells, Sertoli's cells of the testis and a few endothelial cells. Mature granulocytes, however, tend to express very little or no CD99. MIC-2 has also been identified in Lymphoblastic Lymphoma, Rhabdomyosarcoma, Mesenchymal Chondrosarcoma, and Thymoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: LN-1
ISOTYPE: IgM/K
CONTROL: Liver, Prostate, Kidney, Tonsil
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: JCB117
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node, Colon, Spleen
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human, Dog, Cat


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-9
ISOTYPE: IgG1/K
CONTROL: Pancreas, Thymus, Ependyma, Ewing's Sarcoma & Soft Tissue
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2719 | Tinto Predilute | 3.0 ml |
| BSB 2720 | Tinto Predilute | 7.0 ml |
| BSB 2721 | Tinto Predilute | 15.0 ml |
| BSB 2722 | Concentrate | 0.1 ml |
| BSB 2723 | Concentrate | 0.5 ml |
| BSB 2724 | Concentrate | 1.0 ml |
| BSB 2725 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5302 | Tinto Predilute | 3.0 ml |
| BSB 5303 | Tinto Predilute | 7.0 ml |
| BSB 5304 | Tinto Predilute | 15.0 ml |
| BSB 5305 | Concentrate | 0.1 ml |
| BSB 5306 | Concentrate | 0.5 ml |
| BSB 5307 | Concentrate | 1.0 ml |
| BSB 5308 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5309 | Tinto Predilute | 3.0 ml |
| BSB 5310 | Tinto Predilute | 7.0 ml |
| BSB 5311 | Tinto Predilute | 15.0 ml |
| BSB 5312 | Concentrate | 0.1 ml |
| BSB 5313 | Concentrate | 0.5 ml |
| BSB 5314 | Concentrate | 1.0 ml |
| BSB 5315 | Control Slides | 5 |

CD103/ITGAE, RMAb




IHC of CD103 on a FFPE Colon Tissue

CD105, RMAb



IHC of CD105 on a FFPE Testis Tissue

CD117, RMAb



IHC of CD117 on a FFPE GIST Tissue

Integrin, alpha E (ITGAE), also known as CD103, is an integrin protein that in human is encoded by the ITGAE gene. CD103 binds integrin beta 7 ($\beta 7$ -ITGB7) to form the complete heterodimeric integrin molecule $\alpha\beta 7$. CD103 is expressed widely on intraepithelial lymphocyte (IEL) T cells (both $\alpha\beta$ T cells and $\gamma\delta$ T cells) and on some peripheral regulatory T cells (Tregs). It has also been reported on lamina propria T cells. A subset of dendritic cells in the gut mucosa and in mesenteric lymph nodes also expresses this marker and is known as CD103 DCs. The chief ligand for $\alpha\beta 7$ is E-cadherin, an adhesion molecule found on epithelial cells.

CD103 has been found in mononuclear cells in the interfollicular area of lymph nodes and in intraepithelial cells in the overlying mucosa located primarily toward the basal layer of the tonsil. CD103 is useful in identifying Hairy Cell Leukemia, which is positive for this marker in most cases in contrast to other hematologic malignancies which are negative for CD103, with the exception of Splenic Marginal Zone Lymphoma, which rarely expresses CD103. The high sensitivity of anti-CD103 for Hairy Cell Leukemia makes this marker valuable when distinguishing this malignancy from other B-cell neoplasms.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP206

ISOTYPE: IgG

CONTROL: Skin, colon, Tonsil, Thymus, Spleen, Hairy Cell Leukemia

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

CD105/Endoglin is a Type I membrane glycoprotein located on cell surfaces and is part of the TGF beta receptor complex. This protein has been found on endothelial cells, activated macrophages, fibroblasts, and smooth-muscle cells. Endoglin has a role in the development of the cardiovascular system and in vascular remodeling. Its expression is regulated during heart development. In humans, Endoglin may be involved in the autosomal dominant disorder known as Hereditary Hemorrhagic Telangiectasia Type 1.

CD105 is highly expressed in endothelial cells during tumor angiogenesis and inflammation, with weak or negative expression in vascular endothelium of normal tissues. Angiogenesis is controlled by angiogenic factors, mostly secreted by tumor cells. Vascular Endothelial Growth Factor (VEGF) is a potent angiogenic growth factor that stimulates endothelial cell proliferation and induces microvessel permeability. Studies have demonstrated a correlation between VEGF expression and vascular density. Angiogenesis has been proposed as a promising prognostic marker in a variety of tumors. Most studies of angiogenesis have been done with pan-endothelial markers such as CD31 or CD34. Endoglin is a more specific and sensitive marker for tumor angiogenesis than CD31, as it labels only newly-formed blood vessels and may serve as a prognostic marker for Prostate Adenocarcinoma, and cancers of the lung, stomach, breast, and brain. CD105 may serve as a target for anti-angiogenesis therapy.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP274

ISOTYPE: IgG

CONTROL: Spleen, Tonsil, Cervix, Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

CD117 is a tyrosine-kinase receptor for stem cell factor (SCF), also known as "steel factor" or "c-kit ligand". C-kit is a polypeptide that activates bone marrow precursors of a number of blood cells, but its receptor is also present in other cells. C-kit mutations in the interstitial cells of Cajal in the digestive tract are probably the key to Gastrointestinal Stromal Tumors (GISTs).

CD117 is found on interstitial cells of Cajal, germ cells, bone marrow stem cells, melanocytes, breast epithelium and mast cells. This receptor is found on a wide variety of tumor cells (Follicular and Papillary Carcinoma of the Thyroid, Adenocarcinomas from endometrium, lung, ovary, pancreas, breast; Malignant Melanoma, Endodermal Sinus Tumor, Small-cell Carcinoma) but has been particularly useful in differentiating Gastrointestinal Stromal Tumors (GIST) from Kaposi's Sarcoma and tumors of smooth-muscle origin.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP10

ISOTYPE: IgG


CONTROL: Skin, Testis, Breast, GIST, Colon, Brain, Tonsil

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Monkey, Marmoset

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2859 | Tinto Predilute | 3.0 ml | BSB 2866 | Tinto Predilute | 3.0 ml | BSB 5316 | Tinto Predilute | 3.0 ml |
| BSB 2860 | Tinto Predilute | 7.0 ml | BSB 2867 | Tinto Predilute | 7.0 ml | BSB 5317 | Tinto Predilute | 7.0 ml |
| BSB 2861 | Tinto Predilute | 15.0 ml | BSB 2868 | Tinto Predilute | 15.0 ml | BSB 5318 | Tinto Predilute | 15.0 ml |
| BSB 2862 | Concentrate | 0.1 ml | BSB 2869 | Concentrate | 0.1 ml | BSB 5319 | Concentrate | 0.1 ml |
| BSB 2863 | Concentrate | 0.5 ml | BSB 2870 | Concentrate | 0.5 ml | BSB 5320 | Concentrate | 0.5 ml |
| BSB 2864 | Concentrate | 1.0 ml | BSB 2871 | Concentrate | 1.0 ml | BSB 5321 | Concentrate | 1.0 ml |
| BSB 2865 | Control Slides | 5 | BSB 2872 | Control Slides | 5 | BSB 5322 | Control Slides | 5 |

CD117, RMAb




IHC of CD117 on a FFPE GIST Tissue

CD117 is a tyrosine-kinase receptor for stem cell factor (SCF), also known as “steel factor” or “c-kit ligand”. C-kit is a polypeptide that activates bone marrow precursors of a number of blood cells, but its receptor is also present in other cells. C-kit mutations in the interstitial cells of Cajal in the digestive tract are probably the key to Gastrointestinal Stromal Tumors (GISTs), and explain the efficacy of imatinib in the management of these rare malignancies.

CD117 is found on interstitial cells of Cajal, germ cells, bone marrow stem cells, melanocytes, breast epithelium and mast cells. This receptor is found on a wide variety of tumor cells (Follicular and Papillary Carcinoma of the Thyroid, Adenocarcinomas from endometrium, lung, ovary, pancreas, breast; Malignant Melanoma, Endodermal Sinus Tumor, Small-cell Carcinoma) but has been particularly useful in differentiating Gastrointestinal Stromal Tumors (GIST) from Kaposi's Sarcoma and tumors of smooth-muscle origin.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM359
ISOTYPE: IgG
CONTROL: Skin, Testis, Breast, GIST, Colon, Brain, Tonsil
LOCALIZATION: Cytoplasmic, Nuclear, Membranous
SPECIES REACTIVITY: Human, Monkey, Predicted: Marmoset

CD123 IL-3Ra, MAb




IHC of CD123 on a FFPE Kikuchi-Fujimoto Disease Tissue

CD123 is a chain of the IL-3 receptor. This 60-70 kDa transmembrane protein, by itself, binds to IL-3 with rather low affinity. However, when associated with CD131 (common β chain), the protein binds to IL-3 with high affinity. The gene coding for the receptor is located in the pseudoautosomal region of the X and Y chromosomes. The receptor belongs to the Type I cytokine-receptor family and is a heterodimer with a unique alpha chain paired with the common beta (beta c or CDw131) subunit.

The CD123 receptor, found on pluripotent progenitor cells, induces tyrosine phosphorylation within the cell and promotes proliferation and differentiation within the hematopoietic cell lines. CD123 is expressed by myeloid precursors, macrophages, dendritic cells, mast cells, basophils, and megakaryocytes.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-59
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node, Kikuchi-Fujimoto
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

CD137/TNFRSF9, MAb



IHC of CD137/TNFRSF9 on a FFPE Transitional Cell Carcinoma Tissue

CD137, or tumor necrosis factor receptor superfamily member 9 (TNFRSF9), also known as 4-1BB, is encoded by the TNFRSF9 gene on chromosome 1p36.23. CD137/TNFRSF9 is expressed on the surface of certain immune cells, such as Follicular Dendritic Cells, Monocytes, Mast Cells and Granulocytes. CD137/TNFRSF9 is barely detected on resting T cells or T-cell lines, and expression of CD137 is therefore activation dependent. It is found on activated CD4+, CD8+ or Natural Killer Cells. Activation of CD137 by its ligand, CD137L, upregulates survival genes, enhances cell division, induces cytokine production, and prevents activation-induced cell death in T cells.


One study demonstrated CD137/TNFRSF9 positive IHC staining in a majority of Classical Hodgkin Lymphoma (CHL) tissues, indicating its usefulness as a diagnostic biomarker for CHL. CD137/TNFRSF9 was also found to be a highly specific immunohistochemical marker of Neoplastic Follicular Dendritic Cells. Further studies indicate the value of CD137/TNFRSF9 as a prognostic marker in Diffuse Large B-cell Lymphomas, tumor cell differentiation in Gastric Cancer as well as the correlation between CD137/TNFRSF9 expression and bone metastasis of Breast Cancer. CD137/TNFRSF9 is considerably up-regulated in human Gliomas when compared with normal brain tissue.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-159
ISOTYPE: IgM
CONTROL: Colon, Stomach, Tonsil, Testis, Transitional Cell Carcinoma, DBC, Hepatocellular Carcinoma, Diffuse Type Gastric Carcinoma
LOCALIZATION: Nuclear, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3758-3 | Tinto Predilute | 3.0 ml |
| BSB-3758-7 | Tinto Predilute | 7.0 ml |
| BSB-3758-15 | Tinto Predilute | 15.0 ml |
| BSB-3758-01 | Concentrate | 0.1 ml |
| BSB-3758-05 | Concentrate | 0.5 ml |
| BSB-3758-1 | Concentrate | 1.0 ml |
| BSB-3758-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5323 | Tinto Predilute | 3.0 ml |
| BSB 5324 | Tinto Predilute | 7.0 ml |
| BSB 5325 | Tinto Predilute | 15.0 ml |
| BSB 5326 | Concentrate | 0.1 ml |
| BSB 5327 | Concentrate | 0.5 ml |
| BSB 5328 | Concentrate | 1.0 ml |
| BSB 5329 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3717-3 | Tinto Predilute | 3.0 ml |
| BSB-3717-7 | Tinto Predilute | 7.0 ml |
| BSB-3717-15 | Tinto Predilute | 15.0 ml |
| BSB-3717-01 | Concentrate | 0.1 ml |
| BSB-3717-05 | Concentrate | 0.5 ml |
| BSB-3717-1 | Concentrate | 1.0 ml |
| BSB-3717-CS | Control Slides | 5 |

CD138, MAb

IHC of CD138 on a FFPE Tonsil Tissue

CD138/Syndecan-1 is a transmembrane heparin-sulphate proteoglycan which is made up of one core protein and five glycosaminoglycans. CD138 is expected to play a role in cell adhesion. It is expressed on the surface of pre B-cells and plasma cells but is absent from mature B-cells.


Anti-CD138/syndecan-1 is a useful marker for labeling normal and neoplastic plasma cells and Plasmacytoid Lymphomas. It is a selective marker for B-cell Lymphoblastic Leukemia and Lymphoplasmacytoid Leukemia. It is lost from the apoptotic myeloma cells, and thus, is a useful marker for viable Myeloma cells. Various forms of Hodgkin's Disease have also shown positive staining with this antibody.

CD138, RMab

IHC of CD138 on a FFPE Tonsil Tissue

CD138/Syndecan-1 is a transmembrane heparin-sulphate proteoglycan consisting of one core protein and five glycosaminoglycans. CD138 is suspected to play a role in cell adhesion. It is expressed on the surface of pre B-cells and plasma cells but is absent from mature B-cells.

Anti-CD138/syndecan-1 is a useful marker for labeling normal and neoplastic plasma cells and Plasmacytoid Lymphomas. It is a selective marker for B-cell Lymphoblastic Leukemia and Lymphoplasmacytoid Leukemia. It is lost from the apoptotic myeloma cells, and thus, is a useful marker for viable Myeloma cells. Various forms of Hodgkin's Disease have also shown positive staining with this antibody.

CD142/TF/Coagulation Factor III, MAb

IHC of CD142/TF/Coagulation Factor III on a FFPE SARS-Cov-2 Infected Lung Tissue

CD142, also known as Tissue Factor (TF) or Coagulation Factor III or Thromboplastin, is encoded by the F3 gene located on chromosome 1p21.3. CD142 /TF is a 46 kDa sized integral membrane glycoprotein. Upon complex formation with coagulation factor VII, extrinsic blood coagulation is activated by a catalytic cascade that involves specific proteolysis.

CD142/TF is not only associated with the progression, but also with the overall survival rate of many cancers, including breast, gastrointestinal, liver, pancreatic, and prostate cancer. A study that investigated CD142 expression in non-small-cell lung cancer (NSCLC) found that immunohistochemical staining indicates the critical role of CD142 in the progression of NSCLC. A study found increased CD142 levels in breast cancer tissue as well as a correlation between poor survival and high levels of CD142 expression in breast cancer patients.

The SARS-CoV-2 virus triggers the synthesis and release of pro-inflammatory cytokines. One of them is a TNF- α which has been implicated in promoting overexpression of tissue factor (TF) in platelets and macrophages. TF may be a critical mediator associated with the development of thrombotic phenomena in COVID-19.


ANTIBODY TYPE: Mouse Monoclonal**CLONE:** B-A38**ISOTYPE:** IgG1**CONTROL:** Tonsil, Liver, Kidney, Breast, Lymph Node, Cervix, Plasmacytoma**LOCALIZATION:** Membranous**SPECIES REACTIVITY:** Human**ANTIBODY TYPE:** Rabbit Monoclonal**CLONE:** EP201**ISOTYPE:** IgG**CONTROL:** Tonsil, Liver, Kidney, Breast, Lymph Node, Cervix, Plasmacytoma, Adrenal, Skin, Colon, Lung**LOCALIZATION:** Membranous**SPECIES REACTIVITY:** Human, Mouse, Rat**ANTIBODY TYPE:** Mouse Monoclonal**CLONE:** BSB-143**ISOTYPE:** IgG2b**CONTROL:** Placenta, Cervix, Colon, Pancreas, Brain, Kidney, Testis, Pancreatic Carcinoma, Colon Adenocarcinoma**LOCALIZATION:** Cytoplasmic, Membranous**SPECIES REACTIVITY:** Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5330 | Tinto Predilute | 3.0 ml |
| BSB 5331 | Tinto Predilute | 7.0 ml |
| BSB 5332 | Tinto Predilute | 15.0 ml |
| BSB 5333 | Concentrate | 0.1 ml |
| BSB 5334 | Concentrate | 0.5 ml |
| BSB 5335 | Concentrate | 1.0 ml |
| BSB 5336 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6527 | Tinto Predilute | 3.0 ml |
| BSB 6528 | Tinto Predilute | 7.0 ml |
| BSB 6529 | Tinto Predilute | 15.0 ml |
| BSB 6530 | Concentrate | 0.1 ml |
| BSB 6531 | Concentrate | 0.5 ml |
| BSB 6532 | Concentrate | 1.0 ml |
| BSB 6533 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3718-3 | Tinto Predilute | 3.0 ml |
| BSB-3718-7 | Tinto Predilute | 7.0 ml |
| BSB-3718-15 | Tinto Predilute | 15.0 ml |
| BSB-3718-01 | Concentrate | 0.1 ml |
| BSB-3718-05 | Concentrate | 0.5 ml |
| BSB-3718-1 | Concentrate | 1.0 ml |
| BSB-3718-CS | Control Slides | 5 |

CD147, MAb



IHC of CD147 on a FFPE SARS-Cov-2 Infected Lung Tissue


CD147 is a transmembrane glycoprotein also known as Basigin or EMMPRIN (Extracellular Matrix Metalloproteinase Inducer), participating in tumor development and viral entry pathways. CD147 is a member of the immunoglobulin superfamily, upregulated in asthmatic, diabetic, and other inflammatory pathways where it induces expression of Matrix Metalloproteinase 1 and 9. CD147 and its effectors are induced by high-glucose concentration in monocytes and promote their migration. CD147 also participates in fibroblast differentiation through TGF- β 1-induced signaling pathways, and is a marker of undifferentiated human embryonic stem cells.

CD147 has been found to be upregulated in cancer stem cells, resulting in increased metabolism through lactic acid export, increased production of hyaluronan, inhibiting apoptosis, and acting as a main upstream stimulator of matrix metalloproteinases. In clinical studies, CD147 has been associated with larger tumors, deeper invasion, and more lymphocytes. CD147 also facilitates the entry of malaria parasite *Plasmodium falciparum* into red blood cells, and the entry of coronaviruses into human cells by binding the spike glycoprotein, functioning as the second receptor for SARS-CoV-2 and potential participant in pulmonary fibrosis.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-137
ISOTYPE: IgG1
CONTROL: Testis, Colon, Kidney, Stomach, Brain
LOCALIZATION: Membranous, Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3704-3 | Tinto Predilute | 3.0 ml |
| BSB-3704-7 | Tinto Predilute | 7.0 ml |
| BSB-3704-15 | Tinto Predilute | 15.0 ml |
| BSB-3704-01 | Concentrate | 0.1 ml |
| BSB-3704-05 | Concentrate | 0.5 ml |
| BSB-3704-1 | Concentrate | 1.0 ml |
| BSB-3704-CS | Control Slides | 5 |

CD163, MAb



IHC of CD163 on a FFPE Placenta Tissue


CD163 is a 130 kDa membrane glycoprotein. CD163 was recently identified as an acute phase-regulated transmembrane protein whose function is to mediate the endocytosis of haptoglobin-hemoglobin complexes. Solubilized in plasma, CD163 functions as an anti-inflammatory signal and has many roles in disease processes that range from autoimmune conditions such as Rheumatoid Arthritis to Atherosclerosis. CD163 is expressed exclusively on the cell surface of human monocytes and macrophages that evolve predominantly in the late phase of inflammation, and is, therefore, very useful for macrophage-phenotyping. This receptor is expressed on the surface of monocytes (low expression) and histiocytes (high expression).

Staining for CD163 has been helpful in distinguishing synovial macrophages from synovial intimal fibroblasts in the setting of Rheumatoid Arthritis, where its specificity for macrophages was found to be superior to that of CD68, which does not discriminate between these cell types. Flow-cytometry studies have confirmed that CD163 expression is limited to Leukemias with monocytic differentiation. Positive staining can be seen in the skin (histiocytes), gut, Kupffer cells, a few aveolar macrophages, the main population of macrophages in the placenta, and in varying degrees in macrophages in inflamed tissue including tumor tissue, depending on the inflammatory stage. Red-pulp, not white-pulp, macrophages in the spleen and cortical macrophages of the thymus are stained by CD163.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 10D6
ISOTYPE: IgG1
CONTROL: Placenta, Tonsil, Lymph Node, Inflamed Tissue, H. Pylori Infected Tissue
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6303 | Tinto Predilute | 3.0 ml |
| BSB 6304 | Tinto Predilute | 7.0 ml |
| BSB 6305 | Tinto Predilute | 15.0 ml |
| BSB 6306 | Concentrate | 0.1 ml |
| BSB 6307 | Concentrate | 0.5 ml |
| BSB 6308 | Concentrate | 1.0 ml |
| BSB 6309 | Control Slides | 5 |

CDK2, RMAb




IHC of CDK2 on a FFPE Lymphoblastic Lymphoma Tissue

CDK2 (also known as Cyclin-dependent kinase 2 or cell division protein kinase 2) is a 33 kDa enzyme that in humans is encoded by the CDK2 gene on Chromosome 12. CDK2 is a catalytic subunit of the cyclin-dependent kinase complex, whose activity is restricted to the G1-S phase of the cell cycle, and is essential for the cellular G1/S transition.

Studies indicate that over expression of CDK2 may cause the abnormal regulation of cell-cycle, which would could directly contribute to hyperproliferation of cancerous cells. HOX genes, which play a key role in cell differentiation and morphogenesis, are regulated by CDK2. Studies have shown that HOXA7 promoted cell proliferation (mediated by cyclin E1/CDK2) is evident in cases of Hepatocellular carcinoma, indicating that CDK2 may prove to be a useful marker in such cases.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-CDK2
ISOTYPE: IgG
CONTROL: Testis, Tonsil, Prostate, Placenta, Skin, Colon, Transitional Cell Carcinoma, Lymphoblastic Lymphoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2636 | Tinto Predilute | 3.0 ml |
| BSB 2637 | Tinto Predilute | 7.0 ml |
| BSB 2638 | Tinto Predilute | 15.0 ml |
| BSB 2639 | Concentrate | 0.1 ml |
| BSB 2640 | Concentrate | 0.5 ml |
| BSB 2641 | Concentrate | 1.0 ml |
| BSB 2642 | Control Slides | 5 |


CDK4, RMAb

IHC of CDK4 on a FFPE Anal Carcinoma Tissue

Cyclin-dependent kinase 4 (CDK4) is a member of the Ser/Thr protein kinase family. It is a catalytic subunit of the protein kinase complex that is important for cell cycle G1 phase progression. The activity of this kinase is restricted to the G1-S phase, which is controlled by the regulatory subunits D-type cyclins and CDK inhibitor p16 (INK4a).

Overexpression of CDK4 has been observed in many tumor types, including oral squamous cell carcinoma and cancers of the pancreatic (endocrine tumors), lung, breast and colon. The expression of CDK4 is associated with tumor progression.


Binh et al. reported a high expression of CDK4 (92%) in atypical lipomatous tumor/well-differentiated liposarcomas (ALT-WDLPS) and dedifferentiated liposarcomas (DDLPS). CDK4 is useful in differentiating ALT-WDLPS from benign adipose tumors and to separate DDLPS from poorly differentiated sarcomas.

CDX2, RMAb

IHC of CDX2 on a FFPE Colon Adenocarcinoma Tissue

CDX2 is a caudal-type homeobox gene that encodes an intestine-specific transcription factor expressed early in intestinal development and that may be involved in the regulation of proliferation and differentiation of intestinal epithelial cells. It is expressed in the nuclei of epithelial cells throughout the intestine, from duodenum to rectum.

The CDX2 protein is expressed in Primary and Metastatic Colorectal Carcinomas and has also been demonstrated in the intestinal metaplasia of the stomach and intestinal-type gastric cancer. It is not expressed in the normal gastric mucosa. Loss of CDX2 protein expression has been correlated with loss of differentiation in colorectal cancers. Anti-CDX2 antibody has been useful in distinguishing the gastrointestinal origin of Metastatic Adenocarcinomas and carcinoids. Studies have shown that CDX2 is a superior marker compared to CK20. A high percentage of Mucinous Carcinomas of the Ovary also stain positively with this antibody, as well as Carcinomas from the upper gastrointestinal tract.

CEA, RMAb

IHC of CEA on a FFPE Colon Adenocarcinoma Tissue

Carcinoembryonic antigen (CEA) is a glycoprotein involved in cell adhesion. It is normally produced during fetal development, but the production of CEA stops before birth. Therefore, it is not usually present in the blood of healthy adults, although levels are raised in heavy smokers. CEA is synthesized during development in the fetal gut, and is re-expressed in increased amounts in Intestinal Carcinomas and several other tumors.

CEA is employed essentially as a tool to assist in the distinction between Adenocarcinoma and Malignant Mesotheliomas of the epithelial type, along with other markers for mucosubstances such as Leu M1 and Ber-EP4. Another suggested use of CEA is the immunophenotyping of various Metastatic Adenocarcinomas as a means of identifying their origin.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP180
ISOTYPE: IgG
CONTROL: Cervical and Colon Cancer
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP25
ISOTYPE: IgG
CONTROL: Colon, Colon Adenocarcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Rabbit


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-13
ISOTYPE: IgG1/K
CONTROL: Colon, Colon Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2461 | Tinto Predilute | 3.0 ml |
| BSB 2462 | Tinto Predilute | 7.0 ml |
| BSB 2463 | Tinto Predilute | 15.0 ml |
| BSB 2464 | Concentrate | 0.1 ml |
| BSB 2465 | Concentrate | 0.5 ml |
| BSB 2466 | Concentrate | 1.0 ml |
| BSB 2467 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6057 | Tinto Predilute | 3.0 ml |
| BSB 6058 | Tinto Predilute | 7.0 ml |
| BSB 6059 | Tinto Predilute | 15.0 ml |
| BSB 6060 | Concentrate | 0.1 ml |
| BSB 6061 | Concentrate | 0.5 ml |
| BSB 6062 | Concentrate | 1.0 ml |
| BSB 6063 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5337 | Tinto Predilute | 3.0 ml |
| BSB 5338 | Tinto Predilute | 7.0 ml |
| BSB 5339 | Tinto Predilute | 15.0 ml |
| BSB 5340 | Concentrate | 0.1 ml |
| BSB 5341 | Concentrate | 0.5 ml |
| BSB 5342 | Concentrate | 1.0 ml |
| BSB 5343 | Control Slides | 5 |

CEA, RPab




IHC of CEA on a FFPE Colon Adenocarcinoma Tissue

Carcinoembryonic antigen (CEA) is a glycoprotein involved in cell adhesion. It is normally produced during fetal development, but the production of CEA stops before birth. Therefore, it is not usually present in the blood of healthy adults, although levels are raised in heavy smokers. CEA is synthesized during development in the fetal gut, and is re-expressed in increased amounts in Intestinal Carcinomas and several other tumors.

CEA is employed essentially as a tool to assist in the distinction between Adenocarcinoma and Malignant Mesotheliomas of the epithelial type, along with other markers for mucosubstances such as Leu M1 and Ber-EP4. Another suggested use of CEA is the immunophenotyping of various Metastatic Adenocarcinomas as a means of identifying their origin.

Chromogranin A, MAb




IHC of Chromogranin A on a FFPE Pancreas Tissue

Chromogranin A is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Examples of cells producing chromogranin A are the adrenal medulla, enterochromaffin-like cells and beta cells of the pancreas. The function of chromogranin A is unknown but it is a precursor to 3 functional peptides: vasostatin, pancreastatin and parastatin. These peptides negatively modulate the neuroendocrine function of the releasing cell (autocrine) or nearby cells (paracrine).

Chromogranin A is an excellent marker for Carcinoid Tumors, Pheochromocytomas, Paragangliomas, and other Neuroendocrine Tumors. Coexpression of chromogranin A and neuron-specific enolase (NSE) is common in neuroendocrine neoplasms. It has been identified in a wide variety of endocrine tissues including the pituitary, pancreas, hypothalamus, thymus, thyroid, intestine and parathyroid. It is generally accepted that the co-expression of certain keratins and chromogranin means neuroendocrine lineage. The presence of strong chromogranin staining and absence of keratin staining should raise the possibility of paraganglioma. Most pituitary adenomas and prolactinomas readily express chromogranin.

Claudin-1, RPab



IHC of Claudin-1 on a FFPE Skin Tissue

Claudin 1 is an integral membrane protein and a component of tight junction strands. Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. These junctions are composed of sets of continuous networking strands in the outwardly facing cytoplasmic leaflet, with complementary grooves in the inwardly facing extracytoplasmic leaflet.

Claudin-1 stains membranes of cells and is found in nearly all carcinomas, with stains much stronger in carcinoma cells than in normal tissue cells.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Colon, Tonsil, Fallopian Tube, Colon Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Mouse Monoclonal
CLONE: LK2H10
ISOTYPE: IgG1/K
CONTROL: Pancreas, Pituitary, Colon, Brain
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat, Rat, Rabbit, Porcine


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Skin, Small Intestine, Colon Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6534 | Tinto Predilute | 3.0 ml |
| BSB 6535 | Tinto Predilute | 7.0 ml |
| BSB 6536 | Tinto Predilute | 15.0 ml |
| BSB 6537 | Concentrate | 0.1 ml |
| BSB 6538 | Concentrate | 0.5 ml |
| BSB 6539 | Concentrate | 1.0 ml |
| BSB 6540 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5344 | Tinto Predilute | 3.0 ml |
| BSB 5345 | Tinto Predilute | 7.0 ml |
| BSB 5346 | Tinto Predilute | 15.0 ml |
| BSB 5347 | Concentrate | 0.1 ml |
| BSB 5348 | Concentrate | 0.5 ml |
| BSB 5349 | Concentrate | 1.0 ml |
| BSB 5350 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6562 | Tinto Predilute | 3.0 ml |
| BSB 6563 | Tinto Predilute | 7.0 ml |
| BSB 6564 | Tinto Predilute | 15.0 ml |
| BSB 6565 | Concentrate | 0.1 ml |
| BSB 6566 | Concentrate | 0.5 ml |
| BSB 6567 | Concentrate | 1.0 ml |
| BSB 6568 | Control Slides | 5 |

Claudin-5, RMab




IHC of Claudin-5 on a FFPE Liver Tissue

Claudin-5 is a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction (TJ) strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets. Claudin-5 is an endothelial cell-specific component of TJ strands. Mutations in Claudin-5 have been found in patients with velocardiofacial syndrome.

Claudin-5 labels endothelial cells and has been used as a marker for endothelial lesions. Claudin-5 is also found in bronchial and lung epithelial cells. In tumors, Claudin-5 expression has been found in lung adenocarcinoma and squamous carcinoma. In serous ovarian adenocarcinoma, increased Claudin-5 expression is associated with aggressive behavior.

Claudin-7, RMab




IHC of Claudin-7 on a FFPE Transitional Cell Carcinoma Tissue

Claudin-7 is encoded by the CLDN7 gene and is a member of the claudin family. Claudin-7 is a component of tight junctions and plays important roles in maintaining cell polarity and tightly connecting the barriers between cells. Abnormal Claudin-7 protein expression is closely related to tumor occurrence, development, and metastasis.

Claudin-7 expression is down-regulated in Colorectal Cancer (CRC) and Lung Cancer tissues and loss of Claudin-7 is associated with the degree of cancer cell differentiation and metastasis. Claudin-7 may serve as a tumor suppressor, as it is involved in pathways of cell proliferation and epithelial barrier formation through Integrin signaling networks. The deficient extracellular matrix in cells with low Claudin-7 expression promotes abnormal growth patterns and metastasis in Ductal and Lobular Breast Carcinomas, and Claudin-7 expression is lost in ~50% of primary tumors in Breast Cancers. While Claudin-7 downregulation was found in Esophageal, Head/Neck, and Prostate Cancers, overexpression of Claudin-7 is found in many Ovarian Cancers and may also lead to increased tumor invasiveness.

Collagen Type IV, MAb



IHC of Collagen IV on a FFPE Skin Tissue

Collagen is the main protein of connective tissue in animals and the most abundant protein in mammals, making up about 25% of the total protein content. Collagen IV is a major constituent of the basement membranes, along with laminins and enactins. It is composed of the alpha 1 IV chain and alpha 2 IV chain in a 2:1 ratio. It can form insoluble fibers with high tensile strength.

Normal tissue stains with this antibody in a manner consistent with the sites of mesenchymal elements and epithelial basal laminae. Antibody to collagen IV is useful in detecting the loss of parts of basement membrane in carcinomas. Collagen IV can also be useful in the classification of soft tissue tumors; Schwannomas, Leiomyomas, and their well-differentiated malignant counterparts usually immunoreact to this antibody. The vascular nature of neoplasms, Hemangiopericytoma, Angiosarcoma and Epithelioid Hemangioendothelioma can be observed with this antibody.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP224

ISOTYPE: IgG

CONTROL: Liver, Vascular Tissue, Placenta, Colon, Kidney, Fallopian Tube

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP399

ISOTYPE: IgG

CONTROL: Breast, Colon, Fallopian Tube, Pancreas, Kidney, Transitional Cell Carcinoma

LOCALIZATION: Membranous, Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: CIV22

ISOTYPE: IgG1/K


CONTROL: Muscle, Lung

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Rat, Dog, Horse


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|-------------|-----------------|---------|----------|-----------------|---------|
| BSB 2398 | Tinto Predilute | 3.0 ml | BSB-3719-3 | Tinto Predilute | 3.0 ml | BSB 5351 | Tinto Predilute | 3.0 ml |
| BSB 2399 | Tinto Predilute | 7.0 ml | BSB-3719-7 | Tinto Predilute | 7.0 ml | BSB 5352 | Tinto Predilute | 7.0 ml |
| BSB 2400 | Tinto Predilute | 15.0 ml | BSB-3719-15 | Tinto Predilute | 15.0 ml | BSB 5353 | Tinto Predilute | 15.0 ml |
| BSB 2401 | Concentrate | 0.1 ml | BSB-3719-01 | Concentrate | 0.1 ml | BSB 5354 | Concentrate | 0.1 ml |
| BSB 2402 | Concentrate | 0.5 ml | BSB-3719-05 | Concentrate | 0.5 ml | BSB 5355 | Concentrate | 0.5 ml |
| BSB 2403 | Concentrate | 1.0 ml | BSB-3719-1 | Concentrate | 1.0 ml | BSB 5356 | Concentrate | 1.0 ml |
| BSB 2404 | Control Slides | 5 | BSB-3719-CS | Control Slides | 5 | BSB 5357 | Control Slides | 5 |

Collagen Type IV, RMAb




IHC of Collagen Type IV on a FFPE Kidney Tissue

COX-2, RMAb



IHC of COX-2 on a FFPE Liver Tissue

COX-2, RMAb



IHC of COX2 on a FFPE Liver Tissue

Collagen is the main protein of connective tissue in animals and the most abundant protein in mammals, making up about 25% of the total protein content. Collagen IV is a major constituent of the basement membranes, along with laminins and enactins. It is composed of the alpha 1 IV chain and alpha 2 IV chain in a 2:1 ratio. It can form insoluble fibers with high tensile strength.

Normal tissue stains with this antibody in a manner consistent with the sites of mesenchymal elements and epithelial basal laminae. Antibody to collagen IV is useful in detecting the loss of parts of basement membrane in carcinomas. Collagen IV can also be useful in the classification of soft tissue tumors; Schwannomas, Leiomyomas, and their well-differentiated malignant counterparts usually immune react to this antibody. The vascular nature of neoplasms, Hemangiopericytoma, Angiosarcoma and Epithelioid Hemangioendothelioma can be observed with this antibody.

Prostaglandin-endoperoxide synthase 2 (prostaglandin G/H synthase and cyclooxygenase), also known as cyclooxygenase-2 or COX-2, is an enzyme that in humans is encoded by the PTGS2 gene and it is involved in the conversion of arachidonic acid to prostaglandin H2, an important precursor of prostacyclin and thromboxane A2, among others. Pharmacological inhibition of COX can provide relief from the symptoms of inflammation and pain; this is the method of action of well-known drugs such as aspirin and ibuprofen. COX-2 inhibition by nonsteroidal anti-inflammatory agents has been shown to decrease angiogenesis and tumor growth, and promote apoptosis.

The expression of COX-2 is upregulated in many cancers. COX-2 overexpression has been associated with increased microvascular density, and VEGF protein expression in head and neck Squamous Cell Carcinomas and is a poor prognostic indicator in this entity as well. COX-2 overexpression has also been suggested as a poor prognostic indicator in Carcinomas of the Colon, Breast, Pancreas, and Adenocarcinomas of the Lung.

Prostaglandin-endoperoxide synthase 2 (prostaglandin G/H synthase and cyclooxygenase), also known as cyclooxygenase-2 or COX-2, is an enzyme that in humans is encoded by the PTGS2 gene and it is involved in the conversion of arachidonic acid to prostaglandin H2, an important precursor of prostacyclin and thromboxane A2, among others. Pharmacological inhibition of COX can provide relief from the symptoms of inflammation and pain; this is the method of action of well-known drugs such as aspirin and ibuprofen. COX-2 inhibition by nonsteroidal anti-inflammatory agents has been shown to decrease angiogenesis and tumor growth, and promote apoptosis.

The expression of COX-2 is upregulated in many cancers. COX-2 overexpression has been associated with increased microvascular density, and VEGF protein expression in head and neck Squamous Cell Carcinomas and is a poor prognostic indicator in this entity as well. COX-2 overexpression has also been suggested as a poor prognostic indicator in Carcinomas of the Colon, Breast, Pancreas, and Adenocarcinomas of the Lung.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-COL4

ISOTYPE: IgG

CONTROL: Muscle, Lung, Breast, Placenta, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma

LOCALIZATION: Membranous, Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-COX2

ISOTYPE: IgG

CONTROL: Colon, Testis, Kidney, Placenta, Liver, Fallopian Tube, Pancreas, Tonsil, Thymus, Breast, Adenocarcinoma of Colon, Bladder TCC

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP293

ISOTYPE: IgG

CONTROL: Colon, Stomach, Pancreas, Breast, Lung, Adenocarcinoma of Colon

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3777-3 | Tinto Predilute | 3.0 ml |
| BSB-3777-7 | Tinto Predilute | 7.0 ml |
| BSB-3777-15 | Tinto Predilute | 15.0 ml |
| BSB-3777-01 | Concentrate | 0.1 ml |
| BSB-3777-05 | Concentrate | 0.5 ml |
| BSB-3777-1 | Concentrate | 1.0 ml |
| BSB-3777-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5358 | Tinto Predilute | 3.0 ml |
| BSB 5359 | Tinto Predilute | 7.0 ml |
| BSB 5360 | Tinto Predilute | 15.0 ml |
| BSB 5361 | Concentrate | 0.1 ml |
| BSB 5362 | Concentrate | 0.5 ml |
| BSB 5363 | Concentrate | 1.0 ml |
| BSB 5364 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2873 | Tinto Predilute | 3.0 ml |
| BSB 2874 | Tinto Predilute | 7.0 ml |
| BSB 2875 | Tinto Predilute | 15.0 ml |
| BSB 2876 | Concentrate | 0.1 ml |
| BSB 2877 | Concentrate | 0.5 ml |
| BSB 2878 | Concentrate | 1.0 ml |
| BSB 2879 | Control Slides | 5 |

CTLA-4/CD152, MAb



IHC of CTLA-4 on a FFPE Tonsil

CTLA4 or CTLA-4 (cytotoxic T-lymphocyte-associated protein 4), classified as CD152, is a protein receptor known to function as an immune checkpoint which downregulates the immune system. CTLA4 is found on the surface of T cells, and acts as an "off" switch when bound to CD80 or CD86 on the surface of antigen-presenting cells. CTLA4 is a member of the immunoglobulin superfamily that is expressed on the surface of Helper T cells and transmits an inhibitory signal to T cells. CTLA4 transmits an inhibitory signal to T cells, whereas CD28 transmits a stimulatory signal. Intracellular CTLA4 is also found in regulatory T cells and may be important to their function. T cell activation through the T cell receptor and CD28 leads to increased expression of CTLA-4, an inhibitory receptor for B7 molecules.

Mutations in this gene have been associated with insulin-dependent diabetes mellitus, Graves' disease, Hashimoto's thyroiditis, celiac disease, systemic lupus erythematosus, thyroid-associated orbitopathy, primary biliary cirrhosis and other autoimmune diseases. Polymorphisms of the CTLA-4 gene are associated with autoimmune diseases such as autoimmune thyroid disease and multiple sclerosis, though this association is often weak.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-88


ISOTYPE: IgG2a/K

CONTROL: Tonsil, Lymph Node, Colon, Thymus

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human, Mouse, Rat, Dog, Horse

CTLA4, RMAb



IHC of CTLA-4 on a FFPE Colon Tissue

CTLA4 or CTLA-4 (cytotoxic T-lymphocyte-associated protein 4), classified as CD152, is a protein receptor known to function as an immune checkpoint which downregulates the immune system. CTLA4 is found on the surface of T cells, and acts as an "off" switch when bound to CD80 or CD86 on the surface of antigen-presenting cells. The CTLA-4 protein is encoded by the *Ctla4* gene in mice and the *CTLA4* gene in humans. CTLA4 is a member of the immunoglobulin superfamily that is expressed on the surface of Helper T cells and transmits an inhibitory signal to T cells.

Mutations in this gene have been associated with insulin-dependent diabetes mellitus, Graves' disease, Hashimoto's thyroiditis, celiac disease, systemic lupus erythematosus, thyroid-associated orbitopathy, primary biliary cirrhosis and other autoimmune diseases. Polymorphisms of the CTLA-4 gene are associated with autoimmune diseases such as autoimmune thyroid disease and multiple sclerosis, though this association is often weak. In Systemic Lupus Erythematosus (SLE), the splice variant sCTLA-4 is found to be aberrantly produced and found in the serum of patients with active SLE. Germline haploinsufficiency of CTLA4 leads to CTLA4 deficiency or CHAI disease (CTLA4 haploinsufficiency with autoimmune infiltration), a rare genetic disorder of the immune system.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-CTLA4


ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Colon, Thymus

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

CXCL12/SDF-1, MAb



IHC of CXCL12/SDF-1 on a FFPE HER2 Negative Breast Carcinoma Tissue

The cytokine C-X-C motif chemokine 12 (CXCL12) is synthesized by metastasis target tissues and has been shown to attract tumor cells that express the receptor, C-X-C chemokine receptor type 4 (CXCR4). CXCL12 is involved in tumor protection and metastasis, where it induces angiogenesis and epithelial-to-mesenchymal transition and inhibits apoptosis and host immune response. CXCL12 is highly expressed in cancer-associated fibroblasts and tumor stroma, where it sequesters T cells and prevents them from infiltrating and attacking the tumor.

Expression of CXCL12 and CXCR4 in breast, pancreatic, esophageal, lung, prostate, and ovarian Cancers increases angiogenesis and is associated with poor prognosis, especially in the presence of HER2-neu which inhibits degradation of the CXCL12/CXCR4 signaling complex. Expression of CXCL12 in organs like lungs, lymph nodes, bone marrow, liver increases likelihood of metastasis to those sites. CXCL12 expression and activity is controlled by hypoxia, ACKR3 (CXCR7) receptor, and post-translational modifications; CXCL12 hypermethylation has been reported in gastric, breast, colon, lung, and prostate cancer. A large study found high CXCL12 expression associated with reduced overall survival in patients with oesophagogastric, pancreatic and lung cancer, whereas in breast cancer patients high CXCL12 expression conferred an overall survival advantage.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-165

ISOTYPE: IgG1

CONTROL: Testis, Breast, Colon, Fallopian Tube, Tonsil, Transitional Cell Carcinoma, T cell Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2880 | Tinto Predilute | 3.0 ml |
| BSB 2881 | Tinto Predilute | 7.0 ml |
| BSB 2882 | Tinto Predilute | 15.0 ml |
| BSB 2883 | Concentrate | 0.1 ml |
| BSB 2884 | Concentrate | 0.5 ml |
| BSB 2885 | Concentrate | 1.0 ml |
| BSB 2886 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3483 | Tinto Predilute | 3.0 ml |
| BSB 3484 | Tinto Predilute | 7.0 ml |
| BSB 3485 | Tinto Predilute | 15.0 ml |
| BSB 3486 | Concentrate | 0.1 ml |
| BSB 3487 | Concentrate | 0.5 ml |
| BSB 3488 | Concentrate | 1.0 ml |
| BSB 3489 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3713-3 | Tinto Predilute | 3.0 ml |
| BSB-3713-7 | Tinto Predilute | 7.0 ml |
| BSB-3713-15 | Tinto Predilute | 15.0 ml |
| BSB-3713-01 | Concentrate | 0.1 ml |
| BSB-3713-05 | Concentrate | 0.5 ml |
| BSB-3713-1 | Concentrate | 1.0 ml |
| BSB-3713-CS | Control Slides | 5 |

CXCR4/CD184/Fusin, RMAb




IHC of CXCR4/CD184/Fusin on a FFPE Lung Neuroendocrine Tissue

CXCR5/CD185, RPAb



IHC of CXCR5/CD185 on a FFPE Lung Adenocarcinoma Tissue

Cyclin B1, RMAb



IHC of Cyclin B1 on a FFPE Lymphoblastic Lymphoma Tissue

Chemokine Receptor 4 (CXCR4, also CD184 or Fusin) is a homeostatic chemokine. Its major function is to regulate hematopoietic cell trafficking and secondary lymphoid tissue architecture. CXCR4 is widely expressed on Hematopoietic cells, including CD34+ HSC, T Lymphocytes, B Lymphocytes, Monocytes, Macrophages, Neutrophils, and Eosinophils, as well as organs including Brain, Heart, Colon, and Kidney, and Endothelial cells, Epithelial cells, and Smooth Muscle progenitors.

With roles in cell homing and tissue repair, CXCR4 is involved in cancer cell physiology, survival, and metastasis in >20 cancer types. CXCR4 induces and recruits stem cells and promotes angiogenesis in tumors, and responds to CXCL12 expression to promote metastasis in cancers of the Bone, Brain, Breast, Lung, Liver, Kidney, and other tissues. The SDF-1/CXCR4 axis plays a crucial role in engraftment, where CXCR4 expression is correlated with inflammation and graft rejection. The most studied role for CXCL12 chemokine and its receptor CXCR4 in Breast Cancer pathogenesis is the metastasis event, although several reports have demonstrated its involvement in other processes, such as angiogenesis and tumor growth. It has been found that CXCR4 is required for Breast Cancer cell migration to other sites such as Lung, Bone, and Lymph nodes, which express high levels of CXCL12 chemokine. Therefore, CXCR4 is being considered a prognostic marker in Breast Cancer.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP394

ISOTYPE: IgG

CONTROL: Fallopian Tube, Adrenal Gland, Stomach, Kidney, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma, Tonsil, Testis

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Mouse

The G-protein coupled receptor (GPCR) or C-X-C motif chemokine receptor 5, CXCR5, also known as or Burkitt lymphoma receptor 1 (BLR1), plays fundamental roles in inflammatory, infectious and immune responses. Current evidence also indicates that the CXCL13:CXCR5 axis orchestrates cell-cell interactions that regulate lymphocyte infiltration within the tumor microenvironment, thereby determining responsiveness to cytotoxic and immune-targeted therapies. CXCR5 is expressed in mature B-cells and Burkitt's lymphoma.

CXCR5 is highly expressed in primary and secondary follicles within gastric lymphomas. non-small cell lung cancer (NSCLC) tissues express CXCR5, which correlates with stage/grade of the disease. Higher CXCR5 expression and migration by NSCLC cells suggest a role in migration and metastasis of primary lung tumors in response to CXCL13. These findings indicate that differential expression patterns of CXCR5 and CXCL13 in two subtypes (squamous cell carcinoma and adenocarcinoma) of NSCLC are associated with differences in their prognosis and survival. It has been proposed that CXCR5/CXCL13, either alone or in combination, could be used as a prognostic biomarker for lung cancer. Other studies have shown that CXCR5 overexpression in breast cancer patients highly correlates with lymph node metastases, and elevated CXCR5 expression may contribute to abnormal cell survival and migration in breast tumors that lack functional p53 protein.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG2b

CONTROL: Placenta, Brain, Lung, Transitional Cell Carcinoma, Testis, Ovarian Serous Carcinoma, Hepatocellular Carcinoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

Cyclin B1 is a regulatory protein involved with mitosis. It complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Cyclin B1 contributes to the switch-like all or none behavior of the cell in deciding to commit mitosis. Its activation is well regulated, and positive feedback loops ensure that once the cyclin B1-Cdk1 complex is activated it is not deactivated. Cyclin B1-Cdk1 is involved in the early events of mitosis, such as chromosome condensation, nuclear envelope breakdown, and spindle pole assembly. Before mitosis almost all cyclin B1 in the cell is located in the cytoplasm, but in late prophase it relocates to the nucleus. At the end of mitosis, cyclin B1 is targeted for degradation by the APC through its APC localization sequence, permitting the cell to exit mitosis.

Cyclin B1 has been shown to be overexpressed in various tumor types.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-B1

ISOTYPE: IgG

CONTROL: Cervix, HSIL Cervix

LOCALIZATION: Cytoplasmic, Nuclear


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3758-3 | Tinto Predilute | 3.0 ml |
| BSB-3758-7 | Tinto Predilute | 7.0 ml |
| BSB-3758-15 | Tinto Predilute | 15.0 ml |
| BSB-3758-01 | Concentrate | 0.1 ml |
| BSB-3758-05 | Concentrate | 0.5 ml |
| BSB-3758-1 | Concentrate | 1.0 ml |
| BSB-3758-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3721-3 | Tinto Predilute | 3.0 ml |
| BSB-3721-7 | Tinto Predilute | 7.0 ml |
| BSB-3721-15 | Tinto Predilute | 15.0 ml |
| BSB-3721-01 | Concentrate | 0.1 ml |
| BSB-3721-05 | Concentrate | 0.5 ml |
| BSB-3721-1 | Concentrate | 1.0 ml |
| BSB-3721-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6548 | Tinto Predilute | 3.0 ml |
| BSB 6549 | Tinto Predilute | 7.0 ml |
| BSB 6550 | Tinto Predilute | 15.0 ml |
| BSB 6551 | Concentrate | 0.1 ml |
| BSB 6552 | Concentrate | 0.5 ml |
| BSB 6553 | Concentrate | 1.0 ml |
| BSB 6554 | Control Slides | 5 |

Cyclin B1, RMAb



IHC of Cyclin B1 on FPPE Anal Carcinoma Tissue


Cyclin B1 is a regulatory protein involved with mitosis. It complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Cyclin B1 contributes to the switch-like all or none behavior of the cell in deciding to commit mitosis. Its activation is well regulated, and positive feedback loops ensure that once the cyclin

B1-Cdk1 complex is activated, it is not deactivated. Cyclin B1-Cdk1 is involved in the early events of mitosis, such as chromosome condensation, nuclear envelope breakdown, and spindle pole assembly. Before mitosis almost all cyclin B1 in the cell is located in the cytoplasm, but in late prophase it relocates to the nucleus. At the end of mitosis, cyclin B1 is targeted for degradation by the APC through its APC localization sequence, permitting the cell to exit mitosis. Cyclin B1 has been shown to be overexpressed in various tumor types.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM281
ISOTYPE: IgG
CONTROL: Cervix, HSIL Cervix
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3760-3 | Tinto Predilute | 3.0 ml |
| BSB-3760-7 | Tinto Predilute | 7.0 ml |
| BSB-3760-15 | Tinto Predilute | 15.0 ml |
| BSB-3760-01 | Concentrate | 0.1 ml |
| BSB-3760-05 | Concentrate | 0.5 ml |
| BSB-3760-1 | Concentrate | 1.0 ml |
| BSB-3760-CS | Control Slides | 5 |

Cyclin D1, RMAb



IHC of Cyclin D1 on a FPPE Mantle Cell Lymphoma Tissue


Cyclins are a family of proteins involved in the progression of cells through the cell cycle. Cyclins form a complex with their partner, cyclin-dependent kinase (Cdk), which activates the latter's protein kinase function. Cyclins are so named because they are produced or degraded as needed in order to drive the cell through the different stages of the cell cycle. When its concentrations in the cell are low, the cyclin detaches from the Cdk, inhibiting the enzyme's activity, probably by causing a protein chain to block the enzymatic site.

Cyclin D1 or PRAD-1 or bcl-1 is one of the key cell-cycle regulators, and functions in association with Cdk4 and/or Cdk6 by phosphorylating the Rb protein. It is a putative proto-oncogene overexpressed in a wide variety of human neoplasms including Mantle Cell Lymphomas (MCL).

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT14
ISOTYPE: IgG
CONTROL: Tonsil, Placenta, Brain, Pituitary, Adrenal, Cervix, Breast, Mantle Cell Lymphoma, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5365 | Tinto Predilute | 3.0 ml |
| BSB 5366 | Tinto Predilute | 7.0 ml |
| BSB 5367 | Tinto Predilute | 15.0 ml |
| BSB 5368 | Concentrate | 0.1 ml |
| BSB 5369 | Concentrate | 0.5 ml |
| BSB 5370 | Concentrate | 1.0 ml |
| BSB 5371 | Control Slides | 5 |

Cyclin D1, RMAb



IHC of Cyclin D1 on FPPE Mantle Cell Carcinoma Tissue


Cyclins are a family of proteins involved in the progression of cells through the cell cycle. Cyclins form a complex with their partner, cyclin-dependent kinase (Cdk), which activates the latter's protein kinase function. Cyclins are so named because they are produced or degraded as needed in order to drive the cell through the different stages of the cell cycle. When its concentrations in the cell are low, the cyclin detaches from the Cdk, inhibiting the enzyme's activity, probably by causing a protein chain to block the enzymatic site.

Cyclin D1 or PRAD-1 or bcl-1 is one of the key cell-cycle regulators, and functions in association with Cdk4 and/or Cdk6 by phosphorylating the Rb protein. It is a putative proto-oncogene overexpressed in a wide variety of human neoplasms including Mantle Cell Lymphomas. Cyclin D1 has been found to be overexpressed in breast carcinoma.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM241
ISOTYPE: IgG
CONTROL: Tonsil, Placenta, Brain, Cervix, Breast, Mantle Cell Lymphoma, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Predicted: Mouse, Rat


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3761-3 | Tinto Predilute | 3.0 ml |
| BSB-3761-7 | Tinto Predilute | 7.0 ml |
| BSB-3761-15 | Tinto Predilute | 15.0 ml |
| BSB-3761-01 | Concentrate | 0.1 ml |
| BSB-3761-05 | Concentrate | 0.5 ml |
| BSB-3761-1 | Concentrate | 1.0 ml |
| BSB-3761-CS | Control Slides | 5 |

Cyclin E1, RMAb




IHC of Cyclin E1 on a FFPE Breast Carcinoma Tissue

Cytokeratin 4, RMAb



IHC of Cytokeratin 4 on a FFPE Oral Mucosa Tissue

Cytokeratin 5, RMAb



IHC of Cytokeratin 5 on a FFPE Mesothelioma Tissue

Cyclin E1 forms a complex with and functions as a regulatory subunit of CDK2, whose activity is required for cell cycle G1/S transition phase of the cell cycle that determines cell division. The Cyclin E/CDK2 complex phosphorylates p27Kip1 (an inhibitor of Cyclin D), tagging it for degradation and thus promoting expression of Cyclin A, allowing progression to the S phase. This protein accumulates at the G1-S phase boundary and is degraded as cells progress through S phase. Apart from the function in cell cycle progression, cyclin E/CDK2 plays a role in the centrosome cycle by phosphorylating nucleophosmin (NPM). NPM is then released from binding to an unduplicated centrosome, thereby triggering duplication. Cyclin E/CDK2 has also been shown to regulate the apoptotic response to DNA damage via phosphorylation of FOXO1.

Overexpression of Cyclin E correlates with tumorigenesis. It is involved in various types of cancers, including breast, colon, bladder, skin, and lung cancer.

Cytokeratin 4 is a type II cytoke­ratin and is specifically found in differentiated layers of the mucosal and esophageal epithelia together with Cytoke­ratin 13. Mutations in the genes encoding this protein (KRT4) have been associated with White Sponge Nevus, characterized by oral, esophageal, and anal leukoplakia.

A decreased expression of CK4 is associated with head and neck squamous carcinoma. It is helpful in the differentiation of squamous cell carcinoma of esophagus origin from that of thyroid origin.

Cytokeratin 5 is a type II cytoke­ratin found in squamous cell epithelium, myoepithelial cells of the breast and the basal cells of the prostate. Both Cytokeratin 5 and its corresponding partner, Cytokeratin 14, are essential for formation of 8-nm filaments.

Cytokeratin 5 is expressed in most Epithelial Mesotheliomas but not by most Pulmonary Adenocarcinomas and can be used to differentiate between the two.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP126
ISOTYPE: IgG
CONTROL: Placenta, Bladder, Colon, Breast Cancer, Ovarian Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP4
ISOTYPE: IgG
CONTROL: Cornea, Anus, Larynx, Pharynx, Tongue, Prostate, Tonsil, Cervix, Squamous Epithelium of Esophagus, Cervical Squamous Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP24
ISOTYPE: IgG
CONTROL: Prostate, Breast, Placenta, Skin, Mesothelioma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6555 | Tinto Predilute | 3.0 ml |
| BSB 6556 | Tinto Predilute | 7.0 ml |
| BSB 6557 | Tinto Predilute | 15.0 ml |
| BSB 6558 | Concentrate | 0.1 ml |
| BSB 6559 | Concentrate | 0.5 ml |
| BSB 6560 | Concentrate | 1.0 ml |
| BSB 6561 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6590 | Tinto Predilute | 3.0 ml |
| BSB 6591 | Tinto Predilute | 7.0 ml |
| BSB 6592 | Tinto Predilute | 15.0 ml |
| BSB 6593 | Concentrate | 0.1 ml |
| BSB 6594 | Concentrate | 0.5 ml |
| BSB 6595 | Concentrate | 1.0 ml |
| BSB 6596 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6597 | Tinto Predilute | 3.0 ml |
| BSB 6598 | Tinto Predilute | 7.0 ml |
| BSB 6599 | Tinto Predilute | 15.0 ml |
| BSB 6600 | Concentrate | 0.1 ml |
| BSB 6601 | Concentrate | 0.5 ml |
| BSB 6602 | Concentrate | 1.0 ml |
| BSB 6603 | Control Slides | 5 |

Cytokeratin 5 & 6, MMab




IHC of Cytokeratin 5 and 6 on a FFPE Prostatic Adenocarcinoma Tissue

Cytokeratin 5 (58 kDa) is a high-molecular weight, basic type of cytokeratin expressed in basal, intermediate and superficial-cell layers of stratified epithelia as well as transitional epithelia, complex epithelia, mesothelial cells and Mesothelioma. Cytokeratin 6 (56 kD) is also a high-molecular weight, basic type cytokeratin expressed by proliferating squamous epithelium often paired with Cytokeratin 16.

CK 5 and 6 are positively seen in nearly 100% of Malignant Mesotheliomas and is rarely seen in Lung Adenocarcinomas. CK 5 and 6 can positively be seen in undifferentiated Large-cell Carcinoma as well as Squamous Carcinoma. Fewer than 10% of Carcinomas of the breast, colon, and prostate stain positively for this marker. CK 5 and 6 have also been used successfully as a myoepithelial cell marker in the prostate to determine malignancy.

Cytokeratin 5 & 6, RMAb




IHC and IF of Cytokeratin 4 & 6 on an FFPE Meso-othelioma Tissue (IHC) an FFPE Colon Carcinoma Tissue (IF)

Cytokeratin 5 (58 kDa) is a high-molecular weight, basic type of cytokeratin expressed in basal, intermediate and superficial-cell layers of stratified epithelia as well as transitional epithelia, complex epithelia, mesothelial cells and Mesothelioma. Cytokeratin 6 (56 kD) is also a high-molecular weight, basic type cytokeratin expressed by proliferating squamous epithelium often paired with Cytokeratin 16.

CK 5 and 6 are positively seen in nearly 100% of Malignant Mesotheliomas and are rarely seen in Lung Adenocarcinomas. CK 5 and 6 can positively be seen in undifferentiated Large-cell Carcinoma as well as Squamous Carcinoma. Fewer than 10% of Carcinomas of the breast, colon, and prostate stain positively for this marker. CK 5 and 6 have also been used successfully as a myoepithelial cell marker in the prostate to determine malignancy.

Cytokeratin 6, RMAb



IHC of Cytokeratin 6 on a FFPE Cervical Cancer Tissue

Cytokeratin 6 is a type II cytokeratin known for its strong induction in stratified epithelia that features an enhanced cell proliferation rate or abnormal differentiation during wound healing in several diseases (such as psoriasis, actinic keratosis) and in cancer. It can be found on stratified epithelia including oral mucosa, esophagus, basal layer of epidermis, the outer root sheath of hair follicles, and in glandular epithelia.

Together, CK 5 & 6 can be used to differentiate Mesothelioma (positive) from Lung Carcinoma (negative) or metastatic carcinoma (negative), and can also be used to distinguish between Ductal Hyperplasia of the breast (positive) from Solid Papillary DCIS (negative).

ANTIBODY TYPE: Mouse Monoclonal
CLONE: D5/16 B4
ISOTYPE: IgG1
CONTROL: Prostate, Mesothelioma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP24 & EP67
ISOTYPE: IgG
CONTROL: Prostate, Mesothelioma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP67
ISOTYPE: IgG
CONTROL: Oral Mucosa, Esophagus, Skin, Glandular Epithelia
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5400 | Tinto Predilute | 3.0 ml |
| BSB 5401 | Tinto Predilute | 7.0 ml |
| BSB 5402 | Tinto Predilute | 15.0 ml |
| BSB 5403 | Concentrate | 0.1 ml |
| BSB 5404 | Concentrate | 0.5 ml |
| BSB 5405 | Concentrate | 1.0 ml |
| BSB 5406 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6604 | Tinto Predilute | 3.0 ml |
| BSB 6605 | Tinto Predilute | 7.0 ml |
| BSB 6606 | Tinto Predilute | 15.0 ml |
| BSB 6607 | Concentrate | 0.1 ml |
| BSB 6608 | Concentrate | 0.5 ml |
| BSB 6609 | Concentrate | 1.0 ml |
| BSB 6610 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6639 | Tinto Predilute | 3.0 ml |
| BSB 6640 | Tinto Predilute | 7.0 ml |
| BSB 6641 | Tinto Predilute | 15.0 ml |
| BSB 6642 | Concentrate | 0.1 ml |
| BSB 6643 | Concentrate | 0.5 ml |
| BSB 6644 | Concentrate | 1.0 ml |
| BSB 6645 | Control Slides | 5 |

Cytokeratin 7, MMab




IHC of Cytokeratin 7 on a FFPE Lung Adenocarcinoma Tissue

Cytokeratin 7, RMAb



IHC of Cytokeratin 7 on FFPE Lung Carcinoma Tissue

Cytokeratin 8, RMAb



IHC of Cytokeratin 8 on a FFPE Colon Tissue

Cytokeratin 7 (CK7) reacts with proteins that are found in most ductal, glandular and transitional epithelium of the urinary tract and bile duct epithelial cells. CK 7 distinguishes between lung and breast epithelium that stain positive, and colon and prostate epithelial cells that are negative.

This antibody also reacts with many benign and malignant epithelial lesions (e.g., Adenocarcinomas of the ovary, breast and lung). Further, in frozen sections, the antibody has been shown to label the rete epithelium in the testis, epididymis epithelium, and the surface epithelium of the stomach and duodenum. Transitional-cell Carcinomas are positive and Prostate Cancers are negative. This antibody does not recognize intermediate filament proteins, nor does it recognize non-epithelial tissues such as blood vessels, connective tissue, etc.

Cytokeratin 7 (CK7) reacts with proteins that are found in most ductal, glandular and transitional epithelium of the urinary tract and bile duct epithelial cells. CK7 distinguishes between lung and breast epithelium that stain positive, and colon and prostate epithelial cells that are negative.

This antibody also reacts with many benign and malignant epithelial lesions (e.g., Adenocarcinomas of the ovary, breast and lung). Further, in frozen sections, the antibody has been shown to label the rete epithelium in the testis, epididymis epithelium, and the surface epithelium of the stomach and duodenum. Transitional-cell Carcinomas are positive and Prostate Cancers are negative. This antibody does not recognize intermediate filament proteins, nor does it recognize non-epithelial tissues such as blood vessels, connective tissue, etc.

Cytokeratin 8, also known as type II cytoskeletal 8, is a protein that is often paired with Cytokeratin 18. They are perhaps the most commonly found products of the intermediate filament gene family, and are expressed in single-layer epithelial tissues of the body. Cytokeratin 8 is an intermediate filament protein produced early in embryogenesis.

Anti-Cytokeratin 8 can be used to detect Adenocarcinomas with simple epithelium origin. It can be used to distinguish between Duct (peripheral staining) from Lobular (perinuclear staining) Breast Carcinoma.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: OV-TL 13/30

ISOTYPE: IgG1/K

CONTROL: Salivary Gland, Placenta, Breast, Thyroid, Cervix, Pancreas, Fallopian Tube, Transitional Cell Carcinoma, Lung Adenocarcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM284

ISOTYPE: IgG

CONTROL: Salivary Gland, Placenta, Breast, Thyroid, Cervix, Lung Adenocarcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP17

ISOTYPE: IgG

CONTROL: Colon, Prostate, Kidney, Liver, Colon Carcinoma

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5407 | Tinto Predilute | 3.0 ml |
| BSB 5408 | Tinto Predilute | 7.0 ml |
| BSB 5409 | Tinto Predilute | 15.0 ml |
| BSB 5410 | Concentrate | 0.1 ml |
| BSB 5411 | Concentrate | 0.5 ml |
| BSB 5412 | Concentrate | 1.0 ml |
| BSB 5413 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3764-3 | Tinto Predilute | 3.0 ml |
| BSB-3764-7 | Tinto Predilute | 7.0 ml |
| BSB-3764-15 | Tinto Predilute | 15.0 ml |
| BSB-3764-01 | Concentrate | 0.1 ml |
| BSB-3764-05 | Concentrate | 0.5 ml |
| BSB-3764-1 | Concentrate | 1.0 ml |
| BSB-3764-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6660 | Tinto Predilute | 3.0 ml |
| BSB 6661 | Tinto Predilute | 7.0 ml |
| BSB 6662 | Tinto Predilute | 15.0 ml |
| BSB 6663 | Concentrate | 0.1 ml |
| BSB 6664 | Concentrate | 0.5 ml |
| BSB 6665 | Concentrate | 1.0 ml |
| BSB 6666 | Control Slides | 5 |

Cytokeratin 8 & 18, MAb



IHC of Cytokeratin 8 and 18 on a FFPE Colon Carcinoma Tissue


Cytokeratin 8 belongs to the Type II (basic) subfamily of high molecular-weight keratins and exists in combination with Cytokeratin 18 (Type I [acidic] subfamily of low molecular weight keratins). They are perhaps the most commonly found products of the intermediate filament gene family, and are expressed in single-layer epithelial tissues of the body.

Cytokeratins 8 and 18 can be found in most simple epithelium (e.g., thyroid, female breast, gastrointestinal tract, and respiratory tract). Adenocarcinomas and most Non-keratinizing Squamous Carcinomas will stain, but Keratinizing Squamous Carcinomas will not. This antibody is used when attempting to demonstrate the presence of Paget cells; there is very little keratin 18 in the normal epidermis so only Paget cells will stain. This approach facilitates the interpretation using immunostains and is more sensitive than mucin histochemistry.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: B22.1/B23.1
ISOTYPE: IgG1
CONTROL: Breast, Ovary, GI, Prostate, Pancreas, Salivary Gland
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5414 | Tinto Predilute | 3.0 ml |
| BSB 5415 | Tinto Predilute | 7.0 ml |
| BSB 5416 | Tinto Predilute | 15.0 ml |
| BSB 5417 | Concentrate | 0.1 ml |
| BSB 5418 | Concentrate | 0.5 ml |
| BSB 5419 | Concentrate | 1.0 ml |
| BSB 5420 | Control Slides | 5 |

Cytokeratin 10, RMAb



IHC of Cytokeratin 10 on a FFPE Cervical Carcinoma Tissue


Cytokeratin 10 is a type I cytokeratin, which belongs to the superfamily of intermediate filament (IF) proteins. It is expressed in the subprbasal cell layers of certain stratified epithelia, notably epidermis, and is typically associated with Cytokeratin 1.

Anti-Cytokeratin 10 is helpful in identification of more differentiated squamous cell carcinomas.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP97
ISOTYPE: IgG
CONTROL: Squamous Cell Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6667 | Tinto Predilute | 3.0 ml |
| BSB 6668 | Tinto Predilute | 7.0 ml |
| BSB 6669 | Tinto Predilute | 15.0 ml |
| BSB 6670 | Concentrate | 0.1 ml |
| BSB 6671 | Concentrate | 0.5 ml |
| BSB 6672 | Concentrate | 1.0 ml |
| BSB 6673 | Control Slides | 5 |

Cytokeratin 14, MAb



IHC of Cytokeratin 14 on a FFPE Cervix Tissue


Cytokeratin 14 is a Type I polypeptide found in basal cells of squamous epithelia, some glandular epithelia, myoepithelium, and mesothelial cells. It is usually found as a heterotetramer with two cytokeratin 5 molecules, and a Type II keratin. Together, they form the cytoskeleton of epithelial cells. Mutations in the genes for these cytokeratins are associated with Epidermolysis Bullosa Simplex.

Cytokeratin 14 has been studied as a prognostic marker in Breast Cancer. This antibody labels the basal layer of stratifying squamous and non-squamous epithelia. The staining pattern is cytoplasmic. It recognizes Basal Cell Carcinomas and Squamous Cell Carcinomas. Anti-CK 14 has been demonstrated to be useful in differentiating Squamous Cell Carcinomas from other epithelial tumors. This antibody has also been useful in separating oncocytic tumors of the kidney from renal mimics, as well as in determining metaplastic Carcinomas of the Breast.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: LL002
ISOTYPE: IgG3
CONTROL: Squamous Mucosa, Prostate, Breast, Tonsil, Salivary Gland, Skin, Cervix Carcinoma, Squamous Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6219 | Tinto Predilute | 3.0 ml |
| BSB 6220 | Tinto Predilute | 7.0 ml |
| BSB 6221 | Tinto Predilute | 15.0 ml |
| BSB 6222 | Concentrate | 0.1 ml |
| BSB 6223 | Concentrate | 0.5 ml |
| BSB 6224 | Concentrate | 1.0 ml |
| BSB 6225 | Control Slides | 5 |

Cytokeratin 14, RMAb



IHC of Cytokeratin 14 on FPPE Cervical Cancer Tissue

Cytokeratin 14 is a Type I polypeptide found in basal cells of squamous epithelia, some glandular epithelia, myoepithelium, and mesothelial cells. It is usually found as a heterotetramer with two cytokeratin 5 molecules, and a Type II keratin. Together, they form the cytoskeleton of epithelial cells. Mutations in the genes for these cytokeratins are associated with Epidermolysis Bullosa Simplex.

Cytokeratin 14 has been studied as a prognostic marker in Breast Cancer. This antibody labels the basal layer of stratified squamous and non-squamous epithelia. The staining pattern is cytoplasmic and recognizes Basal Cell Carcinomas and Squamous Cell Carcinomas. Anti-CK14 has been demonstrated to be useful in differentiating Squamous Cell Carcinomas from other epithelial tumors. This antibody has also been useful in separating oncocytic tumors of the kidney from renal mimics, as well as in determining metaplastic Carcinomas of the Breast.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM328


ISOTYPE: IgG

CONTROL: Squamous Mucosa, Prostate, Breast, Tonsil, Salivary Gland, Skin, Cervix Carcinoma, Squamous Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Cytokeratin 16, RMAb



IHC of Cytokeratin 16 on a FPPE Anal Carcinoma Tissue

Cytokeratin 16 (CK16) is a protein that in humans is encoded by the KRT16 gene. Cytokeratin 16 is a type I cytokeratin. It is paired with cytokeratin 6 in a number of epithelial tissues, including nail bed, esophagus, tongue, and hair follicles. Mutations in the gene encoding this protein are associated with the genetic skin disorders pachyonychia congenita, non-epidermolytic palmoplantar keratoderma and unilateral palmoplantar verrucous nevus.

Studies have proposed a modulatory role of CK16 in cell proliferation, suggesting its utility as a marker for proliferation. Rapid induction of CK16 expression near the edge of wounds, upregulation in response to epidermal growth factor stimulus, and overexpression in hyperproliferative disorders, including psoriasis and chronic contact dermatitis, support this assertion. In psoriasis, the severity of disease is correlated with the amount of CK16. Additionally, CK16 expression has been described in neoplasms of multiple tissues. Progressive CK16 abundance and intensity were observed with increased grade of severity of cervical intraepithelial neoplasia lesions. Furthermore, 10% of invasive carcinomas were diffusely or focally positive. In keratocystic odontogenic tumors, CK16 was observed in 79% of cases. These observations support CK16 as a marker of hyperproliferation.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP27


ISOTYPE: IgG

CONTROL: Skin, Prostate, Breast, Cervix, Salivary Gland, SCC

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Cytokeratin 17, MMAb



IHC of Cytokeratin 17 on a FPPE Cervical Carcinoma Tissue

Cytokeratin 17 (CK 17) is a Type I cytokeratin with a MW of 46 kD found sometimes in association with Cytokeratin 7. Cytokeratin 17 is found in nail beds, hair follicles, sebaceous glands, and other epidermal appendages. Mutations in the gene encoding this protein lead to Jackson-Lawler type Pachyonychia Congenita and Steatocystoma Multiplex.

Cytokeratin 17 antibody has been used to distinguish immature Cervical Squamous Metaplasia from High Grade Cervical Intraepithelial Neoplasia (CIN III). Anti-CK 17 also labels myoepithelial cells in the benign breast tissue. CK 17 can be useful when included in a panel of antibodies against TTF-1, napsin A, CK 5&6, p63, and SOX-2 for diagnostic differentiation between lung adenocarcinoma (LADC) and lung squamous cell carcinoma (SCLC), especially for poorly-differentiated lung carcinoma. CK 17 is expressed in SCLC much higher than in LADC. In breast carcinomas, approximately 20% of patients show no expression of ER, PR and Her2, which are defined as triple negative tumor. Eighty-five percent of the triple negative breast carcinomas immunoreact with basal cytokeratins including anti-CK 17. The histologic differentiation of ampullary cancer, intestinal vs. pancreaticobiliary, is very important for treatment. Usually anti-CK 17 and anti-MUC1 immunoreactivity represents pancreaticobiliary subtype whereas anti-MUC2 and anti-CDX-2 positivity defines intestinal subtype.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-33

ISOTYPE: IgG2a/K

CONTROL: Skin, Testis, Breast, Cervix, Cervical Carcinoma, Bladder TCC

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Rat, Goat, Pig

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3762-3 | Tinto Predilute | 3.0 ml |
| BSB-3762-7 | Tinto Predilute | 7.0 ml |
| BSB-3762-15 | Tinto Predilute | 15.0 ml |
| BSB-3762-01 | Concentrate | 0.1 ml |
| BSB-3762-05 | Concentrate | 0.5 ml |
| BSB-3762-1 | Concentrate | 1.0 ml |
| BSB-3762-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2887 | Tinto Predilute | 3.0 ml |
| BSB 2888 | Tinto Predilute | 7.0 ml |
| BSB 2889 | Tinto Predilute | 15.0 ml |
| BSB 2890 | Concentrate | 0.1 ml |
| BSB 2891 | Concentrate | 0.5 ml |
| BSB 2892 | Concentrate | 1.0 ml |
| BSB 2893 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2726 | Tinto Predilute | 3.0 ml |
| BSB 2727 | Tinto Predilute | 7.0 ml |
| BSB 2728 | Tinto Predilute | 15.0 ml |
| BSB 2729 | Concentrate | 0.1 ml |
| BSB 2730 | Concentrate | 0.5 ml |
| BSB 2731 | Concentrate | 1.0 ml |
| BSB 2732 | Control Slides | 5 |

Cytokeratin 17, RMAb




IHC of Cytokeratin 17 on a FFPE Cervical Cancer Tissue

Cytokeratin 17 is a Type I cytokeratin with a MW of 46 kD found sometimes in association with Cytokeratin 7. It is found in nail beds, hair follicles, sebaceous glands, and other epidermal appendages. Mutations in the gene encoding this protein lead to Jackson-Lawler type Pachyonychia Congenita and Steatocystoma Multiplex.

Cytokeratin 17 antibody has been used to distinguish immature Cervical Squamous Metaplasia from high grade Cervical Intraepithelial Neoplasia (CIN III). Anti-CK 17 also labels myoepithelial cells in the benign breast tissue. CK 17 labeling of Breast Carcinoma cells (so-called basal phenotype) has been associated with a poor prognosis.

Cytokeratin 18, MAb




IHC of Cytokeratin 18 on a FFPE Colon Tissue

Cytokeratin 18 is a type I cytokeratin and is typically partnered with Cytokeratin 8. They are expressed in simple and glandular and transitional epithelial cells but not in stratified epithelial cells.

Cytokeratin 18 antibody stains positively in Adenocarcinomas originating from simple and glandular epithelium, and also in poorly differentiated tumor cells of Squamous Carcinoma.

Cytokeratin 19, MAb



IHC of Cytokeratin 19 on a FFPE Thyroid Carcinoma Tissue

Cytokeratin 19 is a Type I cytokeratin. Unlike its related family members, this smallest-known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically found in the periderm, the transiently superficial layer that envelops the developing epidermis.

Anti-Cytokeratin 19 reacts with a wide variety of epithelium and epithelial malignancies including Adenocarcinomas of the colon, stomach, pancreas, biliary tract, liver and breast. Perhaps the most useful application is the identification of Thyroid Carcinoma of the papillary type, although Follicular Carcinoma is also labeled by this antibody approximately 50-60% of the time. Cytokeratin 19 is not expressed in hepatocytes; therefore, this antibody is useful in the identification of liver metastasis. The degree of Cytokeratin 19 positivity in Breast Cancer distinguishes malignant from benign tumors. Cytokeratin 19 is often coexpressed with Cytokeratin 7.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: EP98

ISOTYPE: IgG

CONTROL: Skin, Testis, Breast, Cervix, Colon, Salivary Gland, Cervical Carcinoma, Bladder TCC, Transitional Cell Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Rat, Mouse

ANTIBODY TYPE: Mouse Monoclonal

CLONE: EP30

ISOTYPE: IgG

CONTROL: Liver, Kidney, Breast, GI, Prostate, Fallopian Tube

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-34

ISOTYPE: IgG1/K

CONTROL: Colon, Thyroid Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6184 | Tinto Predilute | 3.0 ml |
| BSB 6185 | Tinto Predilute | 7.0 ml |
| BSB 6186 | Tinto Predilute | 15.0 ml |
| BSB 6187 | Concentrate | 0.1 ml |
| BSB 6188 | Concentrate | 0.5 ml |
| BSB 6189 | Concentrate | 1.0 ml |
| BSB 6190 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6695 | Tinto Predilute | 3.0 ml |
| BSB 6696 | Tinto Predilute | 7.0 ml |
| BSB 6697 | Tinto Predilute | 15.0 ml |
| BSB 6698 | Concentrate | 0.1 ml |
| BSB 6699 | Concentrate | 0.5 ml |
| BSB 6700 | Concentrate | 1.0 ml |
| BSB 6701 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2733 | Tinto Predilute | 3.0 ml |
| BSB 2734 | Tinto Predilute | 7.0 ml |
| BSB 2735 | Tinto Predilute | 15.0 ml |
| BSB 2736 | Concentrate | 0.1 ml |
| BSB 2737 | Concentrate | 0.5 ml |
| BSB 2738 | Concentrate | 1.0 ml |
| BSB 2739 | Control Slides | 5 |

Cytokeratin 19, RMab




IHC of Cytokeratin 19 on a FFPE Tonsil Tissue

Cytokeratin 19 is a Type I cytokeratin. Unlike its related family members, this smallest-known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically found in the periderm, the transiently-superficial layer that envelops the developing epidermis.

Anti-Cytokeratin 19 reacts with a wide variety of epithelium and epithelial malignancies including Adenocarcinomas of the colon, stomach, pancreas, biliary tract, liver and breast. Perhaps the most useful application is the identification of Thyroid Carcinoma of the papillary type, although Follicular Carcinoma is also labeled by this antibody approximately 50-60% of the time. Cytokeratin 19 is not expressed in hepatocytes; therefore, this antibody is useful in the identification of liver metastasis. The degree of Cytokeratin 19 positivity in Breast Cancer distinguishes malignant from benign tumors. Cytokeratin 19 is often coexpressed with Cytokeratin 7.

Cytokeratin 19, RMab




IHC of Cytokeratin 19 on FFPE Thyroid Carcinoma Tissue

Cytokeratin 19 is a Type I cytokeratin. Unlike its related family members, this smallest-known acidic cytokeratin is not paired with a basic cytokeratin in epithelial cells. It is specifically found in the periderm, the transiently superficial layer that envelops the developing epidermis.

Anti-Cytokeratin 19 reacts with a wide variety of epithelium and epithelial malignancies including Adenocarcinomas of the colon, stomach, pancreas, biliary tract, liver and breast. Perhaps the most useful application is the identification of Thyroid Carcinoma of the papillary type, although Follicular Carcinoma is also labeled by this antibody approximately 50-60% of the time. Cytokeratin 19 is not expressed in hepatocytes; therefore, this antibody is useful in the identification of liver metastasis. The degree of Cytokeratin 19 positivity in Breast Cancer distinguishes malignant from benign tumors. Cytokeratin 19 is often coexpressed with Cytokeratin 7.

Cytokeratin 20, MMab



IHC of Cytokeratin 20 on a FFPE Colon Adenocarcinoma Tissue

Cytokeratin 20 (CK 20) is a 46 kDa intermediate filament protein whose expression is restricted primarily to gastric and intestinal epithelium, urothelium, and Merkel cells. Cytokeratin 20 is a Type I cytokeratin. It is a major cellular protein of mature enterocytes and goblet cells found in the gastric and intestinal mucosa.

CK 20 is expressed in Adenocarcinomas of the colon, stomach, pancreas and biliary system. It is also expressed in Mucinous Ovarian Tumors, Transitional-cell Carcinomas of the urinary tract, and Merkel-cell Carcinomas. Cytokeratin 20 is useful in the differentiation of specific types of simple epithelial cells of the urinary tract and normal and malignant-transformed epithelia. This antibody is essentially non-reactive in Squamous Cell Carcinomas and Adenocarcinomas of the Breast, Lung, and Endometrium, Non-mucinous Tumors of the Ovary, and Small-cell Carcinomas. This antibody is often used in conjunction with CK 7 and other antibodies to distinguish Colon Carcinomas (CK20+) from Ovarian, Pulmonary, and Breast Carcinomas.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP72
ISOTYPE: IgG
CONTROL: Colon, Bladder, Thyroid Carcinoma, Colon Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM364
ISOTYPE: IgG
CONTROL: Colon, Bladder, Thyroid Carcinoma, Colon Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Predicted: Mouse


ANTIBODY TYPE: Mouse Monoclonal
CLONE: Ks20.8
ISOTYPE: IgG2a/K
CONTROL: Colon Mucosa, Bladder, Colon Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5379 | Tinto Predilute | 3.0 ml |
| BSB 5380 | Tinto Predilute | 7.0 ml |
| BSB 5381 | Tinto Predilute | 15.0 ml |
| BSB 5382 | Concentrate | 0.1 ml |
| BSB 5383 | Concentrate | 0.5 ml |
| BSB 5384 | Concentrate | 1.0 ml |
| BSB 5385 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3763-3 | Tinto Predilute | 3.0 ml |
| BSB-3763-7 | Tinto Predilute | 7.0 ml |
| BSB-3763-15 | Tinto Predilute | 15.0 ml |
| BSB-3763-01 | Concentrate | 0.1 ml |
| BSB-3763-05 | Concentrate | 0.5 ml |
| BSB-3763-1 | Concentrate | 1.0 ml |
| BSB-3763-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5386 | Tinto Predilute | 3.0 ml |
| BSB 5387 | Tinto Predilute | 7.0 ml |
| BSB 5388 | Tinto Predilute | 15.0 ml |
| BSB 5389 | Concentrate | 0.1 ml |
| BSB 5390 | Concentrate | 0.5 ml |
| BSB 5391 | Concentrate | 1.0 ml |
| BSB 5392 | Control Slides | 5 |

Cytokeratin 20, RMAb



IHC of Cytokeratin 20 on a FFPE Colon Cancer Metastasis to Lung Tissue

Cytokeratin 20 (CK20) is a 46 kDa intermediate filament protein whose expression is restricted primarily to gastric and intestinal epithelium, urothelium, and Merkel cells. Cytokeratin 20 is a Type I cytokeratin. It is a major cellular protein of mature enterocytes and goblet cells found in the gastric and intestinal mucosa.

CK 20 is expressed in Adenocarcinomas of the colon, stomach, pancreas and biliary system. It is also expressed in Mucinous Ovarian Tumors, Transitional-cell Carcinomas of the urinary tract, and Merkel-cell Carcinomas. Cytokeratin 20 is useful in the differentiation of specific types of simple epithelial cells of the urinary tract and normal and malignantly-transformed epithelia. This antibody is essentially non-reactive in Squamous Cell Carcinomas and Adenocarcinomas of the Breast, Lung, and Endometrium, Non-mucinous Tumors of the Ovary, and Small-cell Carcinomas. This antibody is often used in conjunction with CK7 and other antibodies to distinguish Colon Carcinomas (CK20+) from Ovarian, Pulmonary, and Breast Carcinomas.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP23


ISOTYPE: IgG

CONTROL: Colon Mucosa, Bladder, Colon Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, at, Goat, Pig, Marmoset

Cytokeratin 8/35βH11, MAb



IHC of Cytokeratin 8/35βH11 on a FFPE Prostatic Adenocarcinoma Tissue

Cytokeratin 8 belongs to the Type II (basic) subfamily of high molecular-weight keratins and exists in combination with cytokeratin 18. Cytokeratin 8 is primarily found in the non-squamous epithelia and is present in the majority of Adenocarcinomas and Ductal Carcinomas. It is absent in Squamous Cell Carcinomas. Hepatocellular Carcinomas are defined by the use of antibodies that recognize only cytokeratin polypeptides 8 and 18.

Anti-Cytokeratin 8/35βH11 stains most Non-Squamous Epithelial tumors; Squamous tumors are negative for this antibody as a rule. This antibody stains Adenocarcinomas of the breast, ovary, gastrointestinal tract, thyroid, pancreas, bile duct, and salivary glands. This antibody does not react with skeletal muscle or nerve cells.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 35betaH11


ISOTYPE: IgM

CONTROL: Prostate, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Cytokeratin Cocktail AE1 & AE3, MAb



IHC of Cytokeratin AE1 & AE3 on a FFPE Colon Tissue

Cytokeratins are intermediate-filament keratins found in the intracytoplasmic cytoskeleton of epithelial tissue. There are two types of cytokeratins: the low-weight, acidic Type I cytokeratins and the high-weight, basic or neutral Type II cytokeratins. Cytokeratins are usually found in pairs comprising a Type I cytokeratin and a Type II cytokeratin. Expression of these cytokeratins is frequently organ or tissue-specific.

Cytokeratin cocktail AE1/AE3 is well suited to distinguish Epithelial Carcinoma from Non-epithelial malignancies and is used to aid Epithelial Tumor classification. This antibody has been used to characterize the source of various neoplasms and to study the distribution of keratin-containing cells in epithelia during normal development and during the development of epithelial neoplasms. This antibody stains cytokeratins present in normal and abnormal human tissues. This antibody has shown high sensitivity and specificity in recognizing epithelial cells of neoplastic origin.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: AE1/AE3

ISOTYPE: Ig1


CONTROL: Prostate, Colon, Skin, Stomach

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat, Mouse, Rat, Monkey, Rabbit, Chicken, Horse

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6702 | Tinto Predilute | 3.0 ml | BSB 5421 | Tinto Predilute | 3.0 ml | BSB 5428 | Tinto Predilute | 3.0 ml |
| BSB 6703 | Tinto Predilute | 7.0 ml | BSB 5422 | Tinto Predilute | 7.0 ml | BSB 5429 | Tinto Predilute | 7.0 ml |
| BSB 6704 | Tinto Predilute | 15.0 ml | BSB 5423 | Tinto Predilute | 15.0 ml | BSB 5430 | Tinto Predilute | 15.0 ml |
| BSB 6705 | Concentrate | 0.1 ml | BSB 5424 | Concentrate | 0.1 ml | BSB 5431 | Concentrate | 0.1 ml |
| BSB 6706 | Concentrate | 0.5 ml | BSB 5425 | Concentrate | 0.5 ml | BSB 5432 | Concentrate | 0.5 ml |
| BSB 6707 | Concentrate | 1.0 ml | BSB 5426 | Concentrate | 1.0 ml | BSB 5433 | Concentrate | 1.0 ml |
| BSB 6708 | Control Slides | 5 | BSB 5427 | Control Slides | 5 | BSB 5434 | Control Slides | 5 |

Cytokeratin HMW 34βE12, MAb



IHC of Cytokeratin HMW/34βE12 on a FFPE Prostatic Adenocarcinoma Tissue


Cytokeratin 34βE12 is a High Molecular Weight cytokeratin that reacts with all squamous and ductal epithelium and stains carcinomas. This antibody recognizes cytokeratins 1, 5, 10, and 14 that are found in complex epithelia. Cytokeratin 34βE12 shows no reactivity with hepatocytes, pancreatic acinar cells, proximal renal tubules or endometrial glands; there has been no reactivity with cells derived from simple epithelia. Nerve cells, glial cells and mesenchymal tissue such as blood vessels containing only non-keratin types of intermediate filaments are not labelled; however, reactivity with smooth-muscle cells has been occasionally observed.

Mesenchymal Tumors, Lymphomas, Melanomas, Neural Tumors and Neuroendocrine Tumors are unreactive with this antibody. Cytokeratin 34βE12 has been shown to be useful in distinguishing Prostatic Adenocarcinoma from Hyperplasia of the Prostate.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 34BetaE12
ISOTYPE: IgG1/K
CONTROL: Prostate, Cervix
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat, Monkey, Rabbit, Cattle, Horse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5393 | Tinto Predilute | 3.0 ml |
| BSB 5394 | Tinto Predilute | 7.0 ml |
| BSB 5395 | Tinto Predilute | 15.0 ml |
| BSB 5396 | Concentrate | 0.1 ml |
| BSB 5397 | Concentrate | 0.5 ml |
| BSB 5398 | Concentrate | 1.0 ml |
| BSB 5399 | Control Slides | 5 |

Cytokeratin HMW AE3. MAb



IHC of Cytokeratin HMW/AE3 on a FFPE Salivary Gland Tissue


Cytokeratins are intermediate-filament keratins found in the intracytoplasmic cytoskeleton of epithelial tissue. There are two types of cytokeratins: the low-weight, acidic Type I cytokeratins and the high-weight, basic or neutral Type II cytokeratins. Cytokeratins are usually found in pairs comprising a Type I cytokeratin and a Type II cytokeratin. Expression of these cytokeratins is frequently organ or tissue-specific.

Cytokeratin, High Molecular Weight AE3 (HMW, CK 8) is capable of recognizing all basic cytokeratins; therefore, it is a broadly reactive antibody staining most epithelia and their neoplasms. Cytokeratin HMW/AE3 stains normal and neoplastic cells of epithelial origin. CK HMW is primarily found in the non-squamous epithelia and is present in the majority of Adenocarcinomas and Ductal Carcinomas. It is absent in Squamous Cell Carcinomas. Hepatocellular Carcinomas are defined by the use of antibodies that recognize only cytokeratin 8 and 18.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: AE3
ISOTYPE: IgG1
CONTROL: Prostate, Salivary Gland, Bladder
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5442 | Tinto Predilute | 3.0 ml |
| BSB 5443 | Tinto Predilute | 7.0 ml |
| BSB 5444 | Tinto Predilute | 15.0 ml |
| BSB 5445 | Concentrate | 0.1 ml |
| BSB 5446 | Concentrate | 0.5 ml |
| BSB 5447 | Concentrate | 1.0 ml |
| BSB 5448 | Control Slides | 5 |

Cytokeratin LMW/AE1, MAb



IHC of Cytokeratin LMW/AE1 on a FFPE Salivary Gland Tissue


Cytokeratins are intermediate-filament keratins found in the intracytoplasmic cytoskeleton of epithelial tissue. There are two types of cytokeratins: the low-weight, acidic Type I cytokeratins and the high-weight, basic or neutral Type II cytokeratins. Cytokeratins are usually found in pairs comprising a Type I cytokeratin and a Type II cytokeratin. Expression of these cytokeratins is frequently organ or tissue-specific.

Cytokeratin Low Molecular Weight AE1 can recognize most acidic keratins, making it a broadly reactive antibody that stains most epithelia and their neoplasms. Members of the acidic and basic subfamilies are found in pairs. Each epithelium contains at least one acidic and one basic keratin so this antibody can show the distribution of keratin-containing cells in epithelia. Cytokeratin AE1 is particularly suited to distinguish poorly-differentiated Carcinomas from non-epithelial Neoplasms. This marker stains both normal and neoplastic cells of epithelial origin.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: AE1
ISOTYPE: IgG1
CONTROL: Prostate, Salivary Gland, Bladder, Breast, Kidney, Pancreas, Cervix
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat, Monkey, Rabbit, Chicken, Cattle, Shrew, Fish

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5435 | Tinto Predilute | 3.0 ml |
| BSB 5436 | Tinto Predilute | 7.0 ml |
| BSB 5437 | Tinto Predilute | 15.0 ml |
| BSB 5438 | Concentrate | 0.1 ml |
| BSB 5439 | Concentrate | 0.5 ml |
| BSB 5440 | Concentrate | 1.0 ml |
| BSB 5441 | Control Slides | 5 |

Cytokeratin LMW CAM5.2, MMab



IHC of Cytokeratin LMW CAM5.2 on a FFPE Breast Carcinoma Tissue


Anti-Cytokeratin (CAM5.2) antibody has a primary reactivity with human keratin proteins that correspond to Moll's peptides #7 and #8, Mr 48 and 52 kDa, respectively. Cytokeratin 7 and 8 are present in secretory epithelia of normal human tissue but not on stratified squamous epithelium.

Anti-Cytokeratin (CAM5.2) stains most epithelial-derived tissue, including liver, renal tubular epithelium, and hepatocellular and renal cell carcinomas. Anti-Cytokeratin (CAM 5.2) may not react with some squamous cell carcinomas.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: CAM5.2
ISOTYPE: IgG2a/K
CONTROL: Colon, Breast, Ovarian Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2056 | Tinto Predilute | 3.0 ml |
| BSB 2057 | Tinto Predilute | 7.0 ml |
| BSB 2058 | Tinto Predilute | 15.0 ml |
| BSB 2059 | Concentrate | 0.1 ml |
| BSB 2060 | Concentrate | 0.5 ml |
| BSB 2061 | Concentrate | 1.0 ml |
| BSB 2062 | Control Slides | 5 |

Cytokeratin, MNF116, MMab



IHC of Cytokeratin MNF116 on a Frozen Basal Cell Carcinoma Tissue


Cytokeratin MNF116 is a broad-spectrum anti-cytokeratin reacting with intermediate and low-molecular-weight keratins, ranging from 40 through 58 kD, corresponding to cytokeratin 5, 6, 8, 17 and 19. It shows a broad pattern of reactivity with human epithelial tissues from simple glandular epithelial to stratified squamous epithelia, like epidermis, mammary gland ducts, and tracheal epithelium.

Cytokeratin MNF116 is a useful aid for the classification of neoplasms of epithelial origin including Squamous Cell Carcinoma, Small Cell Carcinoma, Sarcomatoid Carcinoma, Spindle Cell Carcinoma, Epithelioid and Spindle Cell component of Malignant Mesothelioma and Adenocarcinoma. A wide range of soft tissue tumors are also positive with cytokeratin MNF116: monophasic and biphasic Synovial Sarcoma, vascular neoplasms including Epithelioid Hemangioendothelioma, Epithelioid Angiosarcoma, Epithelioid Sarcoma. Desmoplastic Small Round Cell Tumors require cytokeratin positivity for diagnosis. Smooth muscle tumors and Plasmacytoma may demonstrate aberrant expression of cytokeratin MNF116

ANTIBODY TYPE: Mouse Monoclonal
CLONE: MNF116
ISOTYPE: IgG1/k
CONTROL: Breast, Cervix, Skin, Colon, Colorectal, Gastric, Breast and Prostatic Carcinomas
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3532 | Tinto Predilute | 3.0 ml |
| BSB 3533 | Tinto Predilute | 7.0 ml |
| BSB 3534 | Tinto Predilute | 15.0 ml |
| BSB 3535 | Concentrate | 0.1 ml |
| BSB 3536 | Concentrate | 0.5 ml |
| BSB 3537 | Concentrate | 1.0 ml |
| BSB 3538 | Control Slides | 5 |

Cytokeratin OSCAR, MMab



IHC of CK OSCAR on a FFPE Colon Tissue


Anti-Cytokeratin OSCAR is well-suited to distinguish Epithelial Carcinoma from Non-epithelial malignancies and is used to aid Epithelial Tumor classification. Anti-Cytokeratin OSCAR identifies a number of bands corresponding to cytokeratins 7, 8, 18 and 19 (additional bands – cytokeratins – may also be detected). This antibody has been used to characterize the source of various neoplasms and to study the distribution of keratin-containing cells in epithelia during normal development and during the development of epithelial neoplasms.

In normal tissues, OSCAR is reactive with most epithelial types tested including bile ducts and hepatocytes in liver, bladder epithelium, breast ducts, bronchial epithelium, endometrium, follicular dendritic cells of lymph node and tonsil, intestinal epithelium of the stomach, duodenum, ileum, colon, rectum, pancreas, ovarian epithelium, pancreatic acini, pituitary acini, pneumocytes, prostate, thyroid and skin. In tumors, OSCAR is reactive with most Carcinomas including Breast, TCC, RCC, Lung, Endometrial CA, Prostate CA, Ovarian CA, HCC, Colorectal CA, Stomach CA and Thyroid CA. It is negative in certain normal tissues including brain, lymphocytes and all cells of hemolymphoid origin, muscle, brain, nerves, endothelium and in certain tumors including Melanoma, Sarcoma, Lymphoma, PNET/Ewing's and GIST. This antibody has shown high sensitivity in recognizing epithelial cells and carcinomas.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: OSCAR
ISOTYPE: IgG2a
CONTROL: Breast, Liver, GI, Prostate, Colon, Skin, Stomach, Pancreas, Tonsil
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6177 | Tinto Predilute | 3.0 ml |
| BSB 6178 | Tinto Predilute | 7.0 ml |
| BSB 6179 | Tinto Predilute | 15.0 ml |
| BSB 6180 | Concentrate | 0.1 ml |
| BSB 6181 | Concentrate | 0.5 ml |
| BSB 6182 | Concentrate | 1.0 ml |
| BSB 6183 | Control Slides | 5 |

Cytomegalovirus, MAb



IHC and IF of CMV on a FFPE Infected Lung Tissue

Cytomegalovirus (CMV) is a virus of the Herpesvirus group; in humans it is commonly known as HCMV or Human Herpesvirus 5 (HHV-5). CMV belongs to the Betaherpesvirinae subfamily of Herpesviridae, which also includes Roseolovirus. CMV especially attacks salivary glands. CMV infection can also be life-threatening for patients who are immunocompromised (e.g., patients with HIV or organ-transplant recipients). CMV viruses are found in many mammal species, but CMV species isolated from animals differ from human CMV in terms of genomic structure, and have not been reported to cause human disease.

This Anti-cytomegalovirus antibody cocktail reacts with two different epitopes. The DDG9 antibody reacts with a 76 kDa protein produced by CMV. CCH2 antibody reacts with the early DNA-binding protein p52. There is no cross-reactivity with other Herpesviruses or Adenoviruses. CMV infection is usually seen in immunocompromised patients and involves the GI tract, lung, heart and liver, as well as other organs.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 8B1.2, 1G5.2 & 2D4
ISOTYPE: IgG2a
CONTROL: CMV Infected Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Desmin, MAb



IHC of Desmin on a FFPE Skeletal Muscle Tissue


Desmin is a type of intermediate filament found near the Z line in sarcomeres. Both vimentin and desmin are characteristics of mesenchymal cells.

Desmin antibody detects a protein that is expressed by cells of normal smooth, skeletal and cardiac muscles. Light microscopy studies of Desmin have suggested that it is primarily located at or near the periphery of Z lines in striated muscle fibrils. In smooth muscle, Desmin interconnects cytoplasmic dense bodies with membrane-bound dense plaques. Desmin antibody reacts with Leiomyomas, Rhabdomyomas, and Perivascular cells of Glomus Tumors of the skin (if they are of myogenic nature). This antibody is used to demonstrate the myogenic components/derivation of tumors.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: D33
ISOTYPE: IgG1/K
CONTROL: Skeletal Muscle
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat, Mouse, Rat, Chicken, Sheep, Hamster

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5449 | Tinto Predilute | 3.0 ml |
| BSB 5450 | Tinto Predilute | 7.0 ml |
| BSB 5451 | Tinto Predilute | 15.0 ml |
| BSB 5452 | Concentrate | 0.1 ml |
| BSB 5453 | Concentrate | 0.5 ml |
| BSB 5454 | Concentrate | 1.0 ml |
| BSB 5455 | Control Slides | 5 |

Desmin, RMAb



IHC of Desmin on a FFPE Skeletal Muscle Tissue


Desmin is a class III intermediate filament found near the Z line in sarcomeres. Both vimentin and desmin are characteristics of mesenchymal cells.

Desmin antibody detects a protein that is expressed by cells of normal smooth, skeletal and cardiac muscles. Light microscopy studies of desmin suggests that it is primarily located at or near the periphery of Z lines in striated muscle fibrils. In smooth muscle, Desmin interconnects cytoplasmic dense bodies with membrane bound dense plaques. Desmin antibody reacts with Leiomyomas, Rhabdomyomas, and Perivascular cells of Glomus Tumors of the skin (if they are of myogenic nature). This antibody is used to demonstrate the myogenic components/derivation of tumors.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP15
ISOTYPE: IgG
CONTROL: Skeletal Muscle, Placenta, Colon, Prostate, Skin, Fallopian Tube
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat, Guinea Pig

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6709 | Tinto Predilute | 3.0 ml |
| BSB 6710 | Tinto Predilute | 7.0 ml |
| BSB 6711 | Tinto Predilute | 15.0 ml |
| BSB 6712 | Concentrate | 0.1 ml |
| BSB 6713 | Concentrate | 0.5 ml |
| BSB 6714 | Concentrate | 1.0 ml |
| BSB 6715 | Control Slides | 5 |

Desmoglein-3, RMAb




IHC of Desmoglein-3 on a FFPE Squamous Cell Carcinoma Tissue

Desmoglein 3 is a protein that is encoded by the DSG3 gene on Chromosome 18. Desmoglein 3 is a calcium-binding transmembrane glycoprotein component of desmosomes in vertebrate epithelial cells. Currently, four desmoglein subfamily members have been identified and all are members of the cadherin cell adhesion molecule superfamily. This protein has been identified as the autoantigen of the autoimmune skin blistering disease pemphigus vulgaris.

Desmoglein 3 has been cited as a superior marker for Lung Squamous Cell Carcinomas, and helps distinguish lung squamous cell carcinoma cases from lung adenocarcinoma. Studies have also shown that a panel consisting of Desmoglein-3 utilized with Napsin A can be a useful immunohistochemical marker for differentiation of lung squamous cell carcinoma and adenocarcinoma from other subtypes. Lung cases that are typically positive for Desmoglein 3 tend to have a poor clinical outcome.

DOG1, RMAb




IHC of DOG-1 on a FFPE GIST Tissue

DOG1 (discovered on GIST 1), a cell-surface protein of unknown function, is expressed strongly on the cell surface of GISTs and is rarely expressed in other soft tissue tumors. Among GIST cases with c-Kit mutations, the DOG1 antibody identified 11% more cases than a c-Kit antibody.

DOG1 identifies the vast majority of both c-Kit negative and PDGFRA mutated GIST cases that may still benefit from imatinib mesylate (Gleevec), an inhibitor of the kit tyrosine kinase. In addition, DOG1 immunoreactivity is seen in fewer cases of mesenchymal and epithelial tumors, and melanomas when compared with c-Kit. The use of this highly-sensitive and specific novel marker should increase the accuracy of GIST diagnosis.

DOG-1, RMAb



IHC of DOG1 on a FFPE GIST Tissue

DOG1 (discovered on GIST 1), a cell-surface protein of unknown function, is expressed strongly on the cell surface of GISTs and is rarely expressed in other soft tissue tumors. Among GIST cases with c-Kit mutations, the DOG1 antibody identified 11% more cases than a c-Kit antibody.

DOG1 identifies the vast majority of both c-Kit negative and PDGFRA mutated GIST cases that may still benefit from imatinib mesylate (Gleevec), an inhibitor of the kit tyrosine kinase. In addition, DOG1 immunoreactivity is seen in fewer cases of mesenchymal and epithelial tumors, and melanomas when compared with c-Kit. The use of this highly-sensitive and specific novel marker should increase the accuracy of GIST diagnosis.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP306

ISOTYPE: IgG

CONTROL: Tonsil, Skin, Cervix, Cervical Carcinoma, Lung Squamous Cell Carcinoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-DOG1

ISOTYPE: IgG

CONTROL: GIST

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP332

ISOTYPE: IgG

CONTROL: Salivary Gland, Breast, GIST

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2894 | Tinto Predilute | 3.0 ml |
| BSB 2895 | Tinto Predilute | 7.0 ml |
| BSB 2896 | Tinto Predilute | 15.0 ml |
| BSB 2897 | Concentrate | 0.1 ml |
| BSB 2898 | Concentrate | 0.5 ml |
| BSB 2899 | Concentrate | 1.0 ml |
| BSB 2900 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2901 | Tinto Predilute | 3.0 ml |
| BSB 2902 | Tinto Predilute | 7.0 ml |
| BSB 2903 | Tinto Predilute | 15.0 ml |
| BSB 2904 | Concentrate | 0.1 ml |
| BSB 2905 | Concentrate | 0.5 ml |
| BSB 2906 | Concentrate | 1.0 ml |
| BSB 2907 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6268 | Tinto Predilute | 3.0 ml |
| BSB 6269 | Tinto Predilute | 7.0 ml |
| BSB 6270 | Tinto Predilute | 15.0 ml |
| BSB 6271 | Concentrate | 0.1 ml |
| BSB 6272 | Concentrate | 0.5 ml |
| BSB 6273 | Concentrate | 1.0 ml |
| BSB 6274 | Control Slides | 5 |

E-Cadherin, RMAb




IHC of E-Cadherin on a FFPE Colon Tissue

EGFR, MMAb



IHC of EGFR on a FFPE Colon Carcinoma Tissue

EGFR Phospho, RMAb



IHC of EGFR Phospho on a FFPE Cervix Tissue

Cadherins are a class of transmembrane proteins. They play an important role in cell adhesion by ensuring cells within tissues are bound together. E-Cadherin is an adhesion protein that is expressed in cells of epithelial lineage. It stains positively in glandular epithelium as well as Adenocarcinomas of the lung and G.I. tract, and ovary. E-Cadherin has been useful in distinguishing Adenocarcinoma from Mesothelioma. It has also been shown to be positive in some Thyroid Carcinomas. It can be used to differentiate Ductal Carcinomas (positive for E-Cadherin) from Lobular Breast Carcinomas.

Loss of E-Cadherin function or expression has been implicated in cancer progression and metastasis. E-Cadherin downregulation decreases the strength of cellular adhesion within a tissue, resulting in an increase in cellular motility. This may then allow cancer cells to cross the basement membrane and invade surrounding tissues. Loss of E-Cadherin expression has been suggested as a poor prognostic sign in Breast Carcinoma and Non-Small Cell Lung Carcinomas.

Epidermal Growth Factor Receptor (EGFR) is the receptor for epidermal growth factor (EGF). It is a member of the ErbB family receptors, a subfamily of four closely related receptor tyrosine kinases: EGFR (ErbB-1), HER-2 neu (ErbB-2), HER-3 (ErbB-3) and HER-4 (ErbB-4).

Epidermal Growth Factor Receptor (EGFR) is the receptor for epidermal growth factor (EGF). It is a member of the ErbB family receptors, a subfamily of four closely related receptor tyrosine kinases: EGFR (ErbB-1), HER-2 neu (ErbB-2), HER-3 (ErbB-3) and HER-4 (ErbB-4).

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP6

ISOTYPE: IgG

CONTROL: Breast, Colon, Cervix, Pancreas, Lung, Ovary, GI Tract Adenocarcinoma, Breast Carcinoma

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 31G7

ISOTYPE: IgG1

CONTROL: Skin, Placenta, Testis, Tonsil, Pancreas, Squamous Cell Carcinoma

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP11

ISOTYPE: IgG

CONTROL: Skin, Placenta, Testis, Tonsil, Pancreas, Squamous Cell Carcinoma

LOCALIZATION: Cell Membrane


SPECIES REACTIVITY: Human, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5463 | Tinto Predilute | 3.0 ml |
| BSB 5464 | Tinto Predilute | 7.0 ml |
| BSB 5465 | Tinto Predilute | 15.0 ml |
| BSB 5466 | Concentrate | 0.1 ml |
| BSB 5467 | Concentrate | 0.5 ml |
| BSB 5468 | Concentrate | 1.0 ml |
| BSB 5469 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5470 | Tinto Predilute | 3.0 ml |
| BSB 5471 | Tinto Predilute | 7.0 ml |
| BSB 5472 | Tinto Predilute | 15.0 ml |
| BSB 5473 | Concentrate | 0.1 ml |
| BSB 5474 | Concentrate | 0.5 ml |
| BSB 5475 | Concentrate | 1.0 ml |
| BSB 5476 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6716 | Tinto Predilute | 3.0 ml |
| BSB 6717 | Tinto Predilute | 7.0 ml |
| BSB 6718 | Tinto Predilute | 15.0 ml |
| BSB 6719 | Concentrate | 0.1 ml |
| BSB 6720 | Concentrate | 0.5 ml |
| BSB 6721 | Concentrate | 1.0 ml |
| BSB 6722 | Control Slides | 5 |

EMA, MAb




IHC of EMA on a FFPE Ovarian Adenocarcinoma Tissue

Epithelial Membrane Antigen (EMA) antibody is a mucin-like glycoprotein, shown to be useful as a pan-epithelial marker for detecting early metastatic loci of carcinoma in the bone marrow or liver. It stains normal and neoplastic cells from various tissues, including mammary epithelium, sweat glands and squamous epithelium.

Hepatocellular Carcinoma, Adrenal Carcinoma and Embryonal Carcinomas are consistently EMA negative, so keratin positivity with negative EMA favors one of these tumors. EMA is frequently positive in meningioma, which can be useful when distinguishing it from other intracranial neoplasms. The absence of EMA can also be of value since negative EMA is characteristic of some tumors including Adrenal Carcinoma, Seminomas, Paraganglioma and Hepatoma.

EpCAM/Epithelial Specific Antigen, MAb




IHC of EpCAM on a FFPE Colon Carcinoma Tissue

Epithelial Cell Adhesion Molecule (EpCAM) or Epithelial Specific Antigen is a 40kD cell surface antigen that is broadly distributed in epithelial cells and displays a highly conserved expression in carcinomas. These glycoproteins are located on the cell membrane surface and in the cytoplasm of virtually all epithelial cells, with the exception of most squamous epithelia, hepatocytes, renal proximal tubular cells, gastric parietal cells and myoepithelial cells. However, focal positivity may be seen in the basal layer of squamous cell epithelium of endoderm (e.g., palatine tonsils) and mesoderm (e.g., uterine cervix).

EpCAM expression has been reported to be a possible marker of early malignancy, with expression being increased in tumor cells, and de novo expression being seen in dysplastic squamous epithelium. Epithelial specific antigen has been known to play an important role as a tumor-cell marker in lymph nodes from patients with esophageal carcinoma. EpCAM can be used to distinguish among Basal Cell, Basosquamous Carcinomas and Squamous Cell Carcinomas of the skin.

EpCAM/Epithelial Specific Antigen, MAb



IHC of EpCAM on a FFPE Breast Cancer Tissue

Epithelial Cell Adhesion Molecule (EpCAM) or Epithelial Specific Antigen is a 40kD cell surface antigen that is broadly distributed in epithelial cells and displays a highly conserved expression in carcinomas. These glycoproteins are located on the cell membrane surface and in the cytoplasm of virtually all epithelial cells, with the exception of most squamous epithelia, hepatocytes, renal proximal tubular cells, gastric parietal cells and myoepithelial cells. However, focal positivity may be seen in the basal layer of squamous cell epithelium of endoderm (e.g., palatine tonsils) and mesoderm (e.g., uterine cervix).

EpCAM expression has been reported to be a possible marker of early malignancy, with expression being increased in tumor cells, and de novo expression being seen in dysplastic squamous epithelium. Epithelial specific antigen has been known to play an important role as a tumor-cell marker in lymph nodes from patients with esophageal carcinoma. EpCAM can be used to distinguish among Basal Cell, Basosquamous Carcinomas and Squamous Cell Carcinomas of the skin.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: E29
ISOTYPE: IgG2a/K
CONTROL: Breast, Skin, Colon, Kidney, Cervix
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: Ber-EP4
ISOTYPE: IgG1/K
CONTROL: Colon, Cervix, Salivary Gland, Pancreas, Breast, Thyroid, Liver, Adenocarcinomas
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: MOC-31
ISOTYPE: IgG1/K
CONTROL: Adenocarcinomas
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5477 | Tinto Predilute | 3.0 ml |
| BSB 5478 | Tinto Predilute | 7.0 ml |
| BSB 5479 | Tinto Predilute | 15.0 ml |
| BSB 5480 | Concentrate | 0.1 ml |
| BSB 5481 | Concentrate | 0.5 ml |
| BSB 5482 | Concentrate | 1.0 ml |
| BSB 5483 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6275 | Tinto Predilute | 3.0 ml |
| BSB 6276 | Tinto Predilute | 7.0 ml |
| BSB 6277 | Tinto Predilute | 15.0 ml |
| BSB 6278 | Concentrate | 0.1 ml |
| BSB 6279 | Concentrate | 0.5 ml |
| BSB 6280 | Concentrate | 1.0 ml |
| BSB 6281 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6723 | Tinto Predilute | 3.0 ml |
| BSB 6724 | Tinto Predilute | 7.0 ml |
| BSB 6725 | Tinto Predilute | 15.0 ml |
| BSB 6726 | Concentrate | 0.1 ml |
| BSB 6727 | Concentrate | 0.5 ml |
| BSB 6728 | Concentrate | 1.0 ml |
| BSB 6729 | Control Slides | 5 |

Epstein Barr Virus LMP-1, MAb




IHC of Epstein Barr Virus on a FFPE Hodgkin's Lymphoma Tissue

ERCC1, RMAb



IHC of ERCC1 on a FFPE Lymphoblastic Lymphoma Tissue

ERG, RMAb



IHC of ERG on a FFPE Prostate Carcinoma Tissue

The Epstein-Barr virus (EBV), also called Human Herpesvirus 4 (HHV-4), is a virus of the Herpes family, and is one of the most common viruses in humans. The virus can execute many distinct programs of gene expression, which can be broadly categorized as being lytic cycle or latent cycle. The lytic cycle, or productive infection, results in staged expression of several viral proteins with the ultimate objective of producing infectious virions. The latent cycle (lysogenic) programs are those that do not result in production of virions. A very limited, distinct set of viral proteins are produced during latent cycle infection. These include Epstein-Barr nuclear antigens EBNA-1, EBNA-2, EBNA-3A, EBNA-3B, EBNA-3C, EBNA-leader protein (EBNA-LP), latent membrane proteins LMP-1, LMP-2A and LMP-2B and the Epstein-Barr encoded RNAs (EBERs). In addition, EBV codes for at least twenty microRNAs which are expressed in latently infected cells.

DNA excision repair protein ERCC1 is encoded by the ERCC1 gene on chromosome 19. Excision Repair Cross Complementing 1 (ERCC1) is a mammalian nucleotide excision repair (NER) enzyme involved in repair of damaged DNA. ERCC1 is a homologous to RAD10 in *Saccharomyces cerevisiae*, which is required in mitotic intrachromosomal recombination and repair.

ERCC1 is required in repair of cisplatin-induced DNA adducts and ultraviolet (UV)-induced DNA damage. High expression of ERCC1 has been linked to tumor progression in a variety of cancers including non-small cell lung cancer (NSCLC), squamous cell carcinoma of the head, ovarian cancer and esophageal cancer. Increased levels of ERCC1 expression may correlate with a lower response to platinum based chemotherapies. Studies have also showed however, that in non-small cell lung carcinoma (NSCLC), surgically removed tumors that receive no further therapy have a better survival if ERCC1-positive than if ERCC1-negative.

ERG belongs to the ETS family that plays important roles in cell development, differentiation, proliferation, apoptosis and tissue remodeling. The aberrant expression of several ETS proteins is involved in tumor development and progression. ERG is linked to normal processes such as mesoderm formation. TMPRSS2-ERG fusion, which occurs on account of translocations and interstitial deletions, is implicated in aggressive forms of prostate cancer.

ERG overexpression is associated with aggressive tumor behavior and patient survival in prostate cancer. ERG antibody labels endothelial cells, lymphocytes, and prostate cancer cells.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: CS1-4
ISOTYPE: IgG1
CONTROL: EBV Infected Tissue, Hodgkin's Lymphoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP219
ISOTYPE: IgG
CONTROL: Tonsil, Testis, Breast, Prostate, Fallopian Tube, Ovarian Carcinoma, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP111
ISOTYPE: IgG
CONTROL: Prostate, Colon, Kidney, Fallopian Tube, Tonsil, Myometrium, Skin, Brain, Breast
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 5484 | Tinto Predilute | 3.0 ml | BSB 2482 | Tinto Predilute | 3.0 ml | BSB 6737 | Tinto Predilute | 3.0 ml |
| BSB 5485 | Tinto Predilute | 7.0 ml | BSB 2483 | Tinto Predilute | 7.0 ml | BSB 6738 | Tinto Predilute | 7.0 ml |
| BSB 5486 | Tinto Predilute | 15.0 ml | BSB 2484 | Tinto Predilute | 15.0 ml | BSB 6739 | Tinto Predilute | 15.0 ml |
| BSB 5487 | Concentrate | 0.1 ml | BSB 2485 | Concentrate | 0.1 ml | BSB 6740 | Concentrate | 0.1 ml |
| BSB 5488 | Concentrate | 0.5 ml | BSB 2486 | Concentrate | 0.5 ml | BSB 6741 | Concentrate | 0.5 ml |
| BSB 5489 | Concentrate | 1.0 ml | BSB 2487 | Concentrate | 1.0 ml | BSB 6742 | Concentrate | 1.0 ml |
| BSB 5490 | Control Slides | 5 | BSB 2488 | Control Slides | 5 | BSB 6743 | Control Slides | 5 |

Estrogen Receptor, MAb




IHC of Estrogen Receptor on a FFPE Fallopian Tube Tissue

Estrogen receptor (ER) is a nuclear receptor for estrogens such as estradiol (the main endogenous human estrogen). The two different estrogen receptor proteins produced from the ESR1 and ESR2 genes are usually called the alpha and beta receptors.

This ER antibody recognizes a protein of 67 kDa, which is identified as estrogen receptor alpha.


Estrogen Receptor, RMAb



IHC of Estrogen Receptor on a FFPE Breast Tissue

Estrogen receptor (ER) is a nuclear receptor for estrogens such as estradiol (the main endogenous human estrogen). The two different estrogen receptor proteins produced from the ESR1 and ESR2 genes are usually called the alpha and beta receptors. This ER antibody recognizes a protein of 67 kDa, which is identified as estrogen receptor (ER) alpha.

Estrogen Receptor, RMAb



IHC of Estrogen Receptor on a FFPE Breast Carcinoma Tissue

Estrogen receptor alpha (ER) is a nuclear receptor for estrogens such as estradiol (the main endogenous human estrogen). The two different estrogen receptor proteins produced from the ESR1 and ESR2 genes are usually called the alpha and beta receptors.

This ER antibody recognizes a protein of 67 kDa, which is identified as estrogen receptor (ER) alpha.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-1
ISOTYPE: IgG1/K
CONTROL: Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT11
ISOTYPE: IgG
CONTROL: Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP1
ISOTYPE: IgG
CONTROL: Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2740 | Tinto Predilute | 3.0 ml |
| BSB 2741 | Tinto Predilute | 7.0 ml |
| BSB 2742 | Tinto Predilute | 15.0 ml |
| BSB 2743 | Concentrate | 0.1 ml |
| BSB 2744 | Concentrate | 0.5 ml |
| BSB 2745 | Concentrate | 1.0 ml |
| BSB 2746 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5491 | Tinto Predilute | 3.0 ml |
| BSB 5492 | Tinto Predilute | 7.0 ml |
| BSB 5493 | Tinto Predilute | 15.0 ml |
| BSB 5494 | Concentrate | 0.1 ml |
| BSB 5495 | Concentrate | 0.5 ml |
| BSB 5496 | Concentrate | 1.0 ml |
| BSB 5497 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2489 | Tinto Predilute | 3.0 ml |
| BSB 2490 | Tinto Predilute | 7.0 ml |
| BSB 2491 | Tinto Predilute | 15.0 ml |
| BSB 2492 | Concentrate | 0.1 ml |
| BSB 2493 | Concentrate | 0.5 ml |
| BSB 2494 | Concentrate | 1.0 ml |
| BSB 2495 | Control Slides | 5 |

Estrogen Receptor, RMAb




IHC of Estrogen Receptor on a FFPE Fallopian Tube Tissue

Estrogen receptor alpha (ER) is a nuclear receptor for estrogens such as estradiol (the main endogenous human estrogen). The two different estrogen receptor proteins produced from the ESR1 and ESR2 genes are usually called the alpha and beta receptors. This ER antibody recognizes a protein of 67 kDa, which is identified as estrogen receptor (ER) alpha.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM292
ISOTYPE: IgG
CONTROL: Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Factor H/Complement Factor H, MAb




IHC of Factor H a FFPE SARS-CoV-2 infected Lung Tissue

Factor H or Complement Factor H (CFH), is the major soluble inhibitor of complement, where its binding to self markers (i.e. particular glycan structures) prevents complement activation and amplification on host surfaces. Mutations and polymorphisms that affect recognition of self markers by Factor H are associated with diseases of complement dysregulation, such as age-related macular degeneration and atypical hemolytic uremic syndrome. In addition, pathogens and cancer cells can hijack Factor H to evade the immune response.

Lung, Ovarian, Glial and Colon Cancer cells show enhanced expression and surface binding of soluble regulators, including Factor H. Factor H has been shown to be expressed by human Breast Cancer cells, which correlates with the presence of immunosuppressive macrophages, Breast Cancer recurrence and severity of the disease. Lung cancer cells may develop a protective mechanism against complement attack by expressing and binding Factor H to their cell membranes. Additionally, it has been demonstrated that Factor H is upregulated by constitutive activation of STAT4, which is accounted for by SOCS silencing in Lung Cancer cells. Several studies have also suggested the importance of Factor H in the protection of tumor cells against complement activation. The importance of Factor H expression for the protection of cancer cells in vivo will help to elucidate the mechanisms used by tumor cells to avoid complement activity and assist in the design of more efficient complement-mediated immunotherapies.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-164
ISOTYPE: IgG1
CONTROL: Testis, Liver, Kidney, Pancreas, Adrenal Gland
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

Factor VIII-Related Antigen, RPAb



IHC of Factor VIII on a FFPE Placenta Tissue

Factor VIII (F VIII) is an essential clotting factor. The lack of normal F VIII causes Hemophilia A, an inherited bleeding disorder. FVIII is a glycoprotein procofactor synthesized and released into the bloodstream by the liver. In the circulating blood, it is mainly bound to von Willebrand factor (vWF, also known as Factor VIII-related antigen) to form a stable complex. Upon activation by thrombin or Factor Xa, it dissociates from the complex to interact with Factor IXa, the coagulation cascade. It is a cofactor to Factor IXa in the activation of Factor X, which, in turn, with its cofactor Factor Va, activates more thrombin. Thrombin cleaves fibrinogen into fibrin which polymerizes and crosslinks (using Factor XIII) into a blood clot.

This antibody reacts with endothelial cells in normal, reactive, and neoplastic blood cells. F VIII antibody has helped to establish the endothelial nature of some lesions of disputed histogenesis, e.g., Kaposi's Sarcoma and Cardiac Myxoma. Not all endothelial cells synthesize (or store) this molecule; therefore, it should not be surprising that not all tumors of endothelial differentiation (benign or malignant) react with this antigen.


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: N/A
CONTROL: Skin, Placenta
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3765-3 | Tinto Predilute | 3.0 ml |
| BSB-3765-7 | Tinto Predilute | 7.0 ml |
| BSB-3765-15 | Tinto Predilute | 15.0 ml |
| BSB-3765-01 | Concentrate | 0.1 ml |
| BSB-3765-05 | Concentrate | 0.5 ml |
| BSB-3765-1 | Concentrate | 1.0 ml |
| BSB-3765-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3723-3 | Tinto Predilute | 3.0 ml |
| BSB-3723-7 | Tinto Predilute | 7.0 ml |
| BSB-3723-15 | Tinto Predilute | 15.0 ml |
| BSB-3723-01 | Concentrate | 0.1 ml |
| BSB-3723-05 | Concentrate | 0.5 ml |
| BSB-3723-1 | Concentrate | 1.0 ml |
| BSB-3723-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5498 | Tinto Predilute | 3.0 ml |
| BSB 5499 | Tinto Predilute | 7.0 ml |
| BSB 5500 | Tinto Predilute | 15.0 ml |
| BSB 5501 | Concentrate | 0.1 ml |
| BSB 5502 | Concentrate | 0.5 ml |
| BSB 5503 | Concentrate | 1.0 ml |
| BSB 5504 | Control Slides | 5 |

Factor XIIIa, RMab



IHC of Factor XIIIa on a FFPE Placenta Tissue


Factor XIII or fibrin stabilizing factor is an enzyme of the blood coagulation system that crosslinks fibrin. When thrombin has converted fibrinogen to fibrin, the latter forms a proteinaceous network in which every E-unit is crosslinked to only one D-unit. Factor XIII is activated by thrombin into Factor XIIIa; its activation into Factor XIIIa requires calcium as a cofactor. Factor XIIIa has been identified in platelets, megakaryocytes, and fibroblast-like mesenchymal or histiocytic cells present in the placenta, uterus, and prostate; it is also present in monocytes, macrophages and dermal dendritic cells.

Anti-Factor XIIIa has been found to be useful in differentiating between Dermatofibroma (90% (+)), Dermatofibrosarcoma Protuberans (25%(+)) and Desmoplastic Malignant Melanoma (0%(+)). Factor XIIIa positivity is also seen in Capillary Hemangioblastoma (100%(+)), Hemangiopericytoma (100%(+)), Xanthogranuloma (100%(+)), Xanthoma (100%(+)), Hepatocellular Carcinoma (93%(+)), Glomus Tumor (80%(+)), and Meningioma 80%(+).

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP292
ISOTYPE: IgG
CONTROL: Placenta, Testis, Tonsil, Skin, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5505 | Tinto Predilute | 3.0 ml |
| BSB 5506 | Tinto Predilute | 7.0 ml |
| BSB 5507 | Tinto Predilute | 15.0 ml |
| BSB 5508 | Concentrate | 0.1 ml |
| BSB 5509 | Concentrate | 0.5 ml |
| BSB 5510 | Concentrate | 1.0 ml |
| BSB 5511 | Control Slides | 5 |

Fascin, MMab



IHC of Fascin on a FFPE Hodgkin's Lymphoma Tissue


Fascin, encoded by the human homolog for sn (hsn) gene, has been localized to microspikes and stress fibers of cultured cells where it is thought to be involved in the formation of microfilament bundles. It is expressed predominantly in dendritic cells. Lymphoid cells, myeloid cells and plasma cells are negative. However, Reed Sternberg cells in Hodgkin's Lymphoma are positive for Fascin staining. Epstein-Barr virus may induce expression of Fascin in B-cells.

Fascin is a very sensitive marker for Reed-Sternberg cells and variants in nodular sclerosis, mixed cellularity, and lymphocyte depletion Hodgkin's Disease. This marker might be helpful in distinguishing between Hodgkin's Disease and Non-Hodgkin Lymphoma in difficult cases. Also, the lack of expression of Fascin in the neoplastic follicles in Follicular Lymphoma can be helpful in distinguishing these lymphomas from reactive Follicular Hyperplasia in which the number of follicular dendritic cells is normal or increased.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-36
ISOTYPE: IgG1a/K
CONTROL: Hodgkin's Lymphoma, Lymph Node, Tonsil
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5512 | Tinto Predilute | 3.0 ml |
| BSB 5513 | Tinto Predilute | 7.0 ml |
| BSB 5514 | Tinto Predilute | 15.0 ml |
| BSB 5515 | Concentrate | 0.1 ml |
| BSB 5516 | Concentrate | 0.5 ml |
| BSB 5517 | Concentrate | 1.0 ml |
| BSB 5518 | Control Slides | 5 |

FGFR-3, MMab



IHC of FGFR-3 on a FFPE Transitional Cell Carcinoma Tissue


The FGFR3 gene encodes the protein Fibroblast Growth Factor Receptor 3 (FGFR3), part of a family of four fibroblast growth factor receptors that share similar structures and functions. These proteins play a role in several important cellular processes, including regulation of cell growth and proliferation, determination of cell type, angiogenesis, wound healing, and embryo development.

Approximately 80% of Non-Muscle Invasive Bladder Cancer (NMIBC) tumors have FGFR3 gene mutations. In addition to high prevalence in Bladder Cancer, somatic mutations in the FGFR3 gene have been associated with Multiple Myeloma and Cervical Cancer. FGFR3 gene aberrations have also been reported in Urothelial, Breast, Head/Neck, Lung, Brain, Gastric, Pancreatic, Colorectal, Kidney, Endometrial, Ovarian, and Cervical Cancers. FGFR3, which is highly expressed in Cutaneous Malignant Melanoma (CMM) tissues, is correlated with increased Breslow thickness and lymph node metastasis.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-150
ISOTYPE: IgG2a
CONTROL: Skin, Liver, Brain, Testis, Transitional Cell Carcinoma, Ductal Breast Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3724-3 | Tinto Predilute | 3.0 ml |
| BSB-3724-7 | Tinto Predilute | 7.0 ml |
| BSB-3724-15 | Tinto Predilute | 15.0 ml |
| BSB-3724-01 | Concentrate | 0.1 ml |
| BSB-3724-05 | Concentrate | 0.5 ml |
| BSB-3724-1 | Concentrate | 1.0 ml |
| BSB-3724-CS | Control Slides | 5 |

Fibrinogen, RPab




IHC and IF of Fibrinogen on a FFPE Kidney Tissue (IHC) and on a Frozen Lichen Planus Tissue (IF)

Fibrinogen (factor I) is a glycoprotein that circulates in the blood of vertebrates. During tissue and vascular injury, it is converted enzymatically by thrombin to fibrin and subsequently to a fibrin-based blood clot. Fibrinogen functions primarily to occlude blood vessels and thereby stop excessive bleeding. Fibrin also mediates blood platelet and endothelial cell spreading, tissue fibroblast proliferation, capillary tube formation, and angiogenesis and thereby functions to promote tissue revascularization, wound healing, and tissue repair.

Several disorders (Congenital afibrinogenemia, hypofibrinogenemia, Fibrinogen storage disease, Hereditary fibrinogen A α -Chain amyloidosis, Congenital hypodysfibrinogenemia, Cryofibrinogenemia, acquired hypofibrinogenemia, Chronic Kidney Disease, etc.) in the quantity and/or quality of fibrinogen cause pathological bleeding, pathological blood clotting, and/or the deposition of fibrinogen in the liver, kidneys, and other tissues. Chronic kidney disease (CKD) patients have increased rates of bleeding as well as thrombosis. Fibrinogen and platelets combine to generate a mature clot, but in CKD platelets are dysfunctional.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Breast, Testos, Kidney, Pancreas, Salivary Gland, Skin, Fallopian Tube
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Fli-1, MAb




IHC of Fli-1 on a FFPE PNET Tissue

Fli-1 protein, a member of the ETS family of DNA binding transcription factors, is involved in cellular proliferation and tumorigenesis. Approximately 90% of Ewing's Sarcoma/Primitive Neuroectodermal Tumors (ES/PNET) have a specific translocation, t(11;22)(q24;q12), which results in fusion of EWS to Fli-1, and production of an EWS-Fli-1 fusion protein, which can be detected by this antibody. Among normal tissues only endothelial cells and small lymphocytes express Fli-1. Fli-1 has been found to be expressed in the great majority of vascular tumors including Angiosarcomas, Hemangioendotheliomas, Hemangiomas, and Kaposi's Sarcomas.

It has been reported that the high sensitivity and specificity of Fli-1 is equal to or exceeds that of the established vascular markers, CD31, CD34, and Factor VIII. As the first nuclear marker of endothelium (rather than cytoplasmic or membranous), Fli-1 immunostaining also generally lacks cytoplasmic staining artifacts that are the result of endogenous peroxidases or biotin.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: G146-222
ISOTYPE: IgG2b
CONTROL: Adrenal Gland, Fallopian Tube, Placenta, Cervix, Hemangiomas, PNET, Angiosarcoma & Soft Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

FOXA1/HNF-3A, RPab



IHC of FOXA1 on a FFPE Breast Carcinoma Tissue

FOXA1 is a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. FOXA1 is a downstream target of GATA3 in the mammary gland.

FOXA1 in breast cancer is highly correlated with ER α +, GATA3+, and PR+ protein expression as well as endocrine signaling. FOXA1 absence in ER α + cancer might identify ER α cancers that are resistant to endocrine therapy. In ER α - breast cancer, FOXA1 is highly correlated with improved disease free survival and GATA3. Expression in ER α - cancers may identify a subset of tumors that is responsive to other endocrine therapies such as androgen receptor antagonist treatment.


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Prostate, Breast, Tonsil, Breast Carcinoma, Prostate Carcinoma, Bladder TCC
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3047 | Tinto Predilute | 3.0 ml |
| BSB 3048 | Tinto Predilute | 7.0 ml |
| BSB 3049 | Tinto Predilute | 15.0 ml |
| BSB 3050 | Concentrate | 0.1 ml |
| BSB 3051 | Concentrate | 0.5 ml |
| BSB 3052 | Concentrate | 1.0 ml |
| BSB 3053 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5519 | Tinto Predilute | 3.0 ml |
| BSB 5520 | Tinto Predilute | 7.0 ml |
| BSB 5521 | Tinto Predilute | 15.0 ml |
| BSB 5522 | Concentrate | 0.1 ml |
| BSB 5523 | Concentrate | 0.5 ml |
| BSB 5524 | Concentrate | 1.0 ml |
| BSB 5525 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6744 | Tinto Predilute | 3.0 ml |
| BSB 6745 | Tinto Predilute | 7.0 ml |
| BSB 6746 | Tinto Predilute | 15.0 ml |
| BSB 6747 | Concentrate | 0.1 ml |
| BSB 6748 | Concentrate | 0.5 ml |
| BSB 6749 | Concentrate | 1.0 ml |
| BSB 6750 | Control Slides | 5 |

FOXA1, RMAb




IHC of FOXA1 on a FFPE Colon Tissue

FOXA1 is a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. FOXA1 is a downstream target of GATA3 in the mammary gland.

FOXA1 in breast cancer is highly correlated with ERα+, GATA3+, and PR+ protein expression as well as endocrine signaling. FOXA1 absence in ERα+ cancer might identify ERα cancers that are resistant to endocrine therapy. In ERα- breast cancer, FOXA1 is highly correlated with improved disease free survival and GATA3. Expression in ERα- cancers may identify a subset of tumors that is responsive to other endocrine therapies such as androgen receptor antagonist treatment.

FOXL2, RPAb



IHC of FOXL2 on a FFPE Fallopian Tube Tissue

The Forkhead box L2 (FOXL2) gene encodes FOXL2, a transcription factor involved in ovarian development and function. Before birth and in adulthood, the FOXL2 protein regulates the growth and proliferation of hormone-producing ovarian granulosa cells. FOXL2 is also involved in the breakdown of fats, steroid hormones, and potentially harmful reactive oxygen species in the ovaries.

FOXL2 is a relatively sensitive and highly specific for Sex Cord-Stromal Tumors (SCST). FOXL2 expression is present in almost all SCSTs with a FOXL2 mutation or without a mutation. A specific somatic mutation in the FOXL2 gene has been found in Adult Granulosa Cell Tumor (402C->G / C134W), which is present in 70-95% of Ovarian Adult Granulosa Cell Tumors but not in Ovarian Fibromas. This mutation is also present in 2 of 5 men with Adult Granulosa Cell Tumor but absent in Ovarian Juvenile Granulosa Cell Tumors. Studies have demonstrated that FOXL2 is expressed in Cervical Squamous Cancer. FOXL2 suppresses the proliferation and promotes apoptosis of Cervical Cancer cells mainly through decreasing Ki-67 expression and increasing Fas ligand expression. FOXL2 has also been found to restrain the invasiveness of Cervical Cancer cells; hence, FOXL2 might be a novel tumor suppressor in Cervical Cancer. A study has demonstrated that FOXL2 is expressed in Breast Cancer and influences clinical outcome with improved recurrence-free survival in cases with nuclear expression.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Fallopian Tube, Cervix, Pancreas, Placenta, Extra Marginal Zone Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Mouse, Pig, Bovine

FOXO1, RMAb



IHC of FOXO1 on a FFPE Lymphoblastic Lymphoma Tissue

The Forkhead transcription factor FOXO1, an important downstream target of phosphatidylinositol-3-kinase (PI3K)/AKT signaling pathway, regulates cellular homeostasis by maintaining cell proliferation, apoptosis and viability in normal cells. FOXO1 is a transcription factor that plays important roles in regulation of gluconeogenesis and glycogenolysis by insulin signaling, and is also central to the decision for a preadipocyte to commit to adipogenesis.

FOXO1 is broadly expressed in different types of cells with high level of expression in lymphoid cells and non-Hodgkin's lymphomas. In contrast, in most of classical Hodgkin lymphoma (cHL), Reed-Sternberg cells are FOXO1 negative. Androgens and the androgen receptor (AR) are essential for growth and differentiation of the normal prostate gland as well as proliferation and survival of Prostate Cancer. FOXA1, functions as a pioneer factor to facilitate AR transactivation and Prostate Cancer growth. In contrast, the O-class of FOX proteins such as FOXO1 and FOXO3, which are downstream effectors of the PTEN tumor suppressor, inhibits the transcriptional activity of either full-length AR or constitutively active splice variants of AR in a direct or indirect manner in Prostate Cancer. Translocation of FOXO1 gene with PAX3 has been associated with alveolar rhabdomyosarcoma.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP290

ISOTYPE: IgG

CONTROL: Testis, Thyroid, Tonsil, Lymph Node, Spleen, Lung, Lymphomas

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP277

ISOTYPE: IgG


CONTROL: Prostate, Breast, Tonsil, Colon, Fallopian Tube, Breast Carcinoma, Prostate Carcinoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|-------------|-----------------|---------|----------|-----------------|---------|
| BSB 2908 | Tinto Predilute | 3.0 ml | BSB-3725-3 | Tinto Predilute | 3.0 ml | BSB 2915 | Tinto Predilute | 3.0 ml |
| BSB 2909 | Tinto Predilute | 7.0 ml | BSB-3725-7 | Tinto Predilute | 7.0 ml | BSB 2916 | Tinto Predilute | 7.0 ml |
| BSB 2910 | Tinto Predilute | 15.0 ml | BSB-3725-15 | Tinto Predilute | 15.0 ml | BSB 2917 | Tinto Predilute | 15.0 ml |
| BSB 2911 | Concentrate | 0.1 ml | BSB-3725-01 | Concentrate | 0.1 ml | BSB 2918 | Concentrate | 0.1 ml |
| BSB 2912 | Concentrate | 0.5 ml | BSB-3725-05 | Concentrate | 0.5 ml | BSB 2919 | Concentrate | 0.5 ml |
| BSB 2913 | Concentrate | 1.0 ml | BSB-3725-1 | Concentrate | 1.0 ml | BSB 2920 | Concentrate | 1.0 ml |
| BSB 2914 | Control Slides | 5 | BSB-3725-CS | Control Slides | 5 | BSB 2921 | Control Slides | 5 |

FOXP1, RMAb




IHC of FOXP1 on a FFPE Breast Tissue

FOXP3, RPAb



IHC of FOXP3 on a FFPE Colon Tissue

FOXP3, RMAb



IHC of FOXP3 on a FFPE Thymus Tissue

FOXP1 is part of the forkhead box (FOX) transcription factor family. Forkhead box transcription factors play important roles in regulation of tissue- and cell-type specific gene transcription during both development and adulthood. The FOXP1 protein contains both DNA-binding- and protein-protein binding-domains. FOXP1 is a transcriptional repressor and is responsible for regulating a variety of important aspects of development including tissue development of the lungs, brain, thymus and heart. It is also important in muscle development of the esophagus and esophageal epithelium and for regulating lung airway morphogenesis.

Strong expression of FOXP1 is associated with poor disease-free survival and transformation to Diffuse Large B-cell Lymphomas. Recently, studies suggested a role of FOXP1 in the regulation of ER expression. FOXP1 expression is correlated with ER expression and improved survival in breast cancer patients. Nuclear expression of FOXP1 is associated with ER expression, while cytoplasmic expression of FOXP1 is linked to myometrial invasion in endometrial cancer.

FOXP3, also known as scurf, is a protein involved in immune system responses. A member of the forkhead box protein family, FOXP3 appears to function as a transcription factor in the development and function of regulatory T cells. In regulatory T cell model systems, the FOXP3 transcription factor occupies the promoters of many important for regulatory T-cell function, and may repress transcription of key genes following stimulation of T cell receptors.

Alterations in numbers of regulatory T-cells — in particular those that express FOXP3 — are found in a number of disease states. Patients with tumors have a local relative excess of FOXP3 positive T cells which inhibits the body's ability to suppress the formation of cancerous cells.

FOXP3, also known as scurf, is a protein involved in immune system responses. A member of the forkhead box protein family, FOXP3 appears to function as a transcription factor in the development and function of regulatory T cells. In regulatory T cell model systems, the FOXP3 transcription factor occupies the promoters of many important for regulatory T-cell function, and may repress transcription of key genes following stimulation of T cell receptors.


Alterations in numbers of regulatory T-cells - in particular those that express FOXP3 - are found in a number of disease states. Patients with tumors have a local relative excess of FOXP3 positive T cells which inhibits the body's ability to suppress the formation of cancerous cells.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP137
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Breast
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Colon
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP340
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Thymus
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6751 | Tinto Predilute | 3.0 ml | BSB 6758 | Tinto Predilute | 3.0 ml | BSB 2922 | Tinto Predilute | 3.0 ml |
| BSB 6752 | Tinto Predilute | 7.0 ml | BSB 6759 | Tinto Predilute | 7.0 ml | BSB 2923 | Tinto Predilute | 7.0 ml |
| BSB 6753 | Tinto Predilute | 15.0 ml | BSB 6760 | Tinto Predilute | 15.0 ml | BSB 2924 | Tinto Predilute | 15.0 ml |
| BSB 6754 | Concentrate | 0.1 ml | BSB 6761 | Concentrate | 0.1 ml | BSB 2925 | Concentrate | 0.1 ml |
| BSB 6755 | Concentrate | 0.5 ml | BSB 6762 | Concentrate | 0.5 ml | BSB 2926 | Concentrate | 0.5 ml |
| BSB 6756 | Concentrate | 1.0 ml | BSB 6763 | Concentrate | 1.0 ml | BSB 2927 | Concentrate | 1.0 ml |
| BSB 6757 | Control Slides | 5 | BSB 6764 | Control Slides | 5 | BSB 2928 | Control Slides | 5 |

FSH, MMab


IHC of FSH on a FFPE Pituitary Tissue

Follicle stimulating hormone (FSH) is a hormone synthesized and secreted by gonadotropes in the anterior pituitary gland. In the ovary, FSH stimulates the growth of immature Graafian follicles to maturation. As the follicle grows, it releases inhibin, which deactivates the FSH production. In men, FSH enhances the production of androgen-binding protein by the Sertoli cells of the testis and is critical for spermatogenesis. FSH and LH act synergistically in reproduction.

FSH is a useful marker in the classification of pituitary tumors and the study of pituitary disease. It reacts with FSH-producing cells.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-55
ISOTYPE: IgG1/K
CONTROL: Normal Pituitary
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5533 | Tinto Predilute | 3.0 ml |
| BSB 5534 | Tinto Predilute | 7.0 ml |
| BSB 5535 | Tinto Predilute | 15.0 ml |
| BSB 5536 | Concentrate | 0.1 ml |
| BSB 5537 | Concentrate | 0.5 ml |
| BSB 5538 | Concentrate | 1.0 ml |
| BSB 5539 | Control Slides | 5 |

FSH, RMAb


IHC of FSH on a FFPE Pituitary Adenoma Tissue

Follicle stimulating hormone (FSH) is a hormone synthesized and secreted by gonadotropes in the anterior pituitary gland. In the ovary, FSH stimulates the growth of immature Graaf follicles to maturation. As the follicle grows, it releases inhibin, which deactivates the FSH production. In men, FSH enhances the production of androgen-binding protein by the Sertoli cells of the testis and is critical for spermatogenesis. FSH and LH act synergistically in reproduction.

FSH is a useful marker in the classification of pituitary tumors and the study of pituitary disease. It reacts with FSH-producing cells.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP257
ISOTYPE: IgG
CONTROL: Pituitary
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3539 | Tinto Predilute | 3.0 ml |
| BSB 3540 | Tinto Predilute | 7.0 ml |
| BSB 3541 | Tinto Predilute | 15.0 ml |
| BSB 3542 | Concentrate | 0.1 ml |
| BSB 3543 | Concentrate | 0.5 ml |
| BSB 3544 | Concentrate | 1.0 ml |
| BSB 3545 | Control Slides | 5 |

Fumarate Hydratase/Fumarase, MMab

IHC of Fumarate Hydratase on a FFPE Kidney Tissue


Fumarate hydratase (FH), or fumarase, is encoded by the FH gene, which is an enzymatic component of the tricarboxylic acid (TCA) cycle, or Krebs cycle, where it catalyzes the formation of L-malate from fumarate. Mutations in the FH gene can cause FH deficiency and lead to progressive encephalopathy. It was discovered that succinate dehydrogenase and FH are tumour suppressors and they are associated with metabolic dysfunction and tumorigenesis, providing biochemical evidence to explain enhanced glycolysis in tumours.

Hereditary Leiomyomatosis and Renal Cell Cancer (HLRCC) is an autosomal dominant heritable syndrome with predisposition to development of Renal Cell Carcinoma and Smooth Muscle Tumors of the skin and uterus. Cells of individuals with HLRCC had lower FH enzyme activity than cells from normal controls, making FH enzyme activity testing a useful method for diagnosis and screening. Loss-of-function mutations of FH predisposes individuals to the autosomal dominant syndrome of Multiple Cutaneous and Uterine Leiomyomatosis (MCUL). Biallelic inactivation/mutations of FH are seen in 85% of Hereditary Leiomyomatosis and Renal Cell Carcinoma cases, 100% of Renal Cell Carcinoma with germline FH mutations, 19% of Papillary Renal Cell Carcinoma (Type II) have FH deficiency and 90% of FH deficiency RCC have FH mutations, 1% of unselected Leiomyomas, 2.6% of Leiomyomas in patients < 40 years old and 37- 52% of Leiomyoma with bizarre nuclei.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-151
ISOTYPE: IgG1
CONTROL: Placenta, Breast, Fallopian Tube, Colon, Kidney, Testis, Colon Adenocarcinoma, HER2 Negative Breast Cancer
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3726-3 | Tinto Predilute | 3.0 ml |
| BSB-3726-7 | Tinto Predilute | 7.0 ml |
| BSB-3726-15 | Tinto Predilute | 15.0 ml |
| BSB-3726-01 | Concentrate | 0.1 ml |
| BSB-3726-05 | Concentrate | 0.5 ml |
| BSB-3726-1 | Concentrate | 1.0 ml |
| BSB-3726-CS | Control Slides | 5 |

GAB1/GRB2-Associated-Binding Protein 1, MAb




IHC of GAB1/GRB2 - Associated-Binding Protein 1 on a FFPE Testis Tissue

Galectin-3, MAb



IHC of Galectin-3 on a FFPE Kidney Tissue

Gastrin, RPab



IHC of Gastrin on a FFPE Stomach Tissue

GRB2-associated-binding protein 1 (GAB1) is encoded by the GAB1 gene. GAB1 is a member of the IRS1-like multisubstrate docking protein family that transduces signals from various tyrosine kinases, such as Met, FGFR1, and EGFR. GAB1 plays a central role in cellular growth response, transformation, and apoptosis, as well as inflammatory responses. GAB1 is involved in the amplification of IL-6-induced MAPK pathway and promotes inflammation.

Somatic mutations of the GAB1 gene have been detected in Breast and Colorectal Cancers and studies have shown that elevated expressions of GAB1 have been associated with Breast Cancer metastasis by dissociating the polarity-associated partitioning defective (PAR) complex and promoting epithelial-to-mesenchymal transition. Upregulation of GAB1 is an indication of unfavorable prognosis for Hepatocellular Carcinoma and Epithelial Ovarian Cancer. GAB1 overexpression has also been seen in Adult Acute Lymphoblastic Leukemia and Medulloblastomas. GAB1 can be used with a panel of immunohistochemical markers in the classification of Medulloblastomas into SHH (sonic hedgehog), WNT (wingless-type murine mammary tumor), or non-SHH/WNT subgroups.

Galectin-3 is a 31 kDa beta-galactosidase binding lectin. It has been associated with binding to the basement membrane glycoprotein laminin. Galectin-3 is normally distributed in epithelia of many organs and various inflammatory cells, including macrophages, as well as dendritic cells and Kupffer cells. The expression of this lectin is up-regulated during inflammation, cell proliferation, cell differentiation and through trans-activation by viral proteins.

Anti-Galectin-3 has been demonstrated to be valuable in differentiating between benign and malignant thyroid neoplasms in both histologic sections and Fine Needle Aspiration Biopsy material. Anti-Galectin-3 antibody has also been useful in identifying Anaplastic Large Cell Lymphoma.

Gastrin is a linear peptide hormone produced by G-cells of the duodenum and in the pyloric antrum of the stomach. It is secreted into the bloodstream.

Gastrin antibody gives positive staining of G-cells of human antral/pyloric mucosa and cells producing gastrin or a structural gastrin analogue as is seen in the stomach. No staining of other cells or tissue types has been observed. This antibody may react with sulfated and non-sulfated forms of gastrin. The antibody cross-reacts with more than 50% of the present cholecystokinin octapeptide.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-155

ISOTYPE: IgG2a

CONTROL: Breast, Prostate, Testis, Tonsil, Stomach, Transitional Cell Carcinoma

LOCALIZATION: Cytoplasmic, Membranous, Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 9C4

ISOTYPE: IgG1

CONTROL: Kidney, Testis, Salivary Gland, Breast, Tonsil, Colon, Papillary & Follicular Carcinoma of Thyroid

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Stomach


LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB-3727-3 | Tinto Predilute | 3.0 ml | BSB 5540 | Tinto Predilute | 3.0 ml | BSB 5547 | Tinto Predilute | 3.0 ml |
| BSB-3727-7 | Tinto Predilute | 7.0 ml | BSB 5541 | Tinto Predilute | 7.0 ml | BSB 5548 | Tinto Predilute | 7.0 ml |
| BSB-3727-15 | Tinto Predilute | 15.0 ml | BSB 5542 | Tinto Predilute | 15.0 ml | BSB 5549 | Tinto Predilute | 15.0 ml |
| BSB-3727-01 | Concentrate | 0.1 ml | BSB 5543 | Concentrate | 0.1 ml | BSB 5550 | Concentrate | 0.1 ml |
| BSB-3727-05 | Concentrate | 0.5 ml | BSB 5544 | Concentrate | 0.5 ml | BSB 5551 | Concentrate | 0.5 ml |
| BSB-3727-1 | Concentrate | 1.0 ml | BSB 5545 | Concentrate | 1.0 ml | BSB 5552 | Concentrate | 1.0 ml |
| BSB-3727-CS | Control Slides | 5 | BSB 5546 | Control Slides | 5 | BSB 5553 | Control Slides | 5 |

GATA3, MAb

IHC of GATA-3 on a FFPE Breast Carcinoma Tissue

GATA3, RMab

IHC of GATA-3 on a FFPE Breast Carcinoma Tissue

GCDFP-15, MAb

IHC of GCDFP-15 on a FFPE Breast Carcinoma Tissue

Trans-acting T-cell-specific transcription factor, GATA-3 is, a protein that in humans is encoded by the GATA3 gene. GATA-3 b regulates luminal epithelial cell differentiation in the mammary gland, is an important regulator of T cell development and plays an important role in endothelial cell biology.

GATA-3 is one of the three genes mutated in >10% of breast cancers. Nuclear expression of GATA-3 in breast cancer is considered a marker of luminal cancer in ER+ cancer and luminal androgen responsive cancer in ER-/AR+ tumors. It is highly coexpressed with FOXA1 and serves as negative predictor of basal subtype and HER-2 and is also considered a strong predictor of taxane and platin salts insensitivity.

GATA3 expression is found in urothelial carcinoma, especially in invasive and high grade tumors. Therefore, anti-GATA3 can be used in a panel of antibodies for diagnosis of unknown primary carcinoma, when carcinomas of the breast or bladder are a possibility. Studies have also shown the utility of GATA-3 in differentiating urothelial carcinoma from prostate adenocarcinoma and squamous cell carcinomas of the uterine, cervix, anus and lung.

Trans-acting T-cell-specific transcription factor, GATA-3 is, a protein that in humans is encoded by the GATA3 gene. GATA-3 b regulates luminal epithelial cell differentiation in the mammary gland, is an important regulator of T cell development and plays an important role in endothelial cell biology.

GATA-3 is one of the three genes mutated in >10% of breast cancers. Nuclear expression of GATA-3 in breast cancer is considered a marker of luminal cancer in ER+ cancer and luminal androgen responsive cancer in ER-/AR+ tumors. It is highly coexpressed with FOXA1 and serves as negative predictor of basal subtype and HER-2 and is also considered a strong predictor of taxane and platin salts insensitivity.

GATA3 expression is found in urothelial carcinoma, especially in invasive and high grade tumors. Therefore, anti-GATA3 can be used in a panel of antibodies for diagnosis of unknown primary carcinoma, when carcinomas of the breast or bladder are a possibility. Studies have also shown the utility of GATA-3 in differentiating urothelial carcinoma from prostate adenocarcinoma and squamous cell carcinomas of the uterine, cervix, anus and lung.

Gross Cystic Disease is a common premenopausal disorder in which gross cysts are the predominant pathologic lesion. It is characterized by production of a fluid secretion which accumulates in the breast cysts. Gross Cystic Disease fluid is a pathologic secretion from breast composed of several glycoproteins, including a unique 15 kDa monomer protein, GCDFP-15. The cells within the body that produce GCDFP-15 appear to be restricted primarily to those with apocrine function such as breast cysts and in apocrine glands in the axilla, vulva, eyelid, and ear canal.

Studies have found GCDFP-15 to be a highly specific and sensitive marker for breast cancer. Approximately 70% of breast carcinomas stain positive with antibody to GCDFP-15. In contrast, Colorectal Carcinomas, as well as Mesotheliomas, do not stain with this antibody. Lung Adenocarcinomas rarely stain with this antibody.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: L50-823

ISOTYPE: IgG1/K

CONTROL: Breast Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP368

ISOTYPE: IgG

CONTROL: Breast, Skin, Cervix, Breast Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 23A3

ISOTYPE: IgG2a

CONTROL: Breast, Salivary Gland, Sweat Glands in Skin, Breast Carcinoma

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2670 | Tinto Predilute | 3.0 ml |
| BSB 2671 | Tinto Predilute | 7.0 ml |
| BSB 2672 | Tinto Predilute | 15.0 ml |
| BSB 2673 | Concentrate | 0.1 ml |
| BSB 2674 | Concentrate | 0.5 ml |
| BSB 2675 | Concentrate | 1.0 ml |
| BSB 2676 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3328 | Tinto Predilute | 3.0 ml |
| BSB 3329 | Tinto Predilute | 7.0 ml |
| BSB 3330 | Tinto Predilute | 15.0 ml |
| BSB 3331 | Concentrate | 0.1 ml |
| BSB 3332 | Concentrate | 0.5 ml |
| BSB 3333 | Concentrate | 1.0 ml |
| BSB 3334 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5554 | Tinto Predilute | 3.0 ml |
| BSB 5555 | Tinto Predilute | 7.0 ml |
| BSB 5556 | Tinto Predilute | 15.0 ml |
| BSB 5557 | Concentrate | 0.1 ml |
| BSB 5558 | Concentrate | 0.5 ml |
| BSB 5559 | Concentrate | 1.0 ml |
| BSB 5560 | Control Slides | 5 |

GFAP, MMab




IHC of GFAP on an FFPE on a Brain Tissue

GFAP, RMAb



IHC of GFAP on FFPE Brain Tissue

GH, RPAb



IHC of GH on an FFPE on a Pituitary Tissue

Glial fibrillary acidic protein or GFAP is a Type III protein of the intermediate filaments principally found in astrocytes in the central nervous system, but can also be found in neurons, hepatic stellate cells, kidney mesangial cells, pancreatic stellate cells, and Leydig cells. It has a role in the cytoskeleton of the astrocyte and possibly many other stellate-shaped cells.

Antibodies to GFAP are very useful as markers of astrocytic cells. In addition, many types of brain tumors, presumably derived from astrocytic cells, heavily express GFAP. This marker is mainly used to distinguish neoplasms of astrocytic origin from other neoplasms in the central nervous system.

Glial fibrillary acidic protein or GFAP is a Type III protein of the intermediate filaments principally found in astrocytes in the central nervous system, but can also be found in neurons, hepatic stellate cells, kidney mesangial cells, pancreatic stellate cells, and Leydig cells. It has a role in the cytoskeleton of the astrocyte and possibly many other stellate-shaped cells.

Antibodies to GFAP are very useful as markers of astrocytic cells. In addition, many types of brain tumors, presumably derived from astrocytic cells, heavily express GFAP. This marker is mainly used to distinguish neoplasms of astrocytic origin from other neoplasms in the central nervous system.

Growth hormone (GH or somatotropin) is a 191 amino acid, single-chain polypeptide hormone which is synthesized, stored and secreted by the somatotroph cells within the lateral wings of the anterior pituitary gland, which stimulates growth and cell reproduction in humans and other animals.

GH is a useful marker in classification of pituitary tumors and the study of pituitary disease (acromegaly). It reacts with GH-producing cells.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: G-A-5

ISOTYPE: IgG1

CONTROL: Brain

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat, Rat, Mouse, Rabbit

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM246

ISOTYPE: IgG

CONTROL: Brain

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Predicted: Rat

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Normal Pituitary


LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5561 | Tinto Predilute | 3.0 ml |
| BSB 5562 | Tinto Predilute | 7.0 ml |
| BSB 5563 | Tinto Predilute | 15.0 ml |
| BSB 5564 | Concentrate | 0.1 ml |
| BSB 5565 | Concentrate | 0.5 ml |
| BSB 5566 | Concentrate | 1.0 ml |
| BSB 5567 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3766-3 | Tinto Predilute | 3.0 ml |
| BSB-3766-7 | Tinto Predilute | 7.0 ml |
| BSB-3766-15 | Tinto Predilute | 15.0 ml |
| BSB-3766-01 | Concentrate | 0.1 ml |
| BSB-3766-05 | Concentrate | 0.5 ml |
| BSB-3766-1 | Concentrate | 1.0 ml |
| BSB-3766-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5568 | Tinto Predilute | 3.0 ml |
| BSB 5569 | Tinto Predilute | 7.0 ml |
| BSB 5570 | Tinto Predilute | 15.0 ml |
| BSB 5571 | Concentrate | 0.1 ml |
| BSB 5572 | Concentrate | 0.5 ml |
| BSB 5573 | Concentrate | 1.0 ml |
| BSB 5574 | Control Slides | 5 |

GH, MMab

IHC of Growth Hormone on an FFPE on a Pituitary Tissue

Growth hormone (GH or somatotropin) is a 191 amino acid, single-chain polypeptide hormone which is synthesized, stored and secreted by the somatotroph cells within the lateral wings of the anterior pituitary gland, which stimulates growth and cell reproduction in humans and other animals.


GH is a useful marker in classification of pituitary tumors and the study of pituitary disease (acromegaly). It reacts with GH-producing cells.

GH, RMab

IHC of Growth Hormone on a FFPE on a Pituitary Tissue

Growth hormone (also known as GH or somatotropin) is a 191 amino acid, single-chain polypeptide hormone which is synthesized, stored and secreted by the somatotroph cells within the lateral wings of the anterior pituitary gland, which stimulates growth and cell reproduction in humans and other animals.

GH is a useful marker in classification of pituitary tumors and the study of pituitary disease (acromegaly). It reacts with GH-producing cells.

Glucagon, MMab

IHC of Glucagon on a FFPE Pancreas Tissue

Glucagon is a 29-amino acid polypeptide acting as an important hormone in carbohydrate metabolism. The hormone is synthesized and secreted from alpha cells of the islets of Langerhans, which are located in the endocrine portion of the pancreas. Abnormally-elevated levels of glucagon may be caused by pancreatic tumors such as glucagonoma, symptoms of which include necrolytic migratory erythema (NME), elevated amino acids and hyperglycemia. It may occur alone or in the context of Multiple Endocrine Neoplasia Type 1.

Glucagon antibody detects glucagon-secreting cells and tumors such as glucagonomas. Studies show that approximately 80% of glucagonomas are malignant and these patients have a syndrome most often initially recognized by dermatologists. Symptoms include necrolytic migratory erythema as well as diabetes, anemia, stomatitis, weight loss, frequent venous thromboses, and in some instances, diarrhea and psychiatric disturbances. The diagnosis may be readily confirmed by the demonstration of elevated plasma glucagon concentration.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-99

ISOTYPE: IgG1/K

CONTROL: Normal Pituitary

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP267

ISOTYPE: IgG

CONTROL: Normal Pituitary

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-111

ISOTYPE: IgG1/K

CONTROL: Pancreas, Colon

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3335 | Tinto Predilute | 3.0 ml |
| BSB 3336 | Tinto Predilute | 7.0 ml |
| BSB 3337 | Tinto Predilute | 15.0 ml |
| BSB 3338 | Concentrate | 0.1 ml |
| BSB 3339 | Concentrate | 0.5 ml |
| BSB 3340 | Concentrate | 1.0 ml |
| BSB 3341 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2496 | Tinto Predilute | 3.0 ml |
| BSB 2497 | Tinto Predilute | 7.0 ml |
| BSB 2498 | Tinto Predilute | 15.0 ml |
| BSB 2499 | Concentrate | 0.1 ml |
| BSB 2500 | Concentrate | 0.5 ml |
| BSB 2501 | Concentrate | 1.0 ml |
| BSB 2502 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3342 | Tinto Predilute | 3.0 ml |
| BSB 3343 | Tinto Predilute | 7.0 ml |
| BSB 3344 | Tinto Predilute | 15.0 ml |
| BSB 3345 | Concentrate | 0.1 ml |
| BSB 3346 | Concentrate | 0.5 ml |
| BSB 3347 | Concentrate | 1.0 ml |
| BSB 3348 | Control Slides | 5 |

Glucagon, RMab




IHC of Glucagon on a FFPE Pancreas Tissue

GLUT1, RMab



IHC of GLUT1 on a FFPE Cervix Tissue

Glutamine Synthetase, MAb



IHC of GS on a FFPE Hepatocellular Carcinoma Tissue

Glucagon is a 29-amino acid polypeptide acting as an important hormone in carbohydrate metabolism. The hormone is synthesized and secreted from alpha cells of the islets of Langerhans, which are located in the endocrine portion of the pancreas. Abnormally-elevated levels of glucagon may be caused by pancreatic tumors such as glucagonoma, symptoms of which include necrolytic migratory erythema (NME), elevated amino acids and hyperglycemia. It may occur alone or in the context of Multiple Endocrine Neoplasia Type 1.

Glucagon antibody detects glucagon-secreting cells and tumors such as glucagonomas. Studies show that approximately 80% of glucagonomas are malignant and these patients have a syndrome most often initially recognized by dermatologists. Symptoms include necrolytic migratory erythema as well as diabetes, anemia, stomatitis, weight loss, frequent venous thrombosis, and in some instances, diarrhea and psychiatric disturbances. The diagnosis may be readily confirmed by the demonstration of elevated plasma glucagon concentration.

GLUT1 facilitates the transport of glucose across the plasma membranes of mammalian cells. Energy-yielding metabolism in erythrocytes depends on a constant supply of glucose from blood plasma. Glucose enters the erythrocyte by facilitated diffusion via a specific glucose transporter, at a rate about 50,000 times greater than uncatalyzed transmembrane diffusion. GLUT1 is a type III integral protein with 12 hydrophobic segments, each of which is believed to form a membrane-spanning helix. It is responsible for the low-level of basal glucose uptake required to sustain respiration in all cells. GLUT1 is also a major receptor for the uptake of Vitamin C as well as glucose.

GLUT1 is expressed at variable levels in many human tissues. Overexpression of GLUT1 has been linked to tumor progression or poor survival of patients with carcinomas of the colon, breast, cervix, lung, bladder and mesothelioma. It can be used to distinguish between malignant mesothelioma (positive) from reactive mesothelium (negative).

Glutamine synthetase (GS) is an enzyme that plays an essential role in the metabolism of nitrogen by catalyzing the condensation of glutamate and ammonia to form glutamine. GS is present predominantly in the brain, kidneys, and liver. GS in the brain participates in the metabolic regulation of glutamate, the detoxification of brain ammonia, the assimilation of ammonia, recycling of neurotransmitters, and termination of neurotransmitter signals. In normal liver, GS expression is seen in pericentral hepatocytes, but not in mid-zonal or periportal hepatocytes.

GS positive tumor cells are believed to be derived from GS positive hepatocytes. GS immunoreactivity has been seen in a majority of hepatocellular carcinoma (HCC), including cases of early HCC (70%) and for low grade HCC (59%). A panel composed of antibodies against HSP70, GPC3, and GS has been proposed to be very useful in distinguishing between dysplastic and early malignant hepatocellular nodules arising in cirrhosis. Staining of hepatocellular lesions with anti-GS antibody have been useful in the differential diagnosis of focal nodular hyperplasia (FNH), hepatic adenoma (HCA), and dysplastic nodules, and low grade hepatocellular carcinoma. In the case of FNH, GS stains in a characteristic "map-like" pattern, thus differentiating it from HCA, in which GS staining is usually absent, but may occasionally be present at the border of the lesion or around the veins inside the tumor.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP74
ISOTYPE: IgG
CONTROL: Pancreas, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP141
ISOTYPE: IgG
CONTROL: Placenta, Colon, Prostate, Skin, Kidney, Brain, Tonsil, Breast, Mesothelioma, Colon Carcinoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human, Mouse, Rat


ANTIBODY TYPE: Mouse Monoclonal
CLONE: GS-6
ISOTYPE: IgG2a
CONTROL: Liver, Tonsil, Testis, Prostate, Hepatocellular Carcinoma, Bladder TCC
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2328 | Tinto Predilute | 3.0 ml |
| BSB 2329 | Tinto Predilute | 7.0 ml |
| BSB 2330 | Tinto Predilute | 15.0 ml |
| BSB 2331 | Concentrate | 0.1 ml |
| BSB 2332 | Concentrate | 0.5 ml |
| BSB 2333 | Concentrate | 1.0 ml |
| BSB 2334 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6779 | Tinto Predilute | 3.0 ml |
| BSB 6780 | Tinto Predilute | 7.0 ml |
| BSB 6781 | Tinto Predilute | 15.0 ml |
| BSB 6782 | Concentrate | 0.1 ml |
| BSB 6783 | Concentrate | 0.5 ml |
| BSB 6784 | Concentrate | 1.0 ml |
| BSB 6785 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2929 | Tinto Predilute | 3.0 ml |
| BSB 2930 | Tinto Predilute | 7.0 ml |
| BSB 2931 | Tinto Predilute | 15.0 ml |
| BSB 2932 | Concentrate | 0.1 ml |
| BSB 2933 | Concentrate | 0.5 ml |
| BSB 2934 | Concentrate | 1.0 ml |
| BSB 2935 | Control Slides | 5 |

Glycophorin A, MMab




IHC of Glycophorin A on a FFPE Bone Marrow Tissue

Glycophorins A (GPA) and B (GPB) are single pass membrane sialoglycoproteins. GPA is the carrier of blood group M and N specificities, while GPB accounts for S and U specificities. GPA and GPB provide the cells with a large mucin-like surface and it has been suggested this provides a barrier to cell fusion, thus minimizing aggregation between red blood cells in the circulation.

Anti-Glycophorin A has been used to characterize erythroid cell development and in the diagnosis of Erythroid Leukemias.

Glypican-3, MMab




IHC of Glypican-3 on a FFPE Ovarian Carcinoma Tissue

Glypican 3, also known as GPC3, is a human gene. The protein encoded by this gene is a member of the glypican family. Cell surface heparan sulfate proteoglycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypican-related integral membrane proteoglycan family (GRIPs) contain a core protein anchored to the cytoplasmic membrane via a glycosyl-phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation.

GPC3 has been identified to be a useful tumor marker for the diagnosis of Hepatocellular Carcinoma (HCC), Hepatoblastoma, Melanoma, Testicular Germ Cell Tumors, and Wilms Tumor. In patients with HCC, GPC3 was overexpressed in neoplastic liver tissue and elevated in serum but was undetectable in normal liver, benign liver, and the serum of healthy donors. GPC3 expression was also found to be higher in HCC liver tissue than in cirrhotic liver or liver with focal lesions such as dysplastic nodules and areas of hepatic adenoma (HA) with malignant transformation. In the context of Testicular Germ Cell Tumors, GPC3 expression is up-regulated in certain histologic subtypes, specifically Yolk Sac Tumors and Choriocarcinoma. A high level of GPC3 expression has also been found in some types of embryonal tumors, such as Wilms Tumor and Hepatoblastoma.

Granzyme B, RPA



IHC of Granzyme B on a FFPE Tonsil Tissue

Granzymes are exogenous serine proteases that are released by cytoplasmic granules within cytotoxic T-cells and natural killer cells. Their purpose is to induce apoptosis within virus-infected cells, thus destroying them.

Anti-Granzyme B antibodies have been useful in diagnosing Natural Killer/T-cell Lymphoma, as well as Anaplastic Large Cell Lymphoma. High percentages of cytotoxic T-cells have been shown to be an unfavorable prognostic indicator in Hodgkin's Disease.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: GA-R2

ISOTYPE: IgG2b/K

CONTROL: TBone Marrow, Placenta, Tonsil, Liver, Prostate, Adrenal, Spleen, Colon, Pancreas, Fallopian Tube

LOCALIZATION: Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 1G12

ISOTYPE: IgG1

CONTROL: Bone Marrow, Placenta, Tonsil, Liver, Prostate, Adrenal, Spleen, Colon, Pancreas, Fallopian

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Liver, Testis, Cervix, Tonsil, Lymph Node, Spleen

LOCALIZATION: Cytoplasmic (Granular)


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5582 | Tinto Predilute | 3.0 ml |
| BSB 5583 | Tinto Predilute | 7.0 ml |
| BSB 5584 | Tinto Predilute | 15.0 ml |
| BSB 5585 | Concentrate | 0.1 ml |
| BSB 5586 | Concentrate | 0.5 ml |
| BSB 5587 | Concentrate | 1.0 ml |
| BSB 5588 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6240 | Tinto Predilute | 3.0 ml |
| BSB 6241 | Tinto Predilute | 7.0 ml |
| BSB 6242 | Tinto Predilute | 15.0 ml |
| BSB 6243 | Concentrate | 0.1 ml |
| BSB 6244 | Concentrate | 0.5 ml |
| BSB 6245 | Concentrate | 1.0 ml |
| BSB 6246 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5589 | Tinto Predilute | 3.0 ml |
| BSB 5590 | Tinto Predilute | 7.0 ml |
| BSB 5591 | Tinto Predilute | 15.0 ml |
| BSB 5592 | Concentrate | 0.1 ml |
| BSB 5593 | Concentrate | 0.5 ml |
| BSB 5594 | Concentrate | 1.0 ml |
| BSB 5595 | Control Slides | 5 |

Granzyme B, RMAb




IHC of Granzyme B on a FFPE Spleen Tissue

Granzymes are exogenous serine proteases that are released by cytoplasmic granules within cytotoxic T-cells and natural killer cells. Their purpose is to induce apoptosis within virus-infected cells, thus destroying them.

Anti-Granzyme B antibodies have been useful in diagnosing Natural Killer/T-cell Lymphoma, as well as Anaplastic Large Cell Lymphoma. High percentages of cytotoxic T-cells have been shown to be an unfavorable prognostic indicator in Hodgkin's Disease.

hCG, MAb




IHC of hCG on a FFPE Placenta Tissue

Human chorionic gonadotropin (hCG) is a peptide hormone produced in pregnancy, made by the embryo soon after conception and later by the syncytiotrophoblast. Its role is to prevent the disintegration of the corpus luteum of the ovary and thereby maintain progesterone production that is critical for a pregnancy in humans. hCG may have additional functions; for instance, it is thought to affect the immune tolerance of the pregnancy. Early pregnancy testing generally is based on the detection or measurement of hCG.

hCG antibody detects cells and tumors of trophoblastic origin such as Choriocarcinomas. Large Cell Carcinoma and Adenocarcinoma of the Lung demonstrate hCG positivity in 90% and 60% of cases respectively. 20% of Squamous Cell Lung Carcinomas are positive for hCG. hCG expression by non-trophoblastic tumors may indicate aggressive behavior since it has been observed that hCG may play a role in the host response to a given tumor.

HE4, RMAb



IHC of HE4 on a FFPE Ovarian Cancer Tissue

Human epididymis protein 4 (HE4), is a protein that in humans is encoded by the WFDC2 gene. This gene encodes a protein that is a member of the WFDC domain family. The WFDC domain, or WAP Signature motif, functions as a protease inhibitor in many family members. The encoded protein is a small secretory protein, which may be involved in sperm maturation. HE4 gene expression has been found the highest in normal human trachea and salivary gland, and to a lesser extent, lung, prostate, pituitary gland, thyroid, and kidney. HE4 immunoreactivity has also been found in normal glandular epithelium of the female genital tract and breast, the epididymis and vas deferens, respiratory epithelium, distal renal tubules, colonic mucosa, and salivary glands which is consistent with the HE4 gene expression.

HE4 is a recognized biomarker in ovarian and endometrial cancer and over-expressed in pancreatic adenocarcinoma. In a series of 175 human adult tumors, gene expression was found to be the highest in ovarian serous carcinomas. However, adenocarcinomas of the lung, and occasional breast, transitional cell and pancreatic carcinomas had moderate or high levels of HE4 expression. IHC studies have shown that HE4 is significantly higher expressed in human pancreatic carcinoma tissues than in both normal and adjacent non-tumoral pancreatic tissues, and the staining intensity is inversely correlated with the clinical stage. HE4 is also highly expressed in early stage pancreatic adenocarcinoma.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP230

ISOTYPE: IgG

CONTROL: Liver, Testis, Cervix, Tonsil, Lymph Node, Spleen

LOCALIZATION: Cytoplasmic (Granular)

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-38

ISOTYPE: IgG1/K

CONTROL: Normal Pituitary

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP370

ISOTYPE: IgG


CONTROL: Thyroid, Salivary Gland, Breast, Ovarian Carcinoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2405 | Tinto Predilute | 3.0 ml | BSB 5596 | Tinto Predilute | 3.0 ml | BSB 3349 | Tinto Predilute | 3.0 ml |
| BSB 2406 | Tinto Predilute | 7.0 ml | BSB 5597 | Tinto Predilute | 7.0 ml | BSB 3350 | Tinto Predilute | 7.0 ml |
| BSB 2407 | Tinto Predilute | 15.0 ml | BSB 5598 | Tinto Predilute | 15.0 ml | BSB 3351 | Tinto Predilute | 15.0 ml |
| BSB 2408 | Concentrate | 0.1 ml | BSB 5599 | Concentrate | 0.1 ml | BSB 3352 | Concentrate | 0.1 ml |
| BSB 2409 | Concentrate | 0.5 ml | BSB 5600 | Concentrate | 0.5 ml | BSB 3353 | Concentrate | 0.5 ml |
| BSB 2410 | Concentrate | 1.0 ml | BSB 5601 | Concentrate | 1.0 ml | BSB 3354 | Concentrate | 1.0 ml |
| BSB 2411 | Control Slides | 5 | BSB 5602 | Control Slides | 5 | BSB 3355 | Control Slides | 5 |

HEG1,MMab



IHC of HEG1 on a FFPE Methotelioma Tissue

Sialylated protein HEG homolog 1 or Heart development protein with EGF like domains 1 (HEG1) is a novel calcium binding receptor protein involved in the Cerebral Cavensous Malformations (CCM) pathway and plays a regulatory role in heart and vessel formation by stabilizing endothelial cell junctions. HEG1 is a heavily glycosylated protein that participates in endothelial cell associations that develop the vascular system, and therefore plays a role in angiogenesis and the microenvironment of cellular adhesion.

The exact role of HEG1 in hepatocellular carcinoma is unclear, however, studies have shown HEG1 indicates poor prognosis and plays important roles in hepatocellular carcinoma invasion, metastasis, and epithelial mesenchymal transition by activating Wnt/ β -catenin signaling. HEG1 is potentially a useful prognostic biomarker and therapeutic target for malignant mesothelioma and hepatocellular, thyroid, and ovarian carcinomas.

HEG1 has been found to be a highly specific and sensitive marker of epithelioid Malignant Mesothelioma and offers sensitivity comparable to conventional markers for Epithelioid Mesotheliomas, but provides considerably better specificity, such that the diagnosis of epithelioid mesothelioma versus NSCLC and potentially could be confirmed with a combination of HEG1 and a suitable broad spectrum carcinoma marker such as Claudin-4.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: SKM9-2

ISOTYPE: IgG


CONTROL: Breast, Fallopian Tube, Colon Prostate, Testis, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma, Colon Adenocarcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3728-3 | Tinto Predilute | 3.0 ml |
| BSB-3728-7 | Tinto Predilute | 7.0 ml |
| BSB-3728-15 | Tinto Predilute | 15.0 ml |
| BSB-3728-01 | Concentrate | 0.1 ml |
| BSB-3728-05 | Concentrate | 0.5 ml |
| BSB-3728-1 | Concentrate | 1.0 ml |
| BSB-3728-CS | Control Slides | 5 |

Helicobacter pylori, RPaB



IHC of Helicobacter pylori on a FFPE Infected Stomach Tissue

Helicobacter pylori is a Gram-negative, microaerophilic bacterium found in the stomach, and may be present in other parts of the body, such as the eye. H. pylori is a helix-shaped Gram-negative bacterium about 3 μ m long with a diameter of about 0.5 μ m. It is microaerophilic; that is, it requires oxygen, but at lower concentration than is found in the atmosphere.

H. pylori could be present in patients with chronic gastritis and gastric ulcers and is also linked to the development of duodenal ulcers and stomach cancer. However, over 80% of individuals infected with the bacterium are asymptomatic and it may play an important role in the natural stomach ecology. The strain of H. pylori to which a person is exposed may influence the risk of developing gastric cancer. Strains of H. pylori that produce high levels of two proteins, vacuolating toxin A (VacA) and the cytotoxin-associated gene A (CagA), appear to cause greater tissue damage than those that produce lower levels or that lack those genes completely. These proteins are directly toxic to cells lining the stomach and signal strongly to the immune system that an invasion is underway.

This antibody reacts with Helicobacter pylori on the surface of epithelial cells of infected specimens.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG


CONTROL: H.Pylori Infected Stomach Mucosa

LOCALIZATION: Cell Wall

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5603 | Tinto Predilute | 3.0 ml |
| BSB 5604 | Tinto Predilute | 7.0 ml |
| BSB 5605 | Tinto Predilute | 15.0 ml |
| BSB 5606 | Concentrate | 0.1 ml |
| BSB 5607 | Concentrate | 0.5 ml |
| BSB 5608 | Concentrate | 1.0 ml |
| BSB 5609 | Control Slides | 5 |

Helicobacter pylori, MAb



IHC of Helicobacter pylori on Infected Stomach Tissue

Helicobacter pylori is a helix-shaped Gram-negative bacterium about 3 μ m long with a diameter of about 0.5 μ m. It is microaerophilic; that is, it requires oxygen, but at lower concentration than is found in the atmosphere. It contains a hydrogenase which can be used to obtain energy by oxidizing molecular hydrogen (H₂) produced by intestinal bacteria. It produces oxidase, catalase, and urease. H. pylori has four to six lophotrichous flagella; all gastric and enterohepatic Helicobacter species are highly motile owing to flagella. H. pylori's helical shape (from which the genus name is derived) is thought to have evolved to penetrate the mucoid lining of the stomach. Strains of H. pylori that produce high levels of two proteins, vacuolating toxin A (VacA) and the cytotoxin-associated gene A (CagA), appear to cause greater tissue damage than those that produce lower levels or that lack those genes completely.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-37

ISOTYPE: IgG1/K


CONTROL: H.Pylori Infected Stomach Mucosa

LOCALIZATION: Cell Wall

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2677 | Tinto Predilute | 3.0 ml |
| BSB 2678 | Tinto Predilute | 7.0 ml |
| BSB 2679 | Tinto Predilute | 15.0 ml |
| BSB 2680 | Concentrate | 0.1 ml |
| BSB 2681 | Concentrate | 0.5 ml |
| BSB 2682 | Concentrate | 1.0 ml |
| BSB 2683 | Control Slides | 5 |


Helicobacter pylori, RMAb



IHC of *Helicobacter pylori* on Infected Stomach Tissue

Helicobacter pylori is a helix-shaped Gram-negative bacterium about 3 µm long with a diameter of about 0.5 µm. It is microaerophilic; that is, it requires oxygen, but at lower concentration than is found in the atmosphere. It contains a hydrogenase which can be used to obtain energy by oxidizing molecular hydrogen (H₂) produced by intestinal bacteria. It produces oxidase, catalase, and urease. *H. pylori* has four to six lophotrichous flagella; all gastric and enterohepatic *Helicobacter* species are highly motile owing to flagella. *H. pylori*'s helical shape (from which the genus name is derived) is thought to have evolved to penetrate the mucoid lining of the stomach. Strains of *H. pylori* that produce high levels of two proteins, vacuolating toxin A (VacA) and the cytotoxin-associated gene A (CagA), appear to cause greater tissue damage than those that produce lower levels or that lack those genes completely.

Hemoglobin A, RMAb




IHC of Hemoglobin A on a FFPE Spleen Tissue

Hemoglobin alpha chain belongs to the globin family and is involved in oxygen transport from the lung to various peripheral tissues. It is a heterotetramer of two alpha chains and two beta chains in adult hemoglobin A (HbA); two alpha chains and two delta chains in adult hemoglobin A2 (HbA2). Hemoglobin alpha chain is expressed in red blood cells, and defects in HBA1/HBA2 can lead to alpha thalassemia, the most common of monogenic diseases.

Hemoglobin alpha chain is a useful marker for erythroid cells. An antibody to Hemoglobin alpha has been used for the identification of erythroid cells in myeloproliferative disease.

Hepatitis B Virus Core Antigen, RPAb



IHC of HBcAg on a FFPE Infected Liver Tissue

Hepatitis B virus is spherical in shape with a diameter of 42 nm. It contains a 27 nm partially double-stranded DNA core enclosed within a lipoprotein coat. The antigenic activity of the nucleocapsid core is designated as Hepatitis B core antigen.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP279
ISOTYPE: IgG
CONTROL: *H. Pylori* Infected Stomach Mucosa
LOCALIZATION: Cell Wall
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP124
ISOTYPE: IgG
CONTROL: Bone Marrow, Placenta, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Hepatitis B Infected Liver
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 3279 | Tinto Predilute | 3.0 ml | BSB 6786 | Tinto Predilute | 3.0 ml | BSB 5610 | Tinto Predilute | 3.0 ml |
| BSB 3280 | Tinto Predilute | 7.0 ml | BSB 6787 | Tinto Predilute | 7.0 ml | BSB 5611 | Tinto Predilute | 7.0 ml |
| BSB 3281 | Tinto Predilute | 15.0 ml | BSB 6788 | Tinto Predilute | 15.0 ml | BSB 5612 | Tinto Predilute | 15.0 ml |
| BSB 3282 | Concentrate | 0.1 ml | BSB 6789 | Concentrate | 0.1 ml | BSB 5613 | Concentrate | 0.1 ml |
| BSB 3283 | Concentrate | 0.5 ml | BSB 6790 | Concentrate | 0.5 ml | BSB 5614 | Concentrate | 0.5 ml |
| BSB 3284 | Concentrate | 1.0 ml | BSB 6791 | Concentrate | 1.0 ml | BSB 5615 | Concentrate | 1.0 ml |
| BSB 3285 | Control Slides | 5 | BSB 6792 | Control Slides | 5 | BSB 5616 | Control Slides | 5 |


Hepatitis B Virus Surface Antigen, MAb



IHC of HBsAg on a FFPE Infected Liver Tissue

Hepatitis B virus is spherical in shape with a diameter of 42 nm. It contains a 27 nm partially double-stranded DNA core enclosed within a lipoprotein coat.

Hepatocyte Specific Antigen/ Hep-Par1, MAb




IHC of Hepatocyte Specific Antigen on a FFPE Liver Tissue

Hepatocyte Specific Antigen (HSA or Hep Par 1) has been demonstrated consistently in the vast majority of Hepatocellular Carcinomas. Studies have shown the utility of HSA in the differential diagnosis of Hepatocellular Carcinoma, Cholangiocarcinoma and Hepatoblastomas.

HSA recognizes both benign and malignant liver derived tissues including such tumors as Hepatoblastoma, Hepatocellular Carcinoma, and Hepatic Adenoma. It recognizes both normal adult and fetal liver tissue. The typical pattern is a granular cytoplasmic staining. This antibody is useful in differentiating Hepatocellular Carcinomas with adenoid features from Adenocarcinomas, either primary in the liver or metastatic lesions to the liver. In recognizing Hepatoblastoma, it is useful in differentiating this entity from other small round cell tumors.

HER-2 neu, MAb



IHC of HER-2 neu on a FFPE Breast Carcinoma Tissue

HER-2 neu (also known as c-erbB-2) is a member of the epidermal growth factor receptor (EGFR) family. It is a cell membrane surface-bound tyrosine kinase and is normally involved in the signal transduction pathways leading to cell growth and differentiation.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: A10F1
ISOTYPE: IgG2b/K
CONTROL: Hepatitis B Infected Liver
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: OCH1E5
ISOTYPE: IgG1/K
CONTROL: Liver, Liver Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-3
ISOTYPE: IgG1
CONTROL: Breast Carcinoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5617 | Tinto Predilute | 3.0 ml |
| BSB 5618 | Tinto Predilute | 7.0 ml |
| BSB 5619 | Tinto Predilute | 15.0 ml |
| BSB 5620 | Concentrate | 0.1 ml |
| BSB 5621 | Concentrate | 0.5 ml |
| BSB 5622 | Concentrate | 1.0 ml |
| BSB 5623 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5624 | Tinto Predilute | 3.0 ml |
| BSB 5625 | Tinto Predilute | 7.0 ml |
| BSB 5626 | Tinto Predilute | 15.0 ml |
| BSB 5627 | Concentrate | 0.1 ml |
| BSB 5628 | Concentrate | 0.5 ml |
| BSB 5629 | Concentrate | 1.0 ml |
| BSB 5630 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5631 | Tinto Predilute | 3.0 ml |
| BSB 5632 | Tinto Predilute | 7.0 ml |
| BSB 5633 | Tinto Predilute | 15.0 ml |
| BSB 5634 | Concentrate | 0.1 ml |
| BSB 5635 | Concentrate | 0.5 ml |
| BSB 5636 | Concentrate | 1.0 ml |
| BSB 5637 | Control Slides | 5 |


HER-2 neu, R Mab



IHC of HER-2 neu on a FFPE Breast Carcinoma Tissue

HER-2 neu (also known as c-erbB-2) is a member of the epidermal growth factor receptor (EGFR) family and is notable for its role in the pathogenesis of breast cancer and as a target of treatment. It is a cell membrane surface-bound tyrosine kinase and is normally involved in the signal transduction pathways leading to cell growth and differentiation. HER-2 neu is a proto-oncogene located at the long arm of human chromosome 17(17q11.2-q12).

HER2, R Mab




IHC of HER-2 on a FFPE Breast Carcinoma Tissue

HER2 /ErbB2 is one of four members of the ErbB receptor family of transmembrane receptor-like tyrosine kinases .The kinase activity of ErbB2 can be activated without a ligand if it is overexpressed, and by association with other ErbB proteins.

HER2 is overexpressed in 25–30% of all breast cancers, including primary as well as metastatic breast tumors. HER2 has been widely investigated as a prognostic indicator.

HER2 neu Phospho, R Mab



IHC of HER-2 neu Phospho on a FFPE Breast Carcinoma Tissue

HER2 / ErbB2 is one of four members of the ErbB receptor family of transmembrane receptor-like tyrosine kinases .The kinase activity of ErbB2 can be activated without a ligand if it is overexpressed, and by association with other ErbB proteins. Autophosphorylation did not modulate receptor turnover. A Tyr→Phe substitution of ErbB2 Tyr-877 homologous to pp60c-src Tyr-416 did not alter ErbB2 biological and biochemical properties, thus excluding the possibility that phosphorylation of this residue, located in the kinase domain, modulates ErbB2 gp185 catalytic function.

HER2 is overexpressed in 25–30% of all breast cancers, including primary as well as metastatic breast tumors. HER2 has been widely investigated as a prognostic indicator. Detection of phosphorylated HER2 may aid in predicting breast cancer progression.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-HER2
ISOTYPE: IgG
CONTROL: Breast Carcinoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP3
ISOTYPE: IgG
CONTROL: Breast Carcinoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP123
ISOTYPE: IgG
CONTROL: Breast Carcinoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2035 | Tinto Predilute | 3.0 ml |
| BSB 2036 | Tinto Predilute | 7.0 ml |
| BSB 2037 | Tinto Predilute | 15.0 ml |
| BSB 2038 | Concentrate | 0.1 ml |
| BSB 2039 | Concentrate | 0.5 ml |
| BSB 2040 | Concentrate | 1.0 ml |
| BSB 2041 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2510 | Tinto Predilute | 3.0 ml |
| BSB 2511 | Tinto Predilute | 7.0 ml |
| BSB 2512 | Tinto Predilute | 15.0 ml |
| BSB 2513 | Concentrate | 0.1 ml |
| BSB 2514 | Concentrate | 0.5 ml |
| BSB 2515 | Concentrate | 1.0 ml |
| BSB 2516 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2503 | Tinto Predilute | 3.0 ml |
| BSB 2504 | Tinto Predilute | 7.0 ml |
| BSB 2505 | Tinto Predilute | 15.0 ml |
| BSB 2506 | Concentrate | 0.1 ml |
| BSB 2507 | Concentrate | 0.5 ml |
| BSB 2508 | Concentrate | 1.0 ml |
| BSB 2509 | Control Slides | 5 |

HER-3/c-erbB-3, RMAb



IHC of HER-3 on a FFPE Colon Tissue


Receptor tyrosine-protein kinase erbB-3, also known as HER-3 (Human Epidermal Growth Factor Receptor 3), is a membrane bound protein that in humans is encoded by the ERBB3 gene. HER-3 has been shown to bind the ligands heregulin and NRG-2. Ligand binding causes a change in conformation that allows for dimerization, phosphorylation, and activation of signal transduction. HER-3 can heterodimerize with any of the other three ErbB family members.

During human development, HER-3 is expressed in skin, bone, muscle, nervous system, heart, lungs, and intestinal epithelium and is expressed in normal adult human gastrointestinal tract, reproductive system, skin, nervous system, urinary tract, and endocrine system. HER-3 is overexpressed in a variety of tumors including breast, stomach, pancreas, and colon. While no evidence has been found that HER-3 overexpression, constitutive activation, or mutation alone is oncogenic, the protein as a heterodimerization partner, most critically with HER-2, is implicated in growth, proliferation, resistance, and the promotion of invasion and metastasis.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-HER3
ISOTYPE: IgG
CONTROL: Colon, Cervix, Bladder TCC
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2629 | Tinto Predilute | 3.0 ml |
| BSB 2630 | Tinto Predilute | 7.0 ml |
| BSB 2631 | Tinto Predilute | 15.0 ml |
| BSB 2632 | Concentrate | 0.1 ml |
| BSB 2633 | Concentrate | 0.5 ml |
| BSB 2634 | Concentrate | 1.0 ml |
| BSB 2635 | Control Slides | 5 |

Herpes Simplex Virus I, RPAb




IHC of HSV I on a FFPE Infected Tissue

Herpes Simplex Virus I usually infects the non-genital mucosal surfaces, and may also affect skin or internal organs such as brain, lung, liver, adrenal gland, or GI tract of immunocompromised individuals.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: HSV I Infected Tissues
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2426 | Tinto Predilute | 3.0 ml |
| BSB 2427 | Tinto Predilute | 7.0 ml |
| BSB 2428 | Tinto Predilute | 15.0 ml |
| BSB 2429 | Concentrate | 0.1 ml |
| BSB 2430 | Concentrate | 0.5 ml |
| BSB 2431 | Concentrate | 1.0 ml |
| BSB 2432 | Control Slides | 5 |

Herpes Simplex Virus I, MAb



IHC of HSV I on a FFPE Infected Tissue


Herpes simplex virus I (HSV-I) is a strain of the Herpes virus family, Herpesviridae, which cause infections in humans. The double stranded DNA genome is contained within an icosahedral capsid embedded in a proteinaceous layer (tegument) and surrounded by a lipid envelope, derived from the nuclear membrane of the last host, which is decorated with virus-specific glycoproteins spikes.

HSV I causes a contagious disease, also known as cold sore, night fever, or fever blister. After an initial, or primary, infection, HSV establishes latency, during which the virus is present in the cell bodies of nerves which innervate the area of original outbreak. Herpes symptoms may periodically recur in the form of outbreaks of herpetic sores near the site of original infection. HSV I usually infects the non-genital mucosal surfaces, and Type II typically involves the genitalia. Either type may affect the skin or internal organs (typically brain, lung, liver, adrenal gland, or GI tract) of immunocompromised individuals.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 10A3
ISOTYPE: IgG1
CONTROL: HSV I Infected Tissues
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2684 | Tinto Predilute | 3.0 ml |
| BSB 2685 | Tinto Predilute | 7.0 ml |
| BSB 2686 | Tinto Predilute | 15.0 ml |
| BSB 2687 | Concentrate | 0.1 ml |
| BSB 2688 | Concentrate | 0.5 ml |
| BSB 2689 | Concentrate | 1.0 ml |
| BSB 2690 | Control Slides | 5 |


Herpes Simplex Virus II, RPab



IHC of HSV II on a FFPE Infected Tissue

Herpes Simplex Virus II typically involves the genitalia, and may also affect skin or internal organs such as brain, lung, liver, adrenal gland, or GI tract of immunocompromised individuals.


Herpes Simplex Virus II, MAb



IHC of HSV II on a FFPE Infected Tissue

Herpes Simplex Virus II typically involves the genitalia, and may also affect skin or internal organs such as brain, lung, liver, adrenal gland, or GI tract of immunocompromised individuals

Herpes Simplex Virus I & II, RPab



IHC of HSV I & II on a FFPE Infected Cell Line

The antibody used reacts with antigens common to HSV types 1 and 2. The antibody reacts with all the major glycoproteins present in the viral envelope and at least one core protein as determined by crossed immunoelectrophoresis.

Contaminating antibodies to human and bovine serum have been removed by solid-phase absorption. The antibody shows no reaction with human and bovine plasma when tested by ELISA. The antibody shows no cross-reactivity to cytomegalovirus and Epstein-Barr virus.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: HSV II Infected Tissues
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-116
ISOTYPE: IgG2A
CONTROL: HSV II Infected Tissue
LOCALIZATION: Nuclear, Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: HSV I & II Infected Tissues
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2335 | Tinto Predilute | 3.0 ml |
| BSB 2336 | Tinto Predilute | 7.0 ml |
| BSB 2337 | Tinto Predilute | 15.0 ml |
| BSB 2338 | Concentrate | 0.1 ml |
| BSB 2339 | Concentrate | 0.5 ml |
| BSB 2340 | Concentrate | 1.0 ml |
| BSB 2341 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3462 | Tinto Predilute | 3.0 ml |
| BSB 3463 | Tinto Predilute | 7.0 ml |
| BSB 3464 | Tinto Predilute | 15.0 ml |
| BSB 3465 | Concentrate | 0.1 ml |
| BSB 3466 | Concentrate | 0.5 ml |
| BSB 3467 | Concentrate | 1.0 ml |
| BSB 3468 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5638 | Tinto Predilute | 3.0 ml |
| BSB 5639 | Tinto Predilute | 7.0 ml |
| BSB 5640 | Tinto Predilute | 15.0 ml |
| BSB 5641 | Concentrate | 0.1 ml |
| BSB 5642 | Concentrate | 0.5 ml |
| BSB 5643 | Concentrate | 1.0 ml |
| BSB 5644 | Control Slides | 5 |

Herpes Simplex Virus I & II, MAb




IHC of HSV I & II on a FFPE Infected Tissue

hGAL/GCET2, RMAb



IHC of HGAL on a FFPE Tonsil Tissue

HHV-8, MAb



IHC of HHV-8 on a FFPE Pleura Tissue

The antibody reacts with antigens common to HSV types 1 and 2. The antibody reacts with all major glycoproteins present in the viral envelope and at least one core protein as determined by crossed immunoelectrophoresis.

Contaminating antibodies to human and bovine serum have been removed by solid-phase absorption. The antibody shows no reaction with human and bovine plasma when tested by ELISA. The antibody shows no cross-reactivity to cytomegalovirus and Epstein-Barr virus.

GAL, also known as Germinal Center-associated Lymphoma Protein or GCET, is an important prognostic marker for staging Lymphomas derived from germinal centers. Analysis of the GCET2 protein sequence indicated that it may be involved in signal transduction in the cytoplasm. Two newly characterized germinal center B-cell-associated genes, GCET1 and GCET2, have differential expression in normal and neoplastic B cells. Expression of the HGAL gene is specifically induced in B cells by interleukin-4 (IL-4).

The HGAL protein has been shown to be expressed in the cytoplasm of germinal center B lymphocytes and in B cell lymphomas of germinal center derivation. HGAL is absent in Mantle and Marginal zone B-cells and in the Interfollicular and Paracortical regions in normal Tonsils and Lymph Nodes. HGAL is an ideal marker for the detection of Germinal Center-derived B-cell Lymphomas and has the highest overall sensitivity of detecting Follicular Lymphoma. HGAL has been identified in gene-expression profiling studies of Diffuse Large B-Cell Lymphoma (DLBCL). Among 727 Lymphomas tested by immunohistochemistry on tissue microarrays, HGAL staining was found in Follicular Lymphomas (103 of 107), Burkitt Lymphomas (40 of 40), Mediastinal Large B Lymphomas (7 of 8), and in DLBCLs (103 of 151). Expression of HGAL protein identifies a subset of classic Hodgkin Lymphoma of Germinal Center derivation with improved survival.

Kaposi's Sarcoma-associated herpes virus is the eighth human herpes virus; its formal name according to the International Committee on Taxonomy of Viruses is HHV-8. Anti-HHV-8 labels the latent nuclear antigen protein via immunohistochemistry.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 10A3/BSB-116
ISOTYPE: IgG2a
CONTROL: HSV I & II Infected Tissue
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP316
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Thymus, Fallopian Tube, Germinal Center B-cell Type Diffuse Large B-cell Lymphoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: 13B10
ISOTYPE: IgG1
CONTROL: Kaposi's Sarcoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3469 | Tinto Predilute | 3.0 ml |
| BSB 3470 | Tinto Predilute | 7.0 ml |
| BSB 3471 | Tinto Predilute | 15.0 ml |
| BSB 3572 | Concentrate | 0.1 ml |
| BSB 3573 | Concentrate | 0.5 ml |
| BSB 3574 | Concentrate | 1.0 ml |
| BSB 3575 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2936 | Tinto Predilute | 3.0 ml |
| BSB 2937 | Tinto Predilute | 7.0 ml |
| BSB 2938 | Tinto Predilute | 15.0 ml |
| BSB 2939 | Concentrate | 0.1 ml |
| BSB 2940 | Concentrate | 0.5 ml |
| BSB 2941 | Concentrate | 1.0 ml |
| BSB 2942 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5645 | Tinto Predilute | 3.0 ml |
| BSB 5646 | Tinto Predilute | 7.0 ml |
| BSB 5647 | Tinto Predilute | 15.0 ml |
| BSB 5648 | Concentrate | 0.1 ml |
| BSB 5649 | Concentrate | 0.5 ml |
| BSB 5650 | Concentrate | 1.0 ml |
| BSB 5651 | Control Slides | 5 |

HHV-8, RMAb




IHC of HHV-8 on a FFPE Kaposi's Sarcoma Tissue

HIF-1alpha, RMAb



IHC of HIF-1 alpha on a FFPE Anal Carcinoma Tissue

Histone H3 Phospho, RMAb



IHC of Histone H3 Phospho on a FFPE NH Lymphoma Tissue

Kaposi's sarcoma-associated herpesvirus is the eighth human herpes virus; its formal name according to the International Committee on Taxonomy of Viruses is HHV-8. It causes development of cutaneous lesions such as macules, plaques, and nodules, and may invade tissues such as lymph nodes.

It is associated with Kaposi's sarcoma (KS), primary effusion lymphoma, and multicentric Castleman's disease. HHV-8 detects the latent nuclear antigen (LNA) protein in all stages and pathological presentations of KS, which can differentiate KS from HIV-induced tumors or other neoplasia of similar morphology. Immunohistochemical detection of HHV-8 LNA is particularly useful in early detection, where inflammation but not significant tumor formation has occurred.

Hypoxia-inducible factor 1-alpha, also known as HIF-1-alpha, is a protein that in humans is encoded by the HIF1A gene. Two alternative transcripts encoding different isoforms have been identified, the alpha and beta subunits. HIF-1 plays an essential role in cellular and systemic responses to hypoxia.

HIF-1 is a critical mediator of the hypoxic response that upregulates expression of proteins that promote angiogenesis, anaerobic metabolism, and many other survival pathways. HIF-1alpha is expressed in many types of tumors. The expression of HIF-1alpha is correlated with tumor angiogenesis, cancer progression and clinical outcome in various solid tumors including breast cancer, type 1 endometrial carcinoma, sarcoma, head and neck tumor and brain tumor. Importantly, hypoxia has been clinically demonstrated to predict an adverse treatment outcome in the radiotherapeutic management of cancer of the head and neck, uterine cervix and soft tissue sarcomas. HIF-1 alpha may be of value in analyzing the cancer cell response to therapy.

Phosphohistone-H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Featuring a main globular domain and a long N-terminal tail, H3 is involved with the structure of the nucleosomes of the "beads on a string" structure. In mammalian cells, phosphohistone H3 is negligible during interphase but reaches a maximum for chromatin condensation during mitosis.

Phosphohistone-H3 can serve as a mitotic marker to separate mitotic figures from apoptotic bodies and karyorrhectic debris, which may be a very useful tool in diagnosis of tumor grades, especially in CNS, skin, Gyn., soft tissue, and GIST.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-HHV8
ISOTYPE: IgG
CONTROL: Kaposi's Sarcoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP118
ISOTYPE: IgG
CONTROL: Tonsil, Cervix, Liver, Testis, Kidney, Breast, Thymus, Spleen, Colon, Cervical Cancer
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP233
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB-3778-3 | Tinto Predilute | 3.0 ml | BSB 2517 | Tinto Predilute | 3.0 ml | BSB 2524 | Tinto Predilute | 3.0 ml |
| BSB-3778-7 | Tinto Predilute | 7.0 ml | BSB 2518 | Tinto Predilute | 7.0 ml | BSB 2525 | Tinto Predilute | 7.0 ml |
| BSB-3778-15 | Tinto Predilute | 15.0 ml | BSB 2519 | Tinto Predilute | 15.0 ml | BSB 2526 | Tinto Predilute | 15.0 ml |
| BSB-3778-01 | Concentrate | 0.1 ml | BSB 2520 | Concentrate | 0.1 ml | BSB 2527 | Concentrate | 0.1 ml |
| BSB-3778-05 | Concentrate | 0.5 ml | BSB 2521 | Concentrate | 0.5 ml | BSB 2528 | Concentrate | 0.5 ml |
| BSB-3778-1 | Concentrate | 1.0 ml | BSB 2522 | Concentrate | 1.0 ml | BSB 2529 | Concentrate | 1.0 ml |
| BSB-3778-CS | Control Slides | 5 | BSB 2523 | Control Slides | 5 | BSB 2530 | Control Slides | 5 |

HLA-DR alpha chain, RMAb




IHC of HLA-DR alpha-chain on a FFPE Skin Tissue

HLA class II histocompatibility antigen, DR alpha chain (HLA-DR alpha chain) is a polymorphic cell surface glycoprotein that is crucial for the cellular interaction in the immune system. Class II molecules have limited tissue distribution and are predominantly expressed on B lymphocytes and macrophage; these class II molecules present peptides derived from extracellular proteins to T cells.

HLA-DR alpha chain antibody labels B cells, dendritic cells, monocytes and macrophages. It is also reported to react with tumor cells in many types of cancers including breast, liver, lung and ovary cancer.

HMGA2, RMAb




IHC of HMGA2 on a FFPE Lung Squamous Cell Carcinoma Tissue

High-mobility group AT-hook 2 (HMGA2) belongs to the architectural transcription factor HMGA family and is encoded by the HMGA2 gene. HMGA2 plays a role in chromosomal organization and transcriptional regulation. HMGA2 has three basic DNA-binding domains (AT-hooks) that bind to AT rich regions of nuclear DNA and alter the structure of DNA to promote the assembly of protein complexes that regulate transcription. With few exceptions, HMGA2 is expressed in humans only during early development, and is reduced to undetectable levels of transcription in adult tissues.

Elevated expression of HMGA2 is found in a variety of human cancers, correlates with metastasis and poor prognosis for patients. High HMGA2 expression have been reported in Pituitary Adenoma, Thyroid Carcinoma, Triple-Negative Breast Carcinoma, Breast Carcinoma, Lung Adenocarcinoma, Colorectal Carcinoma, Hepatoblastoma, Pancreatic Adenocarcinoma, Conventional and Intramuscular Lipoma, Liposarcoma, Gastric, and Ovarian Tumors and other conditions. HMGA2 is expressed in most conventional and intramuscular lipomas and can aid in differentiating between Lipomas from dedifferentiated Liposarcomas and distinguishing areas of tumor from normal adipose tissue. In Mesenchymal tumors, HMGA2 is expressed in benign Fibrous Histiocytoma, Nodular Fasciitis, and Vulvovaginal Angiomyxoma. In Thyroid Carcinomas, upregulation of HMGA2 can distinguish between benign and malignant Follicular neoplasias.

hPL, RPAb



IHC of hPL on a FFPE Placenta Tissue

Human placenta lactogen (hPL), also known as human chorionic somatomammotropin (HCS), is a polypeptide placental hormone. Its structure and function are similar to that of human growth hormone. It modifies the metabolic state of the mother during pregnancy to facilitate the energy supply of the fetus. It is first detectable in the maternal serum in the fifth week of gestation and reaches a plateau by the thirty-fourth week.

hPL is expressed in the syncytiotrophoblastic cells of choriocarcinoma. A rare variant of trophoblastic tumor has been reported in the testis with resemblance to uterine placental site trophoblastic tumor. It consisted purely of intermediate trophoblasts, which was diffusely positive for hPL and focally for B-hCG.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP96

ISOTYPE: IgG

CONTROL: Tonsil, Spleen, Liver, Kidney, Adrenal, Colon, Lymph Node, Thymus, Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP398

ISOTYPE: IgG

CONTROL: K Cervix, Lung Squamous Cell Carcinoma, Papillary Thyroid Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Placenta

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6793 | Tinto Predilute | 3.0 ml |
| BSB 6794 | Tinto Predilute | 7.0 ml |
| BSB 6795 | Tinto Predilute | 15.0 ml |
| BSB 6796 | Concentrate | 0.1 ml |
| BSB 6797 | Concentrate | 0.5 ml |
| BSB 6798 | Concentrate | 1.0 ml |
| BSB 6799 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3729-3 | Tinto Predilute | 3.0 ml |
| BSB-3729-7 | Tinto Predilute | 7.0 ml |
| BSB-3729-15 | Tinto Predilute | 15.0 ml |
| BSB-3729-01 | Concentrate | 0.1 ml |
| BSB-3729-05 | Concentrate | 0.5 ml |
| BSB-3729-1 | Concentrate | 1.0 ml |
| BSB-3729-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6800 | Tinto Predilute | 3.0 ml |
| BSB 6801 | Tinto Predilute | 7.0 ml |
| BSB 6802 | Tinto Predilute | 15.0 ml |
| BSB 6803 | Concentrate | 0.1 ml |
| BSB 6804 | Concentrate | 0.5 ml |
| BSB 6805 | Concentrate | 1.0 ml |
| BSB 6806 | Control Slides | 5 |

hPL, RMAb




IHC of hPL on a FFPE Placenta Tissue

HPV, MAb



IHC of HPV on a FFPE LSIL of the Cervix Tissue

HPV16, MAb



IHC of HPV16 on a FFPE Infected Cervix Tissue

Human placenta lactogen (hPL), also known as human chorionic somatomammotropin (HCS), is a polypeptide placental hormone. Its structure and function are similar to that of human growth hormone. It modifies the metabolic state of the mother during pregnancy to facilitate the energy supply of the fetus. It is first detectable in the maternal serum in the fifth week of gestation and reaches a plateau by the thirty-fourth week.

hPL is expressed in the syncytiotrophoblastic cells of choriocarcinoma. A rare variant of trophoblastic tumor has been reported in the testis with resemblance to uterine placental site trophoblastic tumor. It consisted purely of intermediate trophoblasts, which was diffusely positive for hPL and focally for B-hCG.

Papillomaviruses are a diverse group of DNA-based viruses. More than 100 different human papillomavirus (HPV) types have been characterized. Some HPV types cause benign skin warts, or papillomas, for which the virus family is named.

Anti-human papillomavirus, clone BSB-66 (SB24) reacts with an epitope of a major capsid protein of HPV, which is broadly expressed among different HPV subtypes

Papillomaviridae belongs to a taxonomic family of non-enveloped DNA virus, collectively known as papillomavirus. Several hundred types of papillomavirus, have been identified infecting mammals and also other amniotes such as birds, snakes and turtles. Infection by most papillomavirus types, depending on the type, is either asymptomatic (e.g. most Beta-PVs) or causes small benign tumors, known as papillomas or warts (e.g. HPV1, HPV6 or HPV11). Papillomas caused by some types, however, such as HPV16 and HPV18, carry a risk of becoming cancerous.

Papillomaviruses are usually considered as highly host- and tissue-tropic, and are thought to rarely be transmitted between species. Papillomaviruses replicate exclusively in the basal layer of the body surface tissues. All known papillomavirus types infect a particular body surface, typically the skin or mucosal epithelium of the genitals, anus, mouth, or airways. Some papillomavirus types can cause cancer in the epithelial tissues they inhabit, cancer is not a typical outcome of infection. The development of papillomavirus-induced cancers typically occurs over the course of many years. Papillomaviruses have been associated with the development of cervical cancer, penile cancer and oral cancers. An association with vulvar cancer and urothelial carcinoma with squamous differentiation in patients with neurogenic bladder has also been reported.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP241
ISOTYPE: IgG
CONTROL: Placenta
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-66
ISOTYPE: IgG1/K
CONTROL: HPV Infected Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: CAMVIR-1
ISOTYPE: IgG2a
CONTROL: HPV16 Infected Tissues
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2412 | Tinto Predilute | 3.0 ml |
| BSB 2413 | Tinto Predilute | 7.0 ml |
| BSB 2414 | Tinto Predilute | 15.0 ml |
| BSB 2415 | Concentrate | 0.1 ml |
| BSB 2416 | Concentrate | 0.5 ml |
| BSB 2417 | Concentrate | 1.0 ml |
| BSB 2418 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5652 | Tinto Predilute | 3.0 ml |
| BSB 5653 | Tinto Predilute | 7.0 ml |
| BSB 5654 | Tinto Predilute | 15.0 ml |
| BSB 5655 | Concentrate | 0.1 ml |
| BSB 5656 | Concentrate | 0.5 ml |
| BSB 5657 | Concentrate | 1.0 ml |
| BSB 5658 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2943 | Tinto Predilute | 3.0 ml |
| BSB 2944 | Tinto Predilute | 7.0 ml |
| BSB 2945 | Tinto Predilute | 15.0 ml |
| BSB 2946 | Concentrate | 0.1 ml |
| BSB 2947 | Concentrate | 0.5 ml |
| BSB 2948 | Concentrate | 1.0 ml |
| BSB 2949 | Control Slides | 5 |

HSP-27, MAb



IHC of HSP-27 on a FFPE Lymphoblastic Lymphoma Tissue

Heat shock protein 27 (HSP27) also known as heat shock protein beta-1 (HSPB1) is a protein that in humans is encoded by the HSPB1 gene. The common functions of sHsps are chaperone activity, thermotolerance, inhibition of apoptosis, regulation of cell development, and cell differentiation. HSP27 appears in many cell types, especially all types of muscle cells. It is overexpressed during different stages of cell differentiation and development. This suggests an essential role for HSP27 in the differentiation of tissues.

An affinity of high expression levels of different phosphorylated HSP27 species and muscle/neurodegenerative diseases and various cancers has been observed. High levels of HSP27 were also found in sera of breast cancer patients; therefore HSP27 could be a potential diagnostic marker. Phosphorylated Hsp27 has been shown to increase in human prostate cancer (PCa) cell invasion, enhancing cell proliferation, and suppression of Fas-induced apoptosis in human PCa cells.

HSP27 immunohistochemistry is a useful tool for the identification of CIN and cervical squamous cell carcinoma, and is a good complement to p16. HSP27 has been demonstrated to be overexpressed in cervical intraepithelial neoplasia (CIN) and squamous cell carcinoma of the cervix using immunohistochemistry. HSP27 expression has been demonstrated in 47% of CIN1, 75% of CIN2, 92% of CIN3, and 100% of cervical squamous cell carcinomas (SCC).

ANTIBODY TYPE: Mouse Monoclonal

CLONE: G3.1


ISOTYPE: IgG1

CONTROL: Tonsil, Cervix, Prostate, Breast Carcinoma, Cervical Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Rat, Dog, Mouse, Primate

HSP70, RMAb



IHC of HSP70 on a FFPE Colon Adenocarcinoma Tissue

70 kDa heat shock proteins (HSP70) are found ubiquitously in virtually all living organisms, facilitating protein folding and protecting cells from heat stress and toxic chemicals. HSP70 proteins have 3 functional domains: N-terminal ATPase domain, substrate binding domain, and a C-terminal domain that serves as a "lid" for the substrate binding domain. HSP70 binds tightly to partially synthesized peptides and prevents them from aggregating and rendering nonfunctional.

HSP70 is shown to be overexpressed in malignant Melanoma and underexpressed in Renal Cell Carcinoma. A variety of tumor cells can express HSP70 with seemingly contradictory functions. Intracellular HSP70 has a cytoprotective function via suppression of apoptosis and lysosomal cell death (LCD) and extracellular HSP70 can promote tumorigenesis and angiogenesis. Other evidence showed intracellular HSP70 can promote apoptosis and membrane-associated/extracellular HSP70 can elicit antitumor innate and adaptive immune responses.

One study evaluated the expression of HSP70, Estrogen Receptor (ER) and Ki-67 and assessed the relationship between them in Cervical Squamous Cell Neoplasia. It found that HSP70 may play an important role in tumor cell proliferation and is more related with invasive Squamous Cell Carcinoma than Cervical Intraepithelial Neoplasia, but ER may be not related with tumor cell proliferation and differentiation. Therefore, HSP70 may be a useful prognostic factor in Cervical Dysplasia and Cancer.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM432


ISOTYPE: IgG

CONTROL: Breast, Fallopian Tube, Skin, Prostate, Testis, Transitional Cell Carcinoma

LOCALIZATION: Nuclear, Cytoplasmic

SPECIES REACTIVITY: Human

ICOS/CD278, RMAb



IHC of ICOS/CD278 on a FFPE Tonsil Tissue

CD278 is also known as inducible T-cell costimulatory molecule (ICOS). CD278 is homologous to the CD28/CTLA-4 proteins and is expressed on activated T cells and unstimulated thymocytes. CD278 plays a major role in regulation of cell-cell signaling, adaptive immune response, and cell proliferation. Interaction of CD278 and its ligand ICOS-L results in increased production of interleukins, which promote differentiation of Tfh and Tregs and development of Th1, Th2, and Th17 cells.

CD278/ICOS-L interaction has been shown to

promote either antitumor T cell response or pro-tumor responses when triggered in Tregs (such as in Multiple Myeloma, Acute Myeloid Leukemia, and some invasive Breast Carcinomas). Therefore, both agonistic and antagonistic monoclonal antibodies targeting this pathway can be potential cancer immunotherapy.

ICOS is primarily expressed on activated CD4+ and CD8+ T cells where it regulates immune responses and plays a role in the regulation of T Follicular helper cells. ICOS is a sensitive marker for identifying T cell Lymphomas of Follicular Helper T cell origin, especially certain patterns of Angioimmunoblastic T-cell lymphoma (AITL) and Peripheral T-cell lymphomas with T-follicular Helper Phenotype (PTCL-TFH). Immunohistochemical analysis revealed that ICOS is widely expressed by malignant cells in skin biopsy specimens from patients with Mycosis Fungoides and Sézary syndrome (SS), as well as involved in node biopsy specimens from patients with SS.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM417

ISOTYPE: IgG

CONTROL: Colon, Tonsil, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma

LOCALIZATION: Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2950 | Tinto Predilute | 3.0 ml |
| BSB 2951 | Tinto Predilute | 7.0 ml |
| BSB 2952 | Tinto Predilute | 15.0 ml |
| BSB 2953 | Concentrate | 0.1 ml |
| BSB 2954 | Concentrate | 0.5 ml |
| BSB 2955 | Concentrate | 1.0 ml |
| BSB 2956 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3730-3 | Tinto Predilute | 3.0 ml |
| BSB-3730-7 | Tinto Predilute | 7.0 ml |
| BSB-3730-15 | Tinto Predilute | 15.0 ml |
| BSB-3730-01 | Concentrate | 0.1 ml |
| BSB-3730-05 | Concentrate | 0.5 ml |
| BSB-3730-1 | Concentrate | 1.0 ml |
| BSB-3730-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3731-3 | Tinto Predilute | 3.0 ml |
| BSB-3731-7 | Tinto Predilute | 7.0 ml |
| BSB-3731-15 | Tinto Predilute | 15.0 ml |
| BSB-3731-01 | Concentrate | 0.1 ml |
| BSB-3731-05 | Concentrate | 0.5 ml |
| BSB-3731-1 | Concentrate | 1.0 ml |
| BSB-3731-CS | Control Slides | 5 |

IDH1 R132H, MAb




IHC of IDH1 R132H on a FFPE Glioma Tissue

IDH1 R132H, RMAb



IHC of IDH1 R132H on a FFPE Glioblastoma Tissue

IFN-Alpha, MAb



IHC of IFN - Alpha on a FFPE Papillary Thyroid Carcinoma Tissue

The IDH1 gene on chromosome 2q33.3 encodes for isocitrate dehydrogenase 1 (IDH1), located in the cytoplasm and the peroxisomes. This enzyme catalyzes NADPH production via oxidative decarboxylation of isocitrate to alpha-ketoglutarate in the Krebs citric acid cycle. Studies have shown that IDH1 mutation is an early step in gliomagenesis and has been reported to occur in grades II and III astrocytomas, oligodendrogliomas (OG), oligoastrocytomas (OA), and secondary GBM.

Mutations involving IDH1 occur in a high proportion of diffuse gliomas, with implications on diagnosis. About 90% involve exon 4 at codon 132, replacing amino acid arginine with histidine (R132H). Preliminary studies comparing Immunohistochemistry (IHC) with IDH1-R132H mutation-specific antibodies have shown concordance with DNA sequencing and no cross-reactivity with wild-type IDH1 or other mutant proteins.

IDH1 R132H is an isocitrate dehydrogenase isozyme and encoded by the gene IDH1. IDH1 is involved in the citric acid cycle during glucose metabolism and catalyzes the oxidation of isocitrate to alpha-ketoglutarate and reduction of NADP+ to NADPH. Both alpha-ketoglutarate and NADPH play a role in protecting cells from oxidative stress and mitigating oxidative damage. Mutation in residue 132 of IDH1 results in loss of enzymatic function, buildup of 2-hydroxyglutarate, and changes in histone and DNA methylation.

Mutation of IDH1 R132H is implicated in metaphyseal chondromatosis with aciduria as well as diffused gliomas and a number of neoplasms such as acute myeloid leukemia, acute lymphocytic leukemia, myelofibrosis, intrahepatic cholangiocarcinoma, melanoma, chondroid tumors, and certain rare forms of colonic and prostate carcinomas. However studies have shown that IDH1 mutation is not a direct trigger of oncogenesis, but strongly associated with other tumor-promoting mutations. Screening for IDH1 R132H mutation can provide valuable information on diagnosis and prognosis of glioma. Glioma with IDH1 R132H mutation tend to be less aggressive than glioma without IDH1 R132H mutation.

Type I interferons are a large subgroup of interferon proteins that help regulate the immune system by binding to IFN- alpha receptors. IFN- alpha is produced mainly by plasmacytoid dendritic cells and involved in innate immunity against viral infections. Binding of IFN-alpha to its receptor leads to downstream signalling and expression of numerous different IFN-stimulated genes. These genes encode antiviral proteins that directly inhibit viral replication as well as modulate immune function.

Studies have shown IFN- α strongly modulates innate and adaptive immune responses in the host by enhancing the proliferation, cytotoxicity and IFN- γ secretion of NK cells, as well as acting as a pyrogenic factor by altering the activity of thermosensitive neurons in the hypothalamus which causes fever. It does this by binding to opioid receptors and eliciting the release of prostaglandin-E2. IFN- α can also interact with the μ -opioid receptor to act as an analgesic. Additionally, recent studies have shown that Type 1 IFNs stimulate secretion of IP-10 (CXCL10) which is a critical chemokine to recruit effector T cells to the tumor microenvironment and IP-10 knockout mice exhibit a phenotype with compromised effector T cell generation and trafficking. Type 1 IFNs also induce MHC class 1 upregulation on tumor cells which can enhance antitumor CD8 T cell effector response in the tumor microenvironment.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: IHC132
ISOTYPE: IgG1
CONTROL: Glioma, Glioblastoma, Astrocytoma with IDH1 R132H Mutation
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-IDH1
ISOTYPE: IgG
CONTROL: Colon, Tonsil, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-158
ISOTYPE: IgG2b
CONTROL: Placenta, Fallopian Tube, Stomach, Prostate, Testis, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma
LOCALIZATION: Membranous, Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3546 | Tinto Predilute | 3.0 ml |
| BSB 3547 | Tinto Predilute | 7.0 ml |
| BSB 3548 | Tinto Predilute | 15.0 ml |
| BSB 3549 | Concentrate | 0.1 ml |
| BSB 3550 | Concentrate | 0.5 ml |
| BSB 3551 | Concentrate | 1.0 ml |
| BSB 3552 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3732-3 | Tinto Predilute | 3.0 ml |
| BSB-3732-7 | Tinto Predilute | 7.0 ml |
| BSB-3732-15 | Tinto Predilute | 15.0 ml |
| BSB-3732-01 | Concentrate | 0.1 ml |
| BSB-3732-05 | Concentrate | 0.5 ml |
| BSB-3732-1 | Concentrate | 1.0 ml |
| BSB-3732-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3733-3 | Tinto Predilute | 3.0 ml |
| BSB-3733-7 | Tinto Predilute | 7.0 ml |
| BSB-3733-15 | Tinto Predilute | 15.0 ml |
| BSB-3733-01 | Concentrate | 0.1 ml |
| BSB-3733-05 | Concentrate | 0.5 ml |
| BSB-3733-1 | Concentrate | 1.0 ml |
| BSB-3733-CS | Control Slides | 5 |

IFN-Gamma, MMab



IHC of IFN-Gamma on a FFPE SARS-CoV2 Infected Lung Tissue

IFN- γ is a dimerized soluble cytokine and the only member of the type II interferon. In humans, IFN- γ is encoded by the IFNG gene and is critical for innate and adaptive immunity against viral, some bacterial and protozoal infections. IFN- γ is produced predominantly by natural killer and natural killer T cells as part of the innate immune response, and by CD4 Th1 and CD8 cytotoxic T lymphocytes effector T cells once antigen-specific immunity develops as part of the adaptive immune response. IFN- γ is also produced by non-cytotoxic innate lymphoid cells.

IFN- γ interacts with its receptors interferon gamma receptor 1 and interferon gamma receptor 2 and activates them. IFN- γ binding to the receptors activates the JAK-STAT pathway. IFN- γ also binds to the glycosaminoglycan heparan sulfate and inhibits its biological activity. IFN- γ is used to treat Chronic Granulomatous Disease and Osteopetrosis, also has potential as cancer immunotherapy to improve survival in Bladder Carcinoma, Melanoma, and Ovarian Carcinoma.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-161


ISOTYPE: IgG2a

CONTROL: Stomach, Fallopian Tube, Colon, Lung Adenocarcinoma, Ductal Breast Carcinoma, Pancreatic Adenocarcinoma

LOCALIZATION: Membranous, Cytoplasmic

SPECIES REACTIVITY: Human

IgA, RPab



IHC and IF of IgA on a FFPE Kidney Tissue (IHC) and on a Frozen Bullous Dermatitis Tissue (IF)

Immunoglobulin A (IgA) is the main immunoglobulin in mucous secretions, including tears, saliva, and colostrum, as well as respiratory, intestinal, prostatic, and vaginal secretions. It is also found in small amounts in blood. Because it is resistant to degradation by enzymes, secretory IgA provides protection against microbes proliferating in body secretions, especially those of the digestive and respiratory tracts.

IgA antibody reacts with surface immunoglobulin IgA alpha chains. It is extremely useful when identifying Acute Leukemias, IgA Myelomas, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. However, due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Tonsil, Spleen, Lymph Node, Kidney, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

IgA, MMab



IHC and IF of IgA on a FFPE Tonsil Tissue

Immunoglobulin A (IgA) is the main immunoglobulin in mucous secretions, including tears, saliva, and colostrum, as well as respiratory, intestinal, prostatic, and vaginal secretions. It is also found in small amounts in blood. Because it is resistant to degradation by enzymes, secretory IgA provides protection against microbes proliferating in body secretions, especially those of the digestive and respiratory tracts.

IgA antibody reacts with surface immunoglobulin IgA alpha chains. It is extremely useful when identifying Acute Leukemias, IgA Myelomas, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. However, due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-39

ISOTYPE: IgG1/K


CONTROL: Tonsil, Spleen, Lymph Node, Kidney, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB-3734-3 | Tinto Predilute | 3.0 ml | BSB 3054 | Tinto Predilute | 3.0 ml | BSB 5659 | Tinto Predilute | 3.0 ml |
| BSB-3734-7 | Tinto Predilute | 7.0 ml | BSB 3055 | Tinto Predilute | 7.0 ml | BSB 5660 | Tinto Predilute | 7.0 ml |
| BSB-3734-15 | Tinto Predilute | 15.0 ml | BSB 3056 | Tinto Predilute | 15.0 ml | BSB 5661 | Tinto Predilute | 15.0 ml |
| BSB-3734-01 | Concentrate | 0.1 ml | BSB 3057 | Concentrate | 0.1 ml | BSB 5662 | Concentrate | 0.1 ml |
| BSB-3734-05 | Concentrate | 0.5 ml | BSB 3058 | Concentrate | 0.5 ml | BSB 5663 | Concentrate | 0.5 ml |
| BSB-3734-1 | Concentrate | 1.0 ml | BSB 3059 | Concentrate | 1.0 ml | BSB 5664 | Concentrate | 1.0 ml |
| BSB-3734-CS | Control Slides | 5 | BSB 3685 | Control Slides | 5 | BSB 5665 | Control Slides | 5 |

IgD, RPaB




Inset: IHC and IF of IgD on a FFPE Colon Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

IgD makes up about 1% of proteins in the plasma membranes of immature B-lymphocytes (coexpressed with IgM) and is also found in serum in very small amounts. It is monomeric and incorporates the alpha-heavy chain in its structure. IgD's function is currently unknown, as mice lacking IgD seem to retain normal immune responses (implying redundancy if not lack of function), and IgD ceases to be expressed in activated B-lymphocytes. It may function as a regulatory antigen receptor. IgD is the major antigen receptor isotype co-expressed with IgM on the surface of most peripheral B cells in mice and humans.

The IgD antibody reacts with surface immunoglobulin IgD delta chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived from Lymphomas, specifically Marginal Zone Lymphoma. Renal involvement in systemic lupus erythematosus (SLE) is associated with production of antibodies to double stranded DNA, deposition of immune complexes and organ damage. Lupus nephritis patients were characterized by increased percentage of immature/early-transitional B-cells (CD27-IgD+CD21-), higher frequency of activated switched memory (SM, CD27+IgD-CD21-) and exhausted memory B-cells (CD27-IgD-), and decrease in non-switched memory (NSM, CD27+IgD+) B-cells

IgD, RMab




Inset: IHC and IF of IgD on a FFPE Tonsil Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

IgD makes up about 1% of proteins in the plasma membranes of immature B-lymphocytes (coexpressed with IgM) and is also found in serum in very small amounts. It is monomeric and incorporates the alpha-heavy chain in its structure. IgD's function is currently unknown, as mice lacking IgD seem to retain normal immune responses (implying redundancy if not lack of function), and IgD ceases to be expressed in activated B-lymphocytes. It may function as a regulatory antigen receptor.

IgD antibody reacts with surface immunoglobulin IgD delta chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived from Lymphomas, specifically Marginal Zone Lymphoma.

IgE, RPaB



Inset: IHC and IF of IgE on a FFPE Tonsil Tissue (IHC) and on a Frozen Colon Tissue (IF)

IgE, Immunoglobulin E, is an isotype of antibody only found in mammals. IgE is synthesized by plasma cells. Monomers of IgE consist of two heavy chains (ε chain) and two light chains, with the ε chain containing 4 Ig-like constant domains (Cε1-Cε4). IgE's main function is immunity to parasites such as helminths like Schistosoma mansoni, Trichinella spiralis, and Fasciola hepatica. IgE is utilized during immune defense against certain protozoan parasites such as Plasmodium falciparum.

IgE also has an essential role in type I hypersensitivity, which manifests in various allergic diseases, such as allergic asthma, most types of sinusitis, allergic rhinitis, food allergies, and specific types of chronic urticaria and atopic dermatitis. IgE also plays a pivotal role in responses to allergens, such as: anaphylactic drugs, bee stings, and antigen preparations used in desensitization immunotherapy. IgE is known to be elevated in various autoimmune disorders such as Lupus(SLE), Rheumatoid Arthritis(RA) & psoriasis, and is theorized to be of pathogenetic importance in RA and SLE by eliciting a hypersensitivity reaction.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP173
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Thymus, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5666 | Tinto Predilute | 3.0 ml |
| BSB 5667 | Tinto Predilute | 7.0 ml |
| BSB 5668 | Tinto Predilute | 15.0 ml |
| BSB 5669 | Concentrate | 0.1 ml |
| BSB 5670 | Concentrate | 0.5 ml |
| BSB 5671 | Concentrate | 1.0 ml |
| BSB 5672 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2957 | Tinto Predilute | 3.0 ml |
| BSB 2958 | Tinto Predilute | 7.0 ml |
| BSB 2959 | Tinto Predilute | 15.0 ml |
| BSB 2960 | Concentrate | 0.1 ml |
| BSB 2961 | Concentrate | 0.5 ml |
| BSB 2962 | Concentrate | 1.0 ml |
| BSB 2963 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3067 | Tinto Predilute | 3.0 ml |
| BSB 3068 | Tinto Predilute | 7.0 ml |
| BSB 3069 | Tinto Predilute | 15.0 ml |
| BSB 3070 | Concentrate | 0.1 ml |
| BSB 3071 | Concentrate | 0.5 ml |
| BSB 3072 | Concentrate | 1.0 ml |
| BSB 3073 | Control Slides | 5 |

IgG, RPab




IHC and IF of IgG on a FFPE Tonsil Tissue (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

IgG, MAb



IHC and IF of IgG on a FFPE Tonsil Tissue (IHC) and on a Frozen Tonsil (IF)

IgG4, MAb



IHC and IF of IgG4 on a FFPE Tonsil Tissue (IHC) and on a Frozen Tonsil (IF)

IgG is a monomeric immunoglobulin, comprised of two heavy chains and two light chains. This is the most abundant immunoglobulin and is approximately equally distributed in blood and tissue liquids, constituting 75% of serum immunoglobulins in humans. This is the only isotype that can pass through the placenta and bind to many kinds of pathogens. IgG protects the body against them by complement activation (classic pathway), opsonization for phagocytosis and neutralization of their toxins. There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%).

IgG antibody reacts with surface immunoglobulin IgG gamma chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q. Clinically, hematuria and proteinuria are present, with or without nephrotic syndromes. Mesangial IgG glomerulonephritis has been recently recognized as a distinct type of glomerulonephritis. The morphologic criteria detected in these patients included mesangial dense deposits by ultrastructural studies, which were predominantly positive for IgG by immunofluorescence.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Kidney, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

IgG is a monomeric immunoglobulin, composed of two heavy chains and two light chains. This is the most abundant immunoglobulin and is approximately equally distributed in blood and tissue liquids, constituting 75% of serum immunoglobulins in humans. This is the only isotype that can pass through the placenta and bind to many kinds of pathogens. IgG protects the body against them by complement activation (classic pathway), opsonization for phagocytosis and neutralization of their toxins. There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%).

IgG antibody reacts with surface immunoglobulin IgG gamma chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-40
ISOTYPE: IgG2a/K
CONTROL: Tonsil, Lymph Node, Kidney, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

IgG4-related sclerosing disease has been recognized as a systemic disease entity characterized by an elevated serum IgG4 level, sclerosing fibrosis and diffuse lymphoplasmacytic infiltration with the presence of many IgG4-positive plasma cells. As these patients tend to respond favorably to steroid treatment, it is important to recognize this entity and differentiate it from such mimics as lymphoma.

Clinical manifestations are apparent in the pancreas, bile duct, gallbladder, lacrimal gland, salivary gland, retroperitoneum, kidney, lung, breast, thyroid, and prostate. Immunohistochemical analyses in the case of IgG4-related sclerosing disease not only exhibits significantly more IgG4-positive plasma cells in affected tissues but also significantly higher IgG4/ IgG ratios (typically > 30%).


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-96
ISOTYPE: IgG2a/K
CONTROL: Tonsil, Spleen, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3074 | Tinto Predilute | 3.0 ml |
| BSB 3075 | Tinto Predilute | 7.0 ml |
| BSB 3076 | Tinto Predilute | 15.0 ml |
| BSB 3077 | Concentrate | 0.1 ml |
| BSB 3078 | Concentrate | 0.5 ml |
| BSB 3079 | Concentrate | 1.0 ml |
| BSB 5679 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5673 | Tinto Predilute | 3.0 ml |
| BSB 5674 | Tinto Predilute | 7.0 ml |
| BSB 5675 | Tinto Predilute | 15.0 ml |
| BSB 5676 | Concentrate | 0.1 ml |
| BSB 5677 | Concentrate | 0.5 ml |
| BSB 5678 | Concentrate | 1.0 ml |
| BSB 5679 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6807 | Tinto Predilute | 3.0 ml |
| BSB 6808 | Tinto Predilute | 7.0 ml |
| BSB 6809 | Tinto Predilute | 15.0 ml |
| BSB 6810 | Concentrate | 0.1 ml |
| BSB 6811 | Concentrate | 0.5 ml |
| BSB 6812 | Concentrate | 1.0 ml |
| BSB 6813 | Control Slides | 5 |

IgG4, RMab




IHC and IF of IgG4 on a FFPE Tonsil Tissue

IgG4-related sclerosing disease has been recognized as a systemic disease entity characterized by an elevated serum IgG4 level, sclerosing fibrosis and diffuse lymphoplasmacytic infiltration with the presence of many IgG4-positive plasma cells. As these patients tend to respond favorably to steroid treatment, it is important to recognize this entity and differentiate it from such mimics as lymphoma.

Clinical manifestations are apparent in the pancreas, bile duct, gallbladder, lacrimal gland, salivary gland, retroperitoneum, kidney, lung, breast, thyroid, and prostate. Immunohistochemical analyses in the case of IgG4-related sclerosing disease not only exhibits significantly more IgG4-positive plasma cells in affected tissues but also significantly higher IgG4/ IgG ratios (typically 30%).

IgM, RPab




IHC and IF of IgM on a FFPE Kidney Tissue (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

IgM forms polymers where multiple immunoglobulins are covalently linked together with disulfide bonds, normally as a pentamer or occasionally as a hexamer. It has a large molecular mass of approximately 900 kDa (in its pentamer form). In germline cells, the gene segment encoding the constant region of the heavy chain is positioned first among other constant region gene segments. For this reason, IgM is the first immunoglobulin expressed by mature B-cells.

IgM antibody reacts with surface immunoglobulin IgM mu chains. IgM is one of the predominant surface immunoglobulins on B-lymphocytes, and is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q. Clinically, hematuria and proteinuria are present, with or without nephrotic syndromes. Immunoglobulin M (IgM) nephropathy is an uncommon glomerular disease characterized by IgM deposits in the mesangium.

IgM, MAb



IHC and IF of IgM on a FFPE Kidney Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

IgM forms polymers where multiple immunoglobulins are covalently linked together with disulfide bonds, normally as a pentamer or occasionally as a hexamer. It has a large molecular mass of approximately 900 kDa (in its pentamer form). In germline cells, the gene segment encoding the constant region of the heavy chain is positioned first among other constant region gene segments. For this reason, IgM is the first immunoglobulin expressed by mature B-cells.

IgM antibody reacts with surface immunoglobulin IgM mu chains. IgM is one of the predominant surface immunoglobulins on B-lymphocytes, and is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP138

ISOTYPE: IgG

CONTROL: Tonsil, Spleen, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-41

ISOTYPE: IgG1/K

CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Colon

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6814 | Tinto Predilute | 3.0 ml |
| BSB 6815 | Tinto Predilute | 7.0 ml |
| BSB 6816 | Tinto Predilute | 15.0 ml |
| BSB 6817 | Concentrate | 0.1 ml |
| BSB 6818 | Concentrate | 0.5 ml |
| BSB 6819 | Concentrate | 1.0 ml |
| BSB 6820 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3080 | Tinto Predilute | 3.0 ml |
| BSB 3081 | Tinto Predilute | 7.0 ml |
| BSB 3082 | Tinto Predilute | 15.0 ml |
| BSB 3083 | Concentrate | 0.1 ml |
| BSB 3084 | Concentrate | 0.5 ml |
| BSB 3085 | Concentrate | 1.0 ml |
| BSB 5686 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5680 | Tinto Predilute | 3.0 ml |
| BSB 5681 | Tinto Predilute | 7.0 ml |
| BSB 5682 | Tinto Predilute | 15.0 ml |
| BSB 5683 | Concentrate | 0.1 ml |
| BSB 5684 | Concentrate | 0.5 ml |
| BSB 5685 | Concentrate | 1.0 ml |
| BSB 5686 | Control Slides | 5 |

IL-1a, MAb



IHC of IL-1a on a FFPE SARS-CoV2 Infected Lung Tissue

IL-1a is a dual-function cytokine, acting as a transcription factor and a cell-signaling molecule. The 31 kDa precursor form is biologically active in the cytosol and as a membrane protein. IL-1a can be cleaved by calpain into the biologically active 17 kDa mature form that regulates the expression of NFκB, and IFN-γ and its effectors. Cancer stem cells use IL-1a to maintain an inflammatory tumor microenvironment, however IL-1a is also used in recruiting and proliferating immune cells as a necrosis alarmin in damaged epithelia. IL-1a induces NFκB and STAT3 signaling that supports cancer stem cell growth. IL-1a is a major cytokine in innate inflammation signaling pathways, and is thus involved in many inflammatory diseases and reactions to oxidative stress.

IL-1a is involved in cancer cell cross-talk and metastasis, and the secretion of cell growth, angiogenesis, and inflammation factors as a Damage-Associated Molecular Pattern. Higher expression of IL-1a has been found in gastric cancer and squamous carcinomas of the head and neck that exhibit distant metastasis. IL-1a has also been associated with tumor dedifferentiation and metastasis in breast cancer, and with tumor proliferation and angiogenesis in pancreatic cancer.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-138


ISOTYPE: IgG2b

CONTROL: Kidney, Colon, Adrenal gland, Testis, Lung

LOCALIZATION: Membranous, Cytoplasmic

SPECIES REACTIVITY: Human

IL-1b, MAb



IHC of IL-1b on a FFPE SARS-CoV2 Infected Lung Tissue

IL-1b has multiple functions as a major pro-inflammatory cytokine regulated by NFκB. The 31 kDa precursor is biologically inactive in the cytosol, and is processed to its active 17 kDa mature form by Caspase 1 in special lysosomes. The beta-barrel cytokine is secreted to participate in pathways of angiogenesis, antigen presentation, adhesion molecule expression, inflammatory cell activity, and expression of matrix degrading enzymes.

IL-1b is produced by tumor-infiltrating myeloid cells, which then produce VEGF and other angiogenesis-related factors, promoting growth in multiple tumor types. IL-1b has a role in tumor metastasis and helps maintain an inflammatory tumor environment, but has also demonstrated an anti-tumorigenic function in recruiting immune cells and driving activation of the Th1 response against B cell myeloma and lymphoma. Pathogen-Associated Molecular Patterns can also initiate expression of IL-1b, helping to accumulate neutrophils and promote inflammation and fibrosis at the infection site.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-139


ISOTYPE: IgG2b

CONTROL: Colon, Pancreas, Liver, Stomach, Brain, Testis, Lung, Transitional Cell Carcinoma, Tonsil

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

IL-6, MAb



IHC of IL-6 on a FFPE SARS-CoV2 Infected Lung Tissue

Interleukin-6 is a cytokine that regulates inflammation pathways of the classic and trans-signaling cascades, often involved in metabolic, autoimmune and inflammatory conditions such as IgA nephropathy, lupus nephritis, and chronic kidney disease. The ligation of IL-6/IL-6R activates the JAK/STAT3 (cell survival and proliferation) and MAPK (cell growth and protein production) signaling pathways.

IL-6 participates in autoimmune and inflammatory conditions such as arthritis, where it triggers the differentiation of CD4+ T cells into Th17 helper cells, suppresses differentiation into T regulatory cells, and prevents T helper apoptosis along with IL-2. IL-6 also induces differentiation of monocytes into macrophages and uses follicular T helper cells to induce B cell activation, upregulating IgG production. IL-6 is often found in immune tissues, podocytes, mesangial cells, endothelial and epithelial cells. IL-6 is important in cancer cell growth and suppression of the immune system in the tumor microenvironment, metastasis and renewal of cancer stem cells.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-140

ISOTYPE: IgG2b

CONTROL: Testis, Lung, Stomach, Kidney, Transitional Cell Carcinoma

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3705-3 | Tinto Predilute | 3.0 ml |
| BSB-3705-7 | Tinto Predilute | 7.0 ml |
| BSB-3705-15 | Tinto Predilute | 15.0 ml |
| BSB-3705-01 | Concentrate | 0.1 ml |
| BSB-3705-05 | Concentrate | 0.5 ml |
| BSB-3705-1 | Concentrate | 1.0 ml |
| BSB-3705-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3706-3 | Tinto Predilute | 3.0 ml |
| BSB-3706-7 | Tinto Predilute | 7.0 ml |
| BSB-3706-15 | Tinto Predilute | 15.0 ml |
| BSB-3706-01 | Concentrate | 0.1 ml |
| BSB-3706-05 | Concentrate | 0.5 ml |
| BSB-3706-1 | Concentrate | 1.0 ml |
| BSB-3706-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3707-3 | Tinto Predilute | 3.0 ml |
| BSB-3707-7 | Tinto Predilute | 7.0 ml |
| BSB-3707-15 | Tinto Predilute | 15.0 ml |
| BSB-3707-01 | Concentrate | 0.1 ml |
| BSB-3707-05 | Concentrate | 0.5 ml |
| BSB-3707-1 | Concentrate | 1.0 ml |
| BSB-3707-CS | Control Slides | 5 |

IMP-3/IGF2BP3, RMAb




IHC of IMP-3 on a FFPE Placenta Tissue

Inhibin alpha, MAb



IHC of Inhibin alpha on a FFPE Testis Tissue

Inhibin alpha, RMAb



IHC of Inhibin alpha on a FFPE Testis Tissue

Insulin-like growth factor 2 mRNA-binding protein 3 is a protein that in humans is encoded by the IGF2BP3 gene. IMP3 is normally expressed in early embryonic tissues. The IHC of IMP3 may help in the classification of Non-small Cell Lung Carcinomas and Pancreatic Adenocarcinomas as well as subtypes of carcinomas from other organs such as Renal Cell Carcinoma, Adenocarcinoma of the Uterine Cervix, Endometrial Carcinoma, Adenocarcinoma of the Esophagus, Malignant Melanoma, Merkel Cell Carcinoma, Urothelial Carcinoma, Neuroendocrine Carcinoma of the Lung, and triple negative breast cancer.

IMP3 expression for the prognostic evaluation of non-small cell lung carcinomas has been found to exhibit mainly cytoplasmic staining pattern in the NSCLC tissues with a positive rate of IMP3 protein expression of 74.7% in the NSCLC tissues, a significantly higher than the rate of 19.9% in the adjacent non-tumor tissues. IMP3 may be a useful diagnostic marker in the assessment of endometrial cancers and their precursor lesions, particularly when the amount of available tissue material is limited and a concern of type II cancer arises. High frequency of IMP3 expression is present in decidualized endometrial stroma of gestational endometrium and chorionic villi in early pregnancy.

Inhibins are peptide hormones produced by the granulosa cells in female follicles and by Sertoli cells in the male seminiferous tubules. They are selectively expressed by cells of sex-cord stromal derivation, and inhibit the secretion of follitropin by the pituitary gland. Inhibin contains an alpha and beta subunit linked by disulfide bonds. Two forms of inhibin differ in their beta subunits (A or B), while their alpha subunits are identical. Inhibin belongs to the transforming growth factor-beta (TGF-beta) family. Anti-Inhibin Alpha has demonstrated utility in differentiation between Adrenal Cortical Tumors and Renal Cell Carcinoma. Sex-Cord Stromal Tumors of the Ovary as well as Trophoblastic Tumors also demonstrate cytoplasmic positivity with this antibody.

Inhibins are peptide hormones produced by the granulosa cells in female follicles and by Sertoli cells in the male seminiferous tubules. They are selectively expressed by cells of sex-cord stromal derivation, and inhibit the secretion of follitropin by the pituitary gland. Inhibin contains an alpha and beta subunit linked by disulfide bonds. Two forms of inhibin differ in their beta subunits (A or B), while their alpha subunits are identical. Inhibin belongs to the transforming growth factor-beta (TGF-beta) family.

Anti-Inhibin Alpha has demonstrated utility in differentiation between Adrenal Cortical Tumors and Renal Cell Carcinoma. Sex-Cord Stromal Tumors of the Ovary as well as Trophoblastic Tumors also demonstrate cytoplasmic positivity with this antibody.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP286
ISOTYPE: IgG
CONTROL: Placenta, Tonsil, Ovarian Carcinoma, TCC
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: R1
ISOTYPE: IgG2a
CONTROL: Adrenal Cortex, Placenta, Testis, Corpus Luteum
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP378
ISOTYPE: IgG
CONTROL: Testis, Seminoma, Testicular Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2964 | Tinto Predilute | 3.0 ml |
| BSB 2965 | Tinto Predilute | 7.0 ml |
| BSB 2966 | Tinto Predilute | 15.0 ml |
| BSB 2967 | Concentrate | 0.1 ml |
| BSB 2968 | Concentrate | 0.5 ml |
| BSB 2969 | Concentrate | 1.0 ml |
| BSB 2970 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5687 | Tinto Predilute | 3.0 ml |
| BSB 5688 | Tinto Predilute | 7.0 ml |
| BSB 5689 | Tinto Predilute | 15.0 ml |
| BSB 5690 | Concentrate | 0.1 ml |
| BSB 5691 | Concentrate | 0.5 ml |
| BSB 5692 | Concentrate | 1.0 ml |
| BSB 5693 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3356 | Tinto Predilute | 3.0 ml |
| BSB 3357 | Tinto Predilute | 7.0 ml |
| BSB 3358 | Tinto Predilute | 15.0 ml |
| BSB 3559 | Concentrate | 0.1 ml |
| BSB 3560 | Concentrate | 0.5 ml |
| BSB 3561 | Concentrate | 1.0 ml |
| BSB 3562 | Control Slides | 5 |

INI-1, RMAb

IHC of INI-1 on a FFPE Ewing's Sarcoma Tissue

INI-1, RMAb

IHC of INI-1 on a FFPE Ewing's Sarcoma Tissue

iNOS, RMAb

IHC of iNOS on a FFPE Liver Tissue

The INI-1 gene, which encodes a functionally uncharacterized protein component of the hSWI/SNF chromatin remodeling complex, is often mutated or deleted in malignant rhabdoid tumor (MRT). Two isoforms of INI-1, that differ by the variable inclusion of amino acids, potentially are produced by differential RNA splicing.

The morphology of MRTs can present challenges in differential diagnosis. The overall survival of MRTs relative to its potential mimics (medulloblastoma, supratentorial primitive neuroectodermal tumors (sPNETs)) is quite low, and thus differentiation from these other tumors is desirable. Lack of nuclear labeling by anti-INI-1 is characteristic of MRT. The majority of medulloblastomas and sPNETs are labeled by anti-INI-1. MRTs also originate from the kidney and soft tissues.

The INI-1 gene, which encodes a functionally uncharacterized protein component of the hSWI/SNF chromatin remodeling complex, is involved in chromatin remodeling and transcriptional regulation and also known as hSNF5, SMARCB1 and BAF47.

INI-1 is often mutated or deleted in malignant rhabdoid tumor (MRT). Two isoforms of INI-1, that differ by the variable inclusion of amino acids, potentially are produced by differential RNA splicing. The morphology of MRTs can present challenges in differential diagnosis. The overall survival of MRTs relative to its potential mimics such as medulloblastoma and supratentorial primitive neuroectodermal tumors (sPNETs) is quite low, and thus differentiation from these other tumors is desirable. Lack of nuclear labeling by INI-1 is characteristic of MRT. The majority of medulloblastomas and sPNETs are labeled by INI-1. MRTs also originate from the kidney and soft tissues. Germline INI-1 mutations are associated with sporadic schwannomatosis and rhabdoid tumors.

INI-1 has been used to distinguish atypical teratoid / rhabdoid tumor (loss of INI1) from choroid plexus carcinoma (positive for INI-1) and to differentiate epithelioid sarcoma (loss of INI1) from epithelioid hemangioendothelioma (positive for INI-1).

Nitric oxide synthases (NOSs) are a family of enzymes catalyzing the production of nitric oxide (NO) from L-arginine. NO is an important cellular signaling molecule. It helps modulate vascular tone, insulin secretion, airway tone, and peristalsis, and is involved in angiogenesis and neural development. The inducible isoform, iNOS, is involved in immune response, binds calmodulin at physiologically relevant concentrations, and produces NO as an immune defense mechanism, as NO is a free radical with an unpaired electron. It is the proximate cause of septic shock and may function in autoimmune disease.

The inducible isoform iNOS produces large amounts of NO as a defense mechanism. It is synthesized by many cell types in response to cytokines and is an important factor in the response of the body to attack by parasites, bacterial infection, and tumor growth. It is also the cause of septic shock and may play a role in many diseases with an autoimmune etiology. Induction of the high-output iNOS usually occurs in an oxidative environment, and thus high levels of NO have the opportunity to react with superoxide leading to peroxynitrite formation and cell toxicity. These properties may define the roles of iNOS in host immunity, enabling its participation in anti-microbial and anti-tumor activities as part of the oxidative burst of macrophages.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: 25

ISOTYPE: IgG2a

CONTROL: Testis, Brain, Breast, Colon, Kidney, Pituitary Adrenal, Prostate, Thyroid, Lung, Pancreas, Cervix, Salivary Gland, Astrocytoma, Lymphoblastic Lymphoma, Transitional Cell Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-INI1

ISOTYPE: IgG

CONTROL: Testis, Brain, Breast, Colon

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-iNOS

ISOTYPE: IgG

CONTROL: Testis, Adrenal, Lung, Prostate, Liver, Placenta, Spleen

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6828 | Tinto Predilute | 3.0 ml |
| BSB 6829 | Tinto Predilute | 7.0 ml |
| BSB 6830 | Tinto Predilute | 15.0 ml |
| BSB 6831 | Concentrate | 0.1 ml |
| BSB 6832 | Concentrate | 0.5 ml |
| BSB 6833 | Concentrate | 1.0 ml |
| BSB 6834 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3779-3 | Tinto Predilute | 3.0 ml |
| BSB-3779-7 | Tinto Predilute | 7.0 ml |
| BSB-3779-15 | Tinto Predilute | 15.0 ml |
| BSB-3779-01 | Concentrate | 0.1 ml |
| BSB-3779-05 | Concentrate | 0.5 ml |
| BSB-3779-1 | Concentrate | 1.0 ml |
| BSB-3779-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2656 | Tinto Predilute | 3.0 ml |
| BSB 2657 | Tinto Predilute | 7.0 ml |
| BSB 2658 | Tinto Predilute | 15.0 ml |
| BSB 2659 | Concentrate | 0.1 ml |
| BSB 2660 | Concentrate | 0.5 ml |
| BSB 2661 | Concentrate | 1.0 ml |
| BSB 2662 | Control Slides | 5 |

INSM1, MAb



IHC of INSM1 on a FFPE Lung Neuroendocrine Carcinoma Tissue

INSM1, RMAb



IHC of INSM1 on a FFPE Lung Neuroendocrine Carcinoma Tissue

Insulin, MAb



IHC of Insulin on a FFPE Pancreas Tissue

Insulinoma-associated 1 (INSM1) gene encodes a protein containing both a zinc finger DNA-binding domain and a putative prohormone domain, originally isolated from a human insulinoma-glucagonoma subtraction library. INSM1 mRNA is abundantly expressed in fetal NE developmental tissues and expressed in normal adult neuroendocrine tissues (adrenal medulla, pineal gland, pituitary gland, gastrointestinal enterochromaffin cells, pancreatic islet cells, thyroid C cells) and developing neurons, however there is also a high occurrence of INSM1 found in NE tumors, such as Small Cell Lung Cancer (SCLC), Pituitary tumors, Medullary Thyroid Carcinoma, Merkel Cell Carcinoma, Olfactory Neuroblastoma and Pheochromocytoma.

It has been reported that INSM1 expresses exclusively in SCLC specimens using immunohistochemistry, and first elucidated that INSM1 regulates the NE differentiation pathway in lung cancer. In addition, it has demonstrated an increased sensitivity and specificity compared to other NE biomarkers (Chromogranin A, Synaptophysin and CD56) in lung cancer specimens. In addition, it's been shown to be involved in NE differentiation in medullary thyroid carcinoma, pheochromocytoma, intestinal NE carcinoma, islet cell tumor, pituitary tumor, and SCLC cell lines.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-123

ISOTYPE: IgG1/k

CONTROL: Pancreas, Colon, Tonsil, Neuroendocrine Lung Cancer, Endometrial & Colon Carcinomas

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

Insulinoma-associated 1 (INSM1) gene encodes a protein containing both a zinc finger DNA-binding domain and a putative prohormone domain, originally isolated from a human insulinoma-glucagonoma subtraction library.

INSM1 is abundantly expressed in fetal neuroendocrine developmental tissues and expressed in normal adult neuroendocrine (NE) tissues (adrenal medulla, pineal gland, pituitary gland, gastrointestinal enterochromaffin cells, pancreatic islet cells, thyroid C cells) and developing neurons. However there is also a high occurrence of INSM1 found in NE tumors, such as small cell lung cancer (SCLC), pituitary tumors, medullary thyroid carcinoma, Merkel cell carcinoma, olfactory neuroblastoma and pheochromocytoma. It has been reported that INSM1 expresses exclusively in SCLC specimens using immunohistochemistry, and first elucidated that INSM1 regulates the NE differentiation pathway in lung cancer. In addition, it has demonstrated an increased sensitivity and specificity compared to other NE biomarkers (Chromogranin A, Synaptophysin and CD56) in lung cancer specimens. In addition, it's been shown to be involved in NE differentiation in medullary thyroid carcinoma, pheochromocytoma, intestinal NE carcinoma, islet cell tumor, pituitary tumor, and SCLC cell lines.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-INSM1

ISOTYPE: IgG

CONTROL: Pancreas, Colon, Tonsil, Neuroendocrine Lung Cancer, Endometrial & Colon Carcinomas

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Insulin is produced in the beta cells of the Islets of Langerhans in the pancreas. It is a polypeptide hormone that regulates carbohydrate metabolism. Apart from being the primary agent in carbohydrate homeostasis, insulin has effects on fat metabolism and changes the liver's ability in storing or releasing glucose and processing blood lipids, and in other tissues such as fat and muscle. The amount of insulin in circulation has extremely widespread effects throughout the body.

The presence of insulin in the cytoplasm of Islet Tumors is the most reliable indication of functional Insulinomas. Defective insulin storage occurs in Insulinomas; therefore, many sections of the tumor should be stained with both insulin and C-peptide.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-42

ISOTYPE: IgG1/K

CONTROL: Pancreas

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Dog, Cat, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3553 | Tinto Predilute | 3.0 ml |
| BSB 3554 | Tinto Predilute | 7.0 ml |
| BSB 3555 | Tinto Predilute | 15.0 ml |
| BSB 3556 | Concentrate | 0.1 ml |
| BSB 3557 | Concentrate | 0.5 ml |
| BSB 3558 | Concentrate | 1.0 ml |
| BSB 3559 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3780-3 | Tinto Predilute | 3.0 ml |
| BSB-3780-7 | Tinto Predilute | 7.0 ml |
| BSB-3780-15 | Tinto Predilute | 15.0 ml |
| BSB-3780-01 | Concentrate | 0.1 ml |
| BSB-3780-05 | Concentrate | 0.5 ml |
| BSB-3780-1 | Concentrate | 1.0 ml |
| BSB-3780-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5694 | Tinto Predilute | 3.0 ml |
| BSB 5695 | Tinto Predilute | 7.0 ml |
| BSB 5696 | Tinto Predilute | 15.0 ml |
| BSB 5697 | Concentrate | 0.1 ml |
| BSB 5698 | Concentrate | 0.5 ml |
| BSB 5699 | Concentrate | 1.0 ml |
| BSB 5700 | Control Slides | 5 |

Insulin, RMAb




IHC of Insulin on a FFPE Pancreas Tissue

Insulin is produced in the beta cells of the Islets of Langerhans in the pancreas. It is a polypeptide hormone that regulates carbohydrate metabolism. Apart from being the primary agent in carbohydrate homeostasis, insulin has effects on fat metabolism and changes the liver's ability in storing or releasing glucose and processing blood lipids, and in other tissues such as fat and muscle. The amount of insulin in circulation has extremely widespread effects throughout the body.

The presence of insulin in the cytoplasm of Islet Tumors is the most reliable indication of functional Insulinomas. Defective insulin storage occurs in Insulinomas; therefore, many sections of the tumor should be stained with both insulin and C-peptide.

Islet 1/ISL1, RMAb




IHC of Islet 1 on a FFPE Pancreas Tissue

Islet-1, Insulin gene enhancer protein ISL-1, is a protein that in humans is encoded by the ISL1 gene. This gene encodes a transcription factor containing two N-terminal LIM domains and one C-terminal homeo domain. The encoded protein plays an important role in the embryogenesis of pancreatic islets of Langerhans. ISL1 has been shown to interact with Estrogen Receptor alpha.

Islet-1 produces a strong nuclear staining in the islets of normal pancreas and tumor cells of the pancreatic neuroendocrine tumors. Islet-1 has been found to be a reliable marker of pancreatic endocrine tumors and metastasis. It shows a comparable sensitivity and specificity as CDX2 as a marker of ileal and appendiceal neuroendocrine tumors and their metastasis. TTF1 is very rarely expressed in well-differentiated gastroenteropancreatic endocrine tumors. Therefore, the panel of Islet-1 CDX2, and TTF1 and they may be useful for examining metastasis of well-differentiated endocrine carcinomas of unknown origin.

Kappa, RPAb



IHC and IF of Kappa on a FFPE Kidney Tissue (IHC) and on a Frozen Colon Tissue (IF with FluoroMounter DAPI)

Kappa detects surface immunoglobulin on normal and neoplastic B-cells. In paraffin-embedded tissue, Kappa exhibits strong staining of kappa-positive plasma cells and cells that have absorbed exogenous immunoglobulin.

When studying B-cell neoplasms, the determination of light-chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either Kappa or Lambda light chains, whereas reactive proliferations display a mixture of Kappa and Lambda-positive cells. If only a single light-chain type is detected, a lympho-proliferative disorder is very likely. Monoclonality is determined by a Kappa-Lambda ratio greater than or equal to 3:1, a Lambda-Kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population.

In IgG-dominant immune complex-mediated glomerulonephritis, there are multiple pathological findings that strongly suggest the diagnosis of Lupus Nephritis including immunofluorescence staining for IgG, IgM, IgA, Kappa or Lambda, C3 and C1.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP125

ISOTYPE: IgG

CONTROL: Pancreas

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP283

ISOTYPE: IgG

CONTROL: Pancreas, Testis, Thyroid, Cervix, Skin, Pancreatic Neuroendocrine Cancer

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6835 | Tinto Predilute | 3.0 ml |
| BSB 6836 | Tinto Predilute | 7.0 ml |
| BSB 6837 | Tinto Predilute | 15.0 ml |
| BSB 6838 | Concentrate | 0.1 ml |
| BSB 6839 | Concentrate | 0.5 ml |
| BSB 6840 | Concentrate | 1.0 ml |
| BSB 6841 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2971 | Tinto Predilute | 3.0 ml |
| BSB 2972 | Tinto Predilute | 7.0 ml |
| BSB 2973 | Tinto Predilute | 15.0 ml |
| BSB 2974 | Concentrate | 0.1 ml |
| BSB 2975 | Concentrate | 0.5 ml |
| BSB 2976 | Concentrate | 1.0 ml |
| BSB 2977 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3086 | Tinto Predilute | 3.0 ml |
| BSB 3087 | Tinto Predilute | 7.0 ml |
| BSB 3088 | Tinto Predilute | 15.0 ml |
| BSB 3089 | Concentrate | 0.1 ml |
| BSB 3090 | Concentrate | 0.5 ml |
| BSB 3091 | Concentrate | 1.0 ml |
| BSB 3688 | Control Slides | 5 |

Kappa Light Chains, MMab




IHC and IF of Kappa on a FFPE Tonsil Tissue (IHC) and on a Frozen Spleen Tissue (IF)

Kappa detects surface immunoglobulin on normal and neoplastic B-cells. In paraffin-embedded tissue, Kappa exhibits strong staining of kappa-positive plasma cells and cells that have absorbed exogenous immunoglobulin.

When studying B-cell neoplasms, the determination of light-chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-58
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

Ki-67, RMAb




IHC of Ki-67 on a FFPE Colon Tissue

The Ki-67 protein is a cellular marker for proliferation. It is strictly associated with cell proliferation. During the interphase, the Ki-67 antigen can be exclusively detected within the cell nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes. Ki-67 protein is present during all active phases of the cell cycle (G1, S, G2, and mitosis), but is absent from resting cells (G0).

Ki-67 is an excellent marker to determine the growth fraction of a given cell population. The fraction of Ki-67-positive tumor cells (the Ki-67 labeling index) is often correlated with the clinical course of cancer. The best-studied examples in this context are Carcinomas of the Prostate and the Breast.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP5
ISOTYPE: IgG
CONTROL: Testis, Tonsil, Bone Marrow, Placenta, Colon, Tonsil, Fallopian Tube, Astrocytoma, Breast Carcinoma, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Ki-67, RMAb



IHC of Ki-67 on FFPE Breast Carcinoma Tissue

The Ki-67 protein is a cellular marker for proliferation. It is strictly associated with cell proliferation. During the interphase, the Ki-67 antigen can be exclusively detected within the cell nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes. Ki-67 protein is present during all active phases of the cell cycle (G1, S, G2, and mitosis), but is absent from resting cells (G0).

Ki-67 is an excellent marker to determine the growth fraction of a given cell population. The fraction of Ki-67-positive tumor cells (the Ki-67 labeling index) is often correlated with the clinical course of cancer. The best-studied examples in this context are Carcinomas of the Prostate and the Breast.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM360
ISOTYPE: IgG
CONTROL: Testis, Tonsil, Bone Marrow, Placenta, Colon, Tonsil, Fallopian Tube, Astrocytoma, Breast Carcinoma, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5701 | Tinto Predilute | 3.0 ml |
| BSB 5702 | Tinto Predilute | 7.0 ml |
| BSB 5703 | Tinto Predilute | 15.0 ml |
| BSB 5704 | Concentrate | 0.1 ml |
| BSB 5705 | Concentrate | 0.5 ml |
| BSB 5706 | Concentrate | 1.0 ml |
| BSB 5707 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5708 | Tinto Predilute | 3.0 ml |
| BSB 5709 | Tinto Predilute | 7.0 ml |
| BSB 5710 | Tinto Predilute | 15.0 ml |
| BSB 5711 | Concentrate | 0.1 ml |
| BSB 5712 | Concentrate | 0.5 ml |
| BSB 5713 | Concentrate | 1.0 ml |
| BSB 5714 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3767-3 | Tinto Predilute | 3.0 ml |
| BSB-3767-7 | Tinto Predilute | 7.0 ml |
| BSB-3767-15 | Tinto Predilute | 15.0 ml |
| BSB-3767-01 | Concentrate | 0.1 ml |
| BSB-3767-05 | Concentrate | 0.5 ml |
| BSB-3767-1 | Concentrate | 1.0 ml |
| BSB-3767-CS | Control Slides | 5 |

Ksp-Cadherin, MMab




IHC of Ksp-Cadherin on a FFPE Kidney Tissue

LAG-3/CD223. RMAb



IHC of LAG-3 / CD223 on a FFPE Testis Tissue

Lambda, RPAb



IHC and IF of Lambda on a FFPE Tonsil Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

Ksp-Cadherin (Kidney-specific Cadherin) is a novel, kidney-specific member of the Cadherin family of cell-adhesion molecules. Within the kidney, Ksp-Cadherin is found exclusively in the basolateral membrane of renal tubular epithelial cells and collecting duct cells, and not in glomeruli, renal interstitial cells, or blood vessels. Different Cadherins, including E-Cadherin, Cadherin-6, and N-Cadherin, have been investigated in Renal Cell Cancers, demonstrating possible correlations of tumor differentiation and the presence of lymph node metastasis with loss of Cadherins.

Ksp-Cadherin has been used to distinguish Chromophobe Renal-Cell Carcinoma from Oncocytoma. Studies have found a membranous pattern of staining in 96% of 30 Chromophobe Carcinomas, and in only 6% of 31 Oncocytomas, leading to conclude that this is a useful antibody in differentiating these two lesions. On the other hand, another study found Ksp-Cadherin positivity in 100% of 13 chromophobe RCCs, and 95% of 20 Oncocytomas.

Lymphocyte-activation gene 3, also known as LAG-3 and designated as CD223, is a cell surface protein with diverse biologic effects on T cell function encoded by the LAG3 gene. The gene for LAG-3 lies adjacent to the gene for CD4 on human chromosome 12 (12p13) and is approximately 20% identical to the CD4 gene. LAG-3 is expressed on activated T cells, natural killer cells, B cells and plasmacytoid dendritic cells. The LAG-3 protein negatively regulates cellular proliferation, activation, and homeostasis of T cells, in a similar fashion to CTLA-4 and PD-1 and has been reported to play a role in Treg, regulatory T cells, suppressive function. LAG-3 also helps maintain CD8+ T cells in a tolerogenic state and, working with PD-1, helps maintain CD8 exhaustion during chronic viral infection. LAG-3 is known to be involved in the maturation and activation of dendritic cells.

Studies have shown that LAG-3 is a prognostic indicator of poor treatment outcomes in chronic lymphocytic leukemia, has been associated with higher risk of multiple myeloma, may play a modulating role in autoimmune diabetes and has been identified in a subset of HIV-specific LAG-3(+) CD8(+) T cells that negatively correlated with plasma viral load. IHC studies have shown LAG-3 expression on lymphocytes scattered in renal cell carcinoma, melanoma and lymphomas.

Lambda detects surface immunoglobulin on normal and neoplastic B-cells. Lambda staining is seen in B-cell follicles of human lymphoid tissue.

When studying B-cell neoplasms, the determination of light chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population. In IgG-dominant immune complex-mediated glomerulonephritis, there are multiple pathological findings that strongly suggest the diagnosis of Lupus Nephritis including immunofluorescence staining for IgG, IgM, IgA, Kappa or Lambda, C3 and C1.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 4H6/F9

ISOTYPE: IgG1

CONTROL: Kidney, Chromophobe Renal Cell Carcinoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP294

ISOTYPE: IgG

CONTROL: Testis, Lymph Node, Spleen, Lymphoblastic Lymphoma, TCC

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6282 | Tinto Predilute | 3.0 ml |
| BSB 6283 | Tinto Predilute | 7.0 ml |
| BSB 6284 | Tinto Predilute | 15.0 ml |
| BSB 6285 | Concentrate | 0.1 ml |
| BSB 6286 | Concentrate | 0.5 ml |
| BSB 6287 | Concentrate | 1.0 ml |
| BSB 6288 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3363 | Tinto Predilute | 3.0 ml |
| BSB 3364 | Tinto Predilute | 7.0 ml |
| BSB 3365 | Tinto Predilute | 15.0 ml |
| BSB 3366 | Concentrate | 0.1 ml |
| BSB 3367 | Concentrate | 0.5 ml |
| BSB 3368 | Concentrate | 1.0 ml |
| BSB 3369 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3092 | Tinto Predilute | 3.0 ml |
| BSB 3093 | Tinto Predilute | 7.0 ml |
| BSB 3094 | Tinto Predilute | 15.0 ml |
| BSB 3095 | Concentrate | 0.1 ml |
| BSB 3096 | Concentrate | 0.5 ml |
| BSB 3097 | Concentrate | 1.0 ml |
| BSB 3089 | Control Slides | 5 |

Lambda, MAb




IHC and IF of Lambda on a FFPE Tonsil Tissue

Lambda detects surface immunoglobulin on normal and neoplastic B-cells. Lambda staining is seen in B-cell follicles of human lymphoid tissue.

When studying B-cell neoplasms, the determination of light chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-16
ISOTYPE: IgG2a
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

Lamin-B1, RMAb




IHC of Lamin-B1 on a FFPE Colon Adenocarcinoma Tissue

Lamin-B1 is an intermediate filament protein that in humans is encoded by the gene LMNB1. B-type Lamin along with A-type Lamin line the inner surface of the nuclear envelope and form a dynamic structure that is disassembled and reassembled each time a cell enters mitosis. Lamin-B1 plays a central role in chromatin organization, gene positioning, DNA replication and repair, cell cycle progression, stress responses, proliferation and differentiation.

Lamin-B1 loss is associated with cellular senescence and a broad range of aging-related diseases such as cardiovascular diseases and cancers. Recent studies showed Lamin-B1 deficiency promotes Lung Cancer development and metastasis by epigenetic derepression of RET, as well as malignancy and poor prognosis in Gastric Cancer. It has been observed that Lamin B1 levels are reduced in Lung Cancer patients compared to normal Lung tissue and that lower expression of Lamin B1 was associated with higher tumor grade. The levels of lamin A, on the other hand, were not altered in human Lung Tumors, indicating that the different types of lamins have distinct functions in lung carcinogenesis. Overexpression of Lamin-B1 can also indicate poor prognosis in Colon and Pancreatic Cancers.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-LMNB1
ISOTYPE: IgG
CONTROL: Colon, Breast, Fallopian Tube, Tonsil, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma, Colon Adenocarcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

Laminin-R/RPSA, MAb



IHC of Laminin-R/RPSA on a FFPE Colon Carcinoma Tissue

Laminin receptor (Ribosomal Protein SA or 67LR) is a 67kDa protein in the extracellular matrix, coded by a 5,833 bp gene on Chromosome 3, containing seven exons and six introns. Laminin-R is involved in cellular adhesion to basement membranes and signal transduction pathways after binding the ligand Integrin-1. The receptor is also involved in pathways of cytoskeletal movement, chromatin structure, cellular migration, stress response, ribosomal function, translation, and rRNA processing. Laminin-R is also a receptor for green tea polyphenol and mediator in the Nitrous Oxide (eNOS) pathways used to force cancer cells into apoptosis. The 37kDa precursor protein, known also as p40, is a highly-conserved protein in the 40S ribosomal subunit. Laminin-R has been found to be overexpressed in Breast, Colorectal, Pancreatic, Prostate, and Cervical Cancer, and in Lymphomas. Laminin-R interacts with cyclin-dependent kinases and inhibitors to control G1/S, S/G2 and G2/M phase cell cycle arrest; down-regulation of the receptor is correlated with arrested cell growth and a reduction in the development of tumors. Laminin-R has also been found to be involved in regulation of Survivin expression, as well as cancer metastasis and drug resistance.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-144
ISOTYPE: IgG1
CONTROL: Placenta, Kidney, Prostate, Tonsil, Spleen, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5715 | Tinto Predilute | 3.0 ml |
| BSB 5716 | Tinto Predilute | 7.0 ml |
| BSB 5717 | Tinto Predilute | 15.0 ml |
| BSB 5718 | Concentrate | 0.1 ml |
| BSB 5719 | Concentrate | 0.5 ml |
| BSB 5720 | Concentrate | 1.0 ml |
| BSB 5721 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3735-3 | Tinto Predilute | 3.0 ml |
| BSB-3735-7 | Tinto Predilute | 7.0 ml |
| BSB-3735-15 | Tinto Predilute | 15.0 ml |
| BSB-3735-01 | Concentrate | 0.1 ml |
| BSB-3735-05 | Concentrate | 0.5 ml |
| BSB-3735-1 | Concentrate | 1.0 ml |
| BSB-3735-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3736-3 | Tinto Predilute | 3.0 ml |
| BSB-3736-7 | Tinto Predilute | 7.0 ml |
| BSB-3736-15 | Tinto Predilute | 15.0 ml |
| BSB-3736-01 | Concentrate | 0.1 ml |
| BSB-3736-05 | Concentrate | 0.5 ml |
| BSB-3736-1 | Concentrate | 1.0 ml |
| BSB-3736-CS | Control Slides | 5 |

Langerin/CD207, MAb




IHC of Langerin on a FFPE Skin Tissue

Langerin/CD207, RMab



IHC of Langerin on a FFPE Prostate Tissue

LEF-1, RMab



IHC of LEF-1 on a FFPE Testicular Carcinoma Tissue

Langerin is a type II transmembrane cell surface receptor produced by Langerhans Cells, which are immature dendritic cells of the epidermis and mucosa. Epidermal LCs possess strong immunohistochemistry capacity and play a central role in the initiation and regulation of immune responses. Langerin is localized in the Birbeck granules, organelles present in the cytoplasm of Langerhans cells and consisting of superimposed and zippered membranes. It is a C-type lectin with mannose binding specificity, and it has been proposed that mannose binding by this protein leads to internalization of antigen into Birbeck granules and provides access to a nonclassical antigen-processing pathway.

Human spleen, lymph node, thymus, liver, lung, and heart express langerin protein. Langerin protein expression has utility in differentiating Langerhans cell histiocytosis from other non-Langerhans cell histiocytic proliferations.

Langerin is a type II transmembrane cell surface receptor produced by Langerhans Cells, which are immature dendritic cells of the epidermis and mucosa. Epidermal LCs possess strong immunohistochemistry capacity and play a central role in the initiation and regulation of immune responses. Langerin is localized in the Birbeck granules, organelles present in the cytoplasm of Langerhans cells and consisting of superimposed and zippered membranes. It is a C-type lectin with mannose binding specificity, and it has been proposed that mannose binding by this protein leads to internalization of antigen into Birbeck granules and providing access to a nonclassical antigen-processing pathway.

Human spleen, lymph node, thymus, liver, lung and heart express langerin protein. Langerin protein expression has utility in differentiating Langerhans cell histiocytosis from other non-Langerhans cell histiocytic proliferations.

LEF1 is highly overexpressed and associated with disease progression and poor prognosis in B-cell chronic lymphocytic leukemia. Strong nuclear expression of LEF1 has been observed in majority of chronic lymphocytic leukemia/small lymphocytic lymphoma cases and LEF1 is not detected in other small B cell lymphomas. Gene expression profiling revealed overexpression of LEF-1 in chronic lymphocytic leukemia (CLL)/small lymphocytic lymphoma (SLL). LEF-1 immunostaining has been detected in all neoplastic cells of CLL/SLL cases. LEF-1 was identified in 50% of high grade follicular lymphoma and 38% of diffuse large B-cell lymphoma, but not in mantle cell lymphoma or marginal zone lymphoma. Recently, high LEF-1 was demonstrated as a favorable prognostic marker in cytogenetically normal acute myeloid leukemia. Due to its high sensitivity, LEF-1 has been proposed to be a suitable immunohistochemical marker for diagnosis and differential diagnosis for CLL/SLL.

Alternately spliced isoforms may play additional roles in regulating cell growth in colon carcinoma, and nuclear LEF-1 immunostaining was detected in 36% of adenocarcinoma brain metastases.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 12D6

ISOTYPE: IgG2b/K

CONTROL: Skin, Breast, Prostate, Cervix, Liver, Salivary Gland, Langerhans Histiocytosis

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP349

ISOTYPE: IgG

CONTROL: Skin, Prostate, Breast, Liver, Cervix, Salivary Gland, Langerhans Histiocytosis

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP310

ISOTYPE: IgG


CONTROL: Breast, Tonsil, Breast Carcinoma, Small Lymphocytic Lymphoma, Langerhans Histiocytosis

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6842 | Tinto Predilute | 3.0 ml | BSB 3370 | Tinto Predilute | 3.0 ml | BSB 3377 | Tinto Predilute | 3.0 ml |
| BSB 6843 | Tinto Predilute | 7.0 ml | BSB 3371 | Tinto Predilute | 7.0 ml | BSB 3378 | Tinto Predilute | 7.0 ml |
| BSB 6844 | Tinto Predilute | 15.0 ml | BSB 3372 | Tinto Predilute | 15.0 ml | BSB 3379 | Tinto Predilute | 15.0 ml |
| BSB 6845 | Concentrate | 0.1 ml | BSB 3373 | Concentrate | 0.1 ml | BSB 3380 | Concentrate | 0.1 ml |
| BSB 6846 | Concentrate | 0.5 ml | BSB 3374 | Concentrate | 0.5 ml | BSB 3381 | Concentrate | 0.5 ml |
| BSB 6847 | Concentrate | 1.0 ml | BSB 3375 | Concentrate | 1.0 ml | BSB 3382 | Concentrate | 1.0 ml |
| BSB 6848 | Control Slides | 5 | BSB 3376 | Control Slides | 5 | BSB 3383 | Control Slides | 5 |

LH, MMab




IHC of LH on a FFPE Pituitary Tissue

LIN28, RMAb



IHC of LIN28 on an FFPE Seminoma Tissue

LMO2, RMAb



IHC of LMO2 on an FFPE Lymphoblastic Lymphoma Tissue

Luteinizing hormone (LH) is a hormone synthesized and secreted by gonadotropes in the anterior lobe of the pituitary gland. In concert with the other pituitary gonadotropin follicle-stimulating hormone (FSH), it is necessary for proper reproductive function. In the female, an acute rise of LH levels triggers ovulation. In the male, where LH has also been called Interstitial Cell-Stimulating Hormone (ICSH), it stimulates Leydig cell production of testosterone.

LH is a useful marker in classification of Pituitary Tumors and the study of pituitary disease. LH antibody reacts with LH-producing cells (gonadotrophs).

LIN28 homolog A is a protein that in humans is encoded by the LIN28 gene. LIN28 is thought to regulate the self-renewal of stem cells, is highly expressed in human embryonic stem cells and can enhance the efficiency of the formation of induced pluripotent stem (iPS) cells from human fibroblasts.

LIN28 has been found to be a highly sensitive marker for testicular intratubular germ cell neoplasias, classic seminomas, embryonal carcinomas, and yolk sac tumors (YST) with relatively high specificity. LIN28 can be used as a diagnostic marker for these tumors and has demonstrated a similar level of diagnostic utility as SALL4. The major advantage of LIN28 over OCT4 is in diagnosing yolk sac tumors (yolk sac tumors negative for OCT4). In another study, LIN28 was found to be a sensitive marker of ovarian primitive germ cell tumors like Gonadoblastomas, Dysgerminomas, Embryonal Carcinomas, and YSTs. LIN28 can be used to distinguish them from non-testicular germ cell tumors. High expression of Lin28 is associated with poor prognosis and high tumor aggressiveness in esophageal cancer and these effects are mediated through increased proliferation and invasiveness of esophageal cancer cells.

LIM domain only 2 (rhombotin-like 1), also known as LIM Domain Only Protein 2 and T-Cell Translocation Protein 2, is a protein which in humans is encoded by the LMO2 gene. LMO2 encodes a cysteine-rich, two LIM domain protein that is required for yolk sac erythropoiesis. The LMO2 protein has a central and crucial role in hematopoietic development and is highly conserved.

HGAL and LMO2 have been found helpful in classifying difficult cases of Follicular Lymphoma (FL) as an adjunct in the identification of FL of the nongastric GI tract. LMO2 expression has been reported to be special feature of GC DLBCL (Diffuse Large B Cell Lymphoma of germinal center subtype) which can be used as a diagnostic marker. LMO2 has shown usefulness as part of an IHC panel of germinal center-associated markers in eliminating cases of Diffuse Follicle Center Lymphoma. This is accomplished by taking into consideration the histologic and immunohistochemical spectrum of Nodal Marginal Zone Lymphoma (NMZL) and the immunohistochemical analysis for CD43, CD23, CD21, BCL6, HGAL, and LMO2 in the diagnosis of NMZL.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-53
ISOTYPE: IgG1/K
CONTROL: Normal Pituitary
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rabbit, Pig

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP150
ISOTYPE: IgG
CONTROL: Testis, Seminoma, Dysgerminoma, Yolk Sac Tumor, Embryonal Carcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-LM02
ISOTYPE: IgG
CONTROL: Tonsil, Spleen, Placenta, Follicular and Lymphoblastic Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5722 | Tinto Predilute | 3.0 ml |
| BSB 5723 | Tinto Predilute | 7.0 ml |
| BSB 5724 | Tinto Predilute | 15.0 ml |
| BSB 5725 | Concentrate | 0.1 ml |
| BSB 5726 | Concentrate | 0.5 ml |
| BSB 5727 | Concentrate | 1.0 ml |
| BSB 5728 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3560 | Tinto Predilute | 3.0 ml |
| BSB 3561 | Tinto Predilute | 7.0 ml |
| BSB 3562 | Tinto Predilute | 15.0 ml |
| BSB 3563 | Concentrate | 0.1 ml |
| BSB 3564 | Concentrate | 0.5 ml |
| BSB 3565 | Concentrate | 1.0 ml |
| BSB 3566 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3574 | Tinto Predilute | 3.0 ml |
| BSB 3575 | Tinto Predilute | 7.0 ml |
| BSB 3576 | Tinto Predilute | 15.0 ml |
| BSB 3577 | Concentrate | 0.1 ml |
| BSB 3578 | Concentrate | 0.5 ml |
| BSB 3579 | Concentrate | 1.0 ml |
| BSB 3580 | Control Slides | 5 |

Lysozyme, RPab




IHC of Lysozyme on a FFPE Tonsil Tissue

Lysozyme, RMAb



IHC of Lysozyme on a FFPE Tonsil Tissue

Macrophage HAM-56, MAb



IHC of Macrophage on a FFPE Liver Tissue

Lysozyme is a 14.4 kDa enzyme, commonly referred to as the "body's own antibiotic" since it kills bacteria. Lysozyme is an enzyme that destroys bacterial cell walls by hydrolyzing the polysaccharide component of the cell wall. It is abundantly present in a number of secretions, including tears. This protein is present in cytoplasmic granules of the polymorphonuclear neutrophils (PMN) and released through mucosal secretions such as tears and saliva. They can also be found in high concentration in egg white.

Lysozyme stains myeloid cells, histiocytes, granulocytes, macrophages, and monocytes in human tonsil, colon and skin. It is an important marker that may demonstrate the myeloid or monocytic nature of Acute Leukemia. The restrictive nature of Lysozyme antibody staining suggests that Lysozyme may be synthesized predominantly in reactive histiocytes rather than in resting, unstimulated phagocytes. It has not been determined whether Lysozyme stains any other cell or tissue type. Lysozyme may aid in the identification of histiocytic neoplasias and large lymphocytes, as well as classifying lymphoproliferative disorders.

Lysozyme is a 14.4 kDa enzyme, commonly referred to as the "body's own antibiotic" since it kills bacteria. Lysozyme is an enzyme that destroys bacterial cell walls by hydrolyzing the polysaccharide component of the cell wall. It is abundantly present in a number of secretions, including tears. This protein is present in cytoplasmic granules of the polymorphonuclear neutrophils (PMN) and released through mucosal secretions such as tears and saliva. They can also be found in high concentration in egg white.

Lysozyme stains myeloid cells, histiocytes, granulocytes, macrophages, and monocytes in human tonsil, colon and skin. It is an important marker that may demonstrate the myeloid or monocytic nature of Acute Leukemia. The restrictive nature of Lysozyme antibody staining suggests that Lysozyme may be synthesized predominantly in reactive histiocytes rather than in resting, unstimulated phagocytes. It has not been determined whether Lysozyme stains any other cell or tissue type. Lysozyme may aid in the identification of histiocytic neoplasias and large lymphocytes, as well as classifying lymphoproliferative disorders.

Macrophages comprise many forms of mononuclear phagocytes found in tissues that derive from hematopoietic stem cells in the bone marrow. Among the functions of macrophages are nonspecific phagocytosis and pinocytosis, killing of ingested microorganisms, and digestion and presentation of antigens to T and B-lymphocytes. Macrophages work to secrete a large number of diverse products such as lysozyme and collagenases, several complement components and coagulation factors, some prostaglandins and leukotrienes, and many regulatory molecules (Interferon, Interleukin 1).

Macrophage HAM-56 reacts with tingible macrophages (found in the germinal centers of lymph nodes), interdigitating macrophages of lymph nodes and tissue macrophages, (e.g., Kupffer cells of the liver and alveolar macrophages of the lung). This antibody also stains a subpopulation of endothelial cells, most prominently those of the capillaries and smaller blood vessels. HAM-56 reacts with monocytes, but is unreactive with B and T-lymphocytes.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Liver, Kidney, Spleen, Lung, Salivary Gland, Cervix, Pancreas, Bone Marrow, Colon, Lymphoblastic Lymphoma, Transitional Cell Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP134

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Liver, Kidney, Spleen, Salivary Gland, Cervix, Pancreas, Bone Marrow, Colon, Lung

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Mouse Monoclonal

CLONE: HAM-56

ISOTYPE: IgM/K

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human, Monkey

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5729 | Tinto Predilute | 3.0 ml |
| BSB 5730 | Tinto Predilute | 7.0 ml |
| BSB 5731 | Tinto Predilute | 15.0 ml |
| BSB 5732 | Concentrate | 0.1 ml |
| BSB 5733 | Concentrate | 0.5 ml |
| BSB 5734 | Concentrate | 1.0 ml |
| BSB 5735 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6856 | Tinto Predilute | 3.0 ml |
| BSB 6857 | Tinto Predilute | 7.0 ml |
| BSB 6858 | Tinto Predilute | 15.0 ml |
| BSB 6859 | Concentrate | 0.1 ml |
| BSB 6860 | Concentrate | 0.5 ml |
| BSB 6861 | Concentrate | 1.0 ml |
| BSB 6862 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5736 | Tinto Predilute | 3.0 ml |
| BSB 5737 | Tinto Predilute | 7.0 ml |
| BSB 5738 | Tinto Predilute | 15.0 ml |
| BSB 5739 | Concentrate | 0.1 ml |
| BSB 5740 | Concentrate | 0.5 ml |
| BSB 5741 | Concentrate | 1.0 ml |
| BSB 5742 | Control Slides | 5 |

Mammaglobin, RMab



IHC of Mammaglobin on a FFPE Breast Tissue

Mammaglobin is a gene that encodes a 10 kDa glycoprotein. In humans, expression of the gene is limited to the adult mammary gland. A correlation between increased expression of the gene and Breast Cancer has been reported. Mammaglobin mRNA is present in high levels in human Breast Cancer cell lines and primary Breast Cancers. High levels of mRNA have been detected in normal human sweat glands as well, but are absent in Sweat Gland Tumors.

Anti-Mammaglobin (EP249) has been shown to be effective in detecting up to 85% of Breast Carcinomas using immunohistochemical techniques. Studies investigating the detection of mRNA by RT PCR from circulating carcinoma cells in the peripheral blood of Breast Cancer patients have shown that mammaglobin is a highly-specific marker and correlates with several prognostic factors, such as lymph node involvement.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP249


ISOTYPE: IgG

CONTROL: Breast, Skin, Fallopian Tube, Breast Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Monkey, Rat

Mart-1/Melan-A, MMab



IHC of MART-1 M2-7C10 on a FFPE Melanoma Tissue

MART-1 M2-7C is a putative 18 kDa transmembrane protein consisting of 118 amino acids. It has a single transmembrane domain. MART-1/Melan-A is a protein antigen found on melanocytes. Antibodies against this antigen are used to recognize cells of melanocytic differentiation, useful for the diagnosis of Melanoma. The same name is used to refer to the gene which codes for this antigen.

The MART-1 M2-7C antigen is specific for the melanocyte lineage found in normal skin, retina, and melanocytes, but not in other normal tissues. It is thus useful as a marker for Melanocytic Tumors, with the caveat that it is normally found in benign nevi as well. This antibody is very useful in establishing the diagnosis of Metastatic Melanomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: M2-7C10


ISOTYPE: IgG2b/K

CONTROL: Normal Skin, Melanoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Mart-1/Melan-A, MMab



IHC of MART-1/Melan-A on a FFPE Melanoma Tissue

MART-1/Melan-A is a putative 18 kDa transmembrane protein consisting of 118 amino acids. It has a single transmembrane domain. MART-1/Melan-A is a protein antigen found on melanocytes. Antibodies against this antigen are used to recognize cells of melanocytic differentiation, useful for the diagnosis of Melanoma. The same name is used to refer to the gene which codes for this antigen.

The MART-1/Melan-A antigen is specific for the melanocyte lineage found in normal skin, retina, and melanocytes, but not in other normal tissues. It is thus useful as a marker for Melanocytic Tumors, with the caveat that it is normally found in benign nevi as well. This antibody is very useful in establishing the diagnosis of Metastatic Melanomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: A103

ISOTYPE: IgG1

CONTROL: Normal Skin, Melanoma

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5743 | Tinto Predilute | 3.0 ml |
| BSB 5744 | Tinto Predilute | 7.0 ml |
| BSB 5745 | Tinto Predilute | 15.0 ml |
| BSB 5746 | Concentrate | 0.1 ml |
| BSB 5747 | Concentrate | 0.5 ml |
| BSB 5748 | Concentrate | 1.0 ml |
| BSB 5749 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5750 | Tinto Predilute | 3.0 ml |
| BSB 5751 | Tinto Predilute | 7.0 ml |
| BSB 5752 | Tinto Predilute | 15.0 ml |
| BSB 5753 | Concentrate | 0.1 ml |
| BSB 5754 | Concentrate | 0.5 ml |
| BSB 5755 | Concentrate | 1.0 ml |
| BSB 5756 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6870 | Tinto Predilute | 3.0 ml |
| BSB 6871 | Tinto Predilute | 7.0 ml |
| BSB 6872 | Tinto Predilute | 15.0 ml |
| BSB 6873 | Concentrate | 0.1 ml |
| BSB 6874 | Concentrate | 0.5 ml |
| BSB 6875 | Concentrate | 1.0 ml |
| BSB 6876 | Control Slides | 5 |

Maspin, MAb




IHC of Maspin on a FFPE Breast Tissue

MCM2, RMAb



IHC of MCM-2 on a FFPE HSIL of Anal Carcinoma

MCM3, RMAb



IHC of MCM3 on a FFPE Colon Carcinoma Tissue

Maspin (mammary serine protease inhibitor) is a protein that in humans is encoded by the SERPINB5 gene. This protein belongs to the serpin (serine protease inhibitor) superfamily. SERPINB5 was originally reported to function as a tumor suppressor gene in epithelial cells, suppressing the ability of cancer cells to invade and metastasize to other tissues.

Maspin has been shown in primary breast cancer to be regulated by wild-type p53, defining a new category of molecular targets of p53 that have the potential to negatively regulate tumor invasion and metastasis. Loss of Maspin expression correlates with increased tumor aggressiveness and poor prognosis in advanced breast and prostate cancer. In contrast, Maspin has been shown to be overexpressed in pancreatic, ovarian, thyroid, gastric, lung, bladder, breast, skin and colon cancer. Several studies have investigated the prognostic significance of Maspin expression in lung cancer. In primary non-small cell lung cancer (NSCLC) was an independent negative prognostic factor for overall survival, whereas strong nuclear Maspin expression was associated with increased overall and disease-free survival in patients with resectable NSCLC.

MCM-2 (mini-chromosome maintenance 2) is a human gene. The protein encoded by this gene is one of the highly conserved mini-chromosome maintenance proteins (MCM) that is involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex, and may be involved in the formation of replication forks and in the recruitment of other DNA replication-related proteins. This protein forms a complex with MCM-4, 6, and 7, and has been shown to regulate the helicase activity of the complex. This protein is phosphorylated, and thus regulated by protein kinases CDC2 and CDC7.

MCM-2 is essential for eukaryotic DNA replication and drives the formation of pre-replicative complexes, which is the key first step during the G1 phase. Therefore, altered MCM-2 expression may be a hallmark of cell-cycle deregulation, which could be the most essential mechanism in the development and progression of human cancers. MCM2 has been identified by DNA microarray and transcriptional profiling as a gene that is over-expressed in Cervical Carcinomas. This protein is over-expressed in Cervical Dysplasia as a result of HPV infection. The over-expression of MCM-2 provides the link between oncogenic HPV infection and the molecular event of Cervical Dysplasia.

MCM3 is one of the highly conserved mini-chromosome maintenance proteins (MCM) that are involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex and may be involved in the formation of replication forks and in the recruitment of other DNA replication proteins. MCM3 is a subunit of the protein complex that consists of MCM2-7. It has been shown to interact directly with MCM5/CDC46. This protein also interacts with, and thus is acetylated by MCM3AP, a chromatin-associated acetyltransferase. The acetylation of this protein inhibits the initiation of DNA replication and cell cycle progression.

Increased expression of MCM3 has been demonstrated in various tumors by immunohistochemistry and is used as a marker for tumor progression.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-92
ISOTYPE: IgG2a
CONTROL: Prostate, Breast, Tonsil
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-MCM2
ISOTYPE: IgG
CONTROL: HSIL, Cervical, Breast Cancer
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP202
ISOTYPE: IgG
CONTROL: Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3265 | Tinto Predilute | 3.0 ml |
| BSB 3266 | Tinto Predilute | 7.0 ml |
| BSB 3267 | Tinto Predilute | 15.0 ml |
| BSB 3268 | Concentrate | 0.1 ml |
| BSB 3269 | Concentrate | 0.5 ml |
| BSB 3270 | Concentrate | 1.0 ml |
| BSB 3271 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6331 | Tinto Predilute | 3.0 ml |
| BSB 6332 | Tinto Predilute | 7.0 ml |
| BSB 6333 | Tinto Predilute | 15.0 ml |
| BSB 6334 | Concentrate | 0.1 ml |
| BSB 6335 | Concentrate | 0.5 ml |
| BSB 6336 | Concentrate | 1.0 ml |
| BSB 6337 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6884 | Tinto Predilute | 3.0 ml |
| BSB 6885 | Tinto Predilute | 7.0 ml |
| BSB 6886 | Tinto Predilute | 15.0 ml |
| BSB 6887 | Concentrate | 0.1 ml |
| BSB 6888 | Concentrate | 0.5 ml |
| BSB 6889 | Concentrate | 1.0 ml |
| BSB 6890 | Control Slides | 5 |

MCM5, RMAb




IHC of MCM5 on a FFPE Cervical Cancer Tissue

MDM2, MAb



IHC of MDM2 on a FFPE Testicular Cancer Tissue

MDR-1, MAb



IHC of MDR-1 on a FFPE Kidney Tissue

DNA replication licensing factor MCM5 is a member of the MCM family and is responsible for regulating DNA replication. It functions as a replicative helicase, the molecular motor that both unwinds duplex DNA and powers fork progression during DNA replication. MCM5 is upregulated in the transition from the G0 to the G1/S phase of the cell cycle and may actively participate in cell cycle regulation.

MCM5 may be a useful marker for skin cancer, colon cancer, and is of prognostic value in colon cancer and ovarian cancer.

MDM2 is a protein that in humans is encoded by the MDM2 gene. MDM2 is an important negative regulator of the p53 tumor suppressor. The human homologue of this protein is sometimes called HDM2. Further supporting the role of MDM2 as an oncogene, several human tumor types have been shown to have increased levels of MDM2, including soft tissue sarcomas and osteosarcomas as well as breast tumors.

Well Differentiated Liposarcomas (WDLPS), Atypical Lipomatous Tumor/Well-Differentiated Liposarcoma (ALT-WDLPS) and Dedifferentiated Liposarcoma (DDLPS) may be difficult to distinguish from benign Adipose Tumors and from Poorly Differentiated Sarcomas, respectively. Genetically, they are characterized by amplification of MDM2 and CDK4 genes on chromosome 12q13-15. MDM2 and CDK4 protein overexpression have also been identified in these tumors. Detection of MDM2/CDK4 protein overexpression by IHC can be used to diagnose WDLPS and DDLPS. Immunohistochemical expression of MDM2 and CDK4 is specific and provides sensitive markers for the diagnosis of Low-grade Osteosarcomas, helping to differentiate them from benign fibrous and fibro-osseous lesions, particularly in cases with atypical radio-clinical presentation and/or limited biopsy samples.

P-glycoprotein 1, also known as multidrug resistance protein 1 (MDR1) or ATP-binding cassette sub-family B member 1 (ABCB1) or cluster of differentiation 243 (CD243), functions as an energy-dependent efflux pump for structurally diverse agents ranging from ions to peptides. It is implicated in the development of the multiple drug resistance phenomenon observed in human cancer cells following prolonged chemotherapy. The classic form of multiple drug resistance is associated with an increase in the MDR protein, but not in all cases. MDR-1 is an apical transmembrane protein that is an integral part of the blood-brain barrier and functions as a drug transport pump that transports a variety of drugs from the brain back into the blood.

MDR-1 is extensively distributed and expressed in the intestinal epithelium, hepatocytes, renal proximal tubular cells, adrenal gland and capillary endothelial cells comprising the blood-brain and blood-testis barrier.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-MCM5

ISOTYPE: IgG

CONTROL: Tonsil, Testis, Breast, Colon, Spleen, Cervical Carcinoma, Colon Carcinoma, Transitional Cell Carcinoma, Lymphoblastic Lymphoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-64

ISOTYPE: IgG1

CONTROL: Testis, Tonsil, Cervix, Placenta, LipoSarcoma & Soft Tissue, Testicular Cancer

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Mouse Monoclonal

CLONE: JSB-1

ISOTYPE: IgG1

CONTROL: Skeletal Muscle, Kidney, Adrenal, Liver

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6891 | Tinto Predilute | 3.0 ml |
| BSB 6892 | Tinto Predilute | 7.0 ml |
| BSB 6893 | Tinto Predilute | 15.0 ml |
| BSB 6894 | Concentrate | 0.1 ml |
| BSB 6895 | Concentrate | 0.5 ml |
| BSB 6896 | Concentrate | 1.0 ml |
| BSB 6897 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2978 | Tinto Predilute | 3.0 ml |
| BSB 2979 | Tinto Predilute | 7.0 ml |
| BSB 2980 | Tinto Predilute | 15.0 ml |
| BSB 2981 | Concentrate | 0.1 ml |
| BSB 2982 | Concentrate | 0.5 ml |
| BSB 2983 | Concentrate | 1.0 ml |
| BSB 2984 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6898 | Tinto Predilute | 3.0 ml |
| BSB 6899 | Tinto Predilute | 7.0 ml |
| BSB 6900 | Tinto Predilute | 15.0 ml |
| BSB 6901 | Concentrate | 0.1 ml |
| BSB 6902 | Concentrate | 0.5 ml |
| BSB 6903 | Concentrate | 1.0 ml |
| BSB 6904 | Control Slides | 5 |

MDR-1, RMAb




IHC of MDR-1 on a FFPE Adrenal Tissue

P-glycoprotein 1, also known as multidrug resistance protein 1 (MDR1) or ATP-binding cassette sub-family B member 1 (ABCB1) or cluster of differentiation 243 (CD243), functions as an energy-dependent efflux pump for structurally diverse agents ranging from ions to peptides. It is implicated in the development of the multiple drug resistance phenomenon observed in human cancer cells following prolonged chemotherapy. The classic form of multiple drug resistance is associated with an increase in the MDR protein, but not in all cases. MDR-1 is an apical transmembrane protein that is an integral part of the blood-brain barrier and functions as a drug transport pump that transports a variety of drugs from the brain back into the blood.

MDR-1 is extensively distributed and expressed in the intestinal epithelium, hepatocytes, renal proximal tubular cells, adrenal gland and capillary endothelial cells comprising the blood-brain and blood-testis barrier.

Melanoma Cocktail: HMB-45, Mart-1 & Tyrosinase, MAb



IHC of Melanoma Cocktail on a FFPE Melanoma Tissue


HMB-45 reacts against an antigen present in immature melanosomes, cutaneous, melanocytes, prenatal and infantile retinal pigment epithelium and melanoma cells. This antibody is very useful to identify Malignant Melanoma.

MART-1/Melan-A is a protein antigen found on melanocytes. Antibodies against this antigen are used to recognize cells of melanocytic differentiation, useful for the diagnosis of Melanoma. The same name is used to refer to the gene which codes for this antigen. The MART-1/Melan-A antigen is specific for the melanocyte lineage found in normal skin, retina, and melanocytes, but not in other normal tissues. It is thus useful as a marker for Melanocytic Tumors, with the caveat that it is normally found in benign nevi as well.

Tyrosinase is a copper-containing enzyme present in plant and animal tissues that catalyzes the production of melanin and other pigments from tyrosine by oxidation. Anti-Tyrosinase has been found to be quite specific for melanotic lesions such as Malignant Melanoma and Melanotic Neurofibroma. Essentially no carcinomas express this marker.

Melanoma cocktail HMB-45, Mart-1 and Tyrosinase are ideally suited to detect melanomas and melanocytic lesions and may prove to be a valuable marker for melanoma metastasis in sentinel lymph nodes.

Melanoma KBA-6.2, MAb



IHC of Melanoma/KBA.62 on a FFPE Melanoma Tissue

KBA.62 is a mouse monoclonal antibody that reacts against an antigen present in melanocytic tumors such as Melanomas. The antibody was generated to an extract of Melanoma. It reacted positively against Melanocytic Tumors but not other tumors, thus demonstrating specificity and sensitivity. Moreover, this antibody reacts positively against junctional nevus cells but not intradermal nevi, and against fetal melanocytes but not normal adult melanocytes.

This antibody is very useful to identify Malignant Melanoma. Metastatic Amelanotic Melanoma can often be confused with a variety of poorly differentiated Carcinomas, Large Cell Lymphomas, Sarcomas, Spindle Cell Carcinomas and various types of mesenchymal neoplasms. A keratin-negative, vimentin-rich neoplasm that immunoreacts with antibody to S-100 protein and with this melanoma antibody, is, with rare exception, a Melanoma.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP271

ISOTYPE: IgG

CONTROL: Skeletal Muscle, Kidney, Adrenal, Liver

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mammalia

ANTIBODY TYPE: Mouse Monoclonal

CLONE: HMB-45, A103 & BSB-6

ISOTYPE: IgG1/K, IgG1 & IgG2a

CONTROL: Skin, Melanoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog

ANTIBODY TYPE: Mouse Monoclonal

CLONE: KBA6.2

ISOTYPE: IgG1

CONTROL: Melanoma

LOCALIZATION: Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2531 | Tinto Predilute | 3.0 ml |
| BSB 2532 | Tinto Predilute | 7.0 ml |
| BSB 2533 | Tinto Predilute | 15.0 ml |
| BSB 2534 | Concentrate | 0.1 ml |
| BSB 2535 | Concentrate | 0.5 ml |
| BSB 2536 | Concentrate | 1.0 ml |
| BSB 2537 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6877 | Tinto Predilute | 3.0 ml |
| BSB 6878 | Tinto Predilute | 7.0 ml |
| BSB 6879 | Tinto Predilute | 15.0 ml |
| BSB 6880 | Concentrate | 0.1 ml |
| BSB 6881 | Concentrate | 0.5 ml |
| BSB 6882 | Concentrate | 1.0 ml |
| BSB 6883 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6905 | Tinto Predilute | 3.0 ml |
| BSB 6906 | Tinto Predilute | 7.0 ml |
| BSB 6907 | Tinto Predilute | 15.0 ml |
| BSB 6908 | Concentrate | 0.1 ml |
| BSB 6909 | Concentrate | 0.5 ml |
| BSB 6910 | Concentrate | 1.0 ml |
| BSB 6911 | Control Slides | 5 |

Melanoma PNL2, MAb



IHC of Melanoma/PNL2 on a FFPE Melanoma Tissue


PNL2 is a mouse monoclonal antibody that reacts against an antigen present in melanocytic tumors such as Melanomas. The antibody was generated to an extract of Melanoma. It reacted positively against Melanocytic Tumors but not other tumors, thus demonstrating specificity and sensitivity. Moreover, this antibody reacts positively against junctional nevus cells but not intradermal nevi, and against fetal melanocytes but not normal adult melanocytes.

This antibody is very useful to identify Malignant Melanoma. Metastatic Amelanotic Melanoma can often be confused with a variety of poorly differentiated Carcinomas, Large Cell Lymphomas, Sarcomas, Spindle Cell Carcinomas and various types of mesenchymal neoplasms. A keratin-negative, vimentin-rich neoplasm that immunoreacts with antibody to S-100 protein and with this melanoma antibody, is, with rare exception, a Melanoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: PNL2
ISOTYPE: IgG1/K
CONTROL: Melanoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6912 | Tinto Predilute | 3.0 ml |
| BSB 6913 | Tinto Predilute | 7.0 ml |
| BSB 6914 | Tinto Predilute | 15.0 ml |
| BSB 6915 | Concentrate | 0.1 ml |
| BSB 6916 | Concentrate | 0.5 ml |
| BSB 6917 | Concentrate | 1.0 ml |
| BSB 6918 | Control Slides | 5 |

Melanosome HMB45, MAb



IHC of Melanosome HMB-45 on a FFPE Malignant Melanoma Tissue


HMB-45 reacts against an antigen present in immature melanosomes, cutaneous, melanocytes, prenatal and infantile retinal pigment epithelium and melanoma cells. This antibody was generated to an extract of Melanoma. It reacted positively against Melanocytic Tumors but not other tumors, thus demonstrating specificity and sensitivity. Moreover, this antibody reacts positively against junctional nevus cells but not intradermal nevi, and against fetal melanocytes but not normal adult melanocytes.

This antibody is very useful to identify Malignant Melanoma. Metastatic Amelanotic Melanoma can often be confused with a variety of poorly differentiated Carcinomas, Large Cell Lymphomas, Sarcomas, Spindle Cell Carcinomas and various types of mesenchymal neoplasms. A keratin-negative, vimentin-rich neoplasm that immunoreacts with antibody to S-100 protein and with this melanoma antibody, is, with rare exception, a Melanoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: HMB-45
ISOTYPE: IgG1/K
CONTROL: Melanoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5757 | Tinto Predilute | 3.0 ml |
| BSB 5758 | Tinto Predilute | 7.0 ml |
| BSB 5759 | Tinto Predilute | 15.0 ml |
| BSB 5760 | Concentrate | 0.1 ml |
| BSB 5761 | Concentrate | 0.5 ml |
| BSB 5762 | Concentrate | 1.0 ml |
| BSB 5763 | Control Slides | 5 |

Mesothelial Cell, MAb



IHC of Mesothelial Cell on a FFPE Mesothelioma Tissue


Mesothelial Cell HBME-1 has shown to label mesothelial cells, both benign and malignant and thus has been used in distinguishing mesothelioma from adenocarcinomas of various origins. HBME-1 has also been used to distinguish Thyroid carcinomas (both Follicular and Papillary) from benign thyroid lesions.

Mesothelial Cell HBME-1 and MOC-31 have been shown to have a diagnostic efficiency for the distinction between carcinoma and mesothelioma in pleura. HBME-1 staining may be useful for differentiating papillary carcinomas from follicular carcinomas; in papillary lesions it tends to be positive. Several immunohistochemical markers have been used to aid in the diagnosis of follicular-derived lesions of the thyroid (FDLT). HBME-1, ERK, and p16 were found to be more specific for malignancy, whereas CK19 and GAL-3 stained benign lesions with a higher frequency and were not specific for malignant FDLT.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: HBME-1
ISOTYPE: IgGM/K
CONTROL: Breast, Tonsil, Lung, Salivary Gland, TCC, Mesothelioma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3455 | Tinto Predilute | 3.0 ml |
| BSB 3456 | Tinto Predilute | 7.0 ml |
| BSB 3457 | Tinto Predilute | 15.0 ml |
| BSB 3458 | Concentrate | 0.1 ml |
| BSB 3459 | Concentrate | 0.5 ml |
| BSB 3460 | Concentrate | 1.0 ml |
| BSB 3461 | Control Slides | 5 |

Mesothelin, RMaB



IHC of Mesothelin on a FFPE Ovarian Tissue


Mesothelin was first identified by its reactivity with monoclonal antibody K1. The mesothelin gene encodes a precursor protein that is processed to yield mesothelin, which is attached to the cell membrane by a glycoposphatidylinositol linkage and a 31-KDa shed fragment named megakaryocyte-potentiating factor (MPF). Its biological function is not known, but recent studies have shown that mesothelin forms a strong and specific complex with MUC16, which has been suggested to be the basis of ovarian cancer metastasis.

Mesothelin is present on normal mesothelial cells lining the pleura, peritoneum, and pericardium. In tumors, overexpression of Mesothelin has been observed in mesotheliomas, and other tumors including ovarian, pancreatic carcinomas, and cholangiocarcinoma.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP140
ISOTYPE: IgG
CONTROL: Cervix, Tonsil, Mesothelioma, Ovarian Cancer
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6919 | Tinto Predilute | 3.0 ml |
| BSB 6920 | Tinto Predilute | 7.0 ml |
| BSB 6921 | Tinto Predilute | 15.0 ml |
| BSB 6922 | Concentrate | 0.1 ml |
| BSB 6923 | Concentrate | 0.5 ml |
| BSB 6924 | Concentrate | 1.0 ml |
| BSB 6925 | Control Slides | 5 |

MGMT/AGAT, RMaB



IHC of MGMT/AGAT on a FFPE Lung Adenocarcinoma Tissue

MGMT is a house-keeping gene expressed in all tissues, and its promoter methylation and resulting down-regulation of the AGT protein varies among tumor types. The MGMT promoter contains many regulatory domains, and methylation can help decrease protein expression and therefore resistance to drugs, improving effectiveness of alkylating treatments.


High MGMT expression has been reported in Glioma, Myeloma, Melanoma, Colon and Pancreatic Cancers, where inhibition or down-regulation make alkylating drug treatments more effective. MGMT expression is lost in some cancer types such as Lymphomas, Non-Small Cell Lung Cancer, Astrocyte Tumors, Breast, and Prostate Cancer, where it may fail to prevent mutations.

A study detected protein expression of MGMT by IHC in a series of newly diagnosed Glioblastomas and found that the MGMT status detected by either IHC or MSP was significantly correlated with the treatment response and survival of Glioblastoma patients. Another study addresses that MGMT in premalignant lesions from smokers and lung adenocarcinomas, their biological effects on gene regulation, and targeting MGMT for therapy and found strong evidence that the An allele of a MGMT promoter-enhancer SNP is a key determinant for MGMT methylation in lung carcinogenesis. Moreover, TMZ treatment may benefit a subset of lung cancer patients methylated for MGMT.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP337
ISOTYPE: IgG
CONTROL: Breast, Cervix, Fallopian Tube, Prostate, Testis, Transitional Cell Carcinoma, Lung Adenocarcinoma, Ductal Breast Carcinoma, Endometrial Carcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3737-3 | Tinto Predilute | 3.0 ml |
| BSB-3737-7 | Tinto Predilute | 7.0 ml |
| BSB-3737-15 | Tinto Predilute | 15.0 ml |
| BSB-3737-01 | Concentrate | 0.1 ml |
| BSB-3737-05 | Concentrate | 0.5 ml |
| BSB-3737-1 | Concentrate | 1.0 ml |
| BSB-3737-CS | Control Slides | 5 |

MiTF, MMaB



IHC of MiTF on a FFPE Melanoma Tissue


Microphthalmia-associated Transcription Factor (MiTF) is a basic helix-loop-helix leucine zipper transcription factor involved in melanocyte and osteoclast development. Mutations in MiTF cause auditory pigmentary syndromes, such as Waardenburg Syndrome Type II, Type IIa and Tietz Syndrome in humans. There are two known isoforms of MiTF differing by 66 amino acids at the NH2 terminus. Shorter forms are expressed in melanocytes and run as two bands at 52 kDa and 56 kDa, while the longer Mi form runs as a cluster of bands at 60-70 kDa in osteoclasts and in B16 Melanoma cells (but not other Melanoma cell lines), as well as mast cells and heart cells. MiTF plays a critical role in the differentiation of various cell types such as neural crest-derived melanocytes, mast cells, osteoclasts and optic cup-derived retinal pigment epithelium.

This antibody recognizes serine phosphorylated and non-phosphorylated melanocytic isoforms of microphthalmia. It is useful in identifying Malignant Melanoma, and distinguishing mast cell lesions from lesions of myeloid derivation. A relatively rare class of tumors known as PEComas (tumors showing perivascular epithelioid cell differentiation) express MiTF in a high percentage of cases (~90%).

ANTIBODY TYPE: Mouse Monoclonal
CLONE: C5/D5
ISOTYPE: IgG1/K
CONTROL: Skin, Melanoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Dog, Cat


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6247 | Tinto Predilute | 3.0 ml |
| BSB 6248 | Tinto Predilute | 7.0 ml |
| BSB 6249 | Tinto Predilute | 15.0 ml |
| BSB 6250 | Concentrate | 0.1 ml |
| BSB 6251 | Concentrate | 0.5 ml |
| BSB 6252 | Concentrate | 1.0 ml |
| BSB 6253 | Control Slides | 5 |

MLH1, MMab




IHC of MLH1 on a FFPE Colon Carcinoma Tissue

MMP-9, RMab



IHC of MMP-9 on a FFPE Lymphoblastic Lymphoma Tissue

MNDA, MMab



IHC of MNDA on a FFPE Diffuse Large B Cell Lymphoma Tissue

MLH1 is a mismatch repair gene of around 87 kDa, commonly associated with Hereditary Non-Polyposis Colorectal Cancer (HNPCC). This gene was identified as a locus frequently mutated in HNPCC. It is a human homolog of the E. coli DNA mismatch repair gene mutL, consistent with the characteristic alterations in microsatellite sequences (RER+ phenotype) found in HNPCC. Alternatively spliced transcript variants encoding different isoforms have been described, but their full-length natures have not been determined.

In a high proportion of patients with microsatellite instability (MSI-H), the MLH1 protein is typically deficient. This protein deficiency is linked to the autosomal dominant condition of Hereditary Non-Polyposis Colon Cancer. The anti-MLH1 antibody is useful in screening patients and families for this condition. Colon cancers that are microsatellite-unstable have a better prognosis than their microsatellite stable counterparts.

The matrix metalloproteinases (MMPs) are responsible for degradation of the extracellular matrix. The MMPs and their specific tissue inhibitor metalloproteinases (TIMP) have been associated with tumor cell invasion and metastasis in a number of adult tumors. MMP-9, also designated as 92-kDa Type IV Collagenase or gelatinase B, is a member of MMPs, which is produced as a 92- kDa pro-enzyme by neutrophils, macrophages, mast cells and stromal cells, as a normal constituent and released into the extracellular environment after activation in inflammatory tissues.

MMP-9 may be involved in the development of several human malignancies, as degradation of collagen IV in basement membrane and extracellular matrix facilitates tumor progression, including invasion, metastasis, growth and angiogenesis. The expression levels of MMP-9 in tumors are elevated compared with the corresponding normal tissues in a variety of cancer types, including breast, colon, gastric and nasopharyngeal cancers. MMP-9 may play an important role in angiogenesis and neovascularization. For example, MMP9 appears to be involved in the remodeling associated with malignant glioma neovascularization. Increased expression has been seen in a metastatic mammary cancer cell line.

Myeloid cell Nuclear Differentiation Antigen (MNDA) is present in granulocytes and monocytes (myeloid and B cells), found in the marginal zones of germinal centers in tonsil, lymph node, spleen. MNDA is a highly basic protein from 1q22, and is part of the family of nuclear proteins expressed in reaction to interferons. MNDA contains a Pyrin 'death domain' used in self-association and is suggested to induce or help prevent programmed cell death pathways in hematopoietic cells.

MNDA expression is studied mainly in leukemias and lymphomas. MNDA is a marker to distinguish Nodal Marginal Zone Lymphoma from Follicular Lymphoma, and has been shown to be downregulated in Myelodysplastic Syndrome, a precursor for Leukemia. MNDA overexpression in Osteosarcoma induces apoptosis and protects patients from metastasis, and expression in Chronic Lymphocytic Leukemia indicates the degradation of anti-apoptotic mRNA, allowing cells suffering genotoxic stress to undergo apoptosis.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: G168-728
ISOTYPE: IgG2a
CONTROL: Testis, Tonsil, Colon, Kidney, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Rat, Mouse


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP127
ISOTYPE: IgG
CONTROL: Tonsil, Spleen, Liver, Pituitary, Colon, Lymphoblastic Lymphoma, Transitional Cell Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-157
ISOTYPE: IgG1
CONTROL: Breast, Colon, Fallopian Tube, Brain, Tonsil, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5764 | Tinto Predilute | 3.0 ml |
| BSB 5765 | Tinto Predilute | 7.0 ml |
| BSB 5766 | Tinto Predilute | 15.0 ml |
| BSB 5767 | Concentrate | 0.1 ml |
| BSB 5768 | Concentrate | 0.5 ml |
| BSB 5769 | Concentrate | 1.0 ml |
| BSB 5770 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2538 | Tinto Predilute | 3.0 ml |
| BSB 2539 | Tinto Predilute | 7.0 ml |
| BSB 2540 | Tinto Predilute | 15.0 ml |
| BSB 2541 | Concentrate | 0.1 ml |
| BSB 2542 | Concentrate | 0.5 ml |
| BSB 2543 | Concentrate | 1.0 ml |
| BSB 2544 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3738-3 | Tinto Predilute | 3.0 ml |
| BSB-3738-7 | Tinto Predilute | 7.0 ml |
| BSB-3738-15 | Tinto Predilute | 15.0 ml |
| BSB-3738-01 | Concentrate | 0.1 ml |
| BSB-3738-05 | Concentrate | 0.5 ml |
| BSB-3738-1 | Concentrate | 1.0 ml |
| BSB-3738-CS | Control Slides | 5 |

MSH2, MAb

IHC of MSH2 on FFPE Colon Carcinoma Tissue

MSH2, RMAb

IHC of MSH2 on a FFPE Colon Carcinoma Tissue

MSH6, MAb

IHC of MSH6 on a FFPE Colon Carcinoma Tissue

MSH2 is a mismatch repair gene commonly associated with Hereditary Non-Polyposis Colorectal Cancer (HNPCC). This gene was identified as a locus frequently mutated in HNPCC. When cloned, it is a human homolog of the E. coli DNA mismatch repair gene mutS, consistent with the characteristic alterations in microsatellite sequences (RER+ phenotype) found in HNPCC.

MSH2 is abnormally deficient in a high proportion of patients with microsatellite instability (MSI-H). This finding is associated with the autosomal dominant condition found in Hereditary Non-Polyposis Colon Cancer. This anti-MSH2 antibody (along with MLH1 antibody) is useful in screening patients and families for this rare condition. Colon cancers that are microsatellite unstable have a better prognosis than their microsatellite stable counterparts.

MSH2, also known as mutS protein homolog 2, is a mismatch repair gene commonly associated with Hereditary Non-Polyposis Colorectal Cancer (HNPCC). This gene was identified as a locus frequently mutated in HNPCC. When cloned, it is a human homolog of the E. coli DNA mismatch repair gene mutS, consistent with the characteristic alterations in microsatellite sequences (RER+ phenotype) found in HNPCC.

MSH2 is abnormally deficient in a high proportion of patients with microsatellite instability (MSI-H). This finding is associated with the autosomal dominant condition found in Hereditary Non-Polyposis Colon Cancer. This anti-MSH2 antibody (along with MLH1 antibody) is useful in screening patients and families for this rare condition. Colon cancers that are microsatellite unstable have a better prognosis than their microsatellite stable counterparts.

MSH6 is also known as mutS homolog 6, a gene commonly associated with Hereditary Non-Polyposis Colorectal Cancer (HNPCC). HNPCC is an autosomal, dominantly inherited disease associated with marked increase in cancer susceptibility. It is characterized by a familial predisposition to early onset Colorectal Carcinoma and extra-colonic cancers of the gastrointestinal, urological and female reproductive tracts. HNPCC is reported to be the most common form of inherited Colorectal Cancer in the western world. MSH6 is a mismatch repair gene which is deficient in a high proportion of patients with microsatellite instability (MSI-H).

The anti-MSH6 antibody is useful in screening patients and families for HNPCC. Colon cancers that are microsatellite-unstable have a better prognosis than their microsatellite-stable counterparts.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-147

ISOTYPE: IgG1

CONTROL: Colon, Skin, Breast, Tonsil, Fallopian Tube, Colon Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-MSH2

ISOTYPE: IgG

CONTROL: Colon Carcinoma, Colon, Skin, Breast, Tonsil, Fallopian Tube

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 44

ISOTYPE: IgG1


CONTROL: Colon Mucosa, Colon Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Rat, Mouse, Dog


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB-3722-3 | Tinto Predilute | 3.0 ml | BSB 6926 | Tinto Predilute | 3.0 ml | BSB 6142 | Tinto Predilute | 3.0 ml |
| BSB-3722-7 | Tinto Predilute | 7.0 ml | BSB 6927 | Tinto Predilute | 7.0 ml | BSB 6143 | Tinto Predilute | 7.0 ml |
| BSB-3722-15 | Tinto Predilute | 15.0 ml | BSB 6928 | Tinto Predilute | 15.0 ml | BSB 6144 | Tinto Predilute | 15.0 ml |
| BSB-3722-01 | Concentrate | 0.1 ml | BSB 6929 | Concentrate | 0.1 ml | BSB 6145 | Concentrate | 0.1 ml |
| BSB-3722-05 | Concentrate | 0.5 ml | BSB 6930 | Concentrate | 0.5 ml | BSB 6146 | Concentrate | 0.5 ml |
| BSB-3722-1 | Concentrate | 1.0 ml | BSB 6931 | Concentrate | 1.0 ml | BSB 6147 | Concentrate | 1.0 ml |
| BSB-3722-CS | Control Slides | 5 | BSB 5777 | Control Slides | 5 | BSB 6148 | Control Slides | 5 |

MSH6, RMAb



IHC of MSH6 on a FFPE Colon Carcinoma Tissue

MTAP, RMAb



IHC of MTAP on a FFPE Melanoma Tissue

MUC1, MAb



IHC of MUC1 on a FFPE Colon Tissue

MSH6, also known as mutS homolog 6, is a gene commonly associated with Hereditary Non-Polyposis Colorectal Cancer (HNPCC). HNPCC is an autosomal, dominantly inherited disease associated with marked increase in cancer susceptibility. It is characterized by a familial predisposition to early onset Colorectal Carcinoma and extra-colonic cancers of the gastrointestinal, urological and female reproductive tracts. HNPCC is reported to be the most common form of inherited Colorectal Cancer in the western world. MSH6 is a mismatch repair gene which is deficient in a high proportion of patients with microsatellite instability (MSI-H).

The anti-MSH6 antibody is useful in screening patients and families for HNPCC. Colon cancers that are microsatellite-unstable have a better prognosis than their microsatellite-stable counterparts.

The S-methyl-5'-thioadenosine phosphorylase (MTAP) metabolizes methylthioadenosine (MTA), a byproduct of polyamine synthesis and regulator of protein methylation, to salvage adenine and methionine residues for reuse in many pathways affecting cell proliferation, signaling, and apoptosis. The MTAP gene, located at 9p12.3, is co-deleted with CDKN2A (encodes p16 tumor-suppressor) in several types of cancer and is abundant in all normal tissues but is deficient in various tumors.

IHC has been found to be an accurate and useful diagnostic method for detecting MTAP deficiency in NSCLC, and the frequency of MTAP deficiency was found to be relatively high. Accumulated MTA can upregulate transcription factors, activate ERK pathways in Melanoma and MAPK in Liver Cancer causing invasive and metastasizing tumors, and disrupt function of secretory cells in tissues like the prostate. Disrupted DNA methylation affects epigenetic factors and can lead to dedifferentiated stem-like cancer cells in Glioblastoma, although these cells can suffer purine starvation without a functioning salvage pathway. Overexpression of MTAP is also associated with increased proliferation and epithelial-to-mesenchymal transition in Colorectal Carcinoma. A study has found that a combination of MTAP or BAP1 loss detected by IHC can likely detect malignant pleural mesothelioma (MPM) and serve as a useful ancillary IHC for discriminating MPM from reactive mesothelial hyperplasia (RMH).

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-MTAP

ISOTYPE: IgG

CONTROL: Breast, Colon, Prostate, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

Mucin 1, also known as MUC1, is a human gene. This gene is a member of the mucin family and encodes a membrane-bound, glycosylated phosphoprotein. The protein is anchored to the apical surface of many epithelia by a transmembrane domain, the degree of glycosylation varying with cell type. Mucins are high molecular-weight glycoproteins which constitute the major component of the mucus layer that protects the gastric epithelium from chemical and mechanical aggressions. The MUC1 protein serves a protective function by binding to pathogens and also functions in a cell-signaling capacity. Overexpression, aberrant intracellular localization, and changes in glycosylation of this protein have been associated with carcinomas. Multiple alternatively-spliced transcript variants that encode different isoforms of this gene have been reported, but the full-length nature of only some has been determined.

MUC1 is a large cell, surface-mucin glycoprotein expressed by most glandular and ductal epithelial cells and some hematopoietic cell lineages. It is expressed on most secretory epithelium, including mammary gland and some hematopoietic cells. It is expressed in lactating mammary glands and overexpressed in more than 90% Breast Carcinomas and metastases. Transgenic MUC1 has been shown to associate with all four ceBB receptors and localize with erbB1 (EGFR) in lactating glands.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-44

ISOTYPE: IgG1

CONTROL: Breast, Colon, Kidney, Fallopian Tube, Tonsil, Colon Adenocarcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP49

ISOTYPE: IgG

CONTROL: Colon Mucosa, Colon Carcinoma


LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6932 | Tinto Predilute | 3.0 ml |
| BSB 6933 | Tinto Predilute | 7.0 ml |
| BSB 6934 | Tinto Predilute | 15.0 ml |
| BSB 6935 | Concentrate | 0.1 ml |
| BSB 6936 | Concentrate | 0.5 ml |
| BSB 6937 | Concentrate | 1.0 ml |
| BSB 6938 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3739-3 | Tinto Predilute | 3.0 ml |
| BSB-3739-7 | Tinto Predilute | 7.0 ml |
| BSB-3739-15 | Tinto Predilute | 15.0 ml |
| BSB-3739-01 | Concentrate | 0.1 ml |
| BSB-3739-05 | Concentrate | 0.5 ml |
| BSB-3739-1 | Concentrate | 1.0 ml |
| BSB-3739-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6149 | Tinto Predilute | 3.0 ml |
| BSB 6150 | Tinto Predilute | 7.0 ml |
| BSB 6151 | Tinto Predilute | 15.0 ml |
| BSB 6152 | Concentrate | 0.1 ml |
| BSB 6153 | Concentrate | 0.5 ml |
| BSB 6154 | Concentrate | 1.0 ml |
| BSB 6155 | Control Slides | 5 |

MUC1, RMAb

IHC of MUC1 on a FFPE Breast Carcinoma Tissue

MUC2, MMAb

IHC of MUC2 on a FFPE Colon Tissue

MUC2, RMAb

IHC of MUC2 on a FFPE Colon Tissue

Mucin 1, also known as MUC1, is a member of the mucin family and encodes a membrane-bound, glycosylated phosphoprotein. The protein is anchored to the apical surface of many epithelia by a transmembrane domain, the degree of glycosylation varying with cell type. Mucins are high molecular-weight glycoproteins which constitute the major component of the mucus layer that protects the gastric epithelium from chemical and mechanical aggressions. The MUC1 protein serves a protective function by binding to pathogens and also functions in a cell-signaling capacity. Overexpression, aberrant intracellular localization, and changes in glycosylation of this protein have been associated with carcinomas. Multiple alternatively-spliced transcript variants that encode different isoforms of this gene have been reported, but the full-length nature of only some has been determined.

MUC1 is a large cell, surface-mucin glycoprotein expressed by most glandular and ductal epithelial cells and some hematopoietic cell lineages. It is expressed on most secretory epithelium, including mammary gland and some hematopoietic cells. It is expressed in lactating mammary glands and overexpressed in more than 90% Breast Carcinomas and metastases. Transgenic MUC1 has been shown to associate with all four cemb receptors and localize with erbB1 (EGFR) in lactating glands.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP85

ISOTYPE: IgG

CONTROL: Breast, Colon, Kidney, Fallopian Tube, Tonsil, Colon Adenocarcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-45

ISOTYPE: IgG1/K

CONTROL: Small Intestine, Colon, Colon Adenocarcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP187

ISOTYPE: IgG

CONTROL: Small Intestine, Colon, Colon Adenocarcinoma

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6939 | Tinto Predilute | 3.0 ml |
| BSB 6940 | Tinto Predilute | 7.0 ml |
| BSB 6941 | Tinto Predilute | 15.0 ml |
| BSB 6942 | Concentrate | 0.1 ml |
| BSB 6943 | Concentrate | 0.5 ml |
| BSB 6944 | Concentrate | 1.0 ml |
| BSB 6945 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6156 | Tinto Predilute | 3.0 ml |
| BSB 6157 | Tinto Predilute | 7.0 ml |
| BSB 6158 | Tinto Predilute | 15.0 ml |
| BSB 6159 | Concentrate | 0.1 ml |
| BSB 6160 | Concentrate | 0.5 ml |
| BSB 6161 | Concentrate | 1.0 ml |
| BSB 6162 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6946 | Tinto Predilute | 3.0 ml |
| BSB 6947 | Tinto Predilute | 7.0 ml |
| BSB 6948 | Tinto Predilute | 15.0 ml |
| BSB 6949 | Concentrate | 0.1 ml |
| BSB 6950 | Concentrate | 0.5 ml |
| BSB 6951 | Concentrate | 1.0 ml |
| BSB 6952 | Control Slides | 5 |

MUC4, MAb




IHC of MUC4 on a FFPE Colon Tissue

MUC-4, RMab



IHC of MUC-4 on a FFPE Colon Tissue

MUC5AC, MAb



IHC of MUC5AC on a FFPE Stomach Tissue

Mucin 4 (MUC4) is a mucin protein that in humans is encoded by the MUC4 gene. Like other mucins, MUC4 is a high-molecular weight glycoprotein. MUC4 belongs to the human mucin family that is membrane-anchored and can range in molecular weight from 550 to 930 kDa for the actual protein. MUC4 antibody labels normal epithelial cells in the trachea, GI tract and prostate, but not in the pancreas.

MUC-4 has been found to play various roles in the progression of cancer, particularly due to its signaling and anti-adhesive properties which contribute to tumor development and metastasis. It is also found to play roles in other diseases such as endometriosis and inflammatory bowel disease. An abnormal expression of MUC4 has been reported in various carcinomas of the colon, pancreas, breast, and ovaries. Increased expression of MUC4 has been observed in pancreatic carcinoma and cervical squamous carcinoma. MUC4 is helpful in differentiating lung adenocarcinoma (positive) from malignant mesothelioma (negative). Additionally, MUC4 is useful in the identification of low-grade fibromyxoid sarcoma (LGFMS), and sclerosing epithelioid fibrosarcoma. MUC4 expression is also detected in the glandular component of biphasic synovial sarcomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 8G7

ISOTYPE: IgG1

CONTROL: Colon, Prostate, Cervix, Placenta, Salivary Gland, Pancreatic Carcinoma, Cervical Squamous Carcinoma, Bladder TCC

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Mucin 4 (MUC4) is a mucin protein that in humans is encoded by the MUC4 gene. Like other mucins, MUC4 is a high-molecular weight glycoprotein. MUC4 belongs to the human mucin family that is membrane-anchored and can range in molecular weight from 550 to 930 kDa for the actual protein. MUC4 antibody labels normal epithelial cells in the trachea, GI tract and prostate, but not in the pancreas.

MUC-4 has been found to play various roles in the progression of cancer, particularly due to its signaling and anti-adhesive properties which contribute to tumor development and metastasis. It is also found to play roles in other diseases such as endometriosis and inflammatory bowel disease. An abnormal expression of MUC4 has been reported in various carcinomas of the colon, pancreas, breast, and ovaries. Increased expression of MUC4 has been observed in pancreatic carcinoma and cervical squamous carcinoma. MUC4 is helpful in differentiating lung adenocarcinoma (positive) from malignant mesothelioma (negative). Additionally, MUC4 is useful in the identification of low-grade fibromyxoid sarcoma (LGFMS), and sclerosing epithelioid fibrosarcoma. MUC4 expression is also detected in the glandular component of biphasic synovial sarcomas.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP256

ISOTYPE: IgG

CONTROL: Colon, Liver, Kidney, Pancreatic Carcinoma, Cervical Squamous Carcinoma, Bladder TCC

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Mucin 5AC, also known as MUC5AC, is a human gene. The Mucin 5AC antigen is found in columnar mucous cells of surface gastric epithelium and in goblet cells of the fetal and precancerous colon but not in normal colon cells. Mucin genes are expressed in a regulated cell- and tissue-specific manner. MUC1 is detected in mucous cells of the surface epithelium and neck region of the gastric antrum, as well as in pyloric glands and oxyntic glands of the body region. MUC5AC is highly expressed in foveolar epithelium of both body and antrum, whereas MUC6 protein expression is limited to mucous neck cells of the body and pyloric glands of the antrum.

The mucin expression pattern of Gastric Carcinoma is heterogeneous. It includes mucins normally expressed in gastric mucosa (MUC1, MUC5AC and MUC6) and de novo expression of the intestinal mucin MUC2. The heterogeneous pattern of mucin expression, including the expression of the intestinal mucin MUC2, may provide new insights into the differentiation pathways of Gastric Carcinoma. It has been shown that in Gastric Carcinomas evaluated for expression of several mucins (MUC1, MUC2, MUC5AC and MUC6), mucin expression is associated with tumor type (MUC5AC with Diffuse and Infiltrative Carcinomas and MUC2 with Mucinoid Carcinomas) but not with the clinico-biological behavior of the tumors. Mucin expression is associated with tumor location (MUC5AC with Antrum Carcinomas and MUC2 with Cardia Carcinomas), indirectly reflecting differences in tumor differentiation according to tumor location.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: CLH2

ISOTYPE: IgG1

CONTROL: Stomach, Colon, Kidney


LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2985 | Tinto Predilute | 3.0 ml |
| BSB 2986 | Tinto Predilute | 7.0 ml |
| BSB 2987 | Tinto Predilute | 15.0 ml |
| BSB 2988 | Concentrate | 0.1 ml |
| BSB 2989 | Concentrate | 0.5 ml |
| BSB 2990 | Concentrate | 1.0 ml |
| BSB 2991 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2552 | Tinto Predilute | 3.0 ml |
| BSB 2553 | Tinto Predilute | 7.0 ml |
| BSB 2554 | Tinto Predilute | 15.0 ml |
| BSB 2555 | Concentrate | 0.1 ml |
| BSB 2556 | Concentrate | 0.5 ml |
| BSB 2557 | Concentrate | 1.0 ml |
| BSB 2558 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6163 | Tinto Predilute | 3.0 ml |
| BSB 6164 | Tinto Predilute | 7.0 ml |
| BSB 6165 | Tinto Predilute | 15.0 ml |
| BSB 6166 | Concentrate | 0.1 ml |
| BSB 6167 | Concentrate | 0.5 ml |
| BSB 6168 | Concentrate | 1.0 ml |
| BSB 6169 | Control Slides | 5 |

MUC6, MAb

IHC of MUC6 on a FFPE Stomach Tissue

MUM1, RMab

IHC of MUM1 on a FFPE Kidney Tissue

Musashi 2, RMab

IHC of Musashi 2 on a FFPE Lung Neuroendocrine Carcinoma Tissue

Mucin 6, also known as MUC6, is a human gene. Mucin is a high M.W. (>1,000 kDa) glycoprotein, expressed by mucous cells of the gastric epithelium and by goblet cells of the fetal, precancerous and cancerous colon, but not by those of the normal colon. It also appears in other epithelial tissues, which are embryologically derived from the foregut (epigastric and bronchial epithelium) and in Müller ducts (mucous cells of the endocervix and urethral epithelium near the prostatic utricle).

MUC6 antibody works well with ethanol-fixed, cultured epithelial cells and ethanol- or formalin-fixed, paraffin-embedded tissue sections. It stains the surface gastric epithelium of normal human gastrointestinal tract and reacts with fetal, precancerous and cancerous colonic mucosa, but not with normal colon.

MUM1 (multiple myeloma oncogene-1) also known as interferon regulatory factor 4 (IRF4) is a 50 kDa protein and is a member of the interferon regulatory factor family of transcription factors. It is induced by antigen receptor mediated stimuli and plays an important role in cell proliferation, differentiation and survival. MUM1 is expressed in the nuclei and cytoplasm of plasma cells and a small percentage of germinal center (GC) B-cells committed to plasmacytic or memory cell differentiation in the "light zone".

MUM1 is useful for subclassification of lymphoid malignancies and is an excellent marker for Hodgkin's and Reed-Sternberg cells of classic Hodgkin's disease.

Musashi-2 (MSI2) is an RNA-binding protein present in the cytoplasm of hematopoietic, neuronal progenitor, and stem cells. Musashi proteins have an RNA-recognition motif used to bind and stabilize mRNA in stem cell populations. MSI2 is also present in blood cells, astrocytes, and the development of reproductive cells in spermatogenesis and oogenesis.

Musashi 2 can indicate the presence of stem cells in tumors of Colorectal, Lung, and Pancreatic Cancers, and in Glioblastoma, Leukemias, and xenografts, where it supports proliferation and prevents apoptosis. Its RNA-binding ability is increased in Leukemia cells over normal hematopoietic stem cells, and its expression in Acute Myeloid Leukemia is associated with poor prognosis. MSI2-HOXA9 and other fusions/ transformations are rare but MSI2 over-expression is found in almost all hematological disorders, and can be up-regulated further resulting in worse prognosis. Studies comparing normal cervical tissues have found that the expression of MSI2 is increased in cervical cancer tissues and may act as a prognostic biomarker in patients with cervical cancer. MSI2 has been found to be significantly upregulated in bladder cancer cells and tissues compared with normal bladder urothelial cells and tissues. High expression of MSI2 in bladder cancer specimens demonstrated that MSI2 can induce bladder cancer cell migration and invasion by activating the JAK2/STAT3 pathway, and may be a valuable prognostic biomarker for bladder cancer patients.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: CLH5

ISOTYPE: IgG1

CONTROL: Stomach

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP190

ISOTYPE: IgG

CONTROL: Tonsil, Spleen, Colon, Kidney, Breast, Lymph Node, Plasmacytoma, Hodgkin's Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RM422

ISOTYPE: IgG


CONTROL: Testis, Kidney, Colon, Transitional Cell Carcinoma, Breast Cancer

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|-------------|-----------------|---------|
| BSB 6170 | Tinto Predilute | 3.0 ml | BSB 6953 | Tinto Predilute | 3.0 ml | BSB-3740-3 | Tinto Predilute | 3.0 ml |
| BSB 6171 | Tinto Predilute | 7.0 ml | BSB 6954 | Tinto Predilute | 7.0 ml | BSB-3740-7 | Tinto Predilute | 7.0 ml |
| BSB 6172 | Tinto Predilute | 15.0 ml | BSB 6955 | Tinto Predilute | 15.0 ml | BSB-3740-15 | Tinto Predilute | 15.0 ml |
| BSB 6173 | Concentrate | 0.1 ml | BSB 6956 | Concentrate | 0.1 ml | BSB-3740-01 | Concentrate | 0.1 ml |
| BSB 6174 | Concentrate | 0.5 ml | BSB 6957 | Concentrate | 0.5 ml | BSB-3740-05 | Concentrate | 0.5 ml |
| BSB 6175 | Concentrate | 1.0 ml | BSB 6958 | Concentrate | 1.0 ml | BSB-3740-1 | Concentrate | 1.0 ml |
| BSB 6176 | Control Slides | 5 | BSB 6959 | Control Slides | 5 | BSB-3740-CS | Control Slides | 5 |

Mycobacterium tuberculosis, RPab




IHC of Mycobacterium Tuberculosis on a FFPE Infected Lung Tissue

Mycobacterium tuberculosis is a pathogenic bacterial species of the Mycobacteriaceae family and the causative agent of most cases of tuberculosis. M. tuberculosis has an unusual, waxy coating on its cell surface (primarily due to the presence of mycolic acid), which makes the cells impervious to Gram staining; M. tuberculosis can appear Gram negative and Gram positive in clinical settings. Antibiotic resistant strains of mycobacterium tuberculosis have developed resistance to more than one TB drug, due to mutations in their genes.

M. tuberculosis is characterized by caseating granulomas containing Langhans giant cells, which have a "horseshoe" pattern of nuclei. Cells are often seen wrapped together, due to the presence of fatty acids in the cell wall that stick together. This appearance is referred to as chording, like strands of chord that make up a rope. The clinical and histological criteria used to diagnose lymphadenitis caused by Mycobacterium tuberculosis complex organisms have poor specificity. Acid-fast staining and culture have low sensitivity and specificity. The diagnosis of tuberculosis by immunohistochemistry can be used to detect the mycobacterial antigen on formalin-fixed tissue biopsies and it's consider fast, sensitive, and a highly specific method for establishing the etiological diagnosis of tuberculosis in histologic specimens.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Infected Tissue
LOCALIZATION: Cell Wall
SPECIES REACTIVITY: Human

Myelin Basic Protein, RPab




IHC of Myelin Basic Protein on a FFPE Brain Tissue

Myelin Basic Protein (MBP) is a protein believed to be important in the process of myelination of nerves in the central nervous system (CNS). The pool of MBP in the central nervous system is very diverse, with several splice variants being expressed and a large number of post-translational modifications on the protein, which include phosphorylation, methylation, deamidation and citrullination.

MBP has been demonstrated in Neuromas, Neurofibromas, and Neurogenic Sarcomas. However, other spindle-cell neoplasms do not stain with this antibody. Immunoreactivity for MBP in Granular-cell Tumors strengthens the concept of a Schwann-cell derivation of these lesions. Unlike other nervous system proteins such as GFAP and S-100, MBP has not been demonstrated in melanocytes or tumors derived from them.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Brain, Neuroblastoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog

Myelin Basic Protein, RMAb



IHC of Myelin Basic Protein on a FFPE Brain Tissue

Myelin Basic Protein (MBP) is a protein believed to be important in the process of myelination of nerves in the central nervous system (CNS). The pool of MBP in the central nervous system is very diverse, with several splice variants being expressed and a large number of post-translational modifications on the protein, which include phosphorylation, methylation, deamidation and citrullination.

MBP has been demonstrated in Neuromas, Neurofibromas, and Neurogenic Sarcomas. However, other spindle-cell neoplasms do not stain with this antibody. Immunoreactivity for MBP in Granular-cell Tumors strengthens the concept of a Schwann-cell derivation of these lesions. Unlike other nervous system proteins such as GFAP and S-100, MBP has not been demonstrated in melanocytes or tumors derived from them.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP207
ISOTYPE: IgG
CONTROL: Brain, Neuroblastoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2992 | Tinto Predilute | 3.0 ml |
| BSB 2993 | Tinto Predilute | 7.0 ml |
| BSB 2994 | Tinto Predilute | 15.0 ml |
| BSB 2995 | Concentrate | 0.1 ml |
| BSB 2996 | Concentrate | 0.5 ml |
| BSB 2997 | Concentrate | 1.0 ml |
| BSB 2998 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5778 | Tinto Predilute | 3.0 ml |
| BSB 5779 | Tinto Predilute | 7.0 ml |
| BSB 5780 | Tinto Predilute | 15.0 ml |
| BSB 5781 | Concentrate | 0.1 ml |
| BSB 5782 | Concentrate | 0.5 ml |
| BSB 5783 | Concentrate | 1.0 ml |
| BSB 5784 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6960 | Tinto Predilute | 3.0 ml |
| BSB 6961 | Tinto Predilute | 7.0 ml |
| BSB 6962 | Tinto Predilute | 15.0 ml |
| BSB 6963 | Concentrate | 0.1 ml |
| BSB 6964 | Concentrate | 0.5 ml |
| BSB 6965 | Concentrate | 1.0 ml |
| BSB 6966 | Control Slides | 5 |

Myeloperoxidase, RPaB




IHC of Myeloperoxidase on a FFPE Bone Marrow Tissue

Myeloperoxidase (MPO) is a peroxidase enzyme most abundantly present in neutrophil granulocytes. It is a lysosomal protein stored in azurophilic granules of the neutrophil. MPO has a heme pigment, which causes its green color in secretions rich in neutrophils, such as pus and some forms of mucus. Historically, immunohistochemical staining for myeloperoxidase was used in the diagnosis of Acute Myeloid Leukemia to demonstrate that the leukemic cells were derived from the myeloid lineage. Myeloperoxidase staining is still important in the diagnosis of Extramedullary Leukemia or Chloroma.

Myeloperoxidase detects granulocytes and monocytes in blood and precursors of granulocytes in the bone marrow. This antibody can detect myeloid cell populations of the bone marrow as well as in other sites.

Myeloperoxidase, RMab




IHC of Myeloperoxidase on a FFPE Bone Marrow Tissue

Myeloperoxidase (MPO) is a peroxidase enzyme most abundantly present in neutrophil granulocytes. It is a lysosomal protein stored in azurophilic granules of the neutrophil. MPO has a heme pigment, which causes its green color in secretions rich in neutrophils, such as pus and some forms of mucus. Historically, immunohistochemical staining for myeloperoxidase was used in the diagnosis of Acute Myeloid Leukemia to demonstrate that the leukemic cells were derived from the myeloid lineage. Myeloperoxidase staining is still important in the diagnosis of Extramedullary Leukemia or Chloroma.

Myeloperoxidase detects granulocytes and monocytes in blood and precursors of granulocytes in the bone marrow. This antibody can detect myeloid cell populations of the bone marrow as well as in other sites.

MyoD1, RMab



IHC of MyoD1 on a FFPE Rhabdomyosarcoma Tissue

MyoD1 belongs to a family of proteins known as myogenic regulatory factors (MRFs) and has a key role in regulating muscle differentiation. These bHLH (basic helix loop helix) transcription factors act sequentially in myogenic differentiation. MyoD1 is expressed in activated satellite cells, but not in quiescent satellite cells. In development, MyoD1 commits mesoderm cells to a skeletal lineage, and then regulates that process. It may also play a role in muscle repair.

In abnormal tissues, MyoD1 labels tumor cells in Rhabdomyosarcoma and is one of the earliest markers of myogenic commitment.


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Bone Marrow
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP151
ISOTYPE: IgG
CONTROL: Bone Marrow
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP212
ISOTYPE: IgG
CONTROL: Fetal Muscle, RhabdomyoSarcoma & Soft Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 5785 | Tinto Predilute | 3.0 ml | BSB 6967 | Tinto Predilute | 3.0 ml | BSB 6974 | Tinto Predilute | 3.0 ml |
| BSB 5786 | Tinto Predilute | 7.0 ml | BSB 6968 | Tinto Predilute | 7.0 ml | BSB 6975 | Tinto Predilute | 7.0 ml |
| BSB 5787 | Tinto Predilute | 15.0 ml | BSB 6969 | Tinto Predilute | 15.0 ml | BSB 6976 | Tinto Predilute | 15.0 ml |
| BSB 5788 | Concentrate | 0.1 ml | BSB 6970 | Concentrate | 0.1 ml | BSB 6977 | Concentrate | 0.1 ml |
| BSB 5789 | Concentrate | 0.5 ml | BSB 6971 | Concentrate | 0.5 ml | BSB 6978 | Concentrate | 0.5 ml |
| BSB 5790 | Concentrate | 1.0 ml | BSB 6972 | Concentrate | 1.0 ml | BSB 6979 | Concentrate | 1.0 ml |
| BSB 5791 | Control Slides | 5 | BSB 6973 | Control Slides | 5 | BSB 6980 | Control Slides | 5 |

Myogenin, MAb




IHC of Myogenin on a FFPE Rhabdomyosarcoma Tissue

Myogenin is a transcription factor active in muscles. In particular, it is a myogenic regulatory factor. Myogenin is a member of a family of myogenic regulatory genes, which includes MyoD, myf5 and MRF4. These genes encode a set of transcription factors which are essential for muscle development. Expression of myogenin is restricted to cells of skeletal-muscle origin. It is, therefore, a useful marker for tumors of the muscle lineage, being strongly expressed in Alveolar Rhabdomyosarcomas.

Anti-myogenin labels the nuclei of myoblasts in developing muscle tissue, and is expressed in tumor cell nuclei of Rhabdomyosarcoma and some Leiomyosarcomas. Positive nuclear staining may occur in Wilm's Tumor.

Myoglobin, RMAb




IHC of Myoglobin on a FFPE Rhabdomyosarcoma Tissue

Myoglobin is a single-chain globular protein of 153 amino acids, containing a heme (iron-containing porphyrin) prosthetic group in the center around which the remaining apoprotein folds. With a molecular weight of 16.7 kDa, Myoglobin is the primary oxygen-carrying pigment of muscle tissues.

Immunostaining with Myoglobin provides a specific, sensitive and practical procedure for the identification of Rhabdomyosarcoma. Since myoglobin is found exclusively in skeletal and cardiac muscle, it may be used to distinguish Rhabdomyosarcoma from other soft-tissue tumors. Myoglobin staining is also useful when demonstrating rhabdomyoblastic differentiation in other tumors, e.g., Neurogenic Sarcomas and Malignant Mixed Mesodermal Tumors of the uterus and ovary.

Myoglobin, RMAb



IHC of Myoglobin on a FFPE Rhabdomyosarcoma Tissue

Myoglobin is a single-chain globular protein of 153 amino acids, containing a heme (iron-containing porphyrin) prosthetic group in the center around which the remaining apoprotein folds. With a molecular weight of 16.7 kDa, Myoglobin is the primary oxygen-carrying pigment of muscle tissues.

Immunostaining with Myoglobin provides a specific, sensitive and practical procedure for the identification of Rhabdomyosarcoma. Since myoglobin is found exclusively in skeletal and cardiac muscle, it may be used to distinguish Rhabdomyosarcoma from other soft-tissue tumors. Myoglobin staining is also useful when demonstrating rhabdomyoblastic differentiation in other tumors, e.g., Neurogenic Sarcomas and Malignant Mixed Mesodermal Tumors of the uterus and ovary.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: F5D
ISOTYPE: IgG1/K
CONTROL: Rhabdomyosarcoma & Soft Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Dog, Rat, Mouse

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: BSB-104
ISOTYPE: IgG1
CONTROL: Skeletal Muscle Tissue
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP87
ISOTYPE: IgG
CONTROL: Skeletal Muscle Tissue
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5792 | Tinto Predilute | 3.0 ml |
| BSB 5793 | Tinto Predilute | 7.0 ml |
| BSB 5794 | Tinto Predilute | 15.0 ml |
| BSB 5795 | Concentrate | 0.1 ml |
| BSB 5796 | Concentrate | 0.5 ml |
| BSB 5797 | Concentrate | 1.0 ml |
| BSB 5798 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3384 | Tinto Predilute | 3.0 ml |
| BSB 3385 | Tinto Predilute | 7.0 ml |
| BSB 3386 | Tinto Predilute | 15.0 ml |
| BSB 3387 | Concentrate | 0.1 ml |
| BSB 3388 | Concentrate | 0.5 ml |
| BSB 3389 | Concentrate | 1.0 ml |
| BSB 3390 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6981 | Tinto Predilute | 3.0 ml |
| BSB 6982 | Tinto Predilute | 7.0 ml |
| BSB 6983 | Tinto Predilute | 15.0 ml |
| BSB 6984 | Concentrate | 0.1 ml |
| BSB 6985 | Concentrate | 0.5 ml |
| BSB 6986 | Concentrate | 1.0 ml |
| BSB 6987 | Control Slides | 5 |

Myosin Smooth Muscle, MMab




IHC of Myosin Smooth Muscle Heavy Chain on a FFPE Appendix Tissue

Nanog, RMab



IHC of Nanog on an FFPE Seminoma Tissue

Napsin A, MMab



IHC of Napsin A on a FFPE Lung Carcinoma Tissue

Myosins are a large family of motor proteins found in eukaryotic tissues. They are responsible for actin-based motility. Smooth Muscle Myosin, Heavy Chain is a cytoplasmic structural protein that is a major component of the contractile apparatus of the smooth muscle cells, as well as a myoepithelium-associated protein.

SMM-H24 is a mouse monoclonal antibody to Smooth Muscle Myosin, Heavy Chain that reacts with human visceral and vascular smooth muscle cells. The antibody also reacts with human myoepithelial cells. It is very helpful in distinguishing between benign sclerosing breast lesions and infiltrating Carcinomas in difficult cases, since it strongly stains the myoepithelial layer in the benign lesions while it is negative in the infiltrating Carcinomas.

Homeobox protein Nanog is a transcriptional factor that helps embryonic stem cells (ESCs) maintain pluripotency by suppressing cell determination factors. In humans, this protein is encoded by the NANOG gene. Nanog is thought to function in concert with other factors such as Oct-4 and SOX2 to establish ESC identity. Nanog is highly expressed in cancer stem cells and thus may function as an oncogene to promote carcinogenesis. High expression of Nanog correlates with poor survival in cancer patients.

Nanog is highly and specifically expressed in carcinoma in situ (CIS), embryonal carcinomas, and seminomas, but not in teratomas and yolk sac tumors. Additionally, it has been reported that Human embryonic stem cell genes OCT4, NANOG, STELLAR, and GDF3 are expressed in both seminoma and breast carcinoma. Positive Nanog expression is significantly associated with high-grade ovarian serous carcinoma and is absent in benign, borderline, and low-grade serous lesions. A study suggests the expression of Nanog exhibiting cellular shuttling behavior and increasing stromal distribution during the progression of cervical cancer. Recently it was suggested that Nanog overexpression, a hazard factor of differentiation, lymph node metastasis, and tumor size, may predicate decreased overall survival (OS) and disease-free survival (DFS) for lung cancer.

The activation peptides of aspartic proteinases play a role as inhibitors of the active site. These peptide segments, or pro-parts, are deemed important for correct folding, targeting, and control of the activation of aspartic proteinase zymogens. The pronapsin A gene is expressed predominantly in lung and kidney. Its translation product is predicted to be a fully functional glycosylated aspartic proteinase precursor containing an RGD motif and an additional 18 residues at its C-terminus.

In normal tissue, anti-Napsin A labels type II pneumocytes in adult lung and epithelial cells in kidney tissues. In abnormal tissues, Napsin A is a useful marker for lung adenocarcinoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-17
ISOTYPE: IgG1/K
CONTROL: Intestine, Breast, Appendix
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP225
ISOTYPE: IgG
CONTROL: Testis, Cervix, Seminoma, Embryonal Carcinoma, TCC
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-112
ISOTYPE: IgG1/K
CONTROL: Kidney, Lung, Lung Carcinoma, Renal Cell Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5924 | Tinto Predilute | 3.0 ml |
| BSB 5925 | Tinto Predilute | 7.0 ml |
| BSB 5926 | Tinto Predilute | 15.0 ml |
| BSB 5927 | Concentrate | 0.1 ml |
| BSB 5928 | Concentrate | 0.5 ml |
| BSB 5929 | Concentrate | 1.0 ml |
| BSB 5930 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3581 | Tinto Predilute | 3.0 ml |
| BSB 3582 | Tinto Predilute | 7.0 ml |
| BSB 3583 | Tinto Predilute | 15.0 ml |
| BSB 3584 | Concentrate | 0.1 ml |
| BSB 3585 | Concentrate | 0.5 ml |
| BSB 3586 | Concentrate | 1.0 ml |
| BSB 3587 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3392 | Tinto Predilute | 3.0 ml |
| BSB 3393 | Tinto Predilute | 7.0 ml |
| BSB 3394 | Tinto Predilute | 15.0 ml |
| BSB 3395 | Concentrate | 0.1 ml |
| BSB 3396 | Concentrate | 0.5 ml |
| BSB 3397 | Concentrate | 1.0 ml |
| BSB 3398 | Control Slides | 5 |

Napsin A, RMab




IHC of Napsin A on a FFPE Lung Adenocarcinoma Tissue

Nestin, RMab



IHC of Nestin on a FFPE Kidney Tissue

NeuN, MAb



IHC of NeuN on a FFPE Brain Tissue

The activation peptides of aspartic proteinases play a role as inhibitors of the active site. These peptide segments, or pro-parts, are deemed important for correct folding, targeting, and control of the activation of aspartic proteinase zymogens. The pronapsin A gene is expressed predominantly in lung and kidney. Its translation product is predicted to be a fully functional glycosylated aspartic proteinase precursor containing an RGD motif and an additional 18 residues at its C-terminus.

In normal tissue, anti-Napsin A labels type II pneumocytes in adult lung and epithelial cells in kidney tissues. In abnormal tissues, Napsin A is a useful marker for lung adenocarcinoma.

Nestin is a type VI intermediate filament protein; they are expressed mostly in nerve cells where they are implicated in the radial growth of the axon. Nestin is expressed in dividing cells during the early stages of development in the Central Nervous System (CNS), Peripheral Nervous System (PNS) and in myogenic and other tissues. Nestin is expressed by many types of cells during development, although its expression is usually transient and does not persist into adulthood. Nestin is however expressed in the neuronal precursor cells of the subgranular zone in adult organisms. Its expression is also reinduced in the adult during pathological situations, such as the formation of the glial scar after CNS injury and during regeneration of injured muscle tissue.

It has been reported that Nestin expression is significantly increased in melanoma and correlated with more advanced stages of the disease. It has also been reported in tumors of the CNS, including astrocytoma, ependymoma, oligodendroglioma, glioblastoma, and primitive neuroectodermal tumors, as well as in carcinomas such as prostatic adenocarcinoma, pancreatic ductal carcinoma, thyroid carcinoma, and in mesenchymal tumors. In breast carcinoma subtypes, Nestin is highly expressed in basal breast cancer but not in the HER2 subtype or luminal epithelial phenotype. In normal skin, Nestin is expressed in endothelial cells and the bulge area of hair follicles.

NeuN (Feminizing Locus on X-3, Fox-3, or Hexaribonucleotide Binding Protein-3) is a neuron-specific protein that is present in most Central Nervous System (CNS) and Peripheral Nervous System (PNS) neuronal cell types. NeuN protein distributions are restricted to neuronal nuclei, perikarya and some proximal neuronal processes in both fetal and adult brain. However, some neurons fail to be recognized by NeuN at all ages, such as INL retinal cells, Cajal-Retzius cells, Purkinje cells, inferior olivary and dentate nucleus neurons, and sympathetic ganglion cells.

NeuN is widely used to label neurons since the vast majority of neurons are strongly positive. NeuN immunoreactivity becomes obvious as neurons mature, typically after they have downregulated expression of Doublecortin, a marker seen in the earliest stages of neuronal development.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP205

ISOTYPE: IgG

CONTROL: Kidney, Lung, Lung Carcinoma, Renal Cell Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP287

ISOTYPE: IgG

CONTROL: Kidney, Breast, Adrenal, Myometrium, Liver Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: A60

ISOTYPE: IgG1

CONTROL: Brain

LOCALIZATION: Nuclear


SPECIES REACTIVITY: Human, Avian, Chicken, Ferret, Mouse, Pig, Rat, Salamander

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6988 | Tinto Predilute | 3.0 ml |
| BSB 6989 | Tinto Predilute | 7.0 ml |
| BSB 6990 | Tinto Predilute | 15.0 ml |
| BSB 6991 | Concentrate | 0.1 ml |
| BSB 6992 | Concentrate | 0.5 ml |
| BSB 6993 | Concentrate | 1.0 ml |
| BSB 6994 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2000 | Tinto Predilute | 3.0 ml |
| BSB 2001 | Tinto Predilute | 7.0 ml |
| BSB 2002 | Tinto Predilute | 15.0 ml |
| BSB 2003 | Concentrate | 0.1 ml |
| BSB 2004 | Concentrate | 0.5 ml |
| BSB 2005 | Concentrate | 1.0 ml |
| BSB 2006 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2007 | Tinto Predilute | 3.0 ml |
| BSB 2008 | Tinto Predilute | 7.0 ml |
| BSB 2009 | Tinto Predilute | 15.0 ml |
| BSB 2010 | Concentrate | 0.1 ml |
| BSB 2011 | Concentrate | 0.5 ml |
| BSB 2012 | Concentrate | 1.0 ml |
| BSB 2013 | Control Slides | 5 |

NeuN, RMAb



IHC of NeuN on a FFPE Brain Tissue

NeuN (Feminizing Locus on X-3, Fox-3, or Hexaribonucleotide Binding Protein-3) is a neuron-specific protein that is present in most central nervous system (CNS) and peripheral nervous system (PNS) neuronal cell types. NeuN protein distributions are restricted to neuronal nuclei, perikarya and some proximal neuronal processes in both fetal and adult brains. However, some neurons fail to be recognized by NeuN at all ages, such as inner nuclear layer retinal cells, Cajal-Retzius cells, Purkinje cells, inferior olivary and dentate nucleus neurons, and sympathetic ganglion cells.

NeuN is widely used to label neurons since the vast majority of neurons are strongly positive. NeuN immunoreactivity becomes obvious as neurons mature, typically after they have downregulated expression of doublecortin, a marker seen in the earliest stages of neuronal development. NeuN was detected in most of the amyloid bodies, is considered a marker of neuronal differentiation in brain tumor and has been found in all major subtypes except pilocytic astrocytoma. NeuN IHC have demonstrated NeuN immunoreactivity in 56% of epithelial neuroendocrine carcinomas (ENEC) (19/34): 4 of 7 (57%) grade 1 ENECs (Carcinoid), 4 of 5 (90%) grade 2 ENECs (atypical carcinoid), and 11 of 22 (50%) grade 3 ENECs (small and large cell neuroendocrine carcinoma).

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-NeuN


ISOTYPE: IgG

CONTROL: Brain

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

Neurofilament, MMAb



IHC of Neurofilament on a FFPE Brain Tissue

Neurofilaments are the Type IV family of intermediate filaments that are found in high concentrations along the axons of vertebrate neurons.

Neurofilament antibody stains an antigen localized in a number of neural, neuroendocrine and endocrine tumors. Neuromas, Ganglioneuromas, Gangliogliomas, Ganglioneuroblastomas and Neuroblastomas stain positively for neurofilament. Neurofilaments are also present in Paragangliomas and Adrenal and Extra-Adrenal Pheochromocytomas. Carcinoids, Neuroendocrine Carcinomas of the Skin, and Oat Cell Carcinomas of the Lung also express neurofilament.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 2F11


ISOTYPE: IgG1/K

CONTROL: Brain

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Dog, Rat, Mouse, Horse

Neurofilament, RMAb



IHC of Neurofilament on a FFPE Brain Tissue

Neurofilaments are the Type IV family of intermediate filaments that are found in high concentrations along the axons of vertebrate neurons.

Neurofilament antibody identifies an antigen localized in a number of neural, neuroendocrine and endocrine tumors. Neuromas, Ganglioneuromas, Gangliogliomas, Ganglioneuroblastomas and Neuroblastomas stain positively for neurofilament. Neurofilaments are also present in Paragangliomas and Adrenal and Extra-Adrenal Pheochromocytomas. Carcinoids, Neuroendocrine Carcinomas of the Skin, and Oat Cell Carcinomas of the Lung also express neurofilament.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP79

ISOTYPE: IgG

CONTROL: Brain

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3781-3 | Tinto Predilute | 3.0 ml |
| BSB-3781-7 | Tinto Predilute | 7.0 ml |
| BSB-3781-15 | Tinto Predilute | 15.0 ml |
| BSB-3781-01 | Concentrate | 0.1 ml |
| BSB-3781-05 | Concentrate | 0.5 ml |
| BSB-3781-1 | Concentrate | 1.0 ml |
| BSB-3781-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5813 | Tinto Predilute | 3.0 ml |
| BSB 5814 | Tinto Predilute | 7.0 ml |
| BSB 5815 | Tinto Predilute | 15.0 ml |
| BSB 5816 | Concentrate | 0.1 ml |
| BSB 5817 | Concentrate | 0.5 ml |
| BSB 5818 | Concentrate | 1.0 ml |
| BSB 5819 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2014 | Tinto Predilute | 3.0 ml |
| BSB 2015 | Tinto Predilute | 7.0 ml |
| BSB 2016 | Tinto Predilute | 15.0 ml |
| BSB 2017 | Concentrate | 0.1 ml |
| BSB 2018 | Concentrate | 0.5 ml |
| BSB 2019 | Concentrate | 1.0 ml |
| BSB 2020 | Control Slides | 5 |

NGFR, MAb




IHC of NGFR on a FFPE Breast Fibroadenoma Tissue

NGFR (Nerve Growth Factor Receptor), also termed p75 or CD271, is the low-affinity NGFR (LNGFR) which binds NGF and other neurotrophins, including BDNF, NT3 and NT4/5 with similar low-affinity. NGFR p75 is a 75 kD transmembrane glycoprotein that is mainly expressed in Schwann cells and neurons and in a variety of non-neuronal cells. NGFR p75 is necessary for regulating neuronal growth, migration, differentiation and cell death during development of the central and peripheral nervous system. NGFR p75 plays a central role in the regulation of cell number by apoptosis in the developing CNS. During early development, activation of NGFR p75 by NGF induces apoptotic cell death in some neuronal cells, probably through activation of the sphingomyelinase/ceramide pathway, the ICE-like proteases and the JNK pathway. CD271 has recently been described as being expressed in mesenchymal stem cells (bone marrow stromal cells).

NGFR is expressed not only in sympathetic and sensory neurons, but also in various neural crest cell or tumor derivatives such as melanocytes, Melanomas, Neuroblastomas, Pheochromocytomas, Neurofibromas, and neurotized nevi (Type C melanocytes). It is now apparent that expression of NGFR is ubiquitous and not limited to the nervous system, being expressed in mature non-neural cells such as perivascular cells, dental pulp cells, lymphoid follicular dendritic cells, basal epithelium of oral mucosa and hair follicles, prostate basal cells and myoepithelial cells. Studies in Prostate and Urothelial Cancer suggest that NGFR may act as a tumor suppressor, negatively regulating cell growth and proliferation. NGFR labels the myoepithelial cells of breast ducts and intralobular fibroblasts of breast ducts and, thus, aids in the diagnosis of malignancy in the breast.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-18
ISOTYPE: IgG1
CONTROL: Brain, Breast, Prostate, Neuroblastoma, CNS Tumor
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

NKX2.2, R Mab




IHC of NKX2.2 on a FFPE Pancreas Tissue

Homeobox protein NKX2.2 is a protein encoded by the NKX2-2 gene. NKX2.2 is a homeodomain-containing transcription factor that plays a critical role in neuroendocrine/glial differentiation. NKX2.2 is expressed in the developing forebrain and spinal cord. Functionally, the transcription factor is thought to be involved with neuronal developing, patterning, and fate specification of neurons and oligodendrocytes.

NKX2.2 was recently reported as a valuable marker for Ewing's sarcoma. The vast majority of Ewing's sarcomas (85%) harbor a chromosomal translocation, most commonly t(11;22)(q24;q12) encoding an aberrant EWS-FLI transcription factor. NKX2.2 expression is tightly correlated with EWS-FLI expression, a critical downstream target that is required for the cancerous behavior of Ewing's sarcoma. While CD99 is the classical marker for Ewing's sarcoma diagnosis, it is relatively nonspecific. In addition to Ewing's sarcoma, CD99 is also expressed on lymphocytes, hematopoietic cells, endothelial cells and a variety of tumors. In contrast, NKX2.2 labels 93% of Ewing's sarcoma and only a small subset (14/130) of non-Ewing tumors, demonstrating a sensitivity of 93% and specificity of 89%. Staining with NKX2.2 can aid in the differential diagnosis of small round cell tumors.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP336
ISOTYPE: IgG
CONTROL: Pancreas, Brain, Pituitary, Colon, Ewing's Sarcoma & Soft Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

NKX3.1, R Mab



IHC of NKX3.1 on a FFPE Prostatic Carcinoma Tissue

Homeobox protein NKX3.1, also known as BAPX2 and NKX3A is a protein that in humans is encoded by the NKX3.1 gene located on chromosome 8p. NKX3.1 is a prostatic tumor suppressor gene, which is an androgen-regulated, prostate-specific homeobox gene whose expression is predominantly localized in the prostate epithelium. It is a negative regulator of epithelial cell growth in prostate tissue. Loss of NKX3A protein expression is a common finding in human prostate carcinomas and prostatic intraepithelial neoplasia. NKX3-1 expression is seen in prostate epithelium, testis, ureter, and pulmonary bronchial mucous glands.


NKX3-1 has been established as a marker for identifying metastatic tumors. NKX3.1-positive prostate carcinoma cells exhibit nuclear staining. Additionally, most cases of Urothelial Carcinoma have been found to be negative for NKX3.1 and may be helpful to distinguish between high grade Prostate Adenocarcinoma and high grade Infiltrating Urothelial Carcinoma. NKX3.1 has also been found to be expressed in Invasive Ductal Carcinomas (IDC) and Invasive Lobular Carcinomas (ILC) of the breast. NKX3.1 expression is limited to ER, PR, and AR positive carcinomas and is more frequently expressed in ILC than IDC. NKX3.1 has a high specificity and sensitivity for prostate adenocarcinomas and can be used to help distinguish between Prostate Carcinoma and Urothelial Carcinomas.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP356
ISOTYPE: IgG
CONTROL: Prostate, Prostate Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6289 | Tinto Predilute | 3.0 ml |
| BSB 6290 | Tinto Predilute | 7.0 ml |
| BSB 6291 | Tinto Predilute | 15.0 ml |
| BSB 6292 | Concentrate | 0.1 ml |
| BSB 6293 | Concentrate | 0.5 ml |
| BSB 6294 | Concentrate | 1.0 ml |
| BSB 6295 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3106 | Tinto Predilute | 3.0 ml |
| BSB 3107 | Tinto Predilute | 7.0 ml |
| BSB 3108 | Tinto Predilute | 15.0 ml |
| BSB 3109 | Concentrate | 0.1 ml |
| BSB 3110 | Concentrate | 0.5 ml |
| BSB 3111 | Concentrate | 1.0 ml |
| BSB 3112 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3113 | Tinto Predilute | 3.0 ml |
| BSB 3114 | Tinto Predilute | 7.0 ml |
| BSB 3115 | Tinto Predilute | 15.0 ml |
| BSB 3116 | Concentrate | 0.1 ml |
| BSB 3117 | Concentrate | 0.5 ml |
| BSB 3118 | Concentrate | 1.0 ml |
| BSB 3119 | Control Slides | 5 |

NKX3.1, RMAb

IHC of NKX3.1 on a FFPE Prostatic Carcinoma Tissue

NPM1/B23, MMAb

IHC of NPM1/B23 on an FFPE Renal Cell Carcinoma Tissue

NRAS, RMAb

IHC of NRAS on a FFPE Melanoma Tissue

Homeobox protein NKX3.1, also known as BAPX2 and NKX3A is a protein that in humans is encoded by the NKX3.1 gene located on chromosome 8p. NKX3.1 is a prostatic tumor suppressor gene, which is an androgen-regulated, prostate-specific homeobox gene whose expression is predominantly localized in the prostate epithelium. It is a negative regulator of epithelial cell growth in prostate tissue. Loss of NKX3A protein expression is a common finding in human prostate carcinomas and prostatic intraepithelial neoplasia. NKX3-1 expression is seen in prostate epithelium, testis, ureter, and pulmonary bronchial mucous glands.

NKX3-1 has been established as a marker for identifying metastatic tumors. In a study the sensitivity for identifying metastatic prostatic adenocarcinomas was 98.6% for NKX3.1, 94.2% for prostate specific antigen and 98.6% for prostatic acid phosphatase and a specificity of 99.7% for NKX3.1. NKX3.1-positive prostate carcinoma cells exhibit nuclear staining. Additionally, most cases of urothelial carcinoma have been found to be negative for NKX3.1 and may be helpful to distinguish between high grade prostate adenocarcinoma and high grade Infiltrating urothelial carcinoma. NKX3.1 expression is limited to ER, PR, and AR positive carcinomas and is more frequently expressed in ILC than IDC. NKX3.1 has a high specificity and sensitivity for prostate adenocarcinomas and can be used to help distinguish between prostate carcinoma and urothelial carcinomas.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM430
ISOTYPE: IgG
CONTROL: Prostate, Prostate Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Nucleophosmin 1 (NPM1), also known as nucleolar phosphoprotein B23, is a protein that in humans is encoded by the NPM1 gene. The nuclear protein B23 (also referred to as nucleophosmin) is involved in ribosomal assembly, rRNA transport, centrosome duplication, maintenance of genomic stability, and embryonic development.

NPM1 gene is up-regulated, mutated and chromosomally translocated in many tumor types. Chromosomal aberrations involving NPM1 were found in patients with non-Hodgkin lymphoma, Acute Promyelocytic Leukemia, Myelodysplastic Syndrome, and Acute Myelogenous Leukemia. NPM1 is a nuclear protein. In approximately 50% to 60% of cytogenetically normal Acute Myeloid Leukemia (AML), NPM1 is mutated and localized in the cytoplasm. Both wild type and mutant NPM1 can be detected by immunohistochemistry (IHC). The expression of NPM1 is heterogeneous in gastric tumors. NPM1 down-regulation may have a role in gastric carcinogenesis and may help in the selection of anticancer treatment strategies. NPM1 has a critical role in the regulation of colon cancer cells migration and invasion and it may serve as a potential marker for the prognosis of colon cancer patients.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-124
ISOTYPE: IgG1/k
CONTROL: Breast, Cervix, Testis, Pancreas, RCC, TCC
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

NRAS is an enzyme that in humans is encoded by the NRAS gene. The NRAS oncogene is a member of the Ras gene family. It is mapped on chromosome 1. The Ras genes have GTP/GDP binding and GTPase activity, and their normal function may be as G-like regulatory proteins involved in the normal control of cell growth. Mutations which change amino acid residues 12, 13 or 61 activate the potential of N-Ras to transform cultured cells and are implicated in a variety of human tumors.

The Q61R is the most common NRAS mutation found in melanoma that is thought to occur due to UV and radiation exposure. The Q61R mutation results in an amino acid substitution at position 61 in NRAS, from a glutamine (Q) to an arginine (R). It is a commonly acquired mutation found in primary melanomas from melanoma-prone families that have known germline CDKN2A mutations, and it is associated with loss of the wild-type NRAS allele. It is also detected in patients with congenital melanocytic nevi (CMN), which can develop into malignant melanoma. The NRAS Q61R is the most common oncogenic event in encapsulated thyroid tumors of follicular cell origin with high-grade features. 20% of adrenocortical tumors solely carry this RAS alteration; however, adrenocortical hyperplasia does not, suggesting NRAS Q61R is an oncogenic event driving adrenocortical tumor development.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-NRAS
ISOTYPE: IgG
CONTROL: Small Intestine, Melanoma, ASTMA2, Lung Carcinoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3785-3 | Tinto Predilute | 3.0 ml |
| BSB-3785-7 | Tinto Predilute | 7.0 ml |
| BSB-3785-15 | Tinto Predilute | 15.0 ml |
| BSB-3785-01 | Concentrate | 0.1 ml |
| BSB-3785-05 | Concentrate | 0.5 ml |
| BSB-3785-1 | Concentrate | 1.0 ml |
| BSB-3785-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3588 | Tinto Predilute | 3.0 ml |
| BSB 3589 | Tinto Predilute | 7.0 ml |
| BSB 3590 | Tinto Predilute | 15.0 ml |
| BSB 3591 | Concentrate | 0.1 ml |
| BSB 3592 | Concentrate | 0.5 ml |
| BSB 3593 | Concentrate | 1.0 ml |
| BSB 3594 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2663 | Tinto Predilute | 3.0 ml |
| BSB 2664 | Tinto Predilute | 7.0 ml |
| BSB 2665 | Tinto Predilute | 15.0 ml |
| BSB 2666 | Concentrate | 0.1 ml |
| BSB 2667 | Concentrate | 0.5 ml |
| BSB 2668 | Concentrate | 1.0 ml |
| BSB 2669 | Control Slides | 5 |

NSE, MMab




IHC of NSE on a FFPE Pancreas Tissue

OCT-2, RMAb



IHC of OCT-2 on a FFPE Lymphoma Tissue

OCT-4, RMAb



IHC of OCT-4 on a FFPE Seminoma Tissue

Neuron-Specific Enolase (NSE, Enolase 2) is a human gene. It makes a phosphopyruvate hydratase. This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.

NSE is present in high concentration in neurons and in central and peripheral neuroendocrine cells; therefore, NSE reacts with cells of neural and neuroendocrine lineage. If neoplastic cells coexpress keratins and NSE, neuroendocrine differentiation is probable. However, neural tumors that do not express keratin, and show no staining with NSE, would not exclude neural or neuroendocrine differentiation. Thus, detection of neural and neuroendocrine lineage requires the use of panels which include NSE and other markers such as keratin, chromogranin, synaptophysin and neurofilament.

Octamer transcription factor-2 (OCT-2) possesses a leucine zipper domain and belongs to the POU family of transcription factors. It binds to the octamer motif (5-ATTTCAT-3), activates immunoglobulin gene expression and regulates transcription in a number of tissues. OCT-2 is important for the expression of B cell specific genes, such as CD20 and CRISP-3. OCT-2 is expressed in mature B cells, predominantly germinal center B cells.

The OCT-2 antibody labels various B cell lymphomas with strong expression in germinal center-derived lymphomas.

OCT-4 (octamer-binding transcription factor 4) also known as POU5F1 (POU domain, class 5, transcription factor 1) is a protein that in humans is homeodomain transcription factor of the POU family. This protein is critically involved in the self-renewal of undifferentiated embryonic stem cells. Clear cell carcinoma may enter the differential diagnosis of dysgerminoma as both may grow in nests or tubules, contain clear cells, and have a prominent inflammatory infiltrate (lymphocytes in dysgerminoma and plasma cells in clear cell carcinoma).

Expression of the OCT-4 antibody is potentially correlated with tumorigenesis and can affect some aspects of tumor behavior such as tumor recurrence or resistance to therapies. OCT-4 is expressed in undifferentiated pluripotency cells, germ cells in ovary and testes. OCT-4 is a sensitive and specific marker for germ cell tumors. It is consistently detected in carcinoma in situ/ gonadoblastoma, seminomas, germinoma, dysgerminoma, and embryonal carcinoma but not in the differentiated components of nonseminomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-94

ISOTYPE: IgG1/K

CONTROL: Pancreas, Brain, Pituitary, Adrenal, Thyroid

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP115

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP143

ISOTYPE: IgG

CONTROL: Seminoma, Dysgerminoma, Testis Carcinoma


LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5820 | Tinto Predilute | 3.0 ml |
| BSB 5821 | Tinto Predilute | 7.0 ml |
| BSB 5822 | Tinto Predilute | 15.0 ml |
| BSB 5823 | Concentrate | 0.1 ml |
| BSB 5824 | Concentrate | 0.5 ml |
| BSB 5825 | Concentrate | 1.0 ml |
| BSB 5826 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2021 | Tinto Predilute | 3.0 ml |
| BSB 2022 | Tinto Predilute | 7.0 ml |
| BSB 2023 | Tinto Predilute | 15.0 ml |
| BSB 2024 | Concentrate | 0.1 ml |
| BSB 2025 | Concentrate | 0.5 ml |
| BSB 2026 | Concentrate | 1.0 ml |
| BSB 2027 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2028 | Tinto Predilute | 3.0 ml |
| BSB 2029 | Tinto Predilute | 7.0 ml |
| BSB 2030 | Tinto Predilute | 15.0 ml |
| BSB 2031 | Concentrate | 0.1 ml |
| BSB 2032 | Concentrate | 0.5 ml |
| BSB 2033 | Concentrate | 1.0 ml |
| BSB 2034 | Control Slides | 5 |

OLIG2, RMAb


IHC of OLIG2 on a FFPE Brain Tissue

Oligodendrocyte transcription factor (OLIG2) is a basic helix-loop-helix (bHLH) transcription factor encoded by the Olig2 gene. The protein is of 329 amino acids in length, 32kDa in size and contains 1 basic helix-loop-helix DNA-binding domain. It is one of the three members of the bHLH family. The other two members are OLIG1 and OLIG3. The expression of OLIG2 is mostly restricted in central nervous system, where it acts as both an anti-neurigenic and a neurigenic factor at different stages of development. OLIG2 is well known for determining motor neuron and oligodendrocyte differentiation, as well as its role in sustaining replication in early development.

No OLIG2 expression has been found in the non-glial tumors including neuroepithelial tumors, ependymomas, subependymomas, medulloblastomas, and non-neuroepithelial tumors, such as CNS lymphomas, meningiomas, schwannomas, atypical teratoid/rhabdoid tumor, and haemangioblastomas. Compared to the strong staining seen in glioma samples, a weak expression is observed in non-tumoral brain tissue (gliosis). OLIG2 is universally expressed in glioblastoma and other diffuse gliomas (astrocytomas, oligodendrogliomas and oligoastrocytomas), and is a useful positive diagnostic marker of these brain tumors.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP112
ISOTYPE: IgG
CONTROL: Tonsil, Colon, Brain, Astrocytoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2559 | Tinto Predilute | 3.0 ml |
| BSB 2560 | Tinto Predilute | 7.0 ml |
| BSB 2561 | Tinto Predilute | 15.0 ml |
| BSB 2562 | Concentrate | 0.1 ml |
| BSB 2563 | Concentrate | 0.5 ml |
| BSB 2564 | Concentrate | 1.0 ml |
| BSB 2565 | Control Slides | 5 |

Osteonectin/SPARC, MMab


IHC of Osteonectin / SPARC on a FFPE Bladder TCC Tissue

Osteonectin also known as secreted protein acidic and rich in cysteine (SPARC) or basement-membrane protein 40 (BM-40) is a protein that in humans is encoded by the SPARC gene. Osteonectin is a glycoprotein in the bone that binds calcium. Fibroblasts, including periodontal fibroblasts, synthesize Osteonectin. This protein is synthesized by macrophages at sites of wound repair and platelet degranulation, so it may play an important role in wound healing.

Osteonectin also increases the production and activity of matrix metalloproteinases, a function important to invading cancer cells within bone. Overexpression of Osteonectin is reported in many human cancers such as breast, prostate and colon. A correlation between Osteonectin overexpression and ampullary cancers and chronic pancreatitis has been reported. By immunohistochemistry, faint immunoreactivity was detected in the normal pancreas. In contrast, strong staining of the cancer cells was observed in addition to extensive Osteonectin immunoreactivity in surrounding fibroblasts and in the extracellular matrix. In metastatic tissues, strong immunoreactivity was observed in fibroblasts and in extracellular matrix surrounding metastatic cancer cells, whereas the signal was absent in most tumor cells.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-93
ISOTYPE: IgG1
CONTROL: Adrenal, Testis, Placenta, TCC, Testicular Cancer, Cervical Cancer
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3258 | Tinto Predilute | 3.0 ml |
| BSB 3259 | Tinto Predilute | 7.0 ml |
| BSB 3260 | Tinto Predilute | 15.0 ml |
| BSB 3261 | Concentrate | 0.1 ml |
| BSB 3262 | Concentrate | 0.5 ml |
| BSB 3263 | Concentrate | 1.0 ml |
| BSB 3264 | Control Slides | 5 |

OX-40/CD134, MMab

IHC of OX-40/CD134 on a FFPE Lymph Node Tissue


OX-40 also known as CD134 and Tumor Necrosis Factor Receptor Superfamily Member 4 (TNFRSF4), is a member of the TNFR-superfamily of receptors. OX40 and OX40L also regulate cytokine production from T cells, antigen-presenting cells, NK cells, and NKT cells, and modulate cytokine receptor signaling. In line with these important modulatory functions, OX40/OX40L interactions have been found to play a central role in the development of multiple inflammatory and autoimmune diseases. In addition, recent genome-wide association studies have identified single-nucleotide polymorphisms of the OX40L and OX40 genes that are related to cardiovascular diseases and SLE, providing direct evidence for the involvement of the OX40-OX40L interaction in human diseases.

OX40 is a potent costimulatory receptor that can potentiate T-cell receptor signaling on the surface of T lymphocytes, leading to their activation by a specifically recognized antigen. In particular, OX40 engagement by ligands present on dendritic cells dramatically increases the proliferation, effector function, and survival of T cells. Preclinical studies have shown that OX40 agonists increase antitumor immunity and improve tumor-free survival by increasing T and B cell responses to reporter antigen immunizations, led to preferential upregulation of OX40 on CD4(+) FoxP3(+) regulatory T cells in tumor-infiltrating lymphocytes, and increased the antitumor reactivity of T and B cells in patients with melanoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-90
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node, Thymus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3120 | Tinto Predilute | 3.0 ml |
| BSB 3121 | Tinto Predilute | 7.0 ml |
| BSB 3122 | Tinto Predilute | 15.0 ml |
| BSB 3123 | Concentrate | 0.1 ml |
| BSB 3124 | Concentrate | 0.5 ml |
| BSB 3125 | Concentrate | 1.0 ml |
| BSB 3126 | Control Slides | 5 |

p14 ARF/ CDKN2A, RMab




IHC of p14 ARF on an FFPE Anal Carcinoma Tissue

p14 ARF (also called ARF tumor suppressor, ARF, p14ARF) encoded by the p16 tumor suppressor gene is an alternate reading frame protein product of the CDKN2A locus (i.e. INK4a/ARF locus). p14ARF accumulates mainly in the nucleolus where it forms stable complexes with NPM or MDM2. Both p16INK4a and p14ARF are involved in cell cycle regulation. p14ARF inhibits MDM2, thus promoting p53, which promotes p21 activation, which then binds and inactivates certain cyclin-CDK complexes, which would otherwise promote transcription of genes that would carry the cell through the G1/S checkpoint of the cell cycle. Loss of p14ARF by a homozygous mutation in the CDKN2A (INK4A) gene will lead to elevated levels in MDM2 and, therefore, loss of p53 function and cell cycle control.

p14 ARF, has been reported to be associated with the clinicopathological features of different cancers. Very commonly, cancer is associated with a loss of function of INK4a, ARF, Rb, or p53. Without ARF, MDM2 can inappropriately inhibit p53, leading to increased cell survival. The INK4a/ARF locus is found to be deleted or silenced in many kinds of tumors. Homozygous deletions and other mutations of CDK2NA (ARF) have been found to be associated with Glioblastoma. p14ARF expression has been found to be significantly associated with the risk of lung cancer.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-p14
ISOTYPE: IgG
CONTROL: Cervical, Anal and Ovarian Carcinomas
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

p16, MMab




IHC of p16 on an FFPE Transitional Cell Carcinoma Tissue

p16 is a tumour suppressor gene. p16 is an important gene in regulating the cell cycle. p16INK4a regulates the cell cycle by binding and deactivating various cyclin-Cdk complexes. p16 is a G1/S-cell cycle regulator that is involved in the pathway that converges in the tumor suppressor protein Rb.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 16P04,JC2
ISOTYPE: IgG1
CONTROL: Testis, NSCLC, TCC
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

p16, RMab



IHC of p16 on a FFPE Mucoepidermoid Carcinoma Tissue


p16 is a tumor suppressor gene. p16 is an important gene in regulating the cell cycle. p16INK4a regulates the cell cycle by binding and deactivating various cyclin-CDK complexes. p16 is a G1/S-cell cycle regulator that is involved in the pathway that converges in the tumor suppressor protein Rb.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM267
ISOTYPE: IgG
CONTROL: Testis, NSCLC, Transitional Cell Carcinoma
LOCALIZATION: Nuclear, Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3595 | Tinto Predilute | 3.0 ml |
| BSB 3596 | Tinto Predilute | 7.0 ml |
| BSB 3567 | Tinto Predilute | 15.0 ml |
| BSB 3598 | Concentrate | 0.1 ml |
| BSB 3599 | Concentrate | 0.5 ml |
| BSB 3600 | Concentrate | 1.0 ml |
| BSB 3601 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5827 | Tinto Predilute | 3.0 ml |
| BSB 5828 | Tinto Predilute | 7.0 ml |
| BSB 5828 | Tinto Predilute | 15.0 ml |
| BSB 5830 | Concentrate | 0.1 ml |
| BSB 5831 | Concentrate | 0.5 ml |
| BSB 5832 | Concentrate | 1.0 ml |
| BSB 5833 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3768-3 | Tinto Predilute | 3.0 ml |
| BSB-3768-7 | Tinto Predilute | 7.0 ml |
| BSB-3768-15 | Tinto Predilute | 15.0 ml |
| BSB-3768-01 | Concentrate | 0.1 ml |
| BSB-3768-05 | Concentrate | 0.5 ml |
| BSB-3768-1 | Concentrate | 1.0 ml |
| BSB-3768-CS | Control Slides | 5 |

p16, RMab

IHC of p16 on an FFPE Lung Squamous Cell Carcinoma Tissue


p16 is a tumor suppressor gene. p16 is an important gene in regulating the cell cycle. p16INK4a regulates the cell cycle by binding and deactivating various cyclin-CDK complexes. p16 is a G1/S-cell cycle regulator that is involved in the pathway that converges in the tumor suppressor protein Rb.

p21, MMab

IHC of p21 on a FFPE Colon Carcinoma Tissue

p21 is a potent cyclin-dependent kinase inhibitor. The p21 protein binds to and inhibits the activity of cyclin-CDK2 or -CDK1 complexes, and thus functions as a regulator of cell cycle progression at G1. The expression of this gene is tightly controlled by the tumor suppressor protein p53, through which this protein mediates the p53-dependent cell cycle G1 phase arrest in response to a variety of stimuli. In addition to growth arrest, p21 can mediate cellular senescence. p21 can also interact with proliferating cell nuclear antigen (PCNA) and plays a regulatory role in S phase DNA replication and DNA damage repair.

Normal cells typically display a rather intense nuclear p21 expression. Loss of p21 expression has been associated with poor prognosis in several carcinomas including Gastric Carcinoma, Non-Small Cell Lung Carcinoma and Thyroid Carcinoma.

p27, MMab

IHC of p27 on a FFPE Breast Carcinoma Tissue

p27KIP1 is a cell cycle regulatory mitotic inhibitor of Cdk activity. p27KIP1 is a candidate-tumor suppressor gene, and has been proposed to function as a possible mediator of TGF beta induced G1 arrest. p27 is up-regulated in response to antimetabolic stimuli. The increased protein expression of p27 results in cellular arrest by binding to cyclin/Cdk complexes such as cyclin D1/Cdk4.

Low p27 expression has been associated with unfavorable prognosis in Renal-cell Carcinoma, Colon Carcinoma, Breast Carcinomas, Non-small-cell Lung Carcinoma, Hepatocellular Carcinoma, Multiple Myeloma, lymph node metastases in Papillary Carcinoma of the Thyroid, and a more aggressive phenotype of Carcinoma in the Cervix.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-p16
ISOTYPE: IgG
CONTROL: Testis, NSCLC, TCC
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: DCS-60.2
ISOTYPE: IgG2a
CONTROL: Tonsil, Colon, Fallopian Tube, Breast Cancer, Colon Carcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: SX53G8
ISOTYPE: IgG1/K
CONTROL: Testis, Breast, Adrenal, Prostate, Tonsil, Lung, Colon, Non-Small Cell Lung Carcinoma, Colon Adenocarcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3476 | Tinto Predilute | 3.0 ml |
| BSB 3477 | Tinto Predilute | 7.0 ml |
| BSB 3478 | Tinto Predilute | 15.0 ml |
| BSB 3479 | Concentrate | 0.1 ml |
| BSB 3480 | Concentrate | 0.5 ml |
| BSB 3481 | Concentrate | 1.0 ml |
| BSB 3684 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2063 | Tinto Predilute | 3.0 ml |
| BSB 2064 | Tinto Predilute | 7.0 ml |
| BSB 2065 | Tinto Predilute | 15.0 ml |
| BSB 2066 | Concentrate | 0.1 ml |
| BSB 2067 | Concentrate | 0.5 ml |
| BSB 2068 | Concentrate | 1.0 ml |
| BSB 2069 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5834 | Tinto Predilute | 3.0 ml |
| BSB 5835 | Tinto Predilute | 7.0 ml |
| BSB 5836 | Tinto Predilute | 15.0 ml |
| BSB 5837 | Concentrate | 0.1 ml |
| BSB 5838 | Concentrate | 0.5 ml |
| BSB 5839 | Concentrate | 1.0 ml |
| BSB 5840 | Control Slides | 5 |

p40, RMAb




IHC and IF of p40 on an FFPE Prostate Tissue (IHC) and an FFPE Tonsil Tissue (IF)

p53, MMAb



IHC of p53 on a FFPE Breast Carcinoma Tissue

p57, MMAb



IHC of p57 on a FFPE Placenta Tissue

p40 is an antibody that recognizes ΔNp63-a p63 isoform and it is highly specific for squamous/basal cells. It may be a valuable marker in detecting Squamous Cell Carcinoma where p63 is currently used. It recognizes the shortest variant of p53. p40 is superior in specificity to p63 because it does not label lung adenocarcinomas like p63 does, which eliminates the potential of misinterpreting a positive adenocarcinoma as a squamous cell carcinoma.

p53 (also known as tumor protein 53 [TP53]) is a transcription factor that regulates the cell cycle and, hence, functions as a tumor suppressor. p53 has been described as "the guardian of the genome", referring to its role in conserving stability by preventing genome mutation. p53 has many anti-cancer mechanisms. It can activate DNA repair proteins when DNA has sustained damage; it can also hold the cell cycle at the G1/S regulation point on DNA damage recognition. It can initiate apoptosis, programmed cell death, if DNA damage proves to be irreparable. p53 is central to many of the cell's anti-cancer mechanisms. It can induce growth arrest, apoptosis and cell senescence.

Mutations involving p53 have been found in a wide variety of malignant tumors, including Breast, Ovarian, Bladder, Colon, Lung, and Melanoma.

p57 or p57KIP2 is a tumor-suppressor human gene that belongs to the cip/kip gene family. It encodes a cell cycle inhibitor that binds to G1 cyclin-CDK complexes. Thus, p57 causes arrest of the cell-cycle in G1 phase. A mutation of this gene may lead to loss of control over the cell-cycle leading to uncontrolled cellular proliferation. The gene encoding human p57KIP2 is located on chromosome 11p15.5, a region implicated in sporadic cancers, Wilm's Tumor and Beckwith Wiedemann Syndrome (BWS is characterized by increased risk of tumor formation in childhood), making it a tumor suppressor candidate.

Anti-p57 has been used to aide in discriminating Complete Hydatidiform Mole (CHM) (no nuclear labeling of cytotrophoblasts) from Partial Hydatidiform Mole (PHM) and hydropic abortion. In normal placenta, many cytotrophoblast nuclei and stromal cells are labelled with this antibody. Similar findings apply to PHM and hydropic abortion tissues. Intervillous Trophoblastic Islands (IVTIs) demonstrate nuclear labeling in all three entities and serve as an internal control.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: ZR8
ISOTYPE: IgG
CONTROL: Normal Prostate, Breast, Skin
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: DO7
ISOTYPE: IgG2b/K
CONTROL: Lung, Breast, Ovarian Carcinoma, Prostatic Carcinoma, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Cat, Horse, Sheep, Bovine, Monkey


ANTIBODY TYPE: Mouse Monoclonal
CLONE: Kp10
ISOTYPE: IgG1b/K
CONTROL: Placenta, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2070 | Tinto Predilute | 3.0 ml |
| BSB 2071 | Tinto Predilute | 7.0 ml |
| BSB 2072 | Tinto Predilute | 15.0 ml |
| BSB 2073 | Concentrate | 0.1 ml |
| BSB 2074 | Concentrate | 0.5 ml |
| BSB 2075 | Concentrate | 1.0 ml |
| BSB 2076 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5841 | Tinto Predilute | 3.0 ml |
| BSB 5842 | Tinto Predilute | 7.0 ml |
| BSB 5843 | Tinto Predilute | 15.0 ml |
| BSB 5844 | Concentrate | 0.1 ml |
| BSB 5845 | Concentrate | 0.5 ml |
| BSB 5846 | Concentrate | 1.0 ml |
| BSB 5847 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6191 | Tinto Predilute | 3.0 ml |
| BSB 6192 | Tinto Predilute | 7.0 ml |
| BSB 6193 | Tinto Predilute | 15.0 ml |
| BSB 6194 | Concentrate | 0.1 ml |
| BSB 6195 | Concentrate | 0.5 ml |
| BSB 6196 | Concentrate | 1.0 ml |
| BSB 6197 | Control Slides | 5 |

p63, MMab



IHC and IF of p63 on an FFPE Basal Cell Carcinoma Tissue (IHC) an FFPE Tonsil Tissue (IF)


In addition to p53, mammalian cells contain two homologous genes, p63 and p73. These genes give rise to the expression of proteins that are highly similar to p53 in structure and function. In particular, p63 and p73 proteins can induce p53-responsive genes and elicit programmed cell death. p73 and p63 are important during development and differentiation. In particular, p63 appears to be primarily implicated in epithelial development.

Anti-p63 to human p63 protein labels an epitope common to all six p63 isotypes (TAp63 α , TAp63 β , TAp63 γ , Δ Np63 α , Δ Np63 β , Δ Np63 γ). p63 labels the nuclei of myoepithelial cells in the prostate gland as well as breast tissue, making it useful in differentiating benign vs. malignant prostate lesions and breast lesions.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 4A4
ISOTYPE: IgG2a/K
CONTROL: Prostate, Breast, Skin, Salivary Gland
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3602 | Tinto Predilute | 3.0 ml |
| BSB 3603 | Tinto Predilute | 7.0 ml |
| BSB 3604 | Tinto Predilute | 15.0 ml |
| BSB 3605 | Concentrate | 0.1 ml |
| BSB 3606 | Concentrate | 0.5 ml |
| BSB 3607 | Concentrate | 1.0 ml |
| BSB 3608 | Control Slides | 5 |

p63, RMAb



IHC and IF of p63 on a FFPE Prostate Tissue (IHC) and a FFPE Skin Tissue (IF)


In addition to p53, mammalian cells contain two homologous genes, p63 and p73. These genes give rise to the expression of proteins that are highly similar to p53 in structure and function. In particular, p63 and p73 proteins can induce p53-responsive genes and elicit programmed cell death. p73 and p63 are more important during development and differentiation. In particular, p63 appears to be primarily implicated in epithelial development.

Anti-p63 to human p63 protein labels an epitope common to all six p63 isotypes (TAp63 α , TAp63 β , TAp63 γ , Δ Np63 α , Δ Np63 β , Δ Np63 γ). p63 labels the nuclei of myoepithelial cells in the prostate gland as well as breast tissue, making it useful in differentiating benign vs. malignant prostate lesions and breast lesions.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP174
ISOTYPE: IgG
CONTROL: Prostate, Breast, Skin, Salivary Gland
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat, Rabbit

| CAT. # | PRESENTATION | VOL/QTY |
|------------|-----------------|---------|
| BSB 5848 | Tinto Predilute | 3.0 ml |
| BSB 5849 | Tinto Predilute | 7.0 ml |
| BSB 5850 | Tinto Predilute | 15.0 ml |
| BSB 5851 | Concentrate | 0.1 ml |
| BSB 5852 | Concentrate | 0.5 ml |
| BSB 5853 | Concentrate | 1.0 ml |
| BSB 5854-1 | Control Slides | 5 |

p120 Catenin, RMAb



IHC of p120 Catenin on a FFPE Breast Carcinoma Tissue


p120 Catenin is a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. The association of catenins to cadherins produces a complex which is linked to the actin filament network, and which seems to be important for cadherins cell-adhesion properties. Cytoplasmic accumulation of p120 Catenin has been observed in lung cancer, pancreatic cancer, gastric cancer and colon cancers and is associated with poor prognosis in colon cancer patients.

In breast lobular neoplasia, anti-p120 Catenin shows a diffuse cytoplasmic immunostaining pattern, while breast ductal neoplasia retains the membrane immunostaining pattern. p120 Catenin can be useful in differentiating between lobular carcinoma and ductal carcinoma of the breast, and in identifying early lesions of lobular neoplasia.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP66
ISOTYPE: IgG
CONTROL: Breast, Testis, Kidney, Prostate, Pancreas, Tonsil, Salivary Gland, Skin, Cervix, Colon, Malignant, Melanoma, Transitional Cell Carcinoma, Breast Lobular Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse, Rat


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2077 | Tinto Predilute | 3.0 ml |
| BSB 2078 | Tinto Predilute | 7.0 ml |
| BSB 2079 | Tinto Predilute | 15.0 ml |
| BSB 2080 | Concentrate | 0.1 ml |
| BSB 2081 | Concentrate | 0.5 ml |
| BSB 2082 | Concentrate | 1.0 ml |
| BSB 2083 | Control Slides | 5 |

pan TRK, RMAb




IHC panTRK on a Papillary Thyroid Carcinoma Tissue

pan-TRK, RMAb



IHC of pan-TRK on a FFPE Diffuse Gastric Carcinoma Tissue

Parafibromin, MMAb



IHC of Parafibromin on a FFPE Parathyroid Carcinoma Tissue

Pan-TRK IHC has shown to be positive in most cases with NTRK fusion transcripts confirmed by Archer. One study established the Pan-TRK IHC sensitivity and specificity for transcribed NTRK fusions to be 95.2% and 100%, respectively. All positive IHC cases had cytoplasmic staining while the following fusion partner-specific patterns were discovered: all LMNA-NTRK1 fusions displayed nuclear membrane accentuation, all TPM3/4 fusions displayed cellular membrane accentuation, and half of ETV6-NTRK3 fusions displayed nuclear staining. NTRK gene fusions have been identified in both pediatric and adult primary central nervous system tumors, including Glioblastoma Multiforme, Pediatric Gliomas and Astrocytomas. Various translocations involving NTRK1 or NTRK3 have been reported in Spitzoid melanocytic neoplasms, as well as in compound Spitz Nevi. TRK fusions have also been reported in Intrahepatic Cholangiocarcinomas, Breast Cancer, quadruple wild-type (ETV6-NTRK3) Gastrointestinal Stromal Tumors, Gallbladder Adenocarcinomas, Pancreatic Carcinomas, Sinus-Nasal Low-Grade Non-Intestinal-type Adenocarcinomas and Neuroendocrine Tumors of the small bowel. In addition to being present in solid tumors, NTRK gene fusions have been also detected in Acute Lymphoblastic Leukemia and Acute Myeloid Leukemia.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-TRK
ISOTYPE: IgG
CONTROL: Brain, Lung, Papillary Thyroid Carcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

Neurotrophic tyrosine kinase (NTRK) proto-oncogene family codes for proteins Trk A, Trk B, and Trk C, which participate in pathways of neuron cell growth, differentiation, signaling, and survival. The transmembrane neurotrophic receptors are activated by neurotrophins (Nerve Growth Factor, Brain-Derived Growth Factor, and neurotrophins 3/4/5) and in turn the TRKs activate MAPK, AKT, and Phospholipase C pathways. Non-fusion alterations in NTRK have been found in 14% of tested cancers. panTRK staining can be cytoplasmic, nuclear (as in NTRK-ETV6 and NTRK-LMNA fusions), or membranous (in NTRK-TMP3/4) depending on the fusion pair.

NTRK fusions result in NTRK3-ETV6 pairs in 90% of cases, most commonly found in Carcinomas and Sarcomas of the Mammary and Salivary Secretory Glands. NTRK gene fusions are also found in Brain primary tumors and metastases, Lung, Breast, Papillary Thyroid Carcinoma, Colorectal and Pancreatic cancer. NTRK mesenchymal tumors have multiple morphological features and coexpression of S100, CD34, and panTRK.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM423
ISOTYPE: IgG
CONTROL: Brain, Lung Neuroendocrine Cancer
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

Cell division cycle 73, Paf1/RNA polymerase II complex component, homolog (S. cerevisiae), also known as CDC73 and parafibromin, is a protein which in humans is encoded by the CDC73 gene. Mutations in the CDC73 gene are associated with hyperparathyroidism-jaw tumor syndrome (HPT-JT) and parathyroid carcinomas.

The proposed role of HRPT2 as a tumor suppressor was supported by mutation screening in parathyroid adenomas with cystic features, which identified three somatic inactivating mutations, all located in exon 1. None of these mutations were detected in normal controls, and all were predicted to cause deficient or impaired protein function. Sporadic parathyroid carcinomas frequently have HRPT2 mutations that are likely to be of pathogenetic importance. Certain patients with apparently sporadic parathyroid carcinoma carry germ-line mutations in HRPT2 and may have the HPT-JT syndrome or a phenotypic variant.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-50
ISOTYPE: IgG1
CONTROL: Parathyroid, Colon, Testis, Adrenal, Breast, Cervix, Kidney, Pituitary, Brain, Pancreas, Salivary Gland, Lymphoblastic Lymphoma, Transitional Cell Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-2376-3 | Tinto Predilute | 3.0 ml |
| BSB-2376-7 | Tinto Predilute | 7.0 ml |
| BSB-2376-15 | Tinto Predilute | 15.0 ml |
| BSB-2376-01 | Concentrate | 0.1 ml |
| BSB-2376-05 | Concentrate | 0.5 ml |
| BSB-2376-1 | Concentrate | 1.0 ml |
| BSB-2376-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3742-3 | Tinto Predilute | 3.0 ml |
| BSB-3742-7 | Tinto Predilute | 7.0 ml |
| BSB-3742-15 | Tinto Predilute | 15.0 ml |
| BSB-3742-01 | Concentrate | 0.1 ml |
| BSB-3742-05 | Concentrate | 0.5 ml |
| BSB-3742-1 | Concentrate | 1.0 ml |
| BSB-3742-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2747 | Tinto Predilute | 3.0 ml |
| BSB 2748 | Tinto Predilute | 7.0 ml |
| BSB 2749 | Tinto Predilute | 15.0 ml |
| BSB 2750 | Concentrate | 0.1 ml |
| BSB 2751 | Concentrate | 0.5 ml |
| BSB 2752 | Concentrate | 1.0 ml |
| BSB 2753 | Control Slides | 5 |

Parathyroid Hormone/PTH, MAb



IHC of Parathyroid Hormone on a FFPE Parathyroid Tissue


Parathyroid hormone (PTH), parathormone or parathyrin, is secreted by the chief cells of the parathyroid glands as a polypeptide containing 84 amino acids. It acts to increase the concentration of calcium (Ca²⁺) in the blood, whereas Calcitonin (a hormone produced by the parafollicular cells (C cells) of the thyroid gland) acts to decrease calcium concentration. PTH acts to increase the concentration of calcium in the blood by acting upon the parathyroid hormone 1 receptor (high levels in bone and kidney) and the parathyroid hormone 2 receptor (high levels in the central nervous system, pancreas, testis, and placenta). PTH half-life is approximately 4 minutes.

Anti-PTH antibody is also useful to distinguish parathyroid hyperplasia/neoplasms from thyroid and metastatic neoplasms. If the patient's hypercalcemic status is not known, PTH immunohistochemistry is helpful, especially if clear cell parathyroid carcinomas are nonsecretory and there is no abnormality in mineral metabolism. The other instance in which anti-PTH antibodies are useful is in the consideration of parathyroid carcinomas located primarily in the anterior mediastinum. In this situation distinction from primary thymic metastatic carcinomas, non-Hodgkin's lymphoma and germ cell tumors is necessary.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-24
ISOTYPE: IgG1
CONTROL: Parathyroid
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2084 | Tinto Predilute | 3.0 ml |
| BSB 2085 | Tinto Predilute | 7.0 ml |
| BSB 2086 | Tinto Predilute | 15.0 ml |
| BSB 2087 | Concentrate | 0.1 ml |
| BSB 2088 | Concentrate | 0.5 ml |
| BSB 2089 | Concentrate | 1.0 ml |
| BSB 2090 | Control Slides | 5 |

Parvalbumin, RMAb



IHC of Parvalbumin on a FFPE Kidney Tissue

Parvalbumin is a calcium-binding albumin protein with low molecular weight (typically 9-11 kDa), structurally related to calmodulin and troponin C and a stable protein involved in calcium signaling. Parvalbumin plays a role in many physiological processes, namely cell-cycle regulation, second messenger production, muscle contraction, organization of microtubules and phototransduction. Parvalbumin is localized in fast-contracting muscles, where its levels are highest, as well as in the brain and some endocrine tissues. In normal kidney, parvalbumin has been shown to be limited to the distal tubular and collecting duct cells (the intercalated cells).


Studies have demonstrated that parvalbumin may be a suitable IHC marker for distinguishing primary and metastatic chromophobe carcinoma from conventional (clear cell) and papillary renal cell carcinoma. Parvalbumin is strongly expressed in almost all primary, as well as metastatic, chromophobe renal cell carcinoma (RCC) (100%) and oncocytoma (69%), but is essentially negative in other types of RCCs, such as clear cell RCC and papillary RCC.

Additionally, alterations in the function of parvalbumin-expressing neurons have been implicated in various areas of clinical interest such as Alzheimer's disease and age-related cognitive defects.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP300
ISOTYPE: IgG
CONTROL: Brain, Kidney, Tonsil, Lymph Node, Chromophobe RCC
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3406 | Tinto Predilute | 3.0 ml |
| BSB 3407 | Tinto Predilute | 7.0 ml |
| BSB 3408 | Tinto Predilute | 15.0 ml |
| BSB 3409 | Concentrate | 0.1 ml |
| BSB 3410 | Concentrate | 0.5 ml |
| BSB 3411 | Concentrate | 1.0 ml |
| BSB 3412 | Control Slides | 5 |

Parvovirus, MAb




IHC of Parvovirus B19 on a FFPE Placenta Tissue

Parvovirus B19 belongs to the Parvoviridae family of small DNA viruses. It is classified as Erythrovirus because of its capability to invade red blood cell precursors in the bone marrow. Anti-Parvovirus antibody targets the capsid proteins VP1 and VP2 on Human Parvovirus.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: R92F6
ISOTYPE: IgG1
CONTROL: Parvovirus Infected Tissue
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5854 | Tinto Predilute | 3.0 ml |
| BSB 5855 | Tinto Predilute | 7.0 ml |
| BSB 5856 | Tinto Predilute | 15.0 ml |
| BSB 5857 | Concentrate | 0.1 ml |
| BSB 5858 | Concentrate | 0.5 ml |
| BSB 5859 | Concentrate | 1.0 ml |
| BSB 5860 | Control Slides | 5 |

PAX-2, RMAb




IHC of PAX-2 on a FFPE Kidney Tissue

PAX-5, RMAb



IHC of PAX-5 on a FFPE Tonsil Tissue

PAX-6, RMAb



IHC of PAX-6 on a FFPE Pancreas Tissue

PAX-2 is a homeogene strongly expressed during kidney development. PAX-2 gene is expressed in the metanephric mesenchyma after ureter bud induction and is a key factor for the mesenchyma-epithelium conversion. Animals transgenic for PAX-2 have severe renal abnormalities and cysts but no solid tumoral features.

Anti-PAX-2 can be used to distinguish Ovarian Serous Papillary Carcinoma (PAX-2 positive) from Primary Breast Carcinoma (PAX-2 negative). It can also be used to distinguish Clear Cell Renal Carcinoma (positive) from Hepatocellular Carcinoma (negative).

The PAX proteins are important regulators in early development, and alterations in the expression of their genes are thought to contribute to neoplastic transformation. The PAX-5 gene encodes the B-cell lineage-specific activator protein (BSAP) that is expressed at early, but not late, stages of B-cell differentiation. Its expression has also been detected in developing CNS and testis; therefore, PAX-5 gene product may not only play an important role in B-cell differentiation, but also in neural development and spermatogenesis.

PAX-5 expression is not only continuously required for B-cell lineage commitment during early B-cell development but also for B-cell lineage maintenance. PAX-5 is found in most cases of mature and precursor B-cell Non-Hodgkin's Lymphomas/Leukemias. PAX-5 is not detected in Multiple Myeloma and solitary Plasmacytoma, making it useful for such differentiation. Diffuse Large B-cell Lymphomas do express PAX-5, except for those with terminal B-cell differentiation. T-cell neoplasms do not stain with anti-PAX-5; however, there is a strong association with CD20 expression.

Paired box protein PAX-6 also known as aniridia type II protein (AN2) is a protein that in humans is encoded by the PAX6 gene. PAX6 is a transcription factor present during embryonic development of sensory organs (including eye, nasal and olfactory tissues), central nervous and endocrine system. As a transcription factor, PAX6 activates and/or deactivates gene expression patterns to ensure for proper development tissues. Mutations of the PAX6 gene are known to cause various disorders of the eyes. Two common disorders associated with a mutation are: aniridia, the absence of the iris, and Peter's anomaly, thinning and clouding of the cornea.

PAX6 labels neuroendocrine cells and derived tumor cells and is helpful in identification of neuroendocrine tumors. A recent study showed that PAX6 and PAX8 were positive in the majority of neuroendocrine tumors originated from pancreas, duodenum, and colon. Additionally, Neuroendocrine tumors of the lung (NELC), which account for 25% of all lung cancer cases, and transcription factors may drive dedifferentiation of these tumors. SOX4 (p = 0.0002), SOX11 (p < 0.0001) and PAX6 (p = 0.0002) have been found to be significant for tumor type and elevated PAX6 and SOX11 expression correlates with poor outcome in large cell neuroendocrine carcinomas and small cell lung cancer (p < 0.0001 and p = 0.0232, respectively) based on survival data of 34 patients (57%). Therefore, aggressiveness of NELC correlated with increasing expression of transcription factors.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP235

ISOTYPE: IgG

CONTROL: Kidney, Fallopian Tube, Clear Cell Renal Carcinoma, Ovarian Serous Papillary Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-PAX5

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Spleen, Thymus, Colon, Liver, Lymphoblastic Lymphoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP341

ISOTYPE: IgG

CONTROL: Pancreas, Pituitary, Neuroendocrine Tumors


LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human, predicted Mouse, predicted Rat


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2566 | Tinto Predilute | 3.0 ml |
| BSB 2567 | Tinto Predilute | 7.0 ml |
| BSB 2568 | Tinto Predilute | 15.0 ml |
| BSB 2569 | Concentrate | 0.1 ml |
| BSB 2570 | Concentrate | 0.5 ml |
| BSB 2571 | Concentrate | 1.0 ml |
| BSB 2572 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5861 | Tinto Predilute | 3.0 ml |
| BSB 5862 | Tinto Predilute | 7.0 ml |
| BSB 5863 | Tinto Predilute | 15.0 ml |
| BSB 5864 | Concentrate | 0.1 ml |
| BSB 5865 | Concentrate | 0.5 ml |
| BSB 5866 | Concentrate | 1.0 ml |
| BSB 5867 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3127 | Tinto Predilute | 3.0 ml |
| BSB 3128 | Tinto Predilute | 7.0 ml |
| BSB 3129 | Tinto Predilute | 15.0 ml |
| BSB 3130 | Concentrate | 0.1 ml |
| BSB 3131 | Concentrate | 0.5 ml |
| BSB 3132 | Concentrate | 1.0 ml |
| BSB 3133 | Control Slides | 5 |

PAX-7, MAb

IHC of PAX-7 on a FFPE Renal Cell Carcinoma Tissue

PAX-8, RMab

IHC of PAX-8 on a FFPE Ovarian Carcinoma Tissue

PD-1/CD279, MAb

IHC of PD-1 on a FFPE Tonsil Tissue

PAX-7 or Paired Box Gene 7 is a transcription factor coded by a gene on locus 1p36, which can fuse with Forkhead Domain Region (FKHR). PAX-7 protein is involved in developmental pathways in neural tube and mesencephalon formation, muscle cell development, and oxidative stress sensitivity. PAX-7 regulates neural cell adhesion molecules and is expressed in proliferating myoblasts, but it is down-regulated after playing its role in cell differentiation toward muscle-derived specification.

PAX-7 is restricted to muscle satellite and myogenic precursor cells in adults, although it is upregulated in embryonal (ERMS) and alveolar rhabdomyosarcoma (ARMS) and could suggest a de-differentiated cell type when overexpressed in these tumors. Recent studies suggested that PAX-7 is a novel marker, because it was expressed consistently in Ewing sarcoma, in addition to rhabdomyosarcoma and synovial sarcoma. PAX7 was found expressed in 90% of Ewing sarcomas (90%), mainly in a diffuse and strong manner. Although NKX2-2 showed similar sensitivity, PAX7 showed more extensive and strong reactivity. ARMS differs from ERMS by virtue of its occurrence in older patients, distinctive pseudoalveolar pattern, usual absence of strap cells, and strong myogenin rather than MyoD1 expression. Identification of a PAX3 or PAX7/FKHR fusion gene may be necessary for the confident distinction of ARMS from the most primitive forms of ERMS.

PAX-8 is expressed in the thyroid (and associated carcinomas), non-ciliated mucosal cells of the fallopian tubes and simple ovarian inclusion cysts, but not normal ovarian surface epithelial cells. PAX-8 is expressed in a high percentage of ovarian serous, endometrioid, and clear cell carcinomas, but only rarely in primary ovarian mucinous adenocarcinomas. Studies have also found PAX-8 expression in renal tubules as well as renal carcinoma, nephroblastoma and seminoma. Normal lung and lung carcinomas do not express PAX-8. Similarly, the absence of expression of PAX-8 in breast and other non-GYN carcinomas other than those primary to the thyroid indicates that PAX-8 is an important new marker of ovarian cancer and a useful marker for the differential diagnoses in lung and neck tumors, or tumors at distant sites where primary lung carcinoma or thyroid carcinoma are possibilities.

PAX-8, combined with organ system-specific markers such as uroplakin, mammaglobin, and TTF-1 can be a very useful panel to determine the primary site of invasive micropapillary carcinomas of ovary from bladder, lung, and breast.

Programmed Death 1, (PD-1 or CD279), is a Type I membrane protein comprised of 268 amino acids. PD-1 is a member of the extended CD28/CTLA-4 family of T-cell regulators. PD-1 is expressed on the surface of activated T-cells, B-cells, and macrophages. In comparison to CTLA-4, PD-1 more broadly negatively regulates immune responses.

New data suggests that expression of PD-L1 on tumor cells inhibits anti-tumor activity through engagement of PD-1 on effector T-cells. Expression of PD-L1 on tumors is correlated with reduced survival in esophageal, pancreatic and other types of cancers, highlighting the relevance of exploring the PD-1 pathway as a target for immunotherapy. Studies have found that PD-1 is expressed on most T-cells and a small subset of B-cells in the light zone of germinal centers, but not elsewhere in the tonsil. On that basis, it was postulated that PD-1 may play a role in the process of clonal selection of centrocytes, which occurs in this subanatomic site in germinal centers. PD-1 is a new marker of Angioimmunoblastic Lymphoma and suggests a unique cell of origin for this neoplasm. Unlike CD10 and bcl-6, PD-1 is expressed by few B-cells, so it may be a more specific and useful diagnostic marker in Angioimmunoblastic Lymphoma. It also seems to stain a greater percentage of CD3-positive neoplastic cells in Angioimmunoblastic Lymphoma than either CD10 or bcl-6.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-145

ISOTYPE: IgG2a

CONTROL: Placenta, Brain, Testis, Prostate, Papillary Thyroid Carcinoma, Transitional Cell Carcinoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: ZR-1

ISOTYPE: IgG

CONTROL: Ovary, Thyroid

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: NAT-105

ISOTYPE: IgG1


CONTROL: Tonsil, Lymph Node, Thymus, Spleen

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB-3743-3 | Tinto Predilute | 3.0 ml | BSB 2098 | Tinto Predilute | 3.0 ml | BSB 6212 | Tinto Predilute | 3.0 ml |
| BSB-3743-7 | Tinto Predilute | 7.0 ml | BSB 2099 | Tinto Predilute | 7.0 ml | BSB 6213 | Tinto Predilute | 7.0 ml |
| BSB-3743-15 | Tinto Predilute | 15.0 ml | BSB 2100 | Tinto Predilute | 15.0 ml | BSB 6214 | Tinto Predilute | 15.0 ml |
| BSB-3743-01 | Concentrate | 0.1 ml | BSB 2101 | Concentrate | 0.1 ml | BSB 6215 | Concentrate | 0.1 ml |
| BSB-3743-05 | Concentrate | 0.5 ml | BSB 2102 | Concentrate | 0.5 ml | BSB 6216 | Concentrate | 0.5 ml |
| BSB-3743-1 | Concentrate | 1.0 ml | BSB 2103 | Concentrate | 1.0 ml | BSB 6217 | Concentrate | 1.0 ml |
| BSB-3743-CS | Control Slides | 5 | BSB 2104 | Control Slides | 5 | BSB 6218 | Control Slides | 5 |


PD-1/CD279, RMAb



IVD

IHC of PD1 on a FFPE Tonsil Tissue


PDGFR-B, RMAb



IVD

IHC of PDGFR-B on a FFPE Breast Tissue

PD-L1/CD274, RMAb



ASR

IHC PD-L1/CD274 on a FFPE CSL/SLL Lymphoma Tissue

Programmed Death 1 (PD-1 or CD279) is a Type I membrane protein comprised of 268 amino acids. PD-1 is a member of the extended CD28/CTLA-4 family of T-cell regulators. PD-1 is expressed on the surface of activated T-cells, B-cells, and macrophages. In comparison to CTLA-4, PD-1 more broadly negatively regulates immune responses.

New data suggests that expression of PD-L1 on tumor cells inhibits anti-tumor activity through engagement of PD-1 on effector T-cells. Expression of PD-L1 on tumors is correlated with reduced survival in esophageal, pancreatic and other types of cancers, highlighting the relevance of exploring the PD-1 pathway as a target for immunotherapy. Studies have found that PD-1 is expressed on most T-cells and a small subset of B-cells in the light zone of germinal centers, but not elsewhere in the tonsil. On that basis, it was postulated that PD-1 may play a role in the process of clonal selection of centrocytes, which occurs in this subanatomic site in germinal centers. PD-1 is a new marker of Angioimmunoblastic Lymphoma and suggests a unique cell of origin for this neoplasm. Unlike CD10 and bcl-6, PD-1 is expressed by few B-cells, so it may be a more specific and useful diagnostic marker in Angioimmunoblastic Lymphoma. It also seems to stain a greater percentage of CD3-positive neoplastic cells in Angioimmunoblastic Lymphoma than either CD10 or bcl-6.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP239
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Thymus, Spleen
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

Beta-type Platelet-Derived Growth Factor Receptor is a protein that in humans is encoded by the PDGFRB gene. This gene encodes a cell surface tyrosine kinase receptor for members of the Platelet-Derived Growth Factor family (PDGF), which is a mitogen for mesenchyme- and glia-derived cells. PDGF consists of two chains, A and B, which dimerize to form functionally distinct isoforms, PGDF-AA, PDGF-AB and PDGF-BB.

Translocation of the PDGFR gene with the Tel gene is linked to Chronic Myelomonocytic Leukemia (CMML), a myelodysplastic syndrome, and demonstrates the oncogenic potential of the PDGF receptors. A translocation between chromosomes 5 and 12, that fuses this gene to that of the translocation, ETV6, leukemia gene, results in chronic myeloproliferative disorder with eosinophilia.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-PDGFRB
ISOTYPE: IgG
CONTROL: Placenta, Breast, Kidney, Cervix, Skin, Fallopian Tube
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse, Rat

Programmed death-ligand 1 (PD-L1) also known as CD274 or B7 homolog 1 (B7-H1) is a protein that in humans is encoded by the CD274 gene. Programmed death-ligand 1 (PD-L1) is a 40 kDa type 1 transmembrane protein that has been speculated to play a major role in suppressing the immune system during particular events such as pregnancy, tissue allografts, autoimmune disease and other diseases.


The upregulation of PD-L1 may allow cancers to evade the host immune system. An analysis of tumor specimens from patients with renal cell carcinoma found that high tumor expression of PD-L1 was associated with increased tumor aggressiveness and a 4.5-fold increased risk of death. Ovarian cancer patients with higher expression of PD-L1 had a significantly poorer prognosis than those with lower expression. PD-L1 expression correlated inversely with intraepithelial CD8+ T-lymphocyte count, suggesting that PD-L1 on tumor cells may suppress antitumor CD8+ T cells. The PD-1/PD-L1 interaction is implicated in autoimmunity from several lines of evidence. In humans, PD-L1 was found to have altered expression in pediatric patients with Systemic lupus erythematosus. Studying isolated PBMC from healthy children, immature myeloid dendritic cells and monocytes expressed little PD-L1 at initial isolation, but spontaneously up-regulated PD-L1 by 24 hours. In contrast, both mDC and monocytes from patients with active SLE failed to upregulate PD-L1 over a 5 day time course, expressing this protein only during disease remissions.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: 28-8
ISOTYPE: IgG
CONTROL: Tonsil, Placenta, Lymphoblastic Lymphoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3148 | Tinto Predilute | 3.0 ml |
| BSB 3149 | Tinto Predilute | 7.0 ml |
| BSB 3150 | Tinto Predilute | 15.0 ml |
| BSB 3151 | Concentrate | 0.1 ml |
| BSB 3152 | Concentrate | 0.5 ml |
| BSB 3153 | Concentrate | 1.0 ml |
| BSB 3154 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2349 | Tinto Predilute | 3.0 ml |
| BSB 2350 | Tinto Predilute | 7.0 ml |
| BSB 2351 | Tinto Predilute | 15.0 ml |
| BSB 2352 | Concentrate | 0.1 ml |
| BSB 2353 | Concentrate | 0.5 ml |
| BSB 2354 | Concentrate | 1.0 ml |
| BSB 2355 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-2371-3 | Tinto Predilute | 3.0 ml |
| BSB-2371-7 | Tinto Predilute | 7.0 ml |
| BSB-2371-15 | Tinto Predilute | 15.0 ml |
| BSB-2371-01 | Concentrate | 0.1 ml |
| BSB-2371-05 | Concentrate | 0.5 ml |
| BSB-2371-1 | Concentrate | 1.0 ml |
| BSB-2371-CS | Control Slides | 5 |

PD-L1/CD274, RMAb

IHC PD-L1/CD274 on a Hodgkin's Lymphoma Tissue

PDX1, RMAb

IHC of PDX1 on a FFPE Colon Tissue

PELP1, RMAb

IHC of PELP1 on a FFPE Colon Adenocarcinoma Tissue

Programmed death-ligand 1 (PD-L1) also known as CD274 or B7 homolog 1 (B7-H1) is a protein that in humans is encoded by the CD274 gene. Programmed death-ligand 1 (PD-L1) is a 40 kDa type 1 transmembrane protein that has been speculated to play a major role in suppressing the immune system during particular events such as pregnancy, tissue allografts, autoimmune disease and other diseases.

The upregulation of PD-L1 may allow cancers to evade the host immune system. An analysis of tumor specimens from patients with renal cell carcinoma found that high tumor expression of PD-L1 was associated with increased tumor aggressiveness and a 4.5-fold increased risk of death. Ovarian cancer patients with higher expression of PD-L1 had a significantly poorer prognosis than those with lower expression. PD-L1 expression correlated inversely with intraepithelial CD8+ T-lymphocyte count, suggesting that PD-L1 on tumor cells may suppress antitumor CD8+ T cells. The PD-1/PD-L1 interaction is implicated in autoimmunity from several lines of evidence. In humans, PD-L1 was found to have altered expression in pediatric patients with Systemic lupus erythematosus. Studying isolated PBMC from healthy children, immature myeloid dendritic cells and monocytes expressed little PD-L1 at initial isolation, but spontaneously up-regulated PD-L1 by 24 hours. In contrast, both mDC and monocytes from patients with active SLE failed to upregulate PD-L1 over a 5 day time course, expressing this protein only during disease remissions.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-PDL1
ISOTYPE: IgG
CONTROL: Tonsil, Placenta, Lymphoblastic Lymphoma, Hodgkin's Lymphoma
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

PDX1 (Pancreatic and duodenal homeobox 1), also known as insulin promoter factor 1, is a transcription factor necessary for pancreatic development, including β -cell maturation, and duodenal differentiation. PDX1 appears to also play a role in the fating of endocrine cells, encoding for insulin and somatostatin, two pancreatic endocrine products, while repressing glucagon. PDX1 is required for β -cell survival. Cells with reduced PDX1 expression have an increased rate of apoptotic programmed cell death. Mutations in the PDX1 gene may be involved in several pancreatic pathologies, including diabetes mellitus and pancreatic cancer.

Among normal pancreatic tissues, PDX1 nuclear protein is expressed in islet cells, cells of the centroacinar cell compartment, ductal epithelium and is selectively expressed in adult Brunner's glands of the duodenum and pyloric endocrine cells of the stomach. PDX1 expression has been identified in Pancreatic Ductal Adenocarcinomas and endocrine neoplasms. No expression of PDX1 is seen in non-neoplastic acinar cells. Among pancreatic neoplasms, PDX1 consistently labeled >50% of the tumor cells. PDX1 expression is variable in invasive ductal adenocarcinoma and precursor lesions of ductal adenocarcinomas. Solid pseudopapillary neoplasms do not express PDX1. Besides increased expression of PDX1 in Pancreatic cancer, it also has also been reported in tumors of the colon and prostate, indicating that PDX1 may serve as a biomarker in patients with these malignancies.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP139
ISOTYPE: IgG
CONTROL: Pancreas, Colon, Liver, Pancreatic Cancer
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Proline-, glutamic acid- and leucine-rich protein 1 (PELP1) gene expresses PELP1, also known as HMX3 protein, a transcription factor that regulates the expression of many signaling proteins. PELP1 is mainly known as a coactivator for Estrogen Receptor activity, where it participates in complexes that facilitate transcription according to ER-mediated signals. PELP1 interacts with a variety of signaling proto-oncogenes such as Src, HER2, and EGFR, which are known to contribute to tumor formation, and with hormonal pathways that can contribute to hormonal therapy resistance. The wide range of signaling functions makes PELP1 susceptible to expression deregulation in hormone-centered cancers, particularly in Breast, Endometrial, Ovarian, and Prostate Tumors, as well as Ductal and Colorectal Carcinomas.

A study has found that PELP1 expression rate was the highest in breast cancers (70.5%) among different cancers. Compared to GATA3, Mammaglobin and GCDFFP-15, PELP1 was less sensitive than GATA3 for luminal cancers, but was the most sensitive for non-luminal cancers. PELP1 has low expression rate (< 20%) in Colorectal Cancers, Gastric Cancers and Renal Cell Carcinomas, but higher in Lung Cancers (49.1%) and Ovarian Cancers (42.3%). PELP1 expression was associated with poor outcome in non-luminal cancers and modified the prognostic effects of AR, suggesting the potential significance of NR co-regulator in prognostication.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-PELP1
ISOTYPE: IgG
CONTROL: Placenta, Breast, Colon, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma, Seminoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2649 | Tinto Predilute | 3.0 ml |
| BSB 2650 | Tinto Predilute | 7.0 ml |
| BSB 2651 | Tinto Predilute | 15.0 ml |
| BSB 2652 | Concentrate | 0.1 ml |
| BSB-2653 | Concentrate | 0.5 ml |
| BSB 2654 | Concentrate | 1.0 ml |
| BSB 2655 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3155 | Tinto Predilute | 3.0 ml |
| BSB 3156 | Tinto Predilute | 7.0 ml |
| BSB 3157 | Tinto Predilute | 15.0 ml |
| BSB 3158 | Concentrate | 0.1 ml |
| BSB 3159 | Concentrate | 0.5 ml |
| BSB 3160 | Concentrate | 1.0 ml |
| BSB 3161 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3744-3 | Tinto Predilute | 3.0 ml |
| BSB-3744-7 | Tinto Predilute | 7.0 ml |
| BSB-3744-15 | Tinto Predilute | 15.0 ml |
| BSB-3744-01 | Concentrate | 0.1 ml |
| BSB-3744-05 | Concentrate | 0.5 ml |
| BSB-3744-1 | Concentrate | 1.0 ml |
| BSB-3744-CS | Control Slides | 5 |

Perforin, MAb




IHC of Perforin on a FFPE Lymphoma Tissue

PGP 9.5, MAb



IHC of PGP 9.5 on a FFPE Testicle Tissue

PHOX2B, RMAb



IHC of PHOX2B on an FFPE Neuroblastoma Tissue

Perforin is a cytolytic protein found in the granules of Cytotoxic T lymphocytes and NK cells. Upon degranulation, perforin inserts itself into the target cell's plasma membrane, forming a pore. It enables granzymes to enter the target cells and activate apoptosis, the cell death program. Although some investigators report a cytolytic potential of CD4+ T cells, it appears more likely that CD8+ T cells are the major effector population in Th1- associated inflammatory skin diseases. The role of perforin-mediated cytotoxicity has been demonstrated in various autoimmune diseases. In vitro and in vivo studies suggest that the cytotoxicity of CTLs may be mediated by cytotoxic granules in certain inflammatory diseases in humans. In addition, it seems that T-cell cytotoxicity against keratinocytes is mediated by perforin in some inflammatory skin diseases.

Other authors suggest that perforin may have a dual role in alloimmune response (organ transplant applications). In one regard, it has a cytolytic function in acute rejection, and, in contrast, it may be responsible for downregulating both CD4- and CD8-mediated alloimmune response.

Protein gene product 9.5 (PGP 9.5), also known as ubiquitin carboxyl-terminal hydrolase-1 (UCH-L1), is a 27-kDa protein originally isolated from whole brain extracts (1). Although PGP9.5 expression in normal tissues was originally felt to be strictly confined to neurons and neuroendocrine cells (2), it has been subsequently documented in distal renal tubular epithelium, spermatogonia, Leydig cells, oocytes, melanocytes, prostatic secretory epithelium, ejaculatory duct cells, epididymis, mammary epithelial cells, Merkel cells, and dermal fibroblasts.

Immunostaining of a plethora of different mesenchymal neoplasms with this antibody has been demonstrated.

Paired-like homeobox 2b (PHOX2B), also known as neuroblastoma Phox (NBPhox), is a protein that in humans is encoded by the PHOX2B gene located on chromosome 4. It is expressed exclusively in the nervous system, in most neurons that control the viscera (cardiovascular, digestive and respiratory systems). It is also required for neuron differentiation. Mutations in human PHOX2B cause a rare disease of the autonomic nervous system (dysautonomia): congenital central hypoventilation syndrome (associated with respiratory arrests during sleep and, occasionally, wakefulness), Hirschsprung's disease (partial agenesis of the enteric nervous system), ROHHAD, and tumors of the sympathetic ganglia.

PHOX2 gene over-expression in Neuroblastoma (NB) tumors and cell lines suggests these genes may be widely involved in Neuroblastoma development through either a direct mechanism of up-regulation or a failure in maintaining proper transcript levels after embryonic development. The PHOX2B expression has been observed in all peripheral neuroblastic tumors, paragangliomas, and pheochromocytomas tested but in no other pediatric tumors among the 388 cases studied by expression microarray and the 109 cases studied by immunohistochemical analysis. PHOX2B and CD57 have been found to be useful markers of Neuroblastoma. PHOX2B is specific for Neuroblastoma in its differential diagnosis with other small round cell tumors, and its nuclear staining may be helpful for accurate bone marrow tumor quantification.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 5B10
ISOTYPE: IgG1
CONTROL: Spleen
LOCALIZATION: Cytoplasmic, Perinuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-46
ISOTYPE: IgG1/K
CONTROL: Brain, Testis, Colon, Pituitary, Nerve Tissue, Bowel Wall
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP312
ISOTYPE: IgG
CONTROL: Adrenal, Neuroblastoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2105 | Tinto Predilute | 3.0 ml |
| BSB 2106 | Tinto Predilute | 7.0 ml |
| BSB 2107 | Tinto Predilute | 15.0 ml |
| BSB 2108 | Concentrate | 0.1 ml |
| BSB 2109 | Concentrate | 0.5 ml |
| BSB 2110 | Concentrate | 1.0 ml |
| BSB 2111 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2112 | Tinto Predilute | 3.0 ml |
| BSB 2113 | Tinto Predilute | 7.0 ml |
| BSB 2114 | Tinto Predilute | 15.0 ml |
| BSB 2115 | Concentrate | 0.1 ml |
| BSB 2116 | Concentrate | 0.5 ml |
| BSB 2117 | Concentrate | 1.0 ml |
| BSB 2118 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3609 | Tinto Predilute | 3.0 ml |
| BSB 3610 | Tinto Predilute | 7.0 ml |
| BSB 3611 | Tinto Predilute | 15.0 ml |
| BSB 3612 | Concentrate | 0.1 ml |
| BSB 3613 | Concentrate | 0.5 ml |
| BSB 3614 | Concentrate | 1.0 ml |
| BSB 3615 | Control Slides | 5 |

PLA2R1, MAb

IHC of PLA2R1 on a FFPE Membranous Glomerulopathy Tissue

PLAP, MAb

IHC of PLAP on a FFPE Placenta Tissue

PLAP, RMAb

IHC of PLAP on a FFPE Placenta Tissue

PLA2R1 gene encodes phospholipase A2 receptor 1 protein, a 180 kDa transmembrane glycoproteins expressed by podocyte. Approximately 70% of patients with idiopathic membranous glomerulopathy have autoantibodies directed against podocyte PLA2R1. PLA2R1 also promotes tumor suppressive responses including senescence, apoptosis, and inhibition of transformation. Known oncogenes such as HIF2 α and c-Myc repress PLA2R1 expression.

PLA2R1 gain or loss of function experiments in vitro and in vivo shows that this receptor promotes several tumor suppressive responses including senescence, apoptosis and inhibition of transformation. Supporting a tumor suppressive role of PLA2R1, its expression decreases in numerous cancers, and known oncogenes such as HIF2 α and c-MYC repress its expression. PLA2R1 promoter methylation, a classical way to repress tumor suppressive gene expression in cancer cells, is observed in leukemia, in kidney and in breast cancer cells. PLA2R1 also promotes accumulation of reactive oxygen species which induce cell death and senescence. This review compiles recent data demonstrating an unexpected tumor suppressive role of PLA2R1 and outlines the future work needed to improve our knowledge of the functions of this gene in cancer

Alkaline phosphatase, placental type also known as placental alkaline phosphatase (PLAP) is an allosteric enzyme that in humans is encoded by the ALPP gene. PLAP is found in trophoblast cells of normal mature human placenta, Seminomas of testis and Ovarian Carcinomas. Detection of alkaline phosphatase in serum is a marker for Ovarian and Testicular Cancer.

This antibody reacts with a membrane-bound isoenzyme of placental alkaline phosphatase occurring in the placenta during the 3rd trimester of gestation.

This antibody immunoreacts with Germ Cell Tumors and can discriminate between these and other neoplasms. Somatic neoplasms (e.g., breast, gastrointestinal, prostatic and urinary cancers) may also immunoreact with antibodies to PLAP. PLAP positivity, in conjunction with keratin negativity, favors Seminoma over Carcinoma. Germ Cell Tumors are usually keratin positive but they regularly fail to stain with EMA, whereas most Carcinomas stain with anti-EMA. This antibody has shown cross-reaction with human intestinal alkaline phosphatase.

Placental Alkaline Phosphatase (PLAP) is found in trophoblast cells of normal mature human placenta, Seminomas of testis and Ovarian Carcinomas. Detection of alkaline phosphatase in serum is a marker for Ovarian and Testicular Cancer. This antibody reacts with a membrane-bound isoenzyme of placental alkaline phosphatase occurring in the placenta during the 3rd trimester of gestation.

This antibody immunoreacts with Germ Cell Tumors and can discriminate between these and other neoplasms. Somatic neoplasms (e.g., breast, gastrointestinal, prostatic and urinary cancers) may also immunoreact with antibodies to PLAP. PLAP positivity, in conjunction with keratin negativity, favors Seminoma over Carcinoma. Germ Cell Tumors are usually keratin positive but they regularly fail to stain with EMA, whereas most Carcinomas stain with anti-EMA. This antibody has shown cross-reaction with human intestinal alkaline phosphatase.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-129

ISOTYPE: IgG1

CONTROL: Brain, Testis, Kidney, Salivary Gland, Gastric GIST

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-47

ISOTYPE: IgG2b/K

CONTROL: Placenta, Testis, Seminomas, Ovarian Carcinomas

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP194

ISOTYPE: IgG

CONTROL: Placenta

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-2372-3 | Tinto Predilute | 3.0 ml |
| BSB-2372-7 | Tinto Predilute | 7.0 ml |
| BSB-2372-15 | Tinto Predilute | 15.0 ml |
| BSB-2372-01 | Concentrate | 0.1 ml |
| BSB-2372-05 | Concentrate | 0.5 ml |
| BSB-2372-1 | Concentrate | 1.0 ml |
| BSB-2372-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2754 | Tinto Predilute | 3.0 ml |
| BSB 2755 | Tinto Predilute | 7.0 ml |
| BSB 2756 | Tinto Predilute | 15.0 ml |
| BSB 2757 | Concentrate | 0.1 ml |
| BSB 2758 | Concentrate | 0.5 ml |
| BSB 2759 | Concentrate | 1.0 ml |
| BSB 2760 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5868 | Tinto Predilute | 3.0 ml |
| BSB 5869 | Tinto Predilute | 7.0 ml |
| BSB 5870 | Tinto Predilute | 15.0 ml |
| BSB 5871 | Concentrate | 0.1 ml |
| BSB 5872 | Concentrate | 0.5 ml |
| BSB 5873 | Concentrate | 1.0 ml |
| BSB 5874 | Control Slides | 5 |

PMS2, RMAb




IHC of PMS2 on a FFPE Colon Carcinoma Tissue

PMS2 is a gene that encodes for DNA repair proteins involved in mismatch repair. Carriers of the mismatch repair gene mutations have a high lifetime risk of developing Hereditary Non-Polyposis Colon Cancer (HNPCC) and several other cancers including endometrial cancer due to microsatellite instability (MSI) caused by accumulation of DNA replication errors in proliferating cells.

Along with MLH1, MSH2 and MSH6, PMS2 is helpful in diagnosing MSI. Tumors with low-level MSI show unfavorable pathological characteristics compared to tumors with none and tumors with high-level MSI.


Pneumocystis jirovecii, MAb



IHC of Pneumocystis on a FFPE Lung Tissue

Anti Pneumocystis carinii antibody reacts with an epitope on the yeast-like fungal microorganism, Pneumocystis carinii, that is resistant to formalin, picric acid, paraffin, as well as alcohol and xylene. No cross-reactivity has been demonstrated with other fungi or parasitic organisms.

Podoplanin/D2-40, MAb



IHC of Podoplanin/D2-40 on a FFPE Tonsil Tissue

Podoplanin is a transmembrane mucoprotein (38 kDa) recognized by the D2-40 monoclonal antibody. Podoplanin is specifically expressed in the endothelium of lymphatic capillaries but not in the blood vasculature. In normal skin and kidney, podoplanin is co-localized with VEGFR3/FLT4, another marker for lymphatic endothelial cells.

Podoplanin is selectively expressed in lymphatic endothelium as well as Lymphangiomas, Kaposi's Sarcomas and in subset Angiosarcomas with probable lymphatic differentiation. Podoplanin has also been shown to be expressed in Epithelioid Mesotheliomas, Hemangioblastomas and Seminomas.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP51
ISOTYPE: IgG
CONTROL: Colon Musoca, Colon Carcinoma, Breast, Skin, Cervix, Tonsil, TCC
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 3F6
ISOTYPE: IgM/K
CONTROL: Pneumocystis jirovecii Infected Tissue
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: D2-40
ISOTYPE: IgG1
CONTROL: Placenta, Breast, Lung, Cervix, Tonsil, Lymph Node, Lymphangioma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat, Mouse

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2119 | Tinto Predilute | 3.0 ml | BSB 5875 | Tinto Predilute | 3.0 ml | BSB 6064 | Tinto Predilute | 3.0 ml |
| BSB 2120 | Tinto Predilute | 7.0 ml | BSB 5876 | Tinto Predilute | 7.0 ml | BSB 6065 | Tinto Predilute | 7.0 ml |
| BSB 2121 | Tinto Predilute | 15.0 ml | BSB 5877 | Tinto Predilute | 15.0 ml | BSB 6066 | Tinto Predilute | 15.0 ml |
| BSB 2122 | Concentrate | 0.1 ml | BSB 5878 | Concentrate | 0.1 ml | BSB 6067 | Concentrate | 0.1 ml |
| BSB 2123 | Concentrate | 0.5 ml | BSB 5879 | Concentrate | 0.5 ml | BSB 6068 | Concentrate | 0.5 ml |
| BSB 2124 | Concentrate | 1.0 ml | BSB 5880 | Concentrate | 1.0 ml | BSB 6069 | Concentrate | 1.0 ml |
| BSB 2125 | Control Slides | 5 | BSB 5881 | Control Slides | 5 | BSB 6070 | Control Slides | 5 |

PRAME, RMAb



IHC of PRAME on a FFPE Melanoma Tissue

Melanoma antigen preferentially expressed in tumors is a protein that in humans is encoded by the PRAME gene. This gene encodes an antigen that is predominantly expressed in human melanomas and that is recognized by cytolytic T lymphocytes. It is not expressed in normal tissues, except in testis. This expression pattern is like that of other CT antigens, such as MAGE, BAGE and GAGE. However, unlike these other CT antigens, this gene is also expressed in acute leukemias. PRAME overexpression in triple negative breast cancer has also been found to promote cancer cell motility through induction of the epithelial-to-mesenchymal transition.

PRAME mRNA expression is well documented in cutaneous and ocular melanomas. One study concluded that diffuse nuclear immunoreactivity for PRAME was found in 87% of metastatic and 83.2% of primary melanomas. Among melanoma subtypes, PRAME was diffusely expressed in 94.4% of acral melanomas, 92.5% of superficial spreading melanomas, 90% of nodular melanomas, 88.6% of lentigo maligna melanomas, and 35% of desmoplastic melanomas. Most Melanocytic nevi (86.4%), were completely negative for PRAME. Immunoreactivity for PRAME was seen, albeit usually only in a minor subpopulation of lesional melanocytes, in 13.6% of cutaneous nevi, including dysplastic nevi, common acquired nevi, traumatized/recurrent nevi, and Spitz nevi. This study suggests that immunohistochemical analysis for PRAME expression may be useful for diagnostic purposes to support a suspected diagnosis of melanoma.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-PRAME


ISOTYPE: IgG

CONTROL: Testis, Seminoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

Prealbumin/Transthyretin, MAb



IHC of Prealbumin/Transthyretin on a FFPE Hepatocellular Carcinoma Tissue

TTR (TTR, transports thyroxine and retinol, or TBPA) encodes one of the three prealbumins including alpha-1-antitrypsin, transthyretin and orosomucoid. Transthyretin is a carrier protein; it transports thyroid hormones in the plasma and cerebrospinal fluid, and also transports retinol (vitamin A) in the plasma. The liver secretes Transthyretin into the blood, and the choroid plexus secretes TTR into the cerebrospinal fluid.

TTR mutations are associated with amyloid deposition, predominantly affecting peripheral nerves or the heart. TTR misfolding and aggregation is known to be associated with the Amyloid Diseases, Senile Systemic Amyloidosis (SSA), Familial Amyloid Polyneuropathy (FAP), and Familial Amyloid Cardiomyopathy (FAC). Transthyretin amyloidosis is a slowly progressive condition characterized by the buildup of abnormal deposits amyloid (Amyloidosis) in the body's organs and tissues.

It has been reported that TTR can be used as an immunohistochemical marker for choroid plexus papillomas, as well as carcinomas. The TTR gene has been found to be suppressed in Hepatic Carcinoma, where the TTR gene was found to be defective in its gene structure, which may have a relevance in its pathogenesis.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-125


ISOTYPE: IgG2b/K

CONTROL: Liver, Pancreas, Kidney, Leydig Cells, Hepatocellular Carcinoma, Seminoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Mouse, Rat

Progesterone Receptor, MAb



IHC of Progesterone Receptor on a FFPE Breast Carcinoma Tissue

The progesterone receptor (PR) also known as NR3C3 (nuclear receptor subfamily 3, group C, member 3), is an intracellular steroid receptor that specifically binds progesterone. PR is encoded by a single gene PGR residing on chromosome 11q22; it has two main forms, A and B, which differ in their molecular weight. Like all steroid receptors, the progesterone receptor has an amino and a carboxyl terminal, and between them the regulatory domain, a DNA binding domain, the hinge section, and the hormone binding domain.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB2

ISOTYPE: IgG1

CONTROL: Breast, Myometrium, Cervix, Breast Carcinoma

LOCALIZATION: Nuclear


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-2374-3 | Tinto Predilute | 3.0 ml |
| BSB-2374-7 | Tinto Predilute | 7.0 ml |
| BSB-2374-15 | Tinto Predilute | 15.0 ml |
| BSB-2374-01 | Concentrate | 0.1 ml |
| BSB-2374-05 | Concentrate | 0.5 ml |
| BSB-2374-1 | Concentrate | 1.0 ml |
| BSB-2374-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3616 | Tinto Predilute | 3.0 ml |
| BSB 3617 | Tinto Predilute | 7.0 ml |
| BSB 3618 | Tinto Predilute | 15.0 ml |
| BSB 3619 | Concentrate | 0.1 ml |
| BSB 3620 | Concentrate | 0.5 ml |
| BSB 3621 | Concentrate | 1.0 ml |
| BSB 3622 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2126 | Tinto Predilute | 3.0 ml |
| BSB 2127 | Tinto Predilute | 7.0 ml |
| BSB 2128 | Tinto Predilute | 15.0 ml |
| BSB 2129 | Concentrate | 0.1 ml |
| BSB 2130 | Concentrate | 0.5 ml |
| BSB 2131 | Concentrate | 1.0 ml |
| BSB 2132 | Control Slides | 5 |


Progesterone Receptor, RMAb



IHC of Progesterone Receptor on a FFPE Breast Carcinoma Tissue

The progesterone receptor (PR) also known as NR3C3 (nuclear receptor subfamily 3, group C, member 3), is an intracellular steroid receptor that specifically binds progesterone. PR is encoded by a single gene PGR residing on chromosome 11q22; it has two main forms, A and B, which differ in their molecular weight. Like all steroid receptors, the progesterone receptor has an amino and a carboxyl terminal, and between them the regulatory domain, a DNA binding domain, the hinge section, and the hormone binding domain.


Progesterone Receptor, RMAb



IHC of Progesterone Receptor on a FFPE Breast Carcinoma Tissue

The progesterone receptor (PR) also known as NR3C3 (nuclear receptor subfamily 3, group C, member 3), is an intracellular steroid receptor that specifically binds progesterone. PR is encoded by a single gene PGR residing on chromosome 11q22; it has two main forms, A and B, which differ in their molecular weight. Like all steroid receptors, the progesterone receptor has an amino and a carboxyl terminal, and between them the regulatory domain, a DNA binding domain, the hinge section, and the hormone binding domain.

Prolactin, MMab



IHC of Prolactin on a FFPE Pituitary Tissue

Prolactin is a peptide hormone primarily associated with lactation. It is synthesized and secreted by lactotrope cells in the adenohypophysis (anterior pituitary gland). It is also produced in other tissues including the breast and the decidua. Pituitary prolactin secretion is regulated by neuroendocrine neurons in the hypothalamus, most importantly by neurosecretory dopamine neurons of the arcuate nucleus, which inhibit prolactin secretion.

Prolactin is a useful marker in classification of pituitary tumors and the study of pituitary disease. It reacts with lactotrope cells.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT22

ISOTYPE: IgG

CONTROL: Breast, Myometrium, Cervix, Breast Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP2

ISOTYPE: IgG

CONTROL: Breast, Myometrium, Cervix, Breast Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: PRL02

ISOTYPE: IgG1/K

CONTROL: Normal Pituitary

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5882 | Tinto Predilute | 3.0 ml |
| BSB 5883 | Tinto Predilute | 7.0 ml |
| BSB 5884 | Tinto Predilute | 15.0 ml |
| BSB 5885 | Concentrate | 0.1 ml |
| BSB 5886 | Concentrate | 0.5 ml |
| BSB 5887 | Concentrate | 1.0 ml |
| BSB 5888 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2573 | Tinto Predilute | 3.0 ml |
| BSB 2574 | Tinto Predilute | 7.0 ml |
| BSB 2575 | Tinto Predilute | 15.0 ml |
| BSB 2576 | Concentrate | 0.1 ml |
| BSB 2577 | Concentrate | 0.5 ml |
| BSB 2578 | Concentrate | 1.0 ml |
| BSB 2579 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5889 | Tinto Predilute | 3.0 ml |
| BSB 5890 | Tinto Predilute | 7.0 ml |
| BSB 5891 | Tinto Predilute | 15.0 ml |
| BSB 5892 | Concentrate | 0.1 ml |
| BSB 5893 | Concentrate | 0.5 ml |
| BSB 5894 | Concentrate | 1.0 ml |
| BSB 5895 | Control Slides | 5 |

Prolactin, RMab



IHC of Prolactin on a FFPE Pituitary Tissue


Prolactin is a peptide hormone primarily associated with lactation. It is synthesized and secreted by lactotrope cells in the adenohypophysis (anterior pituitary gland). It is also produced in other tissues including the breast and the decidua. Pituitary prolactin secretion is regulated by neuroendocrine neurons in the hypothalamus, most importantly by neurosecretory dopamine neurons of the arcuate nucleus, which inhibit prolactin secretion.

Prolactin is a useful marker in classification of pituitary tumors and the study of pituitary disease. It reacts with lactotrope cells.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP193
ISOTYPE: IgG
CONTROL: Normal Pituitary
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2133 | Tinto Predilute | 3.0 ml |
| BSB 2134 | Tinto Predilute | 7.0 ml |
| BSB 2135 | Tinto Predilute | 15.0 ml |
| BSB 2136 | Concentrate | 0.1 ml |
| BSB 2137 | Concentrate | 0.5 ml |
| BSB 2138 | Concentrate | 1.0 ml |
| BSB 2139 | Control Slides | 5 |

Prostein/P501S, RMab



IHC of Prostein-P501S on a FFPE Prostate Tissue


Prostein/ P501S also known as Solute Carrier Family 45 member 3 (SLC45A3) and prostate cancer-associated protein 6, is a protein that in humans is encoded by the SLC45A3 gene. Prostein is expressed in prostate-specific normal tissues and at a significantly lower level in prostate tumor cell lines, but not in any other normal or malignant tissue examined to date. Prostein is considered to be a good marker to demonstrate prostatic origin in metastatic prostate cancer.

Prostein/P501S stain yields a perinuclear cytoplasmic (Golgi) distribution even in poorly differentiated tumors and metastases. The Immunohistochemistry for P501S is a sensitive and highly specific marker for identifying prostate metastases. The large majority of metastatic prostatic adenocarcinomas are Prostein/P501S positive (99%). A small subset of metastatic prostatic adenocarcinoma shows significant differences in staining intensity and extent for PSA and Prostein/P501S and, therefore, the combined use of these markers may result in increased sensitivity for detecting prostatic origin.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP381
ISOTYPE: IgG
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3162 | Tinto Predilute | 3.0 ml |
| BSB 3163 | Tinto Predilute | 7.0 ml |
| BSB 3164 | Tinto Predilute | 15.0 ml |
| BSB 3165 | Concentrate | 0.1 ml |
| BSB 3166 | Concentrate | 0.5 ml |
| BSB 3167 | Concentrate | 1.0 ml |
| BSB 3168 | Control Slides | 5 |

PSA, MMab



IHC of Prostate-Specific Antigen on a FFPE Prostatic Adenocarcinoma Tissue


Prostate-specific antigen (PSA) is a protein produced by the cells of the prostate gland. PSA is present in small quantities in the serum of normal men, and is often elevated in the presence of prostate cancer and in other prostate disorders. Higher than normal levels of PSA are associated with both localized and metastatic prostate cancer.

The PSA antibody recognizes primary and metastatic prostatic neoplasms but not tumors of nonprostatic origin. The antigen is a 33-34 kDa glycoprotein that is restricted to cells of prostatic origin. An immunohistochemical study showed more than 95% of prostatic carcinomas stained with PSA. PSA is demonstrable in the cytoplasm of acinar and ductal cells of normal or malignant prostate tissue.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-7
ISOTYPE: IgG1/K
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5896 | Tinto Predilute | 3.0 ml |
| BSB 5897 | Tinto Predilute | 7.0 ml |
| BSB 5898 | Tinto Predilute | 15.0 ml |
| BSB 5899 | Concentrate | 0.1 ml |
| BSB 5900 | Concentrate | 0.5 ml |
| BSB 5901 | Concentrate | 1.0 ml |
| BSB 5902 | Control Slides | 5 |

PSA, RMAb




IHC of Prostate-Specific Antigen on a FFPE Prostatic Adenocarcinoma Tissue

Prostate-specific antigen (PSA) is a protein produced by the cells of the prostate gland. PSA is present in small quantities in the serum of normal men, and is often elevated in the presence of prostate cancer and in other prostate disorders. Higher than normal levels of PSA are associated with both localized and metastatic prostate cancer.

The PSA antibody recognizes primary and metastatic prostatic neoplasms but not tumors of nonprostatic origin. The antigen is a 33-34 kDa glycoprotein that is restricted to cells of prostatic origin. An immunohistochemical study showed more than 95% of prostatic carcinomas stained with PSA. PSA is demonstrable in the cytoplasm of acinar and ductal cells of normal or malignant prostate tissue.

PSAP, MMAb




IHC of PSAP on a FFPE Prostate Tissue

Prostatic specific acid phosphatase (PSAP) is an enzyme produced by the prostate. It may be found in increased amounts in men who have prostate cancer or other diseases. The highest levels of acid phosphatase are found in metastasized prostate cancer. Diseases of the bone, such as Paget's disease or hyperparathyroidism, diseases of blood cells, (such as Sickle-Cell Disease), Multiple Myeloma or Lysosomal Storage Diseases, (such as Gaucher's disease), will show moderately increased levels. Certain medications can cause temporary increases or decreases in acid phosphatase levels. Manipulation of the prostate gland through massage, biopsy or rectal exam before a test may increase the levels of PSAP.

This antibody reacts with prostatic specific acid phosphatase in the glandular epithelium of the normal and Hyperplastic Prostate, Carcinoma of the prostate and metastatic cells of Prostatic Carcinoma. This marker may be helpful in pinpointing the site of origin in cases of Metastatic Carcinoma of the prostate, and is considered a more sensitive marker than PSA. However, it also offers less specificity.

PSAP, RMAb



IHC of PSAP on a FFPE Prostate Tissue

Prostatic specific acid phosphatase (PSAP) is an enzyme produced by the prostate. It may be found in increased amounts in men who have prostate cancer or other diseases. The highest levels of acid phosphatase are found in metastasized prostate cancer. Diseases of the bone, such as Paget's disease or hyperparathyroidism, diseases of blood cells, (such as Sickle-Cell Disease), Multiple Myeloma or Lysosomal Storage Diseases, (such as Gaucher's disease), will show moderately increased levels. Certain medications can cause temporary increases or decreases in acid phosphatase levels. Manipulation of the prostate gland through massage, biopsy or rectal exam before a test may increase the levels of PSAP.

This antibody reacts with prostatic specific acid phosphatase in the glandular epithelium of the normal and Hyperplastic Prostate, Carcinoma of the prostate and metastatic cells of Prostatic Carcinoma. This marker may be helpful in pinpointing the site of origin in cases of Metastatic Carcinoma of the prostate, and is considered a more sensitive marker than PSA. However, it also offers less specificity.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-PSA
ISOTYPE: IgG
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: PASE/4LJ
ISOTYPE: IgG1
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP53
ISOTYPE: IgG
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2140 | Tinto Predilute | 3.0 ml |
| BSB 2141 | Tinto Predilute | 7.0 ml |
| BSB 2142 | Tinto Predilute | 15.0 ml |
| BSB 2143 | Concentrate | 0.1 ml |
| BSB 2144 | Concentrate | 0.5 ml |
| BSB 2145 | Concentrate | 1.0 ml |
| BSB 2146 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5903 | Tinto Predilute | 3.0 ml |
| BSB 5904 | Tinto Predilute | 7.0 ml |
| BSB 5905 | Tinto Predilute | 15.0 ml |
| BSB 5906 | Concentrate | 0.1 ml |
| BSB 5907 | Concentrate | 0.5 ml |
| BSB 5908 | Concentrate | 1.0 ml |
| BSB 5909 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2147 | Tinto Predilute | 3.0 ml |
| BSB 2148 | Tinto Predilute | 7.0 ml |
| BSB 2149 | Tinto Predilute | 15.0 ml |
| BSB 2150 | Concentrate | 0.1 ml |
| BSB 2151 | Concentrate | 0.5 ml |
| BSB 2152 | Concentrate | 1.0 ml |
| BSB 2153 | Control Slides | 5 |

PSMA, RMAb

IHC of PSMA on a FFPE Prostate Adenocarcinoma Tissue

PSP94/MSMB, RMAb

IHC of PSP94/MSMB on a FFPE Prostate Tissue

PTEN, RMAb

IHC of PTEN on a FFPE Breast Carcinoma Tissue

PSMA, prostate specific membrane antigen, is a Type 2 integral membrane glycoprotein found in prostate and a few other tissues. Three functionally-distinct proteins are encoded, including folylpoly-gamma-glutamate carboxypeptidase in the intestine, N-acetylated alpha-linked acidic dipeptidase 1 in the brain and prostate-specific membrane antigen in the prostate. A mutation in the intestinal form may be associated with impaired intestinal absorption of dietary folates, resulting in low blood folate levels and consequent hyperhomocysteinemia. The form expressed in the brain may be involved in a number of pathological conditions associated with glutamate cytotoxicity. The prostate form is up-regulated in cancerous cells and is used as an effective diagnostic and prognostic indicator of prostate cancer. This gene likely arose from a duplication event of a nearby chromosomal region. Alternative splicing gives rise to multiple transcript variants.

Although PSMA expression is highest in the prostate, detectable levels of protein are also found in the small intestine and the brain. PSMA is expressed in prostate cancer cells as a noncovalently associated homodimer. Using a secreted form of the protein, it has been demonstrated that the extracellular domain is sufficient for dimerization and that dimerization is required for enzymatic activity. When used as an immunogen, dimeric (but not monomeric) PSMA is capable of efficiently eliciting antibodies that recognize PSMA-expressing tumor cells. It is a possible therapeutic target for prostate cancer and it is being used (with radioactive antibodies) to image prostate tissue.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP192
ISOTYPE: IgG
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Beta-microseminoprotein (MSMB) also called prostate secretory protein of 94 amino acids (PSP94), is one of the three predominant proteins secreted by the prostate gland and found in human seminal fluid along with prostate-specific antigen and prostatic specific acid phosphatase. Using exogenous MSMB, in vitro and in vivo studies indicate that MSMB may have several anti-tumor effects on prostate tumor cells.

MSMB expression is high in normal prostate epithelial cells, but is decreased in prostate cancer cells. Studies have shown that MSMB is a strong independent factor indicating favorable outcome after radical prostatectomy for localized prostate cancer.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP203
ISOTYPE: IgG
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Phosphatase and tensin homolog (PTEN) is a protein that, in humans, is encoded by the PTEN gene. Unlike most of the protein tyrosine phosphatases, this protein preferentially dephosphorylates phosphoinositide substrates. It negatively regulates intracellular levels of phosphatidylinositol-3,4,5-trisphosphate in cells and functions as a tumor suppressor by negatively regulating Akt/PKB signaling pathway.

PTEN is one of the most commonly lost tumor suppressors in human cancer; in fact, up to 70% of men with prostate cancer are estimated to have lost a copy of the PTEN gene at the time of diagnosis. During tumor development, mutations and deletions of PTEN occur that inactivate its enzymatic activity leading to increased cell proliferation and reduced cell death. Frequent genetic inactivation of PTEN occurs in glioblastoma, endometrial cancer, and prostate cancer; and reduced expression is found in many other tumor types such as lung and breast cancer. In breast and prostate cancer, loss of PTEN expression has been shown to correlate positively with advanced stage. Furthermore, PTEN mutation also causes a variety of inherited predispositions to cancer.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-PTEN
ISOTYPE: IgG
CONTROL: Breast, Prostate, Breast & Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6345 | Tinto Predilute | 3.0 ml |
| BSB 6346 | Tinto Predilute | 7.0 ml |
| BSB 6347 | Tinto Predilute | 15.0 ml |
| BSB 6348 | Concentrate | 0.1 ml |
| BSB 6349 | Concentrate | 0.5 ml |
| BSB 6350 | Concentrate | 1.0 ml |
| BSB 6351 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2419 | Tinto Predilute | 3.0 ml |
| BSB 2420 | Tinto Predilute | 7.0 ml |
| BSB 2421 | Tinto Predilute | 15.0 ml |
| BSB 2422 | Concentrate | 0.1 ml |
| BSB 2423 | Concentrate | 0.5 ml |
| BSB 2424 | Concentrate | 1.0 ml |
| BSB 2425 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2643 | Tinto Predilute | 3.0 ml |
| BSB 2644 | Tinto Predilute | 7.0 ml |
| BSB 2645 | Tinto Predilute | 15.0 ml |
| BSB 2646 | Concentrate | 0.1 ml |
| BSB 2647 | Concentrate | 0.5 ml |
| BSB 2648 | Concentrate | 1.0 ml |
| BSB 2649 | Control Slides | 5 |

PTEN, RMab




IHC of PTEN on FFPE Prostate Adenocarcinoma Tissue

Phosphatase and tensin homolog (PTEN) is a protein that, in humans, is encoded by the PTEN gene. Unlike most of the protein tyrosine phosphatases, this protein preferentially dephosphorylates phosphoinositide substrates. It negatively regulates intracellular levels of phosphatidylinositol-3,4,5-trisphosphate in cells and functions as a tumor suppressor by negatively regulating Akt/PKB signaling pathway.

PTEN is one of the most commonly lost tumor suppressors in human cancer; in fact, up to 70% of men with prostate cancer are estimated to have lost a copy of the PTEN gene at the time of diagnosis. During tumor development, mutations and deletions of PTEN occur that inactivate its enzymatic activity leading to increased cell proliferation and reduced cell death. Frequent genetic inactivation of PTEN occurs in glioblastoma, endometrial cancer, and prostate cancer; and reduced expression is found in many other tumor types such as lung and breast cancer. In breast and prostate cancer, loss of PTEN expression has been shown to correlate positively with advanced stage. Furthermore, PTEN mutation also causes a variety of inherited predispositions to cancer.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM265
ISOTYPE: IgG
CONTROL: Colon, Thymus, Skin, Kidney, Breast, Prostate, Breast Carcinoma, Prostatic Carcinoma
LOCALIZATION: Nuclear, Cytoplasmic
SPECIES REACTIVITY: Human

PU.1, RMab




IHC of PU.1 on a FFPE Lymphoma Tissue

PU.1 is a member of the Ets family of transcription factors and is required for the development of multiple hematopoietic lineages. It plays a pivotal role in normal myeloid differentiation and regulates the expression of immunoglobulin and other genes that are important for B-cell development. It is expressed in the myeloid lineage and in immature as well as mature B lymphocytes, with the exception of plasma cells.

PU.1 is expressed in germinal center B-cells and mantle B-cells. The antibody is positive in various lymphomas including B-Chronic Lymphocytic Leukemia, Mantle Cell Lymphoma, Follicular Lymphoma, Marginal Zone Lymphoma, Burkitt Lymphoma, Diffuse Large Cell Lymphoma, Diffuse Large B-cell Lymphoma, T-cell rich B-cell Lymphoma, Nodular Lymphocyte Predominant Hodgkin Lymphoma. It has been demonstrated that a high level of expression of GC antigens (including PU.1) has a positive association with longer overall survival and progression free survival in the case of Follicular Lymphoma.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP18
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

Pygopus 2/Pygo 2, MAb



IHC of Pygopus 2/Pygo 2 on a FFPE Lung Neuroendocrine Tissue

Pygopus Family PHD Finger 2 (PYGO2) is a protein-coding gene associated with G-protein coupled receptor (GPCR) and Wnt signaling pathways, chromatin binding and histone acetyltransferase regulator activity. Pygo2 as a Wnt signaling pathway component has been detected in multiple cancer types. One study found that abnormal Pygo2 expression was associated with poor differentiation and a high tumor, node, and metastases stage and poor prognosis in non-small cell lung cancer (NSCLC) patients, therefore abnormal Pygo2 protein expression may be a marker for advanced NSCLC. Another study found that 59% of the patient tumor specimens exhibited positive Pygo2 immunohistochemistry staining and increased intensity with the grade of malignancy, especially for WHO grade III and IV.

High-level expression has been reported in gliomas compared with normal brain tissues, which suggest an important role of Pygo2 in brain tumor progression. In a Colorectal Cancer study, the expression pattern of Pygo2 was evaluated by immunohistochemistry in tumor tissues and their normal margins, and found the expression of Pygo2 protein was detected in all tumor tissues. Furthermore, this expression was significantly higher in Colorectal Cancer samples than in normal tissues. Pygo2 is also associated with Oligospermia.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-156
ISOTYPE: IgG2a
CONTROL: Fallopian Tube, Adrenal Gland, Kidney, Transitional Cell Carcinoma, Lung Adenocarcinoma, Lung Neuroendocrine Cancer, Papillary Thyroid Carcinoma
LOCALIZATION: Cytoplasmic, Membranous, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3769-3 | Tinto Predilute | 3.0 ml |
| BSB-3769-7 | Tinto Predilute | 7.0 ml |
| BSB-3769-15 | Tinto Predilute | 15.0 ml |
| BSB-3769-01 | Concentrate | 0.1 ml |
| BSB-3769-05 | Concentrate | 0.5 ml |
| BSB-3769-1 | Concentrate | 1.0 ml |
| BSB-3769-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2154 | Tinto Predilute | 3.0 ml |
| BSB 2155 | Tinto Predilute | 7.0 ml |
| BSB 2156 | Tinto Predilute | 15.0 ml |
| BSB 2157 | Concentrate | 0.1 ml |
| BSB 2158 | Concentrate | 0.5 ml |
| BSB 2159 | Concentrate | 1.0 ml |
| BSB 2160 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3745-3 | Tinto Predilute | 3.0 ml |
| BSB-3745-7 | Tinto Predilute | 7.0 ml |
| BSB-3745-15 | Tinto Predilute | 15.0 ml |
| BSB-3745-01 | Concentrate | 0.1 ml |
| BSB-3745-05 | Concentrate | 0.5 ml |
| BSB-3745-1 | Concentrate | 1.0 ml |
| BSB-3745-CS | Control Slides | 5 |

Renal Cell Carcinoma, MAb



IHC of Renal Cell Carcinoma on a FFPE Kidney Tissue

Renal Cell Carcinoma, also known as a Gurnistial Tumor, is the most common form of kidney cancer arising from the renal tubule. It is also the most common type of kidney cancer in adults. Initial treatment is surgery because it is notoriously resistant to radiation therapy and chemotherapy, although some cases respond to immunotherapy.

Renal Cell Carcinoma antibody recognizes a 200 kDa glycoprotein localized in the brush border of the proximal renal tubule. This antibody immunoreacts with approximately 90% of Primary Renal Cell Carcinomas and approximately 85% of Metastatic Renal Cell Carcinomas. Other tumors that may react with this antibody are Parathyroid Adenoma, an occasional Breast Carcinoma, Nephroblastoma, Oncocytoma, Mesoblastic Nephroma, Transitional Cell Carcinoma, and Angiomyolipoma are not labeled with this antibody.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: PN-15


ISOTYPE: IgG1/K

CONTROL: Kidney, Breast, Thyroid, Renal Cell Carcinoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Rat

Retinoblastoma/Rb, MAb



IHC of Retinoblastoma on a FFPE Colon Carcinoma Tissue

The retinoblastoma protein (Rb) is a tumor-suppressor protein that is dysfunctional in many types of cancer. One highly studied function of Rb is to prevent excessive cell growth by inhibiting cell-cycle progression until a cell is ready to divide. Rb prevents the cell from replicating damaged DNA by preventing its progression along the cell cycle through G1 into S.

Should an oncogenic protein (such as that produced by cells infected with high-risk types of human papillomaviruses, SV40 or Adenoviruses) bind and inactivate Rb, this can lead to cancer. Rb protein may act by regulating transcription; loss of its function leads to uncontrolled cell growth. Aberrations in the Rb gene have been implicated in cancers of breast, colon, prostate, kidney, nasopharynx, and Leukemia.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 1F8


ISOTYPE: IgG1

CONTROL: Colon, Breast, Skin, Fallopian Tube, Tonsil, Colon Carcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ROS-1, RMAb



IHC of ROS-1 on a FFPE Non-Small Cell Lung Cancer Tissue

Repressor of Silencing 1 (ROS1) is a receptor tyrosine kinase that undergoes genetic rearrangements in various human cancers and in humans is encoded by the ROS1 gene. The protein encoded by this gene is a type I integral membrane protein with tyrosine kinase activity with structural similarity to the anaplastic lymphoma kinase (ALK) protein. The protein may function as a growth or differentiation factor receptor. ROS1 expression is limited in normal tissues to occasional staining cerebellum, stomach, small intestine, colon and kidney.

Gene rearrangements involving the ROS1 gene were first detected in glioblastoma tumors and cell lines. ROS1 fusion partners include CD74, SLC34A2 and SDC4, leading to oncogenic transformation. ROS1 rearrangement was identified in a cell line derived from a lung adenocarcinoma patient and multiple studies have demonstrated its incidence in lung cancers. While ROS1 is undetectable in the normal lung, studies have described ROS1 rearrangements in 1-2% of NSCLC by FISH. Recent reports have demonstrated strong correlation between ROS1 IHC with FISH positivity. ROS1 fusions have been detected in multiple other tumors, including glioblastoma, non-small cell lung cancer (NSCLC), cholangiocarcinoma, ovarian cancer, gastric adenocarcinoma, colorectal cancer, inflammatory myofibroblastic tumor, angiosarcoma, and epithelioid hemangioendothelioma.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP282

ISOTYPE: IgG

CONTROL: Placenta, Lung, SiHa Cells, NSCL ROS1 +

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5910 | Tinto Predilute | 3.0 ml |
| BSB 5911 | Tinto Predilute | 7.0 ml |
| BSB 5912 | Tinto Predilute | 15.0 ml |
| BSB 5913 | Concentrate | 0.1 ml |
| BSB 5914 | Concentrate | 0.5 ml |
| BSB 5915 | Concentrate | 1.0 ml |
| BSB 5916 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6127 | Tinto Predilute | 3.0 ml |
| BSB 6128 | Tinto Predilute | 7.0 ml |
| BSB 6129 | Tinto Predilute | 15.0 ml |
| BSB 6130 | Concentrate | 0.1 ml |
| BSB 6131 | Concentrate | 0.5 ml |
| BSB 6132 | Concentrate | 1.0 ml |
| BSB 6133 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3623 | Tinto Predilute | 3.0 ml |
| BSB 3624 | Tinto Predilute | 7.0 ml |
| BSB 3625 | Tinto Predilute | 15.0 ml |
| BSB 3626 | Concentrate | 0.1 ml |
| BSB 3627 | Concentrate | 0.5 ml |
| BSB 3628 | Concentrate | 1.0 ml |
| BSB 3629 | Control Slides | 5 |

S-100, MAb




IHC of S-100 on a FFPE Malignant Melanoma Tissue

S100A1, RMAb



IHC of S100A1 on a FFPE Melanoma Tissue

S100A6, RMAb



IHC of S100A6 on a FFPE Ovarian Carcinoma Tissue

S-100 protein is a type of low-molecular weight protein found in vertebrates, characterized by two calcium-binding sites of the helix-loop-helix conformation. S-100 is normally present in cells derived from the neural crest (Schwann cells, melanocytes and glial cells), chondrocytes, adipocytes, myoepithelial cells, macrophages, Langerhans cells, dendritic cells, and keratinocytes. It may be present in some breast epithelial cells. Several members of the S-100 protein family are useful as markers for certain tumors and epidermal differentiation. The S-100 protein can be found in melanomas, 50% of Malignant Peripheral Nerve Sheath Tumors, and Clear Cell Sarcomas.

Almost all Malignant Melanomas and cases of Histiocytosis X are positive for S-100 protein. Despite the fact that S-100 protein is a ubiquitous substance, its demonstration is of great value in the identification of several neoplasms, particularly Melanomas.

The S100A1 protein is a member of the S100 family of proteins containing 4 EF-hand calcium-binding motifs in its dimerized form. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100A1 may function in stimulation of Ca²⁺-induced Ca²⁺ release, inhibition of microtubule assembly and inhibition of protein kinase C-mediated phosphorylation.

In normal tissues, anti-S100A1 is expressed in cardiac muscle, skeletal muscle and neuronal cells. Reduced expression of S100A1 has been implicated in cardiomyopathies. It can also be useful in distinguishing between Renal Oncocytomas and Clear Cell Renal Cell Carcinomas (positive) and Papillary Renal Cell Carcinomas (positive) from Chromophobe Renal Cell Carcinomas (negative). It is a specific and sensitive marker for Nephrogenic Adenoma.

S100 calcium-binding protein A6 (S100A6) is a protein that in humans is encoded by the S100A6 gene. The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21. This protein may function in stimulation of Ca²⁺-dependent insulin release, stimulation of prolactin secretion, and exocytosis. Chromosomal rearrangements and altered expression of this gene have been implicated in melanoma.

S100A6 is a cytoplasmic and nuclear protein abundantly expressed in fibroblasts and epithelial cells. It is also detected in some neurons, glial cells, smooth muscle, myocytes, and lymphocytes. In pancreatic cancer, elevation of S100A6 RNA and protein has been observed in malignant cells. Nuclear S100A6 expression is associated with reduced survival time in pancreatic cancer patients. S100A6 expression is also significantly correlated with melanoma metastases and survival times. It has been shown to be upregulated in 50% of gastric cancers, and highly expressed at the invasive margins of colorectal carcinoma. S100A6 has also been reported as possible diagnostic marker of papillary thyroid carcinoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 4C4.9
ISOTYPE: IgG2a
CONTROL: Melanoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Dog, Cat, Mouse, Rat, Cattle


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP184
ISOTYPE: IgG
CONTROL: Brain, Breast, Thyroid, Tonsil, Pancreas, Salivary Gland, Renal Oncocytomas, Clear Cell Carcinoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP313
ISOTYPE: IgG
CONTROL: Testis, Kidney, Breast, Tonsil, Colon, Carcinomas
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5917 | Tinto Predilute | 3.0 ml |
| BSB 5918 | Tinto Predilute | 7.0 ml |
| BSB 5919 | Tinto Predilute | 15.0 ml |
| BSB 5920 | Concentrate | 0.1 ml |
| BSB 5921 | Concentrate | 0.5 ml |
| BSB 5922 | Concentrate | 1.0 ml |
| BSB 5923 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2161 | Tinto Predilute | 3.0 ml |
| BSB 2162 | Tinto Predilute | 7.0 ml |
| BSB 2163 | Tinto Predilute | 15.0 ml |
| BSB 2164 | Concentrate | 0.1 ml |
| BSB 2165 | Concentrate | 0.5 ml |
| BSB 2166 | Concentrate | 1.0 ml |
| BSB 2167 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3176 | Tinto Predilute | 3.0 ml |
| BSB 3177 | Tinto Predilute | 7.0 ml |
| BSB 3178 | Tinto Predilute | 15.0 ml |
| BSB 3179 | Concentrate | 0.1 ml |
| BSB 3180 | Concentrate | 0.5 ml |
| BSB 3181 | Concentrate | 1.0 ml |
| BSB 3182 | Control Slides | 5 |

S100A8/MRP8, R Mab

IHC of S100A8/MRP8 on a FFPE Kidney Transplant Tissue

S100A9, R Mab

IHC of S100A9 on a FFPE Melanoma Tissue

S100 Beta, R Mab

IHC of S100 Beta on a FFPE Melanoma Tissue

Myeloid Related Protein 8 (MRP8), also known as S100A8, is a calcium binding protein that belongs to the S100 family. By a Ca²⁺ dependent manner, S100A8/A9 forms Calprotectin, a heterodimeric inflammatory mediator of inflammation found in the cytoplasm of neutrophils and expressed on the membrane of monocytes. S100A8 is expressed during myeloid differentiation and chronic inflammations, and it is expressed constitutively or induced in epithelial cells during dermatose.

S100A8 is expressed in cells with myeloid origin, including granulocytes, monocytes and macrophages, and it is observed in blood granulocytes and monocytes. It is also expressed in infiltrate macrophages during inflammatory reactions, but not in normal tissue macrophages. S100A8 also reacts with activated microglial cells in human cerebral malaria. In tumors, positive staining of S100A8 has been observed in various cancers including pancreatic cancer, and it has been linked to inflammation-associated cancers.

S100A9, also known as migration inhibitory factor-related protein 14 (MRP-14) or calgranulin-B is a member of the S100 family of proteins containing 2 EF hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100A9 forms a heterodimer, Calprotectin, with S100A8 in a calcium-dependent manner. S100A9 may function in the inhibition of casein kinase.

S100A9 is expressed in granulocytes, monocytes in peripheral blood and in infiltrating macrophages in inflammatory sites, but not in normal tissue macrophages. Elevated plasma levels of S100A9 has been observed in inflammatory disorders such as chronic bronchitis, cystic fibrosis and rheumatoid arthritis. S100A9 is also detected in tumor cells in carcinomas of the liver, lung, breast and thyroid. It is correlated with tumor differentiation.

S100 calcium binding protein B or S100 Beta is a member of the S100 family. S100 proteins are localized in the cytoplasm and nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 Beta is abundant in glial cells of the central and peripheral nervous system, in melanocytes, chondrocytes, and adipocytes.

Anti-S100 Beta labels langerhans cells, histiocytes, epithelial, myoepithelial cells and integrating reticular cells of lymphoid tissue, and tumors originated from these cells. It is a useful marker for diagnosis of melanoma and tumors of the nervous system.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP90

ISOTYPE: IgG

CONTROL: Tonsil, Breast, Liver Fallopian Tube, Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP185

ISOTYPE: IgG

CONTROL: Tonsil, Liver, Lung, Breast, Cervix, Bone Marrow, Bladder TCC, Thyroid Carcinomas

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP32

ISOTYPE: IgG


CONTROL: Melanoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: predicted Human, predicted Mouse, predicted Rat, predicted Goat, predicted Zebrafish, predicted Macaque Monkey


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2356 | Tinto Predilute | 3.0 ml | BSB 2168 | Tinto Predilute | 3.0 ml | BSB 2175 | Tinto Predilute | 3.0 ml |
| BSB 2357 | Tinto Predilute | 7.0 ml | BSB 2169 | Tinto Predilute | 7.0 ml | BSB 2176 | Tinto Predilute | 7.0 ml |
| BSB 2358 | Tinto Predilute | 15.0 ml | BSB 2170 | Tinto Predilute | 15.0 ml | BSB 2177 | Tinto Predilute | 15.0 ml |
| BSB 2359 | Concentrate | 0.1 ml | BSB 2171 | Concentrate | 0.1 ml | BSB 2178 | Concentrate | 0.1 ml |
| BSB 2360 | Concentrate | 0.5 ml | BSB 2172 | Concentrate | 0.5 ml | BSB 2179 | Concentrate | 0.5 ml |
| BSB 2361 | Concentrate | 1.0 ml | BSB 2173 | Concentrate | 1.0 ml | BSB 2180 | Concentrate | 1.0 ml |
| BSB 2362 | Control Slides | 5 | BSB 2174 | Control Slides | 5 | BSB 2181 | Control Slides | 5 |

S100P, RMAb




IHC of S100P on a FFPE Melanoma Tissue

SALL4, MAb



IHC of SALL4 on a FFPE Testicular Cancer Metastasis to Liver Tissue

SALL4, RMAb



IHC of SALL4 on a FFPE Testicular Cancer Metastasis to Liver Tissue

S100P is a member of the S100 family of proteins containing 2 EF-hand calcium binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100P is expressed in various normal tissues including placenta, bladder, spleen, gastric and intestinal mucosa. Overexpression of S100P has been detected in several cancers such as colon, prostate, pancreatic and lung carcinomas. It has been functionally implicated in carcinogenic processes.

S100P is an early development marker of pancreatic carcinogenesis and can be used as a marker for pancreatic ductal adenocarcinoma. It may also serve as a predictor of distant metastasis and poor survival in non-small cell lung carcinomas.

Sal-like protein 4 (SALL4) is a transcription factor encoded by a member of the Spalt-like (SALL) gene family, SALL4. SALL4 is reactivated and misregulated in various cancer, such as acute myeloid leukemia (AML), B-cell acute lymphocytic leukemia (B-ALL), germ cell tumors, gastric cancer, breast cancer, hepatocellular carcinoma (HCC), lung cancer, and glioma. In many of these cancers, SALL4 expression has been compared in tumor cells to the normal tissue counterpart, e.g. it is expressed in nearly half of primary human endometrial cancer samples, but not in normal or hyperplastic endometrial tissue samples.

Often, SALL4 expression is correlated with worse survival and poor prognosis such as in HCC, or with metastasis such as in endometrial cancer, colorectal carcinoma, and esophageal squamous cell carcinoma. In solid tumors such as germ cell tumors, SALL4 protein expression has become a standard diagnostic biomarker.

SALL4 demonstrates 100% sensitivity and stains more than 90% tumor cells in all intratubular germ cell neoplasia, seminomas, dysgerminomas, embryonal carcinomas, and yolk sac tumor (YST) (both pediatric and postpubertal). SALL4 is also positive in most cases of teratoma and the mononucleated trophoblastic cells in choriocarcinomas. Most non-testicular tumors from various organs and sites are negative for SALL4, though an occasional carcinoma or sarcoma may show weak SALL4 staining in less than 25% of tumor cells.

Sal-like protein 4 (SALL4) is a transcription factor encoded by a member of the Spalt-like (SALL) gene family, SALL4. SALL4 is reactivated and misregulated in various cancer, such as acute myeloid leukemia (AML), B-cell acute lymphocytic leukemia (B-ALL), germ cell tumors, gastric cancer, breast cancer, hepatocellular carcinoma (HCC), lung cancer, and glioma. In many of these cancers, SALL4 expression has been compared in tumor cells to the normal tissue counterpart, e.g. it is expressed in nearly half of primary human endometrial cancer samples, but not in normal or hyperplastic endometrial tissue samples.

Often, SALL4 expression is correlated with worse survival and poor prognosis such as in HCC, or with metastasis such as in endometrial cancer, colorectal carcinoma, and esophageal squamous cell carcinoma. In solid tumors such as germ cell tumors, SALL4 protein expression has become a standard diagnostic biomarker.

SALL4 demonstrates 100% sensitivity and stains more than 90% tumor cells in all intratubular germ cell neoplasia, seminomas, dysgerminomas, embryonal carcinomas, and yolk sac tumor (YST) (both pediatric and postpubertal). SALL4 is also positive in most cases of teratoma and the mononucleated trophoblastic cells in choriocarcinomas. Most non-testicular tumors from various organs and sites are negative for SALL4, though an occasional carcinoma or sarcoma may show weak SALL4 staining in less than 25% of tumor cells.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP186
ISOTYPE: IgG
CONTROL: Colon, Prostate, Pancreatic and Lung Carcinomas
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: 6E3
ISOTYPE: IgG1/K
CONTROL: Testis, Seminoma, Yolk Sac Tumor
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP299
ISOTYPE: IgG1/K
CONTROL: Testis, Seminoma, Yolk Sac Tumor
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2182 | Tinto Predilute | 3.0 ml |
| BSB 2183 | Tinto Predilute | 7.0 ml |
| BSB 2184 | Tinto Predilute | 15.0 ml |
| BSB 2185 | Concentrate | 0.1 ml |
| BSB 2186 | Concentrate | 0.5 ml |
| BSB 2187 | Concentrate | 1.0 ml |
| BSB 2188 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3183 | Tinto Predilute | 3.0 ml |
| BSB 3184 | Tinto Predilute | 7.0 ml |
| BSB 3185 | Tinto Predilute | 15.0 ml |
| BSB 3186 | Concentrate | 0.1 ml |
| BSB 3187 | Concentrate | 0.5 ml |
| BSB 3188 | Concentrate | 1.0 ml |
| BSB 3189 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3190 | Tinto Predilute | 3.0 ml |
| BSB 3191 | Tinto Predilute | 7.0 ml |
| BSB 3192 | Tinto Predilute | 15.0 ml |
| BSB 3193 | Concentrate | 0.1 ml |
| BSB 3194 | Concentrate | 0.5 ml |
| BSB 3195 | Concentrate | 1.0 ml |
| BSB 3196 | Control Slides | 5 |

SARS-CoV-2, MMab


IHC of SARS-CoV-2 on a FFPE Infected Lung Tissue

The Severe Acute Respiratory Syndrome 2 virus (SARS-CoV-2) is a betacoronavirus first isolated in Wuhan, China, in late 2019. The virus has a 29.8 kbp genome, encoding the membrane, envelope, nucleocapsid, and spike glycoprotein. The spike proteins are cleaved by TMPRSS2 serine protease, then the Receptor Binding Domain of the spike protein binds with ACE2 or CD147 to enter the cell.

The SARS-CoV-2 virus has been shown to infect the tracheal and lung epithelium, GI tract, and olfactory neuron, brain, bone marrow and possibly other organs. Cough, fever, and trouble breathing are the main symptoms, although GI distress, fatigue, and neurological distress are also common. Severe symptoms are more likely to appear in patients with advanced age and/or preexisting cardiovascular disease or diabetes. The virus has a 2-11 day incubation period and mortality rate around 2.5%. Severe symptoms include diffuse alveolar damage in the lungs, hyaline membrane formation, microthrombi in the lungs, heart, and brain, and extreme inflammation as "cytokine storms" that flood the body with cytokines (elevated 1L1, IL-6, IL-8, and TNFa among others) and immune cells (especially CD4+ and CD8+ T cells and CD68+ and CD163+ Macrophages). This antibody cross-reacts with SARS-CoV nucleocapsid protein, but not with MERS, 229E, or OC43 coronaviruses.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-134
ISOTYPE: IgG2b
CONTROL: SARS-CoV-2 Infected Tissues
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3701-3 | Tinto Predilute | 3.0 ml |
| BSB-3701-7 | Tinto Predilute | 7.0 ml |
| BSB-3701-15 | Tinto Predilute | 15.0 ml |
| BSB-3701-01 | Concentrate | 0.1 ml |
| BSB-3701-05 | Concentrate | 0.5 ml |
| BSB-3701-1 | Concentrate | 1.0 ml |
| BSB-3701-CS | Control Slides | 5 |

SATB2, RMAb


IHC of SATB2 on a FFPE Colon Carcinoma Tissue

Special AT-rich sequence-binding protein 2 (SATB2) also known as DNA-binding protein SATB2 is a protein that in humans is encoded by the SATB2 gene. SATB2 specifically binds nuclear matrix attachment regions and is involved in transcriptional regulation and chromatin remodeling. SATB2 has been implicated as causative in the cleft or high palate of individuals with 2q32q33 microdeletion syndrome.

SATB2 has been identified as a tissue-specific protein when screening protein expression patterns in human and cancerous tissues, with expression restricted to the lower gastrointestinal tract. SATB2 in combination with CK20 and Cadherin 17 could identify almost all colorectal carcinomas, including poorly differentiated colorectal carcinomas. Upper gastrointestinal (GI) carcinomas and pancreatic ductal carcinomas are usually negative for SATB2, and ovarian carcinomas, lung adenocarcinomas, and adenocarcinomas from other origin are rarely positive for SATB2. Therefore, SATB2 is a good marker for identifying a carcinoma of colorectal origin when working on a tumor of unknown primary. Another potential utility of SATB2 is to identify neuroendocrine neoplasms/carcinomas of the colon and rectum because SATB2 is usually negative in other neuroendocrine neoplasms of the GI tract, pancreas, and lung. SATB2 has been also shown to be a sensitive marker of osteoblastic differentiation in benign and malignant mesenchymal tumors.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP281
ISOTYPE: IgG
CONTROL: Colon, Brain, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3197 | Tinto Predilute | 3.0 ml |
| BSB 3198 | Tinto Predilute | 7.0 ml |
| BSB 3199 | Tinto Predilute | 15.0 ml |
| BSB 3200 | Concentrate | 0.1 ml |
| BSB 3201 | Concentrate | 0.5 ml |
| BSB 3202 | Concentrate | 1.0 ml |
| BSB 3203 | Control Slides | 5 |

SDHB, MMab


IHC SDHB on a FFPE Hepatocellular Carcinoma Tissue

Mutations in the tumor suppressor genes SDHB, SDHC, and SDHD (or collectively SDHx) cause the inherited paraganglioma syndromes, characterized by pheochromocytomas and paragangliomas. The IHC for SDHB is negative in all SDH mutated paragangliomas regardless of whether the B, C or D subunit is involved. However, other tumors have been associated with SDHx mutations, such as Gastrointestinal Stromal Tumors, specifically in the context of Carney-Stratakis syndrome. It has been shown that SDHB immunohistochemistry is a reliable technique for the identification of pheochromocytomas and paragangliomas caused by SDHx mutations. It's been shown that Carney-Stratakis syndrome- and Carney-triad-associated GISTs are negative by immunohistochemistry for SDHB in contrast to KIT- or PDGFRA-mutated GISTs and a majority of sporadic GISTs, and it has been suggested that GISTs of epithelioid cell morphology are tested for SDHB immunohistochemically.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-131
ISOTYPE: IgG1/K
CONTROL: Breast, Adrenal, Prostate, Kidney, Spleen, Tonsil, Breast Carcinoma, Hepatocellular Carcinoma, Lung Adeno Carcinoma, Prostate Carcinoma, Papillary Thyroid Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|--------------|-----------------|---------|
| BSB-2375-3 | Tinto Predilute | 3.0 ml |
| BSB-2375-7 | Tinto Predilute | 7.0 ml |
| BSB-2375-15 | Tinto Predilute | 15.0 ml |
| BSB-2375-01 | Concentrate | 0.1 ml |
| BSB-2375-05 | Concentrate | 0.5 ml |
| BSB-2375-1 | Concentrate | 1.0 ml |
| BSB 2375- CS | Control Slides | 5 |

SF-1/Steroidogenic Factor 1, MMab



IHC of SF-1/Steroidogenic Factor 1 on a FFPE Colon Adenocarcinoma Tissue

The Steroidogenic Factor 1 (SF-1) protein is a transcription factor involved in sex determination by controlling activity of genes related to the reproductive glands or gonads and adrenal glands. This protein is encoded by the NR5A1 gene. SF-1 expression is localized to adult steroidogenic tissues correlating with known expression profiles of steroid hydroxylases. Using in situ hybridization with SF-1 cRNA-specific probes detected gene transcripts in adrenocortical cells, Leydig cells, and ovarian theca and granulosa cells. SF-1 specific antibody studies confirmed the expression profile of SF-1 in rats and humans corresponding to sites of transcript detection.

SF-1 has been found to be a highly valuable IHC marker to determine the adrenocortical origin of an adrenal mass with high sensitivity and specificity. In addition, SF-1 expression is of stage-independent prognostic value in patients with adrenocortical carcinoma. Other SF-1 pathologies include adrenal failure (mutations in the SF-1 DNA-binding interface), adrenal or ovarian insufficiency and gonadal dysgenesis (heterozygous mutations), endometriosis (promoter hypomethylation), and male infertility (mutations in the hinge region of the protein). For the differential diagnosis with endometrioid tumors and carcinoid of the ovary, SF-1 is a sensitive and specific IHC marker for Sertoli cell tumor and that SF-1 is diagnostically comparable with other good sex cord-stromal markers.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-149


ISOTYPE: IgG2a

CONTROL: Breast, Fallopian Tube, Colon, Bone Marrow, Testis, Transitional Cell Carcinoma, Lung Adenocarcinoma, Papillary Thyroid Carcinoma, Prostate Adenocarcinoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

SMAD4/DPC4, MMab



SMAD4 / DPC4 on a FFPE Pancreatic Cancer

SMAD 4, also known as DPC4 or SMAD family member n°4, is a protein involved in cell signaling in mammals. SMAD 4 forms with SMAD 3 a complex which can bind to DNA and modify the expression of several genes related to cellular activities such as proliferation or differentiation. The abbreviation coin co-SMAD stands for common mediator. SMAD 4 is also defined as a signal transducer.

SMAD4, is often found mutated in many cancers. The mutation can be inherited or acquired during an individual's lifetime. If inherited, the mutation affects both somatic and sexual cells. The protein is present in skin, pancreatic, colon, uterus and epithelial cells. It is also produced by fibroblasts. The functional SMAD 4 participates in the regulation of the TGF-β signal transduction pathway, which negatively regulates growth of epithelial cells and the extracellular matrix (ECM). When the structure of SMAD 4 is altered, expression of the genes involved in cell growth is no longer regulated and cell proliferation can go on without any inhibition. The important number of cell divisions leads to the forming of tumors and then to multiploid colorectal cancer and pancreatic carcinoma. It is found inactivated in at least 50% of pancreatic cancers. SMAD 4 is also found mutated in the autosomal dominant disease juvenile polyposis syndrome (JPS). JPS is characterized by hamartomatous polyps in the gastrointestinal (GI) tract.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-63


ISOTYPE: IgG1

CONTROL: Pancreas, Thyroid, Placenta, Cervix, Transitional Cell Carcinoma, Lymphoblastic Lymphoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

SMAD4/DPC4, RMAb



IHC of SMAD4 / DPC4 on a FFPE Pancreatic Cancer

SMAD 4, also known as DPC4 or SMAD family member n°4, is a protein involved in cell signaling in mammals. SMAD 4 forms with SMAD 3 a complex which can bind to DNA and modify the expression of several genes related to cellular activities such as proliferation or differentiation. The abbreviation coin co-SMAD stands for common mediator. SMAD 4 is also defined as a signal transducer.

SMAD4, is often found mutated in many cancers. The mutation can be inherited or acquired during an individual's lifetime. If inherited, the mutation affects both somatic and sexual cells. The protein is present in skin, pancreatic, colon, uterus and epithelial cells. It is also produced by fibroblasts. The functional SMAD 4 participates in the regulation of the TGF-β signal transduction pathway, which negatively regulates growth of epithelial cells and the extracellular matrix (ECM). When the structure of SMAD 4 is altered, expression of the genes involved in cell growth is no longer regulated and cell proliferation can go on without any inhibition. The important number of cell divisions leads to the forming of tumors and then to multiploid colorectal cancer and pancreatic carcinoma. It is found inactivated in at least 50% of pancreatic cancers. SMAD 4 is also found mutated in the autosomal dominant disease juvenile polyposis syndrome (JPS). JPS is characterized by hamartomatous polyps in the gastrointestinal (GI) tract.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-SMAD4

ISOTYPE: IgG

CONTROL: Pancreas, Testis, Breast, Bone Marrow, Colon

LOCALIZATION: Cytoplasmic, Nuclear


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3746-3 | Tinto Predilute | 3.0 ml |
| BSB-3746-7 | Tinto Predilute | 7.0 ml |
| BSB-3746-15 | Tinto Predilute | 15.0 ml |
| BSB-3746-01 | Concentrate | 0.1 ml |
| BSB-3746-05 | Concentrate | 0.5 ml |
| BSB-3746-1 | Concentrate | 1.0 ml |
| BSB-3746-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3204 | Tinto Predilute | 3.0 ml |
| BSB 3205 | Tinto Predilute | 7.0 ml |
| BSB 3206 | Tinto Predilute | 15.0 ml |
| BSB 3207 | Concentrate | 0.1 ml |
| BSB 3208 | Concentrate | 0.5 ml |
| BSB 3209 | Concentrate | 1.0 ml |
| BSB 3210 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3399 | Tinto Predilute | 3.0 ml |
| BSB 3400 | Tinto Predilute | 7.0 ml |
| BSB 3401 | Tinto Predilute | 15.0 ml |
| BSB 3402 | Concentrate | 0.1 ml |
| BSB 3403 | Concentrate | 0.5 ml |
| BSB 3404 | Concentrate | 1.0 ml |
| BSB 3405 | Control Slides | 5 |

Smoothelin, MAb




IHC of Smoothelin on a FFPE Uterus Tissue

Smoothelin is a constituent of the smooth muscle cell cytoskeleton protein exclusively found in differentiated smooth muscle cells (SMC). Cells with SMC-like characteristics, such as myofibroblasts and myoepithelial cells, as well as skeletal and cardiac muscle do not contain smoothelin. To distinguish bladder muscularis mucosae (MM) from muscularis propria (MP) muscle bundles is crucial for accurate staging of bladder carcinoma.

Strong smoothelin expression is nearly exclusively observed in muscularis propria. Therefore, the staining pattern of MP (strongly positive) and MM (negative or weakly positive) makes this technique an attractive diagnostic tool for the sometimes difficult task of staging bladder urothelial carcinoma such as in transurethral resection specimens of urinary bladder tumors. Anti-Smoothelin can also be useful in differentiating between benign (+) and malignant smooth muscle tumors (-).

Somatostatin, MAb




IHC of Somatostatin on a FFPE Pancreas Tissue

Somatostatin is a peptide hormone that regulates the endocrine system and affects neurotransmission and cell proliferation via interaction with G-protein-coupled somatostatin receptors and by inhibition of the release of numerous secondary hormones. Somatostatin has two active forms produced by alternative cleavage of a single preproprotein: one of 14 amino acids; the other of 28 amino acids. Somatostatin is secreted not only by cells of the hypothalamus but also by the stomach, intestine, and delta cells of the pancreas. It binds to somatostatin receptors.

Somatostatin is a useful marker of D-cells of pancreatic islet cells. D-cells are used to identify hyperplasia of the pancreatic islets. Most of these tumors are malignant, giving rise to Somatostatinomas. Somatostatin suppresses gastric acid secretion, gallbladder contractions and pancreatic enzyme secretion.

Somatostatin, RMAb



IHC of Somatostatin on a FFPE Pancreas Tissue

Somatostatin is a peptide hormone that regulates the endocrine system and affects neurotransmission and cell proliferation via interaction with G-protein-coupled somatostatin receptors and by inhibition of the release of numerous secondary hormones. Somatostatin has two active forms produced by alternative cleavage of a single preproprotein: one of 14 amino acids; the other of 28 amino acids. Somatostatin is secreted not only by cells of the hypothalamus but also by the stomach, intestine, and delta cells of the pancreas. It binds to somatostatin receptors.

Somatostatin is a useful marker of D-cells of pancreatic islet cells. D-cells are used to identify hyperplasia of the pancreatic islets. Most of these tumors are malignant, giving rise to Somatostatinomas. Somatostatin suppresses gastric acid secretion, gallbladder contractions and pancreatic enzyme secretion.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: R4A

ISOTYPE: IgG1

CONTROL: Leiomyoma, Bladder, Prostate, Myometrium, Fallopian Tube, Colon, Trantision Cell Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-113

ISOTYPE: IgG2a/K

CONTROL: Pancreas, Colon, Adrenal, Brain

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP130

ISOTYPE: IgG

CONTROL: Pancreas, Brain, Pituitary

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2189 | Tinto Predilute | 3.0 ml |
| BSB 2190 | Tinto Predilute | 7.0 ml |
| BSB 2191 | Tinto Predilute | 15.0 ml |
| BSB 2192 | Concentrate | 0.1 ml |
| BSB 2193 | Concentrate | 0.5 ml |
| BSB 2194 | Concentrate | 1.0 ml |
| BSB 2195 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3413 | Tinto Predilute | 3.0 ml |
| BSB 3414 | Tinto Predilute | 7.0 ml |
| BSB 3415 | Tinto Predilute | 15.0 ml |
| BSB 3416 | Concentrate | 0.1 ml |
| BSB 3417 | Concentrate | 0.5 ml |
| BSB 3418 | Concentrate | 1.0 ml |
| BSB 3419 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2196 | Tinto Predilute | 3.0 ml |
| BSB 2197 | Tinto Predilute | 7.0 ml |
| BSB 2198 | Tinto Predilute | 15.0 ml |
| BSB 2199 | Concentrate | 0.1 ml |
| BSB 2200 | Concentrate | 0.5 ml |
| BSB 2201 | Concentrate | 1.0 ml |
| BSB 5937 | Control Slides | 5 |

Somatostatin Receptor 2/ SSTR2, RMaB




IHC of Somatostatin Receptor 2/SSTR2 on a FFPE Lung Neuroendocrine Carcinoma Tissue

SOX-2, RMaB



IHC of SOX-2 on a FFPE Brain Tissue

SOX-2, RMaB



IHC of SOX-2 on FFPE Brain Tissue

Somatostatin Receptor 2 (SSTR2) is one of five subtypes of the somatostatin receptors, which belong to the superfamily of G protein-coupled receptors (GPCRs). Somatostatin receptor 2 is encoded by the SSTR2 gene, located on chromosome 17q25.1. SSTR2 becomes activated by the hormone somatostatin (SST), which is an inhibitor of hormone secretion and gastrointestinal function. SST and its receptor subtypes also prevent angiogenesis and have anti-proliferative effects on healthy and cancerous cells.

Somatostatin Receptors have been reported to be highly expressed in a wide variety of human tumors. High SSTR2 expression was found via IHC in neuroblastomas, medulloblastomas, paragangliomas, small cell lung cancers, meningiomas, and breast cancers. One study reported a highly specific and increased expression of SSTR2 in a large series of neural and neuroendocrine tumors. Another study investigating patients with gastroenteropancreatic neuroendocrine neoplasm (GEP-NEN) demonstrated the correlation between decreased immunohistochemical staining and advanced stage of tumors. Additionally, an improved survival of patients with high SSTR2 expression was found, indicating the usefulness of SSTR2 as a potential prognostic marker for GEP-NEN.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP149
ISOTYPE: IgG
CONTROL: Placenta, Brain, Testis, Prostate, Papillary Thyroid Carcinoma, Transitional Cell Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human, Rat, Human

SRY (sex determining region Y)-box 2, also known as SOX2, is a transcription factor that is essential for maintaining self-renewal, or pluripotency, of undifferentiated embryonic stem cells. It is required for stem cell maintenance in the central nervous system, and it also regulates gene expression in the stomach.

SOX2 is expressed in fetal brain and is used as a marker for multipotential neural stem cells. In tumors, SOX2 expression is observed in teratoma of the central nervous system, melanoma, testicular germ cell tumor, cervical carcinoma, lung cancer, breast cancer with basal cell phenotype, and squamous cell carcinoma of the gastrointestinal tract. SOX2 may be useful in the identification of embryonal carcinoma. In stage I lung adenocarcinomas, SOX2 seems to be an independent predictor of poor outcome and may help stratify patients at increased risk for recurrence.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP103
ISOTYPE: IgG
CONTROL: Brain, Oligodendrogloma, Squamous Cell Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

SRY (sex determining region Y)-box 2, also known as SOX2, is a transcription factor that is essential for maintaining self-renewal, or pluripotency, of undifferentiated embryonic stem cells. It is required for stem cell maintenance in the central nervous system, and it also regulates gene expression in the stomach.

SOX2 is expressed in fetal brain and is used as a marker for multipotential neural stem cells. In tumors, SOX2 expression is observed in teratoma of the central nervous system, melanoma, testicular germ cell tumor, cervical carcinoma, lung cancer, breast cancer with basal cell phenotype, and squamous cell carcinoma of the gastrointestinal tract. SOX2 may be useful in the identification of embryonal carcinoma. In stage I lung adenocarcinomas, SOX2 seems to be an independent predictor of poor outcome and may help stratify patients at increased risk for recurrence.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM427
ISOTYPE: IgG
CONTROL: Brain, Oligodendrogloma, Squamous Cell Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Predicted: Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3748-3 | Tinto Predilute | 3.0 ml |
| BSB-3748-7 | Tinto Predilute | 7.0 ml |
| BSB-3748-15 | Tinto Predilute | 15.0 ml |
| BSB-3748-01 | Concentrate | 0.1 ml |
| BSB-3748-05 | Concentrate | 0.5 ml |
| BSB-3748-1 | Concentrate | 1.0 ml |
| BSB-3748-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2202 | Tinto Predilute | 3.0 ml |
| BSB 2203 | Tinto Predilute | 7.0 ml |
| BSB 2204 | Tinto Predilute | 15.0 ml |
| BSB 2205 | Concentrate | 0.1 ml |
| BSB 2206 | Concentrate | 0.5 ml |
| BSB 2207 | Concentrate | 1.0 ml |
| BSB 2208 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3770-3 | Tinto Predilute | 3.0 ml |
| BSB-3770-7 | Tinto Predilute | 7.0 ml |
| BSB-3770-15 | Tinto Predilute | 15.0 ml |
| BSB-3770-01 | Concentrate | 0.1 ml |
| BSB-3770-05 | Concentrate | 0.5 ml |
| BSB-3770-1 | Concentrate | 1.0 ml |
| BSB-3770-CS | Control Slides | 5 |

SOX-9, RMAb




IHC of SOX-9 on a FFPE Lymph Node Tissue

SOX-10, MMAb



IHC of Cytokeratin SOX-10 on a FFPE Melanoma Tissue

SOX10, RMAb



IHC of SOX-10 on a FFPE Melanoma Tissue

Transcription factor SOX-9 is a protein that in humans is encoded by the SOX9 gene. SOX9 acts during chondrocyte differentiation and regulates transcription of the anti-Müllerian hormone (AMH) gene. It is expressed during embryogenesis, in the cartilage, neural crest, kidney, and pancreas. SOX-9 plays a pivotal role in male sexual development, interacts with a few other genes to promote the development of male sexual organs and its activity is also required for development, differentiation, and lineage commitment in various tissues including the intestinal epithelium.

SOX9 exhibits several pro-oncogenic properties, including the ability to promote proliferation, inhibit senescence, and collaborate with other oncogenes in neoplastic transformation. In normal colorectal mucosa, SOX9 expression is found predominantly to the lower part of crypts, the proliferative compartment and putative site of stem cells, suggesting SOX9 as a putative stem or progenitor cell biomarker. Although staining is predominantly nuclear, cytoplasmic SOX9 may serve as a valuable prognostic marker for Invasive Ductal Carcinomas and Metastatic Breast Cancer. Additionally, SOX9 upregulation has been associated with higher tumor stage and grade, and overexpression has been recognized as an independent prognostic marker for decreased survival in Colorectal Cancer, NSCLC and HCC patients. In Pancreatic Cancer, SOX9 has been found to regulate the EGFR pathway throughout pancreatic tumorigenesis.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP317

ISOTYPE: IgG

CONTROL: Colon, Prostate, Skin, Breast, Tonsil, Lymph Node, Colon Carcinoma, SCC

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human predicted, predicted Mouse, predicted Rat

Transcription factor SOX-10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. This protein acts as a nucleocytoplasmic shuttle protein and is important for neural crest and peripheral nervous system development. Mutations in this gene are associated with Waardenburg-Shah and Waardenburg-Hirschsprung disease. Anti-SOX-10 has been recently shown to be a sensitive marker of melanoma, including conventional, spindled, and desmoplastic subtypes.

SOX-10 is expressed by metastatic melanomas and nodal capsular nevus in sentinel lymph nodes, but not by other lymph node components such as dendritic cells which usually express S100 protein. In scar specimens, immature fibroblasts, epithelioid granulomas, and histiocytic proliferations can histopathologically mimic residual melanoma and even be positive for MiTF and S100. However, SOX-10 is less likely to be expressed by fibroblasts or histiocytes, especially compared to MiTF and S100. Anti-SOX-10 produces a nuclear stain that provides a clean signal that is much sharper and darker in staining quality when compared to the use of antibodies against MiTF and S100.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-62

ISOTYPE: IgG2b/K

CONTROL: Skin, Salivary Gland, Breast, Melanoma, Schwannoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

Transcription factor SOX-10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. This protein acts as a nucleocytoplasmic shuttle protein and is important for neural crest and peripheral nervous system development. Mutations in this gene are associated with Waardenburg-Shah and Waardenburg-Hirschsprung disease. Anti-SOX-10 has been recently shown to be a sensitive marker of melanoma, including conventional, spindled, and desmoplastic subtypes.

SOX-10 is expressed by metastatic melanomas and nodal capsular nevus in sentinel lymph nodes, but not by other lymph node components such as dendritic cells which usually express S100 protein. In scar specimens, immature fibroblasts, epithelioid granulomas, and histiocytic proliferations can histopathologically mimic residual melanoma and even be positive for MiTF and S100. However, SOX-10 is less likely to be expressed by fibroblasts or histiocytes, especially compared to MiTF and S100. Anti-SOX-10 produces a nuclear stain that provides a clean signal that is much sharper and darker in staining quality when compared to the use of antibodies against MiTF and S100.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP268

ISOTYPE: IgG

CONTROL: Skin, Salivary Gland, Breast, Melanoma, Schwannoma

LOCALIZATION: Nuclear


SPECIES REACTIVITY: Human predicted, predicted Mouse, predicted Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3211 | Tinto Predilute | 3.0 ml |
| BSB 3212 | Tinto Predilute | 7.0 ml |
| BSB 3213 | Tinto Predilute | 15.0 ml |
| BSB 3214 | Concentrate | 0.1 ml |
| BSB 3215 | Concentrate | 0.5 ml |
| BSB 3216 | Concentrate | 1.0 ml |
| BSB 3217 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2209 | Tinto Predilute | 3.0 ml |
| BSB 2210 | Tinto Predilute | 7.0 ml |
| BSB 2211 | Tinto Predilute | 15.0 ml |
| BSB 2212 | Concentrate | 0.1 ml |
| BSB 2213 | Concentrate | 0.5 ml |
| BSB 2214 | Concentrate | 1.0 ml |
| BSB 2215 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2580 | Tinto Predilute | 3.0 ml |
| BSB 2581 | Tinto Predilute | 7.0 ml |
| BSB 2582 | Tinto Predilute | 15.0 ml |
| BSB 2583 | Concentrate | 0.1 ml |
| BSB 2584 | Concentrate | 0.5 ml |
| BSB 2585 | Concentrate | 1.0 ml |
| BSB 2586 | Control Slides | 5 |

SOX-11, MAb




IHC of SOX-11 on a FFPE Mantle Cell Lymphoma Tissue

SOX-11, MAb



IHC of SOX-11 on a FFPE Mantle Cell Lymphoma Tissue

Spectrin, MAb



IHC of Spectrin on a FFPE Bone Marrow Tissue

Transcription factor SOX-11 is a member of the group C SOX (SRY-related HMG-box) transcription factor family involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The protein may function in the developing nervous system and play a role in tumorigenesis and adult neurogenesis. SOX-11 is normally expressed in the developing human central nervous system, Medulloblastoma, and Glioma.

Anti-SOX-11 nuclear protein expression is highly associated with both Cyclin D1-positive and negative Mantle Cell Lymphomas, with a stronger and more homogeneous Immunohistochemistry staining than Cyclin D1.

Transcription factor SOX-11 is a member of the group C SOX (SRY-related HMG-box) transcription factor family involved in the regulation of embryonic development and in the determination of cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. The protein may function in the developing nervous system and plays a role in tumorigenesis and adult neurogenesis. SOX-11 is normally expressed in the developing human central nervous system, medulloblastoma, and glioma.

SOX-11 nuclear protein expression is highly associated with both cyclin D1- positive and negative mantle cell lymphomas, with a stronger and more homogeneous Immunostain than cyclin D1. Absence of SOX-11 expression in mantle cell lymphoma (MCL) may be a characteristic of indolent MCL. Sensitivity and specificity for combined SOX-11 and cyclin D1 immunohistostaining in the diagnosis of cutaneous mantle cell lymphoma are both 100%. SOX-11 has also been detected in rare cases of Burkitt lymphoma, lymphoblastic lymphoma and T cell prolymphocytic leukemia and overexpression can be seen in malignant gliomas.

Spectrin is a cytoskeletal protein that lines the intracellular side of the plasma membrane of many cell types in pentagonal or hexagonal arrangements, forming scaffolds and playing an important role in maintenance of plasma-membrane integrity and cytoskeletal structure. The hexagonal arrangements are formed by tetramers of spectrin associating with short actin filaments at either end of the tetramer. These short actin filaments act as junctional complexes, allowing the formation of the hexagonal mesh.

Spectrin is found in the intracellular side of the plasma membrane of many cell types found in muscles, red blood cells and red cell precursors. Anti-Spectrin antibody is useful in the diagnosis of Erythroid Leukemias.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: CL0142
ISOTYPE: IgG2a
CONTROL: Mantle Cell Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-167
ISOTYPE: IgG1
CONTROL: Mantle Cell Lymphoma, Lung Neuroendocrine Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: RBC2/3D5
ISOTYPE: IgG2b/K
CONTROL: Bone Marrow, Spleen
LOCALIZATION: Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2216 | Tinto Predilute | 3.0 ml |
| BSB 2217 | Tinto Predilute | 7.0 ml |
| BSB 2218 | Tinto Predilute | 15.0 ml |
| BSB 2219 | Concentrate | 0.1 ml |
| BSB 2220 | Concentrate | 0.5 ml |
| BSB 2221 | Concentrate | 1.0 ml |
| BSB 2222 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3783-3 | Tinto Predilute | 3.0 ml |
| BSB-3783-7 | Tinto Predilute | 7.0 ml |
| BSB-3783-15 | Tinto Predilute | 15.0 ml |
| BSB-3783-01 | Concentrate | 0.1 ml |
| BSB-3783-05 | Concentrate | 0.5 ml |
| BSB-3783-1 | Concentrate | 1.0 ml |
| BSB-3783-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5938 | Tinto Predilute | 3.0 ml |
| BSB 5939 | Tinto Predilute | 7.0 ml |
| BSB 5940 | Tinto Predilute | 15.0 ml |
| BSB 5941 | Concentrate | 0.1 ml |
| BSB 5942 | Concentrate | 0.5 ml |
| BSB 5943 | Concentrate | 1.0 ml |
| BSB 5944 | Control Slides | 5 |

STAR, RMAb

IHC of STAR on a FFPE Adrenal Tissue

STAT6, RMAb

IHC of STAT-6 on a FFPE Solitary Fibrous Tumor Tissue

Stathmin, RMAb

IHC of Stathmin on a FFPE Anal Carcinoma Tissue

Steroidogenic acute regulatory protein (STAR) is a protein that in humans is encoded by the STAR gene. The protein encoded by this gene plays a key role in the acute regulation of steroid hormone synthesis by enhancing the conversion of cholesterol into pregnenolone. This protein permits the cleavage of cholesterol into pregnenolone by mediating the transport of cholesterol from the outer mitochondrial membrane to the inner mitochondrial membrane.

STAR is primarily present in steroid-producing cells, including Leydig cells in the testis, theca cells and luteal cells in the ovary and adrenal cells in the adrenal cortex. Low level of STAR expression in other tissues that produce steroid hormones for local use have been reported. STAR is a sensitive and specific marker for Leydig cell tumor. It is useful for differential diagnosis of sex-cord stromal tumor (SCST). Mutations in this gene are a cause of congenital lipid adrenal hyperplasia (CLAH), also called lipid CAH.

STAT6 is a human gene. The protein encoded by this gene is a member of the STAT family of transcription factors. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. This protein plays a central role in exerting IL4 mediated biological responses. It is found to induce the expression of BCL2L1/BCL-X(L), which is responsible for the anti-apoptotic activity of IL4. STAT6 protein expression can be identified by IHC in the cytoplasm and nucleus of several tissues.

Recurrent somatic fusions of the NGFI-A-binding protein 2 (NAB2) gene and STAT6 gene, located at chromosomal region 12q13, have been identified in Solitary Fibrous Tumors (SFT). STAT6 is a highly sensitive and specific immunohistochemical marker for SFT and can be helpful to distinguish this tumor type from histologic mimics. STAT6 is amplified in a subset of dedifferentiated Liposarcoma, resulting in STAT6 protein expression that can be detected by immunohistochemistry and may be a potential pitfall in the differential diagnosis of dedifferentiated Liposarcoma and Solitary Fibrous Tumor. These findings suggest a role for STAT6-mediated transcriptional activity in some cases of dedifferentiated Liposarcoma and highlight the genomic complexity and heterogeneity of dedifferentiated Liposarcomas.

Stathmin 1/oncoprotein 18, also known as STMN1, is a highly conserved 17 kDa protein. Stathmin performs an important function in regulating rapid microtubule remodeling of the cytoskeleton in response to the cell's needs. Regulation of stathmin is cell cycle dependent and controlled by the cell's protein kinases in response to specific cell signals. Stathmin can cause uncontrolled cell proliferation when mutated and not functioning properly.

Overexpression of Stathmin has been associated with tumor progression in endometrial carcinomas, ovarian cancer and oral squamous-cell carcinoma. Stathmin has been found to be positive in 29% of CINs with differential expression based on the grade of the lesion as 9% being CIN1, 45% CIN2, and 93% CIN3; whereas, p16 staining of the same cases was positive in 80% of CINs with 71% CIN1, 100% CIN2, and 94% CIN3. Stathmin shows similar sensitivity for CIN3 to anti-p16 (93% vs 94%) although it drops off for CIN2 (73% vs 96%). The specificity of anti-stathmin for both CIN2-3 (94%) and CIN3 (89%) is higher than that of anti-p16 (44% and 39%, respectively). Therefore, Stathmin has major potential as a diagnostic marker in CIN classification over anti-p16 and it is valuable to distinguish CIN3 from the majority of low-grade precursors and negative/reactive cervical biopsies.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP226

ISOTYPE: IgG

CONTROL: Leydig Cells of Testis, Adrenal & Leydig Cell Tumors

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Mouse, Rat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP325

ISOTYPE: IgG

CONTROL: Solitary Fibrous Tumor

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP247

ISOTYPE: IgG

CONTROL: Testis, Tonsil, HSIL Cervical Carcinoma, Lymphoblastic Lymphoma, Bladder TCC

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human predicted, Mouse predicted, Rat predicted

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3630 | Tinto Predilute | 3.0 ml |
| BSB 3631 | Tinto Predilute | 7.0 ml |
| BSB 3632 | Tinto Predilute | 15.0 ml |
| BSB 3633 | Concentrate | 0.1 ml |
| BSB 3634 | Concentrate | 0.5 ml |
| BSB 3635 | Concentrate | 1.0 ml |
| BSB 3636 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3420 | Tinto Predilute | 3.0 ml |
| BSB 3421 | Tinto Predilute | 7.0 ml |
| BSB 3422 | Tinto Predilute | 15.0 ml |
| BSB 3423 | Concentrate | 0.1 ml |
| BSB 3424 | Concentrate | 0.5 ml |
| BSB 3425 | Concentrate | 1.0 ml |
| BSB 3426 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2594 | Tinto Predilute | 3.0 ml |
| BSB 2595 | Tinto Predilute | 7.0 ml |
| BSB 2596 | Tinto Predilute | 15.0 ml |
| BSB 2597 | Concentrate | 0.1 ml |
| BSB 2598 | Concentrate | 0.5 ml |
| BSB 2599 | Concentrate | 1.0 ml |
| BSB 2600 | Control Slides | 5 |

Surfactant protein D/SP-D, MMab



IHC of Surfactant protein D/SP-D on a FFPE Lung Squamous Cell Carcinoma Tissue

Surfactant protein D, also known as SFTPD or SP-D, is a protein encoded by the SFTPD gene. Pulmonary surfactants are essential to proper respiratory structure and function, consisting of ~90% lipids (mostly phospholipids) and 8–10% surfactant-associated proteins. SP-D is a pattern-recognition molecule in the collectin (collagen-containing C-type lectin) family. In the lungs, it is an opsonin that can agglutinate a range of microbes and enhance their clearance via phagocytosis and super-oxidative burst. It can interfere with allergen-IgE interaction and suppress basophil and mast cell activation.

Other surfactant proteins like SP-B and SP-C showed strong immunohistochemical expression in Lung Hyperplasias and Adenomas, suggesting that SP-B and SP-C are related to lung tumorigenesis. SP-D is likely an innate immune surveillance molecule against tumor development. SP-D may induce apoptosis in Pancreatic Adenocarcinoma via Fas-mediated pathway in a p53-independent manner. Studies Have found low expression of SP-D in Lung, Gastric, and Breast cancers and high expression in different stages and grades of Ovarian cancer.

The purified Spike protein of SARS-CoV-2 bound to Vero but not 293T cells and was itself recognized by SP-D, in the lung alveoli. It suggests that SARS-CoV interacts with innate immune mechanisms in the lung through its S-protein and regulates pulmonary inflammation.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-162


ISOTYPE: IgG1

CONTROL: Placenta, Lung, Pancreas, Adrenal Gland, Lung Squamous Cell Carcinoma, Lung Adenocarcinoma, Gastric GIST

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

Survivin, RMab



IHC of Survivin on a FFPE Colon Tissue

Survivin, also called baculoviral inhibitor of apoptosis repeat-containing 5 or BIRC5, is a member of the inhibitor of apoptosis (IAP) family. The survivin protein functions to inhibit caspase activation, thereby leading to negative regulation of apoptosis or programmed cell death. The survivin protein is expressed highly in most human tumors and fetal tissue, but is completely absent in terminally differentiated cells. Survivin expression is also highly regulated by the cell cycle and is only expressed in the G2-M phase. It is known that survivin localizes to the mitotic spindle by interaction with tubulin during mitosis and may play a contributing role in regulating mitosis.

The association of survivin expression with tumor progression, but not overall patient survival, has been observed in a variety of malignancies including renal cell carcinoma, ovary carcinoma, hepatocellular carcinoma, prostate carcinoma and breast carcinoma.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP119


ISOTYPE: IgG

CONTROL: Colon, Placenta, Testis, Tonsil, Bone Marrow, Colon Carcinoma, Lymphoblastic Lymphoma, Bladder TCC

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

SV40, MMab



IHC and IF of SV40 on a FFPE Infected Kidney Tissue

SV40 is an abbreviation for Simian vacuolating virus 40 or Simian virus 40, a polyomavirus that is found in both monkeys and humans. Like other polyomaviruses, SV40 is a DNA virus that has the potential to cause tumors, but most often persists as a latent infection.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: Pab101

ISOTYPE: IgG2a

CONTROL: SV40 Infected Tissue

LOCALIZATION: Nuclear


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3747-3 | Tinto Predilute | 3.0 ml |
| BSB-3747-7 | Tinto Predilute | 7.0 ml |
| BSB-3747-15 | Tinto Predilute | 15.0 ml |
| BSB-3747-01 | Concentrate | 0.1 ml |
| BSB-3747-05 | Concentrate | 0.5 ml |
| BSB-3747-1 | Concentrate | 1.0 ml |
| BSB-3747-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2223 | Tinto Predilute | 3.0 ml |
| BSB 2224 | Tinto Predilute | 7.0 ml |
| BSB 2225 | Tinto Predilute | 15.0 ml |
| BSB 2226 | Concentrate | 0.1 ml |
| BSB 2227 | Concentrate | 0.5 ml |
| BSB 2228 | Concentrate | 1.0 ml |
| BSB 2229 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2230 | Tinto Predilute | 3.0 ml |
| BSB 2231 | Tinto Predilute | 7.0 ml |
| BSB 2232 | Tinto Predilute | 15.0 ml |
| BSB 2233 | Concentrate | 0.1 ml |
| BSB 2234 | Concentrate | 0.5 ml |
| BSB 2235 | Concentrate | 1.0 ml |
| BSB 2236 | Control Slides | 5 |

Synaptophysin, RPAb




IHC of Synaptophysin on a FFPE Neuroendocrine Tumor

Synaptophysin is a synaptic vesicle glycoprotein weighing 38 kDa. It is present in endocrine cells, the brain, spinal cord, and adrenal glands. It acts as a marker for neuroendocrine cells.

Synaptophysin reacts with neuroendocrine cells of human adrenal medulla, carotid body, skin, pituitary, thyroid, lung, pancreas and gastrointestinal mucosa. Positive staining is seen in neurons of the brain, spinal cord, retina, and Paneth's cells in the gastrointestinal tract and gastric parietal cells. This antibody identifies normal neuroendocrine cells and neuroendocrine neoplasms. Diffuse, finely-granular cytoplasmic staining is observed and probably correlates with the distribution of the antigen within neurosecretory vesicles. The expression of Synaptophysin is independent of the presence of NSE or other neuroendocrine markers. Synaptophysin is an independent broad-range marker of neural and neuroendocrine differentiation.

Synaptophysin, RMAb




IHC of Synaptophysin on a FFPE Pancreas Tissue

Synaptophysin is a synaptic vesicle glycoprotein weighing 38 kDa. It is present in endocrine cells, the brain, spinal cord, and adrenal glands. It acts as a marker for neuroendocrine cells.

Synaptophysin reacts with neuroendocrine cells of human adrenal medulla, carotid body, skin, pituitary, thyroid, lung, pancreas and gastrointestinal mucosa. Positive staining is seen in neurons of the brain, spinal cord, retina, and Paneth's cells in the gastrointestinal tract and gastric parietal cells. This antibody identifies normal neuroendocrine cells and neuroendocrine neoplasms. Diffuse, finely-granular cytoplasmic staining is observed and probably correlates with the distribution of the antigen within neurosecretory vesicles. The expression of Synaptophysin is independent of the presence of NSE or other neuroendocrine markers. Synaptophysin is an independent broad-range marker of neural and neuroendocrine differentiation.

TAG-72, MMAb



IHC of TAG-72 on a FFPE Breast Tissue

Tumor-associated glycoprotein (TAG-72) has been shown to be expressed in a wide variety of epithelial malignant tissues. TAG-72 antigen is a high molecular glycoprotein found on the surface of many cancer cells, including breast, colon and pancreatic cells. It is present in human Adenocarcinomas and in lesser amounts, non-neoplastic tissues. The majority of human Adenocarcinomas including Colorectal, Pancreatic, Gastric, Ovarian, Endometrial, Mammary, and Non-Small Cell Lung Cancer display some cell populations that are positive for TAG-72.

TAG-72 has also been found to be useful for the distinction between Mesothelioma and Adenocarcinoma; however, false positive reactions can occur so results must be interpreted with the utmost caution.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Pancreas, Brain, Pituitary, Adrenal, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP158

ISOTYPE: IgG

CONTROL: Pancreas, Brain, Pituitary, Adrenal, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Rat, Donkey

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-21

ISOTYPE: IgG1/K

CONTROL: Breast Carcinoma

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5945 | Tinto Predilute | 3.0 ml |
| BSB 5946 | Tinto Predilute | 7.0 ml |
| BSB 5947 | Tinto Predilute | 15.0 ml |
| BSB 5948 | Concentrate | 0.1 ml |
| BSB 5949 | Concentrate | 0.5 ml |
| BSB 5950 | Concentrate | 1.0 ml |
| BSB 5951 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2237 | Tinto Predilute | 3.0 ml |
| BSB 2238 | Tinto Predilute | 7.0 ml |
| BSB 2239 | Tinto Predilute | 15.0 ml |
| BSB 2240 | Concentrate | 0.1 ml |
| BSB 2241 | Concentrate | 0.5 ml |
| BSB 2242 | Concentrate | 1.0 ml |
| BSB 2243 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5952 | Tinto Predilute | 3.0 ml |
| BSB 5953 | Tinto Predilute | 7.0 ml |
| BSB 5954 | Tinto Predilute | 15.0 ml |
| BSB 5955 | Concentrate | 0.1 ml |
| BSB 5956 | Concentrate | 0.5 ml |
| BSB 5957 | Concentrate | 1.0 ml |
| BSB 5958 | Control Slides | 5 |

Tau, MMab




IHC of Tau on a FFPE Astrocytoma Tissue

T-Bet/TBX-2, RMab



IHC of T-bet on a FFPE Hairy Cell Leukemia Tissue

TCL1, RMab



IHC of TCL1 on a FFPE Lymphoma Tissue

The tau proteins are the product of alternative splicing from a single gene that in humans is designated MAPT (microtubule-associated protein tau) and is located on chromosome 17. In humans, these proteins are found mostly in neurons compared to non-neuronal cells. One of tau's main functions is to modulate the stability of axonal microtubules. Tau proteins interact with tubulin to stabilize microtubules and promote tubulin assembly into microtubules. Through its isoforms and phosphorylation tau protein interacts with tubulin to stabilize microtubule assembly.

Hyperphosphorylation of the tau protein (tau inclusions, pTau) can result in the self-assembly of tangles of paired helical filaments and straight filaments, which are involved in the pathogenesis of Alzheimer's disease, frontotemporal dementia, and other tauopathies. When misfolded, this otherwise very soluble protein can form extremely insoluble aggregates that contribute to a number of neurodegenerative diseases. Mutations that alter function and isoform expression of tau lead to hyperphosphorylation, which in turn disassembles microtubules and sequesters normal tau, MAP 1, MAP 2, and ubiquitin into neurofibrillary tangles, which are composed of paired helical filaments (PHF). These insoluble structures damage cytoplasmic functions and interferes with axonal transport, which can lead to cell death.

T-box transcription factor TBX21, also known as T-bet, is a T-box transcription factor, is expressed in CD4+ T-lymphocytes committed to T-helper (Th)1 T-cell development from naïve T-helper precursor cells (Thp) and redirects Th2 T cells to Th1 development.

T-bet is expressed in CD4+ T lymphocytes in normal tissues. In lymphoid malignancies, TBX21 has been found in a subset of T-cell lymphomas with Th1 T cell differentiation, a subset of B-cell or Tcells, non-Hodgkin's lymphomas, majority of Hodgkin's lymphomas and precursor B-cell lymphoblastic leukemia/lymphoblastic lymphomas. However, B-cell neoplasms derived from pregerminal center or germinal center B-cells, including mantle cell lymphoma, follicular lymphoma, diffuse large B-cell lymphoma, and Burkitt lymphoma are negative for T-bet. Therefore, anti-T-bet should serve as a useful marker for the diagnosis and subtyping of B-cell and T-cell lymphoproliferative disorders. T-bet is a useful marker for Hodgkin's lymphoma and also helpful in identification of hairy cell leukemia.

T-cell leukemia/lymphoma protein 1 (TCL1, TCL1A, p14TCL1) is involved in T-cell prolymphocytic leukemia (T-PLL) and is normally found in the nucleus and cytoplasm of lymphoid lineage cells during early embryogenesis. Chromosomal translocations may lead to overexpression of TCL1, resulting in T-cell leukemia and B-cell lymphoma. TCL1 is expressed in more differentiated B-cells, under both reactive and neoplastic conditions, from antigen committed B-cells and in germinal center B-cells. TCL1 is down-regulated in the latest stage of B-cell differentiation.

TCL1 is overexpressed in Burkitt Lymphoma, the majority of AIDS-related non-Hodgkin's Lymphoma-designated Immunoblastic Plasmacytoid Lymphoma, Lymphoblastic Lymphoma, Chronic Lymphocytic Leukemia, Mantle Cell Lymphoma, Follicular Lymphoma, Diffuse Large B-cell Lymphoma, and Primary Cutaneous B-cell Lymphoma. Therefore, the most useful application of anti-TCL1 is the discrimination of B-cell lymphomas from T-cell lymphomas, CD30+ Anaplastic Large Cell Lymphomas, Multiple Myeloma, and Marginal Zone B-cell Lymphoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-115
ISOTYPE: IgG1/K
CONTROL: Brain, Kidney, Pituitary, Pancreas, Cervix, Skin, Salivary Gland, Astrocytoma
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat, Dog

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP263
ISOTYPE: IgG
CONTROL: Spleen, Tonsil, Cervix, Liver, Breast, Hairy Cell Leukemia & Histiocytic, Lymphoblastic Lymphoma Bladder TCC
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP105
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3427 | Tinto Predilute | 3.0 ml |
| BSB 3428 | Tinto Predilute | 7.0 ml |
| BSB 3429 | Tinto Predilute | 15.0 ml |
| BSB 3430 | Concentrate | 0.1 ml |
| BSB 3431 | Concentrate | 0.5 ml |
| BSB 3432 | Concentrate | 1.0 ml |
| BSB 3433 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2601 | Tinto Predilute | 3.0 ml |
| BSB 2602 | Tinto Predilute | 7.0 ml |
| BSB 2603 | Tinto Predilute | 15.0 ml |
| BSB 2604 | Concentrate | 0.1 ml |
| BSB 2605 | Concentrate | 0.5 ml |
| BSB 2606 | Concentrate | 1.0 ml |
| BSB 2607 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2251 | Tinto Predilute | 3.0 ml |
| BSB 2252 | Tinto Predilute | 7.0 ml |
| BSB 2253 | Tinto Predilute | 15.0 ml |
| BSB 2254 | Concentrate | 0.1 ml |
| BSB 2255 | Concentrate | 0.5 ml |
| BSB 2256 | Concentrate | 1.0 ml |
| BSB 2257 | Control Slides | 5 |

TCR Alpha, MAb




IHC of TCR Alpha on a FFPE Lung Adenocarcinoma Tissue

TCR Beta, MAb



IHC of TCR Beta on a FFPE Thymus Tissue

TCR Delta, MAb



IHC of TCR Delta on a FFPE Thymus Tissue

T-cell receptor alpha locus is a protein that in humans is encoded by the TRA gene, also known as TCRA or TRAA; It contributes the alpha chain to the larger TCR protein (T-cell receptor). The T cell receptor or TCR is a molecule found on the surface of T lymphocytes (or T cells) that is responsible for recognizing antigens bound to major histocompatibility complex (MHC) molecules. The TCR is composed of two different protein chains (that is, it is a heterodimer). In 95% of T cells, this consists of an alpha (α) and beta (β) chain, whereas in 5% of T cells this consists of gamma and delta (γ/δ) chains. This ratio changes during ontogeny and in cancers, such as Leukemia.

The T cell receptor or TCR is a molecule found on the surface of T lymphocytes (or T cells) that is responsible for recognizing antigens bound to major histocompatibility complex (MHC) molecules. The TCR is composed of two different protein chains (that is, it is a heterodimer). In 95% of T cells, this consists of an alpha (α) and beta (β) chain, whereas in 5% of T cells this consists of gamma and delta (γ/δ) chains. TCR Beta is a member of the immunoglobulin super family and a component of the CD3/TCR complex (along with TCR Alpha).

TCR Beta is expressed by thymocytes and a majority of peripheral (α - β TCR-bearing) T-cells. TCR recognition of self-peptides has been linked to autoimmune disease. Mutant self-peptides have been associated with tumors.

T cell receptor delta locus (symbol TRD), also known as TCRD, is a protein that in humans is encoded by the TRD gene. It contributes the delta (δ) chain to the larger TCR protein (T-cell receptor). The T cell receptor or TCR is a molecule found on the surface of T lymphocytes (or T cells) that is responsible for recognizing antigens bound to major histocompatibility complex (MHC) molecules.

Deletions and mutations of the TRG and TRD gene have been implicated in a variety of cancers. Specifically, $\gamma\delta$ T cells may contribute to the immune response against several tumor types (lymphoma, myeloma, breast, colon, lung, ovary, and others). They act directly through mediation of cytotoxic activity and indirectly through the regulation of other cell types responsible for the anti-tumor response. The presence of $\gamma\delta$ T cells in the tumor microenvironment has been associated with poor prognosis in some cancers. While $\gamma\delta$ T cells have been implicated in T cell lymphomas, there is also a specific subtype known as $\gamma\delta$ T-cell lymphoma, characterized by the proliferation of those cells exclusively. This lymphoma can be quite aggressive with ulcerative plaques and subcutaneous nodules. Apart from carcinomas, TRG has also been correlated with hepatitis B virus (HBV). Specifically, $V\delta 2+$ T cell levels and TCR $\gamma\delta$ T cell cytotoxicity were significantly lower was in patients with chronic HBV infections.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-126

ISOTYPE: IgG2A/k

CONTROL: Thymus, Tonsil, Lymph Node, Spleen, Hodgkin's Lymphoma & BDCM & CCRF-CEM Cell Lines

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-117

ISOTYPE: IgG1/K

CONTROL: Tonsil, Lymph Node

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-127

ISOTYPE: IgG1/K

CONTROL: Thymus, Tonsil, Lymph Node, Spleen

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3637 | Tinto Predilute | 3.0 ml |
| BSB 3638 | Tinto Predilute | 7.0 ml |
| BSB 3639 | Tinto Predilute | 15.0 ml |
| BSB 3640 | Concentrate | 0.1 ml |
| BSB 3641 | Concentrate | 0.5 ml |
| BSB 3642 | Concentrate | 1.0 ml |
| BSB 3643 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2258 | Tinto Predilute | 3.0 ml |
| BSB 2259 | Tinto Predilute | 7.0 ml |
| BSB 2260 | Tinto Predilute | 15.0 ml |
| BSB 2261 | Concentrate | 0.1 ml |
| BSB 2262 | Concentrate | 0.5 ml |
| BSB 2263 | Concentrate | 1.0 ml |
| BSB 2264 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3644 | Tinto Predilute | 3.0 ml |
| BSB 3645 | Tinto Predilute | 7.0 ml |
| BSB 3646 | Tinto Predilute | 15.0 ml |
| BSB 3647 | Concentrate | 0.1 ml |
| BSB 3648 | Concentrate | 0.5 ml |
| BSB 3649 | Concentrate | 1.0 ml |
| BSB 3650 | Control Slides | 5 |

TDP-43/TARDBP, MAb




IHC of TDP-43/TARDBP on a FFPE Lung Adenocarcinoma Tissue

TdT, RPaB



IHC of TdT on a FFPE Thymus Tissue

TdT, RMaB



IHC of TdT on a FFPE Thymus Tissue

Transactive Response DNA binding protein 43 (TDP-43 or TARDBP) is an RNA and DNA binding protein in the heterogeneous nuclear ribonucleoprotein (hnRNP) family, where it regulates RNA associated with glucose and lipid metabolism.

TDP-43/TARDBP is associated with several neurodegenerative conditions amyotrophic lateral sclerosis (ALS) and frontotemporal lobar degeneration (FTLD) and metabolic regulation in cancer cells. TDP-43/TARDBP regulation of miRNAs has both tumor supporting and suppressive roles. Upregulation of TDP-43/TARDBP can lead to continued autophagy for survival and suppression of apoptosis pathways, promoting cancer cell survival under conditions of nutrient stress, as was seen in a study of glioblastomas. TDP-43/TARDBP and its targeted miRNAs have been associated with increased cell growth and migration in hepatocellular and lung cancers, but also with tumor inhibition in rhabdomyosarcomas through miRNA miR-500a-3p. TDP-43/TARDBP may also influence the function and targets of MiR-152, and miRNA associated with hepatocellular carcinoma, endometrial, gastric, and ovarian cancers. TDP-43/TARDBP has been found to be a novel oncogene in melanoma and regulates melanoma proliferation and metastasis potentially through modulation of glucose metabolism. TDP43, has been found to promote triple-negative breast cancer (TNBC) progression.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-166

ISOTYPE: IgG2a

CONTROL: Breast, Fallopian Tube, Testis, Skin, Transitional Cell Carcinoma, Glioblastoma

LOCALIZATION: Nuclear, Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Rat, reactive with Horse, Dog, Bovine

Terminal Deoxynucleotidyl Transferase (also known as TdT and terminal transferase) is a specialized DNA polymerase expressed in immature, pre-B, pre-T lymphoid cells and acute Lymphoblastic Leukemia/Lymphoma cells. TdT catalyzes the addition of nucleotides to the 3' terminus of a DNA molecule. Unlike most DNA polymerases, it does not require a template. The preferred substrate of this enzyme is a protruding 3' overhang, but it can also add nucleotides to blunt or recessed 3' ends.

TdT is normally found in cortical thymocytes and primitive lymphocytes. TdT antibody detects its antigen found in the nucleus of normal hematopoietic cells, normal cortical thymocytes and in the cytoplasm of megakaryocytes of the bone marrow. TdT expression is seen in over 90% of Acute Lymphocytic Leukemia cases with the exception of pre-B-Cell ALL, and normal mature T- or B-lymphocytes. TdT is positive for approximately one third of all cases of Chronic Myeloid Leukemia, making it a good indicator of better response to chemotherapy.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Thymus, Lymphoblastic Lymphoma

LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

Terminal Deoxynucleotidyl Transferase (also known as TdT and terminal transferase) is a specialized DNA polymerase expressed in immature, pre-B, pre-T lymphoid cells and acute Lymphoblastic Leukemia/ Lymphoma cells. TdT catalyzes the addition of nucleotides to the 3' terminus of a DNA molecule. Unlike most DNA polymerases, it does not require a template. The preferred substrate of this enzyme is a protruding 3' overhang, but it can also add nucleotides to blunt or recessed 3' ends.

TdT is normally found in cortical thymocytes and primitive lymphocytes. TdT antibody detects its antigen found in the nucleus of normal hematopoietic cells, normal cortical thymocytes and in the cytoplasm of megakaryocytes of the bone marrow. TdT expression is seen in over 90% of Acute Lymphocytic Leukemia cases with the exception of pre-B-Cell ALL, and normal mature T- or B-lymphocytes. TdT is positive for approximately one third of all cases of Chronic Myeloid Leukemia, making it a good indicator of better response to chemotherapy.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-TdT

ISOTYPE: IgG

CONTROL: Thymus, Lymphoblastic Lymphoma


LOCALIZATION: Nuclear

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3776-3 | Tinto Predilute | 3.0 ml |
| BSB-3776-7 | Tinto Predilute | 7.0 ml |
| BSB-3776-15 | Tinto Predilute | 15.0 ml |
| BSB-3776-01 | Concentrate | 0.1 ml |
| BSB-3776-05 | Concentrate | 0.5 ml |
| BSB-3776-1 | Concentrate | 1.0 ml |
| BSB-3776-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5966 | Tinto Predilute | 3.0 ml |
| BSB 5967 | Tinto Predilute | 7.0 ml |
| BSB 5968 | Tinto Predilute | 15.0 ml |
| BSB 5969 | Concentrate | 0.1 ml |
| BSB 5970 | Concentrate | 0.5 ml |
| BSB 5971 | Concentrate | 1.0 ml |
| BSB 5972 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2265 | Tinto Predilute | 3.0 ml |
| BSB 2266 | Tinto Predilute | 7.0 ml |
| BSB 2267 | Tinto Predilute | 15.0 ml |
| BSB 2268 | Concentrate | 0.1 ml |
| BSB 2269 | Concentrate | 0.5 ml |
| BSB 2270 | Concentrate | 1.0 ml |
| BSB 2271 | Control Slides | 5 |

TdT, RMAb


IHC of TdT on a FFPE Lymphoblastic Lymphoma Tissue

Terminal Deoxynucleotidyl Transferase (also known as TdT and terminal transferase) is a specialized DNA polymerase expressed in immature, pre-B, pre-T lymphoid cells and acute Lymphoblastic Leukemia/Lymphoma cells. TdT catalyzes the addition of nucleotides to the 3' terminus of a DNA molecule. Unlike most DNA polymerases, it does not require a template. The preferred substrate of this enzyme is a protruding 3' overhang, but it can also add nucleotides to blunt or recessed 3' ends.

TdT is normally found in cortical thymocytes and primitive lymphocytes. TdT antibody detects its antigen found in the nucleus of normal hematopoietic cells, normal cortical thymocytes and in the cytoplasm of megakaryocytes of the bone marrow. TdT expression is seen in over 90% of Acute Lymphocytic Leukemia cases with the exception of pre-B-Cell ALL, and normal mature T- or B-lymphocytes. TdT is positive for approximately one third of all cases of Chronic Myeloid Leukemia, making it a good indicator of better response to chemotherapy.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP266
ISOTYPE: IgG
CONTROL: Thymus, Lymphoblastic Lymphoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2608 | Tinto Predilute | 3.0 ml |
| BSB 2609 | Tinto Predilute | 7.0 ml |
| BSB 2610 | Tinto Predilute | 15.0 ml |
| BSB 2611 | Concentrate | 0.1 ml |
| BSB 2612 | Concentrate | 0.5 ml |
| BSB 2613 | Concentrate | 1.0 ml |
| BSB 2614 | Control Slides | 5 |


TFE3, RMAb

IHC of TFE3 on a FFPE Testicular Cancer Tissue

Transcription factor E3 is a protein that in humans is encoded by the TFE3 gene. TFE3, a member of the helix-loop-helix family of transcription factors, binds to the mu-E3 motif of the immunoglobulin heavy-chain enhancer and is expressed in many cell types. A proportion of renal cell carcinomas (RCC) that occur in young patients are associated with translocations involving the TFE3 gene, which results in gene fusions. Subsets of papillary renal cell carcinomas, a t(X;1)(p11;q21) chromosome translocation has been repeatedly reported and is thought to be the cause of this cancer. As a result of the translocation, the transcription factor TFE3 on the X chromosome becomes fused to this gene on chromosome 1. The fused gene results in the fusion of N-terminal proline-rich region of the protein encoded by this gene to the entire TFE3 protein.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP285
ISOTYPE: IgG
CONTROL: Testis, Adrenal, Kidney, Testicular Cancer, RCC with Xp11.2 translocation, Alveolar Soft Part Sarcoma & Soft Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3225 | Tinto Predilute | 3.0 ml |
| BSB 3226 | Tinto Predilute | 7.0 ml |
| BSB 3227 | Tinto Predilute | 15.0 ml |
| BSB 3228 | Concentrate | 0.1 ml |
| BSB 3229 | Concentrate | 0.5 ml |
| BSB 3230 | Concentrate | 1.0 ml |
| BSB 3231 | Control Slides | 5 |

Thrombomodulin/CD141, RMAb

IHC of Thrombomodulin on a FFPE Mesothelioma Tissue


Thrombomodulin, also known as CD141, is an endothelial-specific type 1 membrane receptor that binds thrombin, resulting in the activation of protein C. This causes the degradation of clotting factors Va and VIIIa and reduces the amount of thrombin generated. Defect in Thrombomodulin is a cause of thromboembolic disease, also known as inherited thrombophilia.

Thrombomodulin was initially identified in endothelial cells, but is also found in extra-vascular sites, such as syncytiotrophoblasts in the placenta, epithelial tissues in the gingiva, in skin and in the synovial lining cells. In tumors, Thrombomodulin is expressed in vascular tumors and squamous cell carcinoma in a variety of tissues, including oral mucosa, esophagus, and skin. Thrombomodulin is a useful marker for detecting angiosarcoma, and can also be used to distinguish between mesothelioma (positive) from lung adenocarcinoma (negative).

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP175
ISOTYPE: IgG
CONTROL: Placenta, Liver, Kidney, Tonsil, Cervix, Bladder, Mesothelioma, Transitional Cell Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2279 | Tinto Predilute | 3.0 ml |
| BSB 2280 | Tinto Predilute | 7.0 ml |
| BSB 2281 | Tinto Predilute | 15.0 ml |
| BSB 2282 | Concentrate | 0.1 ml |
| BSB 2283 | Concentrate | 0.5 ml |
| BSB 2284 | Concentrate | 1.0 ml |
| BSB 2285 | Control Slides | 5 |

Thymidylate synthase/TS, MAb



IHC of Thymidylate synthase/TS on a FFPE Colon Tissue

Thymidylate synthase (TS) catalyzes the methylation of deoxyuridine monophosphate to deoxythymidine monophosphate using 5, 10-methylenetetrahydrofolate as a cofactor. This enzyme is critical for DNA repair and replication. TS inhibitors like 5-fluorouracil have been successful in down-regulating tumor progression and increasing immune responses in pancreatic, breast, gastric, ovarian, lung, and head & neck cancer.

A study reported TS protein expression (PE) and gene copy number (GCN) were assayed using IHC and silver in situ hybridization (SISH), respectively, on primary tumors of resected non-small cell lung (NSCLC) patients and concluded that TS PE and GCN vary widely in NSCLC and correlate significantly to each other. TS GCN is higher in SCCs, whereas TS PE does not associate with histological subtypes, clinical features or survival. Variability of TS PE and GCN may indicate potential benefit from pemetrexed therapy in selected SCC patients. A number of studies have investigated the relationship between thymidylate synthase (TS) expression and survival in colorectal cancer (CRC) patients. Most have reported poorer overall and progression-free survival with high TS expression. Another study suggests that immunostaining for TS and p53 protein is useful for pretreatment selection of gastric cancer patients unresponsive to S-1/cisplatin chemotherapy.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-160


ISOTYPE: IgG2a

CONTROL: Bone Marrow, Colon, Tonsil, Testis, T Cell Lymphoblastic Lymphoma, Ductal Breast Carcinoma

LOCALIZATION: Nuclear, Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Mouse, Rat

Thyroglobulin, MAb



IHC of Thyroglobulin on a FFPE Thyroid Tissue

Thyroglobulin (Tg) is a 660 kDa, dimeric protein produced by and used entirely within the thyroid gland. Tg is used by the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3). The active form of thyroxine, triiodothyronine, is produced both within the thyroid gland and on the periphery by 5'-deiodinase, which has been referred to as Tetraiodothyronine-5-deiodinase.

This antibody reacts with human thyroglobulin as demonstrated by a single band of immunoblotting in a lysate of human thyroid tissue. The vast majority of follicular carcinomas of the thyroid will give positive immunoreactivity for thyroglobulin, sometimes only focally. Poorly-differentiated Carcinomas of the Thyroid are frequently thyroglobulin negative. Adenocarcinomas of non-thyroid origin do not react with this antibody.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-49


ISOTYPE: IgG1

CONTROL: Thyroid, Thyroid Carcinoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Thyroglobulin, MAb



IHC of Thyroglobulin on a FFPE Papillary Thyroid Carcinoma Tissue

Thyroglobulin (Tg) is a 660 kDa, dimeric protein produced by and used entirely within the thyroid gland. Tg is used by the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3). The active form of thyroxine, triiodothyronine, is produced both within the thyroid gland and on the periphery by 5'-deiodinase, which has been referred to as Tetraiodothyronine-5-deiodinase.

This antibody reacts with human thyroglobulin as demonstrated by a single band of immunoblotting in a lysate of human thyroid tissue. The vast majority of follicular carcinomas of the thyroid will give positive immunoreactivity for thyroglobulin, sometimes only focally. Poorly-differentiated Carcinomas of the Thyroid are frequently thyroglobulin negative. Adenocarcinomas of non-thyroid origin do not react with this antibody.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 2H11/6E1

ISOTYPE: IgG1&IgG1

CONTROL: Thyroid, Thyroid Carcinoma

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3750-3 | Tinto Predilute | 3.0 ml |
| BSB-3750-7 | Tinto Predilute | 7.0 ml |
| BSB-3750-15 | Tinto Predilute | 15.0 ml |
| BSB-3750-01 | Concentrate | 0.1 ml |
| BSB-3750-05 | Concentrate | 0.5 ml |
| BSB-3750-1 | Concentrate | 1.0 ml |
| BSB-3750-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5973 | Tinto Predilute | 3.0 ml |
| BSB 5974 | Tinto Predilute | 7.0 ml |
| BSB 5975 | Tinto Predilute | 15.0 ml |
| BSB 5976 | Concentrate | 0.1 ml |
| BSB 5977 | Concentrate | 0.5 ml |
| BSB 5978 | Concentrate | 1.0 ml |
| BSB 5979 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2761 | Tinto Predilute | 3.0 ml |
| BSB 2762 | Tinto Predilute | 7.0 ml |
| BSB 2763 | Tinto Predilute | 15.0 ml |
| BSB 2764 | Concentrate | 0.1 ml |
| BSB 2765 | Concentrate | 0.5 ml |
| BSB 2766 | Concentrate | 1.0 ml |
| BSB 2767 | Control Slides | 5 |

Thyroglobulin, RMAb




IHC of Thyroglobulin on a FFPE Papillary Thyroid Carcinoma Tissue

Thyroglobulin (Tg) is a 660 kDa, dimeric protein produced by and used entirely within the thyroid gland. Tg is used by the thyroid gland to produce the thyroid hormones thyroxine (T4) and triiodothyronine (T3). The active form of thyroxine, triiodothyronine, is produced both within the thyroid gland and on the periphery by 5'-deiodinase, which has been referred to as Tetraiodothyronine-5-deiodinase.

This antibody reacts with human thyroglobulin as demonstrated by a single band of immunoblotting in a lysate of human thyroid tissue. The vast majority of follicular carcinomas of the thyroid will give positive immunoreactivity for thyroglobulin, sometimes only focally. Poorly-differentiated Carcinomas of the Thyroid are frequently thyroglobulin negative. Adenocarcinomas of non-thyroid origin do not react with this antibody.

Thyroid Peroxidase, RMAb




IHC of Thyroid Peroxidase on a FFPE Thyroid Adenoma Tissue

Thyroid Peroxidase (TPO) is an enzyme expressed mainly in the thyroid where it is secreted into colloid. Thyroid peroxidase oxidizes iodide ions to form iodine atoms for addition onto tyrosine residues on thyroglobulin for the production of thyroxine (T4) or triiodothyronine (T3), the thyroid hormones. In humans, thyroperoxidase is encoded by the TPO gene. The TPO gene consists of 17 exons, and is located on the short arm of chromosome 2.

Thyroid Peroxidase gene expression is under the regulation of thyroid stimulating hormone. In normal thyroid, expression of TPO described immunohistochemically is reported to produce a diffuse, fine, granular cytoplasmic stain in all follicular cells. TPO labels normal thyroid epithelial cells and thyroid tumor cells. The expression level in thyroid carcinomas is lower than that of normal and benign thyroid tumors. Decreased TPO immunoreactivity is an early event in follicular tumorigenesis, taking place before development of invasiveness in parallel with an acceleration of cell growth and appearance of cell atypia.

TIA-1, MAb



IHC of TIA-1 on a FFPE Spleen Tissue

TIA-1 (T-cell intracytoplasmic antigen) is a 15 kDa cytoplasmic granule-associated protein, expressed in lymphocytes processing cytolytic potential. TIA-1 is a member of an RNA-binding protein family and possesses nucleolytic activity against cytotoxic lymphocyte (CTL) target cells. It has been suggested that this protein may be involved in the induction of apoptosis as it preferentially recognizes poly(A) homopolymers and induces DNA fragmentation in CTL targets. The major granule-associated species is a 15 kDa protein thought to be derived from the carboxyl terminus of the 40 kDa product by proteolytic processing.

The expression of TIA-1 has been studied in Anaplastic Large Cell Lymphomas (ALCL), NK-cell Lymphomas, Peripheral T-cell Lymphomas, T-cell Lymphocytosis, B-cell Lymphomas and Lymphoblastic Leukemia, Hodgkin's, etc. Studies show that 60 to 70% of Anaplastic Large Cell Lymphomas react with TIA-1. TIA-1 reacts with most Large Granular Lymphocytic Leukemias, Hepatosplenic T-cell Lymphomas, intestinal T-cell Lymphomas, NK-like T-cell Lymphomas, NK-cell Lymphomas, nasal T/NK-cell Lymphomas, subcutaneous T-cell Lymphomas and Pulmonary Angiocentric Lymphomas of T- or NK-phenotype. All B-cell Lymphomas, Hodgkin's and Lymphoblastic Leukemias are negative for TIA-1.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP250
ISOTYPE: IgG
CONTROL: Thyroid, Thyroid Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP159
ISOTYPE: IgG
CONTROL: Thyroid, Thyroid Carcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: TIA-1
ISOTYPE: IgG1
CONTROL: Tonsil, Spleen, Anaplastic Large Cell Lymphoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2363 | Tinto Predilute | 3.0 ml |
| BSB 2364 | Tinto Predilute | 7.0 ml |
| BSB 2365 | Tinto Predilute | 15.0 ml |
| BSB 2366 | Concentrate | 0.1 ml |
| BSB 2367 | Concentrate | 0.5 ml |
| BSB 2368 | Concentrate | 1.0 ml |
| BSB 2369 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2615 | Tinto Predilute | 3.0 ml |
| BSB 2616 | Tinto Predilute | 7.0 ml |
| BSB 2617 | Tinto Predilute | 15.0 ml |
| BSB 2618 | Concentrate | 0.1 ml |
| BSB 2619 | Concentrate | 0.5 ml |
| BSB 2620 | Concentrate | 1.0 ml |
| BSB 2621 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6352 | Tinto Predilute | 3.0 ml |
| BSB 6353 | Tinto Predilute | 7.0 ml |
| BSB 6354 | Tinto Predilute | 15.0 ml |
| BSB 6355 | Concentrate | 0.1 ml |
| BSB 6356 | Concentrate | 0.5 ml |
| BSB 6357 | Concentrate | 1.0 ml |
| BSB 6358 | Control Slides | 5 |

TIA-1, RMAb



IHC of TIA-1 on a FFPE CLL/SLL Lymphoma Tissue

TIA-1 (T-cell intracytoplasmic antigen) is a 15 kDa cytoplasmic granule-associated protein, expressed in lymphocytes processing cytolytic potential. TIA-1 is a member of an RNA-binding protein family and possesses nucleolytic activity against cytotoxic lymphocyte (CTL) target cells. It has been suggested that this protein may be involved in the induction of apoptosis as it preferentially recognizes poly(A) homopolymers and induces DNA fragmentation in CTL targets. The major granule-associated species is a 15 kDa protein thought to be derived from the carboxyl terminus of the 40 kDa product by proteolytic processing.

The expression of TIA-1 has been studied in anaplastic large cell lymphomas (ALCL), NK-cell lymphomas, peripheral T-cell lymphomas, T-cell lymphocytosis, B-cell lymphomas and lymphoblastic leukemia, and Hodgkin's lymphoma, etc. Studies show that 60 to 70% of anaplastic large cell lymphomas react with TIA-1. TIA-1 also reacts with most large granular lymphocytic leukemias, hepatosplenic T-cell lymphomas, intestinal T-cell lymphomas, NK-like T-cell lymphomas, NK-cell lymphomas, nasal T/NK-cell lymphomas, subcutaneous T-cell lymphomas and pulmonary angiocentric lymphomas of T-or NK-phenotype. All B-cell lymphomas, Hodgkin's and lymphoblastic leukemias are negative for TIA-1.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: RBT-TIA1


ISOTYPE: IgG

CONTROL: Tonsil, Spleen, Anaplastic Large Cell Lymphoma

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

TIGIT, MMAb



IHC of TIGIT on a FFPE Testis Tissue

TIGIT (T cell immunoreceptor with Ig and ITIM domains) is an immune receptor present on some T cells and Natural Killer Cells (NK). It is also identified as WUCAM and Vstm3. TIGIT could bind to CD155 (PVR) on dendritic cells (DCs), macrophages, and other immune cells with high affinity, and also to CD112 (PVRL2) with lower affinity. Research has shown that TIGIT-Fc fusion protein could interact with PVR on dendritic cells and increase its IL-10 secretion level/decrease its IL-12 secretion level under LPS stimulation, and also inhibit T cell activation in vivo. TIGIT's inhibition of NK cytotoxicity can be blocked by antibodies against its interaction with PVR, and the activity is directed through its ITIM domain.

TIGIT is expressed on regulatory T cells (Tregs) and on activated CD4+ T, CD8+ T, and NK cells. TIGIT and PD-1 has been shown to be over-expressed on tumor antigen-specific (TA-specific) CD8+ T cells and CD8+ tumor infiltrating lymphocytes (TILs) from individuals with melanoma. Blockade of TIGIT and PD-1 led to increased cell proliferation, cytokine production, and degranulation of TA-specific CD8+ T cells and TIL CD8+ T cells, therefore it can be considered an immune checkpoint. Co-blockade of TIGIT and PD-1 pathways elicits tumor rejection in preclinical murine models.

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: BSB-152


ISOTYPE: IgG2b

CONTROL: Testis, Cervix, Fallopian Tube, Skin, Prostate, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma, Lung Adenocarcinoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human, Rat

TIM-3/HAVCR2/CD366, MMAb



IHC of TIM-3/HAVCR2/CD366 on a FFPE Colon Adenocarcinoma Tissue

T-cell Immunoglobulin and Mucin-domain Containing-3 (TIM-3), also known as Hepatitis A virus cellular receptor 2 (HAVCR2), is a protein that in humans is encoded by the HAVCR2 gene. TIM-3 is an immune checkpoint and together with other inhibitory receptors including programmed cell death protein 1 (PD-1) and lymphocyte activation gene 3 protein (LAG3) mediate CD8+ T-cell exhaustion. TIM-3 expression is up-regulated in tumor-infiltrating lymphocytes in Lung, Gastric, Head & Neck Cancers, Schwannoma, Melanoma and Follicular B-cell Non-Hodgkin Lymphoma. The TIM-3 pathway may interact with the PD-1 pathway in the dysfunctional CD8+ T cells and Tregs in cancer. TIM-3 is mainly expressed on activated CD8+ T cells and suppresses macrophage activation following PD-1 inhibition.

Upregulation has been observed in tumors progressing after anti-PD-1 therapy. It has been reported that early breast cancer patients with TIM-3+ iTILs have significantly improved breast cancer-specific survival whereas TIM-3+ sTILs did not reach statistical significance and it was concluded that the presence of TIM-3+ iTILs is an independent favorable prognostic factor in the whole cohort as well as among ER negative patients. In myelogenous leukemia (AML), upregulated TIM-3 during AML could reduce cytokine production. Co-expression of PD-1 and TIM-3 was correlated with AML progression. In glioma patients, TIM-3 was correlated with cancer immune escape and might be a potent target.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-163

ISOTYPE: IgG2c

CONTROL: Colon, Testis, Tonsil, Liver


LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3782-3 | Tinto Predilute | 3.0 ml |
| BSB-3782-7 | Tinto Predilute | 7.0 ml |
| BSB-3782-15 | Tinto Predilute | 15.0 ml |
| BSB-3782-01 | Concentrate | 0.1 ml |
| BSB-3782-05 | Concentrate | 0.5 ml |
| BSB-3782-1 | Concentrate | 1.0 ml |
| BSB-3782-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3751-3 | Tinto Predilute | 3.0 ml |
| BSB-3751-7 | Tinto Predilute | 7.0 ml |
| BSB-3751-15 | Tinto Predilute | 15.0 ml |
| BSB-3751-01 | Concentrate | 0.1 ml |
| BSB-3751-05 | Concentrate | 0.5 ml |
| BSB-3751-1 | Concentrate | 1.0 ml |
| BSB-3751-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3752-3 | Tinto Predilute | 3.0 ml |
| BSB-3752-7 | Tinto Predilute | 7.0 ml |
| BSB-3752-15 | Tinto Predilute | 15.0 ml |
| BSB-3752-01 | Concentrate | 0.1 ml |
| BSB-3752-05 | Concentrate | 0.5 ml |
| BSB-3752-1 | Concentrate | 1.0 ml |
| BSB-3752-CS | Control Slides | 5 |

TLE1, MAb

IHC of TLE1 on a FFPE Synovial Sarcoma Tissue

TLE1, MAb

IHC of TLE1 on FFPE Synovial Sarcoma Tissue

TMPRSS2, MAb

IHC of TMPRSS2 on FFPE SARS-CoV-2 Infected Lung Tissue

The Notch signaling pathway controls cellular interactions important for the specification of a variety of fates in both invertebrates and vertebrates. Key players in the Notch pathway are the TLE genes. TLEs associate with chromatin in live cells and specifically with Histone H3, but not with other core histones. Expression of the TLE genes, TLE1, TLE2, TLE3 and TLE4, correlate with immature epithelial cells that are progressing toward a terminally differentiated state, suggesting a role during epithelial differentiation.

Anti-TLE1 can be used to differentiate synovial sarcoma from other sarcomas, including histologically similar tumors such as malignant peripheral nerve sheath tumor.

The Notch signaling pathway controls cellular interactions important for the specification of a variety of fates in both invertebrates and vertebrates. Key players in the Notch pathway are the TLE genes. TLEs associate with chromatin in live cells and specifically with Histone H3, but not with other core histones. Expression of the TLE genes, TLE1, TLE2, TLE3 and TLE4, correlate with immature epithelial cells that are progressing toward a terminally differentiated state, suggesting a role during epithelial differentiation.

Anti-TLE1 can be used to differentiate synovial sarcoma from other sarcomas, including histologically similar tumors such as malignant peripheral nerve sheath tumor.

Transmembrane Serine Protease 2 is part of the serine protease family, which is active in many physiological and pathological pathways. TMPRSS2 is a 492 amino acid protein, containing a type II transmembrane domain, a receptor class A domain, a scavenger receptor cysteine-rich domain, and a protease domain. It is upregulated by androgen hormones particularly in the prostate, where it may contribute to inflammation by activating PAR2 and prostate cancer through somatic rearrangement. TMPRSS2 is also expressed in the GI tract, stomach, kidney and in the lung epithelium, where it may cleave epithelial sodium channels. The TMPRSS2-ERG fusion pair is a common somatic gene rearrangement occurring in about 50% of primary prostate cancers.

TMPRSS2-ERG fusion-positive tumors may be at higher risk for metastasis and influence from hormones, with a different androgen metabolism and higher insulin signaling than negative tumors. The protease domain also proteolytically cleaves and activates viral envelope glycoproteins, facilitating the cellular entry of human influenza and coronaviruses such as SARS-CoV-2.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: IF5

ISOTYPE: IgG1/K

CONTROL: Synovial Sarcoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-142

ISOTYPE: IgG2a

CONTROL: Synovial Sarcoma

LOCALIZATION: Cytoplasmic, Nuclear

SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-136

ISOTYPE: IgG1

CONTROL: Testis, Colon, Kidney, Brain, Stomach, Pancreas, Prostate

LOCALIZATION: Nuclear, Membranous


SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2314 | Tinto Predilute | 3.0 ml |
| BSB 2315 | Tinto Predilute | 7.0 ml |
| BSB 2316 | Tinto Predilute | 15.0 ml |
| BSB 2317 | Concentrate | 0.1 ml |
| BSB 2318 | Concentrate | 0.5 ml |
| BSB 2319 | Concentrate | 1.0 ml |
| BSB 2320 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3784-3 | Tinto Predilute | 3.0 ml |
| BSB-3784-7 | Tinto Predilute | 7.0 ml |
| BSB-3784-15 | Tinto Predilute | 15.0 ml |
| BSB-3784-01 | Concentrate | 0.1 ml |
| BSB-3784-05 | Concentrate | 0.5 ml |
| BSB-3784-1 | Concentrate | 1.0 ml |
| BSB-3784-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3703-3 | Tinto Predilute | 3.0 ml |
| BSB-3703-7 | Tinto Predilute | 7.0 ml |
| BSB-3703-15 | Tinto Predilute | 15.0 ml |
| BSB-3703-01 | Concentrate | 0.1 ml |
| BSB-3703-05 | Concentrate | 0.5 ml |
| BSB-3703-1 | Concentrate | 1.0 ml |
| BSB-3703-CS | Control Slides | 5 |

TNFa-IP2, MMab




IHC of TNFα-IP2 on FPPE SARS-CoV-2 infected Lung Tissue

Tumor Necrosis Factor alpha Induced-Protein 2 (also known as B94) is a 73-kDa polypeptide involved in pathways of inflammation, metastasis, tumor vasculature, and angiogenesis. The protein has an exocyst complex, Dsl1 complex, conserved oligomeric Golgi (COG) complex and the Golgi-associated retrograde protein (GARP) complex.

TNFαIP2 is found in epithelial cells, and immune cells exposed to Tumor Necrosis Factor alpha (TNFα), IL-1β, LPS, interferon-γ, Retinoic Acid, Latent Membrane Protein 1 (LMP1), and other pro-inflammatory cytokines. TNFαIP can inhibit NFκB to further reduce inflammation in renal dysfunction and septic shock, and can interact with GTPases to regulate breast cancer and HeLa cell actin cytoskeleton and cell structure. TNFαIP2 also participates in T-cell migration as an inflammatory regulator of chemokine secretion, and promotes metastasis and microvessel formation in nasopharyngeal carcinoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-141
ISOTYPE: IgG1
CONTROL: Testis, Tonsil, Lung, Kidney
LOCALIZATION: Membranous, Cytoplasmic
SPECIES REACTIVITY: Human, Mouse, Rat

Topoisomerase IIα, RMAb



IHC of Topoisomerase II alpha on a FPPE HSIL of the Cervix Tissue


DNA Topoisomerase II alpha (Topo IIα) is a nucleic enzyme that affects the topological structure of DNA by interacting with the double-helix DNA, thus playing an important role in DNA replication, transcription, recombination, condensation, and segregation. Type II topoisomerases cut both strands of the DNA helix simultaneously in order to change the linking number of the molecule. Topo IIα is essential in the separation of daughter strands at the end of replication. Failure to separate these strands leads to cell death. In cancers, the Topo IIα is highly expressed in highly-proliferating cells.

Topo IIα has been identified by DNA microarray and transcriptional profiling as a gene that is overexpressed in Cervical Carcinomas. The TOP2A gene is approximately 30 kb in size and encodes a 170 kDa protein. Topo IIα protein is expressed in proliferating cells and in numerous human malignant tumors, including colon, gastric and breast cancers, Lymphomas and others. In certain cancers, such as Peripheral Nerve Sheath Tumors, high expression of this protein is also associated with poor patient survival. Type II topoisomerases are the targets for anticancer drugs, such as topoisomerase II inhibitor therapies like the anthracyclines (Doxorubicin and Epirubicin).

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-topo 2a
ISOTYPE: IgG
CONTROL: Testis, Skin, Colon, Fallopian Tube, Tonsil, Lymph Node, Spleen, HSIL Cervical Cancer, Breast Cancer, Bladder TCC
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6338 | Tinto Predilute | 3.0 ml |
| BSB 6339 | Tinto Predilute | 7.0 ml |
| BSB 6340 | Tinto Predilute | 15.0 ml |
| BSB 6341 | Concentrate | 0.1 ml |
| BSB 6342 | Concentrate | 0.5 ml |
| BSB 6343 | Concentrate | 1.0 ml |
| BSB 6344 | Control Slides | 5 |

Toxoplasma gondii, RPAb



IHC of Toxoplasma gondii on a FPPE Brain Tissue


Toxoplasma gondii is a genus of parasitic protozoa (cats being the definitive host). It can also be carried by the vast majority of warm-blooded animals, including humans. Toxoplasma gondii belongs to the phylum Apicomplexa and is the only known member species of the genus Toxoplasma. The life cycle of T. gondii has two phases. The sexual part of the life cycle (coccidia-like) occurs only in members of the Felidae family (domestic and wild cats), which makes these animals the parasite's primary host. The asexual part of the life cycle can occur in any warm-blooded animal, such as other mammals (including felines) and birds.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Toxoplasma gondii Infected Tissue
LOCALIZATION: Cell Wall
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6043 | Tinto Predilute | 3.0 ml |
| BSB 6044 | Tinto Predilute | 7.0 ml |
| BSB 6045 | Tinto Predilute | 15.0 ml |
| BSB 6046 | Concentrate | 0.1 ml |
| BSB 6047 | Concentrate | 0.5 ml |
| BSB 6048 | Concentrate | 1.0 ml |
| BSB 6049 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3708-3 | Tinto Predilute | 3.0 ml |
| BSB-3708-7 | Tinto Predilute | 7.0 ml |
| BSB-3708-15 | Tinto Predilute | 15.0 ml |
| BSB-3708-01 | Concentrate | 0.1 ml |
| BSB-3708-05 | Concentrate | 0.5 ml |
| BSB-3708-1 | Concentrate | 1.0 ml |
| BSB-3708-CS | Control Slides | 5 |

TRAcP, MAb



IHC of TRAcP on a FFPE Hairy Cell Leukemia Tissue

Tartrate-resistant acid phosphatase (TRAcP) is a glycosylated monomeric metallo-enzyme expressed in mammals. It has a molecular weight of approximately 35 kDa, a basic isoelectric point (7.6 - 9.5), and optimal activity in acidic conditions. TRAcP is synthesized as a latent proenzyme and is activated by proteolytic cleavage and reduction. Normally, TRAcP is highly expressed by osteoclasts, activated macrophages, neurons and endometrium during pregnancy. There are also certain pathological conditions whereby expression of TRAcP is increased. These include patients with Leukemic Reticuloendotheliosis (Hairy Cell Leukemia), Gaucher's Disease, HIV-induced Encephalopathy, Osteoclastoma and in osteoporosis and metabolic bone diseases.

Anti-TRAcP antibody labels the cells of Hairy Cell Leukemia (HCL) with a high degree of sensitivity and specificity. Other cells stained with this antibody are tissue macrophages and osteoclasts, which also express abundant TRAcP activity.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 9C5


ISOTYPE: IgG2b

CONTROL: Tonsil, Spleen, Lymph Node, Hairy cell Leukemia

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

Treponema Pallidum, RPAb



IHC of Treponema Pallidum on a FFPE Infected Skin Tissue

reponema pallidum is a spirochaete bacterium. The treponemes have a cytoplasmic and an outer membrane. The shape of *T. pallidum* is flat and wavy, unlike the other spirochetes, which are helical. Using light microscopy, treponemes are only visible using dark field illumination. They are Gram negative, but some regard them too thin to be Gram stained. *T. pallidum* is a motile spirochaete that is generally acquired by close sexual contact, entering the host via breaches in squamous or columnar epithelium. The organism can also be transmitted to a fetus by transplacental passage during the later stages of pregnancy, giving rise to congenital syphilis. The helical structure of *T. pallidum* allows it to move in a corkscrew motion through a viscous medium such as mucus. It gains access to the host's blood and lymph systems through tissue and mucous membranes.

Detection of *Treponema pallidum* can be difficult, and the correct diagnosis of secondary syphilis is critical. Diagnosis of syphilis is usually based on clinical presentation, dark-field microscope analysis, and serological tests. *T. pallidum* can be also evidenced by Immunohistochemistry in up to 90% of the samples with the bacteria located in the epidermis and the upper dermis of formalin-fixed paraffin-embedded tissues.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Treponema Pallidum Infected Tissue

LOCALIZATION: Cell Wall

SPECIES REACTIVITY: Eubacteria

TROP-2/EGP-1, MAb



IHC of TROP-2/EGP-1 on FFPE Pancreatic Adenocarcinoma Tissue

Trop-2 or Tumor-associated calcium signal transducer 2, also known as epithelial glycoprotein-1 antigen (EGP-1), is a protein that in humans is encoded by the TACSTD2 gene. This transmembrane glycoprotein functions in a variety of cell signaling pathways and was first elucidated as a transducer of an intracellular calcium signal. Trop-2 expression has been demonstrated to depend on a large variety of transcription factors. Trop-2 is involved in several cell signaling pathways, of which many are associated with tumorigenesis.

In Thyroid Cancer cell invasion, Trop-2 signal transduction has been seen as a downstream effect of the ERK/JNK pathways, where its signaling enhances stem cell-like properties of cancer cells, as Trop-2 regulates proliferation and self-renewal through beta-catenin signaling. In a study, the majority of Papillary Thyroid Carcinoma (PTC) specimens were positive for Trop-2; however, the pattern of staining differed significantly between the histopathological variants. All Papillary microcarcinomas (mPTC), PTC classic variant (PTC cv), and tall cell variants (PTC tcv) were Trop-2 positive, with mainly diffuse staining. In contrast, less than half of the PTC follicular variant specimens were positive for Trop-2, with only focal immunoreactivity. Trop-2 may play a role in tumor progression given the involvement in several molecular pathways traditionally associated with cancer development.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-148

ISOTYPE: IgG1

CONTROL: Breast, Prostate, Skin, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma, Pancreatic Carcinoma

LOCALIZATION: Membranous


SPECIES REACTIVITY: Human, Mouse, Rat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5980 | Tinto Predilute | 3.0 ml |
| BSB 5981 | Tinto Predilute | 7.0 ml |
| BSB 5982 | Tinto Predilute | 15.0 ml |
| BSB 5983 | Concentrate | 0.1 ml |
| BSB 5984 | Concentrate | 0.5 ml |
| BSB 5985 | Concentrate | 1.0 ml |
| BSB 5986 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3232 | Tinto Predilute | 3.0 ml |
| BSB 3233 | Tinto Predilute | 7.0 ml |
| BSB 3234 | Tinto Predilute | 15.0 ml |
| BSB 3235 | Concentrate | 0.1 ml |
| BSB 3236 | Concentrate | 0.5 ml |
| BSB 3237 | Concentrate | 1.0 ml |
| BSB 3238 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3754-3 | Tinto Predilute | 3.0 ml |
| BSB-3754-7 | Tinto Predilute | 7.0 ml |
| BSB-3754-15 | Tinto Predilute | 15.0 ml |
| BSB-3754-01 | Concentrate | 0.1 ml |
| BSB-3754-05 | Concentrate | 0.5 ml |
| BSB-3754-1 | Concentrate | 1.0 ml |
| BSB-3754-CS | Control Slides | 5 |

Tryptase, MMab




IHC of Tryptase on a FFPE *H. pylori* Infected Stomach Tissue

Tryptase, RMab



IHC of Tryptase on a FFPE Uterus Tissue

TSH, MMab



IHC of TSH on a FFPE Pituitary Tissue

Tryptase is the most abundant secretory granule-derived serine proteinase contained in mast cells and has recently been used as a marker for mast cell activation. It is involved in allergenic response and is suspected to act as a mitogen for fibroblast lines. Elevated levels of serum tryptase occur in both anaphylactic and anaphylactoid reactions, but a negative test does not exclude anaphylaxis. Mast cells contain a number of preformed chemical mediators such as histamine, chymase, carboxypeptidase and proteolytic tryptase.

Human mast cell tryptase is considered to be an important marker of mast cell activation as well as an important mediator of inflammation. Anti-tryptase is a good marker for mast cells, basophils, and their derivatives.

Tryptase is the most abundant secretory granule-derived serine proteinase contained in mast cells and has recently been used as a marker for mast cell activation. It is involved in allergenic response and is suspected to act as a mitogen for fibroblast lines. Elevated levels of serum tryptase occur in both anaphylactic and anaphylactoid reactions, but a negative test does not exclude anaphylaxis. Mast cells contain a number of preformed chemical mediators such as histamine, chymase, carboxypeptidase and proteolytic tryptase.

Human mast cell tryptase is considered to be an important marker of mast cell activation as well as an important mediator of inflammation. Anti-tryptase is a good marker for mast cells, basophils, and their derivatives.

Thyroid-stimulating hormone (TSH or thyrotropin) is a hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland which regulates the endocrine function of the thyroid gland. TSH stimulates the thyroid gland to secrete the hormones thyroxine (T4) and triiodothyronine (T3). TSH production is controlled by a Thyrotropin-Releasing Hormone (TRH), which is manufactured in the hypothalamus and transported to the pituitary gland, where it increases TSH production and release. Somatostatin is also produced by the hypothalamus and has an opposite effect on the pituitary production of TSH, decreasing or inhibiting its release.

TSH is a useful marker in classification of pituitary tumors and the study of pituitary disease. TSH antibody primarily reacts with TSH-producing cells.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: G3
ISOTYPE: IgG1
CONTROL: Liver, Kidney, Tonsil, Uterus Cervix, Skin, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP259
ISOTYPE: IgG
CONTROL: Liver, Kidney, Tonsil, Uterus, Cervix, Skin, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-56
ISOTYPE: IgG1/K
CONTROL: Pituitary
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 5987 | Tinto Predilute | 3.0 ml | BSB 2377 | Tinto Predilute | 3.0 ml | BSB 5994 | Tinto Predilute | 3.0 ml |
| BSB 5988 | Tinto Predilute | 7.0 ml | BSB 2378 | Tinto Predilute | 7.0 ml | BSB 5995 | Tinto Predilute | 7.0 ml |
| BSB 5989 | Tinto Predilute | 15.0 ml | BSB 2379 | Tinto Predilute | 15.0 ml | BSB 5996 | Tinto Predilute | 15.0 ml |
| BSB 5990 | Concentrate | 0.1 ml | BSB 2380 | Concentrate | 0.1 ml | BSB 5997 | Concentrate | 0.1 ml |
| BSB 5991 | Concentrate | 0.5 ml | BSB 2381 | Concentrate | 0.5 ml | BSB 5998 | Concentrate | 0.5 ml |
| BSB 5992 | Concentrate | 1.0 ml | BSB 2382 | Concentrate | 1.0 ml | BSB 5999 | Concentrate | 1.0 ml |
| BSB 5993 | Control Slides | 5 | BSB 2383 | Control Slides | 5 | BSB 6000 | Control Slides | 5 |

TSH, RMab




IHC of TSH on a FFPE Pituitary Tissue

TTF-1, MAb



IHC of TTF-1 on a FFPE Lung Tissue

Tyrosinase, MAb



IHC of Tyrosinase on a FFPE Malignant Melanoma Tissue

Thyroid-stimulating hormone (TSH or thyrotropin) is a hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland which regulates the endocrine function of the thyroid gland. TSH stimulates the thyroid gland to secrete the hormones thyroxine (T4) and triiodothyronine (T3). TSH production is controlled by a Thyrotropin-Releasing Hormone (TRH), which is manufactured in the hypothalamus and transported to the pituitary gland, where it increases TSH production and release. Somatostatin is also produced by the hypothalamus and has an opposite effect on the pituitary production of TSH, decreasing or inhibiting its release.

TSH is a useful marker in classification of pituitary tumors and the study of pituitary disease. TSH antibody primarily reacts with TSH-producing cells.

Thyroid transcription factor-1 (TTF-1) is a protein that regulates transcription of genes specific to the thyroid, lung and diencephalon. It is also known as thyroid-specific enhancer binding protein and NKX-2. It is used as a marker to determine if a tumor originates in the lung or thyroid. TTF-1 positive cells are found in Type II pneumocytes and Clara cells in the lung. In the thyroid, follicular and parafollicular cells are positive.

TTF-1 is useful in differentiating primary Adenocarcinoma of the Lung from Metastatic Carcinomas of the breast and Malignant Mesothelioma. It can also be used to differentiate Small-Cell Lung Carcinoma from lymphoid infiltrates. For lung cancers, Adenocarcinomas are usually positive, while Squamous Cell Carcinomas and Large Cell Carcinomas are rarely positive. Small-Cell Carcinomas (of any primary site) are usually positive.

Tyrosinase is an enzyme that catalyzes the oxidation of phenols (such as tyrosine) and is widespread in plants and animals. Tyrosinase is a copper-containing enzyme present in plant and animal tissues that catalyzes the production of melanin and other pigments from tyrosine by oxidation. The gene for Tyrosinase is regulated by the Microphthalmia-associated transcription factor. A mutation in the tyrosinase gene leads to impaired tyrosinase production resulting in Type I Oculocutaneous Albinism, a hereditary disease that affects 1 in 17,000 people in the U.S.

Anti-Tyrosinase has been found to be quite specific for melanotic lesions such as Malignant Melanoma and Melanotic Neurofibroma. Essentially no carcinomas express this marker.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP254
ISOTYPE: IgG
CONTROL: Pituitary
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, predicted
 Rabbit

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 8G7G3/1
ISOTYPE: IgG1
CONTROL: Lung, Thyroid, Adenocarcinoma of Lung
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Dog

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-6
ISOTYPE: IgG2a
CONTROL: Skin, Malignant, Melanoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 2622 | Tinto Predilute | 3.0 ml | BSB 6001 | Tinto Predilute | 3.0 ml | BSB 6008 | Tinto Predilute | 3.0 ml |
| BSB 2623 | Tinto Predilute | 7.0 ml | BSB 6002 | Tinto Predilute | 7.0 ml | BSB 6009 | Tinto Predilute | 7.0 ml |
| BSB 2624 | Tinto Predilute | 15.0 ml | BSB 6003 | Tinto Predilute | 15.0 ml | BSB 6010 | Tinto Predilute | 15.0 ml |
| BSB 2625 | Concentrate | 0.1 ml | BSB 6004 | Concentrate | 0.1 ml | BSB 6011 | Concentrate | 0.1 ml |
| BSB 2626 | Concentrate | 0.5 ml | BSB 6005 | Concentrate | 0.5 ml | BSB 6012 | Concentrate | 0.5 ml |
| BSB 2627 | Concentrate | 1.0 ml | BSB 6006 | Concentrate | 1.0 ml | BSB 6013 | Concentrate | 1.0 ml |
| BSB 2628 | Control Slides | 5 | BSB 6007 | Control Slides | 5 | BSB 6014 | Control Slides | 5 |

Uroplakin III, RPAb




IHC of Uroplakin III on a FFPE Transitional Cell Carcinoma Tissue

Uroplakin III, RMAb



IHC of Uroplakin III on a FFPE Bladder TCC Tissue

Varicella Zoster Virus, MAb



IHC of Varicella Zoster Virus on a FFPE Infected Tissue

Uroplakins (UPs) are a family of transmembrane proteins (UPs Ia, Ib, II and III) that are specific differentiation products of urothelial cells. In non-neoplastic mammalian urothelium, UPs are expressed in the luminal surface plasmalemma of superficial (umbrella) cells, forming complexes of 16-nm crystalline particles. Uroplakin III is expressed in urothelial carcinomas, whereas many non urothelial carcinomas were UPIII-negative. Recent studies have shown that UP expression might reflect the malignant potential of urothelial cancer cells as well as being cytodifferential markers of urothelial cells.

UPIII expression is strongly associated with lower tumor grades and lack of UPIII expression in urothelial tumors of the upper urinary tract is associated with much higher rates of metastases. Five-year specific survival is much worse for UPIII negative tumors (54%) than for UPIII positive tumors (100%). Apparently UPIII expression is a better indicator for the malignant potential of the tumor than the grade of the tumor.

Uroplakins (UPs) are a family of transmembrane proteins (UPs Ia, Ib, II and III) that are specific differentiation products of urothelial cells. In non-neoplastic mammalian urothelium, UPs are expressed in the luminal surface plasmalemma of superficial (umbrella) cells, forming complexes of 16-nm crystalline particles. Uroplakin III is expressed in urothelial carcinomas, whereas many non urothelial carcinomas were UPIII-negative. Recent studies have shown that UP expression might reflect the malignant potential of urothelial cancer cells as well as being cytodifferential markers of urothelial cells.

UPIII expression is strongly associated with lower tumor grades and lack of UPIII expression in urothelial tumors of the upper urinary tract is associated with much higher rates of metastases. Five-year specific survival is much worse for UPIII negative tumors (54%) than for UPIII positive tumors (100%). Apparently UPIII expression is a better indicator for the malignant potential of the tumor than the grade of the tumor.

Varicella Zoster Virus (VZV) is a member of the human herpes virus family and causes two distinct clinical manifestations: chickenpox and shingles.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Bladder, Bladder Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP321
ISOTYPE: IgG
CONTROL: Bladder, Bladder Carcinoma
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: SG1-1, SG1-4, NCP-1& IE-62
ISOTYPE: Mixed
CONTROL: Varicella Zoster Virus Infected Tissue
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2286 | Tinto Predilute | 3.0 ml |
| BSB 2287 | Tinto Predilute | 7.0 ml |
| BSB 2288 | Tinto Predilute | 15.0 ml |
| BSB 2289 | Concentrate | 0.1 ml |
| BSB 2290 | Concentrate | 0.5 ml |
| BSB 2291 | Concentrate | 1.0 ml |
| BSB 2292 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3239 | Tinto Predilute | 3.0 ml |
| BSB 3240 | Tinto Predilute | 7.0 ml |
| BSB 3241 | Tinto Predilute | 15.0 ml |
| BSB 3242 | Concentrate | 0.1 ml |
| BSB 3243 | Concentrate | 0.5 ml |
| BSB 3244 | Concentrate | 1.0 ml |
| BSB 3245 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2293 | Tinto Predilute | 3.0 ml |
| BSB 2294 | Tinto Predilute | 7.0 ml |
| BSB 2295 | Tinto Predilute | 15.0 ml |
| BSB 2296 | Concentrate | 0.1 ml |
| BSB 2297 | Concentrate | 0.5 ml |
| BSB 2298 | Concentrate | 1.0 ml |
| BSB 2299 | Control Slides | 5 |

VEGF, RMAb

IHC of VEGF on a FFPE Kidney Tissue

Villin, MAb

IHC of Villin on a FFPE Colon Tissue

Villin, RMAb

IHC of Villin on a FFPE Colon Tissue

Vascular endothelial growth factor (VEGF) is an important signaling protein involved in both vasculogenesis (the de novo formation of the embryonic circulatory system) and angiogenesis (the growth of blood vessels from pre-existing vasculature). As its name implies, VEGF activity is restricted mainly to cells of the vascular endothelium, although it has an effect on a limited number of other cell types (e.g., stimulation monocyte/macrophage migration).

VEGF has been implicated with poor prognosis in breast cancer. The overexpression of VEGF may be an early step in the process of metastasis, a step involved in the "angiogenic" switch. Although VEGF has been correlated with poor survival, its exact mechanism of action in the progression of tumors remains unclear. VEGF is also released in rheumatoid arthritis in response to TNF-alpha, increasing endothelial permeability and swelling and also stimulating angiogenesis (formation of capillaries). Once released, VEGF may elicit several responses. It may cause a cell to survive, move, or further differentiate.

Villin is an actin-binding protein that contains gelsolin domains capped by a "headpiece" consisting of a fast and independently-folding three-helix bundle stabilized by hydrophobic interactions. The headpiece domain is a commonly-studied protein in molecular dynamics due to its fast-folding kinetics and short primary sequence. It is a regulator of the actin cytoskeleton and is expressed mainly in the brush border in animals.

Anti-Villin labels the brush border area in the gastrointestinal mucosal epithelium. This antibody has been useful in differentiating Gastrointestinal Adenocarcinoma, Neuroendocrine Carcinomas and Ovarian Adenocarcinomas from Adenocarcinomas of other organs. Also labeled by this antibody are Merkel cells of the skin.

Villin is an actin-binding protein that contains gelsolin domains capped by a "headpiece" consisting of a fast and independently-folding three-helix bundle stabilized by hydrophobic interactions. The headpiece domain is a commonly-studied protein in molecular dynamics due to its fast-folding kinetics and short primary sequence. It is a regulator of the actin cytoskeleton and is expressed mainly in the brush border of the GI mucosal epithelium in animals.

Anti-Villin labels the brush border area in the gastrointestinal mucosal epithelium. This antibody has been useful in differentiating Gastrointestinal Adenocarcinoma, Neuroendocrine Carcinomas and Ovarian Adenocarcinomas from Adenocarcinomas of other organs. Also labeled by this antibody are Merkel cells of the skin.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-VEGF
ISOTYPE: IgG
CONTROL: Placenta, Angioma, Angiosarcoma & Soft Tissue
LOCALIZATION: Cell Surface, Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: CWWB1
ISOTYPE: IgG1
CONTROL: Kidney, Colon, Small Bowel Mucosa, Colonic Mucosa
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP163
ISOTYPE: IgG
CONTROL: Kidney, Colon, Small Bowel Mucosa, Colonic Mucosa
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6050 | Tinto Predilute | 3.0 ml |
| BSB 6051 | Tinto Predilute | 7.0 ml |
| BSB 6052 | Tinto Predilute | 15.0 ml |
| BSB 6053 | Concentrate | 0.1 ml |
| BSB 6054 | Concentrate | 0.5 ml |
| BSB 6055 | Concentrate | 1.0 ml |
| BSB 6056 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6015 | Tinto Predilute | 3.0 ml |
| BSB 6016 | Tinto Predilute | 7.0 ml |
| BSB 6017 | Tinto Predilute | 15.0 ml |
| BSB 6018 | Concentrate | 0.1 ml |
| BSB 6019 | Concentrate | 0.5 ml |
| BSB 6020 | Concentrate | 1.0 ml |
| BSB 6021 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2300 | Tinto Predilute | 3.0 ml |
| BSB 2301 | Tinto Predilute | 7.0 ml |
| BSB 2302 | Tinto Predilute | 15.0 ml |
| BSB 2303 | Concentrate | 0.1 ml |
| BSB 2304 | Concentrate | 0.5 ml |
| BSB 2305 | Concentrate | 1.0 ml |
| BSB 2306 | Control Slides | 5 |

Vimentin, MAb




IHC of Vimentin on a FFPE Colon Tissue

Vimentin, RMAb



IHC of Vimentin on a FFPE Melanoma Tissue

WT1, MAb



IHC of WT1 on a FFPE Testicular Cancer Tissue

Vimentin is a member of the intermediate filament family of proteins. Intermediate filaments are an important structural feature of eukaryotic cells. Together with microtubules and actin microfilaments, they make up the cytoskeleton.

Expression of vimentin, when used in conjunction with keratin, is helpful in distinguishing melanomas from Undifferentiated Carcinomas and Large-Cell Lymphomas. All Melanomas and Schwannomas react strongly with vimentin. This antibody recognizes a 57 kDa intermediate filament. It labels a variety of mesenchymal cells, including melanocytes, lymph cells, endothelial cells and fibroblasts. Non-reactivity of vimentin antibody is often considered more useful than its presence, since there are a few tumors that do not contain vimentin (e.g., Hepatoma and Seminoma).

Vimentin is a member of the intermediate filament family of proteins. Intermediate filaments are an important structural feature of eukaryotic cells. Together with microtubules and actin microfilaments, they make up the cytoskeleton.

Expression of vimentin, when used in conjunction with keratin, is helpful in distinguishing melanomas from Undifferentiated Carcinomas and Large-Cell Lymphomas. All Melanomas and Schwannomas react strongly with vimentin. This antibody recognizes a 57 kDa intermediate filament. It labels a variety of mesenchymal cells, including melanocytes, lymph cells, endothelial cells and fibroblasts. Non-reactivity of vimentin antibody is often considered more useful than its presence, since there are a few tumors that do not contain vimentin (e.g., Hepatoma and Seminoma).

Wilms' Tumor Protein (WT1) is a suppressor gene located on Chromosome 11p13. Mutations of the WT1 gene on Chromosome 11 are observed in approximately 20% of Wilms tumors. At least half of the Wilms tumors with mutations in WT1 also carry mutations in CTNNB1, the gene encoding the proto-oncogene beta-catenin.

Wilms' tumor is a neoplasm of the kidneys that typically occurs in children. It is also known as a Nephroblastoma. WT1 has been identified in proliferative mesothelial cells, Malignant Mesothelioma, Ovarian Cystadenocarcinoma, Gonadoblastoma, Nephroblastoma and Desmoplastic Small Round Cell Tumor. Lung Adenocarcinomas rarely stain positive with this antibody.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: V9

ISOTYPE: IgG/K

CONTROL: Tonsil, Lymph Node, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat, Rat, Rabbit, Cattle, Horse

ANTIBODY TYPE: Rabbit Monoclonal

CLONE: EP21

ISOTYPE: IgG

CONTROL: Tonsil, Lymph Node, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Predicted: Mouse, Rat, Rhesus Monkey

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 6F-H2

ISOTYPE: IgG1/K

CONTROL: Kidney, Testis, Malignant Mesothelioma, Fallopian Tubes

LOCALIZATION: Nuclear


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6022 | Tinto Predilute | 3.0 ml |
| BSB 6023 | Tinto Predilute | 7.0 ml |
| BSB 6024 | Tinto Predilute | 15.0 ml |
| BSB 6025 | Concentrate | 0.1 ml |
| BSB 6026 | Concentrate | 0.5 ml |
| BSB 6027 | Concentrate | 1.0 ml |
| BSB 6028 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2307 | Tinto Predilute | 3.0 ml |
| BSB 2308 | Tinto Predilute | 7.0 ml |
| BSB 2309 | Tinto Predilute | 15.0 ml |
| BSB 2310 | Concentrate | 0.1 ml |
| BSB 2311 | Concentrate | 0.5 ml |
| BSB 2312 | Concentrate | 1.0 ml |
| BSB 2313 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6029 | Tinto Predilute | 3.0 ml |
| BSB 6030 | Tinto Predilute | 7.0 ml |
| BSB 6031 | Tinto Predilute | 15.0 ml |
| BSB 6032 | Concentrate | 0.1 ml |
| BSB 6033 | Concentrate | 0.5 ml |
| BSB 6034 | Concentrate | 1.0 ml |
| BSB 6035 | Control Slides | 5 |

YAP1, MMab



IHC of YAP1 on FFPE Colon Adenocarcinoma Tissue

Yes-Associated Protein 1 (YAP1) is a regulatory protein that binds the Src family tyrosine kinase YES, promotes transcriptional activation of cyclin E and inhibition of caspase-induced apoptosis. YAP1 is found in the cytoplasm, though transported into the nucleus to perform its roles in proliferation and assist in tissue repair and growth. YAP1 is essential for cancer growth, contributing to stem-cell-like features of proliferation and metastasis and inducing these features in neighboring cells. YAP1 may be upregulated or have higher nuclear localization in more malignant and metastatic tumors, as seen in Lung, Breast, Colorectal, Liver, Gastric, Pancreatic, and Brain Cancers.

YAP1 is frequently overexpressed in Mammary carcinoma, Glioblastoma and Squamous Cell Carcinoma, Pancreatic, Oral, Cervical, Ovarian and Lung Cancers. YAP1 has been found to be expressed in fetal and adult brain regions known to harbor neural progenitor cells but there was little YAP1 immunoreactivity in the adult cerebral cortex. YAP1 protein was also readily detected in the nuclei of human brain tumors. In Medulloblastoma, expression varied between histologic subtypes and was most prominent in Nodular/Desmoplastic tumors. In Gliomas it was frequently expressed in infiltrating Astrocytomas and Oligodendrogliomas, but rarely in Pilocytic Astrocytomas.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: BSB-146


ISOTYPE: IgG1

CONTROL: Testis, Placenta, Breast, Fallopian Tube, Transitional Cell Carcinoma, HER2 Negative Breast Cancer

LOCALIZATION: Nuclear, Cytoplasmic

SPECIES REACTIVITY: Human, Mouse, Rat

ZAP-70, MMab



IHC of ZAP-70 on a FFPE Lymphoma Tissue

ZAP-70 is an abbreviation for Zeta-chain-associated protein kinase 70 (70 is the molecular weight in kDa). The protein is a member of the protein-tyrosine kinase family. ZAP-70 is normally expressed in T-cells and natural killer cells and has a critical role in the initiation of T-cell signaling.

ZAP-70 in B-cells is used as a prognostic marker in identifying different forms of Chronic Lymphocytic Leukemia (CLL). ZAP-70 protein is expressed in leukemic cells in approximately 25% of Chronic Lymphocytic Leukemia (CLL) cases as well. ZAP-70 expression is an excellent surrogate marker for the distinction between the Ig-mutated (ZAP-70 negative) and Ig-unmutated (ZAP-70 positive) CLL subtypes and can identify patient groups with divergent clinical courses. The ZAP-70 positive Ig-unmutated CLL cases have a poorer prognosis.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: 2F3.2





ISOTYPE: IgG2a

CONTROL: Tonsil, Lymph Node, Chronic Lymphocytic Leukemia

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|
| BSB-3755-3 | Tinto Predilute | 3.0 ml | BSB 6036 | Tinto Predilute | 3.0 ml |
| BSB-3755-7 | Tinto Predilute | 7.0 ml | BSB 6037 | Tinto Predilute | 7.0 ml |
| BSB-3755-15 | Tinto Predilute | 15.0 ml | BSB 6038 | Tinto Predilute | 15.0 ml |
| BSB-3755-01 | Concentrate | 0.1 ml | BSB 6039 | Concentrate | 0.1 ml |
| BSB-3755-05 | Concentrate | 0.5 ml | BSB 6040 | Concentrate | 0.5 ml |
| BSB-3755-1 | Concentrate | 1.0 ml | BSB 6041 | Concentrate | 1.0 ml |
| BSB-3755-CS | Control Slides | 5 | BSB 6042 | Control Slides | 5 |



Fast Mohs Immunohistochemistry

Quality Antibodies | Sensitive Detection | Less Time

*A selection of antibodies, detection systems and chromogens for
fast IHC cancer detection after Mohs surgery*

FAST 5-20 MINUTE IHC PROTOCOLS

Antibodies For Fast Mohs IHC

The most common reason for recurrence of skin cancers after Mohs micrographic surgery is residual undetected tumor. Frozen tissue Immunohistochemistry (IHC) has demonstrated greater sensitivity than routine H&E stains, which are difficult to interpret. Bio SB has developed a highly sensitive non-biotin monovalent Fab micropolymer IHC detection system for the detection of IVD antibodies for melanoma, BCC, SCC and other Mohs surgery related conditions. Our innovative IHC detection systems and high affinity monoclonal antibodies, have opened the doors for a faster and accurate immunohistochemistry applicable to Mohs surgery.

BCC, SCC and Carcinoma Antibodies

TintoFast Androgen Receptor, MMab

TintoFast CD31, MMab

TintoFast Cytokeratin 5/6, MMab

TintoFast Cytokeratin Cocktail, AE1/AE3, MMab

TintoFast EpCAM BerEP4, MMab

TintoFast EMA, MMab

TintoFast p40, RMab

TintoFast p63, MMab

TintoFast Podoplanin, MMab

TintoFast NGFR, MMab

Dermatofibrosarcoma Protuberans (DFSP) Antibodies

TintoFast CD34, MMab

Melanoma Antibodies

TintoFast Ki-67, RMab

TintoFast MART-1, MMab

TintoFast Melanoma Cocktail (MART-1, Melanosome HMB-45, Tyrosinase), MMab

TintoFast Melanosome HMB-45, MMab

TintoFast MiTF, MMab

TintoFast PRAME, RMab

TintoFast SOX-10, MMab

Neuroendocrine Tumor Antibodies

TintoFast Chromogranin A, MMab

TintoFast Synaptophysin, RPab

Paget's Disease Antibodies

TintoFast Cytokeratin 7, MMab

Primary Cutaneous Adnexal Neoplasm Antibodies

TintoFast Adipophilin, MMab


TintoFast Androgen Receptor, MMab

TintoFast Cytokeratin 5/6, MMab


TintoFast EMA, MMab

TintoFast p63, MMab


TintoFast Podoplanin, MMab



SOX-10 using Mohs PolyDetector Plus HRP Green on a Frozen Mohs Melanoma




PRAME using Mohs PolyDetector Plus HRP Green and MART-1 using ALK Scarlet on a Frozen Melanoma Tissue



CK AE1 / AE3 using Mohs Poly-Detector Plus HRP Green on a Frozen Mohs Melanoma

TintoFast Adipophilin, MAb



IHC of Adipophilin on a Frozen Acetone-fixed Squamous Cell Carcinoma Tissue

Adipose differentiation-related protein, also known as Perilipin 2 (PLIN2), ADRP or Adipophilin, is a protein which in humans is encoded by the ADFP gene. Adipophilin is associated with the globule surface membrane material. This protein is a major constituent of the globule surface. Increase in mRNA levels is one of the earliest indications of adipocyte differentiation. Adipophilin is present in a wide range of cultured cell lines, including fibroblasts, endothelial, and epithelial cells. In tissues, however, expression of Adipophilin is restricted to certain cell types, such as Lactating Mammary epithelial cells, Adrenal Cortex cells, Sertoli and Leydig cells of the male reproductive system, and Steatosis or fatty change hepatocytes in alcoholic liver Cirrhosis.


Adipophilin expresses in various Sebaceous lesions and other Cutaneous Tumors with a clear cell histology that may mimic sebaceous differentiation. In periocular lesions, it is effective in helping to exclude Basal Cell Carcinoma and Squamous Cell Carcinoma when sebaceous carcinoma is under consideration.

Adipophilin is suitable for immunostaining and is helpful in the identification of intracytoplasmic lipids, as seen in Sebaceous lesions. It is especially helpful in identifying intracytoplasmic lipid vesicles in poorly differentiated Sebaceous Carcinomas in challenging cases such as small periocular biopsies specimens.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-91
ISOTYPE: IgG1
CONTROL: Adrenal, SCC, TCC Sebaceous Neoplasms
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|---------------------|---------|
| BSB-3772-3 | TintoFast Predilute | 3.0 ml |
| BSB-3772-7 | TintoFast Predilute | 7.0 ml |
| BSB-3772-15 | TintoFast Predilute | 15.0 ml |

TintoFast Androgen Receptor, MAb



IHC of Androgen Receptor on a Frozen Acetone-Fixed Melanoma Tissue


The androgen receptor (AR) is a type of nuclear receptor which is activated by binding of either of the androgenic hormones testosterone or dihydrotestosterone. The main function of the androgen receptor is as a DNA-binding transcription factor which regulates gene expression. However, the androgen receptor has additional functions independent of DNA binding. The AR signaling pathway plays a key role in development and function of male reproductive organs, including the prostate and epididymis. AR also plays a role in nonreproductive organs, such as muscle, hair follicles, and the brain.

This antibody reacts with the androgen receptor and also with the newly-described A form of the receptor. This antibody does not cross-react with estrogen, progesterone or glucocorticoid receptors. Abnormalities in the AR-signaling pathway have been linked to a number of diseases, including Prostate Cancer, Kennedy's Disease and male infertility.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-4
ISOTYPE: IgG1
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3690-3 | Tinto Predilute | 3.0 ml |
| BSB-3690-7 | Tinto Predilute | 7.0 ml |
| BSB-3690-15 | Tinto Predilute | 15.0 ml |

TintoFast CD31, MAb



IHC of CD31 on a Frozen Acetone-fixed Squamous Cell Carcinoma Tissue


CD31 is also called PECAM-1 for platelet endothelial cell-adhesion molecule. It plays a key role in removing aged neutrophils from the body. CD-31 is normally found on stem cells, endothelial cells, platelets, macrophages and Kupffer cells, granulocytes, T/NK cells, lymphocytes, megakaryocytes, fibroblasts, osteoclasts and neutrophils. CD-31 is also expressed in certain tumors, including Epithelioid Heman-gioendothelioma, Epithelioid Sarcoma like Hemangioendothelioma, other vascular tumors, Histiocytic malignancies, and Plasmacytomas. It is rarely found in some sarcomas and carcinomas. CD-31 and macrophages play a key role in tissue regeneration.

CD31 is widely used to identify the vascular origin of neoplasms, as it is a highly specific and sensitive marker for vascular endothelial cells.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 1A10
ISOTYPE: IgG1/K
CONTROL: Tonsil, Placenta, Appendix, Spleen, Kidney
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3691-3 | Tinto Predilute | 3.0 ml |
| BSB-3691-7 | Tinto Predilute | 7.0 ml |
| BSB-3691-15 | Tinto Predilute | 15.0 ml |

TintoFast CD34, MAb



IHC of CD34 on a FPPE Dermatofibrosarcoma Protuberans Tissue

CD34 functions as a cell-cell adhesion factor and cell-surface glycoprotein. It may also mediate the attachment of stem cells to bone marrow extracellular matrices or directly to stromal cells. Cells expressing CD34 are normally found in the umbilical cord and bone marrow as hematopoietic cells, and in vascular endothelium. In addition to stem cell recognition, CD34 is expressed by vascular endothelium; it appears that proliferating endothelial cells express this molecule in greater amounts than resting cells. In comparison to factor VIII R Antigen, CD34 stains are stronger and appear to be more sensitive in nature.

In tumors, CD34 is found in Alveolar Soft Part Sarcoma, pre B-ALL (positive in 75%), AML(40%), AMLM7 (most), Dermatofibrosarcoma Protuberans, Gastrointestinal Stromal Tumors, Giant Cell Fibroblastoma, Granulocytic Sarcoma, Kaposi's Sarcoma, Liposarcoma, Malignant Fibrous Histiocytoma, Malignant Peripheral Nerve Sheath tumors, Meningeal Hemangiopericytomas, Meningiomas, Neurofibromas, Schwannomas, and Papillary Thyroid Carcinoma. A negative CD34 may exclude Ewing's Sarcoma/PNET, Myofibrosarcoma of the breast, and Inflammatory Myofibroblastic tumors of the stomach.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: QBEnd/10


ISOTYPE: IgG1

CONTROL: Trichoepithelioma, Desmoplastic Trichilemmoma, Dermatofibrosarcoma Protuberans, Laposi's Sarcoma

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

TintoFast Chromogranin A, MAb



IHC of TintoFast Adipophilin on a FPPE Pancreas Tissue

Chromogranin A is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Examples of cells producing chromogranin A are the adrenal medulla, enterochromaffin-like cells and beta cells of the pancreas. The function of chromogranin A is unknown but it is a precursor to 3 functional peptides: vasostatin, pancreastatin and parastatin. These peptides negatively modulate the neuroendocrine function of the releasing cell (autocrine) or nearby cells (paracrine).

Chromogranin A is an excellent marker for Carcinoid Tumors, Pheochromocytomas, Paragangliomas, and other Neuroendocrine Tumors. Coexpression of Chromogranin A and Neuron-Specific Enolase (NSE) is common in Neuroendocrine Neoplasms. It has been identified in a wide variety of endocrine tissues including the Pituitary, Pancreas, Hypothalamus, Thymus, Thyroid, Intestine and Parathyroid. It is generally accepted that the coexpression of certain Keratins and Chromogranin means Neuroendocrine lineage. The presence of strong Chromogranin staining and absence of Keratin staining should raise the possibility of Paraganglioma. Most Pituitary Adenomas and Prolactinomas readily express Chromogranin.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: LK2H10


ISOTYPE: IgG1/K

CONTROL: Pancreas, Pituitary, Colon, Brain

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human

TintoFast CK 5/6, MAb



IHC of CK 5 & 6 on a Frozen Acetone-Fixed Melanoma Tissue

Cytokeratin 5 (58 kDa) is a high-molecular weight, basic type of cytokeratin expressed in basal, intermediate and superficial-cell layers of stratified epithelia as well as transitional epithelia, complex epithelia, mesothelial cells and Mesothelioma. Cytokeratin 6 (56 kD) is also a high-molecular weight, basic type cytokeratin expressed by proliferating squamous epithelium often paired with Cytokeratin 16.

CK 5 and 6 are positively seen in nearly 100% of Malignant Mesotheliomas and is rarely seen in Lung Adenocarcinomas. CK 5 and 6 can positively be seen in undifferentiated Large-cell Carcinoma as well as Squamous Carcinoma. Fewer than 10% of Carcinomas of the breast, colon, and prostate stain positively for this marker. CK 5 and 6 have also been used successfully as a myoepithelial cell marker in the prostate to determine malignancy.

ANTIBODY TYPE: Mouse Monoclonal

CLONE: D5/16B4

ISOTYPE: IgG1

CONTROL: SCC/BCC

LOCALIZATION: Cytoplasmic


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3681 | Tinto Predilute | 3.0 ml |
| BSB 3682 | Tinto Predilute | 7.0 ml |
| BSB 3683 | Tinto Predilute | 15.0 ml |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|---------------------|---------|
| BSB-3773-3 | TintoFast Predilute | 3.0 ml |
| BSB-3773-7 | TintoFast Predilute | 7.0 ml |
| BSB-3773-15 | TintoFast Predilute | 15.0 ml |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3666 | Tinto Predilute | 3.0 ml |
| BSB 3667 | Tinto Predilute | 7.0 ml |
| BSB 3668 | Tinto Predilute | 15.0 ml |

TintoFast CK7, MMab




IHC of Cytokeratin 7 on a Frozen Acetone-Fixed Extramammary Paget's Disease Tissue

Cytokeratin 7 (CK7) reacts with proteins that are found in most ductal, glandular and transitional epithelium of the urinary tract and bile duct epithelial cells. CK7 distinguishes between lung and breast epithelium that stain positive, and colon and prostate epithelial cells that are negative.

This antibody also reacts with many benign and malignant epithelial lesions (e.g., Adenocarcinomas of the ovary, breast and lung). Further, in frozen sections, the antibody has been shown to label the rete epithelium in the testis, epididymis epithelium, and the surface epithelium of the stomach and duodenum. Transitional Cell Carcinomas are positive and Prostate Cancers are negative. This antibody does not recognize intermediate filament proteins, nor does it recognize non-epithelial tissues such as blood vessels, connective tissue, etc.

TintoFast CK AE1/AE3, MMab




IHC of Cytokeratin AE1/AE3 on a Frozen Acetone-Fixed Basal Cell Carcinoma Tissue

Cytokeratins are intermediate-filament keratins found in the intracytoplasmic cytoskeleton of epithelial tissue. There are two types of cytokeratins: the low-weight, acidic Type I cytokeratins and the high-weight, basic or neutral Type II cytokeratins. Cytokeratins are usually found in pairs comprising a Type I cytokeratin and a Type II cytokeratin. Expression of these cytokeratins is frequently organ or tissue-specific.

Cytokeratin cocktail AE1/AE3 is well suited to distinguish Epithelial Carcinoma from Non-epithelial malignancies and is used to aid Epithelial Tumor classification. This antibody has been used to characterize the source of various neoplasms and to study the distribution of keratin-containing cells in epithelia during normal development and during the development of epithelial neoplasms. This antibody stains cytokeratins present in normal and abnormal human tissues. This antibody has shown high sensitivity and specificity in recognizing epithelial cells of neoplastic origin.

TintoFast EMA, MMab



IHC of EMA on a Frozen Acetone-fixed Melanoma Tissue

Epithelial Membrane Antigen (EMA) antibody is a mucin-like glycoprotein, shown to be useful as a pan-epithelial marker for detecting early metastatic loci of carcinoma in the bone marrow or liver. It stains normal and neoplastic cells from various tissues, including mammary epithelium, sweat glands and squamous epithelium. Hepatocellular Carcinoma, Adrenal Carcinoma and Embryonal Carcinomas are consistently EMA negative, so keratin positivity with negative EMA favors one of these tumors.

EMA is frequently positive in meningioma, which can be useful when distinguishing it from other intracranial neoplasms. The absence of EMA can also be of value since negative EMA is characteristic of some tumors including Adrenal Carcinoma, Seminomas, Paraganglioma and Hepatoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: OV-TL 12/30
ISOTYPE: IgG1/K
CONTROL: TCC/EMPD
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: AE1/AE3
ISOTYPE: IgG1
CONTROL: Skin, SCC, BCC
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Mouse Monoclonal
CLONE: E29
ISOTYPE: IgG2a/K
CONTROL: Breast, Skin, Colon, Kidney, Cervix
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3669 | Tinto Predilute | 3.0 ml |
| BSB 3670 | Tinto Predilute | 7.0 ml |
| BSB 3671 | Tinto Predilute | 15.0 ml |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3663 | Tinto Predilute | 3.0 ml |
| BSB 3664 | Tinto Predilute | 7.0 ml |
| BSB 3665 | Tinto Predilute | 15.0 ml |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3692-3 | Tinto Predilute | 3.0 ml |
| BSB-3692-7 | Tinto Predilute | 7.0 ml |
| BSB-3692-15 | Tinto Predilute | 15.0 ml |

TintoFast EpCAM BerEP4, MAb



IHC of Cytokeratin EpCAM BerEP4 on a Frozen Acetone-Fixed Extramammary Paget's Disease Tissue


Epithelial Cell Adhesion Molecule (EpCAM) or Epithelial Specific Antigen is a 40kD cell surface antigen that is broadly distributed in epithelial cells and displays a highly conserved expression in carcinomas. These glycoproteins are located on the cell membrane surface and in the cytoplasm of virtually all epithelial cells, with the exception of most squamous epithelia, hepatocytes, renal proximal tubular cells, gastric parietal cells and myoepithelial cells. However, focal positivity may be seen in the basal layer of squamous cell epithelium of endoderm (e.g., palatine tonsils) and mesoderm (e.g., uterine cervix).

EpCAM expression has been reported to be a possible marker of early malignancy, with expression being increased in tumor cells, and de novo expression being seen in dysplastic squamous epithelium. Epithelial specific antigen has been known to play an important role as a tumor-cell marker in lymph nodes from patients with esophageal carcinoma. EpCAM can be used to distinguish among Basal Cell, Basosquamous Carcinomas and Squamous Cell Carcinomas of the skin.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: Ber-EP4
ISOTYPE: IgG1/K
CONTROL: BCC, EMPD
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3678 | Tinto Predilute | 3.0 ml |
| BSB 3679 | Tinto Predilute | 7.0 ml |
| BSB 3680 | Tinto Predilute | 15.0 ml |

TintoFast HMB45, MAb



IHC of HMB45 on a Frozen Acetone-fixed Melanoma Tissue


HMB-45 reacts against an antigen present in immature melanosomes, cutaneous melanocytes, prenatal and infantile retinal pigment epithelium and melanoma cells. This antibody was generated to an extract of Melanoma. It reacted positively against Melanocytic Tumors but not other tumors, thus demonstrating specificity and sensitivity. Moreover, this antibody reacts positively against junctional nevus cells but not intradermal nevi, and against fetal melanocytes but not normal adult melanocytes.

This antibody is very useful to identify Malignant Melanoma. Metastatic Amelanotic Melanoma can often be confused with a variety of poorly differentiated Carcinomas, Large Cell Lymphomas, Sarcomas, Spindle Cell Carcinomas and various types of mesenchymal neoplasms. A keratin-negative, vimentin-rich neoplasm that immunoreacts with antibody to S-100 protein and with this melanoma antibody, is, with rare exception, a Melanoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: HMB-45
ISOTYPE: IgG1/K
CONTROL: Melanoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3693-3 | Tinto Predilute | 3.0 ml |
| BSB-3693-7 | Tinto Predilute | 7.0 ml |
| BSB-3693-15 | Tinto Predilute | 15.0 ml |

TintoFast Ki-67, RMAb



IHC of TintoFast Ki-67 on FPPE Breast Carcinoma Tissue


The Ki-67 protein is a cellular marker for proliferation. It is strictly associated with cell proliferation. During the interphase, the Ki-67 antigen can be exclusively detected within the cell nucleus, whereas in mitosis most of the protein is relocated to the surface of the chromosomes. Ki-67 protein is present during all active phases of the cell cycle (G1, S, G2, and mitosis), but is absent from resting cells (G0).

Ki-67 is an excellent marker to determine the growth fraction of a given cell population. The fraction of Ki-67-positive tumor cells (the Ki-67 labeling index) is often correlated with the clinical course of cancer. The best-studied examples in this context are carcinomas of the prostate and the breast.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RM360
ISOTYPE: IgG
CONTROL: Testis, Tonsil, Bone Marrow, Placenta, Colon, Tonsil, Fallopian Tube, Astrocytoma, Breast Carcinoma, Colon Carcinoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|---------------------|---------|
| BSB-3771-3 | TintoFast Predilute | 3.0 ml |
| BSB-3771-7 | TintoFast Predilute | 7.0 ml |
| BSB-3771-15 | TintoFast Predilute | 15.0 ml |

TintoFast MART-1, MMab



IHC of MART-1 on a Frozen Acetone-Fixed Melanoma Tissue


MART-1/Melan-A is a putative 18 kDa transmembrane protein consisting of 118 amino acids. It has a single transmembrane domain. MART-1/Melan-A is a protein antigen found on melanocytes. Antibodies against this antigen are used to recognize cells of melanocytic differentiation, useful for the diagnosis of Melanoma. The same name is used to refer to the gene which codes for this antigen.

The MART-1/Melan-A antigen is specific for the melanocyte lineage found in normal skin, retina, and melanocytes, but not in other normal tissues. It is thus useful as a marker for Melanocytic Tumors, with the caveat that it is normally found in benign nevi as well. This antibody is very useful in establishing the diagnosis of Metastatic Melanomas.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: A103
ISOTYPE: IgG1
CONTROL: Skin, Melanoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3651 | Tinto Predilute | 3.0 ml |
| BSB 3652 | Tinto Predilute | 7.0 ml |
| BSB 3653 | Tinto Predilute | 15.0 ml |

TintoFast Melanoma Cocktail: HMB-45 MART-1 & Tyrosinase, MMab



IHC of Melanoma Cocktail on a Frozen Acetone-fixed Lentigo Maligna Carcinoma Tissue


HMB-45 reacts against an antigen present in immature melanosomes, cutaneous, melanocytes, prenatal and infantile retinal pigment epithelium and melanoma cells. This antibody is very useful to identify Malignant Melanoma. MART-1/Melan-A is a protein antigen found on melanocytes. Antibodies against this antigen are used to recognize cells of melanocytic differentiation, useful for the diagnosis of Melanoma. The same name is used to refer to the gene which codes for this antigen. Tyrosinase is a copper-containing enzyme present in plant and animal tissues that catalyzes the production of melanin and other pigments from tyrosine by oxidation.

The MART-1/Melan-A antigen is specific for the melanocyte lineage found in normal skin, retina, and melanocytes, but not in other normal tissues. It is thus useful as a marker for Melanocytic Tumors, with the caveat that it is normally found in benign nevi as well. Anti-Tyrosinase has been found to be quite specific for melanotic lesions such as Malignant Melanoma and Melanotic Neurofibroma. Essentially no carcinomas express this marker. Melanoma cocktail HMB-45, Mart-1 and Tyrosinase are ideally suited to detect melanomas and melanocytic lesions and may prove to be a valuable marker for melanoma metastasis in sentinel lymph nodes.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: HMB-45, A103 & BSB-6
ISOTYPE: IgG1/K, IgG1 & IgG2a
CONTROL: Prostate, Prostatic Adenocarcinoma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3695-3 | Tinto Predilute | 3.0 ml |
| BSB-3695-7 | Tinto Predilute | 7.0 ml |
| BSB-3695-15 | Tinto Predilute | 15.0 ml |

TintoFast MiTF, MMab



IHC of MiTF on a Frozen Acetone-Fixed Melanoma Tissue


Microphthalmia-associated Transcription Factor (MiTF) is a basic helix-loop-helix leucine zipper transcription factor involved in melanocyte and osteoclast development. Mutations in MiTF cause auditory pigmentary syndromes, such as Waardenburg Syndrome Type II, Type IIa and Tietz Syndrome in humans. There are two known isoforms of MiTF differing by 66 amino acids at the NH2 terminus. Shorter forms are expressed in melanocytes and run as two bands at 52 kDa and 56 kDa, while the longer Mi form runs as a cluster of bands at 60-70 kDa in osteoclasts and in B16 Melanoma cells (but not other Melanoma cell lines), as well as mast cells and heart cells. MiTF plays a critical role in the differentiation of various cell types such as neural crest-derived melanocytes, mast cells, osteoclasts and optic cup-derived retinal pigment epithelium.

This antibody recognizes serine phosphorylated and non-phosphorylated melanocytic isoforms of microphthalmia. It is useful in identifying Malignant Melanoma, and distinguishing mast cell lesions from lesions of myeloid derivation. A relatively rare class of tumors known as PEComas (tumors showing perivascular epithelioid cell differentiation) express MiTF in a high percentage of cases (~90%).

ANTIBODY TYPE: Mouse Monoclonal
CLONE: C5/D5
ISOTYPE: IgG1/K
CONTROL: Skin, Melanoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3657 | Tinto Predilute | 3.0 ml |
| BSB 3658 | Tinto Predilute | 7.0 ml |
| BSB 3659 | Tinto Predilute | 15.0 ml |

TintoFast NGFR, MMab




IHC of NGFR on a Frozen Acetone-fixed Squamous Cell Carcinoma Tissue

NGFR (Nerve Growth Factor Receptor), also termed p75 or CD271, is the low-affinity NGFR (LNGFR) which binds NGF and other neurotrophins, including BDNF, NT3 and NT4/5 with similar low-affinity. NGFR p75 is a 75 kD transmembrane glycoprotein that is mainly expressed in Schwann cells and neurons and in a variety of non-neuronal cells. NGFR p75 is necessary for regulating neuronal growth, migration, differentiation and cell death during development of the central and peripheral nervous system. NGFR p75 plays a central role in the regulation of cell number by apoptosis in the developing CNS. During early development, activation of NGFR p75 by NGF induces apoptotic cell death in some neuronal cells, probably through activation of the sphingomyelinase/ceramide pathway, the ICE-like proteases and the JNK pathway. CD271 has recently been described as being expressed in mesenchymal stem cells (bone marrow stromal cells).

NGFR is expressed not only in sympathetic and sensory neurons, but also in various neural crest cell or tumor derivatives such as melanocytes, Melanomas, Neuroblastomas, Pheochromocytomas, Neurofibromas, and neurotized nevi (Type C melanocytes). It is now apparent that expression of NGFR is ubiquitous and not limited to the nervous system, being expressed in mature non-neural cells such as perivascular cells, dental pulp cells, lymphoid follicular dendritic cells, basal epithelium of oral mucosa and hair follicles, prostate basal cells and myoepithelial cells. Studies in Prostate and Urothelial Cancer suggest that NGFR may act as a tumor suppressor, negatively regulating cell growth and proliferation. NGFR labels the myoepithelial cells of breast ducts and intralobular fibroblasts of breast ducts and, thus, aids in the diagnosis of malignancy in the breast.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-18
ISOTYPE: IgG1
CONTROL: Brain, Breast, Prostate, Neuroblastoma, CNS Tumor
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

TintoFast p40, RMAb




IHC of p40 on a Frozen Acetone-fixed Squamous Cell Carcinoma Tissue

p40 is an antibody that recognizes ΔNp63-a p63 isoform and it is highly specific for squamous/basal cells. It may be a valuable marker in detecting Squamous Cell Carcinoma where p63 is currently used. It recognizes the shortest variant of p53. p40 is superior in specificity to p63 because it does not label lung adenocarcinomas like p63 does, which eliminates the potential of misinterpreting a positive adenocarcinoma as a squamous cell carcinoma.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: Zr8
ISOTYPE: IgG
CONTROL: Normal Prostate, Breast, Skin
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

TintoFast p63, MMab



IHC of p63 on a Frozen Acetone-fixed Squamous Cell Carcinoma Tissue

In addition to p53, mammalian cells contain two homologous genes, p63 and p73. These genes give rise to the expression of proteins that are highly similar to p53 in structure and function. In particular, p63 and p73 proteins can induce p53-responsive genes and elicit programmed cell death. p73 and p63 are important during development and differentiation. In particular, p63 appears to be primarily implicated in epithelial development.

Anti-p63 to human p63 protein labels an epitope common to all six p63 isoforms 5"QA 5"QC 5"Qc 3/QA 3/Qc 3/Qc QMBCFMT UIF OVDMFJPG myoepithelial cells in the prostate gland as well as breast tissue, making it useful in differentiating benign vs. malignant prostate lesions and breast lesions.


ANTIBODY TYPE: Mouse Monoclonal
CLONE: 4A4
ISOTYPE: IgG2a/K
CONTROL: Prostate, Breast, Skin, Salivary Gland
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat, Turtle

| CAT. # | PRESENTATION | VOL/QT |
|-------------|-----------------|---------|
| BSB-3696-3 | Tinto Predilute | 3.0 ml |
| BSB-3696-7 | Tinto Predilute | 7.0 ml |
| BSB-3696-15 | Tinto Predilute | 15.0 ml |

| CAT. # | PRESENTATION | VOL/QT |
|-------------|-----------------|---------|
| BSB-3697-3 | Tinto Predilute | 3.0 ml |
| BSB-3697-7 | Tinto Predilute | 7.0 ml |
| BSB-3697-15 | Tinto Predilute | 15.0 ml |

| CAT. # | PRESENTATION | VOL/QT |
|-------------|-----------------|---------|
| BSB-3698-3 | Tinto Predilute | 3.0 ml |
| BSB-3698-7 | Tinto Predilute | 7.0 ml |
| BSB-3698-15 | Tinto Predilute | 15.0 ml |

TintoFast Podoplanin, MAb



IHC of Podoplanin on a Frozen Acetone-fixed Squamous Cell Carcinoma Tissue


Podoplanin is a transmembrane mucoprotein (38 kDa) recognized by the D2-40 monoclonal antibody. Podoplanin is specifically expressed in the endothelium of lymphatic capillaries but not in the blood vasculature. In normal skin and kidney, podoplanin is co-localized with VEGFR3/FLT4, another marker for lymphatic endothelial cells.

Podoplanin is selectively expressed in lymphatic endothelium as well as Lymphangiomas, Kaposi's Sarcomas and in subset Angiosarcomas with probable lymphatic differentiation. Podoplanin has also been shown to be expressed in Epithelioid Mesotheliomas, Hemangioblastomas and Seminomas.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: D2-40
ISOTYPE: IgG1
CONTROL: Tonsil, Lymph Node, Lymphangioma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat, Mouse

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3699-3 | Tinto Predilute | 3.0 ml |
| BSB-3699-7 | Tinto Predilute | 7.0 ml |
| BSB-3699-15 | Tinto Predilute | 15.0 ml |

TintoFast Prame, RMAb



IHC of TintoFast PRAME on Frozen Melanoma Tissue

Melanoma antigen preferentially expressed in tumors is a protein that in humans is encoded by the PRAME gene. This gene encodes an antigen that is predominantly expressed in human Melanomas and is recognized by cytolytic T lymphocytes. This expression pattern is like that of other CT antigens, such as MAGE, BAGE and GAGE. However, unlike these other CT antigens, this antigen is also expressed in Acute Leukemias.


PRAME overexpression in triple negative breast cancer has also been found to promote cancer cell motility through induction of the epithelial-to-mesenchymal transition. PRAME mRNA expression is well documented in Cutaneous and Ocular Melanomas. One study concluded that diffuse nuclear immunoreactivity for PRAME was found in 87% of metastatic and 83.2% of primary Melanomas. Among Melanoma subtypes, PRAME was diffusely expressed in 94.4% of acral Melanomas, 92.5% of superficial spreading Melanomas, 90% of Nodular Melanomas, 88.6% of Lentigo Maligna Melanomas, and 35% of Desmoplastic Melanomas. When in situ and NonDesmoplastic Invasive Melanoma components were present, PRAME expression was seen in both. Most Melanocytic nevi (86.4%), were completely negative for PRAME.

A study suggests that immunohistochemical analysis for PRAME expression may be useful for diagnostic purposes to support a suspected diagnosis of Melanoma.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-PRAME
ISOTYPE: IgG
CONTROL: Testis, Seminoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|---------------------|---------|
| BSB-3774-3 | TintoFast Predilute | 3.0 ml |
| BSB-3774-7 | TintoFast Predilute | 7.0 ml |
| BSB-3774-15 | TintoFast Predilute | 15.0 ml |

TintoFast SOX-10, MAb



IHC of SOX-10 on a Frozen Acetone-Fixed Melanoma Tissue


Transcription factor SOX-10 is a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional activator after forming a protein complex with other proteins. This protein acts as a nucleocytoplasmic shuttle protein and is important for neural crest and peripheral nervous system development. Mutations in this gene are associated with Waardenburg-Shah and Waardenburg-Hirschsprung disease. Anti-SOX-10 has been recently shown to be a sensitive marker of melanoma, including conventional, spindle, and desmoplastic subtypes.

SOX-10 is expressed by metastatic melanomas and nodal capsular nevus in sentinel lymph nodes, but not by other lymph node components such as dendritic cells which usually express S100 protein. In scar specimens, immature fibroblasts, epithelioid granulomas, and histiocytic proliferations can histopathologically mimic residual melanoma and even be positive for MiTF and S100. However, SOX-10 is less likely to be expressed by fibroblasts or histiocytes, especially compared to MiTF and S100. Anti-SOX-10 produces a nuclear stain that provides a clean signal that is much sharper and darker in staining quality when compared to the use of antibodies against MiTF and S100.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-62
ISOTYPE: IgG2b/K
CONTROL: Skin, Melanoma
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3654 | Tinto Predilute | 3.0 ml |
| BSB 3655 | Tinto Predilute | 7.0 ml |
| BSB 3656 | Tinto Predilute | 15.0 ml |

TintoFast Synaptophysin, RPab



*IHC of TintoFast Synaptophysin on Frozen
 Merkel Cell Carcinoma Tissue*

Synaptophysin is a synaptic vesicle glycoprotein weighing 38 kDa. It is present in endocrine cells, the brain, spinal cord, and adrenal glands. It acts as a marker for neuroendocrine cells. Synaptophysin reacts with Neuroendocrine cells of human Adrenal Medulla, Carotid Body, Skin, Pituitary, Thyroid, Lung, Pancreas and Gastrointestinal mucosa. Positive staining is seen in neurons of the Brain, Spinal Cord, Retina, and Paneth's cells in the Gastrointestinal tract and Gastric parietal cells.

This antibody identifies normal Neuroendocrine cells and Neuroendocrine neoplasms. Diffuse, finely-granular cytoplasmic staining is observed and probably correlates with the distribution of the antigen within Neurosecretory vesicles. The expression of Synaptophysin is independent of the presence of NSE or other Neuroendocrine markers. Synaptophysin is an independent broad-range marker of neural and neuroendocrine differentiation.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Pancreas, Brain, Pituitary, Adrenal, Colon

LOCALIZATION: Cytoplasmic

SPECIES REACTIVITY: Human, Dog, Cat

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|---------------------|---------|
| BSB-3775-3 | TintoFast Predilute | 3.0 ml |
| BSB-3775-7 | TintoFast Predilute | 7.0 ml |
| BSB-3775-15 | TintoFast Predilute | 15.0 ml |




Fast Mohs PolyDetector HRP Detection Systems


The Bio SB Mohs PolyDetector is a non-Biotin, Fab Micropolymer detection system that allows for the demonstration of antigens in FFPE tissues, Mohs cryostat sections, blood smears, cytospreads, and cell preparations. The Bio SB Mohs PolyDetector technology was developed with our proprietary Micropolymer backbone, conjugated to Anti-Mouse and Anti-Rabbit Fab fraction of IgG's, plus high quality HRP. The elimination of the Anti-Mouse and Anti-Rabbit immunoglobulin Fc region reduces non-specific reactions and background. This ensures consistent and reproducible immunostaining for all types of nuclear, cytoplasmic and membranous antigens, in different types of tissues. The increased sensitivity of this system allows for rapid staining procedures without compromising the stain quality.

Fast Mohs Fab PolyDetector HRP Overview

- Non-Biotin, Fab Micropolymer HRP 1-Step immunohistochemistry detection technology.
- Fab Micropolymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal.
- Ready-to-use, high sensitivity system especially designed for immunohistochemistry of frozen tissues.
- Universal: detects mouse or rabbit antibodies.
- Developed as DAB HRP Brown and HRP Green detection kits, but can be used with the AEC HRP Red, HRP Blue and HRP Black substrate-chromogens.
- Designed to be used manually, or can be optimized to work with any open, automated system.
- For in vitro diagnostic use. All kits manufactured according to US FDA and ISO 13485 guidelines.



Fast Mohs Fab PolyDetector HRP Protocol



Mohs PolyDetector DAB HRP Brown

| Product Description | Volume | Catalog # |
|--|----------|-----------|
| Mohs Mouse/Rabbit PolyDetector DAB HRP Brown | 5.0 ml | BSB 0307S |
| Mohs Mouse/Rabbit PolyDetector DAB HRP Brown | 15.0 ml | BSB 0307 |
| Mohs Mouse/Rabbit PolyDetector DAB HRP Brown | 50.0 ml | BSB 0308 |
| Mohs Mouse/Rabbit PolyDetector DAB HRP Brown | 100.0 ml | BSB 0309 |

Mohs PolyDetector HRP Green


| Product Description | Volume | Catalog # |
|--|----------|-----------|
| Mohs Mouse/Rabbit PolyDetector HRP Green | 5.0 ml | BSB 0310S |
| Mohs Mouse/Rabbit PolyDetector HRP Green | 15.0 ml | BSB 0310 |
| Mohs Mouse/Rabbit PolyDetector HRP Green | 50.0 ml | BSB 0311 |
| Mohs Mouse/Rabbit PolyDetector HRP Green | 100.0 ml | BSB 0312 |

Fast Mohs PolyDetector Plus HRP Detection Systems


The Bio SB Mohs PolyDetector Plus is a highly sensitive non-Biotin, Fab Micropolymer detection system that allows for the demonstration of antigens in FFPE tissues, Mohs cryostat sections, blood smears, cytospins, and cell preparations. The Bio SB Mohs PolyDetector Plus technology was developed with our proprietary Micropolymer backbone, conjugated to Anti-Link Fab fraction of IgG's, plus high quality HRP. The elimination of the Anti-Link Fc region reduces non-specific reactions and background. The Mohs PolyDetector Plus kit incorporates an immunoglobulin link and a Fab micropolymer label. This multiple component Fab Micropolymer delivers a highly sensitive and more specific signal than the Mohs PolyDetector kit. It is recommended for easy to detect antigens and more difficult to detect antigens like nuclear targets.

Fast Mohs Fab PolyDetector Plus HRP Overview

- Non-Biotin, 2-Step immunohistochemistry detection technology.
- Recommended for easy to detect antigens and more difficult to detect antigens like nuclear targets.
- Anti-mouse and rabbit Fab Micropolymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal.
- Ready-to-use, high sensitivity system especially designed for immunohistochemistry of frozen tissues.
- Developed as DAB HRP Brown and HRP Green detection kits, but can be used with the AEC HRP Red, HRP Blue and HRP Black substrate-chromogens.
- For in vitro diagnostic use. All kits manufactured according to US FDA and ISO 13485 guidelines.



Fast Mohs Fab PolyDetector Plus HRP Protocol



Mohs PolyDetector Plus DAB HRP Brown

| Product Description | Volume | Catalog # |
|---|----------|--------------|
| Mohs Mouse/Rabbit PolyDetector Plus DAB HRP Brown | 5.0 ml | BSB-0355-5 |
| Mohs Mouse/Rabbit PolyDetector Plus DAB HRP Brown | 15.0 ml | BSB-0355-15 |
| Mohs Mouse/Rabbit PolyDetector Plus DAB HRP Brown | 50.0 ml | BSB-0355-50 |
| Mohs Mouse/Rabbit PolyDetector Plus DAB HRP Brown | 100.0 ml | BSB-0355-100 |

Mohs PolyDetector Plus HRP Green

| Product Description | Volume | Catalog # |
|---|----------|--------------|
| Mohs Mouse/Rabbit PolyDetector Plus HRP Green | 5.0 ml | BSB-0359-5 |
| Mohs Mouse/Rabbit PolyDetector Plus HRP Green | 15.0 ml | BSB-0359-15 |
| Mohs Mouse/Rabbit PolyDetector Plus HRP Green | 50.0 ml | BSB-0359-50 |
| Mohs Mouse/Rabbit PolyDetector Plus HRP Green | 100.0 ml | BSB-0359-100 |


Ancillaries For Fast Mohs IHC

Fast ChromoProtector

| Product Description | Volume | Catalog # |
|----------------------|---------|---------------|
| Fast ChromoProtector | 15 ml | BSB-0327-15 |
| Fast ChromoProtector | 50 ml | BSB-0327-50 |
| Fast ChromoProtector | 100 ml | BSB-0327-100 |
| Fast ChromoProtector | 200 ml | BSB-0327-200 |
| Fast ChromoProtector | 500 ml | BSB-0327-500 |
| Fast ChromoProtector | 1000 ml | BSB-0327-1000 |


Mohs ImmunoDigester

| Product Description | Volume | Catalog # |
|---------------------|----------|-----------|
| Mohs ImmunoDigester | 15.0 ml | BSB 0324 |
| Mohs ImmunoDigester | 50.0 ml | BSB 0325 |
| Mohs ImmunoDigester | 100.0 ml | BSB 0326 |



IVD

IHC of EpCAM BerEP4 on a Frozen Acetone-fixed Basal Cell Carcinoma Tissue using the Mohs PolyDetector DAB HRP Brown



IVD


IHC of MART-1 on a Frozen Acetone-fixed Melanoma Tissue using the Mohs PolyDetector HRP Green



Immunofluorescence for Autoimmunity

Quality Antibodies | Sensitive Detection | IF and IHC Applications

A selection of antibodies, detection systems for IF and IHC Autoimmune disease applications.



Direct Immunofluorescence

Direct immunofluorescence uses a single, primary antibody, chemically linked to a fluorophore or fluorochrome. The primary antibody recognizes the target molecule (antigen) and binds to a specific region called the epitope. The attached fluorophore can be detected via fluorescent microscopy. Bio SB's Direct immunofluorescence, although somewhat less common, has notable advantages over the indirect procedure. The direct detection of the antigen by a primary antibody directly labeled with a fluorochrome reduces the number of steps in the procedure, saving time and reducing non-specific background signal.

All Direct Immunofluorescence Antibodies are available in concentrate and convenient Tinto Predilute formats to meet your laboratory needs.



Tinto Prediluted Antibodies

3 mL Tinto Predilute (30 tests)
7 mL Tinto Predilute (70 tests)
15 mL Tinto Predilute (150 tests)


Concentrated Antibodies

0.5 mL Concentrate
1 mL Concentrate

- High affinity purified antibodies conjugated to FITC.
- Highly sensitive and specific.
- Low background.
- Reduction in time compared to Indirect labeling.
- Validated for transplantation and autoimmunity studies.




IF of IgA on a Frozen Kidney Tissue



IF of Fibrinogen on a Frozen Kidney Tissue

Albumin / FITC, RPaB



IF of Albumin on a FFPE Liver Tissue

The albumins are a family of globular proteins, the most common of which are the serum albumins. Albumins are commonly found in blood plasma and differ from other blood proteins in that they are not glycosylated. Albumin functions primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume. Mutations in this gene on chromosome 4 result in various anomalous proteins.


Low albumin (hypoalbuminemia) may be caused by liver disease, nephrotic syndrome, burns, protein-losing enteropathy, malabsorption, malnutrition, late pregnancy, artefact, genetic variations and malignancy. High albumin (hyperalbuminemia) is almost always caused by dehydration. In some cases of retinol (Vitamin A) deficiency, the albumin level can be elevated to high-normal values.

It has been reported in systemic lupus erythematosus (SLE) patients an increased prevalence of IgG autoantibodies against human serum albumin (anti-HSA IgG) that are associated with SLE disease activity.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Salivary Gland, Kidney, Tonsil, Lupus Erythematosus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3000-3 | Tinto Predilute | 3.0 ml |
| BSB-3000-7 | Tinto Predilute | 7.0 ml |
| BSB-3000-15 | Tinto Predilute | 15.0 ml |
| BSB-3000-05 | Concentrate | 0.5 ml |
| BSB-3000-1 | Concentrate | 1.0 ml |
| BSB-3000-CS | Control Slides | 5 |

C1q / FITC, RPaB



IF of C1q on a FFPE Lupus Erythematosus Tissue


The complement component 1q (C1q) is a protein complex involved in the complement system, which is part of the innate immune system. C1q together with C1r and C1s form the C1 complex. Antibodies of the adaptive immune system can bind antigen, forming an antigen-antibody complex. When C1q binds antigen antibody complexes, the C1 complex becomes activated. Activation of the C1 complex initiates the classical complement pathway of the complement system.

C1q nephropathy is a rare glomerular disease with characteristic mesangial C1q deposition noted on IHC or IF microscopy. It is histologically defined and poorly understood. Light microscopic features are heterogeneous and comprise minimal change disease (MCD), focal segmental glomerulosclerosis (FSGS), and proliferative glomerulonephritis. Clinical presentation is also diverse, and ranges from asymptomatic hematuria or proteinuria to frank nephritic or nephrotic syndrome in both children and adults. Hypertension and renal insufficiency at the time of diagnosis are common findings. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Kidney, Cervix, Spleen, Lupus Erythematosus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3001-3 | Tinto Predilute | 3.0 ml |
| BSB-3001-7 | Tinto Predilute | 7.0 ml |
| BSB-3001-15 | Tinto Predilute | 15.0 ml |
| BSB-3001-05 | Concentrate | 0.5 ml |
| BSB-3001-1 | Concentrate | 1.0 ml |
| BSB-3001-CS | Control Slides | 5 |

C3c / FITC, RPaB



IF of C3c on a FFPE Lupus Erythematosus Tissue


Complement component 3, often simply called C3, is a protein of the immune system. It plays a central role in the complement system and contributes to innate immunity.

C3 glomerulopathy was recently coined to describe renal biopsy appearances characterized by the presence of glomerular deposits composed predominantly of C3 in the absence of significant amounts of Ig. The presence of C3 in the absence of Ig suggests activation of complement by antibody-independent pathways, typically the alternative pathway, and many patients with this type of renal lesion have evidence of genetic or acquired alternative pathway dysregulation. C3 glomerulopathy has been further divided into dense deposit disease (DDD) and C3 glomerulonephritis (C3GN) based on electron microscopy (EM) appearances. The underlying genetic defect has been identified in some hereditary forms of C3GN such as CFHR5 nephropathy. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Placenta, Kidney, Fallopian Tube, Lupus Erythematosus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3002-3 | Tinto Predilute | 3.0 ml |
| BSB-3002-7 | Tinto Predilute | 7.0 ml |
| BSB-3002-15 | Tinto Predilute | 15.0 ml |
| BSB-3002-05 | Concentrate | 0.5 ml |
| BSB-3002-1 | Concentrate | 1.0 ml |
| BSB-3002-CS | Control Slides | 5 |

C4c / FITC, RPab




IF of C4c on a Lupus Erythematosus Tissue

Fibrinogen / FITC, RPab



IF of Fibrinogen on a Frozen Kidney Tissue

IgA / FITC, RPab



IF of IgA on Frozen Kidney Tissue

Complement component 4 (C4), in humans, is a protein involved in the intricate complement system, originating from the human leukocyte antigen (HLA) system. It serves several critical functions in immunity, tolerance, and autoimmunity with the other numerous components. Low serum complement activity or low protein concentrations of complement C4 are found on Systemic Lupus Erythematosus (SLE) and it is often associated with Congenital C4 deficiency. Complete deficiencies of complement components are among the strongest genetic risk factors for SLE or lupus-like disease, across HLA haplotypes and racial backgrounds.

Fibrinogen (factor I) is a glycoprotein that circulates in the blood of vertebrates. During tissue and vascular injury, it is converted enzymatically by thrombin to fibrin and subsequently to a fibrin-based blood clot. Fibrinogen functions primarily to occlude blood vessels and thereby stop excessive bleeding. Fibrin also mediates blood platelet and endothelial cell spreading, tissue fibroblast proliferation, capillary tube formation, and angiogenesis and thereby functions to promote tissue revascularization, wound healing, and tissue repair. Several disorders (Congenital afibrinogenemia, hypofibrinogenemia, Fibrinogen storage disease, hereditary fibrinogen Aa-Chain amyloidosis, congenital hypodysfibrinogenemia, Cryofibrinogenemia, acquired hypofibrinogenemia, Chronic Kidney Disease, etc.) in the quantity and/or quality of fibrinogen cause pathological bleeding, pathological blood clotting, and/or the deposition of fibrinogen in the liver, kidneys, and other tissues. Chronic kidney disease (CKD) patients have increased rates of bleeding as well as thrombosis. Fibrinogen and platelets combine to generate a mature clot, but in CKD platelets are dysfunctional.

Immunoglobulin A (IgA) is the main immunoglobulin in mucous secretions, including tears, saliva, and colostrum, as well as respiratory, intestinal, prostatic, and vaginal secretions. It is also found in small amounts in blood. Because it is resistant to degradation by enzymes, secretory IgA provides protection against microbes proliferating in body secretions, especially those of the digestive and respiratory tracts. IgA antibody reacts with surface immunoglobulin IgA alpha chains. It is extremely useful when identifying Acute Leukemias, IgA Myelomas, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. However, due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Testis, Kidney, Pancreas, Salivary Gland, Colon
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Breast, Testis, Kidney, Pancreas, Salivary Gland, Skin, Fallopian Tube
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Spleen, Lymph Node, Kidney, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3003-3 | Tinto Predilute | 3.0 ml |
| BSB-3003-7 | Tinto Predilute | 7.0 ml |
| BSB-3003-15 | Tinto Predilute | 15.0 ml |
| BSB-3003-05 | Concentrate | 0.5 ml |
| BSB-3003-1 | Concentrate | 1.0 ml |
| BSB-3003-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3004-3 | Tinto Predilute | 3.0 ml |
| BSB-3004-7 | Tinto Predilute | 7.0 ml |
| BSB-3004-15 | Tinto Predilute | 15.0 ml |
| BSB-3004-05 | Concentrate | 0.5 ml |
| BSB-3004-1 | Concentrate | 1.0 ml |
| BSB-3004-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3005-3 | Tinto Predilute | 3.0 ml |
| BSB-3005-7 | Tinto Predilute | 7.0 ml |
| BSB-3005-15 | Tinto Predilute | 15.0 ml |
| BSB-3005-05 | Concentrate | 0.5 ml |
| BSB-3005-1 | Concentrate | 1.0 ml |
| BSB-3005-CS | Control Slides | 5 |

IgD / FITC, RPab




IF of IgD on a FFPE Tonsil Tissue

IgE / FITC, RPab



IF of IgE on a FFPE Lupus Positive Tissue

IgG / FITC, RPab



IF of C3c on a FFPE Lupus Erythematosus Tissue

IgD makes up about 1% of proteins in the plasma membranes of immature B lymphocytes (coexpressed with IgM) and is also found in serum in very small amounts. It is monomeric and incorporates the alpha-heavy chain in its structure. IgD's function is currently unknown, as mice lacking IgD seem to retain normal immune responses (implying redundancy if not lack of function), and IgD ceases to be expressed in activated B-lymphocytes. It may function as a regulatory antigen receptor. IgD is the major antigen receptor isotype co-expressed with IgM on the surface of most peripheral B cells in mice and humans.

The IgD antibody reacts with surface immunoglobulin IgD delta chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived from Lymphomas, specifically Marginal Zone Lymphoma. Renal involvement in systemic lupus erythematosus (SLE) is associated with production of antibodies to double stranded DNA, deposition of immune complexes and organ damage. Lupus nephritis patients were characterized by increased percentage of immature/early-transitional B cells (CD27-IgD+CD21-), higher frequency of activated switched memory (SM, CD27+IgD-CD21-) and exhausted memory B-cells (CD27-IgD-), and decrease in nonswitched memory (NSM, CD27+IgD+) B-cells.

IgE, Immunoglobulin E, is an isotype of antibody only found in mammals. IgE synthesized by plasma cells. Monomers of IgE consist of two heavy chains (e chain) and two light chains, with the e chain containing 4 Ig-like constant domains (Ce1-Ce4). IgE's main function is immunity to parasites such as helminths like *Schistosoma mansoni*, *Trichinella spiralis*, and *Fasciola hepatica*. IgE is utilized during immune defense against certain protozoan parasites such as *Plasmodium falciparum*.

IgE also has an essential role in type I hypersensitivity, which manifests in various allergic diseases, such as allergic asthma, most types of sinusitis, allergic rhinitis, food allergies, and specific types of chronic urticaria and atopic dermatitis. IgE also plays a pivotal role in responses to allergens, such as: anaphylactic drugs, bee stings, and antigen preparations used in desensitization immunotherapy. IgE is known to be elevated in various autoimmune disorders such as Lupus (SLE), Rheumatoid Arthritis (RA) & psoriasis, and is theorized to be of pathogenetic importance in RA and SLE by eliciting a hypersensitivity reaction.

IgG is a monomeric immunoglobulin, comprised of two heavy chains and two light chains. This is the most abundant immunoglobulin and is approximately equally distributed in blood and tissue liquids, constituting 75% of serum immunoglobulins in humans. This is the only isotype that can pass through the placenta and bind to many kinds of pathogens. IgG protects the body against them by complement activation (classic pathway), opsonization for phagocytosis and neutralization of their toxins. There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%).

IgG antibody reacts with surface immunoglobulin IgG gamma chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q. Clinically, hematuria and proteinuria are present, with or without nephrotic syndromes. Mesangial IgG glomerulonephritis has been recently recognized as a distinct type of glomerulonephritis. The morphologic criteria detected in these patients included mesangial dense deposits by ultrastructural studies, which were predominantly positive for IgG by immunofluorescence.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Thymus, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Kidney, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3006-3 | Tinto Predilute | 3.0 ml |
| BSB-3006-7 | Tinto Predilute | 7.0 ml |
| BSB-3006-15 | Tinto Predilute | 15.0 ml |
| BSB-3006-05 | Concentrate | 0.5 ml |
| BSB-3006-1 | Concentrate | 1.0 ml |
| BSB-3006-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3011-3 | Tinto Predilute | 3.0 ml |
| BSB-3011-7 | Tinto Predilute | 7.0 ml |
| BSB-3011-15 | Tinto Predilute | 15.0 ml |
| BSB-3011-05 | Concentrate | 0.5 ml |
| BSB-3011-1 | Concentrate | 1.0 ml |
| BSB-3011-CS | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3007-3 | Tinto Predilute | 3.0 ml |
| BSB-3007-7 | Tinto Predilute | 7.0 ml |
| BSB-3007-15 | Tinto Predilute | 15.0 ml |
| BSB-3007-05 | Concentrate | 0.5 ml |
| BSB-3007-1 | Concentrate | 1.0 ml |
| BSB-3007-CS | Control Slides | 5 |

IgM / FITC, RPab




IF of IgM on a FFPE Colon Tissue

Kappa / FITC, RPab



IF of Kappa on a FFPE Tonsil Tissue

Lambda / FITC, RPab



IF of Lambda on a FFPE Tonsil Tissue

IgM forms polymers where multiple immunoglobulins are covalently linked together with disulfide bonds, normally as a pentamer or occasionally as a hexamer. It has a large molecular mass of approximately 900 kDa (in its pentamer form). In germline cells, the gene segment encoding the constant region of the heavy chain is positioned first among other constant region gene segments. For this reason, IgM is the first immunoglobulin expressed by mature B-cells. IgM antibody reacts with surface immunoglobulin IgM mu chains. IgM is one of the predominant surface immunoglobulins on B-lymphocytes, and is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q. Clinically, hematuria and proteinuria are present, with or without nephrotic syndromes. Immunoglobulin M (IgM) nephropathy is an uncommon glomerular disease characterized by IgM deposits in the mesangium.

Kappa detects surface immunoglobulin on normal and neoplastic B-cells. In paraffin embedded tissue, Kappa exhibits strong staining of kappa-positive plasma cells and cells that have absorbed exogenous immunoglobulin. When studying B-cell neoplasms, the determination of light-chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either Kappa or Lambda light chains, whereas reactive proliferations display a mixture of Kappa and Lambda-positive cells. If only a single light-chain type is detected, a lympho-proliferative disorder is very likely. Monoclonality is determined by a Kappa-Lambda ratio greater than or equal to 3:1, a Lambda-Kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population. In IgG-dominant immune complex-mediated glomerulonephritis, there are multiple pathological findings that strongly suggest the diagnosis of Lupus Nephritis including immunofluorescence staining for IgG, IgM, IgA, Kappa or Lambda, C3 and C1.

Lambda detects surface immunoglobulin on normal and neoplastic B-cells. Lambda staining is seen in B-cell follicles of human lymphoid tissue. When studying B-cell neoplasms, the determination of light chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population. In IgG-dominant immune complex-mediated glomerulonephritis, there are multiple pathological findings that strongly suggest the diagnosis of Lupus Nephritis including immunofluorescence staining for IgG, IgM, IgA, Kappa or Lambda, C3 and C1.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3008-3 | Tinto Predilute | 3.0 ml |
| BSB-3008-7 | Tinto Predilute | 7.0 ml |
| BSB-3008-15 | Tinto Predilute | 15.0 ml |
| BSB-3008-05 | Concentrate | 0.5 ml |
| BSB-3008-1 | Concentrate | 1.0 ml |
| BSB-3008-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3009-3 | Tinto Predilute | 3.0 ml |
| BSB-3009-7 | Tinto Predilute | 7.0 ml |
| BSB-3009-15 | Tinto Predilute | 15.0 ml |
| BSB-3009-05 | Concentrate | 0.5 ml |
| BSB-3009-1 | Concentrate | 1.0 ml |
| BSB-3009-CS | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3010-3 | Tinto Predilute | 3.0 ml |
| BSB-3010-7 | Tinto Predilute | 7.0 ml |
| BSB-3010-15 | Tinto Predilute | 15.0 ml |
| BSB-3010-05 | Concentrate | 0.5 ml |
| BSB-3010-1 | Concentrate | 1.0 ml |
| BSB-3010-CS | Control Slides | 5 |

Immunofluorescence Antibodies for Autoimmunity

Bio SB's high affinity, high sensitivity InDirect Immunofluorescence antibodies work with both our AmpliDetector Plus InDirect FITC Detection system or our PolyDetector micropolymer IHC detection system. This allows for greater flexibility without sacrificing price or staining quality.

All InDirect Immunofluorescence Antibodies are available in concentrate and convenient Tinto Predilute formats to meet your laboratory needs.

Tinto Prediluted Antibodies


3 mL Tinto Predilute (30 tests)
7 mL Tinto Predilute (70 tests)
15 mL Tinto Predilute (150 tests)

Concentrated Antibodies

0.1 mL Concentrate
0.5 mL Concentrate
1 mL Concentrate


| InDirect Immunofluorescence Antibodies | |
|--|----------------------|
| Albumin, RPab | Polyclonal |
| Amyloid A, RMab | EP335 |
| C1q, RPab | Polyclonal |
| C3c, RPab | Polyclonal |
| C3d, RPab | Polyclonal |
| C4c, RPab | Polyclonal |
| C4d, RMab | EP272 |
| Collagen IV, RMab | RBT-COL4 |
| Collagen IV, MMab | CIV22 |
| Cytokeratin 5/6, RMab | EP24 & EP67 |
| Cytomegalovirus, MMab | 8B1.2, 1G5.2 & 2D4.2 |
| Fibrinogen, RPab | Polyclonal |
| IgA, MMab | BSB-39 |
| IgA, RPab | Polyclonal |
| IgD, RMab | EP173 |
| IgD, RPab | Polyclonal |

| InDirect Immunofluorescence Antibodies | |
|--|------------|
| IgE, RPab | Polyclonal |
| IgG, MMab | BSB-40 |
| IgG, RPab | Polyclonal |
| IgG4, MMab | BSB-96 |
| IgG4, RMab | EP138 |
| IgM, RPab | Polyclonal |
| IgM, MMab | BSB-41 |
| Kappa, RPab | Polyclonal |
| Kappa, MMab | BSB-58 |
| Lambda, MMab | BSB-16 |
| Lambda, RPab | Polyclonal |
| p40, RMab | ZR8 |
| p63, MMab | 4A4 |
| p63, RMab | EP174 |
| PLA2R1, MMab | BSB-129 |
| SV40, MMab | pab101 |



C1q


IHC on a FFPE Lupus Erythmatosus Tissue
IF on a Frozen Lupus Erythmatosus Tissue



Fibrinogen


IHC on a FFPE Kidney Tissue
IF on a Frozen Lichen Planus Tissue

Albumin, RPA




IHC and IF of Albumin on a FFPE Lupus Erythematosus (IHC) and a FFPE Kidney Tissue (IF)

Amyloid A, RMA



IHC and IF of Amyloid A on a FFPE Kidney Tissue

C1q, RPA



IHC and IF of C1q on a FFPE Lupus Erythematosus (IHC) and a Frozen Lupus Erythematosus Tissue (IF)

The albumins are a family of globular proteins, the most common of which are the serum albumins. Albumins are commonly found in blood plasma and differ from other blood proteins in that they are not glycosylated. Albumin functions primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume. Mutations in this gene on chromosome 4 result in various anomalous proteins.

Low albumin (hypoalbuminemia) may be caused by liver disease, nephrotic syndrome, burns, protein-losing enteropathy, malabsorption, malnutrition, late pregnancy, artefact, genetic variations and malignancy. High albumin (hyperalbuminemia) is almost always caused by dehydration. In some cases of retinol (Vitamin A) deficiency, the albumin level can be elevated to high-normal values.

It has been reported in systemic lupus erythematosus (SLE) patients an increased prevalence of IgG autoantibodies against human serum albumin (anti-HSA IgG) that are associated with SLE disease activity.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Salivary Gland, Kidney, Tonsil, Lupus Erythematosus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

Serum amyloid A (SAA) proteins are a family of lipoproteins associated with high-density lipoprotein (HDL) in plasma. Different isoforms of SAA are expressed constitutively (constitutive SAAs) at different levels or in response to inflammatory stimuli (acute phase SAAs). These proteins are produced predominantly by the liver. The conservation of these proteins throughout invertebrates and vertebrates suggests that SAAs play a highly essential role in all animals. Acute-phase serum amyloid A proteins (A-SAAs) are secreted during the acute phase of inflammation. A-SAAs are implicated in several chronic inflammatory diseases, such as amyloidosis, atherosclerosis, and rheumatoid arthritis.

Amyloidosis is a disease characterized by the abnormal build-up of amyloid, abnormal non-branching fibrillary β -pleated sheet proteins that are insoluble and highly resistant to proteolytic degradation that result in localized or systemic organ dysfunction. AA amyloidosis is associated with a variety of chronic inflammatory conditions and infections, derived from SAA. Immunohistochemical staining using a panel of antibodies including κ and λ Ig light chains, Amyloid A, and Transthyretin can aid in recognizing most forms of amyloid. The Amyloid A immunostaining detects tissue deposition of serum Amyloid A protein, an acute phase reactive protein. It is positive in AA Amyloidosis and familial Mediterranean fever. SAA concentrations have been reported to be a marker of poor prognosis, elevated in patients with advanced stages of cancer and those with malignant disease.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP335
ISOTYPE: IgG
CONTROL: Kidney, Amyloidosis
LOCALIZATION: Cytoplasmic, Extracellular
SPECIES REACTIVITY: Human

The complement component 1q (C1q) is a protein involved in the complement system, which is part of the innate immune system. C1q together with C1r and C1s form the C1 complex. Antibodies of the adaptive immune system can bind antigen, forming an antigen-antibody complex.

When C1q binds antigen-antibody complexes, the C1 complex becomes activated. Activation of the C1 complex initiates the classical complement pathway of the complement system.

C1q nephropathy is a rare glomerular disease with characteristic mesangial C1q deposition noted on IHC or IF microscopy. It is histologically defined and poorly understood. Light microscopic features are heterogeneous and comprise minimal change disease (MCD), focal segmental glomerulosclerosis (FSGS), and proliferative glomerulonephritis. Clinical presentation is also diverse, and ranges from asymptomatic hematuria or proteinuria to frank nephritic or nephrotic syndrome in both children and adults. Hypertension and renal insufficiency at the time of diagnosis are common findings. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Kidney, Cervix, Spleen, Lupus Erythematosus
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3012 | Tinto Predilute | 3.0 ml |
| BSB 3013 | Tinto Predilute | 7.0 ml |
| BSB 3014 | Tinto Predilute | 15.0 ml |
| BSB 3015 | Concentrate | 0.1 ml |
| BSB 3016 | Concentrate | 0.5 ml |
| BSB 3017 | Concentrate | 1.0 ml |
| BSB 3018 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2803 | Tinto Predilute | 3.0 ml |
| BSB 2804 | Tinto Predilute | 7.0 ml |
| BSB 2805 | Tinto Predilute | 15.0 ml |
| BSB 2806 | Concentrate | 0.1 ml |
| BSB 2807 | Concentrate | 0.5 ml |
| BSB 2808 | Concentrate | 1.0 ml |
| BSB 2809 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3019 | Tinto Predilute | 3.0 ml |
| BSB 3020 | Tinto Predilute | 7.0 ml |
| BSB 3021 | Tinto Predilute | 15.0 ml |
| BSB 3022 | Concentrate | 0.1 ml |
| BSB 3023 | Concentrate | 0.5 ml |
| BSB 3024 | Concentrate | 1.0 ml |
| BSB 3025 | Control Slides | 5 |

C3c, RPab



IHC and IF of C3c on a FFPE Lupus Erythematosus (IHC) and Frozen Lupus Erythematosus Tissue (IF)

Complement component 3, often simply called C3, is a protein of the immune system. It plays a central role in the complement system and contributes to innate immunity.

C3 glomerulopathy was recently coined to describe renal biopsy appearances characterized by the presence of glomerular deposits composed predominantly of C3 in the absence of significant amounts of Ig. The presence of C3 in the absence of Ig suggests activation of complement by antibody-independent pathways, typically the alternative pathway, and many patients with this type of renal lesion have evidence of genetic or acquired alternative pathway dysregulation. C3 glomerulopathy has been further divided into dense deposit disease (DDD) and C3 glomerulonephritis (C3GN) based on electron microscopy (EM) appearances. The underlying genetic defect has been identified in some hereditary forms of C3GN such as CFHR5 nephropathy. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Placenta, Kidney, Fallopian Tube, Lupus Erythematosus

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

C3d, RPab



IHC and IF of C3d on a FFPE Rejected Kidney Transplant (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

Complement component 3, or C3, is a protein of the immune system that plays a central role in the complement system and contributes to innate immunity. Its activation is required for both classical and alternative complement activation pathways. C3d deposition in the renal transplant PTCs (peritubular capillaries) is indicative of AR (acute rejection) with subsequent high probability of graft loss.

Anti-C3d combined with anti-C4d can be utilized as a tool for diagnosis of AR and warrant prompt and aggressive anti-rejection treatment. C3d is also a helpful adjunct in the diagnosis of bullous pemphigoid (BP) and perhaps pemphigus vulgaris (PV), especially in the cases in which only formalin-fixed, paraffin embedded tissue is available for analysis.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal


ISOTYPE: IgG

CONTROL: Rejected Kidney Transplant

LOCALIZATION: Cytoplasmic, Membranous

SPECIES REACTIVITY: Human

C4c, RPab



IHC and IF of C4c on a FFPE Kidney (IHC) and Frozen Glomerulonephritis Tissue (IF)

Complement component 4 (C4), in humans, is a protein involved in the intricate complement system, originating from the human leukocyte antigen (HLA) system. It serves several critical functions in immunity, tolerance, and autoimmunity with the other numerous components.

Low serum complement activity or low protein concentrations of complement C4 are found on Systemic Lupus Erythematosus (SLE) and it is often associated with Congenital C4 deficiency. Complete deficiencies of complement components are among the strongest genetic risk factors for SLE or lupus-like disease, across HLA haplotypes and racial backgrounds.

ANTIBODY TYPE: Rabbit Polyclonal

CLONE: Polyclonal

ISOTYPE: IgG

CONTROL: Testis, Kidney, Pancreas, Salivary Gland, Colon

LOCALIZATION: Cytoplasmic, Membranous


SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3026 | Tinto Predilute | 3.0 ml |
| BSB 3027 | Tinto Predilute | 7.0 ml |
| BSB 3028 | Tinto Predilute | 15.0 ml |
| BSB 3029 | Concentrate | 0.1 ml |
| BSB 3030 | Concentrate | 0.5 ml |
| BSB 3031 | Concentrate | 1.0 ml |
| BSB 3032 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6387 | Tinto Predilute | 3.0 ml |
| BSB 6388 | Tinto Predilute | 7.0 ml |
| BSB 6389 | Tinto Predilute | 15.0 ml |
| BSB 6390 | Concentrate | 0.1 ml |
| BSB 6391 | Concentrate | 0.5 ml |
| BSB 6392 | Concentrate | 1.0 ml |
| BSB 6393 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3033 | Tinto Predilute | 3.0 ml |
| BSB 3034 | Tinto Predilute | 7.0 ml |
| BSB 3035 | Tinto Predilute | 15.0 ml |
| BSB 3036 | Concentrate | 0.1 ml |
| BSB 3037 | Concentrate | 0.5 ml |
| BSB 3038 | Concentrate | 1.0 ml |
| BSB 3039 | Control Slides | 5 |

C4d, RMAb




IHC and IF of C4d on a FFPE Kidney Tissue (IHC) and a Frozen Kidney Rejection Tissue (IF)

Complement component 4, or C4, plays a central role in the complement system. C4d is the final proteolytic remnant of deposited C4b on endothelium and remains covalently attached to endothelium for little more than a week. It is easily detectable by Immunohistochemistry.

Anti-C4d combined with anti-C3d can be utilized as a tool for diagnosis of AR (Acute Rejection) and warrant prompt and aggressive anti-rejection treatment. C4d can be detected in peritubular capillaries in both chronic renal allograft rejection as well as hyperacute rejection, acute vascular rejection, acute cellular rejection, and borderline rejection. It has been shown to be a significant predictor of transplant kidney graft survival and is an aid in treating acute rejection.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP272
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Kidney Transplant Rejection
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

Collagen Type IV, MMAb




IHC and IF of Collagen IV on a FFPE Skin Tissue

Collagen is the main protein of connective tissue in animals and the most abundant protein in mammals, making up about 25% of the total protein content. Collagen IV is a major constituent of the basement membranes, along with laminins and enactins. It is composed of the alpha 1 IV chain and alpha 2 IV chain in a 2:1 ratio. It can form insoluble fibers with high tensile strength.

Normal tissue stains with this antibody in a manner consistent with the sites of mesenchymal elements and epithelial basal laminae. Antibody to collagen IV is useful in detecting the loss of parts of basement membrane in carcinomas. Collagen IV can also be useful in the classification of soft tissue tumors; Schwannomas, Leiomyomas, and their well-differentiated malignant counterparts usually immunoreact to this antibody. The vascular nature of neoplasms, Hemangiopericytoma, Angiosarcoma and Epithelioid Hemangioendothelioma can be observed with this antibody.

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: CIV22
ISOTYPE: IgG1/K
CONTROL: Muscle, Lung
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Rat, Dog, Horse

Collagen Type IV, RMAb



IHC and IF of Collagen Type IV on a FFPE Kidney Tissue

Collagen is the main protein of connective tissue in animals and the most abundant protein in mammals, making up about 25% of the total protein content. Collagen IV is a major constituent of the basement membranes, along with laminins and enactins. It is composed of the alpha 1 IV chain and alpha 2 IV chain in a 2:1 ratio. It can form insoluble fibers with high tensile strength.

Normal tissue stains with this antibody in a manner consistent with the sites of mesenchymal elements and epithelial basal laminae. Antibody to collagen IV is useful in detecting the loss of parts of basement membrane in carcinomas. Collagen IV can also be useful in the classification of soft tissue tumors; Schwannomas, Leiomyomas, and their well-differentiated malignant counterparts usually immunoreact to this antibody. The vascular nature of neoplasms, Hemangiopericytoma, Angiosarcoma and Epithelioid Hemangioendothelioma can be observed with this antibody.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: RBT-COL4
ISOTYPE: IgG
CONTROL: Muscle, Lung, Breast, Placenta, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma
LOCALIZATION: Membranous, Cytoplasmic, Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2831 | Tinto Predilute | 3.0 ml |
| BSB 2832 | Tinto Predilute | 7.0 ml |
| BSB 2833 | Tinto Predilute | 15.0 ml |
| BSB 2834 | Concentrate | 0.1 ml |
| BSB 2835 | Concentrate | 0.5 ml |
| BSB 2836 | Concentrate | 1.0 ml |
| BSB 2837 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5351 | Tinto Predilute | 3.0 ml |
| BSB 5352 | Tinto Predilute | 7.0 ml |
| BSB 5353 | Tinto Predilute | 15.0 ml |
| BSB 5354 | Concentrate | 0.1 ml |
| BSB 5355 | Concentrate | 0.5 ml |
| BSB 5356 | Concentrate | 1.0 ml |
| BSB 5357 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|
| BSB-3777-3 | Tinto Predilute | 3.0 ml |
| BSB-3777-7 | Tinto Predilute | 7.0 ml |
| BSB-3777-15 | Tinto Predilute | 15.0 ml |
| BSB-3777-01 | Concentrate | 0.1 ml |
| BSB-3777-05 | Concentrate | 0.5 ml |
| BSB-3777-1 | Concentrate | 1.0 ml |
| BSB-3777-CS | Control Slides | 5 |

Cytokeratin 5 & 6, RMAb




IHC and IF of Cytokeratin 4 & 6 on an FFPE Mesothelioma Tissue (IHC) an FFPE Colon Carcinoma Tissue (IF)

Cytokeratin 5 (58 kDa) is a high-molecular weight, basic type of cytokeratin expressed in basal, intermediate and superficial-cell layers of stratified epithelia as well as transitional epithelia, complex epithelia, mesothelial cells and Mesothelioma. Cytokeratin 6 (56 kD) is also a high-molecular weight, basic type cytokeratin expressed by proliferating squamous epithelium often paired with Cytokeratin 16.

CK 5 and 6 are positively seen in nearly 100% of Malignant Mesotheliomas and are rarely seen in Lung Adenocarcinomas. CK 5 and 6 can positively be seen in undifferentiated Large-cell Carcinoma as well as Squamous Carcinoma. Fewer than 10% of Carcinomas of the breast, colon, and prostate stain positively for this marker. CK 5 and 6 have also been used successfully as a myoepithelial cell marker in the prostate to determine malignancy.

Cytomegalovirus, MMab




IHC and IF of CMV on a FFPE Infected Lung Tissue

Cytomegalovirus (CMV) is a virus of the Herpesvirus group; in humans it is commonly known as HCMV or Human Herpesvirus 5 (HHV-5). CMV belongs to the Betaherpesvirinae subfamily of Herpesviridae, which also includes Roseolovirus. CMV especially attacks salivary glands. CMV infection can also be life-threatening for patients who are immunocompromised (e.g., patients with HIV or organ-transplant recipients). CMV viruses are found in many mammal species, but CMV species isolated from animals differ from human CMV in terms of genomic structure, and have not been reported to cause human disease.

This Anti-cytomegalovirus antibody cocktail reacts with two different epitopes. The DDG9 antibody reacts with a 76 kDa protein produced by CMV. CCH2 antibody reacts with the early DNA-binding protein p52. There is no cross-reactivity with other Herpesviruses or Adenoviruses. CMV infection is usually seen in immunocompromised patients and involves the GI tract, lung, heart and liver, as well as other organs.

Fibrinogen, RPAb



IHC and IF of Fibrinogen on a FFPE Kidney Tissue (IHC) and on a Frozen Lichen Planus Tissue (IF)

Fibrinogen (factor I) is a glycoprotein that circulates in the blood of vertebrates. During tissue and vascular injury, it is converted enzymatically by thrombin to fibrin and subsequently to a fibrin-based blood clot. Fibrinogen functions primarily to occlude blood vessels and thereby stop excessive bleeding. Fibrin also mediates blood platelet and endothelial cell spreading, tissue fibroblast proliferation, capillary tube formation, and angiogenesis and thereby functions to promote tissue revascularization, wound healing, and tissue repair.

Several disorders (Congenital afibrinogenemia, hypofibrinogenemia, Fibrinogen storage disease, Hereditary fibrinogen A α -Chain amyloidosis, Congenital hypodysfibrinogenemia, Cryofibrinogenemia, acquired hypofibrinogenemia, Chronic Kidney Disease, etc.) in the quantity and/or quality of fibrinogen cause pathological bleeding, pathological blood clotting, and/or the deposition of fibrinogen in the liver, kidneys, and other tissues. Chronic kidney disease (CKD) patients have increased rates of bleeding as well as thrombosis. Fibrinogen and platelets combine to generate a mature clot, but in CKD platelets are dysfunctional.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP24 & EP67
ISOTYPE: IgG
CONTROL: Prostate, Mesothelioma
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 8B1.2, 1G5.2 & 2D4
ISOTYPE: IgG2a
CONTROL: CMV Infected Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Breast, Testis, Kidney, Pancreas, Salivary Gland, Skin, Fallopian Tube
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|----------|-----------------|---------|
| BSB 6604 | Tinto Predilute | 3.0 ml | BSB 5449 | Tinto Predilute | 3.0 ml | BSB 3047 | Tinto Predilute | 3.0 ml |
| BSB 6605 | Tinto Predilute | 7.0 ml | BSB 5450 | Tinto Predilute | 7.0 ml | BSB 3048 | Tinto Predilute | 7.0 ml |
| BSB 6606 | Tinto Predilute | 15.0 ml | BSB 5451 | Tinto Predilute | 15.0 ml | BSB 3049 | Tinto Predilute | 15.0 ml |
| BSB 6607 | Concentrate | 0.1 ml | BSB 5452 | Concentrate | 0.1 ml | BSB 3050 | Concentrate | 0.1 ml |
| BSB 6608 | Concentrate | 0.5 ml | BSB 5453 | Concentrate | 0.5 ml | BSB 3051 | Concentrate | 0.5 ml |
| BSB 6609 | Concentrate | 1.0 ml | BSB 5454 | Concentrate | 1.0 ml | BSB 3052 | Concentrate | 1.0 ml |
| BSB 6610 | Control Slides | 5 | BSB 5455 | Control Slides | 5 | BSB 3053 | Control Slides | 5 |

IgA, MAb




IHC and IF of IgA on a FFPE Tonsil Tissue

IgA, RPaB



IHC and IF of Fibrinogen on a FFPE Kidney Tissue (IHC) and on a Frozen Bullous Dermatitis Tissue (IF)

IgD, RMaB



IHC and IF of IgD on a FFPE Tonsil Tissue

Immunoglobulin A (IgA) is the main immunoglobulin in mucous secretions, including tears, saliva, and colostrum, as well as respiratory, intestinal, prostatic, and vaginal secretions. It is also found in small amounts in blood. Because it is resistant to degradation by enzymes, secretory IgA provides protection against microbes proliferating in body secretions, especially those of the digestive and respiratory tracts.

IgA antibody reacts with surface immunoglobulin IgA alpha chains. It is extremely useful when identifying Acute Leukemias, IgA Myelomas, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. However, due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies.

Immunoglobulin A (IgA) is the main immunoglobulin in mucous secretions, including tears, saliva, and colostrum, as well as respiratory, intestinal, prostatic, and vaginal secretions. It is also found in small amounts in blood. Because it is resistant to degradation by enzymes, secretory IgA provides protection against microbes proliferating in body secretions, especially those of the digestive and respiratory tracts.

IgA antibody reacts with surface immunoglobulin IgA alpha chains. It is extremely useful when identifying Acute Leukemias, IgA Myelomas, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. However, due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q.

IgD makes up about 1% of proteins in the plasma membranes of immature B-lymphocytes (coexpressed with IgM) and is also found in serum in very small amounts. It is monomeric and incorporates the alpha-heavy chain in its structure. IgD's function is currently unknown, as mice lacking IgD seem to retain normal immune responses (implying redundancy if not lack of function), and IgD ceases to be expressed in activated B-lymphocytes. It may function as a regulatory antigen receptor.

IgD antibody reacts with surface immunoglobulin IgD delta chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived from Lymphomas, specifically Marginal Zone Lymphoma.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-39
ISOTYPE: IgG1/K
CONTROL: Tonsil, Spleen, Lymph Node, Kidney, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Spleen, Lymph Node, Kidney, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP173
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5659 | Tinto Predilute | 3.0 ml |
| BSB 5660 | Tinto Predilute | 7.0 ml |
| BSB 5661 | Tinto Predilute | 15.0 ml |
| BSB 5662 | Concentrate | 0.1 ml |
| BSB 5663 | Concentrate | 0.5 ml |
| BSB 5664 | Concentrate | 1.0 ml |
| BSB 5665 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3054 | Tinto Predilute | 3.0 ml |
| BSB 3055 | Tinto Predilute | 7.0 ml |
| BSB 3056 | Tinto Predilute | 15.0 ml |
| BSB 3057 | Concentrate | 0.1 ml |
| BSB 3058 | Concentrate | 0.5 ml |
| BSB 3059 | Concentrate | 1.0 ml |
| BSB 3685 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 2957 | Tinto Predilute | 3.0 ml |
| BSB 2958 | Tinto Predilute | 7.0 ml |
| BSB 2959 | Tinto Predilute | 15.0 ml |
| BSB 2960 | Concentrate | 0.1 ml |
| BSB 2961 | Concentrate | 0.5 ml |
| BSB 2962 | Concentrate | 1.0 ml |
| BSB 2963 | Control Slides | 5 |

IgD, RPab




IHC and IF of IgD on a FFPE Colon Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

IgE, RPab



IHC and IF of IgE on a FFPE Tonsil Tissue (IHC) and on a Frozen Colon Tissue (IF)

IgG, RPab



IHC and IF of IgG on a FFPE Tonsil Tissue (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

IgD makes up about 1% of proteins in the plasma membranes of immature B-lymphocytes (coexpressed with IgM) and is also found in serum in very small amounts. It is monomeric and incorporates the alpha-heavy chain in its structure. IgD's function is currently unknown, as mice lacking IgD seem to retain normal immune responses (implying redundancy if not lack of function), and IgD ceases to be expressed in activated B-lymphocytes. It may function as a regulatory antigen receptor. IgD is the major antigen receptor isotype co-expressed with IgM on the surface of most peripheral B cells in mice and humans.

The IgD antibody reacts with surface immunoglobulin IgD delta chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived from Lymphomas, specifically Marginal Zone Lymphoma. Renal involvement in systemic lupus erythematosus (SLE) is associated with production of antibodies to double stranded DNA, deposition of immune complexes and organ damage. Lupus nephritis patients were characterized by increased percentage of immature/early-transitional B-cells (CD27-IgD+CD21-), higher frequency of activated switched memory (SM, CD27+IgD-CD21-) and exhausted memory B-cells (CD27-IgD-), and decrease in non-switched memory (NSM, CD27+IgD+) B-cells.

IgE, Immunoglobulin E, is an isotype of antibody only found in mammals. IgE is synthesized by plasma cells. Monomers of IgE consist of two heavy chains (ϵ chain) and two light chains, with the ϵ chain containing 4 Ig-like constant domains (C ϵ 1-C ϵ 4). IgE's main function is immunity to parasites such as helminths like *Schistosoma mansoni*, *Trichinella spiralis*, and *Fasciola hepatica*. IgE is utilized during immune defense against certain protozoan parasites such as *Plasmodium falciparum*.

IgE also has an essential role in type I hypersensitivity, which manifests in various allergic diseases, such as allergic asthma, most types of sinusitis, allergic rhinitis, food allergies, and specific types of chronic urticaria and atopic dermatitis. IgE also plays a pivotal role in responses to allergens, such as: anaphylactic drugs, bee stings, and antigen preparations used in desensitization immunotherapy. IgE is known to be elevated in various autoimmune disorders such as Lupus(SLE), Rheumatoid Arthritis(RA) & psoriasis, and is theorized to be of pathogenetic importance in RA and SLE by eliciting a hypersensitivity reaction.

IgG is a monomeric immunoglobulin, comprised of two heavy chains and two light chains. This is the most abundant immunoglobulin and is approximately equally distributed in blood and tissue liquids, constituting 75% of serum immunoglobulins in humans. This is the only isotype that can pass through the placenta and bind to many kinds of pathogens. IgG protects the body against them by complement activation (classic pathway), opsonization for phagocytosis and neutralization of their toxins. There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%).

IgG antibody reacts with surface immunoglobulin IgG gamma chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Thymus, Colon
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Kidney, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5666 | Tinto Predilute | 3.0 ml |
| BSB 5667 | Tinto Predilute | 7.0 ml |
| BSB 5668 | Tinto Predilute | 15.0 ml |
| BSB 5669 | Concentrate | 0.1 ml |
| BSB 5670 | Concentrate | 0.5 ml |
| BSB 5671 | Concentrate | 1.0 ml |
| BSB 5672 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3067 | Tinto Predilute | 3.0 ml |
| BSB 3068 | Tinto Predilute | 7.0 ml |
| BSB 3069 | Tinto Predilute | 15.0 ml |
| BSB 3070 | Concentrate | 0.1 ml |
| BSB 3071 | Concentrate | 0.5 ml |
| BSB 3072 | Concentrate | 1.0 ml |
| BSB 3073 | Control Slides | 5 |


| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3074 | Tinto Predilute | 3.0 ml |
| BSB 3075 | Tinto Predilute | 7.0 ml |
| BSB 3076 | Tinto Predilute | 15.0 ml |
| BSB 3076 | Concentrate | 0.1 ml |
| BSB 3077 | Concentrate | 0.5 ml |
| BSB 3078 | Concentrate | 1.0 ml |
| BSB 3079 | Control Slides | 5 |

IgG, MAb




IHC and IF of IgG on a FFPE Tonsil Tissue

IgG4, MAb



IHC and IF of IgG4 on a FFPE Tonsil Tissue

IgG4, RMAb



IHC and IF of IgG4 on a FFPE Tonsil Tissue

IgG is a monomeric immunoglobulin, comprised of two heavy chains and two light chains. This is the most abundant immunoglobulin and is approximately equally distributed in blood and tissue liquids, constituting 75% of serum immunoglobulins in humans. This is the only isotype that can pass through the placenta and bind to many kinds of pathogens. IgG protects the body against them by complement activation (classic pathway), opsonization for phagocytosis and neutralization of their toxins. There are 4 subclasses: IgG1 (66%), IgG2 (23%), IgG3 (7%) and IgG4 (4%).

IgG antibody reacts with surface immunoglobulin IgG gamma chains. This antibody is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q. Clinically, hematuria and proteinuria are present, with or without nephrotic syndromes. Mesangial IgG glomerulonephritis has been recently recognized as a distinct type of glomerulonephritis. The morphologic criteria detected in these patients included mesangial dense deposits by ultrastructural studies, which were predominantly positive for IgG by immunofluorescence.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-40
ISOTYPE: IgG2a/K
CONTROL: Tonsil, Lymph Node, Kidney, Spleen
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

IgG4-related sclerosing disease has been recognized as a systemic disease entity characterized by an elevated serum IgG4 level, sclerosing fibrosis and diffuse lymphoplasmacytic infiltration with the presence of many IgG4-positive plasma cells. As these patients tend to respond favorably to steroid treatment, it is important to recognize this entity and differentiate it from such mimics as lymphoma.

Clinical manifestations are apparent in the pancreas, bile duct, gallbladder, lacrimal gland, salivary gland, retroperitoneum, kidney, lung, breast, thyroid, and prostate. Immunohistochemical analyses in the case of IgG4-related sclerosing disease not only exhibits significantly more IgG4-positive plasma cells in affected tissues but also significantly higher IgG4/ IgG ratios (typically > 30%).

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-96
ISOTYPE: IgG2a/K
CONTROL: Tonsil, Spleen, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

IgG4-related sclerosing disease has been recognized as a systemic disease entity characterized by an elevated serum IgG4 level, sclerosing fibrosis and diffuse lymphoplasmacytic infiltration with the presence of many IgG4-positive plasma cells. As these patients tend to respond favorably to steroid treatment, it is important to recognize this entity and differentiate it from such mimics as lymphoma.

Clinical manifestations are apparent in the pancreas, bile duct, gallbladder, lacrimal gland, salivary gland, retroperitoneum, kidney, lung, breast, thyroid, and prostate. Immunohistochemical analyses in the case of IgG4-related sclerosing disease not only exhibits significantly more IgG4-positive plasma cells in affected tissues but also significantly higher IgG4/ IgG ratios (typically 30%).


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP138
ISOTYPE: IgG
CONTROL: Tonsil, Spleen, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5673 | Tinto Predilute | 3.0 ml |
| BSB 5674 | Tinto Predilute | 7.0 ml |
| BSB 5675 | Tinto Predilute | 15.0 ml |
| BSB 5676 | Concentrate | 0.1 ml |
| BSB 5677 | Concentrate | 0.5 ml |
| BSB 5678 | Concentrate | 1.0 ml |
| BSB 5679 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6807 | Tinto Predilute | 3.0 ml |
| BSB 6808 | Tinto Predilute | 7.0 ml |
| BSB 6809 | Tinto Predilute | 15.0 ml |
| BSB 6810 | Concentrate | 0.1 ml |
| BSB 6811 | Concentrate | 0.5 ml |
| BSB 6812 | Concentrate | 1.0 ml |
| BSB 6823 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 6814 | Tinto Predilute | 3.0 ml |
| BSB 6815 | Tinto Predilute | 7.0 ml |
| BSB 6816 | Tinto Predilute | 15.0 ml |
| BSB 6817 | Concentrate | 0.1 ml |
| BSB 6818 | Concentrate | 0.5 ml |
| BSB 6819 | Concentrate | 1.0 ml |
| BSB 6820 | Control Slides | 5 |

IgM, RPAb




IHC and IF of IgM on a FFPE Kidney Tissue (IHC) and on a Frozen Lupus Erythematosus Tissue (IF)

IgM forms polymers where multiple immunoglobulins are covalently linked together with disulfide bonds, normally as a pentamer or occasionally as a hexamer. It has a large molecular mass of approximately 900 kDa (in its pentamer form). In germline cells, the gene segment encoding the constant region of the heavy chain is positioned first among other constant region gene segments. For this reason, IgM is the first immunoglobulin expressed by mature B-cells.

IgM antibody reacts with surface immunoglobulin IgM mu chains. IgM is one of the predominant surface immunoglobulins on B-lymphocytes, and is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies. Lupus nephritis is an inflammation of the kidneys caused by Systemic Lupus Erythematosus. Immunofluorescence reveals positively for IgG, IgA, IgM, C3, and C1q. Clinically, hematuria and proteinuria are present, with or without nephrotic syndromes. Immunoglobulin M (IgM) nephropathy is an uncommon glomerular disease characterized by IgM deposits in the mesangium.

ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

IgM, MAb




IHC and IF of IgM on a FFPE Kidney Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

IgM forms polymers where multiple immunoglobulins are covalently linked together with disulfide bonds, normally as a pentamer or occasionally as a hexamer. It has a large molecular mass of approximately 900 kDa (in its pentamer form). In germline cells, the gene segment encoding the constant region of the heavy chain is positioned first among other constant region gene segments. For this reason, IgM is the first immunoglobulin expressed by mature B-cells.

IgM antibody reacts with surface immunoglobulin IgM mu chains. IgM is one of the predominant surface immunoglobulins on B-lymphocytes, and is useful when identifying Leukemias, Plasmacytomas, and B-cell lineage derived Hodgkin's Lymphomas. Due to the restricted expression of heavy and light chains in these diseases, demonstration of B-cell Lymphomas is possible with clonal gene-rearrangement studies.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-41
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node, Spleen, Kidney, Colon
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

Kappa, RPAb



IHC and IF of Kappa on a FFPE Kidney Tissue (IHC) and on a Frozen Colon Tissue (IF with FluoroMounter DAPI)

Kappa detects surface immunoglobulin on normal and neoplastic B-cells. In paraffin-embedded tissue, Kappa exhibits strong staining of kappa-positive plasma cells and cells that have absorbed exogenous immunoglobulin.

When studying B-cell neoplasms, the determination of light-chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either Kappa or Lambda light chains, whereas reactive proliferations display a mixture of Kappa and Lambda-positive cells. If only a single light-chain type is detected, a lympho-proliferative disorder is very likely. Monoclonality is determined by a Kappa-Lambda ratio greater than or equal to 3:1, a Lambda-Kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population.

In IgG-dominant immune complex-mediated glomerulonephritis, there are multiple pathological findings that strongly suggest the diagnosis of Lupus Nephritis including immunofluorescence staining for IgG, IgM, IgA, Kappa or Lambda, C3 and C1.


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3080 | Tinto Predilute | 3.0 ml |
| BSB 3081 | Tinto Predilute | 7.0 ml |
| BSB 3082 | Tinto Predilute | 15.0 ml |
| BSB 3083 | Concentrate | 0.1 ml |
| BSB 3084 | Concentrate | 0.5 ml |
| BSB 3085 | Concentrate | 1.0 ml |
| BSB 3086 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5680 | Tinto Predilute | 3.0 ml |
| BSB 5681 | Tinto Predilute | 7.0 ml |
| BSB 5682 | Tinto Predilute | 15.0 ml |
| BSB 5683 | Concentrate | 0.1 ml |
| BSB 5684 | Concentrate | 0.5 ml |
| BSB 5685 | Concentrate | 1.0 ml |
| BSB 5686 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3086 | Tinto Predilute | 3.0 ml |
| BSB 3087 | Tinto Predilute | 7.0 ml |
| BSB 3088 | Tinto Predilute | 15.0 ml |
| BSB 3089 | Concentrate | 0.1 ml |
| BSB 3090 | Concentrate | 0.5 ml |
| BSB 3091 | Concentrate | 1.0 ml |
| BSB 3092 | Control Slides | 5 |

Kappa Light Chains, MMab




IHC and IF of Kappa on a FFPE Tonsil Tissue (IHC) and on a Frozen Spleen Tissue (IF)

Kappa detects surface immunoglobulin on normal and neoplastic B-cells. In paraffin-embedded tissue, Kappa exhibits strong staining of kappa-positive plasma cells and cells that have absorbed exogenous immunoglobulin.

When studying B-cell neoplasms, the determination of light-chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population.

Lambda, MMab




IHC and IF of Lambda on a FFPE Tonsil Tissue

Lambda detects surface immunoglobulin on normal and neoplastic B-cells. Lambda staining is seen in B-cell follicles of human lymphoid tissue.

When studying B-cell neoplasms, the determination of light chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population.

Lambda, RPaB



IHC and IF of Lambda on a FFPE Tonsil Tissue (IHC) and on a Frozen Tonsil Tissue (IF)

Lambda detects surface immunoglobulin on normal and neoplastic B-cells. Lambda staining is seen in B-cell follicles of human lymphoid tissue.

When studying B-cell neoplasms, the determination of light chain ratios remains the centerpiece. This is sound reasoning because most B-cell Lymphomas express either kappa or lambda light chains, whereas reactive proliferations display a mixture of kappa and lambda-positive cells. If only a single light-chain type is detected, a lymphoproliferative disorder is very likely. Monoclonality is determined by a kappa-lambda ratio greater than or equal to 3:1, a lambda-kappa ratio greater than or equal to 2:1, or a monoclonal population of 75% or more of the total population. In IgG-dominant immune complex-mediated glomerulonephritis, there are multiple pathological findings that strongly suggest the diagnosis of Lupus Nephritis including immunofluorescence staining for IgG, IgM, IgA, Kappa or Lambda, C3 and C1.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-58
ISOTYPE: IgG1/K
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-16
ISOTYPE: IgG2a
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human, Dog, Cat


ANTIBODY TYPE: Rabbit Polyclonal
CLONE: Polyclonal
ISOTYPE: IgG
CONTROL: Tonsil, Lymph Node
LOCALIZATION: Cytoplasmic
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5701 | Tinto Predilute | 3.0 ml |
| BSB 5702 | Tinto Predilute | 7.0 ml |
| BSB 5703 | Tinto Predilute | 15.0 ml |
| BSB 5704 | Concentrate | 0.1 ml |
| BSB 5705 | Concentrate | 0.5 ml |
| BSB 5706 | Concentrate | 1.0 ml |
| BSB 5707 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 5715 | Tinto Predilute | 3.0 ml |
| BSB 5716 | Tinto Predilute | 7.0 ml |
| BSB 5717 | Tinto Predilute | 15.0 ml |
| BSB 5718 | Concentrate | 0.1 ml |
| BSB 5719 | Concentrate | 0.5 ml |
| BSB 5720 | Concentrate | 1.0 ml |
| BSB 5721 | Control Slides | 5 |

| CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|
| BSB 3092 | Tinto Predilute | 3.0 ml |
| BSB 3093 | Tinto Predilute | 7.0 ml |
| BSB 3094 | Tinto Predilute | 15.0 ml |
| BSB 3095 | Concentrate | 0.1 ml |
| BSB 3096 | Concentrate | 0.5 ml |
| BSB 3097 | Concentrate | 1.0 ml |
| BSB 3098 | Control Slides | 5 |


p40, RMab



IHC and IF of p40 on an FFPE Prostate Tissue (IHC)
an FFPE Tonsil Tissue (IF)

p40 is an antibody that recognizes Δ Np63-a p63 isoform and it is highly specific for squamous/basal cells. It may be a valuable marker in detecting Squamous Cell Carcinoma where p63 is currently used. It recognizes the shortest variant of p53. p40 is superior in specificity to p63 because it does not label lung adenocarcinomas like p63 does, which eliminates the potential of misinterpreting a positive adenocarcinoma as a squamous cell carcinoma.

p63, MMab




IHC and IF of p63 on an FFPE Basal Cell Carcinoma
Tissue

In addition to p53, mammalian cells contain two homologous genes, p63 and p73. These genes give rise to the expression of proteins that are highly similar to p53 in structure and function. In particular, p63 and p73 proteins can induce p53-responsive genes and elicit programmed cell death. p73 and p63 are important during development and differentiation. In particular, p63 appears to be primarily implicated in epithelial development.

Anti-p63 to human p63 protein labels an epitope common to all six p63 isoforms (TAp63 α , TAp63 β , TAp63 γ , Δ Np63 α , Δ Np63 β , Δ Np63 γ). p63 labels the nuclei of myoepithelial cells in the prostate gland as well as breast tissue, making it useful in differentiating benign vs. malignant prostate lesions and breast lesions.

p63, RMab



IHC and IF of p63 on a FFPE Prostate Tissue (IHC)
and an FFPE Skin Tissue (IF)

In addition to p53, mammalian cells contain two homologous genes, p63 and p73. These genes give rise to the expression of proteins that are highly similar to p53 in structure and function. In particular, p63 and p73 proteins can induce p53-responsive genes and elicit programmed cell death. p73 and p63 are more important during development and differentiation. In particular, p63 appears to be primarily implicated in epithelial development.

Anti-p63 to human p63 protein labels an epitope common to all six p63 isoforms (TAp63 α , TAp63 β , TAp63 γ , Δ Np63 α , Δ Np63 β , Δ Np63 γ). p63 labels the nuclei of myoepithelial cells in the prostate gland as well as breast tissue, making it useful in differentiating benign vs. malignant prostate lesions and breast lesions.


ANTIBODY TYPE: Rabbit Monoclonal
CLONE: ZR8
ISOTYPE: IgG
CONTROL: Normal Prostate, Breast, Skin
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Mouse Monoclonal
CLONE: 4A4
ISOTYPE: IgG2a/K
CONTROL: Prostate, Breast, Skin, Salivary Gland
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

ANTIBODY TYPE: Rabbit Monoclonal
CLONE: EP174
ISOTYPE: IgG
CONTROL: Prostate, Breast, Skin, Salivary Gland
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human, Mouse, Rat, Rabbit

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|----------|-----------------|---------|----------|-----------------|---------|------------|-----------------|---------|
| BSB 2070 | Tinto Predilute | 3.0 ml | BSB 3602 | Tinto Predilute | 3.0 ml | BSB 5848 | Tinto Predilute | 3.0 ml |
| BSB 2071 | Tinto Predilute | 7.0 ml | BSB 3603 | Tinto Predilute | 7.0 ml | BSB 5849 | Tinto Predilute | 7.0 ml |
| BSB 2072 | Tinto Predilute | 15.0 ml | BSB 3604 | Tinto Predilute | 15.0 ml | BSB 5850 | Tinto Predilute | 15.0 ml |
| BSB 2073 | Concentrate | 0.1 ml | BSB 3605 | Concentrate | 0.1 ml | BSB 5851 | Concentrate | 0.1 ml |
| BSB 2074 | Concentrate | 0.5 ml | BSB 3606 | Concentrate | 0.5 ml | BSB 5852 | Concentrate | 0.5 ml |
| BSB 2075 | Concentrate | 1.0 ml | BSB 3607 | Concentrate | 1.0 ml | BSB 5853 | Concentrate | 1.0 ml |
| BSB 2076 | Control Slides | 5 | BSB 3608 | Control Slides | 5 | BSB 5854-1 | Control Slides | 5 |

PLA2R1, MAb




IHC and IF of PLA2R1 on a FFPE Membranous Glomerulopathy Tissue

PLA2R1 gene encodes phospholipase A2 receptor 1 protein, a 180 kDa transmembrane glycoproteins expressed by podocyte. Approximately 70% of patients with idiopathic membranous glomerulopathy have autoantibodies directed against podocyte PLA2R1. PLA2R1 also promotes tumor suppressive responses including senescence, apoptosis, and inhibition of transformation. Known oncogenes such as HIF2α and c-Myc repress PLA2R1 expression.

PLA2R1 gain or loss of function experiments in vitro and in vivo shows that this receptor promotes several tumor suppressive responses including senescence, apoptosis and inhibition of transformation. Supporting a tumor suppressive role of PLA2R1, its expression decreases in numerous cancers, and known oncogenes such as HIF2α and c-MYC repress its expression. PLA2R1 promoter methylation, a classical way to repress tumor suppressive gene expression in cancer cells, is observed in leukemia, in kidney and in breast cancer cells. PLA2R1 also promotes accumulation of reactive oxygen species which induce cell death and senescence. This review compiles recent data demonstrating an unexpected tumor suppressive role of PLA2R1 and outlines the future work needed to improve our knowledge of the functions of this gene in cancer

ANTIBODY TYPE: Mouse Monoclonal
CLONE: BSB-129
ISOTYPE: IgG1
CONTROL: Brain, Testis, Kidney, Salivary Gland, Gastric GIST
LOCALIZATION: Cytoplasmic, Membranous
SPECIES REACTIVITY: Human

SV40, MAb



IHC and IF of SV40 on a FFPE Infected Kidney Tissue

SV40 is an abbreviation for Simian vacuolating virus 40 or Simian virus 40, a polyomavirus that is found in both monkeys and humans. Like other polyomaviruses, SV40 is a DNA virus that has the potential to cause tumors, but most often persists as a latent infection.

ANTIBODY TYPE: Mouse Monoclonal
CLONE: Pab101
ISOTYPE: IgG2a
CONTROL: SV40 Infected Tissue
LOCALIZATION: Nuclear
SPECIES REACTIVITY: Human

| CAT. # | PRESENTATION | VOL/QTY | CAT. # | PRESENTATION | VOL/QTY |
|-------------|-----------------|---------|----------|-----------------|---------|
| BSB-2372-3 | Tinto Predilute | 3.0 ml | BSB 2230 | Tinto Predilute | 3.0 ml |
| BSB-2372-7 | Tinto Predilute | 7.0 ml | BSB 2231 | Tinto Predilute | 7.0 ml |
| BSB-2372-15 | Tinto Predilute | 15.0 ml | BSB 2232 | Tinto Predilute | 15.0 ml |
| BSB-2372-01 | Concentrate | 0.1 ml | BSB 2233 | Concentrate | 0.1 ml |
| BSB-2372-05 | Concentrate | 0.5 ml | BSB 2234 | Concentrate | 0.5 ml |
| BSB-2372-1 | Concentrate | 1.0 ml | BSB 2235 | Concentrate | 1.0 ml |
| BSB-2372-CS | Control Slides | 5 | BSB 2236 | Control Slides | 5 |


The InDirect Immunofluorescence Detection AmpliDetector Plus FITC kit

Bio SB offers a highly sensitive Immunofluorescence detection system for the IVD antibodies related to autoimmune conditions. Our innovative AmpliDetector Plus FITC detection system and high affinity antibodies, have opened the doors for a faster and accurate Immunofluorescence applicable to Autoimmune Disease like Nephropathies and Lupus.

AmpliDetector Plus FITC is intended for use in Immunofluorescence (IF) applications of formalin-fixed paraffin-embedded tissues (FFPE), frozen tissue sections and cell preparations.

AmpliDetector Plus FITC Detection kit Overview

- Non-Biotin, Fab Micropolymer HRP Amplification detection technology.
- Fab Micropolymer HRP Amplification detection technology allows for better cell penetration to deliver a highly specific and sensitive signal.
- Ready-to-use, high sensitivity system especially designed for InDirect IF detection of formalin-fixed paraffin-embedded or frozen tissue sections and cell preparations.
- Universal: detects mouse monoclonal, rabbit monoclonal and polyclonal antibodies.
- For in vitro diagnostic Use. All kits manufactured according to US FDA and ISO 13485 guidelines.



Mouse / Rabbit AmpliDetector Plus FITC kit


| Product Description | Volume | Catalog # |
|--------------------------------|----------|--------------|
| Rabbit AmpliDetector Plus FITC | 15.0 ml | BSB-0359-15 |
| Rabbit AmpliDetector Plus FITC | 50.0 ml | BSB-0359-50 |
| Rabbit AmpliDetector Plus FITC | 100.0 ml | BSB-0359-100 |

Mouse / Rabbit AmpliDetector Plus FITC with FluoroMounter kit

| Product Description | Volume | Catalog # |
|---|----------|--------------|
| Rabbit AmpliDetector Plus FITC with FluoroMounter | 15.0 ml | BSB-0358-15 |
| Rabbit AmpliDetector Plus FITC with FluoroMounter | 50.0 ml | BSB-0358-50 |
| Rabbit AmpliDetector Plus FITC with FluoroMounter | 100.0 ml | BSB-0358-100 |

Mouse / Rabbit AmpliDetector Plus FITC with FluoroMounter with DAPI kit

| Product Description | Volume | Catalog # |
|---|----------|--------------|
| Rabbit AmpliDetector Plus FITC with FluoroMounter with DAPI | 15.0 ml | BSB-0357-15 |
| Rabbit AmpliDetector Plus FITC with FluoroMounter with DAPI | 50.0 ml | BSB-0357-50 |
| Rabbit AmpliDetector Plus FITC with FluoroMounter with DAPI | 100.0 ml | BSB-0357-100 |



C4d IF on a Frozen Kidney Tissue



www.BioSB.com






BIO-SB SAMPLE MICROARRAY
TISSUE STAIN SLIDE 001



Detection Systems for IHC

ImmunoDetector | PolyDetector | PolyDetector Plus


*A wide array of highly sensitive Biotin & Fab Micropolymer based detection systems.
For use in clinical and research applications.*



PolyDetector AP


For the Immunohistochemical detection of antigens in cells and formalin-fixed or frozen tissues

The Mouse/Rabbit PolyDetector Alkaline Phosphatase (AP) Detection System is a one-step polymeric detection system that allows for the demonstration of antigens in formalin-fixed paraffin-embedded tissues, cryostat sections, blood smears, cytospreads, and cell preparations. The increased sensitivity of this system allows for rapid staining procedures without compromising stain quality.



- Non-Biotin, 1-Step Immunohistochemistry Detection Technology
- Anti-mouse and rabbit Fab micro-polymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed or frozen tissues
- Universal: Detects Mouse or Rabbit antibodies
- ALK Red, ALK Scarlet, ALK Blue, ALK Magenta and ALK Brown Configurations
- For *in Vitro* Diagnostic Use

Recommended Protocol





PolyDetector AP Detection Systems

| Product Description | Volume | Catalog # |
|---|-----------|---------------|
| Mouse/Rabbit PolyDetector AP, w/ALK Scarlet | 5 ml | BSB-0351-15 |
| Mouse/Rabbit PolyDetector AP, w/ALK Scarlet | 15.0 ml | BSB-0351-50 |
| Mouse/Rabbit PolyDetector AP, w/ALK Scarlet | 50.0 ml | BSB-0351-100 |
| Mouse/Rabbit PolyDetector AP, w/ALK Scarlet | 100.0 ml | BSB-0351-200 |
| Mouse/Rabbit PolyDetector AP, w/ALK Scarlet | 200.0 ml | BSB-0351-1000 |
| Mouse/Rabbit PolyDetector AP, w/ALK Scarlet | 1000.0 ml | BSB 0009 |

PolyDetector AP Label

| Product Description | Volume | Catalog # |
|------------------------------------|-----------|-----------|
| Mouse/Rabbit PolyDetector AP Label | 15.0 ml | BSB 0287 |
| Mouse/Rabbit PolyDetector AP Label | 50.0 ml | BSB 0288 |
| Mouse/Rabbit PolyDetector AP Label | 100.0 ml | BSB 0289 |
| Mouse/Rabbit PolyDetector AP Label | 200.0 ml | BSB 0290 |
| Mouse/Rabbit PolyDetector AP Label | 1000.0 ml | BSB 0291 |




IVD For *in Vitro* Diagnostic Use

10 tests per ml considering 100µl per tissue

PolyDetector HRP


For the Immunohistochemical detection of antigens in cells and formalin-fixed or frozen tissues

The Bio SB PolyDetector is a non-Biotin, Fab Micropolymer detection system that allows for the demonstration of antigens in formalin-fixed paraffin-embedded tissues, cryostat sections, blood smears, cytospins, and cell preparations. The PolyDetector technology was developed and is manufactured with a proprietary Micropolymer backbone conjugated to Anti-Mouse and Anti-Rabbit Fab Ig's plus high quality HRP or AP enzymes. The elimination of the Anti-Mouse and Anti-Rabbit immunoglobulin Fc region reduces non-specific reactions. This ensures consistent and reproducible immunostaining for all types of nuclear, cytoplasmic and membranous antigens, in different types of tissues. The increased sensitivity of this system allows for rapid staining procedures without compromising stain quality.



- Non-Biotin, 1-Step Immunohistochemistry Detection Technology
- Anti-mouse and rabbit Fab micro-polymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed or frozen tissues
- Universal: Detects Mouse or Rabbit antibodies
- DAB or AEC Configurations. Also works with HRP Green, HRP Blue and HRP Black
- For *in Vitro* Diagnostic Use

Recommended Protocol




PolyDetector HRP Detection Systems

| Product Description | Volume | Catalog # |
|---|-----------|-----------|
| Mouse/Rabbit PolyDetector DAB HRP Brown | 5 ml | BSB 0201S |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 15.0 ml | BSB 0201 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 50.0 ml | BSB 0203 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 100.0 ml | BSB 0205 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 200.0 ml | BSB 0207 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 1000.0 ml | BSB 0207A |
| Mouse/Rabbit PolyDetector AEC HRP Red | 5 ml | BSB 0202S |
| Mouse/Rabbit PolyDetector AEC HRP Red | 15.0 ml | BSB 0202 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 50.0 ml | BSB 0204 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 100.0 ml | BSB 0206 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 200.0 ml | BSB 0208 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 1000.0 ml | BSB 0208A |

PolyDetector HRP Label

| Product Description | Volume | Catalog # |
|-------------------------------------|-----------|------------|
| Mouse/Rabbit PolyDetector HRP Label | 15.0 ml | BSB 0201H |
| Mouse/Rabbit PolyDetector HRP Label | 50.0 ml | BSB 0203H |
| Mouse/Rabbit PolyDetector HRP Label | 100.0 ml | BSB 0205H |
| Mouse/Rabbit PolyDetector HRP Label | 200.0 ml | BSB 0207H |
| Mouse/Rabbit PolyDetector HRP Label | 1000.0 ml | BSB 0207AH |




IHC of S100 on a FFPE Melanoma Tissue

PolyDetector Plus HRP


For the Immunohistochemical detection of antigens in cells and formalin-fixed or frozen tissues

The BioSB PolyDetector Plus system is a highly sensitive non-Biotin, Fab Micropolymer detection system that allows for the demonstration of antigens in formalin-fixed paraffin-embedded tissues, cryostat sections, blood smears, cytospins, and cell preparations. The PolyDetector Plus technology was developed and is manufactured with a proprietary Micropolymer backbone conjugated to Anti-Mouse and Anti-Rabbit Fab Ig's plus high quality HRP or AP enzymes. The elimination of the Anti-Mouse and Anti-Rabbit immunoglobulin Fc region reduces non-specific reactions. The PolyDetector Plus kit incorporates an Immunoglobulin link and a Fab Micropolymer label. This multiple component Fab Micropolymer delivers a highly sensitive and specific signal in a shorter time-frame than the PolyDetector.



- Non-Biotin, 3-Step Immunohistochemistry Detection Technology
- Anti-mouse and rabbit Fab micro-polymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed or frozen tissues
- Universal: Detects Mouse or Rabbit antibodies
- DAB or AEC Configurations
- For *in Vitro* Diagnostic Use

Recommended Protocol




PolyDetector HRP Detection Systems

| Product Description | Volume | Catalog # |
|---|-----------|-----------|
| Mouse/Rabbit PolyDetector DAB HRP Brown | 5 ml | BSB 0201S |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 15.0 ml | BSB 0201 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 50.0 ml | BSB 0203 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 100.0 ml | BSB 0205 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 200.0 ml | BSB 0207 |
| Mouse/Rabbit PolyDetector DAB HRP Brown | 1000.0 ml | BSB 0207A |
| Mouse/Rabbit PolyDetector AEC HRP Red | 5 ml | BSB 0202S |
| Mouse/Rabbit PolyDetector AEC HRP Red | 15.0 ml | BSB 0202 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 50.0 ml | BSB 0204 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 100.0 ml | BSB 0206 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 200.0 ml | BSB 0208 |
| Mouse/Rabbit PolyDetector AEC HRP Red | 1000.0 ml | BSB 0208A |

PolyDetector HRP Label

| Product Description | Volume | Catalog # |
|-------------------------------------|-----------|------------|
| Mouse/Rabbit PolyDetector HRP Label | 15.0 ml | BSB 0201H |
| Mouse/Rabbit PolyDetector HRP Label | 50.0 ml | BSB 0203H |
| Mouse/Rabbit PolyDetector HRP Label | 100.0 ml | BSB 0205H |
| Mouse/Rabbit PolyDetector HRP Label | 200.0 ml | BSB 0207H |
| Mouse/Rabbit PolyDetector HRP Label | 1000.0 ml | BSB 0207AH |




IHC of Melanoma PNL2 on a FFPE Melanoma Tissue

ImmunoDetector AP


For the Immunohistochemical detection of antigens in cells and formalin-fixed or frozen tissues

The Bio SB ImmunoDetector Alkaline Phosphatase (AP) Detection System is a highly sensitive Biotin-Streptavidin-Alkaline Phosphatase Detection System that allows for the demonstration of antigens in formalin-fixed paraffin-embedded tissue, cryostat sections, cytosmears, and cell preparations. The increased sensitivity of ImmunoDetector AP Detection System allows for rapid staining procedures without compromises in the quality of stains.



- Biotin-Streptavidin, 2-Step Immunohistochemistry Detection Technology
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed or frozen tissues
- Universal: Detects Mouse or Rabbit antibodies
- It can be used with ALK Red, ALK Scarlet, ALK Magenta, ALK Brown and ALK Blue substrate-chromogen
- For *in Vitro* Diagnostic Use

Recommended Protocol



ImmunoDetector AP Detection Systems

| Product Description | Volume | Catalog # |
|---|----------|---------------|
| Mouse/Rabbit ImmunoDetector AP, w/ALK Scarlet | 5 ml | BSB-0350-15 |
| Mouse/Rabbit ImmunoDetector AP, w/ALK Scarlet | 15.0 ml | BSB-0350-50 |
| Mouse/Rabbit ImmunoDetector AP, w/ALK Scarlet | 50.0 ml | BSB-0350-100 |
| Mouse/Rabbit ImmunoDetector AP, w/ALK Scarlet | 100.0 ml | BSB-0350-200 |
| Mouse/Rabbit ImmunoDetector AP, w/ALK Scarlet | 200.0 ml | BSB-0350-1000 |


ImmunoDetector AP Link & Label

| Product Description | Volume | Catalog # |
|--------------------------------------|-----------|-----------|
| Mouse/Rabbit ImmunoDetector AP Label | 15.0 ml | BSB 0082A |
| Mouse/Rabbit ImmunoDetector AP Label | 50.0 ml | BSB 0083A |
| Mouse/Rabbit ImmunoDetector AP Label | 100.0 ml | BSB 0084A |
| Mouse/Rabbit ImmunoDetector AP Label | 200.0 ml | BSB 0085A |
| Mouse/Rabbit ImmunoDetector AP Label | 1000.0 ml | BSB 0086A |

ImmunoDetector AP Link & Label

| Product Description | Volume | Catalog # |
|--|-----------|------------|
| Mouse/Rabbit ImmunoDetector Biotin Link & AP Label | 15.0 ml | BSB 0082LA |
| Mouse/Rabbit ImmunoDetector Biotin Link & AP Label | 50.0 ml | BSB 0083LA |
| Mouse/Rabbit ImmunoDetector Biotin Link & AP Label | 100.0 ml | BSB 0084LA |
| Mouse/Rabbit ImmunoDetector Biotin Link & AP Label | 200.0 ml | BSB 0085LA |
| Mouse/Rabbit ImmunoDetector Biotin Link & AP Label | 1000.0 ml | BSB 0086LA |

| Product Description | Volume | Catalog # |
|---|-----------|-----------|
| Mouse/Rabbit ImmunoDetector Biotin Link | 15.0 ml | BSB 0001L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 50.0 ml | BSB 0003L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 100.0 ml | BSB 0005L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 200.0 ml | BSB 0007L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 1000.0 ml | BSB 0009L |




IHC of p63 on a FFPE Skin Tissue

ImmunoDetector HRP


For the Immunohistochemical detection of antigens in cells and formalin-fixed or frozen tissues

The Bio SB ImmunoDetector is a highly sensitive Biotin-Streptavidin detection system that allows for the demonstration of antigens in formalin-fixed paraffin-embedded tissues, cryostat sections, cytospreads, and cell preparations. The ImmunoDetector Detection System is a universal detection system which is suitable for use with mouse or rabbit primary antibodies. The increased sensitivity of the ImmunoDetector Detection System allows for rapid staining procedures without compromising the quality of stains.



- Biotin-Streptavidin, 2-Step Immunohistochemistry Detection Technology
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed or frozen tissues
- Universal: Detects Mouse or Rabbit antibodies
- It can be used with DAB HRP Brown, AEC HRP Red, HRP Black, HRP Green and HRP Blue substrate chromogen.
- For *in Vitro* Diagnostic Use

Recommended Protocol



ImmunoDetector HRP Detection Systems

| Product Description | Volume | Catalog # |
|---|-----------|-----------|
| Mouse/Rabbit ImmunoDetector DAB HRP Brown | 5 ml | BSB 0001S |
| Mouse/Rabbit ImmunoDetector DAB HRP Brown | 15.0 ml | BSB 0001 |
| Mouse/Rabbit ImmunoDetector DAB HRP Brown | 50.0 ml | BSB 0003 |
| Mouse/Rabbit ImmunoDetector DAB HRP Brown | 100.0 ml | BSB 0005 |
| Mouse/Rabbit ImmunoDetector DAB HRP Brown | 200.0 ml | BSB 0007 |
| Mouse/Rabbit ImmunoDetector DAB HRP Brown | 1000.0 ml | BSB 0009 |
| Mouse/Rabbit ImmunoDetector AEC HRP Red | 5 ml | BSB 0002S |
| Mouse/Rabbit ImmunoDetector AEC HRP Red | 15.0 ml | BSB 0002 |
| Mouse/Rabbit ImmunoDetector AEC HRP Red | 50.0 ml | BSB 0004 |
| Mouse/Rabbit ImmunoDetector AEC HRP Red | 100.0 ml | BSB 0006 |
| Mouse/Rabbit ImmunoDetector AEC HRP Red | 200.0 ml | BSB 0008 |
| Mouse/Rabbit ImmunoDetector AEC HRP Red | 1000.0 ml | BSB 0010 |


ImmunoDetector HRP Link & Label

| Product Description | Volume | Catalog # |
|---|-----------|-----------|
| Mouse/Rabbit ImmunoDetector Biotin Link | 15.0 ml | BSB 0001L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 50.0 ml | BSB 0003L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 100.0 ml | BSB 0005L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 200.0 ml | BSB 0007L |
| Mouse/Rabbit ImmunoDetector Biotin Link | 1000.0 ml | BSB 0009L |

| Product Description | Volume | Catalog # |
|---------------------------------------|-----------|-----------|
| Mouse/Rabbit ImmunoDetector HRP Label | 15.0 ml | BSB 0001H |
| Mouse/Rabbit ImmunoDetector HRP Label | 50.0 ml | BSB 0003H |
| Mouse/Rabbit ImmunoDetector HRP Label | 100.0 ml | BSB 0005H |
| Mouse/Rabbit ImmunoDetector HRP Label | 200.0 ml | BSB 0007H |
| Mouse/Rabbit ImmunoDetector HRP Label | 1000.0 ml | BSB 0009H |


ImmunoDetector HRP Link & Label

| Product Description | Volume | Catalog # |
|---|-----------|------------|
| Mouse/Rabbit ImmunoDetector Biotin Link & HRP Label | 15.0 ml | BSB 0001LH |
| Mouse/Rabbit ImmunoDetector Biotin Link & HRP Label | 50.0 ml | BSB 0003LH |
| Mouse/Rabbit ImmunoDetector Biotin Link & HRP Label | 100.0 ml | BSB 0005LH |
| Mouse/Rabbit ImmunoDetector Biotin Link & HRP Label | 200.0 ml | BSB 0007LH |
| Mouse/Rabbit ImmunoDetector Biotin Link & HRP Label | 1000.0 ml | BSB 0009LH |



IHC of CD10 on a FFPE Kidney Tissue








Complete Multiplex Detection Systems for IHC

Simple, Sensitive & Simultaneous Target Detection

*A selection of highly sensitive micropolymer based complete multiplex detection systems.
For use in clinical and research applications.*



Prostate Intraepithelial Neoplasia (PIN) MultiDetector HRP/AP Kit



CK 34BE12, p63 & AMACR

The Prostate Intraepithelial Neoplasia (PIN) MultiDetector HRP/AP Kit is a triple stain designed to detect prostate cancer in situ in the prostate glands. The basal cell indicators (p63 in DAB brown and CK34BE12 in Scarlet) stain non-affected glands and the AMACR in Green shows prostate tissue that may be affected by carcinoma. Prostate Intraepithelial Neoplasia (PIN) is a pre-cancerous condition of the prostate glands with a high predictive value for adenocarcinoma. An estimated third of men over 50 have a latent form of PIN, which could develop into a higher grade and eventually malignant carcinoma. High Grade (≥ 2) PIN (HGPIN) has morphological and genetic similarity to prostate adenocarcinomas, and adenocarcinoma presence and multifocality may be associated with the size and number of HGPIN foci.

- p63 and CK34BE12 are expressed in basal cells of healthy prostate glands, but absent in adenocarcinoma
- Alpha-methylacyl CoA racemase (AMACR) in the glands has been found at significantly higher near prostate carcinoma foci and can help distinguish PIN from benign lesions.
- Non-Biotin, Multidetector Immunohistochemistry Detection Technology
- Micro-polymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed tissues.
- For *in Vitro* Diagnostic Use

PIN MultiDetector HRP/AP Detection Systems

| Product Description | Volume | Catalog # |
|---------------------|-----------|-------------|
| PIN MultiDetector | 70 tests | BSB-0352-7 |
| PIN MultiDetector | 150 tests | BSB-0352-15 |
| PIN MultiDetector | 500 tests | BSB-0352-50 |




| Components | Volume | | |
|--|------------------------------------|---------------------------------------|---------------------------------------|
| | BSB-0352-7 7 ml Kit 70 tests | BSB-0352-15 15 ml Kit 150 tests | BSB-0352-50 50 ml Kit 500 tests |
| 20X ImmunoRetriever with Citrate | 50 mL | 100 mL | 200 mL |
| MultiDetector Peroxidase Blocker | 15 mL | 30 mL | 100 mL |
| MultiDetector AP Blocker | 15 mL | 30 mL | 100 mL |
| p63 MMab | 7mL | 15 mL | 50 mL |
| CK34BE12 MMab & AMACR RMa cocktail | 7mL | 15 mL | 50 mL |
| Mouse MultiDetector HRP Label | 7mL | 15 mL | 50 mL |
| Mouse MultiDetector AP & Rabbit MultiDetector HRP Label cocktail | 7mL | 15 mL | 50 mL |
| MultiDetector DAB Buffer-Substrate | 15 mL | 30 mL | 100 mL |
| MultiDetector DAB Chromogen | 1 mL | 2 mL | 6 mL |
| MultiDetector ALK Scarlet Buffer-Substrate | 15 mL | 30 mL | 100 mL |
| MultiDetector Scarlet Chromogen | 1 mL | 2 mL | 6 mL |
| MultiDetector HRP Green Buffer-Substrate | 15 mL | 30 mL | 100 mL |
| MultiDetector HRP Green Chromogen | 1 mL | 2 mL | 6 mL |
| 20X ImmunoDNA Washer | 50 mL | 100 mL | 200 mL |
| Fast ChromoProtector | 10 mL | 20 mL | 50 mL |
| Control Slides | 5 slides | 10 slides | 30 slides |


10 tests per ml considering 100µl per tissue

Prostate Intraepithelial Neoplasia (PIN) MultiDetector HRP/AP Kit

CK 34BE12, p63 & AMACR



Recommended Protocol



PNI Carcinoma MultiDetector HRP/AP Kit


CK 5/6 & NGFR

The PNI Carcinoma MultiDetector HRP/AP Kit is a dual stain that allows for the simultaneous visualization of skin carcinomas and nerve tissue. In cutaneous squamous and basal cell carcinoma, Perineural Invasion (PNI) is the infiltration of tumor within the perineural space. PNI is an uncommon manifestation of SCC and BCC but can indicate adverse outcomes including recurrence, metastasis, poor prognosis, and death. This kit has been optimized with MultiDetector HRP and AP Labels and contrasting chromogens to clearly differentiate tumor cells expressing High Molecular Weight cytokeratin (CK 5/6 in Scarlet) from nerve cells expressing Nerve Growth Factor Receptor (NGFR in Brown).


- IHC may be more accurate than routine sections in diagnosing PNI, with significant implications for patient staging, prognosis, and management
- Multiplex IHC may further improve PNI detection by dual staining for cytokeratins (CK AE1/AE3, CK 5/6) and nerve markers (NGFR, SOX-10, NF, etc.)
- Detects Squamous and Basal Cell Carcinomas and nerve bundles using FFPE tissues or frozen Mohs tissue sections
- Non-Biotin, Multidetector Immunohistochemistry Detection Technology
- Micro-polymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed tissues
- For *in Vitro* Diagnostic Use

PNI Carcinoma MultiDetector HRP/AP Detection Systems

| Product Description | Volume | Catalog # |
|-----------------------------|-----------|-------------|
| PNI Carcinoma MultiDetector | 70 tests | BSB-0353-7 |
| PNI Carcinoma MultiDetector | 150 tests | BSB-0353-15 |
| PNI Carcinoma MultiDetector | 500 tests | BSB-0353-50 |



PNI Carcinoma MultiDetector HRP/AP Kit




PNI Carcinoma MultiDetector HRP/AP Kit



| Components | Volume | | |
|--|------------------------------------|---------------------------------------|---------------------------------------|
| | BSB-0353-7 7 ml Kit 70 tests | BSB-0353-15 15 ml Kit 150 tests | BSB-0353-50 50 ml Kit 500 tests |
| ImmunoDNA Retriever with Citrate 20X | 50 mL | 100 mL | 200 mL |
| MultiDetector Peroxidase Blocker | 7 mL | 15 mL | 100 mL |
| MultiDetector AP Blocker | 7 mL | 15 mL | 100 mL |
| NGFR MMab & CK 5/6 RMab cocktail | 7 mL | 15 mL | 50 mL |
| Mouse MultiDetector HRP & Rabbit MultiDetector AP Label cocktail | 7 mL | 15 mL | 50 mL |
| MultiDetector DAB Buffer-Substrate | 7 mL | 15 mL | 100 mL |
| MultiDetector DAB Chromogen | 1 mL | 2 mL | 6 mL |
| MultiDetector ALK Scarlet Buffer-Substrate | 7 mL | 15 mL | 100 mL |
| MultiDetector Scarlet Chromogen | 1 mL | 2 mL | 6 mL |
| ImmunoDNA Washer 20X | 50 mL | 100 mL | 200 mL |
| Fast ChromoProtector | 10 mL | 20 mL | 50 mL |
| Control Slides | 5 slides | 10 slides | 30 slides |

PNI Carcinoma MultiDetector HRP/AP Kit

CK 5/6 & NGFR



Recommended Protocol



HPV MultiDetector HRP/AP Kit


p16 & Ki-67

The HPV MultiDetector HRP/AP Kit dual stain includes a cocktail of Ki67 and p16 antibodies to form a sensitive and specific test for significant cervical lesions. Cervical Intraepithelial Neoplasia (CIN) is the precursor to cervical cancer, the second most common malignancy among women. Colocalized Ki67 and p16 help identify CIN lesions for proper diagnosis and prevention of under- or over-treatment. Detection of nuclear Ki67 with DAB chromogen and nuclear and cytoplasmic p16 with ALK Scarlet provides an easy visualization of the colocalized antibodies for efficient diagnosis on biopsies and liquid cytology samples.


- p16 is a tumor-suppressor gene commonly used to detect squamous and glandular neoplasia in HPV-positive patients
- Ki67 is a nuclear cell proliferation marker, applied to diagnose the grade of abnormal cell growth.
- Non-Biotin, Multidetector Immunohistochemistry Detection Technology
- Micro-polymer detection technology allows for better cell penetration to deliver a highly specific and sensitive signal
- Ready-to-Use, High Sensitivity System especially designed for Immunohistochemistry of formalin-fixed tissues.
- For *in Vitro* Diagnostic Use

HPV MultiDetector HRP/AP Detection Systems

| Product Description | Volume | Catalog # |
|---------------------|-----------|-------------|
| HPV MultiDetector | 70 tests | BSB-0354-7 |
| HPV MultiDetector | 150 tests | BSB-0354-15 |
| HPV MultiDetector | 500 tests | BSB-0354-50 |



HPV MultiDetector HRP/AP Kit




HPV MultiDetector HRP/AP Kit



| Components | Volume | | |
|--|------------------------------------|---------------------------------------|---------------------------------------|
| | BSB-0354-7 7 ml Kit 70 tests | BSB-0354-15 15 ml Kit 150 tests | BSB-0354-50 50 ml Kit 500 tests |
| ImmunoDNA Retriever with Citrate 20X | 50 mL | 100 mL | 200 mL |
| MultiDetector Peroxidase Blocker | 7 mL | 15 mL | 50 mL |
| MultiDetector AP Blocker | 7 mL | 15 mL | 50 mL |
| Ki67 RMab & p16 MMab Cocktail | 7mL | 15 mL | 50 mL |
| Mouse MultiDetector AP & Rabbit MultiDetector HRP Label cocktail | 7mL | 15 mL | 50 mL |
| MultiDetector DAB Buffer-Substrate | 15 mL | 2 x 15 mL | 100 mL |
| MultiDetector DAB Chromogen | 1 mL | 2 mL | 6 mL |
| MultiDetector ALK Scarlet Buffer-Substrate | 7 mL | 15 mL | 50 mL |
| MultiDetector Scarlet Chromogen | 1 mL | 2 mL | 6 mL |
| ImmunoDNA Washer 20X | 50 mL | 100 mL | 200 mL |
| Fast ChromoProtector | 10 mL | 20 mL | 50 mL |
| Control Slides: 7-core HPV CLMA | 5 slides | 10 slides | 30 slides |

HPV MultiDetector HRP/AP Kit

p16 & Ki-67



Recommended Protocol






Substrate-Chromogen Systems


11 Colors Available!

*For Use In AP or HRP Based Detection Systems.
For use in clinical and research applications.*




Substrate-Chromogen Systems

for use with HRP Detection Systems




IVD

IHC of S100 Beta on an FFPE Melanoma Tissue with AEC



IVD

IHC of CD20 on an FFPE Colon Tissue HRP Green



IVD

IHC of CK cocktail on an FFPE Skin Tissue with HRP Blue Plus

- For *in Vitro* Diagnostic Use
- High Sensitivity, Low Background

- AEC HRP Red supplied as Ready-to-Use Format
- DAB Brown, HRP Black, HRP Green and HRP Blue supplied as two components

PolyDetector HRP Yellow kit

| | |
|-----------|---------------|
| 15.0 ml | BSB 0365-15 |
| 50.0 ml | BSB 0365-50 |
| 100.0 ml | BSB 0365-100 |
| 200.0 ml | BSB 0365-200 |
| 500.0 ml | BSB 0365-500 |
| 1000.0 ml | BSB 0365-1000 |

PolyDetector Liquid DAB HRP Brown Kit

| | |
|-----------|-----------|
| 15.0 ml | BSB 0015 |
| 50.0 ml | BSB 0016 |
| 100.0 ml | BSB 0017 |
| 200.0 ml | BSB 0018 |
| 500.0 ml | BSB 0018A |
| 1000.0 ml | BSB 0018B |

PolyDetector HRP Black Kit

| | |
|-----------|---------------|
| 15.0 ml | BSB 0361-15 |
| 50.0 ml | BSB 0361-50 |
| 100.0 ml | BSB 0361-100 |
| 200.0 ml | BSB 0361-200 |
| 500.0 ml | BSB 0361-500 |
| 1000.0 ml | BSB 0361-1000 |

PolyDetector HRP Fuchsia kit

| | |
|-----------|---------------|
| 15.0 ml | BSB 0364-15 |
| 50.0 ml | BSB 0364-50 |
| 100.0 ml | BSB 0364-100 |
| 200.0 ml | BSB 0364-200 |
| 500.0 ml | BSB 0364-500 |
| 1000.0 ml | BSB 0364-1000 |

PolyDetector HRP Green Kit

| | |
|-----------|----------|
| 15.0 ml | BSB 0128 |
| 50.0 ml | BSB 0129 |
| 100.0 ml | BSB 0130 |
| 200.0 ml | BSB 0131 |
| 500.0 ml | BSB 0132 |
| 1000.0 ml | BSB 0133 |

PolyDetector Liquid AEC HRP Red Ready-To-Use

| | |
|-----------|-----------|
| 15.0 ml | BSB 0011 |
| 50.0 ml | BSB 0012 |
| 100.0 ml | BSB 0013 |
| 200.0 ml | BSB 0014 |
| 500.0 ml | BSB 0061A |
| 1000.0 ml | BSB 0061 |

PolyDetector HRP Blue Plus Kit


| | |
|-----------|---------------|
| 15.0 ml | BSB 0363-15 |
| 50.0 ml | BSB 0363-50 |
| 100.0 ml | BSB 0363-100 |
| 200.0 ml | BSB 0363-200 |
| 500.0 ml | BSB 0363-500 |
| 1000.0 ml | BSB 0363-1000 |

PolyDetector HRP Mulberry kit

| | |
|-----------|---------------|
| 15.0 ml | BSB 0362-15 |
| 50.0 ml | BSB 0362-50 |
| 100.0 ml | BSB 0362-100 |
| 200.0 ml | BSB 0362-200 |
| 500.0 ml | BSB 0362-500 |
| 1000.0 ml | BSB 0362-1000 |


Substrate-Chromogen Systems

for use with AP Detection Systems




IVD

IHC of CD10 on an
FFPE Kidney Tissue with ALK Magenta



IVD

IHC of Her2 on an
FFPE Breast Carcinoma Tissue with ALK Blue



IVD

IHC of CD20 on an
FFPE Spleen Tissue with ALK Magenta

- For *in Vitro* Diagnostic Use
- High Sensitivity, Low Background
- Environmentally friendly: No solvents used
- ALK Blue and ALK Red supplied as Ready-to-Use Formats
- ALK Scarlet supplied as two components, ALK Magenta supplied as three components

PolyDetector Alk Blue Red Ready-To-Use

| | |
|-----------|----------|
| 15.0 ml | BSB 0062 |
| 50.0 ml | BSB 0063 |
| 100.0 ml | BSB 0064 |
| 200.0 ml | BSB 0065 |
| 1000.0 ml | BSB 0066 |

PolyDetector Alk Magenta kit

| | |
|-----------|----------|
| 15.0 ml | BSB 0077 |
| 50.0 ml | BSB 0078 |
| 100.0 ml | BSB 0079 |
| 200.0 ml | BSB 0080 |
| 1000.0 ml | BSB 0081 |

PolyDetector Alk Red Ready-To-Use

| | |
|-----------|----------|
| 15.0 ml | BSB 0067 |
| 50.0 ml | BSB 0068 |
| 100.0 ml | BSB 0069 |
| 200.0 ml | BSB 0070 |
| 1000.0 ml | BSB 0071 |

PolyDetector Alk Scarlet kit













| | |
|-----------|----------|
| 15.0 ml | BSB 0138 |
| 50.0 ml | BSB 0139 |
| 100.0 ml | BSB 0140 |
| 200.0 ml | BSB 0141 |
| 1000.0 ml | BSB 0142 |

Compatibility Guide

The Bio SB compatibility guide will help ensure that the correct chromogen, counter-stain, and mounting media is used in any test performed.

KEY

+ = compatible
- = incompatible
A = Aqua Mounter
P = Perma Mounter

| |  DAB HRP Brown |  AEC HRP Red |  HRP Black |  HRP Green |  HRP Blue Plus |  HRP Fuschia |  HRP Mulberry |  HRP Yellow |  ALK Magenta |  ALK Scarlet |
|--|---|---|---|---|---|---|--|--|---|---|
|  Methyl Green | + | - | + | - | +/- | - | - | - | + | + |
|  Hematoxylin | + | + | +/- | +/- | - | + | + | + | + | + |
|  Nuclear Fast Red | +/- | - | + | + | + | - | - | - | - | - |
| Mounting Media | P | A | P | P | A | P | A | A/P | A/P | A/P |


IVD For *in Vitro* Diagnostic Use



Deparaffinization & Retrieval
TintoDeparaffinator Citrate & EDTA
TintoRetriever Pressure cooker




Automated IHC & ISH
TintoStainer Plus
TintoDetector



Antibodies & Probes
TintoAntibodies




Tissue & Cell Line Microarrays



Detection Systems
PolyDetector
Micropolymer IgG Fab'



Solvent Free Mounting
ChromoProtector
XyGreen PermaMOUNTER



Environmentally Friendly Ancillaries for Immunohistochemistry

Solvent Free Ancillaries for IHC & ISH

*Safe, efficient and economical alternatives to traditional solvent containing solutions.
Reduce time from deparafinization to mounting.*



TintoDeparaffinator Citrate and EDTA

Safe, efficient and economical alternatives to traditional deparaffinization. These solvent-free products reduce exposure to toxic solvents like xylenes, toluene and alcohols when handling FFPE tissues for molecular pathology.

Non-Toxic and Biodegradable!


Reduction of up to 1 hour in the Deparaffinization and Mounting Time used with Xylenes and Alcohol!



| Product Description | Volume | Catalog # |
|-----------------------------------|-----------------|-----------|
| TintoDeparaffinator Citrate 20X | 100 mL (2 Lt) | BSB 0175 |
| TintoDeparaffinator Citrate 20X | 1000 mL (20 Lt) | BSB 0176 |
| TintoDeparaffinator EDTA 20X | 100 mL (2 Lt) | BSB 0177 |
| TintoDeparaffinator EDTA 20X | 1000 mL (20 Lt) | BSB 0178 |
| TintoDeparaffinator Hot Rinse 20X | 100 mL (2 Lt) | BSB 0179 |
| TintoDeparaffinator Hot Rinse 20X | 1000 mL (20 Lt) | BSB 0180 |

ChromoProtector and XyGreen PermaMounter

The Bio SB ChromoProtector helps to dissolve micro-paraffin leftovers and preserve stains of substrate-chromogens that are dissolved or sensitive to fading, like AEC, FastRed, HRP Blue, HRP Green, ALK Magenta, ALK Scarlet, ALK Red and ALK Brown, when mounted with solvent-based permanent mounting media. The Bio SB ChromoProtector is also suitable to mount substrate-chromogens such as DAB, HRP Black, ALK Blue or others that are resistant to solvents. The XyGreen PermaMounter is a solvent free, non-toxic and biodegradable permanent mounting media that allows permanent mounting of tissues for molecular pathology applications.



IHC of **ALK-1** on a FFPE Large Anaplastic Lymphoma with **AEC** PermaMounted



IHC of **SATB2** on a FFPE Colon Carcinoma with **HRP Blue** PermaMounted

Non-Toxic and Solvent Free!

| Product Description | Volume | Catalog # |
|----------------------|--------|-----------|
| XyGreen PermaMounter | 15 mL | BSB 0169 |
| XyGreen PermaMounter | 50 mL | BSB 0170 |
| XyGreen PermaMounter | 100 mL | BSB 0171 |
| ChromoProtector | 15 mL | BSB 0151 |
| ChromoProtector | 50 mL | BSB 0152 |
| ChromoProtector | 100 mL | BSB 0153 |

Biodegradable!

Featured IHC Ancillaries

Immuno/DNA Retriever with Citrate

The Bio SB Immuno/DNA Retriever with Citrate can be used in IHC, ICC, IF, and ISH, resulting in standardization of the pretreatment procedure to produce consistent & reliable results.

| Product Description | Volume | Catalog # |
|----------------------------------|---------|-----------|
| ImmunoRetriever 20X with Citrate | 50mL | BSB 0023 |
| ImmunoRetriever 20X with Citrate | 200 mL | BSB 0020 |
| ImmunoRetriever 20X with Citrate | 500 mL | BSB 0021 |
| ImmunoRetriever 20X with Citrate | 1000 mL | BSB 0022 |



Immuno/DNA Retriever with EDTA

The Bio SB Immuno/DNA Retriever with EDTA can be used in IHC, ICC, IF, and ISH Applications. The Bio SB Immuno/DNA Retriever with EDTA is recommended for tissues that have lower concentrations of antigens or when using low affinity antibodies.

| Product Description | Volume | Catalog # |
|-------------------------------|---------|-----------|
| ImmunoRetriever 20X with EDTA | 50 mL | BSB 0033 |
| ImmunoRetriever 20X with EDTA | 200 mL | BSB 0030 |
| ImmunoRetriever 20X with EDTA | 500 mL | BSB 0031 |
| ImmunoRetriever 20X with EDTA | 1000 mL | BSB 0032 |



Antibody Diluent/Protein Blocker

The Bio SB Antibody Diluent/Protein Blocker is designed to dilute & stabilize antibodies while minimizing non specific reactions and promoting specific antigen-antibody binding.

| Product Description | Volume | Catalog # |
|---|---------|-----------|
| ImmunoDetector Protein Blocker / Antibody Diluent | 15 mL | BSB 0113 |
| ImmunoDetector Protein Blocker / Antibody Diluent | 50 mL | BSB 0040 |
| ImmunoDetector Protein Blocker / Antibody Diluent | 100 mL | BSB 0041 |
| ImmunoDetector Protein Blocker / Antibody Diluent | 200 mL | BSB 0114 |
| ImmunoDetector Protein Blocker / Antibody Diluent | 1000 mL | BSB 0115 |



Immuno/DNA Washer

The Bio SB Immuno/DNA Washer is a cost effective washing solution used in-between steps common in IHC, ICC, and IF. The Bio SB Immuno/DNA Washer can be incorporated in capillary gap based or automated IHC protocols.

| Product Description | Volume | Catalog # |
|-----------------------|---------|-----------|
| Immuno/DNA Washer 10X | 200 mL | BSB 0029 |
| Immuno/DNA Washer 10X | 1000 mL | BSB 0042 |
| Immuno/DNA Washer 20X | 200 mL | BSB 0149 |
| Immuno/DNA Washer 20X | 1000 mL | BSB 0150 |



IHC Ancillaries

Enzyme Retrieval

| Product Description | Volume | Catalog # |
|----------------------|---------|-----------|
| ImmunoDNA Digestor | 15 mL | BSB 0108 |
| ImmunoDNA Digestor | 50 mL | BSB 0109 |
| ImmunoDNA Digestor | 100 mL | BSB 0110 |
| ImmunoDNA Digestor | 200 mL | BSB 0111 |
| ImmunoDNA Digestor | 1000 mL | BSB 0112 |
| Mohs Immuno Digestor | 15 mL | BSB 0324 |
| Mohs Immuno Digestor | 50 mL | BSB 0325 |
| Mohs Immuno Digestor | 100 mL | BSB 0326 |

Negative Controls

| Product Description | Volume | Catalog # |
|-------------------------|--------|-----------|
| Mouse Negative Control | 3 mL | BSB 0040A |
| Mouse Negative Control | 6 mL | BSB 0040B |
| Mouse Negative Control | 15 mL | BSB 0040C |
| Rabbit Negative Control | 3 mL | BSB 0041A |
| Rabbit Negative Control | 6 mL | BSB 0041B |
| Rabbit Negative Control | 15 mL | BSB 0041C |

Counterstainers

| Product Description | Volume | Catalog # |
|---------------------------------|---------|---------------|
| Hematoxylin Counterstainer | 15 ml | BSB 0024 |
| Hematoxylin Counterstainer | 50 ml | BSB 0025 |
| Hematoxylin Counterstainer | 100 ml | BSB 0026 |
| Hematoxylin Counterstainer | 200 ml | BSB 0027 |
| Hematoxylin Counterstainer | 1000 ml | BSB 0028 |
| Nuclear Fast Red Counterstainer | 15 ml | BSB 0116 |
| Nuclear Fast Red Counterstainer | 50 ml | BSB 0117 |
| Nuclear Fast Red Counterstainer | 100 ml | BSB 0118 |
| Nuclear Fast Red Counterstainer | 200 ml | BSB 0119 |
| Nuclear Fast Red Counterstainer | 1000 ml | BSB 0120 |
| TintoHematoxylin Automation | 15 mL | BSB 0127 |
| TintoHematoxylin Automation | 50 mL | BSB 0182 |
| TintoHematoxylin Automation | 100 mL | BSB 0183 |
| TintoHematoxylin Automation | 200 mL | BSB 0184 |
| TintoHematoxylin Automation | 1000 mL | BSB 0185 |
| Pink Hematoxylin Counterstainer | 15 mL | BSB-0360-15 |
| Pink Hematoxylin Counterstainer | 50 mL | BSB-0360-50 |
| Pink Hematoxylin Counterstainer | 100 mL | BSB-0360-100 |
| Pink Hematoxylin Counterstainer | 200 mL | BSB-0360-200 |
| Pink Hematoxylin Counterstainer | 1000 mL | BSB-0360-1000 |

Other Mounting Solutions

| Product Description | Volume | Catalog # |
|-------------------------|---------|---------------|
| AquaMounter | 15 ml | BSB 0090 |
| AquaMounter | 50 ml | BSB 0091 |
| AquaMounter | 100 ml | BSB 0092 |
| AquaMounter | 500 ml | BSB 0093 |
| PermaMounter | 118 ml | BSB 0097 |
| Fast ChromoProtector | 15 ml | BSB-0327-15 |
| Fast ChromoProtector | 50 ml | BSB-0327-50 |
| Fast ChromoProtector | 100 ml | BSB-0327-100 |
| Fast ChromoProtector | 200 ml | BSB-0327-200 |
| Fast ChromoProtector | 500 ml | BSB-0327-500 |
| Fast ChromoProtector | 1000 ml | BSB-0327-1000 |
| FluoroMounter | 15 mL | BSB 0157 |
| FluoroMounter | 50 mL | BSB 0158 |
| FluoroMounter | 100 mL | BSB 0159 |
| FluoroMounter with DAPI | 15 mL | BSB 0163 |
| FluoroMounter with DAPI | 50 mL | BSB 0164 |
| FluoroMounter with DAPI | 100 mL | BSB 0165 |

IHC Ancillaries

Washers & Detergents

| Product Description | Volume | Catalog # |
|-----------------------|---------|-----------|
| Immuno/DNA Washer 10X | 200 ml | BSB 0029 |
| Immuno/DNA Washer 10X | 1000 ml | BSB 0042 |
| Immuno/DNA Washer 20X | 200 ml | BSB 0149 |
| Immuno/DNA Washer 20X | 1000 ml | BSB 0150 |
| Tween 20 | 100 ml | BSB 0045 |
| Tween 20 | 500 ml | BSB 0046 |

Slide Treatment Solutions

| Product Description | Volume | Catalog # |
|------------------------------------|---------|-----------|
| CytoLayer Slide Treatment Solution | 100 ml | BSB 4026 |
| CytoLayer Slide Treatment Solution | 500 ml | BSB 4027 |
| CytoLayer Slide Treatment Solution | 1000 ml | BSB 4032 |

Other Blocking Solutions


| Product Description | Volume | Catalog # |
|---------------------------------|---------|-----------|
| Background Blocker | 15 mL | BSB 0103 |
| Background Blocker | 50 ml | BSB 0104 |
| Background Blocker | 100 ml | BSB 0105 |
| Background Blocker | 200 ml | BSB 0106 |
| Background Blocker | 1000 ml | BSB 0107 |
| PolyDetector Peroxidase Blocker | 15 ml | BSB 0050 |
| PolyDetector Peroxidase Blocker | 50 ml | BSB 0051 |
| PolyDetector Peroxidase Blocker | 100 ml | BSB 0052 |
| PolyDetector Peroxidase Blocker | 200 ml | BSB 0053 |
| PolyDetector Peroxidase Blocker | 1000 ml | BSB 0054 |
| PolyDetector AP Blocker | 15 ml | BSB 0055 |
| PolyDetector AP Blocker | 50 ml | BSB 0056 |
| PolyDetector AP Blocker | 100 ml | BSB 0057 |
| PolyDetector AP Blocker | 200 ml | BSB 0058 |
| PolyDetector AP Blocker | 1000 ml | BSB 0059 |
| ImmunoDetector Biotin Blocker | 15 ml | BSB 0098 |
| ImmunoDetector Biotin Blocker | 50 ml | BSB 0099 |
| ImmunoDetector Biotin Blocker | 100 ml | BSB 0100 |
| ImmunoDetector Biotin Blocker | 200 ml | BSB 0101 |
| ImmunoDetector Biotin Blocker | 1000 ml | BSB 0102 |





Tissue and Cell Line Microarrays

Tissue Microarrays | Infectious Disease Arrays | IHC & ISH Validated




Tissue Microarrays and Control Slides

As the diagnostics market continues to grow, researchers and clinicians have a greater need for a wide variety of high quality and cost effective control slides. Control slides are invaluable tools utilized by institutions when validating reagents, qualifying new products, testing protocols or performing research which requires multiple tissue types. Bio SB control slides are cost effective, high quality tissues mounted on Hydrophilic Plus slides which are validated for use in immunohistochemical (IHC) and in situ hybridization (ISH) applications.

Tissue Microarray Features

- Easy method of antibody validation.
- Cost effective diagnostic control.
- Test a large number of tissue types on one slide.
- Available in 11 or 23 core format.
- Validated for use with over 650 antibodies used in Immunohistochemistry .



Above: Depiction of Normal Human 11 Core TMA

Normal Human Tissue Microarray (NH-TMA), Cancer Human Tissue Microarray (CH-TMA) and Cancer Human Cell Line Microarray (CH-CLMA)

The Normal Human Tissue Microarray, or NH-TMA, Cancer Human Tissue Microarray, or CH-TMA, and Cancer Human Cell Line Microarray, or CH-CLMA, are an excellent way to test and validate an antibody, ISH probes or other reagent on multiple tissues. The Bio SB NH-TMA's and CH-TMA's are available in both 11 or 23-core configurations. The Bio SB NH-TMA's, CH-TMA's and CH-CLMA's are an excellent way for clinics and research labs to save time and money by allowing multiple tissues to be tested on one slide.



BSB 0297 - 11 Core Normal Human Tissue Microarray (NH-TMA)

| | | | |
|---------------|--------------|-------------|---------------------|
| PL - Placenta | Blank | LV - Liver | TL - Tonsil |
| CL - Colon | SK - Skin | BRN - Brain | BRS - Breast |
| PR - Prostate | TH - Thyroid | KD - Kidney | FT - Fallopian Tube |

BSB 0298 - 23 Core Normal Human Tissue Microarray (NH-TMA)

| | | | | | |
|---------------|---------------|--------------|-----------------|------------------|-----------------------------|
| PL - Placenta | Blank | BRS - Breast | MY - Myometrium | CX - Cervix | FT - Fallopian Tube |
| BRN - Brain | ST - Stomach | AD - Adrenal | PC - Pancreas | SG - Salivary | CL - Colon |
| LV - Liver | KD - Kidney | TH - Thyroid | LG - Lung | SK - Skin | UR - Urothelial Carcinoma |
| TS - Testis | PR - Prostate | SP - Spleen | TL - Tonsil | BM - Bone Marrow | LL - Lymphoblastic Lymphoma |

BSB 0299 - 7 Core Normal Human Lymphoid Tissue Microarray


| | | | |
|-------------|-----------------------------|-----------------|-------------|
| TL - Tonsil | Blank | LN - Lymph Node | SP - Spleen |
| TL - Tonsil | LL - Lymphoblastic Lymphoma | LN - Lymph Node | SP - Spleen |

The maps above outline the various normal tissue types used. Each slide comes with a "blank" core for easy orientation & interpretation.

Tissue Microarrays and Control Slides

Cancer Human Tissue Microarrays (CH-TMA)

The Human Cancer Tissue Microarrays (CH-TMA) consist of 2 mm cores of cancer human formalin-fixed paraffin-embedded tissues which were assembled in array fashion to allow multiplex molecular pathology analysis and validation of reagents, or to be used as tissue controls for Immunohistochemistry and/or in-situ hybridization (CISH and FISH) applications.



BSB 0230 - 11 Core Cancer Human Tissue Microarray


The Human Cancer 11-core TMA contains the following human cancer tissues: Chronic Granulocytic Leukemia, Colon Cancer, Prostate Cancer, Squamous Cell Carcinoma, Papillary Thyroid Carcinoma, Ovarian Cancer, Bladder TCC, Melanoma, Testicular Cancer, Ductal Breast Carcinoma, and Cervical Cancer.

| Cancer Human Tissue Microarrays | Catalog # |
|---------------------------------|-------------|
| 11-Core Human Cancer TMA | BSB 0230 |
| 23-Core Human Cancer TMA | BSB 0231 |
| 2-Core Human PIN TMA | BSB-0333-CS |

* Tissue selection may vary from lot to lot

Cancer Human Cell Line Microarrays (CH-CLMA)

The Cancer Human Cell Line Microarrays (CH-CLM) consist of 2 mm cores of cancer human formalin-fixed paraffin-embedded cell lines which were assembled in array fashion to allow multiplex molecular pathology analysis and validation of reagents, or to be used as tissue controls for Immunohistochemistry, Immunocytochemistry and/or in-situ hybridization (CISH and FISH) applications.




Above: HER-2 Cell Line Microarray with various signal strengths.

| Cell Line Microarrays | Catalog # |
|---------------------------------|-------------|
| 3-Core ALK CLMA | BSB 0296 |
| 9- Core HER-2 CLMA | BSB 0292 |
| 7-Core ER/PR CLMA | BSB 0293 |
| 7-Core EGFR CLMA | BSB 0295 |
| 9-Core PTEN CLMA | BSB 0300 |
| 7-Core PD-L1 CLMA | BSB 0301 |
| 5-Core BRAF V600E CLMA | BSB 0305 |
| 11-Core Breast Cancer CLMA | BSB 0302 |
| 7-Core HPV Cervical Cancer CLMA | BSB 0294 |
| 3-Core Lung Cancer CLMA | BSB 0296 |
| 3-Core GIST CLMA | BSB 0242 |
| 5-Core Melanoma CLMA | BSB 0243 |
| 3-Core Neuroblastoma CLMA | BSB 0303 |
| 11-Core Immunotherapy CLMA | BSB 0304 |
| 31-Core Multi Cancer CLMA | BSB 0244 |
| 3-Core Androgen Receptor CLMA | BSB-0334-CS |
| 3-Core ROS1 CLMA | BSB-0335-CS |
| 3-Core IDH1 R132 CLMA | BSB-0336-CS |
| 4-Core MMR CLMA | BSB-0337-CS |

Tissue Microarrays and Control Slides

Infectious Disease Cell Line Microarray (ID-CLMA)

The Bio SB Infectious Disease cell line Microarray, or ID-CLMA is a simple and cost effective way to test and validate infectious disease markers by immunohistochemical (IHC) or in situ hybridization (ISH). The ID-Microarray is available in a 5-core or Individual virus configuration, and includes two areas for tissue mounting. All TMA's include negative controls to reduce interpretation error.




Above: 5-Core Multi-Infectious Cell Line Microarray (BSB 0307).


ID-CLMA

- Test multiple infectious disease markers.
- Optimized for IHC & ISH Applications.
- Cost effective solution.
- Test multiple samples at once.
- Two control tissue mounting areas for Manual or Automated IHC/ISH users.

| Infectious Disease CLMA's | Catalog # |
|--|-------------|
| Multi ID-Array (5-Core) | BSB 0307 |
| Adenovirus ID-Array (2-Core) | BSB 0310 |
| HHV-8 ID-Array (2-Core) | BSB 0309 |
| SV-40 ID-Array (2-Core) | BSB 0308 |
| Treponema TMA (2-Tissue) | BSB-0338-CS |
| Gram Positive & Negative (3-Tissue) | BSB-0339-CS |
| Fungus: Aspergillus, Candida, Histoplasma (4-Tissue) | BSB-0340-CS |









ZytoLight® FISH
Fluorescent in situ Hybridization




ZytoLight® FlexISH
Fluorescent in situ Hybridization



ZytoFast® Rapid CISH
Chromogenic in situ Hybridization



ZytoDot® CISH
Chromogenic in situ Hybridization



FISH & CISH Products

ZytoLight® FISH | ZytoLight® FlexISH
ZytoDot® CISH | ZytoFast® Rapid CISH

*A diverse selection of innovative, high quality and cost effective ISH probes, kits,
ancillaries and equipment*



ZytoLight® FISH

This image shows a fluorescence microscopy view of cells. The nuclei are stained blue with DAPI, and several bright green spots are visible within the nuclei, indicating the presence of specific DNA sequences.



ZytoLight® FlexISH

This image shows a fluorescence microscopy view of cells. The nuclei are stained blue with DAPI, and several bright green spots are visible within the nuclei, indicating the presence of specific DNA sequences.



ZytoFast® Rapid CISH

This image shows a brightfield microscopy view of cells. The nuclei are stained blue with hematoxylin, and several red spots are visible within the nuclei, indicating the presence of specific DNA sequences.



ZytoDot® CISH

This image shows a brightfield microscopy view of cells. The nuclei are stained purple with hematoxylin, and several brown spots are visible within the nuclei, indicating the presence of specific DNA sequences.

ZytoLight®

Reliable multi-target detection using Fluorescence in situ Hybridization!

ZytoLight® products for Fluorescence in situ Hybridization (FISH) are designed for the identification of genetic aberrations e.g. translocations, deletions, amplifications, and chromosomal aneuploidies associated with tumors and genetic diseases.

INTRODUCTION

ZytoLight® products are designed for the identification of genetic aberrations e.g. translocations, deletions, amplifications, and chromosomal aneuploidies by Fluorescence in situ Hybridization (FISH) in formalin-fixed, paraffin-embedded tissue sections, cell samples, blood or bone marrow smears, and metaphase chromosome spreads.

HIGH SENSITIVITY AND SPECIFICITY

ZytoLight® FISH probes are direct labeled using the unique ZytoLight® Direct Label System II providing improved signal intensity. All ZytoLight® single copy (SPEC™) probes are processed by the unique ZytoLight® Repeat Subtraction Technique resulting in advanced specificity and less background. No further blocking of repetitive sequences is needed!

ZYTOLIGHT® SPEC™ AND CENT™ PROBES

ZytoLight® SPEC™ probes hybridize to specific, single copy DNA sequences of the human genome. ZytoLight® SPEC™ probes are available for the detection of a variety of chromosomal aberrations associated with tumors and genetic diseases.

ZytoLight® CENT™ probes hybridize to highly repetitive human satellite DNA sequences located at the centromeric regions of chromosomes producing sharp, bright signals specific for each individual chromosome.

ZYTOLIGHT® KITS – CONVENIENT SOLUTIONS

For making FISH analysis reliable and user-friendly, all ZytoLight® FISH probes can be combined with the ZytoLight® FISH-Tissue Implementation Kit (Z-2028-20). Also available now, specific for FISH analyses on cytology specimens, the ZytoLight® FISH-Cytology Implementation Kit (Z-2099-20). Both Implementation Kits include all necessary pretreatment solutions, wash buffers and DAPI/Antifade solution and a detailed protocol to perform successful FISH experiments. Additionally, for some major targets, complete kits including probes and all necessary reagents are available.

The ZytoLight® system uses directly labeled FISH probes:

1. Eliminating the need to detect the probes with fluorophore-coupled antibodies.


The probes are detected by fluorescence microscopy using appropriate filter sets

2. Due to an exciter filter

3. Full-spectrum light, emitted by the microscope lamp

4. Is reduced to light of a defined wavelength that specifically excites the fluorophore of the probe. This light is reflected onto the specimen by a dichroic mirror

5. The fluorophore emits light of longer wavelengths that passes the mirror. Finally, a barrier filter 6. reduces the emitted light to a defined wavelength that can be detected.



ZytoLight® FISH - Solid Tumor Probes By Application

The ZytoLight line of products is designed for use in detecting translocations, deletions, amplifications and chromosomal aneuploidies by fluorescent in situ hybridization (FISH). ZytoLight products can be used on formalin-fixed paraffin-embedded tissue (FFPE), cell samples, blood/bone marrow smears and metaphase chromosome spreads.

Brain and Neural Tumor Probes

ZytoLights Glioma 1p/19q Probe Set

ZytoLight SPEC 1p36/1q25 Dual Color Probe

ZytoLight SPEC 19q13/19p13 Dual Color Probe

ZytoLight SPEC C19MC/19p13 Dual Color Probe

ZytoLight CDKN2A/CEN 9 Dual Color Probe

ZytoLight EGFR/CEN 7 Dual Color Probe

ZytoLight MET/CEN 7 Dual Color Probe

ZytoLight MYCN/2q11 Dual Color Probe

ZytoLight NTRK2 Dual Color Break Apart Probe

ZytoLight PTEN/CEN 10 Dual Color Probe

ZytoLight SPEC TERT Dual Color Break Apart Probe

ZytoLight TP53/17q22 Dual Color Probe

Cervical Cancer Probes

ZytoLight SPEC MYC/CEN 8 Dual Color Probe

ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe

ZytoLight SPEC TERC/CEN 3 Dual Color Probe

ZytoLight SPEC TERT/5q31 Dual Color Probe

Lung Cancer Probes

ZytoLight SPEC ALK/EML4 TriCheck Probe

ZytoLight SPEC ALK Dual Color Break Apart Probe

ZytoLight SPEC ALK/2q11 Dual Color Probe

ZytoLight SPEC BRAF/CEN 7 Dual Color Probe

ZytoLight SPEC CARS Dual Color Break Apart Probe

ZytoLight SPEC CD274, PDCD1LG2/CEN 9 Dual Color Probe

ZytoLight SPEC EGFR/CEN 7 Dual Color Probe

ZytoLight SPEC EML4 Dual Color Break Apart Probe

ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe

ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit

ZytoLight CEN 17/SPEC ERBB2 Dual Color Probe

ZytoLight SPEC ERBB2/D17S122 Dual Color Probe

ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe

ZytoLight SPEC FGFR2 Dual Color Break Apart Probe

ZytoLight SPEC FGFR2/CEN 10 Dual Color Probe

ZytoLight SPEC FGFR3 Dual Color Break Apart Probe

ZytoLight SPEC FGFR3/4p11 Dual Color Probe

ZytoLight SPEC KIF5B Dual Color Break Apart Probe

ZytoLight SPEC KRAS/CEN 12 Dual Color Probe

ZytoLight SPEC MET/CEN 7 Dual Color Probe

ZytoLight SPEC NRG1 Dual Color Break Apart Probe

ZytoLight SPEC NRG1/CD74 TriCheck™ Probe

ZytoLight SPEC NTRK1 Dual Color Break Apart Probe

ZytoLight SPEC NTRK2 Dual Color Break Apart Probe

ZytoLight SPEC NTRK3 Dual Color Break Apart Probe

ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe

ZytoLight SPEC RICTOR/5q31.1 Dual Color Probe

Breast Cancer Probes

ZytoLight SPEC BCL2L1/CEN 20 Dual Color Probe

ZytoLight SPEC CCND1/CEN 11 Dual Color Probe

ZytoLight SPEC EGFR/CEN 7 Dual Color Probe

ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe

ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit

ZytoLight CEN 17/SPEC ERBB2 Dual Color Probe

ZytoLight SPEC ERBB2/D17S122 Dual Color Probe

ZytoLight SPEC ERBB2/TOP2A/CEN 17 Triple Color Probe

ZytoLight SPEC ERBB3/CEN 12 Dual Color Probe

ZytoLight SPEC ERBB4/2q11 Dual Color Probe

ZytoLight SPEC ESR1/CEN 6 Dual Color Probe

ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe

ZytoLight SPEC FGFR2/CEN 10 Dual Color Probe

ZytoLight SPEC MCL1/1p12 Dual Color Probe


ZytoLight SPEC MYC/CEN 8 Dual Color Probe

ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe

ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe


ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe

ZytoLight® FISH - Solid Tumor Probes By Application




Lung Cancer Probes

| |
|--|
| ZytoLight SPEC RET Dual Color Break Apart Probe |
| ZytoLight SPEC ROS1 Dual Color Break Apart Probe |
| ZytoLight SPEC ROS1/CEN 6 Dual Color Probe |
| ZytoLight SPEC SOX2/CEN 3 Dual Color Probe |




Prostate Cancer Probes

| |
|---|
| ZytoLight SPEC ERG Dual Color Break Apart Probe |
| ZytoLight SPEC ERG/TMPRSS2 TriCheck™ Probe |
| ZytoLight SPEC PTEN/CEN 10 Dual Color Probe |
| ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe |




Renal Cell Carcinoma Probes

| |
|---|
| ZytoLight SPEC CCND1 SPEC Break Apart/ 2q11/CEN 6 Quadruple Color Probe |
| ZytoLight SPEC FHIT/CEN 3 Dual Color Probe |
| ZytoLight SPEC CDKN2A/ CEN 3/7/17 Quadruple Color Probe |
| ZytoLight SPEC TFE3 Dual Color Break Apart Probe |
| ZytoLight SPEC VHL/CEN 3 Dual Color Probe |
| ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe |




Sarcoma Probes

| |
|--|
| ZytoLight SPEC ALK Dual Color Break Apart Probe |
| ZytoLight SPEC CDK4/CEN 12 Dual Color Probe |
| ZytoLight SPEC CIC Dual Color Break Apart Probe |
| ZytoLight SPEC COL1A1 Dual Color Break Apart Probe |
| ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe |
| ZytoLight SPEC DDIT3 Dual Color Break Apart Probe |
| ZytoLight SPEC ETV6 Dual Color Break Apart Probe |
| ZytoLight SPEC EWSR1 Dual Color Break Apart Probe |
| ZytoLight SPEC EWSR1/FLI1 TriCheck Probe |
| ZytoLight SPEC FOXO1 Dual Color Break Apart Probe |
| ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe |
| ZytoLight SPEC FOXO1/PAX3 TriCheck Probe |
| ZytoLight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe |
| ZytoLight SPEC FUS Dual Color Break Apart Probe |
| ZytoLight SPEC JAZF1 Dual Color Break Apart Probe |




Sarcoma Probes

| |
|--|
| ZytoLight MDM2/CEN 12 Dual Color Probe |
| ZytoLight MDM4/1p12 Dual Color Probe |
| ZytoLight SPEC MYC/CEN 8 Dual Color Probe |
| ZytoLight SPEC NR4A3 Dual Color Break Apart Probe |
| ZytoLight SPEC NTRK3 Dual Color Break Apart Probe |
| ZytoLight SPEC PDGFB Dual Color Break Apart Probe |
| ZytoLight SPEC PHF1 Dual Color Break Apart Probe |
| ZytoLight SPEC SMARCB1/22q12 Dual Color Probe |
| ZytoLight SPEC SS18 Dual Color Break Apart Probe |
| ZytoLight SPEC SS18/SSX1 TriCheck Probe |
| ZytoLight SPEC TFE3 Dual Color Break Apart Probe |
| ZytoLight SPEC USP6 Dual Color Break Apart Probe |
| ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe |
| ZytoLight SPEC SPEC WT1 Dual Color Break Apart Probe |
| ZytoLight SPEC WWT1 Dual Color Break Apart Probe |
| ZytoLight SPEC YWHAE Dual Color Break Apart Probe |



Salivary Gland Sarcoma Probes

| |
|--|
| ZytoLight SPEC ETV6 Dual Color Break Apart Probe |
| ZytoLight SPEC EWSR1 Dual Color Break Apart Probe |
| ZytoLight SPEC MAML2 Dual Color Break Apart Probe |
| ZytoLight SPEC MYB Dual Color Break Apart Probe |
| ZytoLight SPEC NTRK3 Dual Color Break Apart Probe |
| ZytoLight SPEC NTRK3 Dual Color Break Apart Probe |
| ZytoLight SPEC NUTM1 Dual Color Break Apart Probe |
| ZytoLight SPEC SPEC WT1 Dual Color Break Apart Probe |



Gastrointestinal Cancer Probes

| |
|---|
| ZytoLight SPEC BRAF Dual Color Break Apart Probe |
| ZytoLight SPEC CCND1/CEN 11 Dual Color Probe |
| ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe |
| ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit |
| ZytoLight CEN 17/SPEC ERBB2 Dual Color Probe |
| ZytoLight SPEC ERBB2/D17S122 Dual Color Probe |
| ZytoLight SPEC KRAS/CEN 12 Dual Color Probe |
| ZytoLight SPEC MDM2/CEN 12 Dual Color Probe |
| ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe |

ZytoLight® FISH - Hematology Specific Probes By Application

Acute Lymphoblastic Leukemia (ALL) Probes

| |
|--|
| ZytoLights SPEC ABL1 Dual Color Break Apart Probe |
| ZytoLight SPEC ABL2 Dual Color Break Apart Probe |
| ZytoLight SPEC CRLF2 Dual Color Break Apart Probe |
| ZytoLight SPEC CSF1R Dual Color Break Apart Probe |
| ZytoLight SPEC ETV6 Dual Color Break Apart Probe |
| ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe |
| ZytoLight SPEC KMT2A Dual Color Break Apart Probe |
| ZytoLight SPEC MEF2D/BCL9 TriCheck™ Probe |
| ZytoLight SPEC MYB Dual Color Break Apart Probe |
| ZytoLight SPEC NUP214 Dual Color Break Apart Probe |
| ZytoLight SPEC PNUP98 Dual Color Break Apart Probe |
| ZytoLight SPEC PDGFRA/FIP1L1 TriCheck Probe |
| ZytoLight SPEC SPI1 Dual Color Break Apart Probe |
| ZytoLight SPEC ZNF384 Dual Color Break Apart Probe |

Acute Myelogenous Leukemia (AML) Probes


| |
|---|
| ZytoLight CEN 8 Probe |
| ZytoLight SPEC ABL2 Dual Color Break Apart Probe |
| ZytoLight SPEC CFBF Dual Color Break Apart Probe |
| ZytoLight SPEC CREBBP Dual Color Break Apart Probe |
| ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe |
| ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe |
| ZytoLight SPEC EGR1/5p15 Dual Color Probe |
| ZytoLight SPEC EGR1/D5S23, D5S721 Dual Color Probe |
| ZytoLight SPEC FGFR1 Dual Color Break Apart Probe |
| ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe |
| ZytoLight SPEC KMT2A Dual Color Break Apart Probe |
| ZytoLight SPEC NUP98 Dual Color Break Apart Probe |
| ZytoLight SPEC NUP214 Dual Color Break Apart Probe |
| ZytoLight SPEC PDGFRA/FIP1L1 TriCheck Probe |
| ZytoLight SPEC PDGFRB Dual Color Break Apart Probe |
| ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe |
| ZytoLight SPEC PTPRT/20q11 Dual Color Probe |
| ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe |

Chronic Lymphocytic Leukemia (CLL) Probes

| |
|--|
| ZytoLight SPEC BCL2 Dual Color Break Apart Probe |
| ZytoLight SPEC CCND1 Dual Color Break Apart Probe |
| ZytoLight SPEC CCND1/CEN 11 Dual Color Probe |
| ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe |
| ZytoLight SPEC D13S319/13q34 Dual Color Probe |
| ZytoLight SPEC MYB/CEN 6 Dual Color Probe |
| ZytoLight SPEC MYC/CEN 8 Dual Color Probe |
| ZytoLight SPEC RB1/13q12 Dual Color Probe |
| ZytoLight SPEC TP53/ATM Dual Color Probe |
| ZytoLight SPEC TP53/CEN 17 Dual Color Probe |


Chronic Myelogenous Leukemia (CML) Probes

| |
|--|
| ZytoLight CEN 8 Probe |
| ZytoLight SPEC ABL1 Dual Color Break Apart Probe |
| ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe |
| ZytoLight SPEC PDGFRB Dual Color Break Apart Probe |
| ZytoLight SPEC TP53/17q22 Dual Color Probe |



ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe
Bone marrow biopsy tissue section with translocation affecting the BCR/ABL1 loci as indicated by one separate orange signal, one separate green signal and two orange/green fusion signals.

ZytoLight® FISH - Hematology Specific Probes By Application




Non-Hodgkin Lymphoma Probes

| |
|---|
| ZytoLight SPEC 11q gain/loss Triple Color Probe |
| ZytoLight SPEC BCL2 Dual Color Break Apart Probe |
| ZytoLight SPEC BCL2/CEN 18 Dual Color Probe |
| ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe |
| ZytoLight SPEC BCL6 Dual Color Break Apart Probe |
| ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe |
| ZytoLight SPEC CCND1 Dual Color Break Apart Probe |
| ZytoLight SPEC CCND1/CEN 11 Dual Color Probe |
| ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe |
| ZytoLight SPEC FGFR3 Dual Color Break Apart Probe |
| ZytoLight SPEC IGK Dual Color Break Apart Probe |
| ZytoLight SPEC IGL Dual Color Break Apart Probe |
| ZytoLight SPEC IRF4,DUSP22 Dual Color Break Apart Probe |
| ZytoLight SPEC IGH Dual Color Break Apart Probe |
| ZytoLight SPEC MALT1 Dual Color Break Apart Probe |
| ZytoLight SPEC MYC Dual Color Break Apart Probe |
| ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe |



Multiple Myeloma Probes

| |
|---|
| ZytoLight CCND1 Dual Color Break Apart Probe |
| ZytoLight CCND1/CEN 11 Dual Color Probe |
| ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe |
| ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe |
| ZytoLight SPEC FGFR3 Dual Color Break Apart Probe |
| ZytoLight SPEC FGFR3/IGH Dual Color Dual Fusion Probe |
| ZytoLight SPEC IGH Dual Color Break Apart Probe |
| ZytoLight SPEC MAF/IGH Dual Color Dual Fusion Probe |
| ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe |
| ZytoLight SPEC RB1/13q12 Dual Color Probe |
| ZytoLight SPEC TP53/CEN 17 Dual Color Probe |




Myelodysplastic Syndrome (MDS) Probes

| |
|--|
| ZytoLight CEN 8 Probe |
| ZytoLight SPEC CREBBP Dual Color Break Apart Probe |
| ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe |
| ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe |
| ZytoLight SPEC EGR1/5p15 Dual Color Probe |
| ZytoLight SPEC EGR1/D5S23, D5S721 Dual Color Probe |
| ZytoLight SPEC ETV6 Dual Color Break Apart Probe |
| ZytoLight SPEC NUP98 Dual Color Break Apart Probe |
| ZytoLight SPEC NUP214 Dual Color Break Apart Probe |
| ZytoLight SPEC PDGFRB Dual Color Break Apart Probe |
| ZytoLight SPEC PTPRT/20q11 Dual Color Probe |
| ZytoLight SPEC TERT/5q31 Dual Color Probe |




RUO

ZytoLight SPEC BCL2/CEN 18 Dual Color Probe
Probe hybridized to normal interphase cells as indicated by two orange and two green signals in each nucleus.



RUO

ZytoLight SPEC TP53/CEN 17 Dual Color Probe
Probe hybridized to bone marrow tissue section with deletion of the TP53 gene as indicated by one orange signal and two green signals in each nucleus.



RUO

ZytoLight SPEC TERT/5q31 Dual Color Probe
Probe hybridized to melanoma tissue section showing normal cells as indicated by two green and two orange signals in each nucleus and cells with TERT gene amplification as indicated by multiple green signals per nucleus.


ZytoLight® FISH - Genetic Probes By Application

Sex Mismatched Bone-Marrow Transplantant Management Probes

| |
|---------------------------------------|
| ZytoLight CEN X Probe |
| ZytoLight CEN X/Y Dual Color Probe |
| ZytoLight CEN X/Yq12 Dual Color Probe |
| ZytoLight CEN Y (DYZ3) Probe |
| ZytoLight CEN Yq12 Probe |


Prenatal, Postnatal and Preimplantation Genetics Probes

| |
|---|
| ZytoLight Aneuploidy Panel 18/X/Y and 13/21 |
| ZytoLight Aneuploidy Panel X/Y and 13/18/21 |
| ZytoLight SPEC 13q12 Probe |
| ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe |
| ZytoLight SPEC 13/21 Dual Color Probe |
| ZytoLight CEN 18 Probe |
| ZytoLight SPEC 18/CEN X/Y Triple Color Probe |
| ZytoLight SPEC 21q22 Probe |
| ZytoLight SPEC 21/CEN X/Yq12 Triple Color Probe |
| ZytoLight CEN X Probe |
| ZytoLight CEN X/Y Dual Color Probe |
| ZytoLight CEN X/Yq12 Dual Color Probe |
| ZytoLight CEN Y (DYZ3) Probe |
| ZytoLight CEN Yq12 Probe |




RUO

ZytoLight SPEC NUP214 Dual Color Break Apart Probe
Probe hybridized to normal interphase cells as indicated by two orange/green fusion signals per nucleus and to metaphase chromosomes of a normal cell.




RUO

ZytoLight SPEC ERBB2/CEN17 Dual Color Probe
Probe hybridized to a breast carcinoma tissue section, ERBB2 gene cluster (green), CEN 17 (orange).



RUO

ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe
Probe hybridized to human breast cancer cell line with amplification of the PIK3CA gene as indicated by multiple green signals in each nucleus.




RUO


ZytoLight SPEC MDM2/CEN 12 Dual Color Probe
Probe hybridized to liposarcoma tissue section with amplification of the MDM2 gene (green), CEN 12 (orange).

ZytoLight®


Reliable multi-target detection using Fluorescence in situ Hybridization!



FISH of HER-2 neu on an
FFPE Breast Carcinoma



FISH of EGFR on an
FFPE Colon Carcinoma



FISH of MDM2 on an
FFPE Liposarcoma

ZytoLight Single Color Probes

| Catalog # | Product Description | VOLUME |
|------------|------------------------------|--------|
| Z-2001-200 | ZytoLight CEN 3 Probe | 200 µl |
| Z-2002-200 | ZytoLight CEN 6 Probe | 200 µl |
| Z-2003-200 | ZytoLight CEN 7 Probe | 200 µl |
| Z-2004-200 | ZytoLight CEN 8 Probe | 200 µl |
| Z-2004-50 | ZytoLight CEN 8 Probe | 50 µl |
| Z-2005-200 | ZytoLight CEN 11 Probe | 200 µl |
| Z-2006-200 | ZytoLight CEN 17 Probe | 200 µl |
| Z-2007-200 | ZytoLight CEN 18 Probe | 200 µl |
| Z-2008-200 | ZytoLight CEN X Probe | 200 µl |
| Z-2010-200 | ZytoLight CEN Yq12 Probe | 200 µl |
| Z-2049-200 | ZytoLight SPEC 2q11 Probe | 200 µl |
| Z-2050-200 | ZytoLight CEN 12 Probe | 200 µl |
| Z-2067-200 | ZytoLight CEN 9 Probe | 200 µl |
| Z-2079-200 | ZytoLight CEN 10 Probe | 200 µl |
| Z-2083-200 | ZytoLight SPEC 4p11 Probe | 200 µl |
| Z-2085-200 | ZytoLight SPEC 13q12 Probe | 200 µl |
| Z-2086-200 | ZytoLight SPEC 21q22 Probe | 200 µl |
| Z-2101-200 | ZytoLight SPEC 1p12 Probe | 200 µl |
| Z-2123-200 | ZytoLight CEN Y (DYZ3) Probe | 200 µl |

ZytoLight Dual Color Probes

| Catalog # | Product Description | VOLUME |
|------------|---|--------|
| Z-2013-200 | ZytoLight SPEC MDM2/CEN 12 Dual Color Probe | 200 µl |
| Z-2013-50 | ZytoLight SPEC MDM2/CEN 12 Dual Color Probe | 50 µl |
| Z-2014-200 | ZytoLight SPEC MAML2 Dual Color Break Apart Probe | 200 µl |
| Z-2014-50 | ZytoLight SPEC MAML2 Dual Color Break Apart Probe | 50 µl |
| Z-2015-200 | ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe | 200 µl |
| Z-2015-50 | ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe | 50 µl |
| Z-2016-200 | ZytoLight CEN X/Yq12 Dual Color Probe | 200 µl |
| Z-2016-50 | ZytoLight CEN X/Yq12 Dual Color Probe | 50 µl |
| Z-2018-200 | ZytoLight SPEC FOXO1/PAX3 Dual Color Probe | 200 µl |
| Z-2018-50 | ZytoLight SPEC FOXO1/PAX3 Dual Color Probe | 50 µl |
| Z-2019-200 | ZytoLight SPEC FOXO1/PAX7 Dual Color Probe | 200 µl |
| Z-2019-50 | ZytoLight SPEC FOXO1/PAX7 Dual Color Probe | 50 µl |
| Z-2033-200 | ZytoLight SPEC EGFR/CEN 7 Dual Color Probe | 200 µl |
| Z-2033-50 | ZytoLight SPEC EGFR/CEN 7 Dual Color Probe | 50 µl |
| Z-2056-200 | ZytoLight SPEC ERBB3/CEN 12 Dual Color Probe | 200 µl |
| Z-2057-200 | ZytoLight SPEC ERBB4/2q11 Dual Color Probe | 200 µl |
| Z-2062-200 | ZytoLight SPEC FHIT/CEN 3 Dual Color Probe | 200 µl |
| Z-2063-200 | ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe | 200 µl |
| Z-2063-50 | ZytoLight SPEC CDKN2A/CEN 9 Dual Color Probe | 50 µl |
| Z-2069-200 | ZytoLight SPEC ESR1/CEN 6 Dual Color Probe | 200 µl |
| Z-2069-50 | ZytoLight SPEC ESR1/CEN 6 Dual Color Probe | 50 µl |
| Z-2071-200 | ZytoLight SPEC CCND1/CEN 11 Dual Color Probe | 200 µl |
| Z-2071-50 | ZytoLight SPEC CCND1/CEN 11 Dual Color Probe | 50 µl |
| Z-2072-200 | ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe | 200 µl |
| Z-2072-50 | ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe | 50 µl |
| Z-2074-200 | ZytoLight SPEC MYCN/2q11 Dual Color Probe | 200 µl |
| Z-2074-50 | ZytoLight SPEC MYCN/2q11 Dual Color Probe | 50 µl |
| Z-2075-200 | ZytoLight SPEC 1p36/1q25 Dual Color Probe | 200 µl |
| Z-2075-50 | ZytoLight SPEC 1p36/1q25 Dual Color Probe | 50 µl |
| Z-2076-200 | ZytoLight SPEC 19q13/19p13 Dual Color Probe | 200 µl |

ZytoLight Dual Color Probes

| Catalog # | Product Description | VOLUME |
|------------|---|--------|
| Z-2076-50 | ZytoLight SPEC 19q13/19p13 Dual Color Probe | 50 µl |
| Z-2077-200 | ZytoLight SPEC CEN 17/SPEC ERBB2 Dual Color Probe | 200 µl |
| Z-2077-50 | ZytoLight SPEC CEN 17/SPEC ERBB2 Dual Color Probe | 50 µl |
| Z-2078-200 | ZytoLight SPEC PTEN/CEN 10 Dual Color Probe | 200 µl |
| Z-2078-50 | ZytoLight SPEC PTEN/CEN 10 Dual Color Probe | 50 µl |
| Z-2080-200 | ZytoLight SPEC MDM4/1p12 Dual Color Probe | 200 µl |
| Z-2082-200 | ZytoLight SPEC FGFR3/CEN 4 Dual Color Probe | 200 µl |
| Z-2084-200 | ZytoLight SPEC VHL/CEN 3 Dual Color Probe | 200 µl |
| Z-2087-200 | ZytoLight SPEC MET/CEN 7 Dual Color Probe | 200 µl |
| Z-2087-50 | ZytoLight SPEC MET/CEN 7 Dual Color Probe | 50 µl |
| Z-2090-200 | ZytoLight SPEC MYC Dual Color Break Apart Probe | 200 µl |
| Z-2090-50 | ZytoLight SPEC MYC Dual Color Break Apart Probe | 50 µl |
| Z-2091-200 | ZytoLight SPEC TERT/5q31 Dual Color Probe | 200 µl |
| Z-2091-50 | ZytoLight SPEC TERT/5q31 Dual Color Probe | 50 µl |
| Z-2092-200 | ZytoLight SPEC MYC/CEN 8 Dual Color Probe | 200 µl |
| Z-2092-50 | ZytoLight SPEC MYC/CEN 8 Dual Color Probe | 50 µl |
| Z-2096-200 | ZytoLight SPEC EWSR1 Dual Color Break Apart Probe | 200 µl |
| Z-2096-50 | ZytoLight SPEC EWSR1 Dual Color Break Apart Probe | 50 µl |
| Z-2097-200 | ZytoLight SPEC S18 Dual Color Break Apart Probe | 200 µl |
| Z-2097-50 | ZytoLight SPEC S18 Dual Color Break Apart Probe | 50 µl |
| Z-2100-200 | ZytoLight SPEC DDT3 Dual Color Break Apart Probe | 200 µl |
| Z-2100-50 | ZytoLight SPEC DDT3 Dual Color Break Apart Probe | 50 µl |
| Z-2103-200 | ZytoLight SPEC CDK4/CEN 12 Dual Color Probe | 200 µl |
| Z-2103-50 | ZytoLight SPEC CDK4/CEN 12 Dual Color Probe | 50 µl |
| Z-2105-200 | ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe | 200 µl |
| Z-2105-50 | ZytoLight SPEC MYC/IGH Dual Color Dual Fusion Probe | 50 µl |
| Z-2107-200 | ZytoLight SPEC EGR1/5p15 Dual Color Probe | 200 µl |
| Z-2107-50 | ZytoLight SPEC EGR1/5p15 Dual Color Probe | 50 µl |
| Z-2108-200 | ZytoLight SPEC CCND1 Dual Color Break Apart Probe | 200 µl |
| Z-2108-50 | ZytoLight SPEC CCND1 Dual Color Break Apart Probe | 50 µl |
| Z-2109-200 | ZytoLight SPEC TFE3 Dual Color Break Apart Probe | 200 µl |
| Z-2109-50 | ZytoLight SPEC TFE3 Dual Color Break Apart Probe | 50 µl |
| Z-2110-200 | ZytoLight SPEC IGH Dual Color Break Apart Probe | 200 µl |
| Z-2110-50 | ZytoLight SPEC IGH Dual Color Break Apart Probe | 50 µl |
| Z-2111-200 | ZytoLight SPEC BCR/ABL1 Dual Color Dual Fusion Probe | 200 µl |
| Z-2111-50 | ZytoLight SPEC BCR/ABL1 Dual Color Dual fusion Probe | 50 µl |
| Z-2112-200 | ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe | 200 µl |
| Z-2112-50 | ZytoLight SPEC RUNX1/RUNX1T1 Dual Color Dual Fusion Probe | 50 µl |
| Z-2113-200 | ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe | 200 µl |
| Z-2113-50 | ZytoLight SPEC PML/RARA Dual Color Dual Fusion Probe | 50 µl |
| Z-2114-200 | ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe | 200 µl |
| Z-2114-50 | ZytoLight SPEC BCL2/IGH Dual Color Dual Fusion Probe | 50 µl |
| Z-2115-200 | ZytoLight SPEC KRAS/CEN 12 Dual Color Probe | 200 µl |
| Z-2116-200 | ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe | 200 µl |
| Z-2116-50 | ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe | 50 µl |
| Z-2119-200 | ZytoLight SPEC PDGFB Dual Color Break Apart Probe | 200 µl |
| Z-2119-50 | ZytoLight SPEC PDGFB Dual Color Break Apart Probe | 50 µl |
| Z-2120-200 | ZytoLight CEN Y/X Dual Color Probe | 200 µl |
| Z-2121-200 | ZytoLight SPEC COL1A1 Dual Color Break Apart Probe | 200 µl |
| Z-2122-200 | ZytoLight SPEC FGFR2/CEN 10 Dual Color Probe | 200 µl |
| Z-2124-200 | ZytoLight SPEC ALK Dual Color Break Apart Probe | 200 µl |
| Z-2124-50 | ZytoLight SPEC ALK Dual Color Break Apart Probe | 50 µl |
| Z-2125-200 | ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe | 200 µl |
| Z-2125-50 | ZytoLight SPEC CCND1/IGH Dual Color Dual Fusion Probe | 50 µl |

ZytoLight®

Reliable multi-target detection using Fluorescence in situ Hybridization!

ZytoLight Dual Color Probes

| Catalog # | Product Description | VOLUME |
|------------|--|--------|
| Z-2127-200 | ZytoLight SPEC SOX2/CEN 3 Dual Color Probe | 200 µl |
| Z-2130-50 | ZytoLight SPEC FUS Dual Color Break Apart Probe | 50 µl |
| Z-2131-50 | ZytoLight SPEC KIF5B Dual Color Break Apart Probe | 50 µl |
| Z-2132-50 | ZytoLight SPEC JAZF1 Dual Color Break Apart Probe | 50 µl |
| Z-2136-50 | ZytoLight SPEC EML4 Dual Color Break Apart Probe | 50 µl |
| Z-2137-50 | ZytoLight SPEC CARS Dual Color Break Apart Probe | 50 µl |
| Z-2138-200 | ZytoLight SPEC ERG Dual Color Break Apart Probe | 200 µl |
| Z-2139-50 | ZytoLight SPEC FOXO1 Dual Color Break Apart Probe | 50 µl |
| Z-2140-200 | ZytoLight SPEC PIK3CA/CEN 3 Dual Color Probe | 200 µl |
| Z-2142-50 | ZytoLight SPEC WT1 Dual Color Break Apart Probe | 50 µl |
| Z-2143-200 | ZytoLight SPEC MYB Dual Color Break Apart Probe | 200 µl |
| Z-2143-50 | ZytoLight SPEC MYB Dual Color Break Apart Probe | 50 µl |
| Z-2144-200 | ZytoLight SPEC ROS1 Dual Color Break Apart Probe | 200 µl |
| Z-2144-50 | ZytoLight SPEC ROS1 Dual Color Break Apart Probe | 50 µl |
| Z-2145-50 | ZytoLight SPEC NR4A3 Dual color Break Apart Probe | 50 µl |
| Z-2146-200 | ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe | 200 µl |
| Z-2146-50 | ZytoLight SPEC BIRC3/MALT1 Dual Color Dual Fusion Probe | 50 µl |
| Z-2148-200 | ZytoLight SPEC RET Dual Color Break Part Probe | 200 µl |
| Z-2148-50 | ZytoLight SPEC RET Dual Color Break Part Probe | 50 µl |
| Z-2151-50 | ZytoLight SPEC USP6 Dual Color Break Part Probe | 50 µl |
| Z-2153-200 | ZytoLight SPEC TP53/CEN 17 Dual Color Probe | 200 µl |
| Z-2153-50 | ZytoLight SPEC TP53/CEN 17 Dual Color Probe | 50 µl |
| Z-2157-50 | ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe | 50 µl |
| Z-2157-200 | ZytoLight SPEC ETV6/RUNX1 Dual Color Dual Fusion Probe | 200 µl |
| Z-2159-200 | ZytoLight SPEC TP53/ATM Dual Color Probe | 200 µl |
| Z-2159-50 | ZytoLight SPEC TP53/ATM Dual Color Probe | 50 µl |
| Z-2161-200 | ZytoLight SPEC ALK/2q11 Dual Color Probe | 200 µl |
| Z-2162-200 | ZytoLight SPEC ROS1/CEN 6 Dual Color Probe | 200 µl |
| Z-2164-200 | ZytoLight SPEC 13/21 Dual Color Probe | 200 µl |
| Z-2165-200 | ZytoLight SPEC RB1/13q12 Dual Color Probe | 200 µl |
| Z-2165-50 | ZytoLight SPEC RB1/13q12 Dual Color Probe | 50 µl |
| Z-2167-200 | ZytoLight SPEC NTRK1 Dual Color Break Apart Probe | 200 µl |
| Z-2167-50 | ZytoLight SPEC NTRK1 Dual Color Break Apart Probe | 50 µl |
| Z-2168-200 | ZytoLight SPEC FGFR1 Dual Color Break Apart Probe | 200 µl |
| Z-2168-50 | ZytoLight SPEC FGFR1 Dual Color Break Apart Probe | 50 µl |
| Z-2169-200 | ZytoLight SPEC FGFR2 Dual Color Break Apart Probe | 200 µl |
| Z-2170-200 | ZytoLight SPEC FGFR3 Dual Color Break Apart Probe | 200 µl |
| Z-2170-50 | ZytoLight SPEC FGFR3 Dual Color Break Apart Probe | 50 µl |
| Z-2171-200 | ZytoLight SPEC BCL2L1/CEN 20 Dual Color Probe | 200 µl |
| Z-2173-200 | ZytoLight SPEC MCL1/1p12 Dual Color Probe | 200 µl |
| Z-2174-50 | ZytoLight SPEC BCL2/CEN 18 Dual Color Break Apart Probe | 50 µl |
| Z-2175-50 | ZytoLight SPEC YWHAE Dual Color Break Apart Probe | 50 µl |
| Z-2176-200 | ZytoLight SPEC ETV6 Dual Color Break Apart Probe | 200 µl |
| Z-2176-50 | ZytoLight SPEC ETV6 Dual Color Break Apart Probe | 50 µl |
| Z-2177-200 | ZytoLight SPEC BCL6 Dual Color Break Apart Probe | 200 µl |
| Z-2177-50 | ZytoLight SPEC BCL6 Dual Color Break Apart Probe | 50 µl |
| Z-2178-50 | ZytoLight SPEC SMARCB1/22q12 Dual Color Probe | 50 µl |
| Z-2179-200 | ZytoLight SPEC CD274/PDCD1LG2/CEN 9 Dual Color Probe | 200 µl |
| Z-2179-50 | ZytoLight SPEC CD274/PDCD1LG2/CEN 9 Dual Color Probe | 50 µl |
| Z-2181-200 | ZytoLight SPEC NRG1 Dual Color Break Apart Probe | 200 µl |
| Z-2189-200 | ZytoLight SPEC BRAF Dual Color Break Apart Probe | 200 µl |
| Z-2190-50 | ZytoLight SPEC ERBB2/D17S122 Dual Color Probe | 50 µl |
| Z-2190-200 | ZytoLight SPEC ERBB2/D17S122 Dual Color Probe | 200 µl |
| Z-2191-200 | ZytoLight SPEC BRAF/CEN 7 Dual Color Probe | 200 µl |
| Z-2192-200 | ZytoLight SPEC BCL2 Dual Color Break Apart Probe | 200 µl |
| Z-2192-50 | ZytoLight SPEC BCL2 Dual Color Break Apart Probe | 50 µl |
| Z-2193-200 | ZytoLight SPEC KMT2A Dual Color Break Apart Probe | 200 µl |
| Z-2193-50 | ZytoLight SPEC KMT2A Dual Color Break Apart Probe | 50 µl |
| Z-2195-200 | ZytoLight SPEC VEGFA/CEN 6 Dual Color Probe | 200 µl |
| Z-2196-200 | ZytoLight SPEC MALT1 Dual Color Break Apart Probe | 200 µl |
| Z-2196-50 | ZytoLight SPEC MALT1 Dual Color Break Apart Probe | 50 µl |
| Z-2197-50 | ZytoLight SPEC PDGFRB Dual Color Break Apart Probe | 50 µl |
| Z-2198-50 | ZytoLight SPEC TP53/17q22 Dual Color Break Apart Probe | 50 µl |
| Z-2199-50 | ZytoLight SPEC ABL1 Dual Color Break Apart Probe | 50 µl |
| Z-2200-50 | ZytoLight SPEC ABL2 Dual Color Break Apart Probe | 50 µl |
| Z-2201-50 | ZytoLight SPEC CRLF2 Dual Color Break Apart Probe | 50 µl |
| Z-2202-50 | ZytoLight SPEC CSF1R Dual Color Break Apart Probe | 50 µl |
| Z-2205-200 | ZytoLight SPEC NTRK2 Dual Color Break Apart Probe | 200 µl |
| Z-2205-50 | ZytoLight SPEC NTRK2 Dual Color Break Apart Probe | 50 µl |
| Z-2206-200 | ZytoLight SPEC NTRK3 Dual Color Break Apart Probe | 200 µl |
| Z-2206-50 | ZytoLight SPEC NTRK3 Dual Color Break Apart Probe | 50 µl |
| Z-2207-50 | ZytoLight SPEC CBFB Dual Color Break Apart Probe | 50 µl |
| Z-2208-200 | ZytoLight SPEC NUTM1 Dual Color Break Apart Probe | 200 µl |
| Z-2210-50 | ZytoLight SPEC IRF4, DUSP22 Dual Color Break Apart Probe | 50 µl |
| Z-2211-50 | ZytoLight SPEC EGR1/D5S23,D5S721 Dual Color Probe | 50 µl |

ZytoLight Dual Color Probes

| Catalog # | Product Description | VOLUME |
|------------|---|--------|
| Z-2212-50 | ZytoLight SPEC WWTR1 Dual Color Break Apart Probe | 50 µl |
| Z-2213-50 | ZytoLight SPEC PTPRT/20q11 Dual Color Probe | 50 µl |
| Z-2215-50 | ZytoLight SPEC PHF1 Dual Color Break Apart Probe | 50 µl |
| Z-2265-50 | ZytoLight SPEC NUP214 Dual Color Break Apart Probe | 50 µl |
| Z-2266-50 | ZytoLight SPEC NUP98 Dual Color Break Apart Probe | 50 µl |
| Z-2267-50 | ZytoLight SPEC CREBBP Dual Color Break Apart Probe | 50 µl |
| Z-2268-50 | ZytoLight SPEC CSF1R/D5S23,D5S721 Dual Color Probe | 50 µl |
| Z-2270-50 | ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe | 50 µl |
| Z-2271-50 | ZytoLight SPEC MAFB/IGH Dual Color Dual Fusion Probe | 50 µl |
| Z-2273-50 | ZytoLight SPEC TERT Dual Color Break Apart Probe | 50 µl |
| Z-2274-50 | ZytoLight SPEC C19MC/19p13 Dual Color Probe | 50 µl |
| Z-2275-50 | ZytoLight SPEC ZNF384 Dual Color Break Apart Probe | 50 µl |
| Z-2276-50 | ZytoLight SPEC CKS1B/CDKN2C Dual Color Probe | 50 µl |
| Z-2278-200 | ZytoLight SPEC RICTOR/5q31.1 Dual Color Probe | 200 µl |
| Z-2280-50 | ZytoLight SPEC D13S319/13q34 Dual Color Probe | 50 µl |
| Z-2281-50 | ZytoLight SPEC MYB/CEN 6 Dual Color Probe | 50 µl |
| Z-2282-50 | ZytoLight SPEC FGF3/IGH Dual Color Dual Fusion Probe | 50 µl |
| Z-2284-200 | ZytoLight SPEC TERC/CEN 3 Dual Color Probe | 200 µl |
| Z-2285-50 | ZytoLight SPEC CIC Dual Color Break Apart Probe | 50 µl |
| Z-2286-50 | ZytoLight SPEC IGL Dual Color Break Apart Probe | 50 µl |
| Z-2288-50 | ZytoLight SPEC IGK Dual Color Break Apart Probe | 50 µl |
| Z-2287-50 | ZytoLight SPEC GATA2/MECOM Dual Color Dual Fusion Probe | 50 µl |
| Z-2291-50 | ZytoLight SPEC SPI1 Dual Color Break Apart Probe | 50 µl |
| Z-2294-50 | ZytoLight SPEC JAK2 Dual Color Break Apart Probe | 50 µl |
| Z-2296-50 | ZytoLight SPEC ATM/CEN 12 Dual Color Probe | 50 µl |
| Z-2297-50 | ZytoLight SPEC ATM/CEN 11 Dual Color Probe | 50 µl |
| Z-2299-50 | ZytoLight SPEC DiGeorge/Phelan McDerimid Dual Color Probe | 50 µl |
| Z-2302-50 | ZytoLight Williams-Beuren Dual Color Probe | 50 µl |

ZytoLight Triple Color Probes

| Catalog # | Product Description | VOLUME |
|------------|--|--------|
| Z-2093-200 | ZytoLight SPEC HER2/TOP2A/CEN 17 Triple Color Probe | 200 µl |
| Z-2093-50 | ZytoLight SPEC HER2/TOP2A/CEN 17 Triple Color Probe | 50 µl |
| Z-2095-200 | ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe | 200 µl |
| Z-2095-50 | ZytoLight SPEC 13/CEN 18/SPEC 21 Triple Color Probe | 50 µl |
| Z-2117-200 | ZytoLight SPEC ALK/EML4 TriCheck Probe | 200 µl |
| Z-2117-50 | ZytoLight SPEC ALK/EML4 TriCheck Probe | 50 µl |
| Z-2135-200 | ZytoLight SPEC ERG/TMPRSS2 TriCheck Probe | 200 µl |
| Z-2152-200 | ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe | 200 µl |
| Z-2152-50 | ZytoLight SPEC RREB1/MYB/CEN 6 Triple Color Probe | 50 µl |
| Z-2160-200 | ZytoLight SPEC CLL II Probe SPEC D13S319/13q34/CEN 12 Triple Color Probe | 200 µl |
| Z-2160-50 | ZytoLight SPEC D13S319/13q34/CEN 12 Triple Color Probe | 50 µl |
| Z-2163-200 | ZytoLight SPEC 18/CEN X/Y Triple Color Probe | 200 µl |
| Z-2180-200 | ZytoLight SPEC 12/CEN X/Yq12 Triple Color Probe | 200 µl |
| Z-2183-50 | ZytoLight SPEC EWSR1/FL11 TriCheck Probe | 50 µl |
| Z-2184-50 | ZytoLight SPEC S18/SSX1 TriCheck Probe | 50 µl |
| Z-2185-50 | ZytoLight SPEC FOXO1/PAX3 TriCheck Probe | 50 µl |
| Z-2194-200 | ZytoLight SPEC NRG/CD74 TriCheck Color Probe | 200 µl |
| Z-2209-50 | ZytoLight SPEC PDGFRA/FIP1L1 TriCheck Color Probe | 50 µl |
| Z-2214-50 | ZytoLight SPEC CUX1/EZH2/CEN 7 Triple Color Probe | 50 µl |
| Z-2216-50 | ZytoLight SPEC 11q gain/loss Triple Color Probe | 50 µl |
| Z-2277-50 | ZytoLight SPEC MEF2D/BCL9 TriCheck Probe | 50 µl |
| Z-2289-50 | ZytoLight SPEC DiGeorge Triple Color Probe | 50 µl |
| Z-2307-50 | ZytoLight SPEC 4p11/CEN 10/17 Triple Color Probe | 50 µl |

ZytoLight Quad Color Probes


| Catalog # | Product Description | VOLUME |
|------------|---|--------|
| Z-2081-50 | ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe | 50 µl |
| Z-2081-200 | ZytoLight SPEC CDKN2A/CEN 3/7/17 Quadruple Color Probe | 200 µl |
| Z-2102-200 | ZytoLight SPEC VHL/1p12/CEN 7/17 Quadruple Color Probe | 200 µl |
| Z-2118-200 | ZytoLight SPEC CCND1 Break Apart/2q11/CEN 6 Quadruple Color Probe | 200 µl |
| Z-2305-200 | ZytoLight SPEC Bladder Cancer Quadruple Color Probe | 50 µl |
| Z-2305-50 | ZytoLight SPEC Bladder Cancer Quadruple Color Probe | 200 µl |

ZytoLight FISH Kits

| Catalog # | Product Description | VOLUME |
|-----------|--|--------|
| Z-2020-20 | ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit | 20 |
| Z-2020-5 | ZytoLight SPEC ERBB2/CEN 17 Dual Color Probe Kit | 5 |
| Z-2028-20 | ZytoLight FISH-Tissue Implementation Kit | 20 |
| Z-2028-5 | ZytoLight FISH-Tissue Implementation Kit | 5 |
| Z-2099-20 | ZytoLight FISH-Cytology Implementation Kit | 20 |

FlexISH®


FlexISH brings flexibility to your FISH



FlexISH of ALK/ROS1 DistinguISH Probe with ALK Positivity on Lung Cancer



FlexISH of ALK/ROS1 DistinguISH Probe with ROS1 Positivity on Lung Cancer




FlexISH of ERBB2/CEN 17 on Breast Cancer

INTRODUCTION

With a hybridization temperature of 37°C the FlexISH protocol is fully compatible with routine workflows in pathology laboratories. In a recent comparison study, overnight and short hybridization periods showed excellent correlation between the FISH results obtained with FlexISH, ZytoLight and PathVysion ERBB2 probes. Short hybridization time does not negatively affect the performance, specimen quality or diagnostic results.

IMPORTANCE OF ERBB2 TESTING

Breast cancer is the most commonly diagnosed cancer among women worldwide. It is estimated that 1.7 million new cases are diagnosed per year, accounting for 25% of all new cancer cases in women. The proto-oncogene ERBB2 is amplified in approximately 20% of all breast cancers and is correlated with a poor prognosis. Breast cancer patients harboring ERBB2 amplification are addressed for a targeted therapy with trastuzumab (Herceptin), a humanized monoclonal antibody directed against the extracellular portion of the ERBB2 protein. This treatment is associated with significantly prolonged overall survival and time of tumor progression.



FlexISH Tissue Implementation Kit

| Catalog # | Product Description |
|-----------|-----------------------------------|
| Z-2182-5 | FlexISH-Tissue Implementation Kit |
| Z-2182-20 | FlexISH-Tissue Implementation Kit |

FlexISH Probes

| VOLUME | Catalog # | Product Description | VOLUME |
|--------|-----------|-----------------------------------|--------|
| 5 | Z-2293-50 | FlexISH MYC/IGH TriCheck Probe | 50 µl |
| 20 | Z-2295-50 | FlexISH IGK/IGL DistinguISH Probe | 50 µl |

FlexISH Probes

| Catalog # | Product Description | VOLUME |
|------------|---------------------------------------|--------|
| Z-2166-50 | FlexISH ERBB2/CEN 17 Dual Color Probe | 50 µl |
| Z-2166-200 | FlexISH ERBB2/CEN 17 Dual Color Probe | 200 µl |
| Z-2203-50 | FlexISH ALK/ROS1 DistinguISH Probe | 50 µl |
| Z-2203-200 | FlexISH ALK/ROS1 DistinguISH Probe | 200 µl |
| Z-2269-50 | FlexISH RET/KIF5B TriCheck Probe NEW | 50 µl |
| Z-2269-200 | FlexISH RET/KIF5B TriCheck Probe NEW | 200 µl |
| Z-2283-200 | FlexISH BCL2/BCL6 DistinguISH Probe | 200 µl |
| Z-2283-50 | FlexISH BCL2/BCL6 DistinguISH Probe | 50 µl |

FlexISH®

FlexISH brings flexibility to your FISH

ALK/ROS1 DISTINGUISH PROBE: ONE PROBE - TWO TARGETS

The FlexISH ALK/ROS1 DistingISH Probe is designed to simultaneously detect ALK and ROS1 rearrangements. Additionally, this innovative probe design enables the user to discriminate between possible aberrations affecting the chromosomal regions harboring either the ALK ROS1 gene. Rearrangements affecting the ALK or the ROS1 gene locus are frequently found in non-small cell lung cancer (NSCLC). The specific analysis of ROS1 and ALK rearrangements in NSCLC patients is a very effective and reliable tool for the diagnosis and selection of treatment with e.g. the tyrosine kinase inhibitor crizotinib.

RELIABLE RESULTS OBTAINABLE IN 1 DAY

With the use of the FlexISH ERBB@/CEN 17 Dual Color Probe, or ALK ROS1 DistingISH Probe in combination with the FlexISH Tissue Implementation Kit reliable results can be obtained already within 4.5 hours. The FlexISH protocol can also be incorporated into the routine workflow with overnight hybridization providing the highest flexibility.

FlexISH Workflow Schedule:



Dewaxing



Heat Pretreatment



Pepsin Digestion



Probe Application & Denaturation



120 min.

Adapt Hybridization Time to Your Needs

Overnight



Mouting




Strigency Wash




Analysis

ZytoDot®


Reliable and simple detection of genomic alterations using light microscopy!



CISH of HER-2 neu on an FFPE Breast Carcinoma



CISH of FGR1 on an FFPE Prostate Carcinoma



CISH of TOPO2A on an FFPE Breast Carcinoma

INTRODUCTION

The ZytoDot® products are designed for the detection of aneuploidies and gene amplifications by Chromogenic in situ Hybridization (CISH) in formalin-fixed, paraffin-embedded tissue sections, cell samples, blood or bone marrow smears, and metaphase chromosome spreads.

CISH: A RELIABLE ALTERNATIVE TO FISH

High concordance between CISH and FISH ranging from 92-100% has been shown by numerous international studies for HER2 amplification.

ADVANTAGES OF CISH OVER FISH

- Quick and easy interpretation of results comparable to IHC
- Simultaneous observation of tissue morphology and CISH signals
- Storage of slides at room temperature - CISH signals are permanent
- No costly fluorescent microscope needed

HIGH SIGNAL-TO-NOISE RATIO

The ZytoDot® probes are processed by the unique ZytoVision® Repeat Subtraction Technique resulting in advanced specificity and less background. No further blocking of repetitive sequences is needed!


ADVANTAGES OF ZYTODOT® 2C™

- Simultaneous observation of tissue morphology and CISH signals at 40x using light microscopy
- Two targets detected simultaneously
- High contrasting distinct red and green signals
- Quick and easy interpretation of results comparable to IHC
- Standardized and complete kits
- No costly fluorescent microscope needed

The well established ZytoDot® 2C™ system has been carefully revised to optimize its handling and performance. Handling steps were reduced by approx. 25% due to an optimized protocol including less antibody & washing steps and 2-component chromogenic substrates. This means a reduction by approx. 20% on day 1 and by approx. 50% on day 2! Procedure time is reduced to approx. 1¼ h per day due to shortened incubation times and less antibody & washing steps! Re-designed chromogenic substrates and improved antibody-cocktails lead to stronger signals while showing less background staining and occur even in sub-optimal pre-treated tissue sections.


The ZytoDot® system uses Digoxigenin-labeled probes

1. Which are, after blocking
2. Detected using a Mouse-anti-Digoxigenin antibody
3. This antibody is detected by a polymerized HRP-Goat-anti-Mouse antibody
4. The enzymatic reaction of DAB
5. Leads to the formation of strong permanent brown signals that can be visualized by light microscopy using a 40x objective.



The ZytoDot® 2C™ system uses DIG- and DNP-labeled probe cocktails targeting different genomic sections (see diagram)


1. Which are detected using a Mouse-anti-DIG/Rabbitanti-DNP cocktail
2. These antibodies are detected by a unique cocktail of polymerized HRP-Goat-anti-Mouse/AP-Goat-anti-Rabbit antibodies
3. The enzymatic reaction of AP-Red
4. and HRP-Green
5. leads to the formation of strong permanent red respectively green signals that can be visualized by light microscopy using a 40x objective.



ZytoDot® CISH - Solid Tumor Probes By Application


Like ZytoFast, the ZytoDot line is optimized for use on formalin-fixed paraffin-embedded (FFPE) tissues or cell samples. Results can be easily interpreted via light microscopy.

Brain and Neural Tumor Probes



| |
|------------------------------------|
| ZytoDot 2C Glioma 1p/19q Probe Set |
| ZytoDot 2C SPEC 1p36/1q25 Probe |
| ZytoDot 2C SPEC 19q13/19p13 Probe |
| ZytoDot 2C SPEC CDKN2A/CEN 9 Probe |
| ZytoDot SPEC EGFR Probe |
| ZytoDot 2C SPEC EGFR/CEN 7 Probe |
| ZytoDot 2C SPEC MET/CEN 7 Probe |
| ZytoDot SPEC MYCN Probe |
| ZytoDot 2C SPEC PTEN/CEN 10 Probe |

Breast Cancer Probes




| |
|-----------------------------------|
| ZytoDot EGFR Probe |
| ZytoDot 2C EGFR/CEN 7 Probe |
| ZytoDot ERBB2 Probe |
| ZytoDot ERBB2 Probe Kit |
| ZytoDot 2C ERBB2/CEN 17 Probe |
| ZytoDot 2C ERBB2/CEN 17 Probe Kit |
| ZytoDot 2C ERBB2/D17S122 Probe |
| ZytoDot ESR1 Probe |
| ZytoDot 2C FGFR1/CEN 8 Probe |
| ZytoDot 2C FGFR2/CEN 10 Probe |
| ZytoDot MYC Probe |
| ZytoDot 2C TOP2A/CEN 17 Probe |

Cervical Cancer Probes



| |
|------------------------|
| ZytoDot SPEC MYC Probe |
|------------------------|

Gastrointestinal Cancer




| |
|--|
| ZytoDot SPEC ERBB2 Probe |
| ZytoDot SPEC ERBB2 Probe |
| ZytoDot 2C SPEC ERBB2/CEN 17 Probe |
| ZytoDot 2C SPEC ERBB2/CEN 17 Probe Kit |
| ZytoDot 2C SPEC ERBB2/D17S122 Probe |
| ZytoDot 2C MDM2 Probe |
| ZytoDot 2C SPEC MDM2/CEN 12 Probe |

Salivary Gland Carcinoma




| |
|-------------------------------------|
| ZytoDot SPEC EWSR Break Apart Probe |
|-------------------------------------|

Sarcoma Probes




| |
|--|
| ZytoDot 2C SPEC ALK Break Apart Probe |
| ZytoDot 2C SPEC CDK4/CEN 12 Probe |
| ZytoDot 2C DDIT3 Break Apart Probe |
| ZytoDot 2C EWSR1 Break Apart Probe |
| ZytoDot 2C FOXO1 Break Apart Probe |
| ZytoDot 2C FUS Break Apart Probe |
| ZytoDot MDM2 Probe |
| ZytoDot 2C MDM2/CEN 12 Probe |
| ZytoDot SPEC MYC Probe |
| ZytoDot 2C SS18 Break Apart Probe |
| ZytoDot 2C SPEC USP6 Break Apart Probe |

Prostate Cancer Probes




| |
|---------------------------------------|
| ZytoDot 2C SPEC ERG Break Apart Probe |
| ZytoDot 2C SPEC PTEN/CEN 10 Probe |

Lung Cancer Probes




| |
|--|
| ZytoDot 2C SPEC ALK Break Apart Probe |
| ZytoDot 2C SPEC EML4 Break Apart Probe |
| ZytoDot SPEC EGFR Probe |
| ZytoDot 2C SPEC EGFR/CEN 7 Probe |
| ZytoDot SPEC ERBB2 Probe |
| ZytoDot SPEC ERBB2 Probe Kit |
| ZytoDot 2C SPEC ERBB2/CEN 17 Probe |
| ZytoDot 2C SPEC ERBB2/CEN 17 Probe Kit |
| ZytoDot 2C SPEC ERBB2/D17S122 Probe |
| ZytoDot 2C SPEC FGFR1/CEN 8 Probe |
| ZytoDot 2C SPEC FGFR2/CEN 10 Probe |
| ZytoDot 2C SPEC MET/CEN 7 Probe |
| ZytoDot 2C SPEC RET Break Apart Probe |
| ZytoDot 2C SPEC ROS1 Break Apart Probe |

ZytoDot® CISH - Hematology Specific Probes By Application



Acute Lymphoblastic Leukemia (ALL) Probes

ZytoDot 2C CDKN2A/CEN 9 Probe



Acute Myelogenous Leukemia (AML) Probes

ZytoDot CEN 8 Probe




Chronic Lymphocytic Leukemia (CLL) Probes

ZytoDot 2C SPEC BCL2 Break Apart Probe

ZytoDot 2C SPEC CCND1 Probe

ZytoDot SPEC MYC Probe




Chronic Myelogenous Leukemia (CML) Probes

ZytoDot CEN 8 Probe

Multiple Myeloma Probes


ZytoDot SPEC CCND1 Probe

ZytoDot 2C SPEC IGH Break Apart Probe



Myelodysplastic Syndrome (MDS) Probes

ZytoDot CEN 8 Probe



Non-Hodgkin Lymphoma Probes

ZytoDot 2C SPEC BCL2 Break Apart Probe


ZytoDot 2C SPEC BCL6 Break Apart Probe

ZytoDot 2C SPEC CCND1 Break Apart Probe


ZytoDot 2C SPEC IGH Break Apart Probe

ZytoDot 2C SPEC MALT1 Break Apart Probe

ZytoDot 2C MYC Break Apart Probe



ZytoDot® CISH - Genetic Probes By Application



Sex Mismatched Bone-Marrow Transplantant Management Probes

ZytoDot CEN X Probe

ZytoDot CEN Yq12 Probe

ZytoDot 2C CEN X/Y Probe


Prenatal, Postnatal, and Preimplantation Genetics Probes

ZytoDot SPEC 21q22 Probe

ZytoDot CEN X Probe

ZytoDot CEN Yq12 Probe

ZytoDot 2C CEN X/Y Probe



ZytoDot®

Reliable and simple detection of genomic alterations using light microscopy!

ZytoDot® products for Chromogenic in situ Hybridization (CISH) are designed for the detection of a euploidies and gene amplifications associated with tumors and genetic diseases using an IHC-like procedure and light microscopy.

ZytoDot Single Color Probes

| Catalog # | Product Description | VOLUME |
|------------|--------------------------|--------|
| C-3001-400 | ZytoDot SPEC ERBB2 Probe | 400 µl |
| C-3002-400 | ZytoDot CEN 6 Probe | 400 µl |
| C-3006-400 | ZytoDot CEN 17 Probe | 400 µl |
| C-3007-400 | ZytoDot SPEC EGFR Probe | 400 µl |
| C-3008-400 | ZytoDot CEN 7 Probe | 400 µl |
| C-3012-400 | ZytoDot SPEC MDM2 Probe | 400 µl |
| C-3013-400 | ZytoDot SPEC MYC Probe | 400 µl |
| C-3014-400 | ZytoDot CEN 12 Probe | 400 µl |
| C-3016-400 | ZytoDot CEN 8 Probe | 400 µl |
| C-3020-400 | ZytoDot CEN Yq12 Probe | 400 µl |
| C-3025-400 | ZytoDot CEN X Probe | 400 µl |
| C-3026-400 | ZytoDot SPEC 21q22 Probe | 400 µl |
| C-3029-400 | ZytoDot SPEC MYCN Probe | 400 µl |
| C-3035-400 | ZytoDot SPEC 1p12 Probe | 400 µl |
| C-3045-400 | ZytoDot CEN 3 Probe | 400 µl |
| C-3051-400 | ZytoDot SPEC 2q11 Probe | 400 µl |
| C-3052-400 | ZytoDot SPEC 13q12 Probe | 400 µl |


ZytoDot Dual Color Probes

| Catalog # | Product Description | VOLUME |
|------------|---|--------|
| C-3032-100 | ZytoDot 2C SPEC ERBB2/CEN 17 Probe | 100 µl |
| C-3032-400 | ZytoDot 2C SPEC ERBB2/CEN 17 Probe | 400 µl |
| C-3033-100 | ZytoDot 2C SPEC EGFR/CEN 7 Probe | 100 µl |
| C-3033-400 | ZytoDot 2C SPEC EGFR/CEN 7 Probe | 400 µl |
| C-3036-100 | ZytoDot 2C SPEC 1p36/1q25 Probe | 100 µl |
| C-3036-400 | ZytoDot 2C SPEC 1p36/1q25 Probe | 400 µl |
| C-3037-100 | ZytoDot 2C SPEC 19q13/19p13 Probe | 100 µl |
| C-3037-400 | ZytoDot 2C SPEC 19q13/19p13 Probe | 400 µl |
| C-3040-400 | ZytoDot 2C SPEC TOPO2A/CEN 17 Probe | 400 µl |
| C-3043-100 | ZytoDot 2C SPEC EWSR1 Break Apart Probe | 100 µl |
| C-3046-100 | ZytoDot 2C SPEC SS18 Break Apart Probe | 100 µl |
| C-3047-100 | ZytoDot 2C SPEC DDIT3 Break Apart Probe | 100 µl |
| C-3048-400 | ZytoDot 2C CEN X/Y Probe | 400 µl |
| C-3049-100 | ZytoDot 2C SPEC MDM2/CEN12 Probe | 100 µl |
| C-3049-400 | ZytoDot 2C SPEC MDM2/CEN12 Probe | 400 µl |


| | | |
|------------|---|--------|
| C-3050-400 | ZytoDot 2C SPEC FGFR1/CEN8 Probe | 400 µl |
| C-3053-400 | ZytoDot 2C SPEC PTEN/CEN 10 Probe | 400 µl |
| C-3054-100 | ZytoDot 2C SPEC FUS Break Apart Probe | 100 µl |
| C-3055-100 | ZytoDot 2C SPEC ALK Break Apart Probe | 100 µl |
| C-3055-400 | ZytoDot 2C SPEC ALK Break Apart Probe | 400 µl |
| C-3056-400 | ZytoDot 2C SPEC FGFR 2/CEN 10 Probe | 400 µl |
| C-3057-400 | ZytoDot 2C SPEC MET/CEN 7 Probe | 400 µl |
| C-3058-400 | ZytoDot 2C SPEC ERG Break Apart Probe | 400 µl |
| C-3059-400 | ZytoDot 2C SPEC EML4 Break Apart Probe | 400 µl |
| C-3062-400 | ZytoDot 2C SPEC CDK4/CEN 12 Probe | 400 µl |
| C-3063-100 | ZytoDot 2C SPEC ROS1 Break Apart Probe | 100 µl |
| C-3063-400 | ZytoDot 2C SPEC ROS1 Break Apart Probe | 400 µl |
| C-3064-100 | ZytoDot 2C SPEC RET Break Apart Probe | 100 µl |
| C-3064-400 | ZytoDot 2C SPEC RET Break Apart Probe | 400 µl |
| C-3065-100 | ZytoDot 2C SPEC FOXO1 Break Apart Probe | 100 µl |
| C-3066-400 | ZytoDot 2C SPEC MYC Break Apart Probe | 400 µl |
| C-3067-400 | ZytoDot 2C SPEC CDKN2A/CEN 9 Probe | 400 µl |
| C-3068-100 | ZytoDot 2C SPEC ERBB2/D17S122 Probe | 100 µl |
| C-3071-100 | ZytoDot 2C SPEC IGH Break Apart Probe | 100 µl |
| C-3072-100 | ZytoDot 2C SPEC MALT1 Break Apart Probe | 100 µl |
| C-3073-100 | ZytoDot 2C SPEC BCL2 Break Apart Probe | 100 µl |
| C-3074-100 | ZytoDot 2C SPEC BCL6 Break Apart Probe | 100 µl |
| C-3075-100 | ZytoDot 2C SPEC CCND1 Break Apart Probe | 100 µl |
| C-3076-10 | ZytoDot 2C Glioma 1p/19q Probe Set NEW | |
| C-3076-40 | ZytoDot 2C Glioma 1p/19q Probe Set NEW | |
| C-3077-100 | ZytoDot 2C SPEC USP6 Break Apart Probe NEW | 100µl |
| C-3078-100 | ZytoDot 2C SPEC NTRK1 Break Apart Probe NEW | 100µl |

ZytoDot CISH Kits

| Catalog # | Product Description | VOLUME |
|-----------|--|--------|
| C-3003-40 | ZytoDot SPEC ERBB2 Probe Kit | 40 |
| C-3005-40 | ZytoDot CISH Polymer Detection Kit | 40 |
| C-3018-40 | ZytoDot CISH Implementation Kit | 40 |
| C-3022-10 | ZytoDot 2C SPEC ERBB2/CEN 17 Probe Kit | 10 |
| C-3022-40 | ZytoDot 2C SPEC ERBB2/CEN 17 Probe Kit | 40 |
| C-3028-40 | ZytoDot 2C CISH Polymer Detection Kit | 40 |
| C-3044-10 | ZytoDot 2C CISH Implementation Kit | 10 |
| C-3044-40 | ZytoDot 2C CISH Implementation Kit | 40 |




ZytoDot 2C SPEC PTEN/CEN 10 Probe
Probe hybridized to prostate cancer tissue section with deletion of the PTEN gene as indicated by one green signal.




ZytoLight SPEC IGH Break Apart Probe
Probe hybridized to normal interphase cells as indicated by two red/green fusion signals per nucleus.


Achieving Chromogenic in situ Hybridization results in just 4 hours!



CISH of HPV on an
FFPE Cervical Cancer



CISH of EBV on an
FFPE Hodgkin's Lymphoma



CISH of Ig-Kappa/Lambda on an
FFPE Tonsil Tissue

INTRODUCTION

The ZytoFast® products are designed for outstandingly fast detection and discrimination of human pathogen viruses, e.g. HPV, EBV, CMV, and the determination of lymphocyte clonality by detecting Ig-k and Ig-l light chain RNA by Chromogenic in situ Hybridization (CISH) in formalin-fixed, paraffin-embedded tissue sections and cell samples.

ZYTOFAST®: OUTSTANDINGLY FAST CISH


Optimized protocols and faster tissue penetration due to short oligonucleotide probes of the ZytoFast® system make the ZytoFast® CISH procedure outstandingly fast. Single color results can be achieved within just 4 hours, with hands-on time being only about 2 hours!

HIGH SENSITIVITY AND SPECIFICITY

All ZytoFast® probes are tagged using the unique ZytoFast® HighTag System providing improved signal intensity! High specificity without risk of cross-hybridizations is obtained due to optimized oligonucleotide probes.

ADVANTAGES OF CISH

- Simultaneous observation of tissue morphology and CISH signals
- No risk of false positives due to mispriming or contamination as with PCR
- Easy method comparable to IHC
- No costly equipment needed
- Ability to test archival specimens
- High sensitivity and specificity



ZYTOFAST® KITS – USER FRIENDLY SOLUTIONS

CISH analysis has been made reliable and user friendly because for many targets complete kits are available, including all necessary pretreatment solutions, wash buffers, antibodies, chromogenic substrates, positive and negative control probes as well as a detailed protocol to perform successful CISH experiments. Additionally, all ZytoFast® probes are available separately. Thus, if an increased sensitivity is demanded, the Digoxigenin-labeled probes can be combined easily with any ZytoFast® PLUS CISH Implementation Kit.

Thus, if an increased sensitivity is demanded, the Digoxigenin-labeled probes can be combined easily with any ZytoFast® PLUS CISH Implementation Kit.

The ZytoFast® system uses oligonucleotide probes tagged with Biotin or Digoxigenin 1. Which are detected using enzyme-conjugated antibodies or streptavidin targeting the tags 2. The enzymatic reaction of chromogenic substrates 3. e.g. BCIP/NBT or AEC, leads to the formation of strong color precipitates that can be visualized by light microscopy.

ZytoFast®

Achieving Chromogenic in situ Hybridization results in just 4 hours!

ZytoFast allows users to implement CISH technologies into the laboratory with minimal effort and reduced hands-on time. The ZytoFast line of products allows for the rapid detection of human pathogens and lymphocyte clonality on FFPE tissues and cell samples. ZytoFast CISH results can be generated in as little as 4 hours, and viewed under a standard light microscope.

ZytoFast® CISH - Solid Tumor Probes By Application

| Cervical Cancer Probes | |
|------------------------|--|
| Digoxigenin-labeled | |
| | ZytoFast HPV type 6/11 Probe |
| | ZytoFast HPV type 16/18 Probe |
| | ZytoFast HPV type 31/33 Probe |
| | ZytoFast HPV 16/18/31/33/35/39/45/51/52/56/58/59/66/68/82 High-Risk (HR) Types Probe |
| | ZytoFast HPV type 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82 Screening Probe |

ZytoFast® CISH - Hematology Specific Probes By Application

| Lymphoma, Other Probes | |
|----------------------------|--|
| Digoxigenin-labeled | |
| | ZytoFast human Ig-kappa Probe |
| | ZytoFast human Ig-lambda Probe |
| Digoxigenin/Biotin-labeled | |
| | ZytoFast human Ig-kappa/Ig-lambda Probe |
| | ZytoFast human Ig-kappa/Ig-lambda CISH Kit |
| | ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit |

ZytoFast Complete Kits (Digoxigenin/Biotin labeled)

| Catalog # | Product Description | VOLUME |
|-----------|--|--------|
| T-1005-40 | ZytoFast human Ig-kappa/Ig-lambda CISH Kit | 40 |
| T-1105-40 | ZytoFast human Ig-kappa/Ig-lambda Permanent CISH Kit | 40 |

ZytoFast Biotin & Digoxigenin Labeled Probes (Digoxigenin/Biotin labeled)


| Catalog # | Product Description | VOLUME |
|------------|---|--------|
| T-1017-400 | ZytoFast human Ig-kappa/Ig-lambda Probe | 400 µl |

ZytoFast Digoxigenin Labeled Probes (Digoxigenin labeled)

| Catalog # | Product Description | VOLUME |
|------------|--|--------|
| T-1053-400 | ZytoFast DNA (+) Control Probe | 400 µl |
| T-1054-400 | ZytoFast DNA (-) Control Probe | 400 µl |
| T-1055-400 | ZytoFast HPV type 6/11 Probe | 400 µl |
| T-1056-400 | ZytoFast HPV type 16/18 Probe | 400 µl |
| T-1057-400 | ZytoFast HPV type 31/33 Probe | 400 µl |
| T-1113-400 | ZytoFast CMV Probe | 400 µl |
| T-1114-400 | ZytoFast EBV Probe (Digoxigenin labeled) | 400 µl |
| T-1115-400 | ZytoFast human Ig-kappa Probe | 400 µl |
| T-1116-400 | ZytoFast human Ig-lambda Probe | 400 µl |
| T-1119-400 | ZytoFast RNA (-) Control Probe | 400 µl |
| T-1120-400 | ZytoFast RNA (+) Control Probe (Digoxigenin labeled) | 400 µl |
| T-1140-400 | ZytoFast HPV High-Risk (HR) Types Probe | 400 µl |
| T-1144-400 | ZytoFast HPV Screening Probe | 400 µl |
| | 16/18/31/33/35/39/45/51/52/56/58/59/66/68/82 | 400 µl |
| | 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82 | 400 µl |


ZytoFast CISH DETECTIONS SYSTEMS FOR THE DETECTION OF DIGOXIGENIN PROBES

| Catalog # | Product Description | VOLUME |
|-----------|--|--------|
| T-1061-40 | ZytoFast PLUS CISH Implementation Kit AP-NBT/BCIP | 40 |
| T-1063-40 | ZytoFast PLUS CISH Implementation Kit HRP-DAB | 40 |
| T-1151-40 | ZytoFast PLUS CISH Implementation Kit AP-Permanent Red | 40 |



ZytoFast HPV type 6/11 Probe
Probe Hybridized to HPV infected cervix tissue, detected with the ZytoFast PLUS CISH Implementation Kit AP-Permanent Red.

Ancillaries and Equipment for FISH & CISH




RUO

Ancillaries for FISH & CISH


| Catalog # | Product Description | VOLUME |
|--------------|-------------------------------------|------------|
| AB-0001-30 | Mouse-anti-DIG | 30 ml |
| AB-0001-4 | Mouse-anti-DIG | 4 ml |
| AB-0002-4 | Anti-Mouse-HRP-Polymer | 4 ml |
| AB-0011-4 | Rabbit-anti-DIG | 4 ml |
| AB-0013-4 | HRP/AP-Polymer-Mix | 4 ml |
| AB-0014-4 | Anti-DIG/DNP-Mix | 4 ml |
| AB-0015-4 | Anti-Biotin/DIG-Mix | 4 ml |
| BS-0001-4 | Blocking Solution | 4 ml |
| BS-0002-8 | ZyBlack Quenching Solution | 4 ml |
| CS-0001-20 | Mayer's Hematoxylin Solution | 20 ml |
| CS-0002-20 | Nuclear Blue Solution | 20 ml |
| CS-0003-20 | Nuclear Red Solution | 20 ml |
| CS-0004-20 | Nuclear Green Solution | 20 ml |
| E-4005-125 | Fixogum Rubber Cement | 125 g |
| E-4005-50 | Fixogum Rubber Cement | 50 g |
| E-4007-2 | HER2 Control Slide Set | 1 |
| ES-0001-1000 | Pepsin Solution | 1000 ml |
| ES-0001-4 | Pepsin Solution | 4 ml |
| ES-0001-50 | Pepsin Solution | 50 ml |
| ES-0001-8 | Pepsin Solution Set | 2 x 4 ml |
| ES-0002-4 | Cytology Pepsin Solution | 4 ml |
| ES-0002-50 | Cytology Pepsin Solution | 50 ml |
| ES-0007-50 | VisionArray HPV PreCise Master Mix | 50 ml |
| ES-0008-50 | VisionArray MYCO PreCise Master Mix | 50 ml |
| MT-0004-4 | Mounting Solution (alcoholic) | 50 ml |
| MT-0007-0.8 | DAPI/DuraTect™-Solution | .08 ml |
| MT-0008-0.8 | DAPI/DuraTect™-Solution (ultra) | .08 ml |
| PT-0001-1000 | Heat Pretreatment Solution Citric | 2 x 500 ml |
| PT-0002-500 | Heat Pretreatment Solution EDTA | 500 ml |
| PT-0006-100 | Formaldehyde Dilution Buffer Set | 2 x 50ml |
| WB-0001-500 | Wash Buffer SSC | 500 ml |
| WB-0002-50 | 25x Wash Buffer A | 50 ml |
| WB-0003-50 | 20x SSC Solution | 50 ml |
| WB-0004-1000 | PBS/Tween | 1 tabl. |
| WB-0005-50 | 20x Wash Buffer TBS | 50 ml |
| WB-0007-500 | Cytology Stringency Wash Buffer SSC | 500 ml |
| WB-0008-500 | Cytology Wash Buffer SSC | 500 ml |
| WB-0010-500 | 5x FlexISH Wash Buffer | 500 ml |

Equipment

| Catalog # | Product Description | VOLUME |
|-----------|--|--------|
| E-4010-1 | DAPI/ZyGreen/ZyOrange Triple Bandpass Filter Set | 1 |
| E-4012-1 | ZyGreen Single Bandpass Filter Set v2 | 1 |
| E-4013-1 | ZyOrange Single Bandpass Filter Set v2 | 1 |
| E-4016-1 | ZyGreen/ZyOrange Dual Bandpass Filter Set v2 | 1 |
| E-4017-1 | ZyRed Single Bandpass Filter Set v2 | 1 |
| E-4026-1 | ZyBlue Single Bandpass Filter Set v2 | 1 |
| E-4027-1 | ZyGold Single Bandpass Filter Set v2 | 1 |
| E-4028-1 | ZyBlue/ZyGreen/ZyOrange Triple Bandpass Filter Set | 1 |
| E-4030-1 | DAPI Single Bandpass Filter Set v2 | 1 |
| E-4111-1 | ZEISS Fluorescence Filter Holder "FL EC P&C" | 1 |
| E-4113-1 | ZEISS Fluorescence Filter Holder "FL" | 1 |
| E-4121-1 | OLYMPUS Fluorescence Filter Holder "U-MF 2" | 1 |
| E-4122-1 | OLYMPUS Fluorescence Filter Holder "U-FF" | 1 |
| E-4131-1 | LEICA Fluorescence Filter Holder "DM K" | 1 |
| E-4141-1 | NIKON Fluorescence Filter Holder "C-FL" | 1 |



www.BioSB.com



BIO-SB SAMPLE MICROARRAY
TISSUE STAIN SLIDE 001

Bio SB
MICROARRAYS FOR THE WORLD



Equipment & Slides

Semi-Manual and Automated Equipment
for Immunohistochemistry (IHC), Immunocytochemistry (ICC),
Immunofluorescence (IF) and In Situ Hybridization (FISH & CISH)



TintoStainer Plus Automatic IHC Stainer



| Product Description | Catalog # |
|------------------------------------|-----------|
| TintoStainer Plus Automated System | BSB-7400 |
| TintoChamber (pack of 50) | BSB-7450 |
| Tinto Mixing Station Kit | BSB-7453 |
| TintoReagent Container 7mL | BSB-7454 |
| TintoReagent Container 30mL | BSB-7457 |
| TintoSlide Tray | BSB-7464 |
| TintoReagent Rack | BSB-7465 |



Technical Specifications

Slide Capacity:

30 slides

Temperature Control:

Room temperature ~ 100°C

Reagent Container Capacity:

7ml and 30 ml

Number of Reagent Containers:

36

Bulk Reagent Container Capacity:

2.5L

Hazardous Waste Container Capacity:

2.5L

External Bulk Waste Container Capacity:

18L

Operating Voltage:

100~240V 50/60Hz

Power Consumption:

1200VA

Operation Temperature:

5 °C~35 °C

Operating Humidity:

10~80% RH

Weight




200KG

Dimensions (WxDxH):

750x800x1300mm

BIO SB TINTOSTAINER PLUS AUTOMATIC IHC STAINER

Affordable | Efficient | Reliable



Automated System for Anatomical Pathology

TintoStainer Plus Automatic IHC Stainer

- Fully automated immunohistochemical system. Completes the whole process from deparaffinization to counterstain, improving work efficiency in the laboratory.
- For all tissue types: paraffin sections, frozen sections, and cell specimens.
- Multiple applications: Immunohistochemistry (IHC), Mohs IHC (TintoFast IHC), Immunocytochemistry (ICC), Immunofluorescence (IF) on FFPE, frozen tissues and cell preparations.

High Efficiency Procedure

- Patented sweeping cap gap reagent motion results in better sensitivity and reproducibility.
- Three independent slide trays with up to 10 slides for a total capacity of 30 slides.
- Walkaway fully automated staining from deparaffinization to counterstaining.
- Barcoding identification of slides and reagents.
- Each slide is thermally monitored for precise temperature control.
- 650+ IHC, Mohs and IF IVD antibodies, detection systems and ancillaries available.
- Full service professional support network.



Low Reagent Consumption


- Low reagent consumption makes tests more cost effective.

Sensitive Detection System

- Unique monomeric anti-mouse and anti rabbit Fab' and micropolymer conjugated to enzymes allow for easier tissue penetration which increase sensitivity and reproducibility.
- Universal IHC micropolymer detection systems, with accurate and reliable detection suitable for all antibodies.

High Performance

- Optimized staining procedures based on different types of antibody properties.
- Real-time operation status monitoring of instruments and slides.
- Teflon coated titanium allow needle to avoid cross contamination.



Left: Easy to operate user interface.



Patented Technology

- Sweeping Cap Gap Reagent Motion technology results in better sensitivity, reproducibility and homogenous signals.
- High quality and consistent results with minimal reagent use.

TintoDetector — Capillary Gap System for ICC, IHC, CISH & FISH



The Bio SB TintoDetector is a compact molecular pathology system which allows for full implementation of immunohistochemical (IHC), immunofluorescence (IF), chromogenic in-situ hybridization (CISH) and fluorescent in-situ hybridization (FISH) protocols in your diagnostic laboratory.

The Bio SB TintoDetector is intended for laboratories which require flexibility in their diagnostic needs, yet may not require high volume slide output offered by automated systems, such as the Bio SB TintoStainer.

The Bio SB TintoDetector is an open system that allows for the use of any detection system or reagent. The TintoDetector paired with accessories like the Bio SB TintoRetriever pressure cooker and Bio SB TintoDetector slide holder allows laboratories to implement a wide variety of lab protocols at a low cost.




- 8 staining dishes for reagents
- Built in incubator (RT to 110° C)
- 3 temperature presets
- Stainless steel durable construction
- Xylene resistant
- Chromogen resistant
- Small bench footprint
- Supports ICC protocols
- Supports IHC protocols
- Supports CISH protocols
- Supports FISH protocols
- Open System

TintoDetector Systems and Accessories


| Product Description | Catalog # |
|--|-----------|
| TintoDetector Complete System | BSB 7000 |
| TintoDetector Staining System - Without Incubator | BSB 7001 |
| TintoDetector Incubator | BSB 7002 |
| TintoDetector Slide Holder for Cap Gap Slides | BSB 7003 |
| TintoDetector 30-well Reagent Holder (Pack of 5) | BSB 7004 |
| TintoDetector Slide Holder for Hydrophilic Plus Slides | BSB 7035 |
| TintoDetector Absorbent Pads (Box of 10) | BSB 7036 |

TintoDetector — Capillary Gap System for ICC, IHC, CISH & FISH



Sample TintoDetector Lab Applications

FISH




SPEC ALK/EML4 Tricheck™ Triple color FISH Kit

Bio SB Catalog Number
Z-2117-200

SPEC ALK/EML4 Tricheck™ is an innovative Triple color probe used in the detection of ALK/EML4 rearrangements and translocation on human lung tissue.

CISH




HER2/CEN17 Dual color CISH Kit

Bio SB Catalog Number
C-3032-400

The Dual color HER2/CEN17 CISH kit is used in the detection and amplification of human HER2 gene and alpha satellites of Chromosome 17 on human breast tissue.

ICC




Ki67 Rabbit Monoclonal Clone: EP5

Bio SB Catalog Number
BSB 5713

The Bio SB Ki-67 Rabbit monoclonal antibody can be used in both immunocytochemical detection of Ki67 protein. In this example, Ki67 is used to determine Ki67 positivity in cell culture derived Caski cells.

IHC



bcl-6 Rabbit Monoclonal Clone: RBT-bcl6

Bio SB Catalog Number
BSB 5083

The bcl-6 rabbit monoclonal antibody is used in the immunohistochemical detection of bcl-6 transcriptional regulator gene on FFPE and frozen specimens.

Bio SB offers a complete line of IHC, ICC, CISH and FISH technologies for the molecular pathology laboratory. Not sure which technologies fit your needs? Contact Bio SB to help you implement the right technologies.

TintoDetector Features and Applications



The Bio SB TintoDetector system is a capillary gap based system which can be used for Immunohistochemistry (IHC), Immunocytochemistry (ICC), Immunofluorescence (IF), Fluorescent in-situ hybridization (FISH), and Chromogenic in-situ hybridization (CISH) applications. The TintoDetector is an open system, and reagents from any supplier can be used.



The Bio SB TintoDetector is comprised of several components, which can be found in the diagram below (Figure 1.1). The components of the TintoDetector include...

- | | | | |
|---|----------|--|----------|
| 1. TintoDetector System (Qty: 1)..... | BSB 7000 | 5. TintoDetector Slide Holder (Qty: 1, See Figure 1.2) | BSB 7003 |
| 2. TintoDetector Incubator (Qty: 1) | BSB 7002 | 6. Probe On Plus Slides (Box of 72)..... | BSB 7006 |
| 3. TintoDetector 30-Well Reagent Holder (Qty: 5)..... | BSB 7004 | 7. Probe On Slides (Box of 72)..... | BSB 7007 |
| 4. TintoDetector Absorbent Pads (Qty: 10) | BSB 7036 | 8. Plastic Staining Dish (Qty: 8) | BSB 7009 |

Figure 1.1





Figure 1.2



TintoDetector Features and Applications



TintoDetector Components

Incubator

See Figure 2.1

The TintoDetector Incubator is capable of reaching temperatures up to 110°C, and is capable of storing 3 temperature presets. The incubator is used to apply varying temperatures typically used in CISH and FISH protocols that involve heat retrieval, denaturing of nucleic acids, probe hybridization, stringency washes, and other steps related to CISH/FISH implementation.

Figure 2.1



TintoDetector Slide Holder

See Figure 2.2

The TintoDetector Slide Holder is an extremely durable capillary gap slide holder that is capable of holding 20 capillary gap slides, and easily fits into the TintoDetector Incubator. The TintoDetector Slide holder can be used for all IHC, ICC, IF, CISH, and FISH applications.

Figure 2.2



TintoDetector 30 Well Reagent Holder

See Figure 2.3

The TintoDetector 30-Well Reagent Holder allows for the application of up to 200 microliters of reagent to a paired set of slides. The 30-well reagent holders can be used to apply any reagent used in IHC, ICC, IF, CISH and FISH protocols.

Figure 2.3



TintoDetector Features and Applications



TintoDetector Components, Continued

Figure 2.4



Staining Dish Rack

See Figure 2.4

All TintoDetector staining dishes are capable of holding xylene and alcohol as well as any washes, buffers, or special stains for use in IHC, ICC, IF, CISH, and FISH protocols. All staining dishes are capable of holding 200 mL of reagent.

TintoDetector Specifications

| | |
|---|---|
| TintoDetector Dimensions..... | 36in. x 16in. x 11in. (91cm. x 41cm. x 28cm.) |
| TintoDetector Weight..... | 31 pounds (14kg) |
| TintoDetector Incubator Voltage Requirements..... | 110 Volt / 220 Volt |
| Supported Protocols..... | IHC, ICC, IF, CISH and FISH |

Guidelines

Ensure that TintoDetector is placed in a well ventilated area, as the built in incubator must have access to free air-flow to ensure proper operating temperatures are reached. When operating incubator, ensure that humidity chamber (located on back left of unit) is filled with up to 10 mL of distilled water.

Ensure unit is placed on a level surface.

Always use proper safety guidelines when working with toxic and flammable reagents in your laboratory. All protocols listed within the TintoDetector manual are guidelines only, and are meant as a sample application of the TintoDetector. Always reference the supplier protocol before using the TintoDetector.

TintoDetector Features and Applications




TintoDetector Slide Holder

See Figure 3.2

The TintoDetector System uses a specially designed capillary gap slide with raised white triangles at the bottom of the slide (see figure 3.2). When using the TintoDetector Slide Holder, ensure that the following procedures are followed . . .

- Slides are paired face to face.
- If a single slide or odd number must be used, pair with a blank slide.
- Insert slides so that portion of slide with white triangle faces downward when unit is held.
- Ensure that Probe On or Probe On Plus Slides are used.
- Ensure all bottom edges of slides are aligned to ensure proper capillary gap action.



Figure 3.2



TintoDetector Features and Applications

Step 1 - Slider Holder

Load slides in TintoDetector Holder, face to face and properly ordered. Make sure the raised white triangles are facing each other.





Step 1 - Reagent

Apply reagent to TintoDetector 30-well Reagent Holder. Each reagent well can hold about 200 μ l of reagent.

Step 2 - Application


Place TintoDetector slide holder over 30-well reagent holder, ensuring that reagents line up with slides. Press slide holder gently against reagent holder. Capillary gap action will draw reagent. Transfer TintoDetector Slide Holder to Incubator.



Step 3 - Rinse


After the reagent incubation, eliminate the used reagent into an absorbent pad by gently pressing the slides on the absorbent pad to drain the liquid. Draw washing buffer into the capillary space and repeat the washing process 3 to 5 times. After washing proceed to draw the next step reagent before another incubation.

TintoDetector Mini



The Bio SB TintoDetector Mini is a compact molecular pathology system which allows for the full implementation of immunohistochemical (IHC), Immunocytochemistry (ICC) and Immunofluorescence (IF) protocols. Ideal for use in diagnostic or research laboratories looking for a cost-effective, semi-automated and open platform.

The Bio SB TintoDetector system consists of a durable stainless steel housing, TintoDetector Slide Holder, four solution dishes, 5 30-well reagent holders and 1 package of TintoDetector Pads. With a compact footprint, ease of use and versatility, the TintoDetector offers flexible protocol implementation at a competitive price.



TintoDetector Mini Systems and Accessories

| Product Description | Catalog # |
|---|-----------|
| TintoDetector Mini Complete System | BSB 7085 |
| TintoRetriever Pressure Cooker | BSB 7015 |
| TintoDetector Slide Holder for Cap Gap Slides | BSB 7003 |
| TintoDetector 30-well Reagent Holder, Pack of 5 | BSB 7004 |
| TintoDetector Absorbent Pads, Box of 10 | BSB 7036 |
| TintoDetector Cap Gap Plus Slides, Box of 72 | BSB 7006 |

TintoRetriever Pressure Cooker



Affordable, Precision Heat Epitope and Nucleic Acid Retrieval


Bio SB Epitope & Nucleic Acid Retrieval

Bio SB presents a wide array of epitope and nucleic acid retrieval solutions for clinical and research pathology labs worldwide. All of our ancillaries and nucleic acid/epitope retrieval solutions are cost effective, quality reagents that are widely accepted in the molecular pathology industry. Additionally all products are manufactured according to ISO 13485:2016 and FDA cGMP 820 guidelines.

Bio SB TintoRetriever Pressure Cooker

The Bio SB TintoRetriever Pressure Cooker is a precision controlled heat source capable of maintaining a constant and reliable temperature while minimizing the potential for evaporation of the working solution. The TintoRetriever Pressure cooker can be used for both ISH and IHC applications.

Most formalin fixed tissues require an antigen or nucleic acid retrieval step prior to immunochemical or insitu hybridization staining. The retrieval process breaks some of the methylene bridges that cross-link proteins and nucleic acids in formalin fixed tissue and allows antibodies or DNA/RNA probes to bind to unmasked epitopes or nucleic acids. Conventional methods of epitope retrieval such as microwaves or water baths lead to non-uniform heat distribution or use a large amount of retrieval solution. The TintoRetriever pressure cooker is a solution for laboratories looking for a comprehensive, cost effective, and reliable retrieval system.



TintoRetriever Pressure Cooker



- Cost effective
- Fast epitope retrieval
- Uniform and reproducible heat retrieval
- Built in temperature gauge for validation
- Supports CISH/FISH nucleic acid retrieval
- Supports IHC, ICC and IF epitope retrieval
- Accessories and replacement parts available!



Above: TintoRetriever Pressure Cooker Bundle, Includes 4 Plastic Staining Dishes, 4 Slide Holders, 1 Staining Dish Support, 200 ml ImmunoRetriever 20X Citrate, 200 ml ImmunoRetriever 20X EDTA.

TintoRetriever Pressure Cooker Systems

| Product Description | Catalog # |
|--|-----------|
| TintoRetriever Digital Pressure Cooker Bundle | BSB 7015 |
| TintoRetriever Digital Pressure Cooker | BSB 7008 |



TintoRetriever Pressure Cooker Accessories

| Product Description | Catalog # |
|--|-----------|
| Staining Dish Support | BSB 7086 |
| Plastic Staining Dish | BSB 7009 |
| Slide Holder, 24 places | BSB 7010 |
| TintoRetriever Slide Thermometers, Pack of 5 | BSB 7005 |

Left: Staining Dish Support, holds 96 slides total!

TintoRetriever Pressure Cooker



Using the Bio SB TintoRetriever Pressure Cooker in the Laboratory

The TintoRetriever is a rapid epitope and nucleic acid retrieval system which can be used in a variety of applications in the modern pathology lab.

The Bio SB TintoRetriever comes pre-programmed with 6 temperature settings that can be used for deparaffinization, epitope retrieval, or nucleic acid retrieval (ISH).

IHC Epitope Retrieval & ISH Nucleic Acid Retrieval 1

For formalin-fixed paraffin-embedded (FFPE) tissue sections, the High Pressure setting with a time of 15 minutes is recommended for thorough heat permeabilization of tissue. We recommend using Bio SB Immuno/DNA Retriever with Citrate or Immuno/DNA Retriever with EDTA for IHC/ISH protocols.

Tissue Microarrays 2

Tissue Microarrays (or TMA's) typically need gentler epitope retrieval methods than those used with whole tissue sections. At Bio SB, our 11-core NH-TMA (BSB 0297) and 23-core NH-TMA's (BSB 0298) undergo epitope retrieval using the **100°C Setting**. Additionally, the Bio SB TintoRetriever PT module has gentler epitope retrieval settings (BSB 7030), and may be ideal for labs that utilize TMA's.

Deparaffinization, Retrieval and Hydration #1 or #2 3


For non-toxic deparaffinization, retrieval and hydration, we recommend the Bio SB TintoDeparaffinator Citrate or EDTA 20X, as substitutes of solvents like xylenes, toluenes and alcohols.

Tissues prone to detachment

Should tissues detach after using the "High Pressure" setting, it is recommended that tissues are mounted using Bio SB Hydrophilic Plus Slides (BSB 7028). Hydrophilic Plus Slides prevent tissue detachment while promoting reagent dispersion.

Time

| | Low Pressure | High Pressure |
|---|--------------|---------------|
| 1 | 106°-110° C | 114°-121° C |
| 2 | 100° C | 100-112° C |
| 3 | 80-86° C | 58-62° C |



Tissue Detachment on a Hydrophobic Slide (Left) vs. Hydrophilic Slide (Right)



Hydrophilic Plus Slides for Molecular Pathology




The Hydrophilic Plus Slides are novel positively charged hydrophilic slides that prevent tissue detachment after thermal permeabilization and prior to IHC, ICC, CISH or FISH protocols.

- The Hydrophilic Plus glass slides carry approximately three-times the number of positive charges compared to the commercial positively-charged slides
- The Hydrophilic Plus glass slides are strongly hydrophilic
- The Hydrophilic Plus glass slides display improved tissue-adhesion characteristics compared to other commercially available slides
- All commercially available positively-charged microscope slides were found to be hydrophobic.


Characteristics of Hydrophilic Slides

Water Distribution and Surface Contact Angles for Untreated, Hydrophobic and Hydrophilic Slides



A 200ul drop of colored TBS was spread over the total working area of the different microscope slides (estimated at 1100mm²). The TBS spread over 85% of the working surface of untreated microscope slides and over 86.5% of the Bio SB Hydrophilic Plus slides. In contrast, TBS was only able to cover 15% of working surface of the Prob On Plus slides. These results demonstrated the hydrophilic nature of the Hydrophilic Plus slide surface.

WATER TO SLIDE SURFACE CONTACT ANGLE




Hydrophobicity and hydrophilicity were assessed by measuring the contact angles of water drops on microscope slide surfaces; the smaller the contact angle, the greater the hydrophilicity. Image analyses of contact angles produced the following results:

- Color Frost (untreated) contact angles: 24.90 +/- 2.170
- Probe On Plus (positively charged) contact angles: 37.80 +/- 2.330.
- Hydrophilic Plus slides contact angles: 15.80 +/- 0.70

Results showed that the Hydrophilic Plus slides had significantly smaller contact angles.


Density of Positive Charges and Intensity of Amino Groups on Non-charged, Hydrophobic and Hydrophilic Slides Measured by X-Ray Photoelectron Spectroscopy



Hydrophilic Plus Slides for Molecular Pathology




Relative Percentage of Amino Groups on Untreated Hydrophobic and Hydrophilic Slides Measured by X-Ray Photoelectron Spectroscopy



Hydrophilic Plus Slides are prepared by covalent coupling of positively charged amino groups directly to the slide. The relative concentration of amino groups on various slides was measured by X-ray spectroscopy. Results showed that Hydrophilic Plus slides contained approximately 3-times the number of amino groups compared to hydrophobic slides and greater than 30-times the number of amino groups measured on untreated slides. The high density of amino groups confers hydrophilic properties to the coated slides.

The highest density of positive charges on the surface of hydrophilic slides increases their wetting ability when compared to standard hydrophobic positively charged slides. Increased hydrophilic interaction on the slide surface reduces nonspecific adsorption of hydrophobic agents, and improves tissue adhesion to the slides.

Water and Air Bubbles Trapping, Tissue Damage and Wrinkling after Mounting Tissues on Hydrophobic and Hydrophilic Slides



Reduced water to slide surface contact angles for hydrophilic slides (when compared to a standard hydrophobic charged slides) allows water to form a single thin continuous layer rather than a layer of water drops separated by air bubbles. This prevents micro air bubbles from being captured under tissue sections and in turn reduces the development of possible artifacts, such as bubbles or water drops, which can interfere with tissue attachment to the slide. Generally, a more uniform and thinner layer of water under the cut tissue section results in less trapping of bubbles and tissue wrinkling, and facilitates faster drying after mounting on Hydrophilic Plus slides.

Percentage of Average Tissue Retention and IHC Signal with the Different Experimental Slides

| Slide Name | Brand | Tissue Retention | IHC Signal Background |
|------------------|-------------|------------------|-----------------------|
| Hydrophilic Plus | Bio SB, Inc | 90 - 100% | 3- 4 / + |
| Superfrost Plus | Company A | 10- 50% | 3- 3+ / + |
| Probe On Plus | Company A | 10- 50% | 3- 3+ / + |
| Color Mark | Company B | 5 - 40% | 3- 3+ / + |
| Snow Coat X-tra | Company C | 5 - 30% | 3- 3+ / + |
| Silanized | In-House | 5 - 40% | 3- 3+ / + |

Tissue adhesion **after** was measured on various types of microscope slides. The results showed that Hydrophilic Plus slides had the highest percentage of tissue retention (90-100%).


Hydrophilic Plus Slides for Molecular Pathology



Macroscopic Comparison of Tissue Retention and IHC of Ki67 and CD31 Using ColorMark and Hydrophilic Plus Slides



Microscopic Comparison of Tissue Retention and IHC of Ki67 and CD31 Using Probe On Plus and Hydrophilic Plus Slides




Stable and strong adhesion of biological samples to microscope slide surfaces is important in achieving successful sample preparation and staining. Untreated glass surfaces or positively-charged modified surfaces of microscope slides often do not provide strong enough retention for some biological samples. This can lead to a full or partial sample loss or sample deformation. The Hydrophilic Plus Slides with highly positively charged hydrophilic surfaces, have shown to be superior to other slides with traditional hydrophobic positively charged surfaces for tissue retention.

Our studies have shown that the Hydrophilic Plus Slides are suitable for IHC, ICC, CISH and FISH and are superior in their capacity to retain tissues that otherwise tend to detach from slides after thermal antigen retrieval procedures. In general, 90 to 100% of tissues that were damaged or detached from traditional hydrophobic slides were fully retained by the hydrophilic slides without affecting tissue morphology and the quality of IHC signals.

Hydrophilic Plus and Cap Gap Plus Slides for Molecular Pathology



Microscopic Comparison of IHC Signal Intensity of Melanoma HMB-45 on Superfrost Plus and Hydrophilic Plus Slides



An interesting and unexpected observation surfaced during the course of our studies – IHC signals were more intense when using Hydrophilic Plus Slides. Although the reasons for this finding are unknown, we are investigating the possibility that hydrophilic surfaces improve efficiencies of immunological reactions.

Hydrophilic Plus Slides are very effective tools in handling tissues that are likely to sustain damage, loss, or detachment after thermal antigen retrieval procedures.

Hydrophilic Plus Slides have shown to be reliable and superior alternatives to other slides with traditional hydrophobic positively charged surfaces, and are another effective tools when handling tissues that tend to get damaged, lost, or detached after thermal antigen retrieval procedures prior to IHC, ICC, CISH and FISH procedures.

Hydrophilic Plus Slides & TintoDetector Cap Gap Plus Slides

| Product Description | Catalog # |
|---|-----------|
| Hydrophilic Plus Slides (Box of 100) | BSB 7028 |
| TintoDetector Cap Gap Plus Slides (Box of 72) | BSB 7006 |

The Bio SB TintoDetector Cap Gap Plus Slides are novel positively charged hydrophilic slides that prevent tissue detachment after thermal permeabilization and can be used to mount and process IHC, ICC, FISH and CISH samples for use in the TintoDetector Cap Gap System (BSB 7000). Bio SB TintoDetector Cap Gap Plus Slides feature a 75 μ m raised triangle on the bottom corners of the slide that creates a 150 μ m gap when placed face-to-face in the TintoDetector Cap Gap Slide Holder (BSB 7003).

The Bio SB TintoDetector Cap Gap Plus Slides are prepared by covalent coupling of positively charged groups directly to the slide surface. The Bio SB TintoDetector Cap Gap Plus Slides carry approximately three times the number of surface-bound positive charges compared to other commercially available positively charged slides.

This results in significantly improved tissue-adhesion characteristics, facilitates uniform/reproducible staining, reduces staining artifacts/background and allows better capillary flow of aqueous reagents and solutions.





Molecular Pathology Technical Reference



Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|-----------------------|---------|------------------|---------|--|---|
| A-1-Antichymotrypsin | RMab | EP384 | IgG | Purified human a-1-antichymotrypsin | Human |
| A-1-Antichymotrypsin | RPab | Polyclonal | IgG | A synthetic peptide corresponding to residues of human Alpha1-Antichymotrypsin protein | Human |
| A-1-Antitrypsin | RPab | Polyclonal | IgG | Purified human serum Alpha-1-Antitrypsin | Human, Baboon, Equine, Mink |
| ACE2 | MMab | BSB-135 | IgG1 | Synthetic Peptide corresponding the C-terminal of ACE2 of human origin | Human, Mouse, Rat |
| ACTH | MMab | BSB-25 | IgG1/K | Synthetic peptide corresponding to the N-Terminus of human ACTH | Human |
| Actin Muscle Specific | MMab | HHF35 | IgG1/K | SDS extract of human myocardium | Human |
| Actin Smooth Muscle | MMab | BSB-15 (ASM/H12) | IgG2a/K | Synthetic peptide corresponding to the N-terminus of human smooth muscle actin | Human, Mouse, Rabbit, Dog, Cat, Rat, Sheep, Pig |
| Adenovirus | MMab | 20/11 and 2/6 | IgG1/K | Adeno 3 | Human |
| Adipophilin/ADRP | MMab | BSB-91 | IgG1 | Synthetic peptide corresponding to the C-terminus of human Adipophilin | Human |
| Albumin | RPab | polyclonal | IgG | Recombinant protein corresponding to the N-terminus of the human serum albumin protien. | Human |
| ALDH1A1 | RMab | EP168 | IgG | Synthetic peptide corresponding to residues of human Aldh1A1 protein | Human |
| ALK-1/CD246 | RMab | EP302 | IgG | A synthetic peptide corresponding to residues of human NPM-ALK fusion protein | Human, Predicted: Mouse, Rat |
| ALK-1/CD246 | RMab | RBT-ALK1 | IgG | Recombinant protein corresponding to the tyrosine kinase catalytic region of the human ALK protein | Human |
| Alpha Synuclein | MMab | BSB-114 | IgG2a/K | Recombinant protein of the human alpha synuclein | Human |
| Alpha-Fetoprotein | RMab | EP209 | IgG | Synthetic peptide corresponding to residues of the human Alpha-Fetoprotein | Human |
| Alpha-Fetoprotein | MMab | BSB-23 | IgG2a/K | Synthetic peptide corresponding to amino acids 199-161 of human alpha-fetoprotein | Human |
| AMACRacemase/P504S | RMab | 13H4 | IgG | Synthetic peptide corresponding to the human p504s/AMACR protein | Human |
| Amyloid A | RMab | EP335 | IgG | A synthetic peptide corresponding to residues of human Serum Amyloid A protein | Human |
| Amyloid Beta | RMab | RBT-A4 | IgG | Synthetic peptide corresponding to the C-terminus of the human Amyloid Beta protein | Human, Mouse, Rat |
| Androgen Receptor | MMab | BSB-4 (AR-D12) | IgG1 | Synthetic peptide corresponding to residues in the N-terminus of human Androgen Receptor protein | Human |
| Annexin A1 | MMab | BSB-95 | IgG2b | Recombinant Human Annexin A1 | Human, Rat |
| Annexin VII | RMab | EP367 | IgG | Recombinant human Annexin VII protein | Human, Mouse, Rat |
| Arginase-1 | RMab | EP261 | IgG | Synthetic peptide corresponding to residues of the human ARG-1 protein | Human |
| ARID1A | RMab | EP303 | IgG | A synthetic peptide corresponding to residues of human ARID1A (BAF250a) protein | Human, Mouse, Rat |
| ATM | RMab | EP327 | IgG | Synthetic peptide corresponding to residues of the human ATM protein | Human |
| ATRX | MMab | BSB-108 | IgG2a/K | Synthetic peptide corresponding to the C-terminus of the human ATRX protein | Human, Mouse, Rat |
| ATRX | RMab | RBT-ATRX | IgG | Recombinant human ATRX protein | Human, Mouse, Rat |
| Aurora B | RMab | RM278 | IgG | Synthetic peptide corresponding to the N-terminus of human Aurora kinase B | Human |
| B7H3/CD276 | RMab | RBT-B7H3 | IgG | Synthetic peptide conjugated to KLH corresponding to the C-terminal residues of the human B7H3/CD276 protein | Human, Predicted: Mouse |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|----------------------------------|----------------|--|
| Tonsil, Lymph Node, Thymus, Spleen, Liver, Colon | Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic, Liver Cancer |
| Tonsil, Lymph Node, Thymus, Breast | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Liver Cancer |
| Tonsil, Lymph Node, Spleen, Pancreas, Myometrium | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Liver Cancer |
| Kidney, Testis, Brain, Colon, Fallopian tube | Membranous, Cytoplasmic | 1:25-1:100 | Infectious Disease, Breast Cancer, Lung Cancer, Prostate Cancer, Liver Cancer, Kidney & Urotelial Cancer, Colon & G.I. Cancer |
| Normal Pituitary | Cytoplasmic | 1:250-1:1000 | Pituitary, Neural & Neuroendocrine Cancer |
| Skeletal Muscle, Appendix, Prostate | Cytoplasmic | 1:50-1:200 | Sarcoma & Soft Tissue, Undifferentiated Tumor |
| Appendix, Uterus | Cytoplasmic | 1:250-1:1000 | Breast Cancer, Colon & GI Cancer, Sarcoma & Soft Tissue |
| Adenovirus Infected Tissue | Cytoplasmic, Nuclear | 1:25-1:100 | Infectious Diseases |
| Adrenal, SCC, TCC and Sebaceous Neoplasms | Membranous, Cytoplasmic | 1:50-1:200 | Melanoma & Skin Cancer |
| Salivary Gland, Kidney, Tonsil, Lupus Erythematosus | Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity |
| Kidney, Liver, Testis, Colon Cancer, Breast Cancer | Cytoplasmic | 1:100-1:500 | Ovarian Cancer, Prostate Cancer, Breast Cancer, Colon & GI Cancer, Lymphoma, Leukemia & Histiocytic, Sarcoma & Soft Tissue, Endothelial |
| Anaplastic Large Cell Lymphoma | Cytoplasmic, Nuclear | 1:50-1:200 | Lymphomas, Lung Cancer |
| Anaplastic Large Cell Lymphoma | Cytoplasmic, Nuclear | 1:100-1:500 | Lymphomas, Lung Cancer |
| Brain, Breast, Skin, Tonsil, Bone Marrow, Alzheimer's Disease, Parkinson's Disease | Cytoplasmic, Nuclear | 1:25-1:100 | Neural & Neuroendocrine Cancer |
| Fetal Liver, Liver Lesions, Hepatocellular Carcinoma | Cytoplasmic | 1:50-1:200 | Germ Cell Tumor, Liver Cancer, Undifferentiated Tumor |
| Fetal Liver, Liver Lesions, Hepatocellular Carcinoma | Cytoplasmic | 1:100-1:500 | Germ Cell Tumor, Liver Cancer, Undifferentiated Tumor |
| Kidney, Liver, Salivary Gland, Prostate Lesions, Prostatic Adenocarcinoma | Cytoplasmic | 1:25-1:100 | Prostate Cancer, Kidney & Urotelial Cancer |
| Kidney, Amyloidosis | Extracellular, Cytoplasmic | 1:50-1:200 | Kidney & Urotelial Cancer, Rejection & Autoimmunity |
| Testis, Kidney, Pancreas, Salivary Gland, Alzheimer's Disease | Cytoplasmic, Nuclear | 1:50-1:200 | Neural & Neuroendocrine Cancer |
| Prostate, Testis, Breast, Cervix, Fallopian Tube, Prostatic Adenocarcinoma | Nuclear | 1:100-1:500 | Prostate Cancer, Breast Cancer, Head & Neck Cancer |
| Liver, Tonsil, Spleen, Thymus, Lung, Colon, Hairy Cell Leukemia | Cytoplasmic, Membranous | 1:100-1:500 | Leukemia & Histiocytic |
| Breast, Colon, Kidney, Tonsil, Prostate, Testis, Transitional Cell Carcinoma | Nuclear, Cytoplasmic, Membranous | 1:50-1:200 | Breast Cancer, Prostate Cancer, Liver Cancer, Colon & G.I. Cancer |
| Liver, Hepatocellular Carcinoma | Cytoplasmic, Nuclear | 1:50-1:200 | Liver Cancer, Carcinoma of Unknown Primary Site |
| Breast, Fallopian Tube, Prostate, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Nuclear | 1:25-1:100 | Ovarian Cancer, Breast Cancer, Colon & G.I. Cancer |
| Fallopian Tube, Brain, Colon, Breast, Testis, Tonsil, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Breast Cancer, Gastric Cancer, Lung Cancer, Lymphoma, Leukemia & Histiocytic |
| Breast, Tonsil, Testis, Salivary Gland, Placenta, TCC, Astrocytoma | Nuclear | 1:50-1:200 | Neural & Neuroendocrine Cancer |
| Breast, Tonsil, Testis, Salivary Gland, Placenta, TCC, Astrocytoma | Nuclear | 1:50-1:200 | Neural & Neuroendocrine Cancer |
| Tonsil, Colon, Stomach, Skin, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Prostate Cancer, Breast Cancer, Colon and GI Cancer, Thyroid & Parathyroid Cancer, Neural & Neuroendocrine Cancer, Lung Cancer, Liver Cancer |
| Testis, Adrenal, Tonsil, Breast, Fallopian Tube, Breast Carcinoma, Prostate Carcinoma and Ovarian Carcinoma | Membranous | 1:25-1:100 | Melanoma & Skin Cancer, Prostate Cancer, Gall Bladder & Pancreatic Cancer, Breast Cancer, Immunotherapy |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|------------------------|---------|-----------------|---------|--|---|
| BAP1 | MMab | BSB-109 | IgG1 | Synthetic peptide corresponding to the C-terminus of the human BAP1 protein | Human, Mouse, Rat |
| Bax | RMab | E63 | IgG | Synthetic peptide within Human Bax aa 1-100 (N terminal) | Human |
| BCA-225 | MMab | Cu-18 | IgG1/K | BCA 225 protein secreted by the T47D (clone 11) human breast carcinoma cell line | Human |
| bcl-10 | MMab | BSB-22 | IgG1/K | Recombinant human bcl-10 protein | Human |
| bcl-2 | MMab | BSB-5 (BCL2/A4) | IgG1/K | Synthetic peptide corresponding to residues in the N-terminus of the human bcl2 protein | Human |
| bcl-2 | RMab | EP36 | IgG | Synthetic peptide corresponding to residues between BH3 and BH4 of human Bcl-2 protein | Human, Predicted: Mouse |
| bcl-6 | RMab | RBT-bcl6 | IgG | Synthetic peptide corresponding to residues of the C-terminus of the human bcl-6 protein | Human |
| bcl-6 | RMab | EP278 | IgG | Synthetic peptide corresponding to residues in human bcl6 protein | Human, Predicted: Mouse, Rat |
| bcl-6 | MMab | BSB-26 | IgG1 | Synthetic peptide corresponding to residues of the C-terminus of the human bcl-6 protein | Human |
| Bcl-x | RMab | EP94 | IgG | Synthetic peptide corresponding to residues in human Bcl-x protein | Human, Predicted: Mouse, Rat |
| BCOR | MMab | BSB-128 | IgG1/K | Synthetic peptide corresponding to residues of the N-terminus of the human BCoR protein | Human |
| Beta-Catenin | RMab | RM276 | IgG | A peptide corresponding to human Beta-Catenin. | Human, Predicted: Mouse, Rat, Sheep, Hamster, Cow, Macaque Monkey, African Green Monkey |
| Beta-Catenin | MMab | 14 | IgG1 | Synthetic peptide conjugated to KLH derived from within residues 750 to the C-terminus of Human beta catenin | Human, Dog, Mouse, Rat, Chicken |
| BG8 LewisY | MMab | F3 | IgM | SK-LU-3 lung cancer cell line (Human) | Human |
| BOB.1 | RMab | RBT-BOB1 | IgG | Synthetic peptide corresponding to residues of the human BOB.1 protein | Human |
| Brachyury | RMab | RBT-TBXT | IgG | Synthetic peptide corresponding to the C-terminus of the human brachyury protein | Human |
| BRAF V600E | RMab | RM8 | IgG | A peptide corresponding to BRAF V600E mutant | Human |
| BRG-1/SMARCA4 | MMab | BSB-154 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human Brg-1 protein | Human, Mouse, Rat |
| C1q | RPab | Polyclonal | IgG | KLH conjugated synthetic peptide corresponding to the C-terminus region of human C1QA. | Human |
| C3c | RPab | Polyclonal | IgG | Purified C3c protein isolated from normal human serum. | Human |
| C3d | RPab | Polyclonal | IgG | Purified full length native protein corresponding to the human C7 protein | Human |
| C4c | RPab | Polyclonal | IgG | Recombinant protein corresponding to the human Complement Component 4c protein. | Human |
| C4d | RMab | EP272 | IgG | Prokaryotic recombinant protein corresponding to the C4d protein | Human |
| CA-125 | MMab | OC125 | IgG1/K | Partially purified human mucin fraction from a pool of tissues from patients with epithelial ovarian cancer | Human |
| CA-125 | RMab | EP48 | IgG | Synthetic peptide corresponding to residues of the human CA125 protein | Human |
| CA15-3 | MMab | DF3 | IgG1/K | Membrane enriched fractions of human metastatic breast carcinoma | Human |
| CA19-9 | MMab | 121SLE | IgM | Purified mucins from ovarian cyst | Human |
| Cadherin17/LI-Cadherin | RMab | EP86 | IgG | Synthetic peptide corresponding to residues of the human LI-Cadherin protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|----------------------------------|----------------|--|
| Testis, TCC, Mesothelioma | Nuclear | 1:50-1:200 | Mesothelioma, Lung Cancer, Melanoma & Skin Cancer, Kidney & Urotelial Cancer |
| Breast, Tonsil, Cervix, Hodgkin's Lymphoma | Cytoplasmic, Cell Membranous | 1:50-1:200 | Breast Cancer, Ovarian Cancer, Lung Cancer |
| Breast, Lung, Uterus, Cervical Carcinoma | Cytoplasmic | 1:25-1:100 | Breast Cancer, Cytopathology |
| Tonsil, Spleen, Thymus, MALT Lymphomas | Cytoplasmic, Nuclear | 1:100-1:500 | Lymphoma, Gall Bladder & Pancreatic Cancer |
| Tonsil, Lymph Node, Breast, Bone Marrow, Fallopian Tube | Cytoplasmic, Membranous | 1:50-1:200 | Lymphoma, Endometrial & Genital Cancer, Prostate Cancer, Breast Cancer, Lung Cancer |
| Tonsil, Lymph Node, Breast, Placenta, Fallopian Tube | Cytoplasmic, Membranous | 1:50-1:200 | Lymphoma, Endometrial & Genital Cancer, Prostate Cancer, Breast Cancer, Lung Cancer |
| Tonsil, Lymph Node, Thymus, Skin, Breast, Brain, Follicular Lymphoma | Nuclear | 1:50-1:200 | Hodgkin's and NHD Lymphoma, Lymphoma, Gall Bladder and Pancreatic Cancer |
| Tonsil, Follicular Lymphoma | Nuclear | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Lymphoma, Gall Bladder and Pancreatic Cancer |
| Tonsil, Lymph Node, Thymus, Skin, Breast, Brain, Follicular Lymphoma | Nuclear | 1:50-1:200 | Hodgkin's and NHD Lymphoma, Lymphoma, Gall Bladder and Pancreatic Cancer |
| Kidney, Tonsil, Cervix, Hodgkin's Lymphoma | Cytoplasmic, Membranous | 1:250-1:1000 | Breast Cancer |
| Testis, Cervix, Prostate, TCC, Angiosarcoma | Nuclear | 1:10-1:25 | Sarcoma & Soft Tissue |
| Fibromatosis of the Breast & Abdomen. Breast, Abdomen, Colon, Testis, Pancreas | Cytoplasmic, Membranous, Nuclear | 1:50-1:200 | Breast Cancer, Colon & Gastrointestinal Cancer, Liver Cancer, Gall Bladder & Pancreatic Cancer, Sacoma & Soft Tissue |
| Fibromatosis - Breast Ca/Abdomen | Cytoplasmic, Membranous, Nuclear | 1:50-1:200 | Breast Cancer, Colon & GI Cancer, Liver Cancer, Gall Bladder & Pancreatic Cancer, Sacoma & Soft Tissue |
| Placenta, Tonsil, Pancreas, Cervix, Bladder TCC, Lung Adenocarcinoma | Cytoplasmic | 1:10-1:50 | Lung Cancer, Mesothelioma |
| Tonsil, Lymph Node | Cytoplasmic, Nuclear | 1:250-1:1000 | Hodgkin's and NHD Lymphoma |
| Pituitary, Brain, Chordoma | Cytoplasmic | 1:25-1:100 | Neural & Neuroendocrine Cancer, Liver Cancer |
| BRAF 600E Mutated Melanoma, Papillary Thyroid Cancer | Cytoplasmic | 1:50-1:200 | Thyroid & Parathyroid Cancer, Colon & GI Cancer, Melanoma & Skin Cancer, |
| Colon, Kidney, Prostate, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Nuclear | 1:25-1:100 | Lung Cancer, Ovarian Cancer, Breast Cancer, Melanoma & Skin Cancer, Colon & GI Cancer, Prostate Cancer |
| Kidney, Cervix, Spleen, Lupus Erythematosus | Cytoplasmic | 1:25-1:100 | Rejection & Autoimmunity |
| Placenta, Kidney, Fallopian Tube, Lupus Erythematosus | Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity |
| Kidney Transplant Rejection | Cytoplasmic, Membranous | 1:50-1:200 | Kidney & Urotelial Cancer, Rejection & Autoimmunity |
| Teestis, Kidney, Pancreas, Salivary Gland, Colon | Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity |
| Tonsil, Lymph Node, Spleen, Kidney, Kidney Transplant Rejection | Cytoplasmic, Membranous | 1:50-1:200 | Kidney & Urotelial Cancer, Rejection & Autoimmunity |
| Colon, Pancreas, Epithelioid Mesothelioma, Ovarian Carcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Breast Cancer, Cervical Cancer, Ovarian Cancer |
| Colon, Pancreas, Epithelioid Mesothelioma, Ovarian Carcinoma, | Cytoplasmic, Membranous | 1:50-1:200 | Breast Cancer, Cervical Cancer, Ovarian Cancer |
| Kidney, Breast, Lymph Node, Cervix, Salivary Gland, Bladder TCC | Cytoplasmic, Membranous | 1:250-1:1000 | Breast Cancer, Lung Cancer, Sarcomas & Soft Tissue |
| Colon, Cervix, Pancreas, Pancreas Cancer, Breast, Breast Carcinoma, Colon Carcinoma, Transitional Cell Carcinoma, Ovarian Carcinoma, Thyroid Carcinoma | Cytoplasmic | 1:100-1:500 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer |
| Colon Carcinoma | Cytoplasmic | 1:25-1:100 | Ovarian Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Carcinomas of Unknown Primary Site |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|----------------------|---------|-------------------|---------|---|------------------------------------|
| Cadherin-6 | RMab | EP217 | IgG | Synthetic peptide corresponding to residues of the human CDH6 protein | Human, Predicted: Mouse, Rat |
| Calcitonin | RMab | EP92 | IgG | Synthetic peptide corresponding to residues of the human calcitonin protein | Human |
| Calcitonin | RPab | Polyclonal | IgG | Synthetic peptide corresponding to the human calcitonin protein | Human, Dog, Cat |
| Caldesmon | MMab | BSB-19 (CALD-31) | IgG1/K | Synthetic peptide corresponding to the N-terminus of the human caldesmon protein | Human |
| Calponin | MMab | BSB-20 (CALP-A6) | IgG1/K | Synthetic peptide corresponding to residues of the C-terminus of the human calponin protein | Human, Canine, Feline |
| Calretinin | RMab | EP1798 | IgG | Synthetic peptide corresponding to residues in Human Calretinin | Human, Mouse, Rat |
| Calretinin | RMab | RM324 | IgG | A peptide corresponding to N-terminus of human Calretinin. | Human, Mouse, Rat |
| Carbonic Anhydrase 9 | RMab | EP161 | IgG | Synthetic peptide corresponding to residues in the extracellular domain of the human Carbonic Anhydrase 9 protein | Human |
| Caspase-3 | RMab | RM250 | IgG | Synthetic peptide corresponding to residues within the human Caspase-3 subunit p17 | Human |
| Caveolin-1 | RMab | EP353 | IgG | Synthetic peptide corresponding to residues of human Caveolin-1 protein | Human |
| CD10 | MMab | 56C6 | IgG1 | Recombinant external domain of the human CD10 protein | Human, Dog, Cat |
| CD10 | RMab | EP195 | IgG | Recombinant fragment corresponding to residues in human CD10 protein | Human, Predicted: Rat |
| CD103/ITGAE | RMab | EP206 | IgG | A synthetic peptide corresponding to residues of human CD103 protein | Human |
| CD105 | RMab | EP274 | IgG | A protein fragment corresponding to residues on human CD105 protein | Human |
| CD105/Endoglin | RPab | Polyclonal | IgG | Synthetic KLH conjugated peptide between 380-409 amino acids from the central region of human CD105 | Human, Mouse |
| CD117 | RMab | EP10 | IgG | Synthetic peptide corresponding to residues in the C-terminus of human CD117 protein | Human, Monkey, Predicted: Marmoset |
| CD117 | RMab | RM359 | IgG | A peptide corresponding to the C-terminus of human CD117/c-Kit. | Human, Monkey, Predicted: Marmoset |
| CD11b | RMab | EP45 | IgG | Synthetic peptide corresponding to residues in human CD11b protein | Human, Predicted: Mouse |
| CD11c | RMab | EP157 | IgG | Synthetic peptide corresponding to residues of human CD11c/ITGAX protein | Human |
| CD123 IL-3Ra | MMab | BSB-59 (CD123-D3) | IgG1/K | Recombinant human CD123 protein | Human |
| CD13 | MMab | 38C12 | IgG1 | Recombinant protein encoding the C-terminal half of the extracellular domain of human CD13 | Human |
| CD13 | RMab | EP117 | IgG | Synthetic peptide corresponding to residues in human CD13 protein | Human, Predicted: Rat |
| CD137/TNFRSF9 | MMab | BSB-159 | IgM | Synthetic peptide corresponding to the N-terminus of the human CD137 protein | Human |
| CD138 | RMab | EP201 | IgG | Synthetic peptide corresponding to residues of human CD138 protein | Human, Predicted: Mouse, Rat |
| CD138 | MMab | B-A38 | IgG1 | U266 human peripheral blood myeloma cell line | Human |
| CD14 | RMab | EP128 | IgG | Synthetic peptide corresponding to residues of human CD14 protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|-------------------------|----------------|---|
| Kidney, Renal Cell Carcinoma | Membranous | 1:25-1:100 | Kidney & Urotelial Cancer, Thyroid & Parathyroid Cancer |
| Thyroid Carcinoma, Medullary Carcinoma of Thyroid | Cytoplasmic | 1:100-1:500 | Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Head & Neck Cancer, Cytopathology |
| Thyroid Carcinoma, Medullary Carcinoma of Thyroid | Cytoplasmic | 1:100-1:500 | Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Head & Neck Cancer, Cytopathology |
| Appendix, Uterus, Breast Ducts, Leiomyoma, Colon, Prostate, Skin, Kidney, Myometrium | Cytoplasmic | 1:100-1:400 | Breast Cancer, Sarcoma & Soft Tissue |
| Appendix, Uterus, Breast Ducts, Leiomyoma, Prostate, Colon, Breast, Skin | Cytoplasmic | 1:50-1:200 | Breast Cancer, Sarcoma & Soft Tissue, Head & Neck Cancer |
| Brain, Testis, Colon, Benign Mesothelial Cells, Malignant Mesothelioma | Cytoplasmic, Nuclear | 1:100-1:500 | Lung Cancer, Mesothelioma, Ovarian Cancer, Cytopathology |
| Brain, Testis, Colon, Benign Mesothelial Cells, Malignant Mesothelioma | Cytoplasmic, Nuclear | 1:100-1:500 | Mesothelioma, Lung Cancer, Ovarian Cancer, Cytopathology |
| Stomach, Gallbladder, Kidney Carcinoma, Cervix Carcinoma, Lung Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Kidney & Urotelial Cancer, Lung Cancer, Colon Cancer |
| Colon, Tonsil, Testis, Fallopian Tube, Stomach, Transitional Cell Carcinoma | Nuclear, Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic, Prostate Cancer, Breast Cancer, Gastric Cancer, Ovarian Cancer, Endometrial & Genital Cancer, Colon & GI Cancer |
| Placenta, Liver, Kidney, Spleen, Lung, Mesothelioma, Ewing's Sarcoma, RCC | Membranous | 1:50-1:200 | Lung Cancer, Mesothelioma, Sarcoma & Soft Tissue, Cytopathology |
| Kidney, Tonsil, Lymph Node, Liver, Breast, Prostate, Lung, Pancreas, Salivary Gland, Thymus | Cytoplasmic, Membranous | 1:10-1:50 | Hodgkin's & NHD Lymphoma, Lymphoma, Kidney & Urotelial Cancer, Liver Cancer, Gall Bladder & Pancreatic Cancer, Endometrial & Genital Cancer, Breast Cancer |
| Kidney, Tonsil, Lymph Node, Liver, Breast, Prostate, Lung, Pancreas, Salivary Gland, Thymus | Cytoplasmic, Membranous | 1:10-1:50 | Hodgkin's & NHD Lymphoma, Lymphoma, Kidney & Urotelial Cancer, Liver Cancer, Gall Bladder & Pancreatic Cancer, Endometrial & Genital Cancer, Breast Cancer |
| Skin, Colon, Tonsil, Thymus, Spleen, Hairy Cell Leukemia | Cytoplasmic, Membranous | 1:25-1:100 | Leukemia & Histiocytic, Lymphoma, Ovarian Cancer, Colon & GI Cancer |
| Spleen, Tonsil, Cervix, Lymphoblastic Lymphoma | Cytoplasmic | 1:50-1:200 | Endothelial, Prostate Cancer, Breast Cancer, Colon & GI Cancer, Kidney & Urothelial Cancer, Lung Cancer |
| Spleen, Tonsil, Cervix, Lymphoblastic Lymphoma | Cytoplasmic | 1:25-1:100 | Endothelial, Prostate Cancer, Breast Cancer, Colon & GI Cancer, Kidney & Urothelial Cancer, Lung Cancer |
| Skin, Testis, Breast, GIST, Colon, Brain, Tonsil | Cytoplasmic, Membranous | 1:100-1:500 | GIST, Cervical Cancer, Colon & GI Cancer, Germ Cell Tumor, Head & Neck Cancer, Kidney & Urotelial Cancer, Leukemia & Histiocytic, Sarcoma & Soft Tissue, Thyroid & Parathyroid Cancer, Undifferentiated Tumor |
| Skin, Testis, Breast, GIST, Colon, Brain, Tonsil | Cytoplasmic, Membranous | 1:50-1:200 | GIST, Cervical Cancer, Colon & GI Cancer, Germ Cell Tumor, Head & Neck Cancer, Kidney & Urotelial Cancer, Leukemia & Histiocytic, Sarcoma & Soft Tissue, Thyroid & Parathyroid Cancer, Undifferentiated Tumor |
| Spleen, Leukemia | Cytoplasmic | 1:100-1:400 | Leukemia & Histiocytic |
| Bone Marrow, Spleen, Tonsil, Colon, Liver, Hairy Cell Leukemia | Cytoplasmic | 1:25-1:50 | Leukemia & Histiocytic |
| Tonsil, Lymph Node, Kikuchi-Fujimoto Disease | Cytoplasmic, Membranous | 1:25-1:100 | Lymphomas, Leukemia & Histiocytic |
| Spleen, Tonsil, Prostate, Liver | Cytoplasmic, Membranous | 1:25-1:50 | Leukemia & Histiocytic, Sarcoma & Soft Tissue, Liver Cancer |
| Spleen, Tonsil, Prostate, Liver | Cytoplasmic, Membranous | 1:50-1:200 | Leukemia & Histiocytic, Sarcoma & Soft Tissue, Liver Cancer |
| Colon, Stomach, Tonsil, Testis, Transitional Cell Carcinoma, DBC, Hepatocellular Carcinoma, Diffuse Type Gastric Carcinoma | Nuclear, Membranous | 1:25-1:100 | Hodgkin's & NHD Lymphoma, Lymphoma, Leukemia & Histiocytic, Immunotherapy, Gastric Cancer, Breast Cancer, Infectious Diseases |
| Tonsil, Liver, Kidney, Breast, Lymph Node, Cervix, Plasmacytoma, Adrenal, Skin, Colon, Lung | Membranous | 1:25-1:100 | Hematopoietic, Lymphoma, Rejection & Autoimmunity |
| Tonsil, Liver, Kidney, Breast, Lymph Node, Cervix, Plasmacytoma | Membranous | 1:50-1:200 | Hematopoietic, Lymphoma, Rejection & Autoimmunity |
| Placenta, Tonsil, Spleen, Diffuse Large B-cell Lymphoma | Cytoplasmic, Membranous | 1:25-1:100 | Leukemia & Histiocytic, Lymphoma, Lung Cancer, Sarcoma and Soft Tissue |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--|---------|----------|---------|--|---|
| CD14 | MMab | 7 | Ig2a | Recombinant protein encoding the external domain of human CD14 | Human |
| CD142/TF/Coagulation Factor III/Thromboplastin | MMab | BSB-143 | IgG2b | Synthetic peptide corresponding to the human internal region of the human TF protein | Human, Mouse, Rat |
| CD147 | MMab | BSB-137 | IgG1 | Recombinant protein corresponding to the extracellular domain of the human CD147 protein | Human |
| CD15 | RMab | EP273 | IgG | U937 cell line | Human |
| CD15 | MMab | BSB-119 | IgM | The U937 histiocytic cell line was used as the immunogen for the CD15 Leu-M1 antibody | Human |
| CD16 | RMab | EP364 | IgG | Synthetic peptide corresponding to residues of human CD16 protein | Human |
| CD163 | MMab | 10D6 | IgG1 | Recombinant protein encoding the domains 1-4 of the N-terminal region of human CD163 | Human |
| CD19 | MMab | BSB-97 | IgG1 | Recombinant CD19 protein | Human |
| CD1a | RMab | EP80 | IgG | Synthetic peptide corresponding to residues in human CD1a protein | Human |
| CD2 | MMab | AB75 | IgG1/K | Recombinant fragment encoding the external domain of the human CD2 molecule | Human |
| CD20 | MMab | L26 | IgG2a/K | Human tonsil B cells | Human, Canine, Feline |
| CD21 | RMab | EP64 | IgG | Synthetic peptide corresponding to residues on the C-terminus of human CD21 protein | Human, Predicted: Mouse |
| CD23 | RMab | EP75 | IgG | Synthetic peptide corresponding to residues in human CD23 protein | Human |
| CD23 | MMab | 1B12 | IgG | Recombinant external domain of CD23 protein | Human |
| CD25 | MMab | 4C9 | IgG2b | Recombinant protein corresponding to the external domain of the Interleukin-2 Receptor molecule | Human, Mouse |
| CD25 | RMab | RBT-CD25 | IgG | Synthetic peptide corresponding to human the C-terminus of human CD25 | Human |
| CD3 | RMab | RBT-CD3 | IgG | Synthetic peptide corresponding to residues in the cytoplasmic domain of the human CD3 protein | Human |
| CD3 Epsilon | RMab | RBT-CD3e | IgG | Synthetic peptide corresponding to residues of the epsilon chain of the human CD3 protein | Human |
| CD30 | MMab | Ber-H2 | IgG1/K | L428 cell line cells | Human |
| CD31 | MMab | 1A10 | IgG1/K | Recombinant protein corresponding to the extracellular domain downstream of the signal sequence of the CD31 molecule | Human, Dog, Cat, Mouse |
| CD33 | RMab | RBT-CD33 | IgG | Recombinant protein encoding the extracellular domain of human CD33 | Human |
| CD34 | MMab | QBEnd/10 | IgG1 | Human endothelial vesicles | Human |
| CD34 | RMab | EP88 | IgG | Synthetic peptide corresponding to C-terminal of human CD34 protein | Human, Predicted: Mouse, Rat, Sheep, Dog, Pig, Loxodonta Africana |
| CD35 | RMab | EP197 | IgG | Synthetic peptide corresponding to residues of human CD35 protein | Human |
| CD35 | MMab | RLB25 | IgG2b | Prokaryotic recombinant fusion protein corresponding to the first four complement control protein domains of the CD35 molecule | Human |
| CD38 | MMab | SPC32 | IgG1 | Recombinant protein encoding the extracellular domain of human CD38 | Human, Rabbit |
| CD38 | RMab | EP135 | IgG | Synthetic peptide corresponding to residues of human CD38 protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|--|-------------------------|----------------|--|
| | Placenta, Tonsil, Spleen, Diffuse Large B-cell Lymphoma | Cytoplasmic, Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma, Lung Cancer, Sarcoma and Soft Tissue |
| | Placenta, Cervix, Colon, Pancreas, Brain, Kidney, Testis, Pancreatic Carcinoma, Colon Adenocarcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Lung Cancer, Breast Cancer, Colon & Gastrointestinal Cancer, Prostate Cancer, Liver Cancer, Gall Bladder & Pancreatic Cancer, Infectious Diseases |
| | Testis, Colon, Kidney, Stomach, Brain | Membranous, Cytoplasmic | 1:50-1:200 | Infectious Disease, Kidney Cancer & Urothelial, Breast Cancer, Colon & G.I. Tract, Ovarian Cancer, Head & Neck Cancer, Prostate Cancer, Gall Bladder & Pancreatic Cancer, Endometrial Cancer, Liver Cancer |
| | Tonsil, Lymph Node, Hodgkin's Lymphoma | Cytoplasmic, Membranous | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Lung Cancer |
| | Tonsil, Lymph Node, Hodgkin's Lymphoma | Cytoplasmic, Membranous | 1:50-1:200 | Hodgkin's and NHD Lymphoma, Lung Cancer |
| | Placenta, Liver, Breast, Spleen, Thymus, Lung | Cytoplasmic, Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma, Colon & GI Cancer |
| | Placenta, Tonsil, Lymph Node, Inflamed Tissue, H Hylori Infected Tissues | Cytoplasmic, Membranous | 1:25-1:100 | Leukemia & Histiocytic, Sarcoma & Soft Tissue, Melanoma & Skin Cancer |
| | Tonsil, Lymph Node, Spleen, Colon | Membranous | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Leukemia & Histiocytic, Rejection & Autoimmunity |
| | Skin, Thymus, Lymphoblastic Lymphoma | Cytoplasmic, Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma, Colon & GI Cancer |
| | Tonsil, Lymph Node, Colon, Fallopian Tube, Thymus, Spleen | Membranous | 1:25-1:100 | Lymphoma |
| | Tonsil, Lymph Node | Membranous | 1:250-1:1000 | Hodgkin's and NHD Lymphoma, Leukemia & Histiocytic, Rejection & Autoimmunity |
| | Tonsil, Lymph Node, Spleen | Membranous | 1:50-1:200 | Hodgkin's and NHD Lymphoma, Lymphoma, Sacroma |
| | Tonsil, Lymph Node | Membranous | 1:50-1:200 | Hodgkin's and NHD Lymphoma, Leukemia & Histiocytic |
| | Tonsil, Lymph Node | Membranous | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Leukemia & Histiocytic |
| | Tonsil, Small Bowel, Colon, Spleen, Mastocytosis, Hodgkin's Lymphoma | Cytoplasmic, Membranous | 1:10-1:50 | Leukemia & Histiocytic, Lymphoma, Liver Cancer, Melanoma & Skin Cancer |
| | Tonsil, Small Bowel, Colon, Spleen, Mastocytosis, Hodgkin's Lymphoma | Cytoplasmic, Membranous | 1:10-1:50 | Leukemia & Histiocytic, Lymphoma, Liver Cancer, Melanoma & Skin Cancer |
| | Tonsil, Lymph Node, Liver, Testis, Kidney, Colon, Spleen, Thymus, Lymphoblastic Lymphoma | Membranous | 1:50-1:200 | Hodgkin's and NHD Lymphoma, Lymphoma |
| | Tonsil, Lymph Node | Membranous | 1:100-1:500 | Hodgkin's and NHD Lymphoma, Lymphoma, Leukemia & Histiocytic |
| | Tonsil, Lymph Node, Hodgkin's Lymphoma | Membranous | 1:100-1:500 | Hodgkin's and NHD Lymphoma, Lymphoma, Testicular Cancer, Ovarian Cancer |
| | Tonsil, Placenta, Appendix, Spleen, Kidney | Cytoplasmic, Membranous | 1:50-1:200 | Breast Cancer, Endothelial, Hematopoietic, Rejection & Autoimmunity |
| | Placenta, Myometrium, Lung, Colon, Spleen, Lymph Node, Tonsil, Acute Myeloid Leukemia | Membranous | 1:10-1:50 | Placenta, Myometrium, Lung, Colon, Spleen, Lymph Node, Tonsil, Acute Myeloid Leukemia |
| | Tonsil, Placenta, Appendix | Cytoplasmic, Membranous | 1:100-1:500 | Endothelial, Hematopoietic, Leukemia & Histiocytic, Sarcoma & Soft Tissue, Liver Cancer, GIST, Colon & GI Cancer, Undifferentiated Tumor |
| | Tonsil, Placenta, Appendix | Cytoplasmic, Membranous | 1:50-1:200 | Endothelial, Hematopoietic, Leukemia & Histiocytic, Sarcoma & Soft Tissue, Liver Cancer, GIST, Colon & GI Cancer, Undifferentiated Tumor |
| | Tonsil, Lymph Node | Membranous | 1:50-1:200 | Lymphoma, Sarcoma & Soft Tissue |
| | Tonsil, Lymph Node | Membranous | 1:25-1:100 | Lymphoma, Sarcoma & Soft Tissue |
| | Tonsil, Lymph Node, Spleen, Prostate, Salivary Gland | Membranous | 1:25-1:100 | Tonsil, Lymph Node, Spleen, Prostate, Salivary Gland |
| | Synthetic peptide corresponding to residues of human CD38 protein | Membranous | 1:25-1:100 | Leukemia & Histiocytic, Lymphoma, Rejection & Autoimmunity |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--------------------------|---------|------------------|---------|---|---|
| CD4 | RMab | RBT-CD4 | IgG | Synthetic peptide corresponding to residues in the internal region of the human CD4 protein | Human |
| CD41/ Integrin alpha IIb | RMab | EP178 | IgG | Synthetic peptide corresponding to residues of human CD41/ Integrin alpha IIb protein | Human |
| CD42b | RMab | EP409 | IgG | A synthetic peptide corresponding to residues of human CD42b protein. | Human |
| CD43 | MMab | MT1 | IgG1 | Human lymph node cells | Human |
| CD44 | MMab | BSB-12 | IgG2a | Recombinant human CD44 protein | Human |
| CD45 | MMab | 2B11 & PD7/26 | IgG1/K | PD7/26/16: human peripheral blood lymphocytes maintained in T cell growth factor and 2B11: isolated neoplastic cells from T cell lymphoma | Human |
| CD45R | MMab | MB1 | IgG1 | Isolated from membrane and human lymphocytes | Human |
| CD45RA | MMab | 4KB5 | IgG1/K | Hairy cell leukemia cells | Human |
| CD45RO | MMab | UCHL-1 | IgG2a/K | Interleukin-2-dependent human T lymphocytes | Human, Mouse, Rat, Non-human Primate |
| CD5 | RMab | RBT-CD5 | IgG | Synthetic peptide corresponding to residues from the intercellular region of the human CD5 protein | Human |
| CD5 | RMab | RM314 | IgG | A peptide corresponding to the C-terminus of human CD5. | Human |
| CD56 | MMab | 123C3.D5 | IgG1/K | Membrane preparation of a small cell lung carcinoma | Human |
| CD57 | MMab | BSB-10 (CD57/B8) | IgG1/K | Synthetic peptide corresponding to residues of human CD57 protein | Human, Dog, Cat |
| CD6 | MMab | BSB-54 | IgG1 | Recombinant human CD6 protein | Human |
| CD61 | RMab | EP65 | IgG | A synthetic peptide corresponding to residues of human CD61 protein. | Human |
| CD61 | MMab | 2f2 | IgG1/K | Recombinant protein encoding part of the external domain of human CD61 | Human |
| CD63 | MMab | NKI/C3 | IgG1/K | Smooth plasma membrane fraction of MeWo cells | Human |
| CD68 | MMab | BSB-8 (CD68/G2) | IgG1 | Synthetic peptide corresponding to residues from the internal region of the human CD68 protein | Human, Mouse, Rabbit, Rat |
| CD68 | MMab | KP-1 | IgG1/K | Subcellular fraction of human alveolar macrophages | Human, Hamster, Mouse, Non-Human Primate, Porcine, Rabbit, Rat, Cat, Monkey |
| CD7 | MMab | LP15 | IgG2b | Full length of the human CD7 protein | Human |
| CD7 | RMab | EP132 | IgG | Synthetic peptide corresponding to residues of the human CD7 protein | Human |
| CD71 | MMab | 10F11 | IgGb2 | Prokaryotic recombinant protein corresponding to a region of the N-terminal intracellular domain of the human CD71 molecule | Human |
| CD73/NTSE | RMab | RM431 | IgG | Recombinant human CD73 protein | Human |
| CD74 | MMab | LN2 | IgG1 | SU-DHL-4 lymphoma cells | Human |
| CD75 | MMab | LN-1 | IgG1/k | Nuclei from pokeweed mitogen stimulated peripheral blood lymphocytes | Human |
| CD79a | MMab | JCB117 | IgG1/K | Recombinant protein containing part of the extracellular portion of the CD79a glycoprotein | Human |
| CD8 | RMab | EP334 | IgG | A synthetic peptide corresponding to residues of human CD8 α chain | Human |
| CD8 | MMab | C8/144B | IgG/K | A 13 aminoacid synthetic peptide from the C-terminal cytoplasmic domain of alpha chain of the human CD8 protein | Human |
| CD99 | MMab | BSB-9 (CD99/B5) | IgG1/K | Recombinant human CD99 protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|---|-------------------------|----------------|--|
| | Tonsil, Lymph Node | Membranous | 1:25-1:100 | Melanoma & Skin Cancer, Lymphoma |
| | Spleen, Bone Marrow | Cytoplasmic, Membranous | 1:25-1:100 | Hematopoietic, Melanoma & Skin Cancer |
| | Bone Marrow, Spleen | Cytoplasmic | 1:25-1:100 | Leukemia & Histiocytic |
| | Human lymph node cells | Membranous | 1:100-1:500 | Lymphoma |
| | Urothelium, Tonsil, Kidney, Breast, Liver, Skin, Prostate, Thymus, Spleen, Lymph Node, Esophageal Carcinoma | Membranous | 1:250-1:1000 | Kinley & Urothelial Cancer, Breast Cancer, Endometrial & Genital Cancer, Colon & GI Cancer, Melanoma & Skin Cancer |
| | Tonsil, Lymph Node, Thymus, Spleen | Membranous | 1:250-1:1000 | Hodgkin's and NHD Lymphoma, Leukemia & Histiocytic, Undifferentiated Tumor |
| | Tonsil, Lymph Node | Membranous | 1:25-1:100 | Lymphoma |
| | Tonsil, Lymph Node | Membranous | 1:25-1:100 | Colon & GI Cancer, Breast Cancer, Cervical Cancer |
| | Tonsil, Lymph Node | Membranous | 1:250-1:1000 | Lymphoma, Breast Cancer, Colon & GI Cancer, Kidney & Urothelial Cancer |
| | Tonsil, Lymph Node, Spleen, Thymus | Membranous | 1:25-1:100 | Leukemia & Histiocytic, Lymphoma |
| | Tonsil, Lymph Node, Spleen, Thymus | Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma |
| | Pancreas, Neuroblastoma, Brain, Thyroid, Prostate, Colon, Lung SCC | Membranous | 1:250-1:1000 | Leukemia & Histiocytic, Lymphoma, Lung Cancer, Neural & Neuroendocrine Cancer, Undifferentiated Tumor |
| | Tonsil, Lymph Node, Spleen, Prostate, Breast, Brain | Membranous | 1:100-1:500 | Hodgkin's and NHD Lymphoma, Neural & Neuroendocrine Cancer, Kidney & Urothelial Cancer |
| | Tonsil, Lymph Node, Prostate, Colon, Spleen, Bone Marrow | Cytoplasmic | 1:250-1:1000 | Lymphoma, Prostate Cancer |
| | Brain, Kidney, Testis, Bone Marrow | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic |
| | Bone Marrow | Cytoplasmic | 1:100-1:500 | Hematopoietic, Leukemia & Histiocytic, Kidney & Urothelial Cancer |
| | Skin, Malignant Melanoma | Cytoplasmic, Membranous | 1:250-1:1000 | Melanoma & Skin Cancer, Kidney & Urothelial Cancer |
| | Tonsil, Lymph Node | Cytoplasmic, Membranous | 1:250-1:1000 | Leukemia & Histiocytic, Kidney & Urothelial Cancer, Breast Cancer |
| | Tonsil, Lymph Node | Cytoplasmic, Membranous | 1:250-1:1000 | Leukemia & Histiocytic, Kidney & Urothelial Cancer, Breast Cancer |
| | Tonsil, Lymph Node, Colon, Liver, Spleen, Bone Marrow | Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma |
| | Tonsil, Lymph Node, Colon, Liver, Spleen, Bone Marrow, Lymphoblastic Lymphoma | Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma |
| | Bone Marrow, Placenta, Adrenal, Tonsil, Skin | Cytoplasmic, Membranous | 1:25-1:100 | Leukemia & Histiocytic, Breast Cancer, |
| | Placenta, Adrenal Gland, Liver, Testis, Transitional Cell Carcinoma, Ovarian Serous Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Lung Cancer, Colon & G.I. Cancer, Melanoma and Skin Cancer, Breast Cancer, Gall Bladder & Pancreatic Cancer, Immunotherapy |
| | Tonsil, Lymph Node | Cytoplasmic, Membranous | 1:250-1:1000 | Lymphoma, Leukemia & Histiocytic, Lung Cancer |
| | Liver, Prostate, Kidney, Tonsil | Cytoplasmic, Membranous | 1:250-1:1000 | Lymphoma, Kidney & Urothelial Cancer |
| | Recombinant protein containing part of the extracellular portion of the CD79a glycoprotein | Membranous | 1:250-1:1000 | Lymphoma, Leukemia & Histiocytic, Hodgkin's and NHD Lymphoma |
| | Tonsil, Lymph Node, Liver, Colon | Membranous | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma, Melanoma & Skin Cancer, Immunotherapy |
| | Tonsil, Lymph Node, Liver, Fallopian Tube, Breast, Colon, Prostate | Membranous | 1:250-1:1000 | Leukemia & Histiocytic, Lymphoma, Melanoma & Skin Cancer, Immunotherapy |
| | Pancreas, Thymus, Ependyma, Ewing's Sarcoma & Soft Tissue | Cytoplasmic, Membranous | 1:50-1:200 | Ovarian Cancer, Sarcoma & Soft Tissue, Lymphoma, Leukemia & Hystioctic, Undifferentiated Tumor |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|------------------|---------|----------------|---------|--|--|
| CDK2 | RMab | RBT-CDK2 | IgG | Synthetic peptide corresponding to the N-terminus of the human CDK2 protein | Human |
| CDK4 | RMab | EP180 | IgG | Synthetic peptide corresponding to human CDK4 protein | Human |
| CDX2 | RMab | EP25 | IgG | Synthetic peptide corresponding to residues near the C-terminus of the human CDX-2 protein | Human, Predicted: Rabbit |
| CEA | MMab | BSB-13 (CEA31) | IgG1/K | Recombinant human CEA | Human |
| CEA | RPab | Polyclonal | IgG | Recombinant human CEA | Human, Dog, Cat |
| Chromogranin A | MMab | LK2H10 | IgG1/K | Purified human pheochromocytoma | Human, Dog, Cat, Rat, Rabbit, Porcine |
| Claudin-1 | RPab | Polyclonal | IgG | Synthetic peptide derived from the C-terminus of the human Claudin-1 protein | Human |
| Claudin-5 | RMab | EP224 | IgG | Synthetic peptide corresponding to residues of human Claudin-5 protein | Human |
| Claudin-7 | RMab | EP399 | IgG | Synthetic peptide corresponding to residues of the human Claudin-7 protein | Human |
| c-Met/HGFR | RMab | EP1454Y | IgG | Synthetic peptide corresponding to residues near the C-terminus of human Met protein | Human |
| c-Myc | MMab | 9 E10 | IgG1 | Human c-myc gene product, (amino-acid residues 408-432)-KLH conjugate | Human |
| c-Myc | RMab | EP121 | IgG | Synthetic peptide (the amino acid sequence is considered to be commercially sensitive) within Human c-Myc aa 1-100 (N terminal). | Human, Predicted: Mouse, Rat |
| Collagen Type IV | MMab | CIV22 | IgG1/K | Purified human glomeruli | Human, Rat, Canine, Horse |
| Collagen Type IV | RMab | RBT-COL4 | IgG | Recombinant human Collagen IV protein | Human |
| COX-2 | RMab | RBT-COX2 | IgG | Synthetic peptide corresponding to residues near the C-terminus of human COX-2 | Human, Dog, Cat, Mouse, Rat |
| COX-2 | RMab | EP293 | IgG | A synthetic peptide corresponding to residues of human COX-2 protein | Human, Predicted Mouse |
| CTLA-4/CD152 | MMab | BSB-88 | IgG2a/K | Synthetic peptide corresponding to residues of the C-terminus of the human CTLA-4 protein | Human, Mouse, Rat, Predicted: Canine, Equine |
| CTLA-4/CD152 | RMab | RBT-CTLA-4 | IgG | Synthetic peptide corresponding to residues of the C-terminus of the human CTLA-4 protein. | Human |
| Cyclin B1 | RMab | RBT-B1 | IgG | Synthetic peptide corresponding to residues of the C-terminus of the human cyclin B1 protein | Human |
| Cyclin B1 | RMab | RM281 | IgG | A peptide corresponding to Cyclin B1 | Human |
| Cyclin D1 | RMab | RM241 | IgG | A peptide corresponding to Cyclin D1 | Human, Predicted: Mouse, Rat |
| Cyclin D1 | RMab | RBT-14 | IgG | Synthetic peptide corresponding to residues of the C-terminus of the human cyclin D1 protein | Human, Mouse, Rat |
| Cyclin E1 | RMab | EP126 | IgG | Synthetic peptide corresponding to residues of human Cyclin E1 protein | Human |
| Cytokeratin 10 | RMab | EP97 | IgG | Synthetic peptide corresponding to residues in the C-term of human Cytokeratin 10 (CK10) protein | Human, Predicted: Mouse, Rat |
| Cytokeratin 14 | RMab | RM328 | IgG | A peptide corresponding to C-terminus of human Cytokeratin-14 | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|----------------------------------|----------------|---|
| Testis, Tonsil, Prostate, Placenta, Skin, Colon, Transitional Cell Carcinoma, Lymphoblastic Lymphoma | Nuclear | 1:25-1:100 | Breast Cancer, Liver Cancer, Lung Cancer |
| Cervical and Colon Cancer | Nuclear | 1:25-1:100 | Breast Cancer, Sarcoma & Soft Tissue |
| Colon, Colon Adenocarcinoma | Nuclear | 1:50-1:200 | Colon & GI Cancer, Liver Cancer, Lung Cancer, Ovarian Cancer, Gall Bladder & Pancreatic Cancer, Carcinomas of Unknown Primary Site |
| Colon, Colon Adenocarcinoma | Cytoplasmic | 1:250-1:1000 | Colon & GI Cancer, Liver Cancer, Lung Cancer, Ovarian Cancer, Cervical Cancer, Carcinomas of Unknown Primary Site Gall Bladder & Pancreatic Cancer, Mesothelioma, Lung Cancer |
| Colon, Tonsil, Fallopian Tube, Colon Adenocarcinoma | Cytoplasmic | 1:100-1:500 | Colon & GI Cancer, Liver Cancer, Lung Cancer, Ovarian Cancer, Cervical Cancer, Carcinomas of Unknown Primary Site Gall Bladder & Pancreatic Cancer, Mesothelioma, Lung Cancer |
| Pancreas, Pituitary, Colon, Brain | Cytoplasmic | 1:250-1:1000 | Gall Bladder & Pancreatic Cancer, Lung Cancer, Neural & Neuroendocrine Cancer, Colon & GI Cancer, Carcinomas of Unknown Primary Site |
| Skin, Small Intestine, Colon Carcinoma | Membranous | 1:50-1:200 | Breast Cancer, Colon & GI Cancer, |
| Liver, Vascular Tissue, Placenta, Colon, Kidney, Fallopian Tube | Cytoplasmic, Membranous | 1:25-1:100 | Endothelial, Breast Cancer |
| Breast, Colon, Fallopian Tube, Pancreas, Kidney, Transitional Cell Carcinoma | Nuclear, Cytoplasmic, Membranous | 1:50-1:200 | Colon & G.I. Cancer, Lung Cancer, Prostate Cancer, Ovarian Cancer, Breast Cancer |
| Breast, Tonsil, Cervix, Papillary Thyroid Carcinoma, Colon Carcinoma | Cytoplasmic, Membranous | 1:10-1:50 | Kidney & Urothelial Cancer, Liver Cancer, Breast Cancer, Neural & Neuroendocrine Cancer, Colon & GI Cancer, Lung Cancer |
| Burkitt Lymphoma, Lung Cancer, Prostate Cancer | Nuclear, Cytoplasmic | 1:25-1:100 | Leukemia & Histiocytic, Lymphoma, Prostate Cancer |
| | Nuclear, Cytoplasmic | 1:10-1:50 | Leukemia & Histiocytic, Lymphoma, Prostate Cancer |
| Muscle, Lung | Cytoplasmic | 1:50-1:200 | Sarcoma & Soft Tissue, Breast Cancer, Gall Bladder & Pancreatic Cancer |
| Muscle, Lung, Breast, Placenta, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Cytoplasmic, Membranous | 1:50-1:200 | Sarcoma & Soft Tissue, Breast Cancer, Gall Bladder & Pancreatic Cancer |
| Colon, Testis, Kidney, Placenta, Liver, Fallopian Tube, Pancreas, Tonsil, Thymus, Breast, Adenocarcinoma of Colon, Bladder TCC | Cytoplasmic | 1:25-1:100 | Colon & GI Cancer |
| Colon, Stomach, Pancreas, Breast, Lung, Adenocarcinoma of Colon | Cytoplasmic | 1:50-1:200 | Colon & GI Cancer |
| Tonsil, Lymph Node, Colon, Thymus | Membranous | 1:25-1:100 | Rejection & Autoimmunity, Lymphoma; Immunotherapy |
| Tonsil, Lymph Node, Colon, Thymus | Membranous | 1:25-1:100 | Rejection & Autoimmunity, Lymphoma; Immunotherapy |
| Testis, Breast, Colon, Cervix, HSIL Cervix, Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Cervical Cancer, Breast Cancer, Lung Cancer |
| Testis, Breast, Colon, Cervix, HSIL Cervix, Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Cervical Cancer, Breast Cancer, Lung Cancer |
| Tonsil, Placenta, Brain, Pituitary, Adrenal, Cervix, Breast, Mantle Cell Lymphoma, Breast Carcinoma | Nuclear | 1:25-1:100 | Lymphoma, Breast Cancer, Lung Cancer |
| Tonsil, Placenta, Brain, Pituitary, Adrenal, Cervix, Breast, Mantle Cell Lymphoma, Breast Carcinoma | Nuclear | 1:25-1:100 | Lymphoma, Breast Cancer, Lung Cancer |
| Placenta, Bladder, Colon, Breast Cancer, Ovarian Carcinoma | Nuclear | 1:50-1:200 | Cervical Cancer, Breast Cancer, Colon & GI Cancer, Endometrial & Genital Cancer, Hodgkin's and NHD Lymphoma |
| Squamous Cell Carcinoma | Cytoplasmic | 1:25-1:100 | Melanoma & Skin Cancer |
| Squamous Mucosa, Prostate, Breast, Tonsil, Salivary Gland, Skin, Cervix Carcinoma | Cytoplasmic | 1:50-1:200 | Breast Cancer, Melanoma & Skin Cancer, Kidney & Urothelial Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--------------------------------|---------|---------------|---------|--|---|
| Cytokeratin 14 | MMab | LL002 | IgG3 | Synthetic peptide of 15 amino acid residues from the C-terminus of the human cytokeratin 14 | Human, Rats |
| Cytokeratin 16 | RMab | EP27 | IgG | A synthetic peptide corresponding to residues of human Cytokeratin 16 protein | Human |
| Cytokeratin 17 | MMab | BSB-33 | IgG2a/K | Synthetic peptide corresponding to human Cytokeratin 17 protein | Human, Rat, Goat and Pig |
| Cytokeratin 17 | MMab | EP98 | IgG | Synthetic peptide corresponding to residues on the C-terminus of human Cytokeratin 17 (CK17) protein | Human, Rat, Mouse |
| Cytokeratin 18 | RMab | EP30 | IgG | Synthetic peptide corresponding to residues of human Cytokeratin 18 protein | Human |
| Cytokeratin 19 | MMab | BSB-34 | IgG1/K | Synthetic peptide against N-terminus of human cytokeratin 19 protein | Human |
| Cytokeratin 19 | RMab | RM364 | IgG | A peptide corresponding to the C-terminus of human CK-19 (Cytokeratin-19) | Human, Predicted: Mouse |
| Cytokeratin 20 | MMab | Ks20.8 | IgG2a/K | Semi-purified human cytokeratin preparation | Human, Dog, Cat, Mouse |
| Cytokeratin 20 | RMab | EP23 | IgG | Synthetic peptide corresponding to residues near the C-terminus of human CK20 | Human, Predicted: Rat, Goat, Pig, Marmoset |
| Cytokeratin 35BH11 | MMab | 35BetaH11 | IgM | Triton/High salt-insoluble material from the human hepatocellular carcinoma cell line Hep3B | Human |
| Cytokeratin 4 | RMab | EP4 | IgG | Synthetic peptide corresponding to residues on the C-terminus of human CK4 protein | Human |
| Cytokeratin 5 | RMab | EP24 | IgG | Synthetic peptide corresponding to residues near the C-term of human CK-5 protein | Human, Predicted: Mouse |
| Cytokeratin 5 & 6 | RMab | EP24 & EP67 | IgG | "Cytokeratin 5: A synthetic peptide corresponding to residues near the C-terminus of human CK-5 protein. Cytokeratin 6: A synthetic peptide corresponding to residues on the C-terminus of human CK-6 protein." | Human |
| Cytokeratin 5 & 6 | MMab | D5 & 16B4 | IgG1 | Purified human cytokeratin 5 | Human |
| Cytokeratin 6 | RMab | EP67 | IgG | Synthetic peptide corresponding to residues on the C-terminus of human Cytokeratin 6 protein | Human |
| Cytokeratin 7 | RMab | RM284 | IgG | A peptide corresponding to the C-terminus of human Cytokeratin-7 | Human |
| Cytokeratin 7 | MMab | OV-TL 12/30 | IgG1/K | OTN II ovarian carcinoma cell line | Human, Dog, Cat |
| Cytokeratin 8 | RMab | EP17 | IgG | Synthetic peptide corresponding on the C-terminus of human Cytokeratin 8 protein | Human, Predicted: Mouse |
| Cytokeratin 8 & 18 | MMab | B22.1 & B23.1 | IgG1 | Cytoskeleton preparations from HeLa cells and PMC-42 human breast carcinoma cells | Human |
| Cytokeratin Cocktail AE1 & AE3 | MMab | AE1 & AE3 | IgG1 | Purified human epidermal keratin | Human, Dog, Cat, Mouse, Rat, Monkey, Rabbit, Chicken, Horse |
| Cytokeratin HMW 34BE12 | MMab | 34BetaE12 | IgG1/K | Purified keratin from human stratum corneum | Human, Dog, Cat, Monkey, Rabbit, Cattle, Horse |
| Cytokeratin HMW AE3 | MMab | AE3 | IgG1 | Purified human epidermal keratin | Human |
| Cytokeratin LMW | MMab | CAM5.2 | IgG2a/K | Purified human keratin 8 | Human |
| Cytokeratin LMW AE1 | MMab | AE1 | IgG1 | Purified human epidermal keratin | Human, Mouse, Rat, Monkey, Rabbit, Chicken, Cattle, Shrew, Fish |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|--------------|----------------|---|
| Squamous Mucosa, Prostate, Breast, Tonsil, Salivary Gland, Skin, Cervix Carcinoma, Squamous Carcinoma | Cytoplasmic | 1:50-1:200 | Breast Cancer, Melanoma & Skin Cancer, Kidney & Urothelial Cancer |
| Skin, Prostate, Breast, Servix, Salivary Gland, SCC | Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity, Melanoma & Skin Cancer, Cervical Cancer, Head & Neck Cancer |
| Skin, Testis, Breast, Cervix, Cervical Carcinoma, Bladder TCC | Cytoplasmic | 1:25-1:100 | Ovarian Cancer, Cervical Cancer, Lung Cancer, Breast Cancer, Gall Bladder & Pancreatic Cancer, Carcinoma of Unknown Primary Site |
| Skin, Testis, Breast, Cervix, Colon, Salivary Gland, Cervical Carcinoma, Bladder TCC, Transitional Cell Carcinoma | Cytoplasmic | 1:25-1:100 | Ovarian Cancer, Cervical Cancer, Lung Cancer, Breast Cancer, Gall Bladder & Pancreatic Cancer, Carcinoma of Unknown Primary Site |
| Liver, Kidney, Breast, GI, Prostate, Fallopian Tube | Cytoplasmic | 1:25-1:100 | Breast Cancer, Colon & GI Cancer, Neural & Neuroendocrine Cancer |
| Colon, Thyroid Carcinoma | Cytoplasmic | 1:100-1:500 | Colon & GI Cancer, Liver Cancer, Thyroid & Parathyroid Cancer, Breast Cancer |
| Colon, Bladder, Thyroid Carcinoma, Colon Carcinoma | Cytoplasmic | 1:50-1:200 | Colon & GI Cancer, Liver Cancer, Thyroid & Parathyroid Cancer, Breast Cancer |
| Colon Mucosa, Bladder, Colon Carcinoma | Cytoplasmic | 1:100-1:500 | Carcinomas of Unknown Primary Site, Colon & GI Cancer, Kidney & Urothelial Cancer, Lung Cancer, Ovarian Cancer |
| Colon Mucosa, Bladder, Colon Carcinoma | Cytoplasmic | 1:50-1:200 | Carcinomas of Unknown Primary Site, Colon & GI Cancer, Kidney & Urothelial Cancer, Lung Cancer, Ovarian Cancer |
| Prostate, Colon | Cytoplasmic | 1:50-1:200 | Breast Cancer, Ovarian Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer |
| Cornea, Anus, Larynx, Pharynx, Tongue, Prostate, Tonsil, Cervix, Squamous Epithelium of Esophagus, Cervical Squamous Carcinoma | Cytoplasmic | 1:25-1:100 | Head & Neck Tumor, Colon & GI Cancer |
| Prostate, Breast, Placenta, Skin, Mesothelioma | Cytoplasmic | 1:25-1:100 | Mesothelioma, Breast Cancer, Lung Cancer, Prostate Cancer, Carcinoma of Unknown Primary Site |
| Prostate, Mesothelioma | Cytoplasmic | 1:25-1:100 | Mesothelioma, Lung Cancer, Prostate Cancer, Breast Cancer, Melanoma & Skin Cancer |
| Prostate, Mesothelioma | Cytoplasmic | 1:25-1:100 | Mesothelioma, Lung Cancer, Prostate Cancer, Breast Cancer, Melanoma & Skin Cancer |
| Oral Mucosa, Esophagus, Skin, Glandular Epithelia | Cytoplasmic | 1:25-1:100 | Rejection & Autoimmunity, Lung Cancer |
| Salivary Gland, Placenta, Breast, Thyroid, Cervix, Lung Adenocarcinoma | Cytoplasmic | 1:50-1:200 | Carcinomas of Unknown Primary Site, Kidney & Urothelial Cancer, Colon & GI Cancer, Lung Cancer, Breast Cancer, Ovarian Cancer, Mesothelioma, Head & Neck Cancer, Melanoma & Skin Cancer |
| Salivary Gland, Placenta, Breast, Thyroid, Cervix, Pancreas, Fallopian Tube, Transitional Cell Carcinoma, Lung Adenocarcinoma | Cytoplasmic | 1:100-1:500 | Carcinomas of Unknown Primary Site, Kidney & Urothelial Cancer, Colon & GI Cancer, Lung Cancer, Breast Cancer, Ovarian Cancer, Mesothelioma, Head & Neck Cancer, Melanoma & Skin Cancer |
| Colon, Prostate, Kidney, Liver, Colon Carcinoma | Cytoplasmic | 1:25-1:100 | Breast Cancer |
| Breast, Ovary, GI, Prostate, Pancreas, Salivary Gland | Cytoplasmic | 1:250-1:1000 | Colon & GI Cancer, Kidney & Urothelial Cancer, Liver Cancer |
| Prostate, Colon, Skin, Stomach | Cytoplasmic | 1:100-1:500 | Carcinomas of Unknown Primary Site, Undifferentiated Tumor, Breast Cancer, Melanoma & Skin Cancer, Lung Cancer, Germ Cell Tumor, Sarcoma & Soft Tissue, Testicular Cancer |
| Prostate, Cervix | Cytoplasmic | 1:50-1:200 | Prostate Cancer, Liver Cancer, Melanoma & Skin Cancer |
| Prostate, Salivary Gland, Bladder | Cytoplasmic | 1:100-1:500 | Undifferentiated Tumor, Breast Cancer, Melanoma & Skin Cancer, Lung Cancer |
| Colon, Breast, Ovarian Carcinoma | Cytoplasmic | 1:10-1:50 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Liver Cancer, Lung Cancer, Carcinoma Unknown Origin |
| Prostate, Salivary Gland, Bladder, Breast, Kidney, Pancreas, Cervix | Cytoplasmic | 1:100-1:500 | Undifferentiated Tumor, Breast Cancer, Melanoma & Skin Cancer, Lung Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|-----------------------------------|---------|----------------------|---------|---|--|
| Cytokeratin OSCAR | MMab | OSCAR | IgG2a | Crude cyokeratin extract prepared from RT-4 and MCF-7 cells | Human |
| Cytokeratin, MNF116 | MMab | MNF116 | IgG | BALB/C mice immunized with crude extract of splenic cells from a nude mouse engrafted with MCF-7 cells, a human breast carcinoma cell line. | Human |
| Cytomegalovirus | MMab | 8B1.2, 1G5.2 & 2D4.2 | IgG2a | Purified human CMV antigen | Human |
| CXCL12/SDF-1 | MMab | BSB-165 | IgG1 | Recombinant human CXCL12/SDF-1 alpha protein | Human |
| CXCR4/CD184/Fusin | RMab | EP394 | IgG | Synthetic peptide corresponding to residues of the human CXCR4 protein | Human, Mouse |
| CXCR5/CD185 | RPab | Polyclonal | IgG | Recombinant human CXCR5 protein | Human |
| Desmin | MMab | D33 | IgG1/K | Purified desmin from human muscle | Human, Dog, Cat, Mouse, Rat, Chicken, Sheep, Hamster |
| Desmin | RMab | EP15 | IgG | Synthetic peptide corresponding to C-terminus of human Desmin protein | Human, Predicted: Mouse, Rat, Guinea Pig |
| Desmoglein-3 | RMab | EP306 | IgG | A synthetic peptide corresponding to residues of human Desmoglein 3 protein | Human |
| DOG1 | RMab | RBT-DOG1 | IgG | Synthetic peptide corresponding to the N-terminus of human DOG-1 protein | Human |
| DOG-1 | RMab | EP332 | IgG | A synthetic peptide corresponding to residues of human DOG1 protein | Human |
| E-Cadherin | RMab | EP6 | IgG | Synthetic peptide corresponding to residues in the 5th cadherin domain of human E-Cadherin protein | Human |
| EGFR | MMab | 31G7 | IgG1 | EGFR derived from A-431 cells | Human |
| EGFR Phospho | RMab | EP11 | IgG | Synthetic peptide corresponding to residues surrounding Tyr1068 of human EGF receptor protein | Human, Predicted: Mouse |
| EMA | MMab | E29 | IgG2a/K | Purified human milk fat globule membrane preparation | Human |
| EpCAM/Epithelial Specific Antigen | MMab | MOC-31 | IgG1/K | Neuraminidase treated cells from a variant small cell lung carcinoma cell line (GLS-1) | Human |
| EpCAM/Epithelial Specific Antigen | MMab | Ber-EP4 | IgG1/K | MCF-7 human breast carcinoma cell line | Human |
| Epstein Barr Virus LMP-1 | MMab | CS1-4 | IgG1 | Recombinant fusion protein of bacterial beta-galactosidase and EBV LMP-1 | Human |
| ERCC1 | RMab | EP219 | IgG | Synthetic peptide corresponding to residues of the human ERCC1 protein | Human, Predicted: Mouse |
| ERG | RMab | EP111 | IgG | Synthetic peptide corresponding to residues on the C-terminus in human ERG protein | Human, Predicted: Mouse, Rat |
| Estrogen Receptor | MMab | BSB-1 | IgG1/K | A synthetic peptide against the N-terminus of human estrogen receptor alpha. | Human |
| Estrogen Receptor | RMab | EP1 | IgG | Recombinant protein of human estrogen receptor alpha corresponding to aminoacids 1-300 | Human |
| Estrogen Receptor | RMab | RBT-11 | IgG | Synthetic peptide against the C-terminus of human estrogen receptor alpha. | Human |
| Estrogen Receptor | RMab | RM292 | IgG | A peptide corresponding to residues near the N-terminus of human ER-alpha. | Human |
| Factor H/Complement Factor H | MMab | BSB-164 | IgG1 | Recombinant human complement protein factor H | Human |
| Factor VIII-Related Antigen | RPab | Polyclonal | IgG | Purified human Factor VIII | Human, Dog, Cat |
| Factor XIIIa | RMab | EP292 | IgG | Recombinant protein corresponding to A-subunit of coagulation Factor XIII | Human |
| Fascin | MMab | BSB-36 | IgG2a/K | Recombinant human fascin protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|--|-------------------------|----------------|---|
| | Breast, Liver, GI, Prostate, Colon, Skin, Stomach, Pancreas, Tonsil | Cytoplasmic | 1:50-1:200 | Carcinomas of Unknown Primary Site, Undifferentiated Tumor, Breast Cancer, Melanoma & Skin Cancer, Lung Cancer, Germ Cell Tumor, Sarcoma & Soft Tissue, Testicular Cancer |
| | Breast, Cervix, Skin, Colon, Colorectal, Gastric, Breast and Prostatic Carcinomas | Cytoplasmic | 1:25-1:100 | Carcinomas of Unknown Primary Site, Undifferentiated Tumor, Breast Cancer, Melanoma & Skin Cancer, Lung Cancer, Germ Cell Tumor, Sarcoma & Soft Tissue, Testicular Cancer |
| | CMV Infected Tissue | Nuclear | 1:25-1:100 | Infectious Diseases, Cytopathology |
| | Testis, Breast, Colon, Fallopian Tube, Tonsil, Transitional Cell Carcinoma, T cell Lymphoblastic Lymphoma | Cytoplasmic, Membranous | 1:25-1:100 | Breast Cancer, Lung Cancer, Prostate Cancer, Ovarian Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer |
| | Fallopian Tube, Adrenal Gland, Stomach, Kidney, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma, Tonsil, Testis | Membranous, Cytoplasmic | 1:25-1:100 | Brain Cancer, Breast Cancer, Lung Cancer, Liver Cancer, Kidney & Urothelial Cancer |
| | Placenta, Brain, Lung, Transitional Cell Carcinoma, Testis, Ovarian Serous Carcinoma, Hepatocellular Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Lymphomas, Lung Cancer |
| | Skeletal Muscle, Placenta, Colon, Prostate, Skin, Fallopian Tube | Cytoplasmic | 1:25-1:100 | Sarcoma & Soft Tissue, Melanoma & Skin Cancer, GIST, Undifferentiated Tumor |
| | Skeletal Muscle, Placenta, Colon, Prostate, Skin, Fallopian Tube | Cytoplasmic | 1:25-1:100 | Sarcoma & Soft Tissue, Melanoma & Skin Cancer, GIST, Undifferentiated Tumor |
| | Tonsil, Skin, Cervix, Cervical Carcinoma, Lung Squamous Cell Carcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Lung Cancer, Rejection & Autoimmunity, Head & Neck Cancer, Melanoma & Skin Cancer |
| | Salivary Gland, Breast, GIST | Cytoplasmic, Membranous | 1:100-1:500 | GIST, Head & Neck Cancer, Sarcoma & Soft Tissue, Colon & GI Cancer |
| | A synthetic peptide corresponding to residues of human DOG1 protein | Cytoplasmic, Membranous | 1:25-1:100 | GIST, Head & Neck Cancer, Sarcoma & Soft Tissue, Colon & GI Cancer |
| | Breast, Colon, Cervix, Pancreas, Lung, Ovary, GI Tract Adenocarcinoma, Breast Carcinoma | Membranous | 1:100-1:500 | Breast Cancer, Kidney & Urothelial Cancer, Mesothelioma |
| | Skin, Placenta, Testis, Tonsil, Pancreas, Squamous Cell Carcinoma | Cell Membrane | 1:25-1:100 | Breast Cancer, Colon & GI Cancer, Lung Cancer |
| | Skin, Placenta, Testis, Tonsil, Pancreas, Squamous Cell Carcinoma | Cell Membrane | 1:25-1:100 | Breast Cancer, Colon & GI Cancer, Lung Cancer |
| | Breast, Skin, Colon, Kidney, Cervix | Cytoplasmic, Membranous | 1:250-1:1000 | Breast Cancer, Ovarian Cancer, Lung Cancer |
| | Breast, Colon, Prostate, Kidney, Thyroid, Pancreas, Salivary Gland, Adenocarcinomas | Cytoplasmic | 1:50-1:200 | Mesothelioma, Lung Cancer, Colon & GI Cancer, Liver Cancer, Cytopathology |
| | Colon, Cervix, Salivary Gland, Pancreas, Breast, Thyroid, Liver, Adenocarcinomas | Cytoplasmic | 1:25-1:100 | Mesothelioma, Lung Cancer, Colon & GI Cancer, Cytopathology |
| | EBV Infected Tissue, Hodgkin's Lymphoma | Cytoplasmic | 1:25-1:100 | Infectious Diseases, Hodgkin's and NHD Lymphoma |
| | Tonsil, Testis, Breast, Prostate | Nuclear | 1:50-1:200 | Lung Cancer, Ovarian Cancer |
| | Prostate, Colon, Kidney, Fallopian Tube, Tonsil, Myometrium, Skin, Brain, Breast | Nuclear | 1:100-1:500 | Endothelial, Prostate Cancer, Kidney & Urothelial Cancer |
| | Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Carcinomas of Unknown Primary Site, Endometrial Genital Cancer |
| | Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Carcinomas of Unknown Primary Site, Endometrial Genital Cancer |
| | Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Carcinomas of Unknown Primary Site, Endometrial Genital Cancer |
| | Breast, Myometrium, Cervix, Fallopian Tube, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Carcinomas of Unknown Primary Site, Endometrial Genital Cancer |
| | Testis, Liver, Kidney, Pancreas, Adrenal Gland | Cytoplasmic, Membranous | 1:50-1:200 | Lung Cancer, Ovarian Cancer, Brain Cancer, Colon Cancer, Rejection & Autoimmunity, Infectious Diseases |
| | Skin, Placenta | Cytoplasmic | 1:50-1:200 | Endothelial, Sarcoma & Soft Tissue |
| | Placenta, Colon, Adrenal, Testis, Dermatofibroma | Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic, Melanoma & Skin Cancer |
| | Lymph Node, Tonsil, Liver, Kidney, Adrenal, Hodgkin's Lymphoma | Cytoplasmic | 1:250-1:1000 | Hodgkin's and NHD Lymphoma, Lymphoma |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|---|---------|--------------|---------|---|---|
| FGFR-3 | MMab | BSB-150 | IgG2a | Synthetic peptide corresponding to the N-terminus of the human FGFR-3 protein | Human |
| Fibrinogen | RPab | Polyclonal | IgG | KLH conjugated synthetic peptide corresponding to the N-terminus of human FGA. | Human |
| Fli-1 | MMab | G146-222 | IgG2b | Fli-1 ets Domain Fusion Protein | Human |
| FOXA1/HNF-3A | RPab | Polyclonal | IgG | Synthetic peptide sequence (GVYSRPVLNTS) corresponding to the C-terminus amino acids of FOXA1 | Human, Mouse |
| FOXL2 | RPab | Polyclonal | IgG | Synthetic peptide corresponding to the C-terminus portion of the human FOXL2 protein | Human, Mouse, Pig, Bovine |
| FOXO1 | RMab | EP290 | IgG | A synthetic peptide corresponding to residues of human FOXO1 protein | Human |
| FOXP1 | RMab | EP137 | IgG | Synthetic peptide corresponding to residues of human FOXP1 protein | Human, Predicted: Mouse, Rat |
| FOXP3 | RMab | EP340 | IgG | A synthetic peptide corresponding to residues of human FOXP3 protein | Human |
| FOXP3 | RPab | Polyclonal | IgG | Synthetic peptide conjugated to KLH derived from within residues 50 - 150 of Mouse FOXP3 | Human |
| FSH | MMab | BSB-55 | IgG1/K | Synthetic peptide against the N-terminus of human follicle stimulating hormone receptor | Human |
| FSH | RMab | EP257 | IgG | A synthetic peptide corresponding to residues of human FSH (beta subunit) protein. | Human |
| Fumarate Hydratase/Fumarase | MMab | BSB-151 | IgG1 | Synthetic peptide corresponding to the C-terminus of the human fumarate hydratase protein | Human, Mouse, Rat |
| GAB1/ GRB2-associated-binding protein 1 | MMab | BSB-155 | IgG2a | Synthetic peptide corresponding to the N-terminus of the human Gab 1 protein | Human, Mouse, Rat |
| Galectin-3 | MMab | 9C4 | IgG1 | Recombinant protein corresponding to full length Galectin-3 | Human |
| Gastrin | RPab | Polyclonal | IgG | Synthetic peptide derived from N-terminus of gastrin | Human, Dog, Cat |
| GATA3 | MMab | L50-823 | IgG1/K | Synthetic peptide between trans-activation and DNA-binding domains of GATA-3 | Human, Mouse |
| GATA3 | RMab | EP368 | IgG | Synthetic peptide corresponding to residues of human GATA3 protein | Human |
| GCDFP-15 | MMab | 23A3 | IgG2a | Recombinant protein encoding the excreted domain of human GCDFP15 | Human, Rat |
| GFAP | RMab | RM246 | IgG | A peptide corresponding to the N-terminus of human GFAP | Human, Predicted: Rat |
| GFAP | MMab | G-A-5 | IgG1 | Purified GFAP isolated from porcine spinal cord | Human, Canine, Feline, Rat, Mouse, Rabbit |
| GH | MMab | BSB-99 | IgG1/K | Recombinant human growth hormone protein | Human |
| GH | RPab | Polyclonal | IgG | Recombinant human growth hormone protein | Human |
| GH | RMab | EP267 | IgG | Synthetic peptide corresponding to residues of human growth hormone protein | Human |
| Glucagon | MMab | BSB-111 | IgG1/K | Recombinant human glucagon protein | Human |
| Glucagon | RMab | EP74 | IgG | Synthetic peptide corresponding to residues in human glucagon protein | Human, Dog, Mouse, Rat, Horse |
| GLUT1 | RMab | EP141 | IgG | Synthetic peptide corresponding to residues of human Glut-1 protein | Human, Predicted: Mouse, Rat |
| Glutamine Synthetase | MMab | GS-6 | IgG2a | Human Glutamine Synthetase aa. 1-373 | Human, Mouse, Rat |
| Glycophorin A | MMab | GA-R2 (HIR2) | IgG2b/K | GYP A / CD235a / Glycophorin A antibody was raised against human GYP A | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|--|-------------------------|----------------|---|
| | Skin, Liver, Brain, Testis, Transitional Cell Carcinoma, Ductal Breast Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Kidney & Urothelial Cancer, Cervical Cancer, Breast Cancer, Head and Neck Cancer, Lung, Cancer, Brain Cancer, Ovarian Cancer, Colon and GI Cancer, Leukemia & Histiocytic, Melanoma & Skin Cancer |
| | Breast, Testis, Kidney, Pancreas, Salivary Gland, Skin, Fallopian Tube | Cytoplasmic | 1:25-1:100 | Rejection & Autoimmunity |
| | Adrenal Gland, Fallopian Tube, Placenta, Cervix, Hemangiomas, PNET, AngioSarcoma & Soft Tissue | Nuclear | 1:25-1:100 | Endothelial, Hematopoetic, Sarcoma & Soft Tissue |
| | Prostate, Breast, Tonsil, Colon, Fallopian Tube, Breast Carcinoma, Prostate Carcinoma | Cytoplasmic, Nuclear | 1:25-1:100 | Breast Cancer, Prostate Cancer |
| | Fallopian Tube, Cervix, Pancreas, Placenta, Extra Marginal Zone Lymphoma | Nuclear, Cytoplasmic | 1:50-1:200 | Ovarian Cancer, Germ Cell Tumors, Cervical Cancer |
| | Testis, Thyroid, Tonsil, Lymph Node, Spleen, Lung, Lymphomas | Nuclear | 1:25-1:100 | Hodgkin's & NHD Lymphoma, Cervical Cancer, Prostate Cancer, Sarcoma & Soft Tissue |
| | Tonsil, Lymph Node, Breast | Nuclear, Cytoplasmic | 1:25-1:100 | Lymphoma, Breast Cancer, Ovarian Cancer, Testicular Cancer, Germ Cell Tumor |
| | Tonsil, Lymph Node, Thymus, Colon | Nuclear | 1:25-1:100 | Breast Cancer, Ovarian Cancer, Prostate Cancer, Hodgkin's and NHD Lymphoma, Immunotherapy |
| | Tonsil, Lymph Node, Thymus, Colon | Nuclear | 1:25-1:100 | Breast Cancer, Ovarian Cancer, Prostate Cancer, Hodgkin's and NHD Lymphoma, Immunotherapy |
| | Normal Pituitary | Cytoplasmic | 1:100-1:500 | Pituitary, Neural & Neuroendocrine Cancer |
| | Pituitary | Cytoplasmic | 1:50-1:200 | Pituitary, Neural & Neuroendocrine Cancer |
| | Placenta, Breast, Fallopian Tube, Colon, Kidney, Testis, Colon Adenocarcinoma, HER2 Negative Breast Cancer | Cytoplasmic | 1:25-1:100 | Renal and Urothelial Cancer, Sarcoma & Soft Tissue Cancers |
| | Breast Prostate, Testis, Tonsil, Stomach, Transitional Cell Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Neural & Neuroendocrine Cancer, Breast Cancer, Colon & GI Cancer, Liver Cancer, Ovarian Cancer |
| | Kidney, Testis, Salivary Gland, Breast, Tonsil, Colon, Papillary & Follicular Carcinoma of Thyroid | Cytoplasmic | 1:50-1:200 | Thyroid & Parathyroid Cancer, Undifferentiated Tumor, Head & Neck Cancer, Cytopathology |
| | Stomach | Cytoplasmic | 1:250-1:1000 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer |
| | Breast, Skin, Placenta, Colon, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Carcinomas of Unknown Primary Site |
| | Breast, Skin, Cervix, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Carcinomas of Unknown Primary Site |
| | Breast, Sweat Glands in Skin, Breast Carcinoma | Cytoplasmic | 1:100-1:500 | Breast Cancer, Germ Cell Tumor |
| | Brain | Cytoplasmic | 1:25-1:100 | Head & Neck Cancer, Neural & Neuroendocrine Cancer |
| | Brain | Cytoplasmic | 1:250-1:1000 | Head & Neck Cancer, Neural & Neuroendocrine Cancer |
| | Normal Pituitary | Cytoplasmic | 1:250-1:1000 | Pituitary, Neural & Neuroendocrine Cancer |
| | Normal Pituitary | Cytoplasmic | 1:250-1:1000 | Pituitary, Neural & Neuroendocrine Cancer |
| | Normal Pituitary | Cytoplasmic | 1:250-1:1000 | Pituitary, Neural & Neuroendocrine Cancer |
| | Pancreas, Colon | Cytoplasmic | 1:250-1:1000 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer |
| | Pancreas, Colon | Cytoplasmic | 1:50-1:200 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer |
| | Placenta, Colon, Prostate, Skin, Kidney, Brain, Tonsil, Breast, Mesothelioma, Colon Carcinoma | Membranous | 1:50-1:200 | Lung Cancer, Mesothelioma, Endotelial, Cytopathology |
| | Liver, Tonsil, Testis, Prostate, Hepatocellular Carcinoma, Bladder TCC | Cytoplasmic | 1:50-1:200 | Liver Cancer |
| | Bone Marrow, Placenta, Tonsil, Liver, Prostate, Adrenal, Spleen, Colon, Pancreas, Fallopian Tube | Membranous | 1:100-1:500 | Hematopoetic, Leukemia & Histiocytic |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|---------------------------------------|---------|----------------|--------------|---|------------------------------|
| Glypican-3 | MMab | 1G12 | IgG1 | Fragment corresponding to the C-terminal last 70 amino acids of the human Glypican-3 protein | Human, Dog |
| Granzyme B | RMab | EP230 | IgG | Synthetic peptide corresponding to residues of human Granzyme B protein | Human |
| Granzyme B | RPab | Polyclonal | IgG | Synthetic peptide corresponding to residues E(20) I I G G H E A K P H S R P Y M A Y L(38) of human Granzyme B | Human |
| hCG | MMab | BSB-38 | IgG1/K | Recombinant human hCG beta protein | Human |
| HE4 | RMab | EP370 | IgG | Synthetic peptide corresponding to residues of human HE4 protein | Human |
| HEG1 | MMab | SKM9-2 | IgG1 | 893-SKSPSLVSLPT-903 of the human HEG1 protein | Human |
| Helicobacter pylori | RPab | Polyclonal | IgG | Synthetic peptide corresponding to a region within amino acids 341 and 398 of Helicobacter pylori urease B. | Human |
| Helicobacter Pylori | MMab | BSB-37 | IgG1/K | Purified Helicobacter pylori bacteria microdissected from FFPE infected cell lines | Human |
| Helicobacter Pylori | RMab | EP279 | IgG | Helicobacter pylori | Human |
| Hemoglobin A | RMab | EP124 | IgG | Synthetic peptide corresponding to residues of human Hemoglobin alpha chain protein | Human, Predicted: Mouse, Rat |
| Hepatitis B Virus Core Antigen | RPab | Polyclonal | IgG | Purified Hepatitis B virus | Human |
| Hepatitis B Virus Surface Antigen | MMab | A10F1 | IgG2b/K | Purified HbsAg, subtype ayw3. | Human |
| Hepatocyte Specific Antigen/ Hep-Par1 | MMab | OCH1E5 | IgG1/K | Formalin-fixed, failed human allograft liver that was mechanically disrupted | Human, Dog, Cat |
| HER2 neu | RMab | EP3 | IgG | Synthetic peptide corresponding to residues near the C-terminus of human HER2 protein. | Human |
| HER-2 neu | RMab | RBT-HER2 | IgG | Recombinant protein encoding extracellular domain of human HER-2 protein | Human |
| HER-2 neu | MMab | BSB-3 (HER-24) | IgG1 | Recombinant protein encoding extracellular domain of human HER-2 protein | Human |
| HER2 neu Phospho | RMab | EP123 | IgG | Synthetic peptide corresponding to residues surrounding tyrosine 877 of the human HER-2 protein | Human |
| HER-3/c-erbB-3 | RMab | RBT-HER3 | IgG | Synthetic peptide corresponding to human c-erb-3 (HER-3) protein | Human |
| Herpes Simplex Virus I | MMab | 10A3 | IgG1 | ICP8 purified from U-35-VERO cells | Human |
| Herpes Simplex Virus I | RPab | Polyclonal | IgG | Recombinant protein directed against the major glycoproteins present in the viral envelope of HSV I. | Human |
| Herpes Simplex Virus I & II | MMab | 10A3 & BSB-116 | IgG1 & IgG2a | ICP 8 purified from U-35 vero cells and a recombinant protein directed against the major glycoproteins present in the viral envelope of HSV I & II. | Human |
| Herpes Simplex Virus I & II | RPab | Polyclonal | IgG | Recombinant protein directed against the major glycoproteins present in the viral envelope of HSV I & II | Human |
| Herpes Simplex Virus II | RPab | Polyclonal | IgG | Recombinant protein directed against the major glycoproteins present in the viral envelope of HSV II | Human |
| Herpes Simplex Virus II | MMab | BSB-116 | IgG2a | Recombinant protein directed against the major glycoproteins present in the viral envelope of HSV II | Human |
| hGAL/GCET2 | RMab | EP316 | IgG | A synthetic peptide corresponding to residues of human hGAL (GCET2) protein | Human |
| HHV-8 | MMab | 13B10 | IgG1 | Purified HHV-8 virus | Human |
| HHV-8 | RMab | RBT-HHV8 | IgG | Recombinant protein corresponding to the latent nuclear antigen 1 molecule of HHV8 | Human |
| HIF-1alpha | RMab | EP118 | IgG | Synthetic peptide corresponding to residues near the C-terminus of human HIF-1 alpha protein | Human, Predicted: Mouse, Rat |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|---|-------------------------|----------------|---|
| Melanoma, Hepatoblastoma, Hepatocellular Carcinoma | Cytoplasmic | 1:100-1:500 | Liver Cancer, Ovarian Cancer, Testicular Cancer, Melanoma & Skin Cancer, Germ Cell Tumors |
| Liver, Testis, Cervix, Tonsil, Lymph Node, Spleen | Cytoplasmic (Granular) | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Lymphoma, Leukemia & Histiocytic, Melanoma & Skin Cancer, Rejection & Autoimmunity, Immunotherapy |
| Liver, Testis, Cervix, Tonsil, Lymph Node, Spleen | Cytoplasmic (Granular) | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Lymphoma, Leukemia & Histiocytic, Melanoma & Skin Cancer, Rejection & Autoimmunity, Immunotherapy |
| Placenta, Brain, Pituitary, Lung Adenocarcinoma, Choriocarcinoma | Cytoplasmic | 1:250-1:1000 | Ovarian Cancer, Testicular Cancer, Germ Cell Tumor, Undifferentiated Tumor |
| Thyroid, Salivary Gland, Breast, Ovarian Carcinoma | Cytoplasmic, Nuclear | 1:50-1:200 | Ovarian Cancer, Gall Bladder & Pancreatic Cancer, Lung Cancer |
| Testis, Prostate, Mesothelioma, Papillary Thyroid Carcinoma, Ovarian Serous Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Mesothelioma, Liver Cancer, Thyroid and Parathyroid Cancer, Ovarian Cancer |
| Helicobacter pylori Infected Stomach Mucosa | Cell Wall | 1:100-1:500 | Infectious Diseases, Colon & GI Cancer |
| H.Pylori Infected Stomach Mucosa | Cell Wall | 1:25-1:100 | Infectious Diseases, Colon & GI Cancer |
| H.Pylori Infected Stomach Mucosa | Cell Wall | 1:50-1:200 | Infectious Diseases, Colon & GI Cancer |
| Bone Marrow, Placenta, Spleen | Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic |
| Hepatitis B Infected Liver | Nuclear | 1:25-1:100 | Infectious Diseases, Liver Cancer |
| Hepatitis B Infected Liver | Cytoplasmic | 1:50-1:200 | Infectious Diseases, Liver Cancer |
| Liver, Liver Carcinoma | Cytoplasmic | 1:100-1:500 | Liver Cancer, Carcinoma of Unknown Primary Site |
| Breast Carcinoma | Membranous | 1:50-1:200 | Breast Cancer, Colon & GI Cancer |
| Breast Carcinoma | Membranous | 1:100-1:500 | Breast Cancer, Colon & GI Cancer |
| Breast Carcinoma | Membranous | 1:100-1:500 | Breast Cancer, Colon & GI Cancer |
| Breast Carcinoma | Membranous | 1:100-1:500 | Breast Cancer, Colon & GI Cancer |
| Colon, Cervix, Bladder TCC | Cytoplasmic | 1:25-1:100 | Breast Cancer |
| HSV I Infected Tissues | Nuclear, Cytoplasmic | 1:50-1:200 | Infectious Diseases |
| HSV I Infected Tissue | Nuclear, Cytoplasmic | 1:25-1:100 | Infectious Diseases |
| HSV I & II Infected Tissues | Cytoplasmic, Nuclear | 1:50-1:200 | Infectious Diseases |
| HSV I & II Infected Tissues | Nuclear, Cytoplasmic | 1:25-1:100 | Infectious Diseases |
| HSV II Infected Tissues | Nuclear, Cytoplasmic | 1:25-1:100 | Infectious Diseases |
| HSV II Infected Tissues | Cytoplasmic, Nuclear | 1:50-1:200 | Infectious Diseases |
| Tonsil, Lymph Node, Thymus, Fallopian Tube, Germinal Center B-cell Type Diffuse Large B-cell Lymphoma | Cytoplasmic | 1:25-1:100 | Lymphoma, Hodgkin's & NH Lymphoma |
| Kaposi's Sarcoma | Nuclear | 1:25-1:100 | Sarcoma & Soft Tissue, Lymphoma, Infectious Diseases |
| Kaposi's Sarcoma | Nuclear | 1:50-1:200 | Sarcoma & Soft Tissue, Lymphoma, Infectious Diseases |
| Tonsil, Cervix, Liver, Testis, Kidney, Breast, Thymus, Spleen, Colon, Cervical Cancer | Cytoplasmic, Nuclear | 1:250-1:1000 | Cervical Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--------------------|---------|---------------|---------|---|---------------------------------|
| Histone H3 Phospho | RPab | Polyclonal | IgG | Synthetic peptide corresponding to Human Histone H3 aa 1-100 (phospho S28) conjugated to Keyhole Limpet Haemocyanin (KLH) | Human |
| Histone H3 Phospho | RMab | EP233 | IgG | Synthetic phospho-peptide corresponding to residues surrounding Thr3 of Human Histone H3 | Human |
| HLA-DR alpha chain | RMab | EP96 | IgG | Synthetic peptide corresponding to residues in human HLA-DRA protein | Human |
| HMGA2 | RMab | EP398 | IgG | A synthetic peptide corresponding to residues of the human HMGA2 protein | Human |
| hPL | RMab | EP241 | IgG | Synthetic peptide corresponding to residues of human placental lactogen protein | Human |
| hPL | RPab | Polyclonal | IgG | Purified human placental lactogen | Human |
| HPV | MMab | BSB-66 (SB24) | IgG1/K | Recombinant major capsid protein of human papillomavirus type 16. | Human |
| HPV16 | MMab | CAMVIR-1 | IgG2a | Human papilloma virus type 16, major capsid protein L1 | Human |
| HSP-27 | MMab | G3.1 | IgG1 | Partially purified human HSP27 | Human, Rat, Dog, Mouse, Primate |
| HSP70 | RMab | RM432 | IgG | Recombinant human HSP70 protein | Human |
| ICOS/CD278 | RMab | RM417 | IgG | Synthetic peptide corresponding to residues near the C-terminus of the human ICOS protein | Human |
| IDH1 R132H | MMab | IHC132 | IgG1 | Synthetic peptide corresponding to the human IDH1 surrounding residual 132. | Human |
| IDH1 R132H | RMab | RBT-IDH1 | IgG | Synthetic peptide corresponding to IDH1 R132H mutant | Human |
| IgA | MMab | BSB-39 | IgG1/K | Purified human IgA secretory Component protein | Human |
| IgA | RPab | Polyclonal | IgG | Purified human IgA secretory Component protein. | Human |
| IgD | RPab | Polyclonal | IgG | IgD isolated from a pool of normal human plasma | Human |
| IgD | RMab | EP173 | IgG | Synthetic peptide corresponding to residues of human IgD protein. | Human |
| IgE | RPab | Polyclonal | IgG | Purified IgE ϵ -Heavy Chain isolated from a pool of normal human sera. | Human |
| IFN-Alpha | MMab | BSB-158 | IgG2b | Recombinant human IFN- α 1/13 protein | Human |
| IFN-Gamma | MMab | BSB-161 | IgG2a | Recombinant human IFN- γ protein | Human |
| IgG | MMab | BSB-40 | IgG2a/K | Purified human IgG gamma chain | Human, Canine, Feline |
| IgG | RPab | Polyclonal | IgG | Purified IgG gamma chain. | Human |
| IgG4 | RMab | EP138 | IgG | Synthetic peptide corresponding to residues in the hinge region of Human IgG4 | Human |
| IgG4 | MMab | BSB-96 | IgG1 | Purified human IgG4 | Human |
| IgM | RPab | Polyclonal | IgG | Purified human IgM heavy chain. | Human |
| IgM | MMab | BSB-41 | IgG1/K | Purified human IgM heavy chain | Human |
| IL-1a | MMab | BSB-138 | IgG2b | Recombinant protein corresponding to the IL-1 α of human origin | Human |
| IL-1b | MMab | BSB-139 | IgG2b | Recombinant protein corresponding to the IL-1 β of human origin | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|--|-------------------------|----------------|--|
| | Tonsil, Lymph Node, Thymus, Spleen | Nuclear | 1:25-1:100 | Lymphoma, Melanoma & Skin Cancer, Neural & Neuroendocrine Cancer, Sarcoma & Soft Tissue |
| | Tonsil, Lymph Node, Thymus, Spleen | Nuclear | 1:100-1:500 | Lymphoma, Melanoma & Skin Cancer, Neural & Neuroendocrine Cancer, Sarcoma & Soft Tissue |
| | Tonsil, Spleen, Liver, Kidney, Adrenal, Colon, Lymph Node, Thymus, Lymphoblastic Lymphoma | Cytoplasmic | 1:50-1:200 | Breast Cancer, Liver Cancer, Lung Cancer, Ovarian Cancer |
| | Cervix, Lung Squamous Cell Carcinoma, Papillary Thyroid Carcinoma | Nuclear | 1:50-1:200 | Thyroid and Parathyroid Cancer, Pituitary Cancer, Breast Cancer, Lung Cancer, Colon and GI Cancer, Ovarian Cancer, Liver Cancer, Gallbladder and Pancreatic Cancer, Sarcoma and Soft Tissue Cancer |
| | Placenta | Cytoplasmic | 1:100-1:500 | Endometrial & Genital Cancer, Testicular Cancer, Ovarian Cancer |
| | Placenta | Cytoplasmic | 1:100-1:500 | Endometrial & Genital Cancer, Testicular Cancer, Ovarian Cancer |
| | HPV Infected Tissue | Nuclear | 1:100-1:500 | Cervical Cancer, Head & Neck Cancer, Infectious Diseases |
| | HPV16 Infected Tissues | Nuclear | 1:250-1:1000 | Cervical Cancer, Head & Neck Cancer, Infectious Diseases |
| | Tonsil, Cervix, Prostate, Breast Carcinoma, Cervical Carcinoma | Cytoplasmic | 1:50-1:200 | Cervical Cancer, Breast Cancer, Prostate Cancer |
| | Breast, Fallopian Tube, Skin, Prostate, Testis, Transitional Cell Carcinoma | Nuclear, Cytoplasmic | 1:50-1:200 | Cervical Cancer, Melanoma and Skin Cancer, Kidney and Urothelial Cancers |
| | Colon, Tonsil, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Cytoplasmic, Membranous | 1:50-1:200 | Immunotherapy, Leukemia and Histiocytic Cancer, Lymphoma, Breast Cancer, Melanoma and Skin Cancer |
| | Glioma, Glioblastoma, Astrocytoma with IDH1 R132H Mutation | Cytoplasmic | 1:25-100 | Neural & Neuroendocrine Cancer, Leukemia and Histiocytic Cancer, Liver Cancer, Melanoma and Skin Cancer, Colon and GI Cancers, Prostate Cancer |
| | Glioma, Glioblastoma, Astrocytoma with IDH1 R132H Mutation | Cytoplasmic | 1:25-100 | Neural & Neuroendocrine Cancer, Leukemia and Histiocytic Cancer, Liver Cancer, Melanoma and Skin Cancer, Colon and GI Cancers, Prostate Cancer |
| | Tonsil, Spleen, Lymph Node, Kidney, Colon | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Spleen, Lymph Node, Kidney, Colon | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Lymph Node, Spleen | Cytoplasmic | 1:100-1:500 | Lymphomas, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Lymph Node, Spleen | Cytoplasmic | 1:50-1:200 | Lymphomas, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Thymus, Colon | N/A | 1:250-1:1000 | Rejection & Autoimmunity |
| | Placenta, Fallopian Tube, Stomach, Prostate, Testis, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma, | Cytoplasmic, Membranous | 1:25-1:100 | Immunotherapy, Infectious Disease, Leukemia and Histiocytic Cancer, Melanoma and Skin Cancer, Kidney and Urothelial Cancer |
| | Stomach, Fallopian Tube, Colon, Lung Adenocarcinoma, Ductal Breast Carcinoma, Pancreatic Adenocarcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Immunotherapy, Infectious Disease, Kidney and Urothelial Cancer, Melanoma and Skin Cancer, Ovarian Cancer |
| | Tonsil, Lymph Node, Kidney, Spleen | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Lymph Node, Kidney, Spleen | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Spleen, Colon | Cytoplasmic | 1:50-1:200 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Thyroid & Parathyroid Cancer |
| | Tonsil, Spleen, Colon | Cytoplasmic | 1:50-1:200 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Thyroid & Parathyroid Cancer |
| | Tonsil, Lymph Node, Spleen, Kidney, Colon | Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Tonsil, Lymph Node, Spleen, Kidney, Colon | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic, Hodgkin's & NH Lymphoma, Rejection & Autoimmunity |
| | Kidney, Colon, Adrenal gland, Testis, Lung | Membranous, Cytoplasmic | 1:10-1:50 | Rejection & Autoimmunity, Head and Neck Cancer, Gall Bladder & Pancreatic Cancer, Sarcoma & Soft Tissue, Breast Cancer, Colon & G.I. Tract |
| | Colon, Pancreas, Liver, Stomach, Brain, Testis, Lung, Transitional Cell Carcinoma, Tonsil | Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity, Breast Cancer, Colon & G.I. Tract, Lung Cancer, Liver Cancer, Head and Neck Cancers, Melanoma & Skin Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--------------------|---------|-----------------|---------|---|--------------------------------|
| IL-6 | MMab | BSB-140 | IgG2b | Recombinant protein corresponding to the IL-6 of human origin | Human |
| IMP-3/1GF2BP3 | RMab | EP286 | IgG | A synthetic peptide corresponding to the human IMP-3 protein | Human, Predicted Mouse and Rat |
| Inhibin alpha | RMab | EP378 | IgG | Synthetic peptide corresponding to residues of human Inhibin alpha chain protein | Human |
| Inhibin alpha | MMab | R1 | IgG2a | Synthetic peptide: MVLHLLFLITPGGHSCQGLELARELVLAK, corresponding to amino acids 1-32 of human inhibin alpha | Human, Dog, Cat |
| INI-1 | MMab | 25 | IgG2a | Mouse BAF47 aa.257-359 | Human |
| INI-1 | RMab | RBT-INI1 | IgG | Synthetic peptide corresponding to residues of mouse BAF47 protein | Human |
| iNOS | RMab | RBT-iNOS | IgG | Synthetic peptide corresponding to the N-terminus of the human iNOS protein | Human |
| INSM1 | MMab | BSB-123 | IgG/K | Synthetic peptide corresponding to residues from the N-terminal domain of the human INSM1 protein. | Human |
| INSM1 | RMab | RBT-INSM1 | IgG | "Synthetic peptide corresponding to residues from the N-terminal domain of the human INSM1 protein." | Human |
| Insulin | RMab | EP125 | IgG | Synthetic peptide corresponding to residues of the human Insulin protein | Human |
| Insulin | MMab | BSB-42 | IgG1/K | Recombinant full length human insulin protein | Human, Canine, Feline, Rat |
| Islet 1/ISL1 | RMab | EP283 | IgG | A synthetic peptide corresponding to residues of human Islet-1 protein | Human |
| Kappa | RPab | Polyclonal | IgG | Purified Kappa Light chains from human myeloma serum. | Human |
| Kappa Light Chains | MMab | BSB-58 (Kap-56) | IgG1/K | Purified Kappa Light chains from human myeloma serum | Human, Dog, Cat |
| Ki-67 | RMab | EP5 | IgG | Synthetic peptide corresponding to residues of the human Ki-67 protein | Human |
| Ki-67 | RMab | RM360 | IgG | Synthetic peptide corresponding to residues of the human Ki-67 protein | Human |
| Ksp-Cadherin | MMab | 4H6/F9 | IgG1 | Maltose-binding protein fusion protein containing the C-terminal 267 amino acids of the human isoform of Ksp-cadherin | Human |
| LAG-3/CD223 | RMab | EP294 | IgG | Synthetic peptide corresponding to residues of human LAG-3 (CD223) protein | Human |
| Lambda | MMab | BSB-16 (Lamb14) | IgG2a | Purified Kappa Light chains from human myeloma serum | Human, Dog, Cat |
| Lambda | RPab | Polyclonal | IgG | Purified human Lambda Light Chains from human serum. | Human |
| Lamin-B1 | RMab | RBT-LMNB1 | IgG | Recombinant human Lamin B1 protein | Human, Mouse, Rat |
| Laminin-R/RPSA | MMab | BSB-144 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human Laminin-R protein | Human, Mouse, Rat |
| Langerin | MMab | 12D6 | IgG2b | Recombinant protein representing the external domain of the Langerin molecule | Human |
| Langerin/CD207 | RMab | EP349 | IgG | Synthetic peptide corresponding to residues of human Langerin (CD207) protein | Human |
| LEF-1 | RMab | EP310 | IgG | Synthetic peptide corresponding to residues of human LEF-1 protein | Human |
| LH | MMab | BSB-53 | IgG1/K | Synthetic peptide corresponding to N-terminus of the human luteinizing hormone | Human |
| LIN28 | RMab | EP150 | IgG | Synthetic peptide corresponding to residues of the C-terminus of human LIN28A protein. | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|---|-------------------------|----------------|--|
| Testis, Lung, Stomach, Kidney, Transitional Cell Carcinoma | Membranous, Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity, Colon & G.I. Tract, Breast Cancer, Lung Cancer, Prostate Cancer, Ovarian Cancer, all Bladder & Pancreatic Cancer, Kidney Cancer and Urothelial |
| Placenta, Tonsil, Ovarian Carcinoma, TCC | Cytoplasmic, Nuclear | 1:25-1:100 | Gall Bladder & Pancreatic Cancer, Lung Cancer, Kidney & Urothelial Cancer, Melanoma & Skin Cancer, Cervical Cancer, Endometrial & Genital Cancer, Colon & GI Cancer, Cytopathology |
| Testis, Seminoma, Testicular Carcinoma | Cytoplasmic | 1:50-1:200 | Endometrial & Genital Cancer, Ovarian Cancer, Kidney & Urothelial Cancer |
| Placenta, Testis, Corpus Luteum, Adrenal Cortex, Testicular Carcinoma | Cytoplasmic | 1:25-1:100 | Endometrial & Genital Cancer, Ovarian Cancer, Kidney & Urothelial Cancer |
| Testis, Brain, Breast, Colon, Kidney, Pituitary, Adrenal, Prostate, Thyroid, Lung, Pancreas, Cervix, Salivary Gland, Astrocytoma, Lymphoblastic Lymphoma, Transitional Cell Carcinoma | Nuclear | 1:25-1:100 | Sarcoma & Soft Tissue, Neural & Neuroendocrine Cancer, Kidney & Urothelial Cancer |
| Testis, Brain, Breast, Colon, Kidney, Pituitary, Adrenal, Prostate, Thyroid, Lung, Pancreas, Cervix, Salivary Gland, Astrocytoma, Lymphoblastic Lymphoma, Transitional Cell Carcinoma | Nuclear | 1:25-1:100 | Sarcoma & Soft Tissue, Neural & Neuroendocrine Cancer, Kidney & Urothelial Cancer |
| Testis, Adrenal, Lung, Prostate, Liver, Placenta, Spleen | Cytoplasmic | 1:25-1:100 | Rejection & Autoimmunity |
| Pancreas, Colon, Tonsil, Neuroendocrine Lung Cancer, Endometrial & Colon Carcinomas | Nuclear | 1:50-1:200 | Lung Cancer, Neural & Neuroendocrine Cancer, Pituitary, Colon & GI Cancer |
| Pancreas, Colon, Tonsil, Neuroendocrine Lung Cancer, Endometrial & Colon Carcinomas | Nuclear | 1:25-1:100 | Lung Cancer, Neural & Neuroendocrine Cancer, Pituitary, Colon & GI Cancer |
| Pancreas | Cytoplasmic | 1:50-1:200 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer |
| Pancreas | Cytoplasmic | 1:250-1:1000 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer |
| Pancreas, Testis, Thyroid, Cervix, Skin, Pancreatic Neuroendocrine Cancer | Nuclear | 1:25-1:100 | Gall Bladder & Pancreatic Cancer |
| Tonsil, Lymph Node | Cytoplasmic | 1:250-1:1000 | Lymphoma, Rejection & Autoimmunity |
| Tonsil, Lymph Node | Cytoplasmic | 1:250-1:1000 | Lymphoma, Rejection & Autoimmunity |
| Testis, Tonsil, Bone Marrow, Placenta, Colon, Tonsil, Fallopian Tube, Astrocytoma, Breast Carcinoma, Colon Carcinoma | Nuclear | 1:25-1:100 | Breast Cancer, Cervical Cancer, Colon & GI Cancer, Lung Cancer, Lymphoma, Melanoma & Skin Cancer, Ovarian Cancer |
| Testis, Tonsil, Bone Marrow, Placenta, Colon, Tonsil, Fallopian Tube, Astrocytoma, Breast Carcinoma, Colon Carcinoma | Nuclear | 1:25-1:100 | Breast Cancer, Cervical Cancer, Colon & GI Cancer, Lung Cancer, Lymphoma, Melanoma & Skin Cancer, Ovarian Cancer |
| Kidney, Chromophobe Renal Cell Carcinoma | Cytoplasmic, Membranous | 1:250-1:1000 | Kidney & Urothelial Cancer |
| Testis, Lymph Node, Spleen, Lymphoblastic Lymphoma, TCC | Cytoplasmic, Membranous | 1:10-1:50 | Leukemia & Histiocytic, Rejection & Autoimmunity, Hodgkin's and NHD Lymphoma, Lymphoma, Immunotherapy |
| Tonsil, Lymph Node | Cytoplasmic | 1:250-1:1000 | Lymphoma, Rejection & Autoimmunity |
| Tonsil, Lymph Node | Cytoplasmic | 1:250-1:1000 | Lymphoma, Rejection & Autoimmunity |
| Colon, Breast, Fallopian Tube, Tonsil, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma, Colon Adenocarcinoma | Nuclear | 1:50-1:200 | Lung Cancer, Colon and GI Cancer, Gallbladder and Pancreatic Cancer |
| Placenta, Kidney, Prostate, Tonsil, Spleen, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Cytoplasmic, Membranous | 1:50-1:200 | Breast Cancer, Colon and GI Cancer, Gallbladder and Pancreatic Cancer, Prostate Cancer, Cervical Cancer, Lymphoma |
| Skin, Langerhans Histiocytosis | Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic, Sarcoma & Soft Tissues |
| Skin, Prostate, Breast, Liver, Cervix, Salivary Gland, Langerhans Histiocytosis | Cytoplasmic | 1:50-1:200 | Leukemia & Histiocytic, Sarcoma & Soft Tissues |
| Breast, Tonsil, Breast Carcinoma, Small Lymphocytic Lymphoma | Nuclear | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma, Colon & GI Cancer, Brain Cancer |
| Normal Pituitary | Cytoplasmic | 1:100-1:500 | Pituitary, Neural & Neuroendocrine Cancer |
| Testis, Seminoma, Dysgerminoma, Yolk Sac Tumor, Embryonal Carcinoma | Cytoplasmic | 1:25-1:100 | Ovarian Cancer, Testicular Cancer, Liver Cancer, Breast Cancer, Endometrial and Genital Cancer, Colon and GI Cancer, Germ Cell Tumors |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--|---------|----------------------|----------------------|---|--------------------------------------|
| LM02 | RMab | RBT-LM02 | IgG | Synthetic peptide corresponding to the N-terminal of the human LM02 protein. | Human |
| Lysozyme | RMab | EP134 | IgG | Synthetic peptide corresponding to residues of human Lysozyme C protein | Human, Predicted: Mouse |
| Lysozyme | RPab | Polyclonal | IgG | Purified human lysozyme from urine of patients with monocytic leukemia | Human |
| Macrophage HAM-56 | MMab | HAM-56 | IgM/K | Sonicated human alveolar macrophages | Human, Monkey |
| Mammaglobin | RMab | EP249 | IgG | Synthetic peptide corresponding to residues of human Mammaglobin A protein | Human, Mouse, Monkey, Predicted: Rat |
| Mart-1/Melan A | MMab | M2-7C10 | IgG2b/K | Recombinant human MART-1 protein | Human |
| Mart-1/Melan-A | MMab | A103 | IgG1 | Recombinant human Melan-A protein | Human |
| Maspin | MMab | BSB-92 | IgG2a | Synthetic peptide against the N-terminus of the human maspin protein | Human |
| MCM2 | RMab | RBT-MCM2 | IgG | Synthetic peptide corresponding to residues of human MCM2 protein | Human |
| MCM3 | RMab | EP202 | IgG | Synthetic peptide corresponding to residues of human MCM3 protein | Human, Predicted: Mouse, Rat |
| MCM5 | RMab | RBT-MCM5 | IgG | Synthetic peptide corresponding to residues of human MCM5 protein | Human |
| MDM2 | MMab | BSB-64 | IgG1 | Synthetic peptide against the N-terminus of the human MDM2 protein | Human |
| MDR-1 | MMab | JSB-1 | IgG1 | Multidrug-resistant Chinese hamster ovary cell line (ChR5) | Human |
| MDR-1 | RMab | EP271 | IgG | A protein fragment corresponding to human MDR-1 protein | Human, Mammalia |
| Melanoma Cocktail: HMB-45, MART-1 & Tyrosinase | MMab | HMB-45, A103 & BSB-6 | IgG1/K, IgG1 & IgG2a | Pigmented melanoma metastases from lymph nodes (HMB45); Recombinant human Melan-A protein (A103) and Recombinant full length human tyrosinase protein (BSB-6) | Human, Dog |
| Melanoma KBA.62 | MMab | KBA.62 | IgG1 | Human KAL cells derived from lymph node metastasis of malignant melanoma | Human |
| Melanoma PNL2 | MMab | PNL2 | IgG1 | Recombinant C-Terminal of the human SST2 protein | Human, Dog |
| Melanosome HMB45 | MMab | HMB-45 | IgG1/K | Purified pigmented melanoma metastases from lymph nodes (HMB45) | Human |
| Mesothelial Cell | MMab | HBME-1 | IgM/K | Human mesothelioma cells from patients with malignant epithelial mesothelioma | Human |
| Mesothelin | RMab | EP140 | IgG | Synthetic peptide corresponding to residues on C-terminus of human mesothelin protein | Human, Dog, Cat |
| MGMT/AGAT | RMab | EP337 | IgG | A synthetic peptide corresponding to residues of the human MGMT protein | Human |
| MiTF | MMab | C5 & D5 | IgG1/K & IgG1/K | N-terminal fragment of the human MiTF protein | Human, Dog, Cat |
| MLH1 | MMab | G168-728 | IgG2a | Full length recombinant MLH1 protein | Human, Rat, Mouse |
| MMP-9 | RMab | EP127 | IgG | Synthetic peptide corresponding to residues of human MMP-9 protein | Human, Predicted: Rat |
| MNDA | MMab | BSB-157 | IgG1 | Synthetic peptide corresponding to the C-terminus of the human MNDA protein | Human |
| MSH2 | RMab | RBT-MSH2 | IgG | Synthetic peptide corresponding to the N-terminus of the human MSH2 protein | Human |
| MSH2 | MMab | BSB-147 | IgG1 | Recombinant Human MSH2 Protein | Human |
| MSH6 | MMab | 44 | IgG1 | Synthetic peptide human MSH6 | Human, Rat, Mouse, Dog |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|--|-------------------------|----------------|---|
| | Tonsil, Spleen, Placenta, Follicular and Lymphoblastic Lymphoma | Cytoplasmic | 1:50-1:200 | Lymphoma |
| | Tonsil, Lymph Node, Liver, Kidney, Spleen, Salivary Gland, Cervix, Pancreas, Bone Marrow, Colon, Lung | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic |
| | Tonsil, Lymph Node, Liver, Kidney, Spleen, Lung, Salivary Gland, Cervix, Pancreas, Bone Marrow, Colon, Lymphoblastic Lymphoma, Transitional Cell Carcinoma | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic |
| | Tonsil, Lymph Node, Thymus, Spleen, Placenta, Colon | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic |
| | Breast, Skin, Fallopian Tube, Breast Carcinoma | Cytoplasmic | 1:25-1:100 | Breast Cancer, Carcinomas of Unknown Primary Site |
| | Normal Skin, Melanoma | Cytoplasmic | 1:250-1:1000 | Melanoma & Skin Cancer |
| | Normal Skin, Melanoma | Cytoplasmic | 1:250-1:1000 | Melanoma & Skin Cancer |
| | Prostate, Breast, Tonsil | Nuclear, Cytoplasmic | 1:25-1:100 | Breast Cancer, Prostate Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Lung Cancer |
| | HSIL, Cervical, Breast Cancer | Nuclear | 1:100-1:500 | Cervical Cancer |
| | Colon Carcinoma | Nuclear | 1:25-1:100 | Breast Cancer, Cervical Cancer, Head & Neck Cancer |
| | Tonsil, Testis, Breast, Colon, Spleen, Cervical Carcinoma, Colon Carcinoma, Transitional Cell Carcinoma, Lymphoblastic Lymphoma | Nuclear | 1:100-1:300 | Cervical Cancer |
| | Testis, Tonsil, Cervix, Placenta, LipoSarcoma & Soft Tissue, Testicular Cancer | Nuclear | 1:25-1:100 | Sarcoma & Soft Tissues, Breast Cancer |
| | Skeletal Muscle, Kidney, Adrenal, Liver | Cytoplasmic | 1:25-1:100 | Breast Cancer |
| | Skeletal Muscle, Kidney, Adrenal, Liver | Cytoplasmic | 1:50-1:200 | Breast Cancer |
| | Skin, Melanoma | Cytoplasmic | 1:100-1:500 | Melanoma & Skin Cancer |
| | Melanoma | Membranous | 1:50-1:200 | Melanoma & Skin Cancer |
| | Melanoma | Cytoplasmic | 1:100-1:500 | Melanoma & Skin Cancer, Sarcoma & Soft Tissue, Undifferentiated Tumor |
| | Melanoma | Cytoplasmic | 1:250-1:1000 | Melanoma & Skin Cancer |
| | Breast, Tonsil, Lung, Salivary Gland, TCC, Mesothelioma | Cytoplasmic, Membranous | 1:25-1:100 | Mesothelioma, Lung Cancer, Head & Neck Cancer, Gall Bladder & Pancreatic Cancer, Cytopathology |
| | Synthetic peptide corresponding to residues on C-terminus of human mesothelin protein | Cytoplasmic, Membranous | 1:50-1:200 | Mesothelioma, Lung Cancer, Head & Neck Cancer, Gall Bladder & Pancreatic Cancer, Cytopathology |
| | Breast, Cervix, Fallopian Tube, Prostate, Testis, Transitional Cell Carcinoma, Lung Adenocarcinoma, Ductal Breast Carcinoma, Endometrial Carcinoma | Nuclear, Cytoplasmic | 1:50-1:200 | Neural & Neuroendocrine Cancer, Breast Cancer, Leukemia and Histiocytic Cancers, Melanoma and Skin Cancer, Colon and GI Cancer, Gallbladder and Pancreatic Cancer, Lung Cancer, Prostate Cancer, Lymphoma |
| | Skin, Melanoma | Nuclear | 1:250-1:1000 | Skin, Melanoma |
| | Testis, Tonsil, Colon, Kidney, Colon Carcinoma | Nuclear | 1:10-1:50 | Colon & GI Cancer, Melanoma & Skin Cancer |
| | Tonsil, Spleen, Liver, Pituitary, Colon, Lymphoblastic Lymphoma, Transitional Cell Carcinoma | Cytoplasmic | 1:50-1:200 | Breast Cancer, Colon & GI Cancer |
| | Breast, Colon, Fallopian Tube, Brain, Tonsil, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Leukemia and Histiocytic Cancers, Lymphomas, Sarcoma and Soft Tissue Tumors |
| | Colon, Skin, Breast, Tonsil, Fallopian Tube, Colon Carcinoma | Nuclear | 1:50-1:200 | Colon & GI Cancer, Melanoma & Skin Cancer |
| | Colon, Skin, Breast, Tonsil, Fallopian Tube, Colon Carcinoma | Nuclear | 1:100-1:500 | Colon & GI Cancer, Melanoma & Skin Cancer |
| | Colon Mucosa, Colon Carcinoma | Nuclear | 1:50-1:200 | Colon & GI Cancer, Melanoma & Skin Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|----------------------------|---------|------------------|---------|--|--|
| MSH6 | RMab | EP49 | IgG | Synthetic peptide corresponding to Human MSH6 aa 1-100 (N terminal) | Human, Predicted: Mouse, Rat |
| MTAP | RMab | RBT-MTAP | IgG | Recombinant human MTAP protein | Human |
| MUC1 | RMab | EP85 | IgG | A synthetic peptide corresponding to residues on the C-terminus of human MUC1 protein | Human |
| MUC1 | MMab | BSB-44 | IgG1/K | Synthetic peptide corresponding to the N-terminus of human MUC1 | Human |
| MUC2 | MMab | BSB-45 | IgG1/K | Synthetic peptide corresponding to the C-terminus of human MUC2 | Human |
| MUC2 | RMab | EP187 | IgG | Synthetic peptide corresponding to residues of human MUC2 protein | Human |
| MUC4 | MMab | 8G7 | IgG1 | KLH-conjugated linear peptide corresponding to the beta chain tandem repeat region of human MUC4 | Human |
| MUC4 | RMab | EP256 | IgG | Synthetic peptide corresponding to residues of human MUC4 protein | Human |
| MUC5AC | MMab | CLH2 | IgG1 | A synthetic peptide of the human MUC5AC tandem repeat protein | Human |
| MUC6 | MMab | CLH5 | IgG1 | Synthetic peptide of the Gastric Mucin (Muc-6) tandem repeat sequence purified by HPLC | Human |
| MUM1 | RMab | EP190 | IgG | Synthetic peptide corresponding to residues of human IRF4 (MUM1) protein | Human, Canine, Feline |
| Musashi 2 | RMab | RM422 | IgG | Synthetic peptide corresponding to the internal region of human Musashi 2 | Human |
| Mycobacterium tuberculosis | RPab | Polyclonal | IgG | Purified PPD | Human |
| Myelin Basic Protein | RPab | Polyclonal | IgG | Synthetic peptide conjugated to KLH derived from within residues 150 to the C-terminus of Mouse Myelin Basic Protein | Human, Dog |
| Myelin Basic Protein | RMab | EP207 | IgG | Synthetic peptide corresponding to residues of human Myelin Basic Protein | Human |
| Myeloperoxidase | RMab | EP151 | IgG | Synthetic peptide corresponding to residues in human MPO protein | Human |
| Myeloperoxidase | RPab | Polyclonal | IgG | Purified human granulocytic MPO | Human, Dog, Cat |
| MyoD1 | RMab | EP212 | IgG | Synthetic peptide corresponding to residues of human MyoD1 protein | Human |
| Myogenin | MMab | F5D | IgG1/K | Recombinant protein containing rat myogenin aa 30-224 | Human, Canine, Rat, Mouse |
| Myoglobin | RMab | EP87 | IgG | Synthetic peptide corresponding to residues on the C-terminus of human myoglobin protein | Human |
| Myoglobin | MMab | BSB-104 | IgG1 | Recombinant protein representing the full length of the human Myoglobin protein | Human |
| Myosin Smooth Muscle | MMab | BSB-17 (SMM-H24) | IgG1/K | Synthetic peptide corresponding to residues of human myosin heavy chain protein | Human, Dog, Cat |
| Nanog | RMab | EP225 | IgG | A synthetic peptide corresponding to residues of human Nanog protein. | Human |
| Napsin A | RMab | EP205 | IgG | Synthetic peptide corresponding to residues of human Napsin A protein | Human |
| Napsin A | MMab | BSB-112 | IgG1/K | Synthetic peptide corresponding to residues of human Napsin A protein | Human |
| Nestin | RMab | EP287 | IgG | Synthetic peptide corresponding to C-terminus of human Nestin protein | Human |
| NeuN | MMab | A60 | IgG1 | Purified cell nuclei from mouse brain | Human, Avian, Chicken, Ferret, Mouse, Pig, Rat, Salamander |
| NeuN | RMab | RBT-NeuN | IgG | Recombinant NeuN human protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|---|-------------------------|----------------|--|
| Colon Mucosa, Colon Carcinoma | Nuclear | 1:50-1:200 | Colon & GI Cancer, Melanoma & Skin Cancer |
| Breast, Colon, Prostate, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma | Cytoplasmic, Nuclear, | 1:50-1:200 | Mesothelioma, Melanoma and Skin Cancer, Liver Cancer, Colon and GI Cancer, Prostate Cancer |
| Breast, Colon, Kidney, Fallopian Tube, Tonsil, Colon Adenocarcinoma | Cytoplasmic | 1:50-1:200 | Breast, Colon, Kidney, Fallopian Tube, Tonsil, Colon Adenocarcinoma |
| Breast, Colon, Kidney, Fallopian Tube, Tonsil, Colon Adenocarcinoma | Cytoplasmic | 1:250-1:1000 | Breast Cancer, Colon & GI Cancer |
| Small Intestine, Colon, Colon Adenocarcinoma | Cytoplasmic | 1:250-1:1000 | Colon & GI Cancer |
| Small Intestine, Colon, Colon Adenocarcinoma | Cytoplasmic | 1:25-1:100 | Colon & GI Cancer |
| Colon, Prostate, Cervix, Placenta, Salivary Gland, Pancreatic Carcinoma, Cervical Squamous Carcinoma, Bladder TCC | Cytoplasmic | 1:25-1:100 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Ovarian Cancer, Lung Cancer, Endometrial & Genital Cancer, Cervical Cancer, Sarcoma & Soft Tissue |
| Colon, Liver, Kidney, Pancreatic Carcinoma, Cervical Squamous Carcinoma, Bladder TCC | Cytoplasmic | 1:25-1:100 | Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Ovarian Cancer, Lung Cancer, Endometrial & Genital Cancer, Cervical Cancer, Sarcoma & Soft Tissue |
| Stomach, Colon, Kidney | Cytoplasmic | 1:25-1:100 | Colon & GI Cancer |
| Stomach | Cytoplasmic | 1:100-1:400 | Colon & GI Cancer |
| Tonsil, Spleen, Colon, Kidney, Breast, Lymph Node, Plasmacytoma, Hodgkin's Lymphoma | Cytoplasmic, Nuclear | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Lymphoma |
| Testis, Kidney, Colon, Transitional Cell Carcinoma, Breast Cancer | Nuclear, Cytoplasmic | 1:50-1:200 | Lung Cancer, Leukemias and Histiocytic Cancer, Cervical Cancer, Breast Cancer, Neural and Neuroendocrine Cancer, Liver Cancer, Gallbladder and Pancreatic Cancer |
| Infected Tissue | Cell Wall | 1:250-1:1000 | Infectious Diseases, Cytopathology |
| Brain, Neuroblastoma | Cytoplasmic | 1:50-1:200 | Neural & Neuroendocrine Cancer |
| Brain, Neuroblastoma | Cytoplasmic | 1:25-1:100 | Neural & Neuroendocrine Cancer |
| Bone Marrow | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic |
| Bone Marrow | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic |
| Fetal Muscle, RhabdomyoSarcoma & Soft Tissue | Nuclear | 1:50-1:200 | Sarcoma & Soft Tissue |
| RhabdomyoSarcoma & Soft Tissue | Nuclear | 1:50-1:200 | Sarcoma & Soft Tissue, Undifferentiated Tumor |
| Skeletal Muscle Tissue | Cytoplasmic | 1:100-1:500 | Sarcoma & Soft Tissue |
| Skeletal Muscle Tissue | Cytoplasmic, Membranous | 1:100-1:500 | Sarcoma & Soft Tissue |
| Intestine, Breast, Appendix | Cytoplasmic | 1:250-1:1000 | Breast Cancer |
| Testis, Cervix, Seminoma, Embryonal Carcinoma, TCC | Cytoplasmic | 1:50-1:200 | Ovarian Cancer, Testicular Cancer, Germ Cell Tumor, Breast Cancer, Lung Cancer |
| Kidney, Lung, Renal Cell Carcinoma, Lung Carcinoma | Cytoplasmic | 1:10-1:50 | Lung Cancer, Carcinoma of Unknow Primary Site, Cytopathology |
| Kidney, Lung, Renal Cell Carcinoma, Lung Carcinoma | Cytoplasmic | 1:100-1:500 | Lung Cancer, Carcinoma of Unknow Primary Site, Cytopathology |
| Kidney, Breast, Adrenal, Myometrium, Liver Carcinoma | Cytoplasmic | 1:50-1:200 | Melanoma & Skin Cancer, Breast Cancer |
| Brain | Nuclear | 1:25-1:100 | Neural & Neuroendocrine Cancer |
| Brain | Nuclear | 1:100-1:250 | Neural & Neuroendocrine Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|-------------------------|---------|------------|---------|---|--|
| Neurofilament | MMab | 2F11 | IgG1/K | Purified human neurofilament preparation. | Human, Canine, Rat, Mouse, Horse |
| Neurofilament | RMab | EP79 | IgG | Synthetic peptide corresponding to residues in human Neuro-filament M protein | Human |
| NGFR | MMab | BSB-18 | IgG1 | Recombinant human p75 nerve growth factor receptor protein | Human |
| NKX2.2 | RMab | EP336 | IgG | A synthetic peptide corresponding to residues of human NKX2.2 protein | Human, Predicted Mouse and Rat |
| NKX3.1 | RMab | EP356 | IgG | A synthetic peptide corresponding to residues of human NKX3.1 protein | Human |
| NKX3.1 | RMab | RM430 | IgG | A synthetic peptide corresponding to residues of human NKX3.1 protein. | Human |
| NPM1/B23 | MMab | BSB-124 | IgG1/K | Synthetic peptide corresponding to the C-terminals of the human NPM1/B23 protein. | Human |
| NRAS | RMab | RBT-NRAS | IgG | Synthetic peptide of human N-Ras protein containing the (Q61R) point mutation | Human |
| NSE | MMab | BSB-94 | IgG1/K | Purified human neuron-specific enolase (NSE) | Human |
| NUT/NUTM1 | RMab | RBT-NUTM1 | IgG | Recombinant protein corresponding to Human NUTM1 | Human |
| Oct-2 | RMab | EP115 | IgG | Synthetic peptide corresponding to residues of human OCT-2 protein | Human |
| Oct-4 | RMab | EP143 | IgG | Synthetic peptide corresponding to residues of human OCT-4 protein | Human |
| OLIG2 | RMab | EP112 | IgG | Synthetic peptide corresponding to residues in human OLIG2 protein | Human, Predicted: Mouse |
| Osteonectin/SPARC | MMab | BSB-93 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human SPARC protein | Human, Rat |
| OX-40/CD134 | MMab | BSB-90 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human OX40 protein | Human |
| p120 Catenin | RMab | EP66 | IgG | Synthetic peptide corresponding to residues in human p120 Catenin protein | Human, Predicted: Mouse, Rat, Rabbit |
| p14 ARF | RMab | RBT-p14 | IgG | Synthetic peptide corresponding to the C-terminus of the human p14 ARF/CDKN2A. | Human |
| p16 | MMab | 16P04, JC2 | IgG1 | Purified human recombinant full length p16 protein. | Human |
| p16 | RMab | RBT-p16 | IgG | Purified human recombinant full length p16 protein. | Human |
| p16 | RMab | RM267 | IgG | A peptide corresponding to the C-terminus of human p16INK4a (Cyclin-dependent kinase inhibitor 2A). | Human |
| p21 | MMab | DCS-60.2 | IgG2a | Purified human recombinant fusion p21 protein | Human |
| p27 | MMab | SX53G8 | IgG1/K | Purified glutathione S-transferase (GST)-p27Kip1 fusion protein | Human |
| p40 | RMab | ZR8 | IgG | Synthetic peptide corresponding to the N-terminal domain of human p63 | Human |
| p53 | MMab | D07 | IgG2b/K | Recombinant human wild-type p53 protein | Human, Cat, Cattle, Horse, Sheep, Bovine, Monkey |
| p57 | MMab | Kp10 | IgG2b/K | Recombinant human p57Kip2 protein | Human |
| p63 | MMab | 4A4 | IgG2a/K | Recombinant fragment corresponding to Human p63 aa 1-205. | Human |
| p63 | RMab | EP174 | IgG2a/K | Recombinant protein fragment corresponding to amino acids 1-203 of human p63 protein. | Human |
| Parafibromin | MMab | BSB-50 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human parafibromin protein | Human |
| Parathyroid Hormone/PTH | MMab | BSB-24 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human parathyroid hormone protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|---|-------------------------|----------------|---|
| | Brain | Cytoplasmic | 1:250-1:1000 | Neural & Neuroendocrine Cancer |
| | Brain | Cytoplasmic | 1:100-1:500 | Neural & Neuroendocrine Cancer |
| | Brain, Breast, Prostate, Neuroblastoma, CNS Tumor | Cytoplasmic | 1:250-1:1000 | Melanoma & Skin Cancer, Breast Cancer |
| | Pancreas, Brain, Pituitary, Colon, Ewing's Sarcoma & Soft Tissue | Nuclear | 1:50-1:200 | Sarcoma & Soft Tissue, |
| | Prostate, Prostate Carcinoma | Nuclear | 1:50-1:200 | Prostate Cancer, Breast Cancer, Carcinoma of Unknow Primary Site |
| | Prostate, Prostate Carcinoma | Nuclear | 1:50-1:200 | Prostate Cancer, Breast Cancer, Carcinoma of Unknow Primary Site |
| | Breast, Cervix, Testis, Pancreas, RCC, TCC | Cytoplasmic | 1:50-1:200 | Lymphoma, Leukemia & Histoctytic, Colon & GI Cancer. |
| | Small Intestine, Melanoma, ASTMA2, Lung Carcinoma | Membranous | 1:25-1:50 | Melanoma & Skin Cancer |
| | Pancreas, Brain, Pituitary, Adrenal, Thyroid | Cytoplasmic | 1:250-1:1000 | Neural & Neuroendocrine Cancer, Lung Cancer |
| | Testis | Nuclear | 1:50-1:200 | Head and Neck Cancer, Sarcoma and Sof Tissue Tumor, Lymphoma |
| | Tonsil, Lymph Node | Nuclear | 1:100-1:500 | Hodgkin's and NHD Lymphoma, Lymphoma |
| | Seminoma, Dysgerminoma, Testis Carcinoma | Nuclear | 1:50-1:200 | Ovarian Cancer, Testicular Cancer, Germ Cell Tumor |
| | Tonsil, Colon, Brain, Astrocytoma | Nuclear | 1:50-1:200 | Neural & Neuroendocrine Cancer |
| | Adrrenal, Testis, Placenta, TCC, Testicular Cancer, Cervical Cancer | Cytoplasmic | 1:100-1:500 | Breast Cancer, Prostate Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer |
| | Tonsil, Lymph Node, Thymus | Cytoplasmic, Membranous | 1:25-1:100 | Rejection & Autoimmunity, Melanoma & Skin Cancer; Immunotherapy |
| | Breast, Testis, Kidney, Prostate, Pancreas, Tonsil, Salivary Gland, Skin, Cervix, Colon, Malignant Melanoma, Transitional Cell Carcinoma, Breast Lobular Carcinoma | Cytoplasmic, Membranous | 1:100-1:500 | Breast Cancer |
| | Cervical, Anal and Ovarian Carcinomas | Cytoplasmic | 1:25-1:100 | Breast Cancer, Colon & GI Cancer, Cervical Cancer, Neural and Neuroendocrine Cancer |
| | Testis, NSCLC, TCC | Nuclear, Cytoplasmic | 1:25-1:100 | Cervical Cancer, Breast Cancer, Head & Neck Cancer, |
| | Testis, NSCLC, TCC | Cytoplasmic, Nuclear | RTU | Cervical Cancer, Breast Cancer, Head & Neck Cancer |
| | Testis, NSCLC, TCC | Cytoplasmic, Nuclear | 1:25-1:100 | Cervical Cancer, Breast Cancer, Head & Neck Cancer |
| | Tonsil, Colon, Fallopian Tube, Breast Cancer, Colon Carcinoma | Cytoplasmic, Nuclear | 1:50-1:200 | Endometrial & Genital Cancer |
| | Testis, Breast, Adrenal, Prostate, Tonsil, Lung, Colon, Non-Small Cell Lung Carcinoma, Colon Adenocarcinoma | Nuclear | 1:25-1:100 | Breast Cancer |
| | Prostate, Breast, Skin | Nuclear | 1:50-1:200 | Lung Cancer, Prostate Cancer, Breast Cancer, Melanoma & Skin Cancer, Carcinoma of Unknown Primary Site |
| | Lung , Breast, Ovarian Prostate, or Colon Carcinomas | Nuclear | 1:250-1:1000 | Breast Cancer, Colon & GI Cancer, Endometrial & Genital Cancer, Liver Cancer |
| | Placenta, Colon Carcinoma | Nuclear | 1:250-1:1000 | Germ Cell Tumor |
| | Prostate, Breast, Skin, Salivary Gland | Cytoplasmic | 1:50-1:200 | Cervical Cancer, Breast Cancer, Head & Neck Cancer, Kidney & Urotelial Cancer, Lung Cancer, Prostate Cancer, Thyroid & Parathyroid Cancer, Melanoma & Skin Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Breast, Skin, Salivary Gland | Cytoplasmic | 1:50-1:200 | Cervical Cancer, Breast Cancer, Head & Neck Cancer, Kidney & Urotelial Cancer, Lung Cancer, Prostate Cancer, Thyroid & Parathyroid Cancer, Melanoma & Skin Cancer, Carcinomas of Unknown Primary Site |
| | Parathyroid, Colon, Testis, Adrenal, Breast, Cervix, Kidney, Pituitary, Brain, Pancreas, Salivary Gland, Lymphoblastic Lymphoma, Bladder TCC, Transitional Cell Carcinoma | Nuclear | 1:25-1:100 | Thyroid & Parathyroid Cancer |
| | Parathyroid | Cytoplasmic, Membranous | 1:250-1:1000 | Thyroid & Parathyroid Cancer, Head & Neck Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|--------------------------|---------|------------|---------|--|--------------------------------|
| Parvalbumin | RMab | EP300 | IgG | Synthetic peptide corresponding to residues of human Parvalbumin protein | Human |
| Parvovirus | MMab | R92F6 | IgG1 | Native parvovirus B19 purified from human plasma | Human |
| PAX-2 | RMab | EP235 | IgG | Recombinant fragment corresponding to residues of the human PAX-2 protein | Human |
| PAX-5 | RMab | RBT-PAX5 | IgG | Synthetic peptide corresponding to the C-terminus of the human PAX-5 protein | Human, Mouse |
| PAX-6 | RMab | EP341 | IgG | A synthetic peptide corresponding to residues of human PAX6 protein | Human, Predicted Mouse and Rat |
| PAX-7 | MMab | BSB-145 | IgG2a | Recombinant PAX-7 of human origin | Human, Mouse, Rat |
| PAX-8 | RMab | ZR-1 | IgG | Synthetic peptide corresponding to the C-terminus of Human PAX8 protein | Human |
| PD-1/CD279 | RMab | EP239 | IgG | A synthetic peptide corresponding to residues of human PD-1 protein | Human |
| PD-1/CD279 | MMab | NAT-105 | IgG1 | TY cells (human T/NK cell Leukemia) | Human |
| PDGFR-B | RMab | RBT-PDGFRB | IgG | Highly pure (>98%) recombinant hPDGF-B (human Platelet Derived Growth Factor) | Human, Mouse, Rat |
| PD-L1/ CD274 | RMab | 28-8 | IgG | Recombinant full length protein corresponding to Human PD-L1 (extracellular) | Human |
| PD-L1/CD274 | RMab | RBT-PDL1 | IgG | Synthetic peptide corresponding to N-terminus residues of the human PD-L1 protein | Human |
| PDX1 | RMab | EP139 | IgG | A synthetic peptide corresponding to residues of human PDX1 protein | Human |
| PELP1 | RMab | RBT-PELP1 | IgG | Recombinant PELP1 protein | Human, Mouse, Rat |
| Perforin | MMab | 5B10 | IgG1 | Recombinant protein corresponding to C-terminal region of human perforin | Human |
| PGP 9.5 | MMab | BSB-46 | IgG1/K | Recombinant full length human PGP 9.5 protein | Human |
| PHOX2B | RMab | EP312 | IgG | Synthetic peptide corresponding to residues of human PHOX2B protein. | Human |
| PLA2R1 | MMab | BSB-129 | IgG1 | Recombinant fragment corresponding to residues of the human phospholipase A2 receptor 1, protein | Human |
| PLAP | RMab | EP194 | IgG | Synthetic peptide corresponding to residues of human PLAP protein | Human |
| PLAP | MMab | BSB-47 | IgG2b/K | Recombinant full length human PLAP protein | Human, Rat, Mouse |
| PMS2 | RMab | EP51 | IgG | Synthetic peptide corresponding to residues in human PMS2 protein | Human |
| Pygopus 2/Pygo 2 | MMab | BSB-156 | IgG2a | Synthetic peptide corresponding to the C-terminus of the human Pygopus 2 protein | Human, Mouse, Rat |
| Pneumocystis jirovecii | MMab | 3F6 | IgM/K | Pneumocystis jirovecii cysts isolated from human lung | Human |
| Podoplanin/D2-40 | MMab | D2-40 | IgG1 | Resected tissue from dysgerminoma of the ovary | Human, Rat, Mouse |
| PRAME | RMab | RBT-PRAME | IgG | Synthetic peptide corresponding to residues of human PRAME protein | Human |
| Prealbumin/Transthyretin | MMab | BSB-125 | IgG | Synthetic peptide corresponding to the C-terminus of the human Prealbumin protein. | Human |
| Progesterone Receptor | RMab | EP2 | IgG | Synthetic peptide corresponding to residues near the N-terminus of human progesterone receptor protein | Human |
| Progesterone Receptor | RMab | RBT-22 | IgG | Recombinant human Progesterone Receptor protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|--------------------------|----------------|--|
| Brain, Kidney, Tonsil, Lymph Node, Chromophobe RCC | Cytoplasmic | 1:25-1:100 | Kidney & Urothelial Cancer |
| Parvovirus Infected Tissue | Cytoplasmic, Nuclear | 1:100-1:500 | Infectious Diseases |
| Kidney, Fallopian Tube, Clear Cell Renal Carcinoma, Ovarian Serous Papillary Carcinoma | Nuclear | 1:50-1:200 | Kidney & Urothelial Cancer, Ovarian Cancer |
| Tonsil, Lymph Node, Spleen, Thymus, Colon, Liver, Lymphoblastic Lymphoma | Nuclear | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Leukemia & Histiocytic, Colon & GI Cancer |
| Pancreas, Pituitary, Neuroendocrine Tumors | Nuclear | 1:50-1:200 | Gall Bladder & Pancreatic Cancer, Ovarian Cancer, Colon & GI Cancer, Lung Cancer, Carcinomas of Unknown Primary Site |
| Placenta, Brain, Testis, Prostate, Papillary Thyroid Carcinoma, Transitional Cell Carcinoma | Nuclear, Cytoplasmic | 1:25-1:100 | Sarcoma & Soft Tissue, Neural and Neuroendocrine Cancer |
| Ovary, Thyroid | Nuclear | 1:25-1:100 | Kidney & Urothelial Cancer, Ovarian Cancer, Thyroid & Parathyroid Cancer, Carcinomas of Unknown Primary Site |
| Tonsil, Lymph Node, Thymus, Spleen | Cytoplasmic, Membranous | 1:50-1:200 | Lymphoma, Hodgkin's and NHD Lymphoma, Immunotherapy |
| Tonsil, Lymph Node, Thymus, Spleen | Cytoplasmic | 1:100-1:500 | Lymphoma, Hodgkin's and NHD Lymphoma, Immunotherapy |
| Placenta, Breast, Kidney, Cervix, Skin, Fallopian Tube, | Cytoplasmic | 1:10-1:50 | Leukemia & Histiocytic, Sarcoma & Soft Tissue, Endotelial |
| Tonsil, Placenta, Lymphoblastic Lymphoma | Membranous | 1:25-1:100 | Immunotherapy, Lymphoma, Hodgkin's and NHD Lymphoma, Lung Cancer, Melanoma, Kidney & Urothelial Cancer, Ovarian Cancer |
| Tonsil, Placenta, Lymphoblastic Lymphoma | Membranous | 1:25-1:100 | Immunotherapy, Lymphoma, Hodgkin's and NHD Lymphoma, Lung Cancer, Melanoma, Kidney & Urothelial Cancer, Ovarian Cancer |
| Pancreas, Colon, Liver, Pancreatic Cancer | Nuclear | 1:50-1:200 | Gall Bladder & Pancreatic Cancer, Colon & GI Cancer, Prostate Cancer |
| Placenta, Breast, Colon, Testis, Transitional Cell Carcinoma, T Cell Lymphoblastic Lymphoma, Seminoma | Nuclear | 1:50-1:200 | Breast Cancer, Endometrial Cancer, Ovarian Cancer, Prostate Cancer, Colon and GI Cancer |
| Spleen | Cytoplasmic, Perinuclear | 1:50-1:200 | Rejection & Autoimmunity |
| Brain, Testis, Colon, Pituitary, Nerve Tissue, Bowel Wall | Cytoplasmic | 1:100-1:500 | Lung Cancer, Gall Bladder & Pancreatic Cancer, Colon & GI Cancer |
| Adrenal, Neuroblastoma | Cytoplasmic | 1:25-1:100 | Neural & Neuroendocrine Cancer |
| Brain, Testis, Kidney, Salivary Gland, Gastic GIST | Cytoplasmic, Membranous | 1:50-1:200 | Rejection & Autoimmunity , Kidney & Urothelial Cancer |
| Placenta | Cytoplasmic | 1:100-1:500 | Ovarian Cancer, Testicular Cancer, Germ Cell Tumor, Undifferentiated Tumor |
| Placenta | Cytoplasmic | 1:250-1:1000 | Ovarian Cancer, Testicular Cancer, Germ Cell Tumor, Undifferentiated Tumor |
| Synthetic peptide corresponding to residues in human PMS2 protein | Nuclear | 1:25-1:100 | Colon & GI Cancer, Melanoma & Skin Cancer |
| Fallopian Tube, Adrenal Gland, Kidney, Transitional Cell Carcinoma, Lung Adenocarcinoma, Lung Neuroendocrine Cancer, Papillary Thyroid Carcinoma | Nuclear | 1:50-1:200 | Lung Cancer, Neural & Neuroendocrine Cancer, Colon and GI Cancer |
| Pneumocystis jirovecii Infected tissue | Membranous | 1:10-1:50 | Infectious Diseases, Cytopathology |
| Placenta, Breast, Lung, Cervix, Tonsil, Lymph Node, Lymphangioma | Cytoplasmic | 1:100-1:500 | Endothelial, Sarcoma & Soft Tissue, Mesothelioma, Lung Cancer |
| Tonsil, Seminoma | Cytoplasmic, Nuclear | 1:50-1:200 | Melanoma & Skin Cancer, Breast Cancer |
| Liver, Pancreas, Kidney, Leydig Cells, Hepatocellular Carcinoma, Seminoma | Cytoplasmic | 1:50-1:200 | Neural & Neuroendocrine Cancer, Liver Cancer |
| Breast, Myometrium, Cervix, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Endometrial & Genital Cancer, Carcinomas of Unknown Primary Site |
| Breast, Myometrium, Cervix, Breast Carcinoma | Nuclear | 1:50-1:200 | Breast Cancer, Endometrial & Genital Cancer, Carcinomas of Unknown Primary Site |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|-----------------------|---------|----------|---------|---|---|
| Progesterone Receptor | MMab | BSB-2 | IgG1 | Recombinant human Progesterone Receptor protein | Human |
| Prolactin | RMab | EP193 | IgG | Synthetic peptide corresponding to residues of human Prolactin protein | Human |
| Prolactin | MMab | PRL02 | IgG1/K | Purified human prolactin | Human |
| Prostein/P5015 | RMab | EP381 | IgG | Synthetic peptide corresponding to the N-terminus of human prostein | Human |
| PSA | MMab | BSB-7 | IgG1/K | Synthetic peptide corresponding to the C-terminus of human PSA protein | Human |
| PSA | RMab | RBT-PSA | IgG | Synthetic peptide corresponding to the C-terminus of human PSA protein | Human |
| PSAP | MMab | PASE/4LJ | IgG1 | Purified prostatic acid phosphate from human seminal plasma | Human |
| PSAP | RMab | EP53 | IgG | Synthetic peptide corresponding to residues of human PSAP protein | Human |
| PSMA | RMab | EP192 | IgG | Synthetic peptide corresponding to residues of human PSMA protein | Human |
| PSP94/MSMB | RMab | EP203 | IgG | Synthetic peptide corresponding to residues of human MSMB (PSP94) protein | Human |
| PTEN | RMab | RBT-PTEN | IgG | Synthetic peptide corresponding to the N-terminus of the human PTEN protein | Human |
| PTEN | RMab | RM265 | IgG | A peptide corresponding to the C-terminus of human PTEN | Human |
| PU.1 | RMab | EP18 | IgG | Synthetic peptide corresponding to residues near the N terminus of human transcription factor PU.1 protein | Human |
| Renal Cell Carcinoma | MMab | PN-15 | IgG1/K | Microsomal fraction of human renal cortical tissue homogenate | Human, Rat |
| Retinoblastoma/ Rb | MMab | 1F8 | IgG1/K | Recombinant human Rb protein | Human |
| ROS-1 | RMab | EP282 | IgG | Synthetic peptide corresponding to residues of human ROS1 protein. | Human |
| S-100 | MMab | 4C4.9 | IgG2a | Purified bovine brain S100 protein | Human, Dog, Cat, Mouse, Rat, Cattle |
| S100A1 | RMab | EP184 | IgG | Synthetic peptide corresponding to residues of human S100A1 protein | Human, Predicted: Mouse, Rat |
| S100A6 | RMab | EP313 | IgG | A synthetic peptide corresponding to residues of human S100A6 protein | Human |
| S100A8/MRP8 | RMab | EP90 | IgG | Synthetic peptide corresponding to N-terminal region of human MRP8 / S100A8 protein | Human, Predicted: Mouse |
| S100A9 | RMab | EP185 | IgG | Synthetic peptide corresponding to residues of human S100A9 protein | Human |
| S100Beta | RMab | EP32 | IgG | Synthetic peptide corresponding to residues on the C-terminus of human S100 Beta protein | Human, Predicted: Mouse, Rat, Goat, Zebrafish, Macaque Monkey |
| S100P | RMab | EP186 | IgG | Synthetic peptide corresponding to residues of human S100P protein | Human |
| SALL4 | RMab | EP299 | IgG | A synthetic peptide corresponding to residues of human SALL4 protein | Human |
| SALL4 | MMab | 6 E 3 | IgG1/K | A recombinant polypeptide corresponding to aminoacids 954-1054 of human SALL4 expressed as a GST fusion protein | Human, Mouse |
| SARS-CoV-2 | MMab | BSB-134 | IgG2b | Recombinant SARS-CoV-2 Nucleocapsid | SARS-CoV-2 virus |
| SATB2 | RMab | EP281 | IgG | A synthetic peptide corresponding to residues of human SATB2 protein | Human |
| SDHB | MMab | BSB-131 | IgG1/K | A Recombinant full length of the human succinate dehydrogenase iron-sulfur protein | Human, Mouse, Rat |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| | Tissue Control | Localization | Dilution Range | Application |
|--|--|-------------------------|----------------|---|
| | Breast, Myometrium, Cervix, Breast Carcinoma | Nuclear | 1:100-1:500 | Breast Cancer, Endometrial & Genital Cancer, Carcinomas of Unknown Primary Site |
| | Normal Pituitary | Cytoplasmic | 1:100-1:500 | Pituitary, Neural & Neuroendocrine Cancer |
| | Normal Pituitary | Cytoplasmic | 1:100-1:500 | Pituitary, Neural & Neuroendocrine Cancer |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:50-1:200 | Prostate Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:250-1:1000 | Prostate Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:250-1:1000 | Prostate Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:250-1:1000 | Prostate Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:50-1:200 | Prostate Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:50-1:200 | Prostate Cancer, Carcinomas of Unknown Primary Site |
| | Prostate, Prostatic Adenocarcinoma | Cytoplasmic | 1:25-1:100 | Prostate Cancer |
| | Breast, Prostate, Breast & Prostatic Adenocarcinoma | Cytoplasmic, Nuclear | 1:25-1:100 | Breast Cancer, Endometrial & Genital Cancer, Neural & Neuroendocrine Cancer, Prostate Cancer |
| | Colon, Thymus, Skin, Kidney, Breast, Prostate, Breast Carcinoma, Prostatic Carcinoma | Cytoplasmic, Nuclear | 1:50-1:200 | Breast Cancer, Endometrial & Genital Cancer, Neural & Neuroendocrine Cancer, Prostate Cancer |
| | Tonsil, Lymp Node | Nuclear | 1:50-1:200 | Leukemia & Histiocytic, Lymphoma, Hodgkin's & NH Lymphoma |
| | Kidney, Breast, Thyroid, Renal Cell Carcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Kidney & Urotelial Cancer |
| | Colon, Breast, Skin, Fallopian Tube, Tonsil, Colon Carcinoma | Nuclear | 1:25-1:100 | Endometrial & Genital Cancer |
| | Placenta, Lung, SiHa Cells, NSCL ROS1 + | Cytoplasmic | 1:50-1:200 | Lung Cancer, Neural and Neuroendocrine Cancer, Ovarian Cancer, Colon and GI Cancer |
| | Melanoma | Cytoplasmic, Nuclear | 1:100-1:500 | Carcinomas of Unknown Primary Site, Melanoma & Skin Cancer, GIST, Undifferentiated Tumor |
| | Brain, Breast, Thyroid, Tonsil, Pancreas, Salivary Gland, Renal Onocytomas, Clear Cell Carcinoma | Cytoplasmic, Nuclear | 1:25-1:100 | Kidney & Urotelial Cancer |
| | Testis, Kidney, Breast, Tonsil, Colon, Carcinomas | Nuclear, Cytoplasmic | 1:50-1:200 | Kidney & Urotelial Cancer, Gall Bladder & Pancreatic Cancer, Colon & GI Cancer, Thyroid & Parathyroid Cancer |
| | Tonsil, Breast, Liver, Fallopian Tube, Lymphoblastic Lymphoma | Cytoplasmic, Nuclear | 1:10-1:50 | Gall Bladder & Pancreatic Cancer, Breast Cancer |
| | Tonsil, Liver, Lung, Breast, Cervix, Bone Marrow, Bladder TCC, Thyroid Carcinomas | Cytoplasmic, Nuclear | 1:10-1:50 | Rejection & Autoimmunity, Lung Cancer, Breast Cancer, Liver Cancer |
| | Melanoma | Cytoplasmic, Nuclear | 1:100-1:500 | Carcinomas of Unknown Primary Site, Melanoma & Skin Cancer, GIST |
| | Colon, Prostate, Pancreatic and Lung Carcinomas | Cytoplasmic, Nuclear | 1:25-1:100 | Gall Bladder and Pancreatic Cancer, Prostate Cancer, Lung Cancer, Colon & GI Cancer |
| | Testis, Seminoma, Yolk Sac Tumor | Nuclear | 1:10-1:50 | Ovarian Cancer, Testicular Cancer, Liver Cancer, Breast Cancer, Endometrial and Genital Cancer, Colon and GI Cancer, Germ Cell Tumors, Undifferentiated Tumor |
| | Testis, Seminoma, Yolk Sac Tumor | Nuclear | 1:50-1:200 | Ovarian Cancer, Testicular Cancer, Liver Cancer, Breast Cancer, Endometrial and Genital Cancer, Colon and GI Cancer, Germ Cell Tumors, Undifferentiated Tumor |
| | SARS-CoV-2 Infected Tissues | Cytoplasmic | 1:25-1:100 | Infectious Diseases |
| | Colon, Brain, Colon Carcinoma | Nuclear | 1:50-1:200 | Colon, Brain, Colon Carcinoma |
| | Breast, Adrenal, Prostate, Kidney, Spleen, Tonsil, Breast Carcinoma, Hepatocellular Carcinoma, Lung Adeno Carcinoma, Prostate Carcinoma, Papillary Thyroid Carcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Neural & Neuroendocrine Cancer, Pituitary, Kidney & Urotelial Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|-------------------------------|---------|-------------------|---------|---|--------------------------------------|
| SF-1/Steroidogenic Factor 1 | MMab | BSB-149 | IgG2a | Synthetic Peptide corresponding to the internal region of the SF-1 of human origin | Human, Mouse, Rat |
| SMAD4/DPC4 | MMab | BSB-63 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human SMAD4 protein | Human |
| SMAD4/DPC4 | RMab | RBT-SMAD4 | IgG | Synthetic peptide corresponding to the C-terminus of the human SMAD4 protein | Human |
| Smoothelin | MMab | R4A | IgG1 | Cytoskeletal smoothelin extract of chicken gizzard | Human |
| Somatostatin | RPab | Polyclonal | IgG | Synthetic cyclic peptide (1-14) of somatostatin | Human, Dog, Horse, Lizard |
| Somatostatin | MMab | BSB-113 | IgG2a/K | Synthetic peptide corresponding to the N-terminus of the human somatostatin protein | Human |
| Somatostatin | RMab | EP130 | IgG | Synthetic peptide corresponding to residues of human Somatostatin-28 protein | Human |
| Somatostatin Receptor 2/SSTR2 | RMab | EP149 | IgG | Synthetic peptide corresponding to residues on the C-terminus of the human SSTR2 protein | Mouse, Rat, Human |
| SOX-10 | RMab | EP268 | IgG | Recombinant fragment corresponding to residues in human SOX-10 protein | Human, Predicted: Mouse, Rat |
| SOX-10 | MMab | BSB-62 | IgG2b/K | Recombinant fragment corresponding to N-Terminus of human SOX-10 | Human |
| SOX-11 | MMab | CL0142 | IgG2a | Recombinant fragment corresponding to amino acids 56-166 of human SOX-11 | Human |
| SOX-11 | MMab | BSB-167 | IgG1 | Synthetic peptide from human SOX-11 protein | Human |
| SOX-2 | RMab | EP103 | IgG | Synthetic peptide corresponding to residues in human SOX2 protein | Human, Predicted: Mouse |
| SOX-2 | RMab | RM427 | IgG | A peptide corresponding to the C-terminus of human SOX2 | Human, Predicted: Mouse |
| SOX-9 | RMab | EP317 | IgG | A synthetic peptide corresponding to residues of human SOX9 protein | Human, Predicted Mouse and Rat |
| Spectrin | MMab | RBC2/3D5 | IgG2b/K | Native spectrin from human red blood cells membrane | Human |
| STAR | RMab | EP226 | IgG | Synthetic peptide corresponding to residues of human Steroidogenic Acute Regulatory Protein (STAR) protein. | Human |
| STAT6 | RMab | EP325 | IgG | Synthetic peptide corresponding to residues of human STAT6 protein | Human |
| Surfactant protein D/SP-D | MMab | BSB-162 | IgG1 | Synthetic peptide corresponding to the N-terminus of the human SP-D protein | Cytoplasmic, Membranous |
| Stathmin | RMab | EP247 | IgG | Synthetic peptide corresponding to residues of human the stathmin protein | Human, Predicted: Mouse, Rat |
| Survivin | RMab | EP119 | IgG | Synthetic peptide corresponding to residues on the N-terminus of the human Survivin protein | Human |
| SV40 | MMab | Pab101 | IgG2a | SV40 transformed mouse line B4 | Human |
| Synaptophysin | RMab | EP158 | IgG | A synthetic peptide corresponding to residues on the C-terminus (cytoplasmic domain) of human Synaptophysin protein | Human, Predicted: Mouse, Rat, Donkey |
| Synaptophysin | RPab | Polyclonal | IgG | Synthetic peptide conjugated to KLH derived from within residues 1 - 100 of Human Synaptophysin. | Human, Dog, Cat |
| TAG-72 | MMab | BSB-21 (Tag72-22) | IgG1/K | Recombinant TAG-72 human tumor associated glycoprotein | Human |
| Tau | MMab | BSB-115 | IgG1/K | Synthetic peptide conjugated to KLH corresponding to C-terminal residues of the human Tau protein | Human, Mouse, Rat, Canine |
| T-Bet/TBX-2 | RMab | EP263 | IgG | Synthetic peptide corresponding to residues of human TBX21 protein | Human |
| TCL1 | RMab | EP105 | IgG | Synthetic peptide corresponding to residues in human TCL1 protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|----------------------------------|----------------|---|
| Breast, Fallopian Tube, Colon, Bone Marrow, Testis, Transitional Cell Carcinoma, Lung Adenocarcinoma, Papillary Thyroid Carcinoma, Prostate Adenocarcinoma | Nuclear | 1:50-1:200 | Germ Cell Tumor Antibodies, Ovarian Cancer, Endometrial Cancer |
| Pancreas, Thyroid, Placenta, Cervix, TCC & Lymphoblastic Lymphoma | Cytoplasmic, Nuclear | 1:50-1:200 | Gall Bladder and Pancreatic Cancer, Liver Cancer; Colon & GI Cancer |
| Pancreas, Testis, Breast, Bone Marrow, Colon | Nuclear, Cytoplasmic | 1:50-1:200 | Gall Bladder and Pancreatic Cancer, Liver Cancer; Colon & GI Cancer |
| Leiomyoma, Bladder, Prostate, Myometrium, Fallopian Tube, Colon, Transitional Cell Carcinoma | Cytoplasmic | 1:25-1:100 | Kidney & Urothelial Cancer |
| Pancreas | Cytoplasmic | 1:100-1:500 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer, Colon and GI Cancer |
| Pancreas, Colon, Adrenal, Brain | Cytoplasmic | 1:250-1:1000 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer, Colon and GI Cancer |
| Pancreas, Brain, Pituitary | Cytoplasmic | 1:250-1:1000 | Gall Bladder & Pancreatic Cancer, Neural & Neuroendocrine Cancer, Colon and GI Cancer |
| Placenta, Brain, Testis, Prostate, Papillary Thyroid Carcinoma, Transitional Cell Carcinoma | Nuclear, Cytoplasmic, Membranous | 1:25-1:100 | Neural & Neuroendocrine Cancer, Lung Cancers, Breast Cancer, Gallbladder and Pacreatic Cancer |
| Skin, Salivary Gland, Breast, Melanoma, Schwannoma | Nuclear | 1:50-1:200 | Melanoma & Skin Cancer, Head & Neck Cancer, Undifferentiated Tumor |
| Skin, Salivary Gland, Breast, Melanoma, Schwannoma | Nuclear | 1:25-1:100 | Melanoma & Skin Cancer, Head & Neck Cancer, Undifferentiated Tumor |
| Mantle Cell Lymphoma, Lung Neuroendocrine Carcinoma | Nuclear | 1:25-1:100 | Lymphomas, Lung Cancer |
| Mantle Cell Lymphoma, Lung Neuroendocrine Carcinoma | Nuclear | 1:25-1:100 | Lymphomas, Lung Cancer |
| Brain, Oligodendrogloma, Squamous Cell Carcinoma | Nuclear | 1:50-1:200 | Lung Cancer, Germ Cell Tumor |
| Brain, Oligodendrogloma, Squamous Cell Carcinoma | Nuclear | 1:50-1:200 | Lung Cancer, Germ Cell Tumor |
| Colon, Prostate, Skin, Breast, Tonsil, Lymph Node, Colon Carcinoma | Nuclear | 1:50-1:200 | Colon & GI Cancer, Lung Cancer, Prostate Cancer, Gall Bladder & Pancreatic Cancer, Liver Cancer, Ovarian Cancer |
| Bone Marrow, Spleen | Membranous | 1:25-1:100 | Hematopoietic, Leukemia & Histiocytic |
| Leydig cells of Testis, Adrenal & Leydig Cell Tumors | Cytoplasmic | 1:50-1:200 | Pituitary, Ovarian Cancer, Kidney and Urothelial Cancer, Testicular Cancer |
| Solitary Fibrous Tumor | Nuclear | 1:50-1:200 | Sarcoma & Soft Tissues, Lung Cancer |
| Placenta, Lung, Pancreas, Adrenal Gland, Lung Squamous Cell Carcinoma, Lung Adenocarcinoma, Gastric GIST | Human | 1:25-1:100 | Lung Cancer, Ovarian Cancer, Gallbladder and Pancreatic Cancer |
| Testis, Tonsil, HSIL Cervical Carcinoma, Lymphoblastic Lymphoma, Bladder TCC | Cytoplasmic, Membranous | 1:50-1:200 | Cervical Cancer, Endometrial & Genital Cancer, Ovarian Cancer, Cytopathology |
| Colon, Placenta, Testis, Tonsil, Bone Marrow, Colon Carcinoma, Lymphoblastic Lymphoma, Bladder TCC | Nuclear | 1:50-1:200 | Kidney & Urothelial Cancer, Ovarian cancer, Liver Cancer, Prostate Cancer, Breast Cancer |
| SV40 Infected Tissue | Nuclear | 1:25-1:100 | Infectious Diseases |
| Pancreas, Brain, Pituitary, Adrenal, Colon | Cytoplasmic | 1:100-1:500 | Lung Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Sarcoma & Soft Tissue, Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Undifferentiated Tumor, Carcinomas of Unknown Primary Site |
| Pancreas, Brain, Pituitary, Adrenal, Colon | Cytoplasmic | 1:250-1:1000 | Lung Cancer, Colon & GI Cancer, Gall Bladder & Pancreatic Cancer, Sarcoma & Soft Tissue, Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer, Undifferentiated Tumor, Carcinomas of Unknown Primary Site |
| Testis, Spleen, Cervix, Colon, Breast Carcinoma | Cytoplasmic | 1:250-1:1000 | Breast Cancer, Gall Bladder & Pancreatic Cancer, Mesothelioma |
| Brain, Kidney, Pituitary, Pancreas, Cervix, Skin, Salivary Gland, Astrocytoma | Nuclear, Cytoplasmic | 1:100-1:500 | Neural & Neuroendocrine Cancer |
| Spleen, Tonsil, Cervix, Liver, Breast, Hairy Cell Leukemia & Histiocytic, Lymphoblastic lymphoma, Bladder TCC | Nuclear | 1:25-1:100 | Lymphoma, Leukemia & Histiocytic |
| Tonsil, Lymph Node | Nuclear, Cytoplasmic | 1:100-1:500 | Lymphoma, Leukemia & Histiocytic |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|---------------------------|---------|--------------------|-----------------|--|---------------------------------------|
| TCR Alpha | MMab | BSB-126 | IgG1/K | Synthetic peptide corresponding the N-Terminus of the human TCR alpha protein. | Human |
| TCR Beta | MMab | BSB-117 | IgG1/K | Human TCR β chain constant region | Human |
| TCR Delta | MMab | BSB-127 | IgG1/K | Recombinant protein corresponding the human TCR delta protein. | Human |
| TDP-43/TARDBP | MMab | BSB-166 | IgG2a | Synthetic peptide corresponding to the N-terminus of the human TARDBP protein | Human, Mouse, Rat, Horse, Dog, Bovine |
| TdT | RMab | RBT-TdT | IgG | Synthetic peptide conjugated to KLH corresponding to N-terminal residues of the human TdT protein | Human |
| TdT | RMab | EP266 | IgG | Synthetic peptide corresponding to residues of the human TdT protein | Human |
| TdT | RPab | Polyclonal | IgG | Synthetic peptide conjugated to KLH corresponding to N-terminal residues of the human TdT protein | Human |
| TFE3 | RMab | EP285 | IgG | A synthetic peptide corresponding to residues of human TFE3 protein | Human |
| TIGIT | MMab | BSB-152 | IgG2b | Recombinant human TIGIT protein | Human, Rat |
| TIM-3/HAVCR2/CD366 | MMab | BSB-163 | IgG2c | Recombinant extracellular domain of the human TIM3 protein | Human |
| Thrombomodulin/ CD141 | RMab | EP175 | IgG | A synthetic peptide corresponding to residues at the C-terminus of human thrombomodulin protein | Human |
| Thyroglobulin | RMab | EP250 | IgG | Synthetic peptide corresponding to residues of human Thyroglobulin protein | Human |
| Thyroglobulin | MMab | BSB-49 | IgG1 | Recombinant human thyroglobulin protein | Human, Cat |
| Thyroglobulin | MMab | 2H11 & 6E1 | IgG1/K & IgG1/K | Purified human thyroglobulin protein | Human, Dog, Cat |
| Thymidylate synthase/TS | MMab | BSB-160 | IgG2a | Synthetic peptide corresponding to the C-terminus of human TS protein | Human, Mouse, Rat |
| Thyroid Peroxidase | RMab | EP159 | IgG | Synthetic peptide corresponding to residues in the C-terminus of human TPO protein | Human |
| TIA-1 | MMab | TIA-1 | IgG1 | Recombinant fragment corresponding to a region within amino acids 1 and 165 of Human TIA1 | Human |
| TIA-1 | RMab | RBT-TIA1 | IgG | Recombinant fragment corresponding to the N Terminus of the human TIA1. | Human |
| TLE1 | MMab | 1F5 | IgG1/K | Synthetic peptide corresponding to amino acids amino acids 200-350 of TLE1 of mouse origin | Human |
| TLE-1 | MMab | BSB-142 | IgG2a | Recombinant human TLE1 fragment | Human, Mouse |
| TMPRSS2 | MMab | BSB-136 | IgG1 | Synthetic peptide corresponding to residues in the internal region of the human TMPRSS2 protein | Human, Mouse, Rat |
| TNF α -IP2 | MMab | BSB-141 | IgG1 | Recombinant protein corresponding to the TNF α -induced protein 2 of human origin | Human, Mouse, Rat |
| Topoisomerase II α | RMab | RBT-Topo2 α | IgG | Synthetic peptide corresponding to the C-terminal residues of human Topoisomerase II alpha protein | Human |
| Toxoplasma gondii | RPab | Polyclonal | IgG | Toxoplasma gondii | Human |
| TRAcP | MMab | 9C5 | IgG2b | Purified human TRAcP protein | Human |
| Treponema pallidum | RPab | Polyclonal | IgG | Purified treponema pallidum | Eubacteria |
| TRK | RMab | RBT-TRK | IgG | Synthetic peptide corresponding to the C-terminal residues of human pan TRK protein | Human |
| TRK | RMab | RM423 | IgG | Synthetic peptide corresponding to the C-terminus of the human TRK protein | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range





| | Tissue Control | Localization | Dilution Range | Application |
|--|---|----------------------------------|----------------|--|
| | Thymus, Tonsil, Lymph Node, Spleen, Hodgkin's Lymphoma & BDCM & CCRF-CEM Cell Line | Cytoplasmic | 1:25-1:100 | Lymphoma, Leukemia & Histiocytic |
| | Tonsil, Lymph Node | Cytoplasmic | 1:25-1:100 | Lymphoma, Leukemia & Histiocytic |
| | Thymus, Tonsil, Lymph Node, Spleen | Cytoplasmic, Membranous | 1:10-1:50 | Lymphoma, Leukemia & Histiocytic |
| | Breast, Fallopian Tube, Testis, Skin, Transitional Cell Carcinoma, Glioblastoma | Nuclear, Cytoplasmic | 1:50-1:200 | Breast Cancer, Melanoma and Skin Cancer, Liver Cancer, Lung Cancer, Sarcoma and Soft Tissue Cancer, Endometrial Cancer, Gastric Cancer, Ovarian Cancer |
| | Thymus, Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Lymphoma, Leukemia & Histiocytic |
| | Thymus, Lymphoblastic Lymphoma | Nuclear | 1:50-1:200 | Lymphoma, Leukemia & Histiocytic |
| | Thymus, Lymphoblastic Lymphoma | Nuclear | 1:100-1:500 | Lymphoma, Leukemia & Histiocytic |
| | Testis, Adrenal, Kidney, Testicular Cancer, RCC with Xp11.2 translocation, Alveolar Soft Part Sarcoma & Soft Tissue | Nuclear | 1:50-1:200 | Kidney & Urotelial Cancer, Sarcoma & Soft Tissue, Carcinoma of Unknown Primary Site |
| | Cervix, Fallopian Tube, Skin, Prostate, Testis, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma, Lung Adenocarcinoma | Cytoplasmic, Membranous | 1:50-1:200 | Immunotherapy, Melanoma and Skin Cancer |
| | Colon, Testis, Tonsil, Liver | Cytoplasmic | 1:25-1:100 | Immunotherapy, Leukemia and Histiocytic Cancer, Lymphoma, Breast Cancer, Melanoma and Skin Cancer |
| | Placenta, Liver, Kidney, Tonsil, Cervix, Bladder, Mesothelioma, Transitional Cell Carcinoma | Membranous, Cytoplasmic | 1:10-1:50 | Kidney & Urotelial Cancer, Mesothelioma |
| | Thyroid, Thyroid Carcinoma | Cytoplasmic | 1:250-1:1000 | Thyroid & Parathyroid Cancer, Head & Neck Cancer |
| | Thyroid, Thyroid Carcinoma | Cytoplasmic | 1:250-1:1000 | Thyroid & Parathyroid Cancer, Head & Neck Cancer |
| | Thyroid, Thyroid Carcinoma | Cytoplasmic | 1:250-1:1000 | Thyroid & Parathyroid Cancer, Head & Neck Cancer |
| | Bone Marrow, Colon, Tonsil, Testis, T Cell Lymphoblastic Lymphoma, Ductal Breast Carcinoma | Nuclear, Cytoplasmic, Membranous | 1:25-1:100 | Gallbladder and Pancreatic Cancer, Breast Cancer, Ovarian Cancer, Lung Cancer, Head and Neck Cancer |
| | Thyroid, Thyroid Carcinoma | Cytoplasmic | 1:100-1:500 | Thyroid & Parathyroid Cancer, Head & Neck Cancer |
| | Tonsil, Spleen, Liver, Anaplastic Large Cell Lymphoma | Cytoplasmic (Granular) | 1:25-1:100 | Hodgkin's and NHD Lymphoma, Lymphoma, Leukemia & Histiocytic, Melanoma & Skin Cancer |
| | Tonsil, Spleen, Anaplastic Large Cell Lymphoma | Cytoplasmic | 1:50-1:200 | Hodgkin's and Non-Hodgkin Lymphoma, Lymphoma, Leukemia & Histiocytic, Melanoma & Skin Cancer |
| | Synovial Sarcoma | Cytoplasmic, Nuclear | 1:10-1:50 | Sarcoma & Soft Tissue |
| | Synovial Sarcoma | Cytoplasmic, Nuclear | 1:50-1:200 | Sarcoma & Soft Tissue |
| | Testis, Colon, Kidney, Brain, Stomach, Pancreas, Prostate | Nuclear, Membranous | 1:50-1:200 | Infectious Diseases, Prostate Cancer, Colon & GI Cancer, Kidney Cancer, and Lung Cancer |
| | Testis, Tonsil, Lung, Kidney | Membranous, Cytoplasmic | 1:50-1:200 | Rejection & Autoimmunity, Kidney & Urothelial, Breast Cancer, Head and Neck Cancer, Gall Bladder & Pancreatic Cancer, Lymphoma, Ovarian Cancer, Cervical Cancer, Infectious Diseases |
| | Testis, Skin, Colon, Fallopian Tube, Tonsil, Lymph Node, Spleen, HSIL Cervical Cancer, Breast Cancer, Bladder TCC | Nuclear | 1:50-1:200 | Breast Cancer, Cervical Cancer |
| | Toxoplasma gondii Infected Tissue | Nuclear | 1:50-1:200 | Infectious Diseases |
| | Tonsil, Spleen, Lymph Node, Hairy Cell Leukemia | Cytoplasmic | 1:100-1:500 | Leukemia & Histiocytic, Lymphoma |
| | Infected Tissue | Cell Wall | 1:250-1:1000 | Infectious Diseases |
| | Brain, Lung Neuroendocrine | Nuclear, Cytoplasmic | 1:10-1:50 | Lung Cancer, Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer |
| | Brain, Lung Neuroendocrine | Nuclear, Cytoplasmic | 1:25-1:100 | Lung Cancer, Neural & Neuroendocrine Cancer, Thyroid & Parathyroid Cancer |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Antibody Name | Ab Type | Clone | Isotype | Immunogen | Species Reactivity |
|------------------------|---------|---------------|------------|---|---|
| Trop-2/EGP-1 | MMab | BSB-148 | IgG1 | Synthetic peptide corresponding to the extracellular domain of the human TROP-2 protein | Human, Mouse, Rat |
| Tryptase | RMab | EP259 | IgG | Recombinant protein of human Tryptase protein | Human |
| Tryptase | MMab | G3 | IgG1 | Purified human tryptase protein | Human |
| TSH | MMab | BSB-56 | IgG1/K | Recombinant human thyroid stimulating hormone receptor protein | Human |
| TSH | RMab | EP254 | IgG | Synthetic peptide corresponding to residues of human TSH (subunit beta) protein | Human, Predicted: Rabbit |
| TTF-1 | MMab | 8G7G3/1 | IgG1 | Recombinant rat TTF-1 protein | Human, Dog |
| Tyrosinase | MMab | BSB-6 (Ty/G5) | IgG2a | Recombinant full length human tyrosinase protein | Human |
| Uroplakin III | RMab | EP321 | IgG | A synthetic peptide corresponding to residues of human Uroplakin III protein | Human |
| Uroplakin III | RPab | Polyclonal | IgG | Synthetic peptide corresponding to the C-terminus of the human uroplakin III protein | Human |
| Varicella Zoster Virus | MMab | | SG1, SG1-S | Varicella-Zoster virus lysate | Human |
| VEGF | RMab | | IgG | Recombinant human vascular endothelial growth factor protein | Human |
| Villin | RMab | EP163 | IgG | A synthetic peptide corresponding to residues in human Villin-1 protein | Human |
| Villin | MMab | CWWB1 | IgG1 | Human Villin protein | Human |
| Vimentin | RMab | EP21 | IgG | Synthetic acetylated peptide corresponding to the C-term of human Vimentin protein was used | Human, Predicted: Mouse, Rat, Rhesus Monkey |
| Vimentin | MMab | V9 | IgG1/K | Synthetic acetylated peptide corresponding to the C-term of human Vimentin protein was used | Human, Dog, Cat, Rat, Rabbit, Cattle, Horse, Guinea Pig |
| WT1 | MMab | 6F-H2 | IgG1/K | Recombinant protein corresponding to amino acids 1-181 of human WT1 | Human |
| YAP1 | MMab | BSB-146 | IgG1 | Synthetic peptide corresponding to the C-terminus of the human YAP1 protein | Human, Mouse, Rat |
| Zap-70 | MMab | 2F3.2 | IgG2a | GST-fusion to tandem SH2 domains of human ZAP-70 corresponding to residues 1-254 | Human |

Antibody Type, Clone, Isotype, Control, Localization, Dilution Range

| Tissue Control | Localization | Dilution Range | Application |
|--|---------------------------|----------------|--|
| Breast, Prostate, Skin, Transitional Cell Carcinoma, Papillary Thyroid Carcinoma, Pancreatic Carcinoma | Membranous | 1:25-1:100 | Thyroid and Parathyroid Cancer, Gallbladder and Pancreatic Cancer, Cervical Cancer, Gastric Cancer |
| Liver, Kidney, Tonsil, Uterus, Cervix, Skin, Colon | Cytoplasmic | 1:50-1:200 | Hematopoietic |
| Liver, Kidney, Tonsil, Uterus, Cervix, Skin, Colon | CytoPlasmic | 1:100-1:500 | Hematopoietic |
| Pituitary | Cytoplasmic | 1:250-1:1000 | Pituitary, Neural & Neuroendocrine Cancer |
| Pituitary | Cytoplasmic | 1:250-1:1000 | Pituitary, Neural & Neuroendocrine Cancer |
| Lung, Thyroid, Adenocarcinoma of Lung | Nuclear | 1:250-1:1000 | Lung Cancer, Thyroid & Parathyroid, Mesothelioma, Carcinomas of Unknown Primary Site, Liver Cancer |
| Skin, Malignant Melanoma | Cytoplasmic | 1:50-1:200 | Skin, Malignant Melanoma |
| A synthetic peptide corresponding to residues of human Uroplakin III protein | Nuclear | 1:50-1:200 | Kidney & Urotelial Cancer |
| Bladder, Bladder Carcinoma | Cytoplasmic, Membranous | 1:25-1:100 | Kidney & Urotelial Cancer |
| Varicella Zoster Virus Infected Tissue | Cytoplasmic, Membranous | 1:25-1:100 | Infectious Diseases |
| Placenta, Angioma, AngioSarcoma & Soft Tissue | Cytoplasmic, Cell Surface | 1:10-1:50 | Endothelial, Cervical Cancer, Liver Cancer |
| Kidney, Colon, Small Bowel Mucosa, Colonic Mucosa | Cytoplasmic, Membranous | 1:100-1:500 | Colon & GI Cancer, Liver Cancer, Gall Bladder and Pacreatic cancer, Carcinomas of Unknown Primary Site |
| Kidney, Colon, Small Bowel Mucosa, Colonic Mucosa | Cytoplasmic, Membranous | 1:25-1:100 | Colon & GI Cancer, Liver Cancer, Gall Bladder and Pacreatic cancer, Carcinomas of Unknown Primary Site |
| Tonsil, Lymph Node, Colon | Cytoplasmic | 1:50-1:200 | Endometrial & Genital Cancer, Kidney & Urotelial Cancer, Melanoma & Skin Cancer, Sarcoma & Soft Tissue, Undifferentiated Tumor, Liver Cancer |
| Tonsil, Lymph Node, Colon | Cytoplasmic | 1:250-1:1000 | Endometrial & Genital Cancer, Kidney & Urotelial Cancer, Melanoma & Skin Cancer, Sarcoma & Soft Tissue, Undifferentiated Tumor, Liver Cancer |
| Kidney, Testis, Malignant Mesothelioma | Nuclear | 1:100-1:500 | Lung Cancer, Ovarian Cancer, Kidney & Urotelial Cancer, Mesothelioma, Carcinomas of Unknown Primary Site, Sarcoma and Soft Tissue |
| Placenta, Breast, Fallopian Tube, Testis, Transitional Cell Carcinoma, HER2 Negative Breast Cancer | Nuclear, Cytoplasmic | 1:25-1:100 | Lung Cancer, Head and Neck Cancer, Breast Cancer, Colon and GI Cancer, Liver Cancer, Gallbladder and Pancreatic Cancer, Kidney and Urothelial Cancer, Neural and Neuroendocrine Cancer |
| Tonsil, Lymph Node, Thymus, Chronic Lymphocytic Leukemia | Cytoplasmic | 1:250-1:1000 | Leukemia & Histiocytic, Lymphoma |



www.BioSB.com



BIO-SB SAMPLE MICROARRAY
TISSUE STAIN SLIDE 001

Bio SB
BIO SCIENCE FOR THE WORLD

Immunophenotyping of Cancer by Immunohistochemistry



Carcinomas

| CARCINOMAS 1 | CK Cocktail | CK7 | CK20 | CK, LMW | CK, HMW | CK5 | p63 | Vimentin | TTF-1 (cytoplasm.) | Thyroglobulin |
|---------------------------------|-------------|-----|------|---------|---------|-----|-----|----------|--------------------|---------------|
| Hepatocellular Carcinoma | - | - | - | - | - | - | - | - | + | - |
| Renal Cell Carcinoma | + | - | - | + | - | - | - | + | - | - |
| Bladder Adenocarcinoma | + | + | +/- | + | + | - | - | - | - | - |
| Salivary Gland Carcinoma | + | + | - | + | + | + | + | - | - | - |
| Thyroid Carcinoma | + | + | - | + | - | - | - | + | + | + |
| Breast Carcinoma | + | + | - | + | + | - | - | - | - | + |
| Lung Adenocarcinoma | + | + | - | + | + | - | - | - | + | - |
| Colorectal Adenocarcinoma | + | - | + | + | - | - | - | - | - | - |
| Prostatic Adenocarcinoma | + | - | - | + | - | - | - | - | - | - |
| Transitional Cell Carcinoma | + | + | + | + | + | + | + | - | - | - |
| Ovarian Carcinoma, Non Mucinous | + | + | - | + | + | + | - | - | - | - |
| Cervical Carcinoma | + | + | - | - | - | - | - | - | - | - |
| Sweat Gland Carcinoma | + | + | - | - | - | + | + | - | - | - |
| Pancreatic Carcinoma | + | + | - | + | +/- | - | - | - | - | - |
| Gastric Carcinoma | + | + | - | - | - | - | - | - | - | - |
| Squamous Cell Carcinoma | + | - | - | + | + | + | + | - | - | - |
| Endometrial Adenocarcinoma | + | + | - | + | - | - | - | + | - | - |

| CARCINOMAS 1 | GCDFP-15 | ER/PR | CEA (poly) | CDX-2 | Villin | Hep-Par1 | RCC | CD10 | Beta-Catenin | A-1-Antitrypsin |
|---------------------------------|----------|-------|------------|-------|--------|----------|-----|------|--------------|-----------------|
| Hepatocellular Carcinoma | - | - | + | - | - | + | - | + | - | +/- |
| Renal Cell Carcinoma | - | - | - | - | - | - | + | + | - | + |
| Bladder Adenocarcinoma | - | - | - | + | + | - | - | - | + | - |
| Salivary Gland Carcinoma | - | - | + | - | - | - | - | - | - | - |
| Thyroid Carcinoma | - | - | - | - | - | - | - | - | - | - |
| Breast Carcinoma | + | +/- | - | - | - | - | - | + | - | - |
| Lung Adenocarcinoma | - | - | + | - | - | - | - | - | - | - |
| Colorectal Adenocarcinoma | - | - | + | + | + | - | - | + | + | - |
| Prostatic Adenocarcinoma | - | - | - | - | - | - | - | + | - | - |
| Transitional Cell Carcinoma | - | - | - | - | - | - | - | + | - | - |
| Ovarian Carcinoma, Non Mucinous | - | + | - | - | - | - | - | - | - | - |
| Cervical Carcinoma | - | - | + | - | - | - | - | - | - | - |
| Sweat Gland Carcinoma | - | - | + | - | - | - | - | - | - | - |
| Pancreatic Carcinoma | - | - | + | - | - | - | - | +/- | - | - |
| Gastric Carcinoma | - | - | + | + | + | - | - | - | - | - |
| Squamous Cell Carcinoma | - | - | - | - | - | - | - | - | - | - |
| Endometrial Adenocarcinoma | - | + | - | - | - | - | - | - | + | - |

| CARCINOMAS 2 | CK7 | CK20 | CDX-2 | Villin | TTF-1 | Napsin A | PSA/PSAP | Caveolin-1 | Mammaglobin | GCDFP-15 |
|---------------------------------------|-----|------|-------|--------|-------|----------|----------|------------|-------------|----------|
| Breast Carcinoma | + | - | - | - | - | - | - | + | + | + |
| Lung Adenocarcinoma | + | - | - | - | + | + | - | - | - | - |
| Transitional Cell (Bladder) Carcinoma | + | + | - | - | - | - | - | - | - | - |
| Colorectal Adenocarcinoma | - | + | + | + | - | - | - | + | - | - |
| Prostatic Adenocarcinoma | - | - | - | - | - | - | + | +/- | - | - |

| SQUAMOVS VS. TRANSITIONAL CARCINOMAS | CK, 34bE12 | p63 | CK5 | Thrombomodulin | CK7 | CK20 | Uroplakin III | COX-2 |
|--------------------------------------|------------|-----|-----|----------------|-----|------|---------------|-------|
| Squamous Carcinoma | + | + | + | + | - | - | - | - |
| Transitional Cell Carcinoma | + | + | -/+ | + | + | + | + | + |

Carcinomas

| SMALL CELL VS. MERKEL CELL CARCINOMAS | TTF-1 | CEA | CK20 | Chromogranin A | E-Cadherin (nuclear) | Neurofilament | CD117 | Synaptophysin |
|--|-------|-----|------|----------------|----------------------|---------------|-------|---------------|
| Merkel Cell Carcinoma | - | - | + | + | + | + | + | + |
| Small Cell Carcinoma | + | - | - | - | - | - | +/- | + |

| CUTANEOUS NEOPLASMS | CD10 | Androgen Receptor | CK20 | CD34 | Ber-EP4 | bcl-2 | CK19 |
|-------------------------------|------|-------------------|------|------|---------|-------|------|
| Basal Cell Carcinoma | + | + | - | - | + | + | + |
| Trichoepithelioma | - | - | + | + | + | + | + |
| Merkel Cell Carcinoma | - | - | + | - | + | + | + |
| Microcystic Adnexal Carcinoma | +/- | - | - | - | -/+ | + | + |
| Sebaceous Carcinoma | +/- | + | - | - | + | +/- | - |
| Sebaceous Adenoma | - | + | - | - | + | + | - |

| CARCINOMAS: DIFFERENTIAL DIAGNOSIS | Androgen Receptor | BCA-225 | GDDFP-15 | ER/PR | Mammaglobin | Ny-BR-1 | PSA/PSAP | CD44 |
|---|-------------------|---------|----------|-------|-------------|---------|----------|------|
| Salivary Duct Carcinoma | + | + | + | - | - | - | - | - |
| Breast Carcinoma | +(apocrine) | + | + | +/- | + | + | - | + |
| Prostate Carcinoma (nuclear) | + | - | - | - | - | - | + | + |

| PROSTATE LESIONS | PSA/PSAP | P504s | CK, 34bE12 | p63 | CK7 | Thrombomodulin | Uroplakin III | PAX-2 |
|-------------------------|----------|-------|------------|-----|-----|----------------|---------------|-------|
| Prostate Carcinoma | + | + | - | - | - | - | - | - |
| Urothelial Carcinoma | - | - | + | + | + | + | + | - |
| Nephrogenic Adenoma | - | + | +/- | - | + | - | - | + |

| LUNG CARCINOMAS | CD56 | Chromogranin A | Synaptophysin | CK7 | Napsin A | p63 | TTF-1 | BG8 | GLUT-1 |
|-------------------------------------|------|----------------|---------------|-----|----------|-----|-------|-----|--------|
| Small Cell Carcinoma | + | -/+ | +/- | - | - | - | + | + | - |
| Adenocarcinoma | - | - | - | + | + | - | + | + | +/- |
| Squamous Cell Carcinoma | - | - | - | - | - | + | - | + | + |
| Large Cell Neuroendocrine Carcinoma | + | + | + | +/- | - | - | +/- | - | - |

| COLON VS. OVARIAN CARCINOMA | CK7 | CK20 | CEA | CDX-2 | Villin | CA19-9 | Ep-CAM | WT1 | CA-125 | CK5/6 | PAX-8 |
|------------------------------------|-----|------|-----|-------|--------|--------|--------|-----|--------|-------|-------|
| Ovarian Carcinoma, Serous | + | - | + | - | + | + | + | + | + | - | + |
| Ovarian Carcinoma, Mucinous | +/- | +/- | - | + | + | + | + | - | - | - | - |
| Ovarian Endometrioid Carcinoma | + | - | - | - | - | +/- | + | -/+ | + | - | + |
| Colon Carcinoma | - | + | + | + | + | + | + | - | - | - | - |

Carcinomas

| CARCINOMAS: MUCIN EXPRESSION IN NEOPLASMS | MUC1 | MUC2 | MUC5AC | MUC6 |
|--|------|------|--------|------|
| Pancreatic Adenocarcinoma | + | - | + | - |
| Lung Carcinoma, Signet Ring | + | - | + | - |
| Ovarian Adenocarcinoma | + | - | + | - |
| Gastric Carcinoma | - | -/+ | + | -/+ |
| Cervical Adenocarcinoma | + | - | + | - |
| Esophageal Carcinoma | + | - | + | - |
| Paget's Extramammary | + | -/+ | + | - |
| Cholangiocarcinoma | + | - | +/- | - |
| Salivary Duct Adenocarcinoma | - | + | - | + |
| Colon Carcinoma, Signet Ring | - | + | - | - |
| Prostate Carcinoma | - | +/- | - | - |
| Pancreatic Intraductal Papillary Carcinoma | - | + | + | + |
| HCC | - | - | - | - |
| Adrenocortical Carcinoma | - | - | - | - |
| Breast Carcinoma | + | - | - | - |
| Lung Carcinoma | + | - | - | - |
| Kidney Carcinoma | + | - | - | - |
| Bladder Carcinoma | + | - | - | - |
| Endometrial Carcinoma | + | - | - | - |
| Ovarian Carcinoma, Mucinous | + | - | + | - |
| Renal Cell Carcinoma | + | - | - | - |
| Urothelial Carcinoma | + | - | - | - |
| Appendiceal Adenocarcinoma | - | + | + | - |
| Barrett's Esophagus | + | + | + | - |
| Pancreatic Mucinous Cystic | - | - | + | - |
| Breast Colloid Carcinoma | + | + | - | + |

| MUCIN EXPRESSION IN ORGANS | MUC1 | MUC2 | MUC4 | MUC5AC | MUC6 |
|----------------------------|------|------|------|--------|------|
| Stomach | + | - | + | + | + |
| Small Intestine | - | + | - | - | + |
| Colon | - | + | - | - | - |
| Pancreas | + | - | - | - | + |

| AMPULLARY CARCINOMA (ENTERIC VS DUCTAL) | CK17 | Hep-Par1 | CDX-2 |
|---|------|----------|-------|
| Enteric | - | + | - |
| Ductal | + | - | -/+ |

| COLON VS. PROSTATE ADENOCARCINOMA | CDX-2 | CK20 | CEA | CA19-9 | PSA | P504s |
|-----------------------------------|-------|------|-----|--------|-----|-------|
| Colon Adenocarcinoma | + | + | + | + | - | - |
| Prostate Adenocarcinoma | - | - | - | - | + | + |

Organ Specific

| ORGAN SPECIFIC: PANCREATIC TUMORS | Synaptophysin | Chromogranin_A | Insulin | Glucagon | Gastrin | Somatostatin | MUC4 | CD56 |
|--|---------------|----------------|---------|----------|---------|--------------|------|------|
| Pancreatic Adenocarcinoma | - | - | - | - | - | - | + | - |
| Neuroendocrine Tumor | + | + | +/- | +/- | +/- | +/- | - | + |
| Solid Pseudopapillary | + | - | - | - | - | - | - | + |
| Pancreatic Ductal Carcinoma | - | - | - | - | - | - | - | - |
| Acinic Cell Carcinoma | - | - | - | - | - | - | - | - |
| Pancreatoblastoma | - | + | - | - | - | - | - | + |
| Benign Pancreas | + | + | + | + | - | + | - | - |

| ORGAN SPECIFIC: PANCREATIC TUMORS | B-Catenin | PGP 9.5 | CK19 | CA19-9 | bcl-10 | E-Cadherin | CD10 | Maspin |
|--|-----------|---------|------|--------|--------|------------|------|--------|
| Pancreatic Adenocarcinoma | - | - | + | + | - | - | +/- | + |
| Neuroendocrine Tumor | + | + | +/- | +/- | - | - | - | - |
| Solid Pseudopapillary | + | - | - | - | - | - | + | - |
| Pancreatic Ductal Carcinoma | +/- | - | - | + | - | +/- | +/- | + |
| Acinic Cell Carcinoma | + | - | + | -/+ | + | + | +/- | - |
| Pancreatoblastoma | + | - | - | - | + | - | - | + |
| Benign Pancreas | + | - | - | - | - | - | - | - |

| ORGAN SPECIFIC: BRAIN CNS TUMORS | GFAP | Neurofilament | Synaptophysin | S-100 | CK Cocktail | PR | EMA | Vimentin | NGFR | INI-1 | SOX2 |
|---|------|---------------|---------------|-------|-------------|-----|-----|----------|------|-------|------|
| Astrocytoma | +/- | - | - | + | - | - | - | + | + | + | + |
| Glioblastoma | + | - | - | + | - | - | - | + | - | + | + |
| Oligodendroglioma | - | - | - | + | - | - | - | + | - | + | + |
| Ependymoma | +/- | - | - | + | - | - | - | -/+ | + | + | + |
| Choroid Plexus Carcinoma | -/+ | - | + | + | + | - | - | - | - | + | - |
| Central Neurocytoma | -/+ | - | +/- | - | - | - | - | - | + | + | + |
| Neuroblastoma | +/- | + | +/- | +/- | - | - | - | - | + | + | - |
| Pineocytoma | -/+ | - | + | + | - | - | - | - | - | + | - |
| Meningioma | - | - | - | - | - | + | + | + | - | + | - |
| Schwannoma | + | - | - | + | -/+ | - | - | + | + | - | - |
| Rhabdoid Tumors | - | +/- | - | +/- | + | - | + | + | - | - | - |
| Metastatic Carcinoma | - | - | - | - | + | -/+ | + | -/+ | - | - | - |

| ORGAN SPECIFIC: KIDNEY RENAL EPITHELIAL TUMORS | RCC | CD10 | PAX-2 | Vimentin | Ksp-Cadherin | Parvalbumin | CD117 | Ep-CAM | Caveolin-1 | PAX-8 | pVHL |
|---|-----|------|-------|----------|--------------|-------------|-------|--------|------------|-------|------|
| Clear Cell Renal Cell Carcinoma | + | + | + | + | - | - | - | - | + | + | + |
| Chromophobe RCC | - | - | - | - | + | + | + | + | + | + | + |
| Oncocytoma | - | - | - | - | + | + | + | - | +/- | + | + |

| ORGAN SPECIFIC: BREAST CARCINOMAS | CK7 | CK20 | ER/PR | CD44 | CA15-3 | CA19-9 | p63 | CK5 | CD117 |
|--|-----|------|-------|------|--------|--------|-----|-----|-------|
| Infiltrating Ductal Carcinoma | + | - | + | + | + | - | - | - | - |
| Adenoid Cystic Carcinoma | + | - | - | - | + | + | + | + | + |

| ORGAN SPECIFIC: BREAST VS. LUNG VS. PROSTATE CARCINOMAS | GCDFP-15 | Mammaglobin | CA15-3 | Caveolin-1 | PSA | TTF-1 | Napsin A |
|--|----------|-------------|--------|------------|-----|-------|----------|
| Breast Carcinoma | + | + | + | + | - | - | - |
| Lung Carcinoma | - | - | - | - | - | + | + |
| Prostate Carcinoma | - | - | - | +/- | + | - | - |

Organ Specific

| ORGAN SPECIFIC: DIFFERENTIAL DIAGNOSIS OF PARATHYROID TUMORS | Chromogranin A | Synaptophysin | PTH | S-100 | TTF-1 | Calcitonin | PAX-8 |
|---|----------------|---------------|-----|-------|-------|------------|-------|
| Parathyroid Tumors | + | + | + | - | - | - | - |
| Follicular Cell Tumors | - | - | - | +/- | + | - | + |
| Medullary Thyroid Carcinoma | + | + | - | - | + | + | + |

| ORGAN SPECIFIC: SEX CORD STROMAL TUMORS | Calretinin | Inhibin | CD99 | CK7 | EMA | Vimentin | MART-1 |
|--|------------|---------|------|-----|-----|----------|--------|
| Granulosa Cell Tumors | + | + | + | - | - | + | + |
| Sertoli-Leydig Cell Tumors | + | + | -/+ | + | - | + | + |
| Gynandroblastoma | + | + | | | | | |
| Gonadoblastomas | + | + | + | - | - | + | - |

| ORGAN SPECIFIC: UTERUS: TROPHOBLASTIC CELLS | 1st Trimester | | 2nd Trimester | | 3rd Trimester | |
|--|---------------|--------|---------------|--------|---------------|-------|
| | hCG | hPL | hCG | hPL | hCG | hPL |
| Cytotrophoblast | - | - | - | - | - | - |
| Intermediate Trophoblast | 1-24% | 25-49% | -/+ | 50-74% | 1-24% | 1-49% |
| Syncytiotrophoblast | >75% | 1-24% | 25-49% | 50-74% | 1-24% | >75% |

| ORGAN SPECIFIC: UTERUS : TROPHOBLASTIC PROLIFERATIONS | p57 | hCG | PLAP | hPL | CK Cocktail | Vimentin |
|--|-----|-----------------|-----------------|-----------------|-----------------|-----------------|
| Partial Mole | + | weak, diffuse | + | weak, diffuse | strong, diffuse | - |
| Complete Mole | - | strong, diffuse | weak, focal | weak, focal | strong, diffuse | - |
| Choriocarcinoma | - | strong, diffuse | weak, focal | weak, focal | strong, diffuse | - |
| Placental Site Tumor | | strong, focal | strong, diffuse | strong, diffuse | strong, diffuse | strong, diffuse |

| ORGAN SPECIFIC: SKIN: ADNEXAL TUMORS | CK7 | CK20 | S-100 | EMA | GCDFP-15 | CD15 |
|---|-----|------|-------|-----|----------|------|
| Merkel Cell Carcinoma | - | + | - | + | - | - |
| Sebaceous Tumor | + | - | - | - | - | + |
| Apocrine Tumor | | - | - | +/- | + | +/- |
| Eccrine Tumor | | - | + | + | - | - |

| ORGAN SPECIFIC: BREAST LESION | GCDFP-15 | Mammaglobin | B-Catenin | E-Cadherin | CK, 34bE12 | p120 |
|--------------------------------------|----------|-------------|---------------|------------|------------|----------------|
| Lobular | + | + | - | - | + | +(cytoplasmic) |
| Ductal | + | + | +(membranous) | + | - | +(membranous) |

| ORGAN SPECIFIC: MENINGIOMAS FROM HISTOLOGIC MIMICS | Claudin 1 | EMA | S-100 | CD34 | GFAP |
|---|-----------|-----|-------|------|------|
| Meningothelial Meningioma | + | + | - | - | - |
| Atypical Meningioma | + | + | - | + | - |
| Fibrous Meningioma | - | + | + | - | - |
| Solitary Fibrous Tumor | - | - | - | + | - |
| Meningeal Hemangiopericytoma | - | - | - | + | - |
| Schwannoma | +/- | - | + | - | + |

| ORGAN SPECIFIC: RENAL CELL CARCINOMA VS. HEMANGIOBLASTOMA | D2-40 | FLI-1 | CD31 | CK Cocktail | CD10 |
|--|-------|-------|------|-------------|------|
| Metastatic Renal Cell Carcinoma | - | - | - | + | + |
| Hemangioblastoma | + | + | + | - | - |

Organ Specific

| ORGAN SPECIFIC: OVARIAN CARCINOMAS | PAX-8 | WT1 | CA-125 | CEA | pVHL |
|---|--------------------|-----------------|------------------|--------------|-----------------|
| Ovarian Carcinoma, Serous | + | + | + | + | - |
| Ovarian Carcinoma, Mucinous | - | - | - | - | - |
| Ovarian Endometrioid Carcinoma | + | - | + | - | - |
| Ovarian Clear Cell Carcinoma | + | - | + | - | + |
| ORGAN SPECIFIC: SKIN: PAGETOID TUMORS | CK, LMW | CK, HMW | S-100 | CEA | Vimentin |
| Melanoma | - | - | + | - | + |
| Paget's Disease | + | - | -/+ | + | - |
| Bowen's Disease | + | + | - | - | - |
| ORGAN SPECIFIC: SKIN: BASAL VS. SQUAMOUS CELL CARCINOMA | CK Cocktail | Ep-CAM | EMA | bcl-2 | |
| Basal Cell Carcinoma | + | + | - | + | |
| Squamous Cell Carcinoma | + | - | + | - | |
| ORGAN SPECIFIC: PERINEURIOMA VS. NEUROFIBROMA | Claudin 1 | EMA | S-100 | GLUT1 | |
| Perineurioma | + | + | - | + | |
| Neurofibroma | + | + | - | - | |
| ORGAN SPECIFIC: BREAST CARCINOMA IN SITU VS. INFILTRATING BREAST CARCINOMA | | Calponin | SM Myosin | p63 | |
| Breast Carcinoma in-situ | | + | + | + | |
| Infiltrating Breast Carcinoma | | - | - | - | |
| ORGAN SPECIFIC: PANCREAS VS. PANCREATIC ADENOCARCINOMA | pVHL | S-100 | Maspin | | |
| Pancreas | + | - | - | | |
| Pancreatic Intraepithelial Neoplasia | - | + | + | | |
| Pancreatic Adenocarcinoma | - | + | + | | |

Lymphomas

| LYMPHOMAS: B-CELL LYMPHOMAS | CD45 | CD20 | CD79a | bcl-2 | bcl-6 | CD10 | CD23 | Cyclin D1 | PAX-5 | BOB.1 | Oct-2 | PU.1 |
|--|------|------|-------|-------|-------|------|------|-----------|-------|-------|-------|------|
| Follicular | + | + | + | + | + | + | - | - | + | + | + | + |
| CLL/SLL | + | + | + | + | - | - | + | - | + | + | + | + |
| Mantle Cell | + | + | + | + | - | - | - | + | + | + | + | + |
| Marginal Zone BCL | + | + | + | + | - | - | - | - | + | + | + | + |
| Lymphoplasmacytic | + | + | + | + | - | - | - | - | + | + | | |
| Diffuse Large Cell Lymphoma | + | + | + | + | + | - | - | - | + | + | + | + |
| Burkitt Lymphoma | + | + | + | - | + | + | - | - | + | | | |
| Hairy Cell Leukemia | + | + | + | + | - | - | - | - | + | | | |

| LYMPHOMAS: B-CELL LYMPHOMAS | p27/Kip1 | IgD | MUM1 | T-beta | TRAcP | Annexin A1 | CD5 | CD43 | CD3 | ZAP-70 | CD25 | FOXP1 | TCL1 |
|--|----------|-----|------|--------|-------|------------|-----|------|-----|--------|------|-------|------|
| Follicular | + | + | - | - | - | - | - | - | - | - | - | - | + |
| CLL/SLL | + | + | + | + | - | - | + | + | - | +/- | | | + |
| Mantle Cell | + | + | - | - | - | - | + | + | - | - | - | -/+ | + |
| Marginal Zone BCL | + | -/+ | + | + | +/- | - | - | - | - | - | - | | - |
| Lymphoplasmacytic | + | - | + | + | - | - | - | - | - | - | - | | + |
| Diffuse Large Cell Lymphoma | - | - | + | - | - | - | - | - | - | - | | + | -/+ |
| Burkitt Lymphoma | - | - | - | - | - | - | - | + | - | - | | | + |
| Hairy Cell Leukemia | - | - | | + | + | + | - | - | - | | + | | |

| LYMPHOMAS: T-CELL LYMPHOMAS | CD45 | CD2 | CD3 | CD4 | CD5 | CD7 | CD8 | CD25 | CD45RO | PD-1 | CD56 | CD57 | Perforin | Granzyme_B | TCL1 |
|--|------|-----|-----|-----|-----|-----|-----|------|--------|------|------|------|----------|------------|------|
| Angioblastic | + | + | + | + | + | + | +/- | + | + | + | | | | | |
| Lymphoblastic | + | +/- | + | +/- | + | + | +/- | + | + | - | | | | | - |
| Subcutaneous Panniculitic | + | + | + | - | + | + | +/- | - | + | - | - | | + | + | |
| Lennert's | + | + | + | - | + | | - | + | + | - | | | | - | |
| NK-Type | + | + | + | - | - | + | - | + | + | - | +/- | +/- | - | + | + |
| Cutaneous | + | + | + | + | - | + | - | - | - | - | | | + | + | - |
| Peripheral | + | + | + | + | + | - | - | + | + | - | - | - | - | - | - |
| Mycosis Fungoides | + | + | + | + | + | - | - | + | + | - | - | | | | - |

| LYMPHOMAS: ACUTE MYELOID LEUKEMIA | MPO | CD68 | Factor VIII | CD61A | Hemoglobin A | BOB.1 | Oct-2 | Glycophorin A | Spectrin | CD34 | CD43 | CD74 | CD45 | Lysozyme | CD138 |
|--|-----|------|-------------|-------|--------------|-------|-------|---------------|----------|------|------|------|------|----------|-------|
| Acute Myeloid, M0 | - | - | - | - | - | - | - | - | - | + | + | + | + | + | + |
| Myeloblast, M1&2 | + | + | - | - | - | - | - | - | - | - | + | + | + | + | + |
| Promyelocytic, M3 | + | + | - | - | - | + | + | | | - | + | | - | | |
| Myelomonocytic, M4 | + | + | - | - | - | - | + | - | - | + | + | + | + | + | |
| Monoblastic, M5 | + | + | - | - | - | - | + | - | - | - | + | + | + | + | |
| Acute Myeloid, M6 | + | - | - | | - | - | - | + | + | - | | | - | | + |
| Megakaryocytic, M7 | - | - | + | + | - | +/- | - | | | - | | | | | - |
| Megakaryoblast | + | + | + | | + | | | | | | | | | | - |

| LYMPHOMAS: LYMPHOBLASTIC LYMPHOMAS: BCL VS. TCL | TDT | CD10 | PAX-5 | CD20 | CD19 | CD3 | CD5 | CD7 | CD117 | CD74 |
|--|-----|------|-------|------|------|-----|-----|-----|-------|------|
| Lymphoblastic BCL | + | +/- | + | +/- | + | - | - | - | - | + |
| Lymphoblastic TCL | + | + | - | - | - | + | +/- | + | - | - |

| LYMPHOMAS: PLASMA CELLS | CD138 | CD79a | EMA | MUM1 | CD56 | Cyclin D1 | CD43 | CD20 | CD19 |
|--------------------------------|-------|-------|-----|------|------|-----------|------|------|------|
| Plasma Cell Tumor | + | + | + | + | + | - | - | - | - |

Lymphomas

| LYMPHOMAS: LYMPH NODE HISTIOCYTOSIS | CD68 | S-100 | CD1a | Lysozyme | CD21/CD35 | FDC | PD-1 |
|---|-------------|--------------|-------------|-----------------|------------------|------------|-------------|
| Reactive | + | - | - | + | - | - | - |
| Langerhans | + | + | + | + | - | - | - |
| Sinus Histiocytosis with Massive Lymphadenopathy | + | + | - | + | - | - | - |
| Follicular Dendritic Cell Tumor | - | - | - | - | + | + | + |
| Dermatopathic Lymphadenitis | - | + | + | + | - | - | - |

| LYMPHOMAS: HISTIOCYTIC LESIONS | CD45 | CD4 | CD68 | Lysozyme | CD163 | Factor XIIIa | CD20 | CD3 |
|---|-------------|------------|-------------|-----------------|--------------|---------------------|-------------|------------|
| Histiocytic Lesions | + | + | + | + | + | + | - | - |

| LYMPHOMAS: HISTIOCYTIC PROLIFERATION | S-100 | CD68 | Vimentin | Lysozyme | CD1a | Factor XIIIa | HAM56 |
|---|--------------|-------------|-----------------|-----------------|-------------|---------------------|--------------|
| Juvenile Xanthogranuloma | - | + | + | + | - | + | + |
| Langerhans Cell Histiocytosis | + | + | + | + | + | - | + |
| Dermatobroma | - | + | + | - | - | + | + |

| LYMPHOMAS: IMMUNOGLOBULIN: HEAVY AND LIGHT CHAIN | IgA | IgG | IgD | IgM | Kappa | Lambda |
|---|------------|------------|------------|------------|--------------|---------------|
| Secretory Meningioma | + | - | - | + | | |
| Microvillous Lymphoma | - | - | - | + | | |
| Cutaneous Lymphoma | - | - | - | - | | |
| Myeloma | + | + | +/- | - | + | + |
| Diffuse LBCL | - | + | - | + | + | + |
| Marginal Zone Lymphoma | - | - | -/+ | + | + | + |
| SLL/CLL | - | - | + | + | + | + |

| LYMPHOMAS: B-CELL MATURATION | CD20 | CD79a | CD19 |
|---|-------------|--------------|-------------|
| Pro-B | - | + | + |
| Pre-B | - | + | + |
| Late Pre-B | + | + | + |
| Mantle B-cells | + | + | + |
| Centroblastic | + | + | + |
| Centrocytic | + | + | + |
| Plasmablasts | - | + | - |
| Plasma cells | - | + | - |

| LYMPHOMAS: HIGH GRADE LYMPHOMA: BCL VS. TCL | CD20 | CD45R | CD43 | CD45R0 | CD3 |
|--|-------------|--------------|-------------|---------------|------------|
| Low Grade BCL | + | + | - | - | - |
| High Grade BCL | + | + | - | - | - |
| Low Grade TCL | - | - | + | + | + |
| High Grade TCL | - | - | + | + | +/- |


| LYMPHOMAS: FOLLICULAR LYMPHOMA VS. REACTIVE FOLLICLES | bcl-2 | bcl-6 | CD10 | Ki-67 |
|--|--------------|--------------|-------------|--------------|
| Follicular Lymphoma | + | + | + | - |
| Reactive Follicles | - | + | + | + |

Lymphomas

| LYMPHOMAS: MASTOCYTOSIS | Tryptase | CD117 | CD25 | CD163 | CD2 |
|--------------------------------|-----------------|--------------|-------------|--------------|------------|
| Mastocytosis | + | + | + | - | + |
| Reactive Mast Cells | + | + | - | + | - |
| Myelomastocytic Leukemia | + | + | - | | - |

| LYMPHOMAS: SPLENIC HEMATOPOIETIC PROLIFERATIONS IN NEOPLASTIC AND BENIGN DISORDERS | MPO | CD34 | CD117 | CD68 | Hemoglobin A |
|---|------------|-------------|--------------|-------------|---------------------|
| Chronic Myelogenous Leukemia | + | -/+ | +/- | + | - |
| Chronic Idiopathic Myelofibrosis | + | +/- | -/+ | | - |
| Myelodysplastic Syndrome | | + | -/+ | | - |
| Myelodysplastic / Myeloproliferative Disorders | + | - | - | + | - |
| Mastocytosis | + | - | + | | - |
| Erythroid Disorders | +/- | - | - | -/+ | + |
| Splenic Lymphoma | -/+ | - | - | | - |
| Acute Myeloid Leukemia | + | + | + | + | - |
| Polycythemia Vera | | + | + | | + |

ZytoLight Interpretation Guides



HER-2 neu/ERBB2 Interpretation Guide



HER-2 neu/ERBB2 Interpretation Guide







Zytolight SPEC ALK/EML4 TrichCheck Probe - Recommended Evaluation Scheme

Initial Screening using a Green/Orange Dual Bandpass Filter Set



Secondary Screening using a Blue Sing Bandpass Filter Set



Signal patterns other than those shown above may be observed. It is advised to analyze metaphase spreads for the interpretation of those atypical signal patterns.

Also Observable:



HER-2 neu/ERBB2 FISH Signal Interpretation Guide

HER-2 neu/ERBB2 non-amplified Cell



Count: 2 green and 2 orange signals

HER-2 neu/ERBB2 non-amplified Cell




Count: 2 green and 2 orange signals
One green signal is split but 2 signals of the same color separated by a distance of < 1 Signal diameter, are counted as one

Cell with low level of amplification
HER-2 neu/ERBB2



Count: 7 green and 2 orange signals

Cell with high level of amplification
HER-2 neu/ERBB2



Green signals overlapping orange signals. Signal cluster overlapping signal. Check signals in single bandpass filter.

Cell with monosomy of
chromosome 17



Count: 3 green and 1 orange signal


Cell with monosomy of
chromosome 17



Count: 5 green and 5 orange signals


- Artifacts (crush or edge artifacts) that make interpretation difficult should be excluded from counting.
- Do not count if controls are not as expected.
- If > 25% of signals are weak the test cannot be scored.
- The test should be repeated if >10% of signals occur every cytoplasm.

Over-digested cell



Over-digestion can be recognized by dark areas visible inside of the nuclei. Over-digested nuclei - Do not count!

Cell with autofluorescence



Strong autofluorescence hinders signal recognition. Autofluorescence - Do not count!

Selected Lung Cancer Probes of the ZytoLight Portfolio

| Prod. No. | Product | Label | Tests (Volume)* |
|-------------------------|---|-----------|---------------------|
| Z-2124 | ZytoLight SPEC ALK Dual Color Break Apart Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2117 | ZytoLight SPEC ALK/EML4 TriCheck Probe | ● / ● / ● | 5/20 (50 ul/200 ul) |
| Z-2137 | ZytoLight SPEC CARS Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2179 | ZytoLight SPEC CD274, PDCD1LG2/CEN 9 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2033 | ZytoLight SPEC EGFR/CEN 7 Dual Color Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2136 | ZytoLight SPEC EML4 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2072 | ZytoLight SPEC FGFR1/CEN 8 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2122 | ZytoLight SPEC FGFR2/CEN 10 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2082 | ZytoLight SPEC FGFR3/ CEN 4 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2131 | ZytoLight SPEC KIF5B Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2115 | ZytoLight SPEC KRAS/CEN 12 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2087 | ZytoLight SPEC MET/CEN 7 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2167 | ZytoLight SPEC NTRK1 Dual Color Break Apart Probe | ● / ● | 20 (200 ul) |
| Z-2148 | ZytoLight SPEC RET Dual Color Break Apart Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2144 | ZytoLight SPEC ROS1 Dual Color Break Apart Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2127 | ZytoLight SPEC SOX2/CEN 3 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2133 | ZytoLight SPEC TFG Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Related Products | | | |
| Z-2028-5 | ZytoLight FISH-Tissue Implementation Kit Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 150 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml | | 5 |
| Z-2028-20 | ZytoLight FISH-Tissue Implementation Kit Incl. Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 500 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml | | 20 |
| Z-2099-20 | ZytoLight FISH-Cytology Implementation Kit Incl. Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/DuraTect-Solution, 0.8 ml | | 20 |

ZytoLight Fluorochromes

Two Factors that mainly influence FISH analyses:

- **Fluorochromes of the FISH probes**
- **Appropriate filter sets**

| Fluorochrome | Excitation | Emission | Equivalent to |
|--------------|------------|----------|---------------|
| ● ZyBlue | 418 nm | 467 nm | DEAC |
| ● ZyGreen | 503 nm | 528 nm | FITC |
| ● ZyOrange | 547 nm | 572 nm | Rhodamine |

Recommended Filter Sets

All Filter sets have a superior signal-to-noise ratio and need to be assembled in fluorescence filter holders specific for the respective microscope.

| Prod. No. | Product | Detected Fluorochrome |
|-----------|---|-----------------------|
| E-4030-1 | DAPI Single Bandpass Filter Set v2 | DAPI |
| E-4026-1 | ZyBlue Single Bandpass Filter Set v2 | ● |
| E-4012-1 | ZyGreen Single Bandpass Filter Set v2 | ● |
| E-4013-1 | ZyOrange Single Bandpass Filter Set v2 | ● |
| E-4016-1 | ZyGreen/ZyOrange Dual Bandpass Filter Set v2 | ● / ● |
| E-4010-1 | DAPI/ZyGreen/ZyOrange Triple Bandpass Filter Set v2 | DAPI / ● / ● |
| E-4025-1 | ZyBlue/ZyGreen/ZyOrange Triple Bandpass Filter Set v2 | ● / ● / ● |

Selected Sarcoma Probes of the ZytoLight Portfolio

| Prod. No. | Product | Label | Tests (Volume)* |
|-------------------------|---|-------|---------------------|
| Z-2103 | ZytoLight SPEC CDK4/CEN 12 Dual Color Probe | ● / ● | 20 (200 ul) |
| Z-2121 | ZytoLight SPEC COL1A1 Dual Color Break Apart Probe | ● / ● | 20 (200 ul) |
| Z-2116 | ZytoLight SPEC COL1A1/PDGFB Dual Color Dual Fusion Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2100 | ZytoLight SPEC DDIT3 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2096 | ZytoLight SPEC EWSR1 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2139 | ZytoLight SPEC FOXO1 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2018 | ZytoLight SPEC FOXO1/PAX3 Dual Color Single Fusion Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2019 | ZytoLight SPEC FOXO1/PAX7 Dual Color Single Fusion Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2130 | ZytoLight SPEC FUS Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2132 | ZytoLight SPEC JAZF1 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2013 | ZytoLight SPEC MDM2/CEN 12 Dual Color Probe | ● / ● | 5/20 (50 ul/200 ul) |
| Z-2145 | ZytoLight SPEC NR4A3 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2119 | ZytoLight SPEC PDGFB Dual Color Break Apart Probe | ● / ● | 20 (200 ul) |
| Z-2178 | ZytoLight SPEC SMARCB1/22q12 Dual Color Probe | ● / ● | 5 (50 ul) |
| Z-2097 | ZytoLight SPEC SS18 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2109 | ZytoLight SPEC TFE3 Dual Color Break Apart Probe | ● / ● | 20 (200 ul) |
| Z-2151 | ZytoLight SPEC USP6 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2142 | ZytoLight SPEC WT1 Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Z-2175 | ZytoLight SPEC YWHAE Dual Color Break Apart Probe | ● / ● | 5 (50 ul) |
| Related Products | | | |
| Z-2028-5 | ZytoLight FISH-Tissue Implementation Kit Incl. Heat Pretreatment Solution Citric, 150 ml; Pepsin Solution, 1 ml; Wash Buffer SSC, 150 ml; 25x Wash Buffer A, 50 ml; DAPI/DuraTect-Solution, 0.2 ml | | 5 |
| Z-2028-20 | ZytoLight FISH-Tissue Implementation Kit Incl. Heat Pretreatment Solution Citric, 500 ml; Pepsin Solution, 4 ml; Wash Buffer SSC, 500 ml; 25x Wash Buffer A, 100 ml; DAPI/DuraTect-Solution, 0.8 ml | | 20 |
| Z-2099-20 | ZytoLight FISH-Cytology Implementation Kit Incl. Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/DuraTect-Solution, 0.8 ml | | 20 |

ZytoLight Fluorochromes

Two Factors that mainly influence FISH analyses:

- **Fluorochromes of the FISH probes**
- **Appropriate filter sets**

| Fluorochrome | Excitation | Emission | Equivalent to |
|--------------|------------|----------|---------------|
| ● ZyGreen | 503 nm | 528 nm | FITC |
| ● ZyOrange | 547 nm | 572 nm | Rhodamine |

Recommended Filter Sets

All Filter sets have a superior signal-to-noise ratio and need to be assembled in fluorescence filter holders specific for the respective microscope.

| Prod. No | Product | Detected Fluorochrome |
|----------|---|-----------------------|
| E-4030-1 | DAPI Single Bandpass Filter Set v2 | DAPI |
| E-4012-1 | ZyGreen Single Bandpass Filter Set v2 | ● |
| E-4013-1 | ZyOrange Single Bandpass Filter Set v2 | ● |
| E-4016-1 | ZyGreen/ZyOrange Dual Bandpass Filter Set v2 | ● / ● |
| E-4010-1 | DAPI/ZyGreen/ZyOrange Triple Bandpass Filter Set v2 | DAPI / ● / ● |

Ordering Information

BIO SB, INC.

Please call our Customer Service department from 9:00 A.M. to 5:00 P.M. (PST), Monday through Friday, to place an order or to enquire about an existing order.

Telephone: +1-805-692-2768
Fax: +1-805-692-2769
Mail Orders: Bio SB, Inc.
Attention to: Customer Service
5385 Hollister Avenue
Building 8, #108
Santa Barbara, CA 93111
U.S.A

To expedite the order process, please include the following information on your purchase order or correspondence:

- Purchase Order Number
- Name and phone person ordering
- Shipping address (please do not use P.O. Box number)
- Billing address (if different from above)
- Name of product, catalog number and quantity
- Special shipping instructions
- Credit card number and expiration date (for credit card payments)

INTERNATIONAL ORDERS

Please visit <http://www.biosb.com/disributors.html> or E-mail us at sales@biosb.com

NEW AND EXISTING BIO SB CUSTOMER ACCOUNT:

First time orders paid by credit card (see Payment section) will be processed and shipped immediately.
For other payment methods please accept a delivery time of up to five business days for credit verification purposes.
Please note that for ll new accounts, prepayment is required.

CREDIT TERMS

Net 30 days in U.S. Dollars, upon approval. Overdue accounts are subject to a finance charge of 1.5% per month (18% per annum).

CONFIRMING ORDERS

To avoid duplication of your shipment, please mark boldly "confirming order - please do not ship" on your order.

PRICING

All prices are quoted in U.S. dollars, exclusive of state and county sales tax, where applicable. Prices are valid only for shipments within the U.S. and are subject to change without notice. Please inquire about our standing order and quantity discount policies.

SHIPPING

Shipping and handling charges are prepaid and added to the invoice. They vary with the destination, weight and content, and are available upon request at order entry and are indicated on the invoice. Reagent orders received by 2:00 P.M. (PST), Monday through Thursday, will generally be Expedited Shipping for Next Day Delivery. Early A.M. and Saturday deliveries are available upon request.

PAYMENT

All payments must be made in U.S. dollars. We accept MasterCard, VISA, American Express, Discovery and bank transfer. Please note that we do not take responsibility for charges or delays incurred at receiver's end as this can vary depending on the rules of customs in each country.

RETURN POLICY

Reagents are covered by the following Total Quality Assurance policy which states:

If you are not completely satisfied with the quality of our reagents, you may return them to us for a refund, replacement or credit memo, at our option. Bio SB's liability is limited to a refund or replacement, at our option. Please obtain a Return Material Authorization (RMA) number from Customer Service prior to the return of a product. Returns, which are not caused by unsatisfactory product performance, must be made within 30 days of delivery and will be subject to a 30% processing fee. Returns or replacements cannot be accommodated for expired products. All products sent without a RMA number will be returned to sender.

TROUBLESHOOTING

Should you have technical difficulties, please submit a Product troubleshooting and complaint form online using the at the following URL:
<https://www.biosb.com/product-troubleshootingcomplaint-form/>.

COMPLAINTS

Complaints need to be sent with a Prodciut Troubleshooting and Complaint (Form 00020) per FDA/ISO Regulations. To recieve form email techncial@biosb.com. Returns will be processed with a returned authorization form.

Bio SB, Inc. Distributors

ARGENTINA



Bioars, S. A.
Estomba 961
Buenos Aires
Argentina C1427COU
Tel: +5 (11) 4555 4601
Website: bioars.com.ar

AUSTRALIA



Diagnostic Technology Pty Ltd.
Suite 45, 7 Narabang Way
Belrose NSW 2085 – Australia
Tel: +61-2-9986-2011
Fax: +61-2-9986-2022
Website: diagnostictechnology.com.au

AUSTRIA



THP Medical Products GmbH
Shuttleworthstr. 19
Vienna
Austria 1210
Tel.: + 43 1 292 8280
Website: www.thp.at

Biozol Diagnostica GmbH

Leipziger Str. 4,
Eching, Germany 85386
Tel.: +49 89 37996666
Website: www.thp.at

BANGLADESH



Daffodil Medi World
84/3 Galaxy Habul Tower,
Ahmed Nager, Mirpur-1
Dhaka, 1216
Bangladesh
Tel.: +88-02-5805-1249
Website: www.daffodilmediworld.com/

BELARUS



Avako, LLC
Timiryazeva str., 65B, 109-3
Minsk 220035
Belarus
Tel.: +33-6-8313-4101
Website: http://www.avakollc.com/

BOLIVIA



Laboratorio Catedral SRL
Andres Ibanez 115,
Santa Cruz de la Sierra, Bolivia
Tel.: +591 3 3333305
Website: https://www.laboratoriocatedral.com/

BELGIUM



Diagomics
7 Av. Didier Daurat BP30044

31702 Blagnac Cedex – France
Tel: +33-6-8313-4101
Website: diagomics.com

BOSNIA AND HERZEGOVINA



EMGEN D.O.O.
Francuske Revolucije bb
Sarajevo
Bosnia and Herzegovina, 71000
Tel: +387-6-675-7059
Website: http://www.emgen.ba

BOLIVIA



Laboratorio Catedral SRL
Andres Ibanex # 115
Santa Cruz de la sierra, Santa Cruz
Bolivia
Tel: +591-7049-9988
Website: laboratoriocatedral.com

BRAZIL



Nova Analitica Imp. Exp. Ltda
Rua Assungui, 431
Sao Paulo, 04131-000
Brazil
Tel: +55 11 21628080
Website: analiticaweb.com.br

CANADA



FroggaBio Inc.
50 Great Gulf Drive, Unit #41
Concord, Ontario
Canada L45K 5W1
Tel: +1-877-318-7277
Website: http://www.froggabio.com/

CHILE



Polycompany International Ltda.
Poeta Vicente Huidobro #4041
Macul – Santiago – Chile
Tel: +56-2-221-7760
Fax: +56-2-221-8255

Imerlab Importaciones

Segunda Transversal #6164
San Miguel – Santiago – Chile
Tel: +56-2-359-0043, 312-3855 & 9-231-1714
Fax: +56-2-359-0043
Website: imerlab.cl

CHINA



Biohub International Trade Co., Ltd.
Chuansha Rd. # 6619, Pudong
Shanghai, 202100
China
Tel: +86-021-5272-4187
http://www.qfbio.com/

Wason Biotech Inc.

Rm. 402, Bldg 6, Liando U Vally West,
Bo. 10 Jingsheng South Second Street,
Jinqiao Sci & Tech Industrial Zone, Tongzhou District
Beijing 101102
China
+86-10-5637-0298 ext. 8005
http://www.wasonbiotech.com/

COLOMBIA



Immunopat S.A.S.
Calle 50 num. 43-91 Ofc. 1001
Medellin, Antioquia
Colombia, 1506
Tel: +573-17-500-6887

COSTA RICA



Grupo Labiopath S.A.
Villa Flores, 750 mts Oeste de la Bomba del Barrio
San Jose
Cuidados de Alejuela, Costa Rica
Centro América
Tel: +506 2431-3953
Fax: +506 8993-6751
Website: labiopath.com

CYPRUS



C&V Kriticos Suppliers Ltd.
27a Kanaris STR
Nicosia
Cyprus 1080
Tel: +357-22-33-8234

CZECH REPUBLIC



Pragostem
17. Listopadu 198
Pardubice
Czech Republic
Tel: +420-727 824 243

DENMARK



Nuppulinnan Laboratoriopalvelu Oy
Sallintie 1C
Kellokoski, Uusimaa
Finland 04500
Tel: +420-727 824 243
Website: pragostem.cz

ECUADOR



Acces-Lab
Mariano Jimbo N40-90 y Gaspar de Villaroel
Quito – Ecuador
Tel: +593-245-6500, 243-7131 & 243-7139

EGYPT



SNF Medical
54 (B), Bashaer Elsalam Bldgs,
Mariota, Fasal St.

Giza, Egypt
Tel: +20237454550
+201228167210

EL SALVADOR



Rowalt Pharmaceutical S.A. de C.V.
Avenida Libertad No. 419
Colonia Libertad
San Salvador – El Salvador, C.A.
Tel: +503-225-8974
Contact: Walter Ernesto Rodríguez Rivas

ESTONIA



Nuppulinnan Laboratoriopalvelu Oy
Sallintie 1C
Kellokoski, Uusimaa
Finland 04500
Tel: +358-40-5305367
<http://www.nuppulinnanlaboratoriopalvelu.fi/>

FINLAND



Nuppulinnan Laboratoriopalvelu Oy
Sallintie 1 C
Kellokoski, Uusimaa
Finland 04500
Tel: +358 40 5305367
E-mail: nuppulinna@dlc.fi
Website: www.nuppulinnanlaboratoriopalvelu.fi

France, Tunisia, Algeria & Switzerland

Diagomics
7 Av. Didier Daurat BP30044
31702 Blagnac Cedex – France
Tel: +33-6-8313-4101
Website: diagomics.com



GEORGIA



Deltamed Georgia
Tsinandali str. 9
Tbilisi 0144
Georgia
Tel: +955-32-292-0910
E-mail: info@deltamedgeorgia.com
Website: www.deltamedgeorgia.com

GERMANY



Medac GmbH
Theaterstrasse 6
Wedel 22880
Germany
Tel: +49-4103 8006-342
E-mail: contact@medac.de
Website: www.medac.de

Biozol Diagnostica GmbH

Leipziger Straße 4
Eching Germany
Tel: +49-89-3799-666-6
Fax: +49-89-3799-666-99

GREECE



Spanos Andreas Leonidas Bioinnotech
Filopappou 34
Athens, Acropolis
Website: www.bioinnotech.gr
Greece 11741
Tel: +30-211-012-9260

GUATEMALA



SKG International
1141 SE 36 St.
Cape Coral, FL
USA, 33904
Website: www.skginternationalgroup.com
Tel: +786-609-4890

HONDURAS



RC Lab S. de R.L.
Col. Alameda #712
Tegucigalpa – Honduras
Tel: +504-232-6398

SKG International

1141 SE 36 St.
Cape Coral, FL
USA, 33904
Website: www.skginternationalgroup.com
Tel: +786-609-4890

HUNGARY



BioMarker Ltd.
BioMarker Kft
Moór János utca 3-5
Gödöllő Ipari Park
2100 Gödöllő
Hungary
Tel: +36-2-841-9986
Website: biomarker.hu

INDONESIA



PT Biosains Medika Indonesia
Jalan Margasatwa Raya 47A Cilandak
South Jakarta, DKI Jakarta
Indonesia 12450
Tel: +62 21 7590 0664
Website: www.biosm-indonesia.com

PT. MAB Supply

Jl. Rawa Bambu No. 15 B, Pasar Minggu
Jakarta Selatan, DKI Jakarta
Indonesia 12520
Tel: +62-21-7884-9315

INDIA



Sentier
Plot No. 780 Block "A", 2nd Floor, 3rd Avenue,
Besides 5th Avenue Bakery
Above Chai Kahani Defence Colony, Saninikpuri
(Post)
Secunderabad-500094
Tel: +919949861367
+91-9959223383

Website: <http://www.sentierhealth.co.in/>

IRAQ



Nearest Scientific Buearu
Al-Mansur, Mutanabi
Baghdad, Iraq
Tel: +964-770-271-6356

ISRAEL



Biotest Ltd.
6 Hasadna St.
Kfar Saba, 442405
Israel
Tel: +972-9-766-7454
Website: www.biotest.co.il

ITALY



DBA Italia
Via Umbria 10
Segrate
Italy 20090
Tel: +390226922300
Fax: +390226923535
Website: dbaitalia.com

Aurogene SRL

Via Dei Lucani 51,53,55,
Rome, Lazio
Italy, 00185
Tel: +39-069-818-5510
Website: <http://www.aurogene.eu>

JAPAN



BIOSJ
138 Enmado-cho
Moriyama City
Shiga, 524-0042 – Japan
Tel: +81-77-582-4463
Fax: +81-77-582-4463
Website: biosj.com

JORDAN



Medical Business Supply Center
Gardens Street Building No.: 80. 2nd Floor
P.O. Box: 2597 Amman 11953
Jordan
Tel: +962 6 4029697

Genetics Company (El Weratha Company)

Wasfi Al-Tal Street
BLDG 195 7th Floor
Amman
Jordan, 11953
Tel: +952-6-553-6402
Website: www.genetics-jo.com/

KAZAKHSTAN



Konsultant Biotech Limited Liability Partnership
26/118 Masanchy Street
Almaty
Republic of Kazakhstan
Tel: +7 727-3798949

LATVIA

Nuppulinnan Laboratoriopalvelu Oy
Sallintie 1C
Kellokoski, Uusimaa
Finland 04500
Tel: +358 40 530536

KENYA

Caroga Diagnostics
Nairobi County, Kenya
Tel: +254 20 6006793
Website: <http://www.carogadx.com>

LEBANON

Kharif for General Trading
Mazraa
Beirut
Lebanon
Tel: +961-130-2583
Website: www.kharif.info

Amedic

El Kharayeb – Main Road – Ein Street
Tyre
Lebanon
Tel: +91 7645-7480

LITHUANIA

Nuppulinnan Laboratoriopalvelu Oy
Sallintie 1C
Kellokoski, Uusimaa
Finland 04500
Tel: +358 40 5305367

LUXEMBOURG

ITK Diagnostic bv
Joh. Enschedeweg 13 / P.O. 73
1422 DR Uithoorn / 1420 AB Uithoorn
The Netherlands
Tel: +31-(0)297-568893
Website: www.itk.nl

MALAYSIA

Ultimed Scientific SDN BHD (UTMS)
No. B-15-1, Jalan C180/1, Dataran C180
Cheras, Selangor,
Malaysia, 43200
Tel: +603-9547-9944
Website: www.ultimed.com.my

Bitu LifeScience Sdn Bhd

C-23-05, Dataran 32, No 2, Jalan 19/1
Petaling Jaya, Selangor D.E.
Malaysia, 46300
Tel: +603-7-960-8286
Website: www.bitlifesc.com.my

MEXICO

Bio SB Mexico
Cantil 175-4, Colonia Jardines del Pedregal
CP 01900, Alcaldía Álvaro Obregón,
Ciudad de México.
Tel: +52 (55) 5264 1692
Whatsapp: +52-56-1623-9840
Website: biosb.mx

MONGOLIA

Monobiopharm Co. LTD
Sukhbaatar District - 4 khoroo 1-69, First 40000
Ulaanbaatar, Mongolia
Tel: +976 99156493
+97611330234
Fax: +976 11 324420
Website: www.monbiopharm.mn

MOROCCO

Progen Sarl
Rue Ibnou Katir, Redidence Al Mawlid
Imm C, #11, 20380
Casablanca – Morocco
Tel: +212-522-98-5005
Fax: +212-522-98-5006

THE NETHERLANDS

ITK Diagnostics bv
Joh. Enschedeweg 13 / P.O. Box 73
1422 DR Uithoorn / 1420 AB Uithoorn
The Netherlands
Tel: +31-(0)297-568893
Website: www.itk.nl

NICARAGUA

SKG International
1141 SE 36 St.
Cape Coral, FL
USA, 33904
Website: www.skginternationalgroup.com
Tel: +786-609-4890

NIGERIA

BioRapid Diagnostics Nig Ltd.
Emaddo Plaza, Suite B010
Beside Fortisbank, Mararaba
Nigeria
Tel: + 234 803 7259813
Website: www.biorapidiag.com

PAKISTAN

Med Lab Services
Office # 1, 1st Floor, ABC Plaza, Commercial Center
Satellite Town, Rawalpindi – Pakistan
Tel: +92 51-442-3895
Fax: +92 51-445-6995
Website: medlab.com.pk

Lab Diagnostic Systems Pvt Ltd

111, Hali Toad, Westbridge 1
Rawalpindi, Punjab
Pakistan 46000
Tel: +92-51-111-145-236

PERU

BioMaxim, S.R.L.
Av. Arequipa 2450, Of. 1109
Lima 14 – Peru
Tel: +511-221-0672
Fax: +511-221-3493

PHILIPPINES

**3-I Philippines Internation
Industrial Investments**
704 Aurora Boulevard
Quezon City, Metro Manila
Philippines zip 1112
Tel: +63-2-721-5211
Website: www.three-i.net

Mediverse Corp.

8th Mabuhay Tower Cebu IT Park
Cebu City Philippines 6000
Tel: +63-32-341-5566
Website: www.mediversecorp.com

POLAND

Elektro Med Grzegorz Palkowski
Ul. Zebierzowska 11
Niepolomice
Poland 32 005
Tel: +48 12 288 91 40
Fax: +48 12 288 91 42
Website: www.elektromed.pl

PUERTO RICO

Caribbean Biotech and Medical Devices
Villas De Cupey
San Juan, Puerto Rico
Tel: +1-787-249-4466

QATAR

Atrium Medical Care WLL
Office 303, GREDCO Building III, Near Aster Clinic
Doha 5483 Qatar
Tel: +974-4436-9438

ROMANIA

Laboratorium Life Science SRL
Str. Baia no. 16, District 1
Bucharest, Sector 3
Romania 013495
Tel: +40-723-515206

SAUDI ARABIA

Anwa'a Al Techniyah Trade EST
King Fahad Rd., Olyia District,
Gulf Commercial Complex
Riyadh, 12221

Saudi Arabia
Tel: +966-11-216-6160
Website: www.typesoftech.com

SERBIA

Magna Pharmacia d.o.o. Beograd
Milutina Milankovica 7b
Novi Beograd Serbia
Tel: +381-11-3622-500
Website: www.magnapharmacia.rs

SINGAPORE


RAS Pharma & Biochem Pte, Ltd.
18, Kaki Bukit Road 3, #02-02,
Entrepreneur Business Centre,
Singapore 415978
Tel: 68427188
Fax: +65 68427688
Website: <http://eraspharma.com/>

SOUTH AFRICA

SMM Instruments (Pty) Ltd
Tillbury Business Park, Unit B7, 1030 16th Road,
Randjespark
Midrand, Johannesburg, Gauteng 1686
Tel: +27-011-540-6000

SOUTH KOREA

J.One LifeScience (J1LS)
104-1207, 8, Hangang-daero 43-gil,
Yongsan-gu, Seoul, 04376,
Republic of Korea
Tel: +82-2-577-8943
Fax: +82-2-577-8946

SRI LANKA

Biomedite (PVT) LTD.
No.276/2A Hospital Rd. Kalubowila.
Dehiwala Colombo
Sri Lanka 10350
Tel: +94 770055028

SWITZERLAND

ABC Biopply
Gewerbstrasse 8
Cham, Switzerland, 6330
Tel: +41-49-834-95-56
Fax: +41-41-743-25-36
Website: biopply.com

Diagomics
7 Av. Didier Daurat BP30044
31702 Blagnac Cedex – France
Tel: +33-6-8313-4101
Website: diagomics.com

Medite Service AG
Binzstrasse 5 Dietikon,
Zurich. CH-8953 Switzerland
Tel: +41-44-744-3020
Website: <https://www.mediteservice.ch/>

SWEDEN

Nuppulinnan Laboratoriopalvelu Oy
Sallintie 1C
Kellokoski, Uusimaa
Finland 04500
Tel: +358-40-5305367
<http://www.nuppulinnanlaboratoriopalvelu.fi/>

TAIWAN

Hong Jing Co. Ltd.
Chung-shan Rd., Young-Ho City,
234 Taipei, Taiwan, R.O.C.
Tel: +886 2 32338585 # 307
Website: www.hongjing.com.tw

TURKEY

BioGen Medikal Cihazlar Tic. A.Ş.
Oruçreis mah. Tekstilent Cad. 12 B Koza Plaza
Kat 18 No: 18334235 Esenler, Istanbul, Turkey
Tel: +90-212-438-6590
Fax: +90 212-438-6594
<http://www.biogen.com.tr/>

Pars Med
Baris mah. Enver Adakan Cad. West Side 1
Turkey 34520
Tel: +90 533-088-8031
Website: www.parsmed.info/

UKRAINE

Global Biomarketing Group
18-a, Mykhayla Boychuka St.
Kyiv, 01103
Ukraine
Tel: +38094 8272381
Website: <http://www.gbglab.com/>

UNITED ARAB EMIRATES

Phoenicia Drug LLC
Deira, Abdulaziz Al majid Building
Dubai U.A.E.
Tel: +971-4-222-5497
Website: www.pddubai.com


UNITED KINGDOM

Stratech Scientific Ltd.
Unit 7 Acorn Business Centre, Oaks Drive,
Newmarket, Suffolk, United Kingdom, CB8 7SY
Tel: +44-163-878-2600
Website: stratech.co.uk

UNITED STATES

Bio SB, Inc.
5385 Hollister Ave. Building 8, #108
Santa Barbatra, CA 93111 – U.S.A.
Tel: +1-805-692-2768
E-mail: sales@biosb.com
Website: biosb.com

Cancer Diagnostics, Inc
4300 Emperor Blvd. #400
Durham, NC 27703
Tel: +1-877-846-5393 x727
Website: cancerdiagnostics.com

URUGUAY

Visur Productos Diagnósticos Limitada
26 de Marzo 1208, Ap # 501
Montevideo, C.P. 11300, Uruguay
Tel: +598 2708-6099
Fax: +598 2708-4739
Website: visurltda.com

VENEZUELA

Distribuidora Macrosearch
Av. Araure, Centro Comercial Chuao, Piso2, #7B
Caracas, Miranda
Venezuela 1060-1061
Tel: +58-212-992-9841


VIETNAM

Technimex JSC
70 Tran Hung Dao str. Hoan Kiem dist.,
Hanoi, Vietnam
Tel: +844 3822-1504 EXT. 35
Fax: +844 3822-0377
Website: technimexvn.com

YEMEN

Ezo Medical Lab
Al-Amal Bank Building - Taiz Street - SANAA
Tel: +967 1 633461
Website: <https://www.ezomedical.com/>

NOTES:



At Bio SB our passion is providing biomedical laboratories with the tools to improve the diagnosis, prognosis and therapies that benefit patients worldwide. Bio SB manufactures and develops products in accordance with FDA QSR 21 CFR Part 820 cGMP and ISO 13485:2016. These guidelines enable us to produce an IVD product that meets the highest in vitro diagnostic standards.

Bio SB Inc.

5385 Hollister Ave., Building 8, #108, Santa Barbara, CA 93111, USA
Tel (USA): 1-800-561-1145 | Tel (Int): +1-805-692-2768 | Fax: 805-692-2769
E-mail: sales@biosb.com | Website: www.biosb.com

MM 0001 - 033121